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
EN Instructions for use

Pump tubing for ulrich medical CT/MRI injectors

About This Document

- These instructions for use are part of the product and they describe its safe use as intended.
- Read the instructions for use before using the product.
- Keep the instructions for use with the product so that they are accessible at all times.

Highlighting and Symbols

Highlighting	Meaning
WARNING	Hazards to persons. Non-compliance can cause death or serious injuries to persons and damage to the product.
→	Handling instructions.

Symbols on the Product and Packaging

Icon	Description	Icon	Description
	Caution		Do not use if package is damaged
	Manufacturer's CE mark		Do not use tubing if protective cap not properly attached
	Product number		Pyrogen-free nepirogeniška
	Lot code		Keep dry
	Use by date ant pakuotės pažymėtas produkto galiojimo laikas		Temperature limit
	Manufacturer		Relative humidity limit
	Consult the instructions for use		Protect from sunlight
	Sterilized using ethylene oxide sterili (simbolis ant pakuotės)		Do not re-use! Pump tubing intended for single use only for 24 hours. vienkartinė

Intended use

The pump tubing for ulrich medical CT/MRI injectors is used for venous-side administration of contrast agents and physiological saline solutions (NaCl) during a CT/MRI examination. The pump tubing creates the connection from the injector (device side) to the patient tubing. The pump tubing can be used for any number of injections within 24 hours.

Product description

On the patient side, the pump tubing consists of a tube with a protective cap which has a male Luer lock connector (XD8003: with integrated check valve. The check valve is self-closing. It only opens when a medium flows in the defined direction. If back pressure is exerted or if there is no media flow, the valve remains closed). The spikes contain an integrated bacteria filter in the supply air duct.

A particle filter in the pump tubing retains foreign bodies. The pressure sensor unit controls the flow of injection volume in the pump tubing. The hygienic safety has been proven by an independent testing facility.

- A Device side with pump tubing
- B Patient side with patient tubing
- 1 Spike for media supply for contrast agent
- 2 Spike for media supply for NaCl
- 3 Spike for media supply for contrast agent
- 4 W-piece (XD 8003), crosspiece (XD 2020)
- 5 Tubing section for roll pump
- 6 Pressure sensor unit
- 7 Particle filter
- 8 Male Luer lock, with check valve (XD 8003), without check valve (XD 2020)
- 9 Protective caps for tubing ends
- 10 Media flow

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Delivery contents and scope of application

Art. no.	Designation
XD 8003	Pump tubing for CT motion injector
XD 2020	Pump hose for XD200x injectors

Product for single use

The product is EO sterilized and delivered in sterile packaging.

- Do not reuse product.
- Use product for a maximum of 24 h
- Do not resterilize product.
- Do not use product from open or damaged sterile packaging.
- Discontinue using product if the use-by date has expired.
- Do not use product that was not stored under the recommended storage conditions.

Safe handling

→ Before using the product, check its operability and proper condition.

- To avoid damage due to improper use, use product according to these instructions for use and the instructions for use of the injector.
- Product and accessories should only be used by persons who have the required knowledge, training or experience.

Risk of injury!

The use of tubing combinations or systems that have not been approved by ulrich medical jeopardizes the safety of the patient and/or the user as well as the proper functioning of the injector.

- Only use original products approved by ulrich medical.

Risk of injury! Air injection, overdose and flow reduction due to pressure regulation if tubing extensions are used.

Tubing extensions jeopardize the safety of the patient and risk contaminating the tubing system.

- Do not use tubing extensions (e.g. "Heidelberg extension").

Risk of injury! Air injection and spraying of contrast agent. Risk of damage to property.

If screwed in too tightly, the thread of the Luer lock connector may break and lead to open cracks.

- Close the Luer lock connector properly.

Risk of injury due to infections and contaminations!

Improper handling jeopardizes the safety of the patient.

- Handle all sterile components carefully (see chart, pos. 1, 2, 3, 8) to ensure sterility.
- Do not use product if the protective caps are missing or are not attached securely.
- If new patient tubing is not connected immediately after disconnecting the patient tubing from the pump tubing: Place a sterile protective cap on the pump tubing.
- Replace and dispose of the pump tubing after every working day at the latest or after 24 hours at the latest.

Risk of injury due to mixing of contrast agents with the Tandem software option!

- If a type of contrast agent other than the type previously used is loaded onto the injector: Change pump tubing.

Use

Preparation

- Remove the product from the packaging.
- Check product before use for damaged parts.
- Do not use product with damaged parts.

Connect product to injector

- Connect product only as directed in the instructions for use of the injector.

Performance testing

- Conduct performance testing only as directed in the instructions for use of the injector.

Operation

- Operate product only as directed in the instructions for use of the injector.
- Replace the pump tubing after every working day at the latest or after 24 hours at the latest.

Validated processing method

Single-use product

WARNING! Patient/user infection due to re-use and/or processing.

If re-used and/or processed, there is a risk of contamination and material fatigue. Soiling, contamination, impaired function may lead to injury, illness or death.

- Do not process product.

Technical service

ulrich GmbH & Co. KG
Injector Systems Customer Service
Buchbrunnenweg 12
89081 Ulm / Germany
Telephone: +49 (0)731 9654-0
Fax: +49 (0)731 9654-2805
Email: injectors@ulrichmedical.com

Technical Data

Parameter	Value
Pressure tightness	20 bar max.
Material	PUR, PVC (DEHP-free, latex-free)
Max. permissible duty	24 hours, for multiple patients
Filling volume (XD 2020)	17.3 mL
Filling volume (XD 8003)	19.5 mL

Storage

- Store product in sterile packaging in a dry, dark and temperature-controlled room.
- Protect product in sterile packaging from environmental influences such as soil, dust, moisture, heat and sunlight.