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<b>Method of sterilization:</b>	<b>Steam, E.T.O.-gas, Formaldehyde</b> <i>Not for electron-beam, γ-rays-irradiation or hydrogenperoxide-sterilisation !</i>
<b>Packaging:</b>	Pouches bundled per 100, sealed into PE-film-bag, Reels wrapped into PE-film, packed into cardboard carton
<b>Sealing temperature:</b>	170° - 190°C
<b>Shelflife:</b>	5 years from manufacture, min. 4 years from delivery

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**Industriestrasse 6**

**91555 Feuchtwangen**

declare under sole responsibility, that our products

**See-through packaging paper and film**

VP-code/quality

**MM70-1/PP50-1**

Article numbers:

**3FKFB210... / 3FKSB220... / 3FKFS230... / 3FKSS230...**

**See-through packaging non-woven and film**

VP-code/quality

**MM60-11/PP50-1**

Article numbers:

**3FKFB240... / 3FKFS240...**

**See-through packaging Tyvek® and film**

VP-code/quality

**TYO-3/PE62-4**

Article numbers:

**3FKFB250... / 3FKFS250...**

are sterile barrier systems, and therefore accessories to medical devices intended by their manufacturer to be used terminally sterilized. As such, the products named above are

medical devices of **class I** according to annex IX, rule 1 of the directive 93/42/EEC,

and meet the basic requirements of the directive 93/42/EEC and its transposition into national laws.

The conformity assessment procedure is carried out in accordance with annex VII of the directive 93/42/EEC.

Applied standards are

**DIN EN ISO 11607 - part 1 und part 2**

**DIN EN 868 - part 5**

*Any modifications to the device which are not authorized will invalidate this Declaration of Conformity.*

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D-91555 Feuchtwangen · Industriestrasse 6

Fon +49 (0) 9852 901-0 · Fax -174

Heiko Brandl  
Managing Director  
Medical Packaging

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**Industriestrasse 6**

**91555 Feuchtwangen**

herewith we are able to confirm, that the process indicators for

- Steam-sterilization,
- Gas-sterilization,
- Plasma-sterilization and
- Formaldehyd-sterilization

printed on VP packaging material do comply to the performance requirements of

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Thomas Lechner  
Head of Management Systems VP-Group

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Günter Späth  
Head of Quality Assurance VP-Group

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Thomas Lecher

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Gunter Spath

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
## See-through-reels flat and gusseted paper/film

### with STEAM-, EO- and FORM-indicator - ISO 11607, CE

See-through pouches and reels (catalogue page 2-10 to 2-30) are the universal packaging solution for sterilisation applications and perfectly meet the packaging demands of almost all light and medium weight instruments and sets.

The sterilisation proof 70 g/m<sup>2</sup> heavy-duty paper is an effective barrier against germs and bacteria, yet it permits the transmission of air, steam or sterilisation gases. The three indicators for STEAM-, EO- and FORM-sterilisation are imprinted outside the filling space.


#### See-through-reels flat - paper/film - length: 200 m

Article-No.	Size in cm	Pack (reels)
3FKFS230102	5	6
3FKFS230104	7,5	4
3FKFS230106	10	3
3FKFS230108	12,5	2
3FKFS230110	15	2
3FKFS230112	20	2
3FKFS230114	25	2
3FKFS230116	30	1
3FKFS230118	38	1
3FKFS230120	42	1

#### See-through-reels gusseted - paper/film - length: 100 m

Article-No.	Size in cm	Pack (reels)
3FKSS230204	7,5 x 2,5	4
3FKSS230206	10 x 5,0	3
3FKSS230210	15 x 5,0	2
3FKSS230212	20 x 5,0	2
3FKSS230214	25 x 6,5	2
3FKSS230216	30 x 6,5	1
3FKSS230218	38 x 8,0	1

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*Technical data is submitted to the usual variations of operating procedures. Therefore it is not compulsory and does not exempt from validation and operating tests.*

*Tolerances which differ from this product specification including the shipping container and the labelling system have to be confirmed - after negotiation - in our article description which includes also the artwork.*

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*Packaging material for disposable and for reprocessible medical products.*

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*Steam-, ethyleneoxide- and formaldehyde sterilization*

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*The packaging material complies with the DIN EN ISO 11607 part 1 and part 2*

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*Concerning points 3 to 8 of this product specification it is important to know that these stericlin® products are standardized products. Due to technical reasons it is not possible to develop or discuss customized specifications regarding the above mentioned points.*

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*The intended seal strength is a result of the variable sealing parameters temperature, pressure and dwell time resp. sealing speed. As a result of validation of individual sealing processes the values of sealing temperature may be lower or higher.*

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*Protect from daylight or direct sunlight.*

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*Do not open protecting packaging until the point of use of the product.*

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*SH*

## See-through packaging from steriCLIN® – unsurpassed in quality and performance


Sterile barrier systems (SBS) from steriCLIN® are the result of decades of research and development by our safety experts. We have shaped the development of modern sterilisation solutions like no other manufacturer. Today, our SBS products are the benchmark standard for an entire industry. Quality you can trust without reservation. SteriCLIN® – Responsibility is our commitment.

steriCLIN®


Our trademark is found on all see-through packaging products made in Feuchtwangen.

Product ID

Compliance with standards can be recognised from the complete coding: Standard codes, size, REF = Article code, LOT = Batch code, Date of manufacture in MMYYYY format.

Warning


Do not use if pack is damaged.


EO-indicator

Colour change white → yellow  
The process indicator complies with all requirements of the ISO 11140 standard and maintains the required size of 100 mm².

Bags and reels

with gussets are particularly suitable for packaging voluminous medical articles.

FORM Indicator


Colour change purple → green  
The process indicator complies with all requirements of the ISO 11140 standard and maintains the required size of 100 mm². The recesses increase stability of the sealing seam.

Peel Symbol


Indicates the peeling direction when opening the tubing segment. To ensure precise traceability, the number on the peel symbol shows the number of the roll on the production line.

STEAM Indicator

Colour change pink → brown  
The process indicator complies with all requirements of the ISO 11140 standard and maintains the required size of 100 mm².

Our poster on the correct use of see-through packaging products is indispensable for every CSSD. We would be pleased to send you a copy.

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