EVIS EUS Endoscopic Ultrasound Center

EU-ME3

Advancing the Dimensions of Endosonography

Advancing the Dimensions of Endosonography

Focused on Your Expertise

With more functions, better visualization, and enhanced operability, the EU-ME3 expands the dimensions of endosonography.



Improved Ultrasound Imaging

Enhanced Visualization


Enhanced B-mode

The EU-ME3 provides outstanding image quality and functionality – compatible to a high-end ultrasound center – in a compact body. B-mode image quality has been substantially enhanced compared to our conventional processor (EU-ME2).

Improved Elastography


The EU-ME3 features an elastography function which visualizes the amount of strain in the tissue (tissue stiffness) during compression and retraction, making it possible to obtain more information about tissue properties.

Contrast Harmonic Echo (CHE)

Contrast Harmonic Echo (CHE) images harmonic components from ultrasound contrast agents.

The newly added C-THE mode images signals from biological tissue and the contrast.





Tissue Harmonic Echo (THE)

When ultrasound waves are propagated through tissue, distortion is produced and harmonic components are generated. The Tissue Harmonic Echo (THE) mode uses these components to build an image of the targeted area, providing a more detailed granular depiction. Advantages of harmonic imaging include improved resolution, improved signal-to-noise ratio, and fewer artifacts.

Doppler Modes

The EU-ME3 offers three basic Doppler modes to distinguish blood flow more clearly – Color Flow, Power Flow, and Pulsed Wave Doppler (PWD). Doppler modes can be used to support safer procedures, benefitting both the patient and the physician.

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

Designed for Enhanced Usability

Enhancing Functionality

Shear Wave Quantification (SWQ)
SWQ provides an absolute value of tissue stiffness within a region of interest. It performs this quantitative tissue assessment by calculating the propagation velocity of shear waves, generated from a push-pulse.


Elastography (i-ELST)
i-ELST is a new technology incorporated into the EU-ME3 that makes it easier to display elastic images, even when displacement due to pulsation is modest.

Eil. Nr. 1.1.5

s-FOCUS

The EU-ME3 is equipped with an s-FOCUS mode that reduces the change in resolution with distance from the ultrasound transducer surface. s-FOCUS eliminates the need to manually adjust the focal zones during the procedure.

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Enhancing Versatility

Wide Range of Compatibility

Integrating both electronic and mechanical scanning technologies, the EU-ME3 is compatible with echoendoscopes and miniature probes, creating a total endosonography solution for a full range of applications.

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Comparison of Ultrasound Functions

	EU-ME2	EU-ME2 PREMIER	EU-ME2 PREMIER PLUS	EU-ME3
Eil. Nr. 1.1.4 a B-mode	✓	✓	✓	✓
Eil. Nr. 1.1.4 d THE (Tissue Harmonic Echo)	-	✓	✓	✓
Eil. Nr. 1.1.4 b Flow	✓	✓	✓	✓
Eil. Nr. 1.1.4 c PWD (Pulsed Wave Doppler)	✓	✓	✓	✓
Eil. Nr. 1.1.4 e CHE (Contrast Harmonic Echo)	-	✓	✓	✓ (Software Option)
Elastography	-	-	✓	✓ (Software Option)
Eil. Nr. 1.1.4 f SWQ (Shear Wave Quantification)*	-	-	-	✓ (Software Option)

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Specifications

EVIS EUS ENDOSCOPIC ULTRASOUND CENTER OLYMPUS EU-ME3

Power Supply	Voltage		220 – 240 V AC	
	Voltage fluctuation		Within ±10%	
	Frequency		50/60 Hz	
	Frequency fluctuation		Within ±1 Hz	
	Consumption electric power		340 VA	
Size	Dimensions	Main unit	371 (W) × 175 (H) × 480 (D) mm 445 (W) × 184 (H) × 530 (D) mm (max.)	
		Keyboard	392 (W) × 39 (H) × 210 (D) mm	
	Weight	Main unit	21.5 kg (without software option case) 21.8 kg (with software option case)	
		Keyboard	2.5 kg	
Classification	Type of protection against electric shock	Class I		
	Degree of protection against electric shock or applied part	TYPE BF applied part where no classification mark appears, the device is a TYPE BF applied part.		
	Degree of protection against explosion	The Ultrasound Center should be kept away from flammable gases.		
Ultrasound Scanning Format			Mechanical scanning, electronic scanning	
1.1.2.1	Display mode	B-mode	Eil. Nr. 1.1.3.1	
	Scanning	Radial scanning, helical scanning	Eil. Nr. 1.1.3.2	
	Usable frequencies	12 MHz, 20 MHz	Eil. Nr. 1.1.3.3	
Mechanical Scanning	Display range	2, 3, 4, 6, 9, 12 cm		
	Display processing	Rotation	Rotatable	
		Display area	Full circle, bottom sector, top sector, scroll	
		Direction	Normal/Inverse	
	Cine memory	Over 1,500 frames storable depending on the conditions. Cine review function		
3D	3D display, MPR display	Eil. Nr. 1.1.3.4		
1.1.2.2	Measurement	Distance, area, circumstance		
	Display mode	B-mode, FLOW mode, PW mode, CHE mode, ELST mode		
	Scanning	Radial scanning, curved linear array scanning		
	Usable frequencies	5 MHz, 6 MHz, 7.5 MHz, 10MHz, 12MHz	Eil. Nr. 1.1.4.2	
	Display range	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12 cm		
	Display processing	Rotation	Rotatable during radial scanning	
		Display area	Radial: Full circle, bottom sector, top sector, scroll, Curved linear array: Fixed	
		Direction	Normal/Inverse	
	Cine memory	Over 2,000 frames storable depending on the conditions. Cine review function		
	Focus	Auto preset	s-FOCUS, AUTO, MANUAL	Eil. Nr. 1.1.5
Focus settings		Focus location and Focus number adjustable.		
Electronic Scanning	FLOW mode	COLOR-FLOW mode, POWER-FLOW mode, H-FLOW mode		
	PW mode	B+PW, COLOR+PW, POWER+PW, H-FLOW+PW		
	Measurement	Distance, area, circumstance, PW measurement		
	THE mode	THE-P mode, THE-R mode		
	CHE mode (Software Options)	Display pattern	CHE, C-THE	
		Preset (CH agent type)	2 types (Low acoustic pressure, Middle acoustic pressure), selectable	
		Frequency selection	2 types (CHE-P, CHP-R)	
	ELST mode (Software Options)	Pressurization guide	Pressurization bar, Strain graph	
		Strain ratio	Measures strain or ratio of strain of 2 areas.	
	SWQ (Software Options)	Calculates and displays transmission speed and elasticity of shear wave in ROI.		
	Recording Data	Data format	Movie data	AVI Eil. Nr. 1.1.6
		Keyboard	Built-in track pad and touch panel. Eil. Nr. 1.1.7	
		Recording device	DVR	
Ancillary Equipment	Video system center	Monitor display selection	Endoscopic/ultrasound image Eil. Nr. 1.1.8.2	
		Sub screen	Endoscopic image can be displayed in sub screen.	
		Patient data	Patient data can be shared with video system center.	

Trademarks, product names, logos, or trade names used in this document are generally registered trademarks or trademarks of each company.
Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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
EUS Diagnostics and Therapeutics

EUS Diagnostics and Therapeutics


Precise EUS Diagnosis: A Compact Ultrasound Processor Dedicated to EUS Applications, EUS Scopes with Outstanding Maneuverability and a Comprehensive Variety of EUS Needles.

The **EU-ME3** is a high-quality, compact ultrasound processor for use with Olympus endoscopic and endobronchial ultrasound equipment. Various new features help to explore the future of endosonography and bring real clarity to your EUS procedures.


- **Enhanced B-mode:** The improved B-mode image quality supports the localization of tumors and more accurate identification of tissue properties and boundaries.
- **Shear Wave Quantification (SWQ):** SWQ provides an absolute value of tissue stiffness within a region of interest by calculating the propagation velocity of shear waves, generated from a push-pulse.
- **Contrast Harmonic EUS (CHE)*:** Technology designed to depict higher harmonics, the newly added C-THE mode images signals from biological tissue and the contrast.
- **Elastography:** An advanced form of ultrasound, elastography displays the relative stiffness of tissues by taking advantage of the deformation caused by the compression or vibrations generated by the heartbeat or vascular pulsations.
- **s-FOCUS:** Reduces the change in resolution with distance from the ultrasound transducer surface. s-FOCUS eliminates the need to manually adjust the focal zones during the procedure.


B-mode

SWQ

Contrast Harmonic EUS (C-THE)

Elastography

* Regulations and usage of ultrasound contrast agents vary according to the country where they are used and the type of agents. Please use the ultrasound contrast medium according to the instructions provided with the products.

EVIS EUS

Designed for Enhanced Usability

Keyboard Usability

The keyboard was designed with a simple layout in mind and includes a user-friendly built-in touch panel, LED backlit keys and a trackpad for ease of use and cleaning. The larger LCD touch panel allows for a greater range of functions to be displayed at one time.

Ease of Targeting

The position and size of the Doppler region of interest (ROI) can be conveniently adjusted with a trackpad or buttons on the touch panel.

Customizable Features

Software options are available to meet the needs of any facility. Because the functions are optional, you can select and add the necessary functions according to your needs and budget.


Doppler Modes


The EU-ME3 offers three basic Doppler modes to distinguish blood flow more clearly - Color Flow, Power Flow, and Pulsed Wave Doppler (PWD). Doppler modes can be used to support safer procedures, benefitting both the patient and the physician.

In addition to the three basic Doppler modes, the EU-ME3 also features H-Flow. H-Flow is a more sensitive Doppler mode that shows directional blood flow with less blooming. It is especially useful for imaging small vessels around the tip of the echoendoscope.


Color Flow


Power Flow

Pulsed Wave Doppler

H-Flow

EUS Scopes

A Complete EUS Scope Lineup – Each Scope Is a High-Precision Instrument Helping You Find the Right Diagnosis.

The new **GF-UE190*** is an electronic radial-scanning scope with a 360° ultrasound view with advanced new features. It offers impressive and flexible performance with greater safety and enhanced ease of use, complemented by improved ultrasound and endoscopic image quality.

- Carefully redesigned to improve insertion and maneuverability for increased safety and patient comfort.
- Increased penetration depth and higher resolution for improved ultrasound images.
- Unique D-Cut for easy balloon removal and detachable ultrasound cable for ease of handling.

Eil. Nr. 1.2


The **GF-UCT180** delivers high-quality ultrasound images with greater B-mode imaging depth, safe device control through the round transducer design and a remarkably short, rigid distal end.

- Wide 180° ultrasonic view.
- 3.7 mm instrument channel and optimized forceps elevator to ensure reliable control.
- Ideal for EUS-guided FNA.

The **TGF-UC180J** forward-viewing scope delivers exceptional maneuverability and control ensured by the short, rigid portion, straight channel and wide 180° angulation capability.

- Allows a wide variety of therapeutic and FNA applications.
- Extremely short distal end and wide angulation.
- Straight channel port for better control of EndoTherapy devices.

Wide Selection of Ultrasonic Miniature Probes

The MAJ-1720 probe-driving unit has been designed to work with the wide selection of ultrasonic probes offered by Olympus. Available in a variety of frequencies, outer diameters and tip designs, any one of these can be combined with the MAJ-1720 and your endoscope to match the specific requirements of any given procedure, including:

- Intraductal ultrasonography (IDUS).
- Dual-plane reconstruction (DPR).
- GuideSheath procedures.

EUS-FNA and -FNB Needles

Exceptional Flexibility for Smooth Insertion

Multilayered coil sheath and nitinol* needle tip ensure exceptional flexibility for smooth insertion even when the scope is angulated.

Excellent Force Transmission

A coil sheath with multilayered construction facilitates excellent force transmission from the handle to the needle tip for improved ease of puncturing.

No Needle Deformation after Multiple Passes

Nitinol* retains its shape without deformation, even after multiple passes. This ensures the needle follows a predictable and consistent trajectory.

Enhanced Visibility under Ultrasound



An echogenic dimple finish on the needle tip enhances ultrasound visibility and ensures precise and confident targeting of fine needle aspiration.

Smooth Puncturing and Side Hole Option


A Menghini-type tip with sharp edges is designed for tissue core sampling and facilitates smooth puncturing, even at an oblique angle.

Designed for FNB Intended Use

EZShot 3 Plus is available in 19G, 22G and 25G and has received CE certification for fine needle aspiration (FNA) and fine needle biopsy (FNB).

* EZShot 3 Plus 19G and 22G.

Endoscopic Ultrasound

Comparison of Ultrasound Functions				EVIS EUS
	EU-ME2	EU-ME2 Premier	EU-ME2 Premier Plus	EU-ME3
B-mode	✓	✓	✓	✓
Tissue Harmonic Echo (THE)	✓	✓	✓	✓
Flow	✓	✓	✓	✓
PWD (Pulsed Wave Doppler)	✓	✓	✓	✓
CHE (Contrast Harmonic Echo)	-	✓	✓	✓ (Software Option)
Elastography	-	-	✓	✓ (Software Option)
SWQ (Shear Wave Quantification)*	-	-	-	✓ (Software Option)

* For GI. Only compatible with GF-UCT180/260 and GF-UE190/290.

Endoscopic Ultrasonic Probes		Eil. Nr. 1.3.		Eil. Nr. 1.4.
		UM-S20-17S	UM-S20-20R	UM-3R
Display Mode	Skenavimo režimas	B-Mode	B-Mode	B-Mode
Scanning Mode		Mechanical radial scanning	Mechanical radial scanning	Mechanical radial scanning
Working Length (mm)		2150	2050	2050
Total Length (mm)		2225	2140	2140
Scanning Direction		360°	360°	360°
Ultrasonic Frequency	Skenavimo dažnis	20 MHz	20 MHz	20 MHz
Insertion Tube Max. Outer Diameter (mm)	Įvedamo vamzdelio diametras	ø 1.8	ø 2.0	ø 2.5
Insertion Tube Outer Diameter (mm)		ø 1.4 (distal end 1085 mm) ø 1.7 (on connector side)	ø 1.7 (distal end 850 mm) ø 2.0 (on connector side)	ø 2.4
Scope Compatibility	Suderinamas su endoskopais	2.2 or more	2.2 or more	2.8 or more
Contact Method for ROI		Direct contact	Direct contact	Direct contact

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Eil. Nr. 1.2

Ultrasonic Videoscopes for EUS

	GF-UE190	GF-UCT180	TGF-UC180J
Scanning Mode	Electronic, 360° radial array	Electronic, curved linear array	Forward-viewing, electronic, curved linear array
Insertion Tube OD	10.9	12.6	12.6
Channel ID (mm)	2.2	3.7	3.7
Bending (Up/Down, Left/Right)	130°/90°, 90°/90°	130°/90°, 90°/90°	180°/90°, 90°/90°
Working Length (mm)	1250	1250	1245
Ultrasound Field of View	360°	180°	90°
Endoscopic Field of View (Forward/Oblique)	100°/50°	100°/55°	120°/0°
Depth of Field (mm)	3–100	3–100	3–100
Frequency (MHz)	5/6/7.5/10/12*	5/6/7.5/10/12*	5/6/7.5/10/12*

* With EU-ME3. Please contact your Olympus representative for compatibility with other systems.

Single-Use FNA and FNB Needle

EZShot3Plus


Article Name	Order Number	Quantity	Min. Working Channel Ø	Working Length	Max. Needle Length	Needle Ø	Stylet Tip Shape
NA-U200H-8019	N5409330	5	2.8 mm	1400 mm	80 mm	19G	Rounded
NA-U200H-8019S	N5409530	5	2.8 mm	1400 mm	80 mm	19G with side hole	Rounded
NA-U200H-8022	N5409730	5	2.8 mm	1400 mm	80 mm	22G	Rounded
NA-U200H-8022S	N5409930	5	2.8 mm	1400 mm	80 mm	22G with side hole	Rounded
NA-U200H-8025	N5783830	5	2.8 mm	1400 mm	80 mm	25G	Rounded

CE 0197 CE 0344 CE 2797

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.



OLYMPUS EUROPA SE & CO. KG
Wendenstraße 20, 20097 Hamburg, Germany
www.olympus-europa.com

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