
Specifications

	LMD-XH550ST	LMD-XH320ST
Picture Performance		
Panel	TFT Active Matrix LCD	
Picture Size (Diagonal)	1387.8 mm (54 3/4 inches)	800.757 mm (31 5/8 inches)
Effective Picture Size (H x V)	1209.6 x 680.4mm (47 5/8 x 26 7/8 inches)	697.92 x 392.58 mm (27 1/2 x 15 1/2 inches)
Pixel pitch	0.315 x 0.315 mm (0.0124 x 0.0124 inches)	0.18175 x 0.18175 mm (0.00716 x 0.00716 inches)
Resolution (H x V)	840 x 480 pixels	840 x 480 pixels
Aspect	3.1.3	16:9
Pixel Efficiency	99.99%	99.99%
Backlight	LED	LED
Panel Technology	LCD with IPS	LCD with IPS
Luminance (Panel Specification)	550cd/m ² (Typical) 1500cd/m ² (Typical, Peak)	450cd/m ² (Typical) 1750cd/m ² (Typical, Peak)
Contrast Ratio	1,000,000:1 (Typical)	1,000,000:1 (Typical)
Colors	10 bit colors (1,073,741,824)	10 bit colors (1,073,741,824)
Panel Frame Rate	50/60 Hz	50/60 Hz
Viewing Angle (Panel Specification)	89°/89°/89°/89° (Typical) (up/down/left/right, contrast > 10:1)	
Vertical Viewing Angle (3D Mode)	37° at a viewing distance more than 1,000mm, crosstalk ratio less than 7% (Typical)	32° at a viewing distance more than 690mm, crosstalk ratio less than 7% (Typical)
Gamma	Endoscope, 1.8, 2.0, 2.2, 2.4, 2.6, DICOM, HLG	
Input		
HDMI Input	HDMI (x1) (HDCP 2.3 correspondence)	
DVI-D Input	DVI-D (x1) (TMDS single link, HDCP 1.4 correspondence)	
SDI Input	BNC (x1) 3G/HD/SD-SDI	
Display Port	Display Port (x1) (SS1 HDCP 1.3 correspondence)	
Serial Remote (LAN)	D-sub 9-pin (RS-232C) (x1) RJ-45 (x1) (Ethernet, 10BASE-T/100BASE-TX)	
Remote	Stereo mini jack (x1)	
AC Input	AC input connector (x1) AC 100-240 V, 50/60 Hz	
DC Input	DC input connector (x1) DC 26 V	
Output		
SDI Output	BNC (x1) 3G/HD/SD-SDI	
CLONE output	BNC (x1) 3G/HD/SD-SDI	
DC 5 V / 12 V Output	5 V Output (x1) 2.0 A max 12 V Output (x1) 2.5 A max	
General		
Power Requirements	AC IN: 100 V - 240 V, 50/60 Hz, 3.1 A - 1.1 A	AC IN: 100 V - 240 V, 50/60 Hz, 2.4 A - 1.1 A DC IN: 26 V, 8.8 A (max) (Supplied from AC adaptor)
Power Consumption	Approx. 290 W (max.)	
Operating Temperature	0°C to 40°C (32°F to 104°F)	
Operating Humidity	30% to 85% (no condensation)	
Storage/Transport Temperature	-20°C to +60°C (-4°F to +140°F)	
Storage/Transport Humidity	20% to 90%	
Operating/Storage/Transport Pressure	700 hPa to 1060 hPa	
Dimensions (W x H x D)	1249.6 x 763.3 x 92.8 mm (49 1/4 x 30 1/8 x 3 3/4 inches)	753.9 x 476.3 x 84.7 mm (29 3/4 x 18 7/8 x 3 3/8 inches)
Mass	Approx. 27 kg (Approx. 59 lb 8.4 oz)	Approx. 10.5 kg (Approx. 23 lb 2.38 oz)
Mounting (W x H)	200 x 200 mm 400 x 200 mm	100 x 100 mm
IP Code	IP45 (front)/IP32 (others)	
Supplied Accessories		
	AC power cord (1)	
	Plug holder for the supplied AC power cord (2)	
	Before Using This Unit (1)	
	Instructions for Use (1)	
	CD-ROM (including the Instructions for Use) (1)	
	Service Contact List (1)	
	Information for Customers in Europe (1)	
	M6 x 12 mm Screw (4)	M4 x 12 mm Screw (4)
	3D Eye Shield Kit: CFV-E30SK (1)	
	Instructions for Use of the 3D Eye Shield Kit (1)	
Optional Accessories		
	3D Eye Shield Kit: CFV-E30SK	
	2D Eye Shield Kit: CFV-E20SK	
	3D Eye Shield: CFV-E30D	
	2D Eye Shield: CFV-E20D	
	Shield Cover: CFV-505C	
	Shield Frame: CFV-B100	
	Foot Switch FS-24	
	AC Adaptor AC-300MD	
	Monitor Stand SU-600MD	


These products are distributed to US and EU as medical devices. They satisfy product safety standards (e.g. IEC 60601-1). For more details, please contact your nearest sales office or an authorized dealer.

Dimensions

LMD-XH550ST


LMD-XH320ST

Unit: mm (inches)

Optional Accessories

- CFV-E30SK 3D Eye Shield Kit
- CFV-E20SK 2D Eye Shield Kit
- CFV-E30D 3D Eye Shield
- CFV-E20D 2D Eye Shield
- CFV-B100 Shield Frame

Distributed by

N8602544-042022

Innovation by Sony & Olympus

©2022 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features, design, and specifications are subject to change without notice. The values for mass and dimension are approximate. Some images in this document are simulated. "SONY" is a registered trademark of Sony Corporation. "A.I.M.E." is a trademark of Sony Corporation. All other trademarks are the property of their respective owners. The "Innovation by Sony & Olympus" phrase and "Innovation by Sony & Olympus" are makers that indicate advanced medical imaging solutions from Sony and Olympus. This device is not available in some areas.

1 psl.

High Quality 4K UHD 3D/2D 55-inch and 32-inch Medical Monitor

The LMD-XH550ST and LMD-XH320ST medical grade monitors with high brightness, high contrast and wide color gamut can display very high quality 4K UHD color video images in 3D and 2D. High-resolution 4K and 3D imaging provides depth perception and spatial orientation suitable for use in hospital operating theaters, surgical centers and classrooms for education.

Superb Picture Quality

4K Ultra HD Resolutions

4K resolution (3840 x 2160 pixels) on a large 55-inch or 32-inch screen offers high picture quality.

High Brightness and High Contrast

The LMD-XH550ST and LMD-XH320ST feature Sony's unique Local Dimming Technology, which controls the brightness of the LEDs depending on the scene. Thanks to this, the monitors can achieve a peak brightness of more than 1,500 cd/m² (panel specification) and a contrast ratio of 1,000,000:1.

HDR-compatible

The LMD-XH550ST and LMD-XH320ST have the HDR display function. When receiving the HLG (Hybrid Log-Gamma) signals, the monitors can reproduce the images more realistically. Greater detail capture in the highlights and dark areas - no detail remains undiscovered in the dark or lost in overexposed highlights.



SDR

HDR

Wider Color Gamut

The LCD panel and signal processing technology offer a wide color gamut. This conforms to the ITU-R recommendation BT.2020 with a wider color space than the BT.709 color space, allowing for better color reproduction and realistic visualization.

Anti-Reflection Structure

The LMD-XH550ST and LMD-XH320ST are designed to minimize reflected glare in brightly lit ORs. The panel structure achieves clear, high-contrast images by controlling reflection and minimizing light dispersion in the LCD panel.

3D Features

Various Display Formats in 3D

The LMD-XH550ST and LMD-XH320ST can display 3D and 2D signals in HD and 4K. These monitors offer SDI (12G/3G), Display Port, DVI, and HDMI terminals supporting a variety of different 2D and 3D signal formats such as side-by-side, line-by-line, and top & bottom. Furthermore, dual stream is available when accepting an HD 3D signal via SDI. 4K dual stream via two 12G-SDI is also available (XH550ST only). The 3D/2D operating mode can be changed directly with the button on the front panel and the mode is indicated on the display.

*SDI Only

Advanced Features

Advanced Image Multiple Enhancer (A.I.M.E.™)

Sony's unique A.I.M.E. technology enhances the displayed images by adjusting the color, contrast and visibility of shadow areas for more comfortable viewing.

A.I.M.E. OFF

A.I.M.E. ON

Continuous Auto Focus in True 4K to Minimize Distractions

Experience the world's first 4K camera head with Continuous Auto Focus (C-AF), which automatically tracks the target object and keeps it in focus without the need to rely on manual adjustment.

EDOF (Extended Depth Of Field) to broaden the focus area

EDOF expands the range of surgical viewing area without compromising resolution and brightness, promoting accurate diagnosis and precise treatment. The EDOF technology addresses traditional camera head challenges such as out-of-focus and blurry/distorted images, by providing a deeper depth of focused field and clearer view of anatomy (as compared to previous generation camera heads).

True 4K Image Quality for More Precise Clinical Images

Our exclusive Full 4K image sensor from Sony produces Full 4K image quality, enabling the Camera Head to deliver incredibly fine details that conventional full HD images cannot provide.

Cost-effective Reprocessing Options

Reprocessing is easy and cost-effective due to Autoclave compatibility.

Specifications		
Size	Camera head dimensions (unit: mm)	
	Camera head weight	270 g
Observation	Cable	ø 5.1 mm x 3 m
	Pickup system	CMOS image sensor
NBI Observation Mode*	Focal length	f = 19.7 mm
		Available
IR Observation*		Available
Electronic Shutter Function*		Available
Electronic Zoom Function*		Available
Cleaning/Disinfection/Sterilization	Cleaning/disinfection	Immersible in detergent/disinfectant solution, Washer Disinfectant
	Sterilization	Autoclave/V-PRO®/STERRAD®
Classification (Electromedical Equipment)	Degree of protection against electric shock	TYPE CF (the entire camera head)
	Degree of protection against explosion	The product should be kept away from flammable gases.

* For more details, refer to the instruction manual for the OTV-S700.

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f=19.7 mm

3.0 Focusing

Focus can be adjusted when the diopter of the connected rigid endoscope is -3.0 to +0.5 (1/M).

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See the Invisible – IR

Fluorescence-Guided Surgery with IR


- ✓ For better perfusion control and easier identification of biliary structures.
- ✓ Upgradable feature to be activated via software licensing.
- ✓ Three IR modes available with VISERA ELITE III, to suit various requirements.

3.2 4.1

Fluorescence imaging during white-light observation (full color IR white-light overlay mode).

Partial white light and IR light at the same time (IR magenta mode).

Pure black-and-white IR light (IR independently mode).

3.2 4.3

The Power of Fluorescence-Guided Surgery in Detail

Choose the IR gain that suits best and experience three selectable IR modes.

Full Color IR White-Light Overlay Mode

- ✓ providing real-time 4K fluorescence imaging during white-light observation

IR Magenta Mode

- ✓ Providing the classically known magenta IR image with fluorescence

IR Independently Black-and-White-Only Mode

- ✓ allowing to generally check blood perfusion with independent infrared light

IR black and white only

IR white-light overlay

IR magenta


Adjustment of the IR Gain

Adjust the visibility of the fluorescence by controlling IR gain to produce the best visualization regardless of scope, tissue or environment changes. The IR gain can be adjusted in three levels: low, medium, and high, and the higher the level, the stronger the fluorescence is expressed.

3.2 7.2

3.3 1 ENDOEYE 3D 10 MM

Video Telescope

Specifications - Observation

[ZPS10413]

1.0 Outputting videos

1.1 Image output terminals

- The following video output terminals are provided.

Terminal	Purpose	Standard
12G-SDI OUT (4K) output 1	For 4K main monitor	SMPTE ST 2082
12G-SDI OUT (4K) output 2	For 4K sub monitor	SMPTE ST 2082
3G-SDI(4K) A/B/C/D	For 3G-SDI Quad-Link	SMPTE424M
3G/HD-SDI output 1	For HD main monitor	SMPTE424M/SMPTE292M
3G/HD-SDI output 2	For HD sub monitor	SMPTE424M/SMPTE292M
3G/HD-SDI Recorder output	For image recording device	SMPTE424M/SMPTE292M
HD-SDI printer output	For printer	SMPTE292M

1.2 Supported video outputs

- This product supports the following video outputs.

Video	2D/3D *1	Output signal	Transmission type	HDR (BT.2100)	Standard/ Color gamut	Aspect ratio	Resolution			
							[HxV]	[Hz]	i/p	
HD	3D	3G-SDI	3G-SDI LevelB Dual Stream	-	Scope connection : BT.709 Displaying color bars : BT.809 *2	-	1920x1080	59.94/50	i	
		HD-SDI	Side by Side						p	
	2D	3G-SDI	3G-SDI LevelA						i	
		HD-SDI	-						i	
4K	3D	12G-SDI	2-Sample Interleave Division Line by Line	ON/OFF	HD scope connection : BT.709 4K scope connection : BT.2020 Displaying color bars : BT.709 *2,3	16:09	3480x2160	59.94/50	p	
		3G-SDI Quad-Link		-						
	2D	12G-SDI		2-Sample Interleave Division						ON/OFF
		3G-SDI Quad-Link								-

*1: "2D" is for applications where only 2D videos are output.

"3D" is for applications where if 3D output is combined, 3D observation can be turned off.

*2: Scope refers to both camera head and videoscope.

*3: BT.709 is applied during frozen image playback even when a 4K scope is connected.

Important: MAJ-2511 (OTV-S700 Upgrade Pack 3D) is required for enabling the 3D observation function.

Kommentaras:
 "p" - progresivis
 "i" - interlaced

ENDOEYE 3D 10 MM

3D-Image Rotation without Loss of Horizon

- Change the view direction while maintaining horizontal orientation of images, enabling a continuous critical view in 3D

Chip-on-the-Tip Technology

- Bright, clear, and natural 3D depth perception
- Focus-free handling; no manual focusing required

Autoclave Compatibility

- Reduced costs compared with other sterilization methods

Full Control During Surgery

- Three programmable remote control buttons

Plug and Play

- Simple to set up and easy to use

Ordering Information

Order Number	Outer Diameter	Direction of View	Field of View	Working Length	Max. Insertion Portion Width	Cable Length	Weight
WA50080A	10 mm	0°	67° <i>3.3.5</i>	330 mm	10.07 mm <i>3.3.4</i>	2780 mm	≤1500g
WA50082A	10 mm	30° <i>3.3.3</i>	67°	335 mm	10.20 mm	2780 mm	≤1500g

Specifications

3.3.1


Reprocessing Methods	Autoclavable, ETO, STERRAD, washer-disinfector		
Ambient Conditions	Operating conditions	Ambient temperature	10 to 35°C (50 to 95°F)
		Relative humidity	30 to 85%
	Storage conditions	Temperature	10 to 40°C (50 to 104°F)
		Relative humidity	30 to 85%
Transport conditions	Temperature	-47 to 70°C (-53 to 158°F)	
	Relative humidity	10 to 95%	
Classification of ME Equipment and ME Systems	Protection against electric shock: Classification of applied parts	Type BF	
	Protection against harmful ingress of water or particulate matter	IPX7	

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS


OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
 www.olympus-europa.com

Main body

3.3.6


- 1) Remote control button SW1
- 2) Remote control button SW2
- 3) Adjustment ring
- 4) Remote control button SW3

**10mm ENDOEYE 3D**

- ① **WA50080A** Video telescope "ENDOEYE 3D",
10 mm,
0°
▪ Autoclavable
- ② **WA50082A** Video telescope "ENDOEYE 3D",
10 mm,
30°
▪ Autoclavable


Technical Data

Optical system	
Field of view	67° (16:9)
Direction of view WA50080A	0°
Direction of view WA50082A	30°
Depth of field	25 mm – 120 mm
Insertion portion	
Outer diameter	10 mm
Max. Working length WA50080A	330 mm
Max. Working length WA50082A	335 mm
Max. Insertion portion width WA50080A	10.07 mm
Max. Insertion portion width WA50082A	10.20 mm
Safety	
Protection against electric shock	Type BF applied part
IP classification	IPX7
Ambient conditions	
Operation temperature	10–35 °C
relative humidity	3–85 %
Storage temperature	10–40 °C
relative humidity	30–85 %
Transport temperature	–47–70 °C
relative humidity	10–95 %

Video System

- N5426950** "OTV-S300"
Video System Center "VISERA ELITE II"
All-in-One 2D/3D Processor and Light Source

For VISERA ELITE II, see page 489.

Instrument Tray

- WA05938A** Instrument tray,
wire basket,
for ENDOEYE video telescopes

3-3.9

For wire instrument baskets, see page 657.

3-4-1 **4K Multi-Specialty Video Processor with Upgradeable 3D/IR Capability**

All-in-one System Integrating 4K, 3D, and IR

VISERA ELITE III is an evolutions product which integrates 4K, 3D and IR in one system. The combination of those features in one system allows users to standardize for one imaging system across ORs, leading to increased operational efficiency.

Customized Preferences to Meet Clinical Needs and Budgets

The 3D and IR observation modes can be activated via a software upgrade available via Portable memory port. Users can easily select a customized configuration to meet the specific clinical needs and budget of their facilities.

Color Customization to Support a Range of Procedures

More accurate and more nuanced color reproduction can be achieved by making small adjustments to hue, saturation and color tone for Red, Orange, Yellow and Violet respectively. Image settings can be optimized to meet clinical conditions and physician preferences.

Specifications		
Power Supply	Rated voltage	110–240 V AC; within ±10%
	Rated frequency	50/60 Hz; within ±1 Hz
	Rated input	Less than 170 VA
Size	Dimensions (maximum)	390 (W) x 198 (H) x 500 (D) mm
	Weight	12.1 kg
Observation	Digital signal output	12G-SDI (SMPTE ST 2082), 3G-SDI (SMPTE424M), HD-SDI (SMPTE292M)
	Switching the recording output signal	This function allows you to switch the recording method of the video signal output to the recording device. · HD-SDI recording method · 3G-SDI recording method
	Electronic zoom	Electronic magnification of endoscopic images is available. 4K camera head: 6 steps (x1.0, x1.2, x1.4, x1.6, x1.8, x2.0) non 4K camera head: 3 steps (x1.0, x1.2, x1.5) 3.4.3
	Observation mode	The optical-digital observation can be performed. The endoscope compatible with the optical-digital observation is required. NBI observation mode, IR observation mode*, YE observation mode, 3D observation mode*
Documentation	Remote control	The following peripheral device can be controlled (specified models only). · Portable memory · Video recorder
	Recording format	Standard image quality: TIFF Low image quality: JPEG
Classification (Electromedical Equipment)	Degree of protection against electric shock	TYPE CF (the entire camera head)
	Degree of protection against explosion	The product should be kept away from flammable gases.

* IR/3D observation can be available by upgrading.

Innovation by
Sony & Olympus

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS MEDICAL SYSTEMS CORP.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
www.olympus-global.com

Global Professional Education website
www.olympusprofed.com

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
Unique Features

Multi-Input Recording

Directly toggle between two inputs and perform parallel or picture-in-picture recordings.

Multi-Storage/DLM

Recordings can be saved up to two locations in different recording resolutions. Data Life Management optimizes storage space.

MVR Manager

Browse and review all patient files. Zoom, take snapshots, cut video. Add notes, report in PDF, print and archive as you wish.

MVR Secure

Protect your files from hacking or power outages. Includes user authentication and historical logs for data protection.

Upgrade Packages

The MVR Pro Upgrade Packages can be installed at any time by simply using an assigned activation key provided by MediCapture.

DICOM¹

Easily retrieve patient DICOM worklist and store to your PACS.

4K Recording^{1,5}

Records delicate surgeries in stunning 4K (2160 p30/25).

MVR Remote App⁴

- 1 Upgrade package option
- 2 Subject to changes without notice
- 3 Please ask for a quotation for higher memory
- 4 Install the MVR Remote App on your Android smartphone or tablet
- 5 Check camera compatibility with MediCapture

MediCapture, Inc. – North America

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462, USA

+1 888 922 7887

www.medicapture.com
sales@medicapture.com

MediCapture Europe B.V.

Kraanpoort 10
6041 EG Roermond
The Netherlands

+31 475 214 305

www.medicapture.com
sales.europe@medicapture.com

MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,
Xinyi District 110, Taipei City,
Taiwan, Republic Of China

+886 2 2345 3668

www.medicapture.com
sales.asiapacific@medicapture.com

© 2019 MediCapture, Inc.
MVRPRO-DS01-EN-191218-MC

MVR PRO Specs²

Device Type	HD Medical Video Recorder 3.4.5
	Records High-Definition video and still images from medical video sources to: USB, internal storage, network storage. Exports video and images to PACS in DICOM format ¹
Case/Buttons	Stainless steel, fluid resistant
User Interface	Overlay on internal touchscreen monitor, or external monitor
Video Inputs	Multi-Input: 2x HDMI, DVI (with adapter), 4K-SDI with MVC Pro SDI to HDMI (optional) ⁵
Video Outputs	1x HDMI, DVI (with adapter)
Video Input Resolutions	4K on Input 1 only ⁶ : 3840x2160p30/25 Both inputs (1+2): 1920x1200p60, 1080p60/50/30/25, 1080i60/50, 1280x1024p60/50, 720p 60/50
Video Recording Resolutions	3.4.6 3840x2160 ¹ , 1920x1200, 1920x1080, 1280x1024, 1280x720. For 3840x2160 recording, "4K recording" needs to be activated ^{1,5}
Video Recording Bitrates	MP4 / coding H.264 up to 32 Mbit/s (1920x1080 p60), lower resolutions have lower bitrates MP4 / coding HEVC (H.265) reduces bitrate by 50%.
Image Format	JPG, PNG
Sound Recording	Bluetooth Mic, Mic by Stereo Headset (4 poles), USB Mic / Audio embedded in video
Internal Storage	128GB. Extension on request ³
External Storage	USB flash drive, USB hard drive (FAT32, NTFS, exFAT)
USB Support	3.4.5 USB3.0, SuperSpeed, 1 front connector, 2 rear connectors
Network Storage	External Server, 1000/100/10Mbit/s
Footswitch/Trigger	3,5mm mini jack. Compatible with standard footswitches (not included) and closed-contact camera triggers, Bluetooth footswitches
Remote Control	MVR Remote App ⁴ , Administration access via Network (Settings, Samba), IR remote control (included)
Mechanical (touchscreen closed)	Size: 173 x 153 x 50mm / 6,7 x 6,0 x 2" Weight: 1,20 kg / 2lb 13oz
Environmental	Storage: -40°C to +85°C Operating: -20°C to +40°C
Power	AC 100-240V, 50-60Hz, 0.4-0.2A, 3 prong AC jack
Certifications	CE 93/42/EEC ANSI/AAMI ES60601-1 (2005(R)2012 + A:2012, C1:2009(R)2012 + A2:2010(R)2012) (EN 60601-1 ed 3.1) CAN/CSA C22.2 No.60601-1-1:2014-Edition 3 FCC Part 15, Class 1 Equipment according to US-FDA, RoHS3 / ISO 13485: 2016
Order Number	MVR409 MVR Pro recorder, USB flash drive, HDMI cable, Wi-Fi / BT antenna, HDMI to DVI adapter, Split Remote Cable Adapter, Power Cable, IR Remote Control, Rubber Feet (4), User Guide (English)
Order information for optional Upgrade Packages:	AKDICOM : DICOM Upgrade Package AK4K : 4K Recording Upgrade Package

Medical Imaging Made Easy™

1.7 Enable Internal Storage

- Go to: ► **Settings** ► **Archive** ► **Storage**

Internal storage needs to be enabled for PACS export. The recorder will store the study to the internal storage and will export the study later to the PACS.

- Set internal storage to the maximum video resolution you want to allow for PACS export.

Note: If you do not want to allow videos to be sent to PACS, you can disable video upload in PACS Server settings → Only Images
See section 3.3

3.4.6
(2160p → 4K
rezolucija)

1.8 Alternative storages

Besides the internal storage for later PACS export, you can enable a second option for direct storage.

-USB storage


or


-Network storage to a network server.

Note: MVR / MVR Pro stores files onto this second option in standard format "MP4", "JPG", "PNG".


a) USB storage

USB storage is recommended for backup storage, for presentations or to hand-out a copy to the patient.

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Full Color Infrared LED Light Source 3.5.1

Full Color in IR Mode for Better Decision Making

"Fluorescence Imaging Overlay Mode" provides real time 4K fluorescence imaging during white light observation.

Lower Operating Costs with Long-lasting LEDs

LED lights are used in place of a conventional Xenon light source, reducing the frequency of lamp replacement and lowering operational costs. In addition to LEDs for white light observation, an LED dedicated to IR observation is incorporated to facilitate greater clinical versatility.

Specifications		
Power Input	Rated voltage	100-240 V AC; within ±10%
	Rated frequency	50/60 Hz; within ±1 Hz
	Rated input	215 VA
Size	Dimensions (maximum)	390 (W) x 162 (H) x 502 (D) mm
	Weight	13.5 kg
Illumination	Providing maximum light intensity	Less than 5.64 W
	Cooling	Forced-air cooling
	Observation mode	Normal light observation mode, NBI observation mode, IR observation mode*, YE observation mode*, 3D observation mode* * It can be used by performing an upgrade to add the corresponding observation function.
Automatic Brightness Adjustment 3.5.3	Brightness adjustment method	LED drive current control
	Exposure	17 steps
Classification (Electromedical Equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock of applied part	CF type
	Degree of protection against explosion	The light source should be kept away from flammable gases.

3.5.2

Innovation by
Sony & Olympus

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS MEDICAL SYSTEMS CORP.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
www.olympus-global.com

Global Professional Education website
www.olympusprofed.com

N8602406-062022

ULTRA Telescopes


3.6-2

VISERA 4K UHD

Surgery

10 mm IR Laparoscopes


- ① **WAIR100A** Telescope "IR", 10 mm, 0°
- ② **WAIR130A** Telescope "IR", 10 mm, 30°

Please check availability and legal situation of medical dye, which is necessary to perform fluorescence imaging. 3.6.5

Technical Data

Direction of view	
WAIR100A	0°
WAIR130A	30°
Field of view	
WAIR100A, WAIR130A	88°
Working length	
WAIR100A	316.4 mm
WAIR130A	318.6 mm
Maximum insertion portion width	
WAIR100A, WAIR130A	10.2 mm

10 mm IR Telescopes Short


- WA4KT130** Telescope "ULTRA", IR, 10 mm, 30°, short

Technical Data

Direction of view	
WA4KT130	30°
Field of view	
WA4KT130	88°
Working length	
WA4KT130	219.0 mm
Maximum insertion portion width	
WA4KT130	10.2 mm

5 mm IR Telescopes


- WAIR500A** Telescope "ULTRA", IR, 5.4 mm, 0°
- WAIR530A** Telescope "ULTRA", IR, 5.4 mm, 30°

Technical Data

Direction of view	
WAIR500A	0°
WAIR530A	30°
Field of view	
WAIR500A, WAIR530A	84°
Working length	
WAIR500A	315.8 mm
WAIR530A	317.0 mm
Maximum insertion portion width	
WAIR500A, WAIR530A	5.45 mm

Video System

- ① **N5426850** "OTV-S200"
Video System Center "VISERA ELITE II"
All-in-One 2D Processor and Light Source
- N5426950** "OTV-S300"
Video System Center "VISERA ELITE II"
All-in-One 2D/3D Processor and Light Source
- ② **N5427050** Light Source "CLV-S200-IR"

For VISERA ELITE II, see page 489.


Camera Head

- N5766550** CH-S200-XZ-EA Camera head autoclavable
- N5409160** CH-S200-XZ-EB 3CMOS Camera Head


Instrument Tray

WA05970A Instrument tray,
with lid,
outer dimensions: 537 x 139 x 268 mm

Multi-Purpose Insert Trays

① **WA05971A** Insert tray,
for the upper part of WA05970A,
outer dimensions:
482 x 60 x 224 mm


② **WA05973A** Insert tray,
for the bottom part of WA05970A,
outer dimensions:
460 x 39 x 207 mm

Lid for Insert Trays

WA05974A Lid,
for all insert trays,
except for WA05973A


Insert Tray for Telescopes

3.6.7

WA05994A Insert tray,
for 2 telescopes, VISERA camera head, and
video adapter

Insert Trays for Surgery/Laparoscopy

WA05972A Insert tray,
for 5 HIQ+ shafts and jaws inserts,
maximum length of instruments: 330 mm,
outer dimensions: 480 x 47 x 224 mm

27

- Eco-friendly in terms of RoHS compliance
- CF applied part classification for use in cardiac surgery
- 3.4.1 ■ Optimal light transmission thanks to fused optical fibers
- Durability: Optimized tube and fused fiber bundles for extremely high durability against mechanical stress and autoclaving
- High temperature resistance
- Compatible with a wide variety of endoscopes if used with suitable adapters

Light-Guide Cables

Delivery of all light-guide cables includes light-guide adapter WA03320A. For light-guide adapters, see sheet ACC-150.

WA03300A Light-guide cable, 2.8 mm, 3 m, CF type

WA03310A Light-guide cable, 4.25 mm, 3.7.2 3 m, CF type

PDD Light-Guide Cable

A93200A Light-guide cable, plug type, for fluorescence diagnosis, 3.0 mm, fluid, not autoclavable

Colposcope Light-Guide Cable

N1008800 Light-guide cable "OCS5-LG", for colposcope OCS-500

Adapting Olympus Light-Guide Cables to Light Sources

WA03320A Spare light-guide adapter, Olympus light-guide cable to Olympus light source

A3199 Light-guide adapter, Olympus light-guide cable to Winter & Ibe light source


with screw-type connection

A3201 Light-guide adapter, Olympus light-guide cable to Storz light source


A3202 Light-guide adapter, Olympus light-guide cable to Wolf light source

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





The WM-NP3 is supplied with a separation transformer based on one of the following voltage options

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