


# Cardio7e



## Streamlining workflow with Bionet ECG



EASY & Simplifying WORKFLOW with 7" Touchscreen & AUTO button operation



Advanced analysis & Interpretation



Effective data management with 5min. of Full disclosure



Excellent connectivity LAN/Wifi network bidirectional connection to PACS/EMR/HIS

4.2. Maksimali užrašomos ritmo kardiogramos trukmė 5 min.

## EASY Operation

- With 7" Bright Color LCD Touchscreen operation
- Alphanumeric keyboard on screen
- Quick patient ID entry with barcode scanner
- Convenient A4 print-out with 12CH interpretation on/off
- Lead fault detection & notice / Pace-maker detection

## Simplifying workflow

- Easy patient ID entry & worklist orders retrieval
- One touch "AUTO" workflow with acquisition, save, print-out, and transfer to HIS
- Direct DICOM file transfer to PACS
- Convenient various file conversions in PDF, XML, MFER, BMP, JPEG

## Efficient data management & Full disclosure

- 5min. of Full disclosure for Arrhythmia tracing & review
- 200 ECGs internal data memory
- Long-term recording 1min, 3min, 5min, and 10min for HRV analysis
- 10 sec. pre-acquisition for immediate printing of any instant abnormal rhythm at sight

## Advanced Analysis & Report

- VECTOR, GUIDE, TEXT report
- Patented Bionet ST-MAP for enhanced interpretation
- 120 kinds of interpretation based on the Minnesota Code

## Excellent connectivity


- LAN, Wifi network connection to PACS/EMR/HIS
- USB : Barcode scanner, external memory, keyboard

11. 3. Aparatas turi būti pajėgus medicininio EKG tyrimo ataskaitą pdf formatu persiusti į medicininių tyrimų duomenų archyvavimo sistemą MedDream PACS  
4. Aparatas turi būti pajėgus bevieliu tinklu gauti pacientų sąrašą (DICOM Worklist) iš medicininių tyrimų duomenų archyvavimo sistemos MedDream PACS.

## Simple AUTO workflow



## Advanced analysis, Report, & Data Management



### Long-term recording

1, 3, 5, 10 min.

### pre-acquisition

10 sec' real time print-out

### User-friendly UI


7" LCD Touch screen



9. Galimybė kardiogramą peržiūrėti aparato ekrane

### Preview screen

Data preview before print-out




One touch  
AUTO Button

AUTO

### Enhanced data capacity

5 min.  
Full disclosure



200ECGs  
Internal memory

No	ID	Name	Age	Date	Study
1	ujj	ryhj	0	20150708	ECG
2	1234		0		ECG
3	234		0		ECG
4	cdee		0	20150512	MVV
5	ddrrg	rdvgy	0	20150512	FVC
6	ed		0	20150512	FVC
7	gijn		0	20150512	ECG
8	ujj		3	20150512	FVC
9	wrtu		0	20150512	FVC
10	1234		0		ECG
11	234		0		ECG
12	cdee		0	20150512	MVV
13	ddrrg	rdvgy	0	20150512	FVC
14	ed		0	20150512	FVC
15	gijn		0	20150512	ECG
16	ujj		3	20150512	FVC
17	wrtu		0	20150512	FVC
18	1234		0		ECG
19	234		0		ECG

Easy data search  
ID, Name, Age, etc

No	ID	Name	Age	Date	Study
1	11270410	phudg	0	20151127	ECG
2	11275510	ggafb	0	20151127	ECG
3	11275510	ggafb	0	20151127	ECG
4	7890	dvvv	32	20151127	FVC

Patient Search

Search Condition: ID

-ID:

Name:

Day:

Age:

YMMO:

Find    Close    All

# Technical specifications

ECG Leads	Simultaneous 12 Leads Resting ECG	8. Vartotojo pasirenkamos derivacijos
Recording Channel	3CH+1RHY, 3CH+3RHY, 6CH+1RHY, 12CH, 6CH+ST map 1CH Long Time(1min, 3min, 5min, 10min) and Special Beat Report (Text, Guide, Vector, ST map)	7. Keičiamas elektrokardiografo spausdintuvo jautrumas (priklausomai nuo aparato darbo režimo)
Gain	2.5, 5, 10, 20, Auto (I-aVF: 10, V1-V6: 5) mm/mV	
Printing speed	5, 12.5, 25, 50, 100 mm/sec	5. Keičiamas EKG spausdinimo (popieriaus traukimo) greitis
Sampling rate	Analysis Sampling Rate - 500 Digital Sampling Rate - 8,000	
Filters	AC (50/60 Hz, -20dB or better), Muscle (25-35Hz, -3dB or better), Base line drift (0.05Hz, 0.1Hz, 0.2Hz, -3dB or better), Low pass filter(off, 40Hz, 100Hz, 150Hz)	
Display	800 x 480 pixel resolution, 7" Color TFT wide display 12 channels preview, ID, HR, Gain, Speed, Power status	2. Aparato ekranas
User Interface	Touch screen, Rotary push-knob	
Printer Resolution (Speed 25 mm/s)	Thermal print head, Roll Paper Report paper - width : A4/Letter 215mm(8.5") - length : A4 297 mm(11.7") Letter 279mm(11") Resolution : Vertical 8dot/mm(0.125mm pitch), Horizontal 16dot/mm.	6. Spausdinimo raiška, esant 25 mm/s greičiui
Patient data	ID, name, age, gender, height, weight, race, smoke, department, room number, study desc., accession no., referring physician	
Basic Measurement	Heart rate(30-300bpm, ±3bpm), PR/RR int., QRS dur., QT/QTc int., P-R-T axis, SV1/RV5/R+S amp.	
ECG Data Storage	Internal Storage for 200 ECGs : Built-in memory	
Electrical	Internal noise: 20µV(p-p)max Input impedance: ≥ 50MΩ Input voltage range: ±5mV CMRR: > 105dB DC offset voltage: ≥ ±400mV Patient leakage current: < 10µA Frequency response: 0.05 ~ 200 with in -3dB Isolated, defibrillation and ESU protected	12. Aparato vidinė atmintis
Signal Quality Control	Pacemaker pulse detection, Lead fault detection, Signal Saturation Detection	
Line Power	Input: 100-240Vac, 1.5-0.75A, 50-60Hz Output: 15Vdc, 4.2A	14. Aparato naudojami maitinimo šaltiniai
Battery type	Replaceable and Rechargeable, Lithium ion, 10.8v, 3250mA	
Battery Capacity	6 hours of normal use or print 200 ECG pages (12 channel format at 25mm/s and 10mm/mV) Battery recharge to full capacity in 3hours.	15. Reikalavimai aparato akumuliatoriui
Communication	LAN, WIFI(option), USB flash driver, USB barcode scanner	
Safety Conformity	Class I, Type CF defibrillation-proof applied part.	10.3. Beviolio ryšio WLAN (WIFI) sąsaja;
Environmental Operation	Ambient Temperature: 10-40°C Relative Humidity: 30-85%RH Atmospheric Pressure: 70-106KPa	13. Apsauga nuo defibriliacijos išskrovos 10.2. Potencialų išlyginimo jungtis;
Environmental Storage/Ship	Ambient Temperature: -20-60°C Relative Humidity: 10-95%RH Atmospheric Pressure: 50-106KPa	
Dimensions	300(W) x 290(H) x 89.5(D) mm, Approx. 3.5 kg	
Packing dim. & weight	415(W) X 400(H) X 215(D)mm & approx. ≈6kg	
Standard Accessory	Patient cable(1 EA), Limbs electrodes(1 SET), Chest electrodes(1 SET), ECG chart paper(1 EA), Power cord (1 EA), Operation manual & ECG diagnosis guide (1EA), Battery (1EA)	16. Komplektacija:
Options	Rechargeable battery(1 EA), Cart(1 EA), Hanger(1 EA), ECG spray (1 EA)	

\* Specifications can be changed without prior notification.



## Headquarters Bionet Co., Ltd.

5F, Dream Markone 61 Digital-ro 3gil,  
Guro-gu, Seoul, KOREA 08375  
TEL +82 2 6300 6410  
E-Mail. sales@ebionet.com

10. 1.EKG kabelio jungtis;

1. Aparato vežimėlis

Medical, ECG

## Cardio7e

Captures, stores, analyzes, and transmits by **simplifying workflow with user-friendly usability**

### Features

#### Simplifying workflow with "AUTO" operation

- One button "AUTO" operation from acquisition to print-out & data record in 3 sec.

#### Easy ECG file data conversion

- Easy file conversion to PDF, XML, JPG, BMP, MFER, DICOM that transfers with easy user interface.

#### Enhanced sampling rate for data accuracy

- With **8,000Hz per second sampling rate**, it is designed to improve clinical decision for cardiology with accurate ECG interpretation algorithm and **full storage of 200 ECG data per 5 minutes**.

#### Advanced algorithm & analysis


- With **long term recordings of 3,5,10 min**, it features various analysis algorithms of ECG diagnosis.

#### BMS Plus Patient data management in mind


- All **patient data is transferable to BMS Plus** and PC based data management SW with real time remote operation.

#### BMS plus offers caregivers the workflow flexibility for

- Data preview
- Device remote control
- File conversion to PDF, XML, JPG, BMP, MFER, **DICOM** for data transfer **through PACS**.
- Direct print-out, magnifying and comparing ECG trace, measuring ECG complex, adding comments




4.1. Ramybės kardiograma užrašoma pasirinktinai rankiniu arba automatinu režimu;




# Cardio7 Operation Manual

■ Front View



① Printer Cover Switch

■ Rear View



10. 4. 2 USB jungtis.

- ① Protective Ground Terminal
- ② Power Switch
- ③ AC Power Connection Port ( Fuse : 250V , 5A ( High breaking type ) )
- ④ USB Port
- ⑤ Spirometry port / RS232 protocol
- ⑥ RJ45 LAN Port

3) After pulling out the paper, remove the jammed paper by cutting.

- When printing more than 10 minutes with the RHYTHM key, the paper may get blackened.

#### Warning

Make sure to use only the ECG record paper provided by our company. If not, it might negatively influence on the diagnosis when confirmed by a doctor due to poor printing status. Our company is not responsible for all the issues occurring by ignoring this requirement.

#### WARNING

No modification of this equipment is allowed.

Do not modify this equipment without authorization of the manufacturer.

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of equipment.

## Network Connection

-No one other than service personnel can connect this device to the network.

Consult with hospital IT staff in advance. It is recommended to follow IEC 80001-1

(Risk management of IT networks connected to medical devices)

### LAN Network

LAN networks are usually configured through a star topology. Individual devices can be grouped together via a layer-n-switch. Other data traffic is separated by a separate VLAN network. Configure the device's network settings according to this user guide and network specifications.

LAN connection specifications are described in the following standard specifications.

- Wired network: IEEE 802.3
- **Wireless network: IEEE 802.11 (a, b, g, n)**

If the device is used as a layer-2-switch or layer-3-switch, the port settings must be configured

11. 1. Aparatē turi būtī  
īdiegta bevielio ryšio  
sąsaja, atitinkanti  
IEEE 802.11n  
standartą

# Cardio7 Operation Manual

When clicking the 'Format' button, a line window will show up and it will disappear when selecting EKG/FVC/SVC/MVV, MFER(ECG), XML(ECG), JPEG, PDF. 'Format' menu is used to select the user's options when transferring the data.

- EKG/FVC/SVC/MVV: To save the data by type of '\*\*\*.ekg' / '\*\*\*.fvc' / '\*\*\*.svc' / '\*\*\*.mvv'
- MFER: To save the data by type of '\*\*\*.mwf'
- XML: To save the data by type of '\*\*\*.xml'
- JPEG: To save the data by type of '\*\*\*.jpg'
- PDF: To save the data by type of '\*\*\*.pdf'
- DICOM : To save the data by type of '\*\*\*.dcm'

11. 2. Aparatas turi būti pajėgus bevieliu tinklu, naudojant DICOM ECG (modality) standartą, perduoti medicininių tyrimų duomenis į medicininių tyrimų duomenų archyvavimo sistemą MedDream PACS

When clicking the 'Delete' button, a line window will show up and it will disappear when selecting Yes or No. 'Delete' menu is used to select the user's options after the data has transferred.

- Yes: To delete the data after it has transferred
- No: To not delete the data after it has transferred

## NOTE

- If 'Format' menu has only 'EKG/FVC/SVC/MVV' options, it is possible to select 'All' or 'All New' in 'Method' menu; for other formats, 'Selected' is the only option.

## 5) Switching menu

When selecting the 'Main' button at the main page of 'File management', a line window will show up. You can choose ECG, Spiro or Worklist to go to the chosen page.

The patient's information of selected data does not transmit when skipping to 'Worklist Management' screen. But it depends on the setting value of 'PAT Info.' menu when skipping to 'ECG' main screen.

- ECG : Skip to ECG main screen
- Worklist : Skip to 'Worklist Management' screen
- Spiro : Skip to spirometry measurement screen