PHILIPS

Attention - Before using this device, Consult Instructions for Use

PHOENIX ATHERECTOMY SYSTEM

INSTRUCTIONS FOR USE

PHOENIX ATHERECTOMY SYSTEM

INSTRUCTIONS FOR USE

Table of Contents

Pg	Language		Instruction For Use
3	English	English	INSTRUCTIONS FOR USE
11	German	Deutsch	GEBRAUCHSANWEISUNG
19	Spanish	Español	INSTRUCCIONES DE USO
27	French	Français	MODE D'EMPLOI
35	Italian	Italiano	ISTRUZIONI PER L'USO
43	Dutch	Nederlands	GEBRUIKSAANWIJZING
51	Polish	Polski	SPOSÓB UŻYCIA
59	Turkish	Türkçe	KULLANMA TALİMATI
67	Bulgarian	Български	ИНСТРУКЦИИ ЗА УПОТРЕБА
75	Portuguese Brazil	Português BR	INSTRUÇÕES DE USO
83	Czech	Čeština	NÁVOD K POUŽITÍ
91	Danish	Dansk	BRUGSANVISNING
99	Estonian	Eesti	KASUTUSJUHEND
107	Finnish	Suomi	KÄYTTÖOHJEET
115	Greek	Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
124	Croatian	Hrvatski	UPUTE ZA UPORABU
132	Hungarian	Magyar	HASZNÁLATI ÚTMUTATÓ
140	Lithuanian	Lietuviškai	NAUDOJIMO INSTRUKCIJOS
148	Latvian	Latviski	LIETOŠANAS INSTRUKCIJAS
156	Norwegian	Norsk	BRUKSANVISNING
164	Portuguese	Português	INSTRUÇÕES DE UTILIZAÇÃO
172	Romanian	Română	INSTRUCȚIUNI DE UTILIZARE
180	Russian	Русский	ИНСТРУКЦИИ ПО ЭКСПЛУАТАЦИИ
188	Swedish	Svenska	BRUKSANVISNING
196	Serbian	Srpski	UPUTSTVO ZA UPOTREBU
204	Simplified Chinese	简体中文	使用说明
212	Korean	한국어	사용 지침

Attention - Before using this device, Consult Instructions for Use

PHOENIX ATHERECTOMY SYSTEM

INSTRUCTIONS FOR USE

DEVICE NAMES

Phoenix Atherectomy Catheter

Phoenix Atherectomy Handle

Phoenix Wire Support Clip

CAUTION: FEDERAL (US) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON ORDER OF A PHYSICIAN.

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE, NOTING ALL WARNINGS AND PRECAUTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

INDICATIONS

The Phoenix Atherectomy System is intended for use in atherectomy of the peripheral vasculature. The system is not intended for use in the coronary, carotid, iliac or renal vasculature.

CONTRAINDICATIONS

The use of the Phoenix Atherectomy System is contraindicated where introduction of any catheter would constitute an unacceptable risk to the patient. Contraindications include but are not limited to:

- Recent history of myocardial infarction or thrombotic or hemorrhagic stroke
- Any evidence or history of intracranial bleeding or aneurysm
- Known hypercoagulable state or coagulopathy
- Bleeding diathesis, platelet dysfunction, thrombocytopenia, or thrombocytosis
- Patients in whom antiplatelet, anticoagulant, or thrombolytic therapy is contraindicated
- Allergy or intolerance to contrast agents or medications used to perform endovascular procedures
- Do NOT operate the Phoenix Atherectomy System in vessels smaller than the indicated size as perforation, dissection or injury may occur.

SYSTEM DESCRIPTION

The Phoenix Atherectomy System has two main components: the Phoenix Atherectomy Catheter and the Phoenix Atherectomy Handle. The Phoenix Wire Support Clip is an accessory to the Handle to facilitate guidewire management during the procedure. All system components are sterile, single-use devices designed for atherectomy of the peripheral vasculature.

The Phoenix Atherectomy Catheter is a long, flexible double-lumen catheter that contains a torque shaft that is attached to a metal cutting element at the distal Catheter tip. When activated, the torque shaft rotates, causing the cutting element to excise diseased plaque from the artery. The excised diseased material is brought into the Catheter where it is mechanically conveyed from the patient using an Archimedes Screw that is fixed on the outer surface of the torque shaft. The central lumen accommodates a standard 0.014" guidewire. The distal section of the Catheter (including the cutter head at the Catheter tip) is radiopaque for fluoroscopic visualization of the Catheter position and tip.

orientation. The proximal end of the Catheter is composed of an adapter with two ports; a side port for debris conveyance and a central port for guidewire introduction. The Catheter installs into the Handle by snapping the proximal end of the Catheter into the Handle.

The Phoenix Handle includes a battery-operated motor for driving the cutter of the Phoenix Atherectomy Catheter at a nominal speed of 10,000 to 12,000 rpm. The Phoenix System is activated by sliding the ON/OFF slider switch on the top of the Handle.

The Phoenix Wire Support Clip is an accessory with an integrated guidewire torquer, which can be used to help support and secure the guidewire in position during the procedure.

See Figure 1 for the assembled Phoenix Atherectomy System with the Phoenix Atherectomy Catheter, Handle and Wire Support Clip.

Table 1. Phoenix Atherectomy Catheter Key Specifications

Catheter Tip Diameter	Minimum Introducer Size	Crossing Profile	Working Length	Maximum Guide Wire Diameter	Minimum Vessel Diameter ¹	Anatomical Locations
1.5 mm	4Fr (1.5 mm) or larger	1.5 mm	149 cm	0.014" (0.36 mm)	2.0 mm	Femoral, popliteal, or distal arteries located below the knee
1.8 mm	5F (1.8 mm) or larger	1.8 mm	130 cm	0.014" (0.36 mm)	2.5 mm	
2.2 mm	6F (2.2 mm) or larger	2.2 mm	130 cm	0.014" (0.36 mm)	3.0 mm	
1.8 mm	5F (1.8 mm) or larger	1.8 mm	149 cm	0.014" (0.36 mm)	2.5 mm	
2.2 mm	6F (2.2 mm) or larger	2.2 mm	149 cm	0.014" (0.36 mm)	3.0 mm	
2.4 mm	7F (2.4mm) or larger	2.4 mm	130 cm	0.014" (0.36 mm)	3.0 mm	

¹**Warning: Do not use the Phoenix Atherectomy Catheter in vessels smaller than the indicated size or harm to patient (vessel perforation, dissection or injury) could occur.**


Figure 1: Phoenix Atherectomy System: Catheter distal tip (left) and Handle (right), shown with the Catheter inserted into Handle and Wire Support Clip attached.

WARNINGS

- The Phoenix Atherectomy System is intended for one-time use only. Do NOT resterilize and/or reuse. Resterilization or reuse may potentially compromise device performance and safety and may increase the risk of infection.
- Use the Phoenix Atherectomy System products prior to the “Use By” (“Expires”) date specified on each individual package label.
- Do NOT operate the Phoenix Atherectomy System if the guidewire is subintimal as perforation, dissection or injury may occur.
- Do NOT operate the Phoenix Atherectomy System in vessels smaller than the indicated size as perforation, dissection or injury may occur.
- When the Phoenix Atherectomy Catheter is exposed to the vascular system, it should not be advanced or retracted except under direct fluoroscopic observation. If resistance is met during advancement or retraction, determine the cause of the resistance before continuing.
- Do NOT turn ON the Phoenix Atherectomy System except under fluoroscopy.
- Do NOT leave the distal tip of Catheter in a stationary position with System turned ON, or vessel perforation, dissection or injury may occur.
- Do NOT turn ON the Phoenix Atherectomy System unless the Phoenix Atherectomy Catheter is tracked over a guidewire and the distal tip of the guidewire is placed in an appropriate section of the vessel.
- All air must be removed from all catheter lumens by priming the Catheter as described in the Instructions for Procedure section prior to inserting the Catheter into the patient.

- Do not use, or attempt to correct, a Phoenix Atherectomy Catheter if it is bent or kinked or has any evidence of damage as this may result in breakage or compromised performance.
- The Phoenix Catheter only may be used with the Phoenix Handle.
- The Phoenix Atherectomy System is not to be used in the presence of combustible or flammable gases, anesthetics, or cleaners/disinfectants.

PRECAUTIONS

- Prior to use, all Phoenix Atherectomy System components should be examined to verify functionality.
- For use only by physicians trained in the use of the Phoenix Atherectomy System and percutaneous peripheral interventional procedures.
- Test the Phoenix Catheter prior to insertion into the patient as described in the Instructions for Procedure section.
- Close the hemostatic valve tight enough to prevent blood leakage around the Catheter shaft, but still allow axial movement of the Catheter through the valve. Special care should be taken when adjusting the valve tightness if a Tuohy-Borst valve is used.
- Take care not to immerse the Handle because damage to the electrical components of the Handle could occur.
- As in all atherectomy procedures, the duration of the procedure should be determined by the physician based on patient condition.
- The maximum recommended length of run time for the Phoenix Handle during a procedure is 20 minutes.
- Use only listed compatible guidewires with the Phoenix Atherectomy System. Use of any guidewires not listed as compatible may compromise performance or damage the Phoenix Atherectomy System.
- Use of PTFE-coated guidewires is not recommended, unless the guidewire is specifically indicated for atherectomy.
- Verify that the proximal end of the guidewire is secured within the torquer when the Handle is ON. Failure to do so may result in guidewire rotation and/or loss of position within the vessel.
- Ensure the Catheter tracks smoothly and easily over guidewire during Handle operation. If the Catheter does not track easily over the guidewire, immediately turn OFF the System and replace the System.
- The Phoenix System must be OFF while advancing forward through the Introducer Sheath.
- Carefully monitor the Catheter tip and guidewire under fluoroscopy during Handle operation. If rotation or “whipping” of the guidewire is suspected, immediately turn OFF the system and investigate. Do not turn ON the System unless the Catheter is verified to track easily over the guidewire.
- Monitor flow of excised material into the disposal reservoir or drainage tube during operation of the Catheter. If flow of excised material ceases during the procedure, this is a sign that the Catheter drive system (cutter, torque shaft, or Handle motor drive) may not be operating properly.
 - Pull the Catheter back or proximal to the lesion to re-establish blood flow.
 - If flow of excised material is not re-established, turn OFF the Handle. Remove the Catheter from the patient. Otherwise embolization or an inability to remove excised plaque could occur.
- Hold the guidewire firmly with a torque device during the Catheter removal process either by hand or with the Phoenix Wire Support Clip if it is used. Failure to do so may result in guidewire rotation and/or loss of position within the vessel.

ADVERSE EFFECTS

Potential adverse effects associated with the use of this device and other interventional catheters include but are not limited to the following (alphabetical order):

- Access site injury, including pain
- Arterial dissection
- Arterial occlusion
- Arterial perforation or pseudoaneurysm
- Arterial restenosis
- Arterial rupture
- Arterial spasm or abrupt or sub-acute closure
- Bleeding complications which may require transfusion
- Death
- Embolism, including thrombus, plaque, air, device, etc.
- Emergency or non-emergency arterial bypass surgery
- Fracture of any component of the device that may or may not lead to serious injury or surgical intervention
- Hematoma
- Hypotension
- Infection or fever
- Ischemia
- Reaction to contrast media, procedure medications or Catheter materials, including allergic reaction
- Vascular injury which may require surgical repair

The occurrence of these adverse effects may lead to the need for repeat catheterization/angioplasty, emergency bypass surgery, or death.

HOW SUPPLIED, STERILIZATION, AND EXPIRATION

The Phoenix Atherectomy Catheter, Handle and Wire Support Clip are packaged and sterilized individually. The Phoenix Atherectomy Catheter, Handle, and Wire Support Clip are sterilized using ethylene oxide gas. Each device is sterile if the sealed pouch is unopened and undamaged. All are intended for single use only and should not be reused or re-sterilized. Use before the “Use by” dated printed on the package labels. Store products indoors at room temperature in their original packaging away from sunlight and keep dry.

SUPPLIES REQUIRED FOR THE PHOENIX ATHERECTOMY PROCEDURE

NOTE: The Phoenix Atherectomy System is protected against electrical shock (defibrillation proof- type CF)

Item description (Single use items only - Do not re-sterilize or reuse):

1. A Phoenix Atherectomy Catheter and Handle
2. Disposal bag (250 ml or larger) or drainage tube for material conveyed out of the Catheter's disposal outlet
3. Appropriately sized introducer sheath with cross-cut valve per Key Catheter Specifications Table
4. Sterile heparinized (10,000 IU/L) 0.9% normal saline for priming the Catheter prior to use
5. 10cc or larger slip-tip syringe for priming the Catheter prior to use
6. Hemostasis valve
7. Compatible 0.014" sized, exchange length guidewire (260 cm minimum):
 - Abbott Vascular Hi-Torque S'PORT Wire (300 cm)
 - Ev3 Nitrex Guidewire (300 cm)
 - Abbott Hi-Torque Iron Man Wire (300 cm)
 - Abbott Hi-Torque Floppy II Wire (300 cm)
 - CSI ViperWire Guide Wire (300 cm)
 - Abbott Hi-Torque Standard Guide Wire (300 cm)
 - Abbott Hi-Torque Cross-it Guide Wire (300 cm)
 - Asahi Astat XS20 Guide Wire (300 cm)
 - Asahi Regalia XS1 Guide Wire (300 cm)
 - Phoenix Guidewire (300 cm)
 - USE GS14-300A Guidewire (300cm)
8. Torque device
9. Phoenix Wire Support Clip

INSTRUCTIONS FOR PROCEDURE

INSPECTION OF PHOENIX ATHERECTOMY SYSTEM COMPONENTS

Prior to use, all equipment to be used for the procedure should be examined carefully for defects. Examine the packaging for cuts, tears, or other breach of the sterile barrier. Examine the Phoenix Atherectomy Catheter for bends, kinks or other damage. Examine the Handle for sign of damage. Turn ON the Handle and ensure that it can be activated. Do not use any defective equipment.

PREPARATION OF THE PHOENIX ATHERECTOMY CATHETER

Complete the following steps using aseptic technique to prepare the Catheter for use:

1. Withdraw the Catheter from the sterile protective packaging and set the Catheter on the sterile table.
2. Priming the Phoenix Catheter: Attach a 10cc syringe filled with heparinized saline onto the disposal outlet and flush slowly with gentle pressure until saline drips out of the guidewire port. Cover the central guidewire port with a finger and continue to flush until saline drips out of the distal tip of the Catheter.
3. Handle Test: Briefly turn the handle on/off switch on (then off) to confirm that the Handle gear spins. Do not use any Handle found to be broken or non-functional during this preparation.
4. Insert the proximal end of the Catheter into the Handle. Ensure the proximal end of the Catheter is snapped into the Handle completely. Do not use any Catheter found to be broken or non-functional during this preparation.
5. Attach the disposal bag or drainage tube to the disposal outlet.

INSERTION AND OPERATION

I. Insertion

1. Refer to the Key Catheter Specifications Table for the appropriate introducer sheath size requirements. Insert the appropriately-sized sheath with a cross-cut hemostasis valve using standard techniques.
2. Advance a 0.014" guidewire through the sheath beyond the lesion to be treated, taking care to remain intraluminal.
3. Backload the end of the guidewire into the distal tip and out of the proximal end of the Phoenix Atherectomy Catheter guidewire lumen. Insert the distal tip of the Phoenix Atherectomy Catheter into the introducer sheath with the Phoenix System OFF until the tip exits the Introducer Sheath.

II. Guidewire Preparation

1. While holding the guidewire stationary and using fluoroscopic guidance, advance the Phoenix Atherectomy Catheter distal tip over the guidewire to within a few millimeters proximal to the target lesion.
2. Adjust the guidewire position relative to the lesion as needed. Confirm that the guidewire is intraluminal.
3. Confirm that the distal tip of the guidewire is positioned a minimum of 20 cm from the distal tip of the Catheter.
4. Apply torque device tightly to the guidewire approximately 20 cm to 30 cm from central guidewire port.
5. Loop the guidewire back through 180°, creating a backwards 'C' shape and, feed the guidewire into the integrated torquer on the Wire Support Clip to prevent guidewire rotation. This configuration provides wire support and allows the Catheter to track more easily over the guidewire (see Figure 1 for resulting guidewire loop).
6. Use of Wire Support Clip to prevent guidewire rotation:
 - a. Snap the Wire Support Clip onto the Handle as it appears in Figure 1.
 - b. Feed the guidewire through the integrated torquer and form a support loop of about 10-15 cm that spans from the central guidewire exit port (of the Catheter) to the Wire Support Clip.
 - c. Tighten the torquer onto the guidewire.
 - d. Verify that the guidewire does not spin in the Wire Support Clip.

Warning: Do Not Operate the Phoenix Atherectomy System if the guidewire is subintimal as perforation or dissection of the vessel may occur. A second angiographic viewing angle prior to device insertion to confirm wire placement is recommended.

Precaution: Monitor the smooth tracking of the Catheter over the guidewire during operation. If resistance is encountered during the procedure, remove Catheter and flush guidewire lumen by running the Catheter within a heparinized saline bowl where the System is turned ON with the Catheter tip completely submerged in heparinized saline to actively clear blood and/or debris that may be within the catheter lumen.

Precaution: Verify that the proximal end of the guidewire is secured within the torquer when the Handle is ON. Failure to do so may result in guidewire rotation and/or loss of position within the vessel.

III. Crossing the Lesion and Debulking

1. Under fluoroscopic guidance, while proximal to the lesion, turn the Phoenix Atherectomy System ON using the switch on the Handle for 5 to 10 seconds. Make sure that fluid exits out of the disposal outlet. Do not use any Catheter found to be non-functional.
2. Advance the Catheter slowly at a rate of 1 mm/sec and carefully through the lesion. In highly stenosed lesions or lesions ≥ 10 cm in length, periodically pause and withdraw the Catheter slightly to allow improved blood flow and plaque removal during cutting. Continue to advance until the distal tip of Catheter has crossed the lesion.
3. The Phoenix Atherectomy System must remain ON to remove plaque.
4. Withdraw the Catheter until the distal tip is proximal to the lesion. Image the lumen and/or repeat straight cutting through the lesion if desired.
5. Turn OFF the Phoenix Atherectomy System with the switch on the Handle.

IV. Phoenix Atherectomy Catheter Removal

1. To remove the Phoenix Atherectomy Catheter, stabilize the guidewire across the lesion and use a standard over-the-wire guidewire management technique to carefully remove the Phoenix Atherectomy Catheter out of the sheath under fluoroscopic guidance. Hold the guidewire firmly during the Catheter removal process by manually securing the guidewire or securing ("locking") a torque device or Wire Support Clip to the guidewire.
2. Removal of the Catheter can be accomplished by switching ON the Handle ("running the handle") per the following steps:
 - a. Lock a torque device or Wire Support Clip on the proximal end of the guidewire a distance back from the Introducer Sheath.
 - b. Hold the torque device manually or use the Wire Support Clip to prevent guidewire rotation, and turn ON the Phoenix Atherectomy System using the switch on the handle.
 - c. Under fluoroscopy, remove the Catheter by feeding the wire into the handle until the Catheter tip exits the Introducer Sheath.
 - d. Turn OFF the Phoenix Atherectomy System with the switch on the Handle.
 - e. **Note:** The Wire Support Clip should be used as described in **Step II.6** to immobilize the guidewire from rotation during Catheter removal.
3. Perform a post-atherectomy angiogram.
4. If the Phoenix System is used to treat multiple lesions where the Catheter is removed and re-inserted through the Introducer Sheath, the Catheter must be flushed between re-insertions by running the Catheter within a heparinized saline bowl where the System is turned ON with the Catheter tip completely submerged in heparinized saline to actively clear blood and/or debris that may be within the catheter lumen.

V. Post-Phoenix Atherectomy Procedure

1. If adjunctive angioplasty or stenting is judged to be necessary by the physician, insert the device of choice over the guidewire. Perform angioplasty or stenting per standard practice.
2. Perform a post-procedure angiogram.

VI. Removing the Sheath

Leave the sheath in situ until the hemodynamic profile becomes normal. Close the puncture site per routine practice.

CLINICAL DATA:

The EASE clinical study was a prospective, multi-center, non-randomized, single-arm study to evaluate the safety, effectiveness, and performance of the Phoenix Atherectomy System in the percutaneous atherectomy treatment of de novo and restenotic atherosclerotic lesions of infra-inguinal lower extremity arteries. 128 total patients were enrolled across 16 centers for treatment. 105 subjects (with 123 treated lesions) comprising the Per-Protocol group were enrolled in accordance with the vessel sizing and anatomic location requirements corresponding to each Phoenix Catheter device size. The primary safety endpoint was freedom from major adverse events (MAE) at 30 days, as assessed by an independent physician adjudicator. The primary effectiveness endpoint was acute technical success, defined as the ability to achieve a residual diameter stenosis $\leq 50\%$ as assessed by the study investigator at the time of treatment with the Phoenix Atherectomy System (prior to any adjunctive therapy) using quantitative angiography or visual estimate. Patients were monitored for 6 months following baseline atherectomy treatment.

The 30-day freedom from MAE rate was 94.3% (99/105) and technical success was achieved in 95.1% (117/123) of the lesions treated in the per protocol group. Both primary endpoints met the predefined IDE study success criteria. Table 2 summarizes the safety and effectiveness data resulting from the EASE study.

Following treatment of the first 36 subjects enrolled in the EASE Study, a higher than anticipated complication was observed. Of the first 36 subjects treated, there were 19 serious procedure or device related adverse events and 6 adverse events met the definition of major adverse events (MAE). Specific to the MAE, there have been 5 dissections/perforations and one vessel occlusion requiring treatment. AtheroMed initiated a review of available data to assess the underlying causes. The analysis of these early results led to additional sizes of the Phoenix Atherectomy Catheter being added to the EASE Study protocol, as well as updated recommendations regarding minimum vessel diameters and anatomical locations for device use (as reflected in Table 1 of this document).

Table 2: Per Protocol Summary of Safety and Effectiveness in EASE (n=105 patients and 123 lesions)

Parameter	Results (n=105 patients, 123 lesions)
Demographics:	
% Male	60% (63/105)
Age	71.7 yr ± 10.0 years
% CLI	32% (34/105)
Diabetes	53% (56/105)
Target Lesion Locations & Diameters	
SFA	28.5% (35)
Popliteal	19.5% (24)
Tibial-Peroneal Trunk	11.4% (14)
Anterior Tibial	17.1% (21)
Poster Tibial	14.6% (18)
Peroneal	8.9% (11)
RVD - proximal	3.6 ± 0.7 mm
RVD - distal	3.5 ± 0.7 mm
Target Lesion Length	34 ± 30 mm
Pre-Phoenix % Diameter Stenosis	89.1% ± 9.8%
Post-Phoenix % Diameter Stenosis	34.6% ± 16.2%
Final Post-Procedure % Diameter Stenosis	10.5% ± 9.8%
Phoenix Treatment Time	5.9 ± 4.7 minutes
Technical Success (lesions, Phoenix Catheters) ¹	95.1% (117/123)
Procedural Success (lesions) ²	99.2% (122/123)
Clinical Results	
ABI	
Baseline	0.70 ± 0.23
30 day	0.96 ± 0.21
6 mo.	0.84 ± 0.29
Rutherford Score	
Baseline	3.3 ± 1.1
6 mo.	1.2 ± 1.4
Rutherford Class Improvement ≥ 1 at 6 months	80%
MAE through 30 Days (by type):	
Clinically-Driven TLR	0.9% (1/105)
Perforation	1.9% (2/105)
Grade C or Greater Dissection	0.9% (1/105)
Symptomatic Distal Emboli ³	0.9% (1/105)
Unplanned (Toe) Amputation	2.9% (3/105)
Abrupt closure/Occlusion	0% (0/105)
MAE through 6 mo. (by type), Cumulative:	
Clinically Driven TLR	12% (12/101)
Perforation	2% (2/101)
Grade C or Greater Dissection	1% (1/101)
Symptomatic Distal Emboli ³	1% (1/101)
Unplanned (Toe) Amputation	4% (4/101)
Abrupt closure/Occlusion	0% (0/101)
Cardiovascular Related Death	1% (1/101)
Freedom from TLR at 6 mo ⁴	88%
Freedom from TVR at 6 mo ⁴	86%
¹ Residual Stenosis ≤ 50% post-Phoenix	
² Residual Stenosis ≤ 30% final	
³ Clinical signs or symptoms of distal emboli detected in treated limb distal to the treated lesion after the index procedure and requiring mechanical or pharmacological means to improve flow.	
⁴ Kaplan-Meier Estimates	

SPECIAL PATIENT POPULATIONS:

The safety and effectiveness of the Phoenix Atherectomy System has not been established in the following patient populations:

- Patients with known allergies to metals (such as Nitinol, stainless steel or other stent materials)
- Patient unable to take appropriate anti-platelet therapy
- Patients with no distal run-off
- Patients with in-stent restenosis at the peripheral vascular site
- Patients with severe, persistent hypertension (systolic pressure > 180 mm Hg)
- Patients with a history of coagulopathy or hypercoagulable blood disorder
- Patients receiving hemodialysis or with impaired renal function (creatinine is > 2.5mg/dL) at the time of treatment
- Patients with evidence of intracranial or gastric intestinal bleeding within the past three months
- Patients with severe trauma, fracture, major surgery or biopsy of a parenchymal organ within the past 3 months
- Patients who are pregnant or nursing a child
- Patients < 18 years


DEVICE DISPOSAL:

The Phoenix Atherectomy Catheter is disposable per hospital biohazard procedures. The Phoenix Handle is NOT disposable per standard biohazard procedures but should be disposed per standard hazardous waste procedures.

CAUTION: Do not incinerate the Phoenix Handle. It contains Lithium batteries and should be disposed per standard hospital hazardous waste procedures.

The Phoenix Atherectomy System has been tested and found to comply with the EMC limits for the Medical Device Directive 93/42/EEC (IEC/EN 60601-1-2). These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates, uses and can radiate radiofrequency energy and, if not used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular setting.

	Use Before Date		Nonpyrogenic
	Store in a dry, dark, cool place		Content: One (1)
	Do not use open or damaged packages		Do not Re-Sterilize
	Catalogue number		Batch Code
	Sterilized using Ethylene Oxide		Single Use Only
	Defibrillation Proof Type CF Applied Part		Prescription Only

	Volcano AtheroMed, Inc. 1530 O'Brien Dr. Suite A Menlo Park CA 94025 USA (800) 228-4728 (916) 638-8008 (916) 638 8112 fax		Volcano Europe BVBA/SPRL Excelsiorlaan 41 B-1930 Zaventem, Belgium +32.2.679.1076 +32.2.679.1079 fax
---	---	--	---

©2019 Volcano AtheroMed, Inc. All rights reserved.

ATHEROMED, PHOENIX ATHERECTOMY, PHOENIX, and the ATHEROMED and PHOENIX ATHERECTOMY logos are trademarks or registered trademarks and are the property of Koninklijke Philips N.V. in the U.S. and other countries.

The PHOENIX ATHERECTOMY System (including components and/or methods thereof) may be protected under one or more pending United States Patent Applications, as well as corresponding international patent applications.

[PATENT www.philips.com/patent](http://www.philips.com/patent)

REPORTING OF A SERIOUS INCIDENT

If a serious incident has occurred in relation to the device, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. A serious incident means any incident that directly or indirectly led, might have led or, in case of recurrence, could lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's, fetus or other person's state of health, or a serious public health threat.

REACH DECLARATION:

REACH requires Philips Healthcare (PH) to provide chemical content information for "Substances of Very High Concern (SVHC) if they are present above 0.1% of the product weight. The SVHC list is updated on a regular basis. Therefore, refer to the following Philips REACH website for the most up-to-date list of products containing SVHC above the threshold: www.philips.com/REACH.

LIMITED WARRANTY AND DISCLAIMER

Volcano AtheroMed, Inc. warrants that reasonable care has been used in the design and manufacture of this medical device.

Handling, storage and preparation of this medical device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the control of Volcano AtheroMed, Inc. directly affect this medical device and the results obtained from its use. Further, no representation or warranty is made that a Volcano AtheroMed, Inc. product will not fail. Volcano AtheroMed, Inc. disclaims responsibility for any medical complications – including death – directly or indirectly resulting from the use of this product. Except as expressly provided by this limited warranty, VOLCANO ATHEROMED, INC. IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF ITS PRODUCTS WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Volcano AtheroMed, Inc. assume no liability with respect to medical devices reused, reprocessed or resterilized, or if the device is not used by the stated "Use By" or "Expiration" date, or if the packaging is opened or damaged before use and makes no warrants – expressed or implied – including but not limited to merchantability of fitness for a particular purpose with respect to such medical device.

THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, AND VOLCANO ATHEROMED, INC. MAKES NO WARRANTY – EXPRESS OR IMPLIED – INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

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