## ATB™ DENSITOMETER CONTROL KIT

Product for the DENSIMAT quality control

### INTENDED USE

The ATB™ DENSITOMETER CONTROL KIT is designed to test the DENSIMAT.

Regular use of the kit enables the accuracy of the instrument to be monitored and therefore the quality of the inoculum to be controlled.

### PRINCIPLE

The kit consists of several ampules enabling to test the entire instrument reading scale.

Each ampule has a reference (number) and its true value, expressed in [McFarland (McF)], printed on the label.

By simply placing these ampules in the DENSIMAT one at a time, the value indicated by the instrument can be compared to the value printed on the ampule label.

### CONTENT OF THE KIT (5 ampules)

- 1 ampule **0**
- 1 ampule **0.5**
- 1 ampule **3**
- 1 ampule **6**
- 1 ampule **> 7.5 McF**
- 1 package insert provided in the kit or downloadable from [www.biomerieux.com/techlib](http://www.biomerieux.com/techlib)

### COMPOSITION OF THE MEDIUM

Demineralized water	to make required McF values
Gelling agent/Latex/Ca <sup>++</sup>	

### MATERIAL REQUIRED BUT NOT PROVIDED

- Fluffless cloth
- Swabs
- Alcohol (90°)
- Cleaning tool (provided with the DENSIMAT)

### WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.** This test is intended for use by trained laboratory professionals.
- Do not open the ampules.
- Do not shake the ampules.

### STORAGE CONDITIONS

- The kit should be stored at room temperature (**+15°C/+25°C**).
- The kit can be stored until the expiry date indicated on the packaging.
- Do not heat the ampules.
- Do not expose the ampules to freezing temperatures.

### TEST FREQUENCY

It is recommended to test the accuracy of the DENSIMAT once a month.

### INSTRUCTIONS FOR USE

Clean the instrument before testing.

For cleaning instructions, refer to the User Manual provided with the DENSIMAT (see chapter on maintenance).

Prior to use, wipe each ampule clean using a dry fluffless cloth.

### READING AND INTERPRETATION

1. Place the ampule in the instrument and turn it in the well until the lowest stable value is obtained.
2. Record the results in the table at the end of this package insert.
3. Test the instrument by inserting the ampules into the reading well in the following order:
  - Ampule 0

The result should be between 0 and 0.2 McF.

- Ampule 0.5
- Ampule 3

The expected results are the value (McF) indicated on the ampule label  $\pm 0.2$  McF.

- Ampule 6

The expected result is the value (McF) indicated on the ampule label  $\pm 0.5$  McF.

- Ampule > 7.5 McF

When this ampule is inserted, DENSIMAT displays the value > 7.5 McF.

#### LIMITATIONS OF THE METHOD

- ATB™ DENSITOMETER CONTROL KIT ampules containing cracked agar may produce results that are not conform. Do not use the ampules.
- ATB™ DENSITOMETER CONTROL KIT ampules containing liquid agar may produce results that are not conform. Do not use the ampules.

#### EXPECTED RESULTS

If the results obtained are conform to the specifications indicated in the INSTRUCTIONS FOR USE paragraph, the instrument is functioning correctly.

If they are not conform, perform the test again as described in the INSTRUCTIONS FOR USE paragraph.

If the instrument is not functioning correctly, contact bioMérieux SA or your local distributor.

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#### INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Date of manufacture

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\* **True values:** True values, expressed in McF, which are printed on the ampule label.

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Release date	Part Number	Change Type	Change Summary
2020/05	09070H	Administrative	Improvements to match the bioMérieux templates and style guide and comply with the IVDR (EU) 2017/746 regulation.

BIOMERIEUX, the BIOMERIEUX logo and ATB are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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
**REF 20150 / 70640 / 70700**

## API® NaCl 0.85 % Medium API® Suspension Medium

Media for the preparation of organism suspensions

### INTENDED USE

API® NaCl 0.85 % Medium and API® Suspension Medium are designed to be used as suspension media with API®, ID 32 or ATB™ products, in the conditions defined in the package inserts of these products.

### CONTENT OF THE KIT


- 100 ampules
- 1 package insert provided in the kit or downloadable from [www.biomerieux.com/techlib](http://www.biomerieux.com/techlib)

### COMPOSITION

<b>API® NaCl 0.85 % Medium</b> 2 mL (Ref 20070) 3 mL (Ref 20040) 5 mL (Ref 20230)	Sodium chloride	8.5 g
	Demineralized water	1000 mL
<b>API® Suspension Medium</b> 2 mL (Ref 70700) 3 mL (Ref 70640) 5 mL (Ref 20150)	Demineralized water	

### WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.** This test is intended for use by trained laboratory professionals.
- **For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.**
- Do not use media after the expiry date.
- Upon receipt, discard any kits which are visibly deteriorated: for example, damaged packaging, stains.
- The media are for single use only and should not be reused.
- Allow media to come to room temperature before use.
- Open ampules carefully as follows:

- Place the ampule in the ampule protector.
- Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
- Press the cap down as far as possible.
- Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
- Take the ampule out of the ampule protector and put the protector aside for subsequent use.
- Carefully remove the cap.

- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI M29-A, Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories – CDC/NIH – Latest edition", or to the regulations currently in use in each country.

### STORAGE CONDITIONS

- The media should be stored at +2°C/+30°C until the expiry date indicated on the ampule label.
- After opening the kit, check that the ampules are intact. If any ampules are broken, the kit must be disposed of in an appropriate container.

**WASTE DISPOSAL**

Unused reagents may be considered as non hazardous waste and disposed of accordingly.

Dispose of all used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

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**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalog number
	<i>In Vitro</i> Diagnostic Medical Device
	For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

**LIMITED WARRANTY**

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the "System") other than as set forth in the IFU.

**For more detailed information, consult the package insert of the strips with which the media are used.**

**REVISION HISTORY**

Change type categories

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

**Note:** *Minor typographical, grammar, and formatting changes are not included in the revision history*

Release date	Part Number	Change Type	Change Summary
2020/05	11791D	Administrative	Improvements to match the bioMérieux templates and style guide and comply with the IVDR (EU) 2017/746 regulation.

BIOMÉRIEUX, the BIOMÉRIEUX logo, API and ATB are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

For users in the European Union (Regulation (EU) 2017/746) and in countries with similar requirements: Should a serious incident occur during the use of this device or as a result of its use, please report it to the manufacturer and/or their authorized representative as well as to your national authority.

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CLSI yra prekės ženklas, priklausantis „Clinical Laboratory and Standards Institute, Inc“.