# Xpert Check™ Package Insert

**REF** XPERTCHECK-CE-5

This product is sold under license from Molecular Probes, Inc.

Copyright ©2019 Cepheid.

**Trademark Information**

Cepheid®, the Cepheid logo, GenXpert® and Xpert® are registered trademarks of Cepheid.

Xpert Check is a trademark of Cepheid.

Adobe® and Reader® are registered trademarks of Adobe Systems Incorporated.

Windows® is a registered trademark of Microsoft Corporation.

Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089-1189  
USA

Phone: +1.408.541.4191  
Fax: +1.408.541.4192

# Preface

---

## About This Document

The *Xpert Check™ Package Insert* provides instructions on running Xpert Check software for checking module performance.

## Safety Information

Read and understand any safety information presented in this document before you begin operating the instrument. Make sure you follow the precautionary statements presented in this guide:

### Caution

---

**Indicates that damage to the system, loss of data, or invalid results could occur if the user fails to comply with the advice given.**

---

### Important

---

**Highlights information that is critical for the completion of a task or the optimal performance of the system.**

---

### Note

---

Identifies information that applies only in special cases.

---

## Related Documents

For other information outside the scope of this document, see the following publications:

- GeneXpert Dx Operation Manual
- Infinity Operation Manual

## Cepheid Headquarters Locations

Corporate Headquarters	European Headquarters
Cepheid 904 Caribbean Drive Sunnyvale, CA 94089-1189 USA	Cepheid Europe SAS Vira Solelh 81470 Maurens-Scopont France
Telephone: +1 408.541.4191	Telephone: +33.563.825 300
Fax: +1 408.541.4192	Fax: +33.563.825 301
<a href="http://www.cepheid.com">www.cepheid.com</a>	<a href="http://www.cepheidinternational.com">www.cepheidinternational.com</a>

## Technical Assistance

Before contacting Cepheid Technical Support, collect the following information:

- Product name
- Serial number of the instrument
- Error messages (if any)
- Software version and, if applicable, Computer Service Tag number

Region	Telephone	Email
US	+ 1 888 838 3222	techsupport@cepheid.com
Australia and New Zealand	+ 1800 130 821 + 0800 001 028	techsupportANZ@cepheid.com
Belgium, Netherlands and Luxembourg	+ 33 563 825 319	support@cepheideurope.com
Brazil and Latin America	+ 55 11 3524 8373	latamsupport@cepheid.com
China	+ 86 021 5406 5387	techsupportchina@cepheid.com
France	+ 33 563 825 319	support@cepheideurope.com
Germany	+ 49 69 710 480 480	support@cepheideurope.com
India, Bangladesh, Bhutan, Nepal and Sri Lanka	+ 91 11 48353010	techsupportindia@cepheid.com
Italy	+ 39 800 902 567	support@cepheideurope.com
Portugal	+ 351 800 913 174	support@cepheideurope.com
Spain	+ 34 919 90 67 62	support@cepheideurope.com
South Africa	+ 27 861 22 76 35	support@cepheideurope.com
United Kingdom	+ 44 3303 332 533	support@cepheideurope.com
Other European, Middle East and African countries	+ 33 563 825 319 + 971 4 253 3218	support@cepheideurope.com
Other countries not listed above	+ 1 408 400 8495	techsupport@cepheid.com

Contact information for other Cepheid offices is available on our website at [www.cepheid.com](http://www.cepheid.com) or [www.cepheidinternational.com](http://www.cepheidinternational.com) under the **SUPPORT** tab. Select the **Contact Us** option.

## Table of Symbols

Symbol	Meaning
	Catalog number
	<i>In vitro</i> diagnostic medical device
	Batch code
	Do not reuse
	Caution
	Consult instructions for use
	Manufacturer
	Country of manufacture
	Contains sufficient for <n> tests
	Expiration date
	Control
	CE marking - European Conformity
	Authorized representative in the European Community
	Temperature limitation
	Biological risks

Cepheid  
 904 Caribbean Drive  
 Sunnyvale, CA 94089-1189  
 USA

Phone: +1.408.541.4191  
 Fax: +1.408.541.4192

Cepheid Europe SAS  
 Vira Solelh 81470 Maurens-Scopont  
 France

Tel: +33.563.82.53.00  
 Fax: +33.563.82.53.01

Email: support@cepheideurope.com

# Table of Contents

---

<b>1</b>	<b>Introduction</b> . . . . .	<b>1-1</b>
1.1	Proprietary Name . . . . .	1-1
1.2	Common or Usual Name . . . . .	1-1
1.3	Intended Use. . . . .	1-1
1.4	Summary and Explanation . . . . .	1-1
1.5	Reagents and Instruments . . . . .	1-2
1.5.1	Materials Provided . . . . .	1-2
1.5.2	Storage and Handling. . . . .	1-2
1.5.3	Materials Required but Not Provided . . . . .	1-3
1.6	Limitations. . . . .	1-3
1.7	Warnings and Precautions . . . . .	1-3
1.8	Chemical Hazards. . . . .	1-4
1.9	Assistance and Contact Information . . . . .	1-4
1.10	Software Buttons, Icons and Symbols . . . . .	1-5
<b>2</b>	<b>Procedure</b> . . . . .	<b>2-1</b>
2.1	System Preparation . . . . .	2-1
2.1.1	GeneXpert Dx Preparation . . . . .	2-1
2.1.2	Infinity-48 Preparation . . . . .	2-1
2.1.3	Infinity-48s or Infinity-80 Preparation . . . . .	2-2
2.2	GeneXpert and Infinity Optics Cleaning . . . . .	2-2
2.2.1	Lens Cleaning Procedure. . . . .	2-2
2.3	Data Collection Procedure: GeneXpert Dx and Infinity. . . . .	2-5
2.3.1	Xpert Check Completion For Non-internet Connected Users. . . . .	2-24
2.3.2	Obtaining the Xpert Check Code for Non-Internet Connected Users . . . . .	2-33
2.4	Return System to Normal Operation. . . . .	2-40
2.4.1	GeneXpert Dx. . . . .	2-40
2.4.2	Infinity-48 . . . . .	2-40
2.4.3	Infinity-48s or Infinity-80 . . . . .	2-40
2.5	Information Key Screen . . . . .	2-41
2.5.1	Reasons to Repeat Xpert Check with a New Cartridge . . . . .	2-41
2.5.2	Reasons to Repeat Xpert Check with the Same Cartridge . . . . .	2-41
2.5.3	Application of Xpert Check Code . . . . .	2-41

# List of Figures

---

Figure 2-1	Lens Cleaning Brush (300-8330) . . . . .	2-2
Figure 2-2	Inserting the Cleaning Brush into the I-CORE Slit . . . . .	2-3
Figure 2-3	Terms of Service Screen . . . . .	2-6
Figure 2-4	Xpert Check Login Screen . . . . .	2-7
Figure 2-5	Login Error Screen . . . . .	2-7
Figure 2-6	Home Screen . . . . .	2-8
Figure 2-7	Contact Information Screen - Page 1 . . . . .	2-9
Figure 2-8	Contact Information Screen - Page 2 . . . . .	2-9
Figure 2-9	Open Module Doors Screen . . . . .	2-10
Figure 2-10	Module preparation Screen . . . . .	2-10
Figure 2-11	Select modules for data collection Screen . . . . .	2-11
Figure 2-12	Confirm module selection Screen . . . . .	2-12
Figure 2-13	Error Screen Examples . . . . .	2-12
Figure 2-14	Scanning the Cartridge Barcode Screen . . . . .	2-13
Figure 2-15	Confirm skip Screen . . . . .	2-14
Figure 2-16	Venting the Cartridge by Opening and Closing the Cartridge Lid - Animated Screen . . . . .	2-15
Figure 2-17	Loading the Cartridge into the Module . . . . .	2-16
Figure 2-18	Data Collection Screen . . . . .	2-17
Figure 2-19	Test Completion Screen - Successful . . . . .	2-18
Figure 2-20	Test Completion Screen - Unsuccessful Module Checking Example . . . . .	2-18
Figure 2-21	Confirm Continue Screen Overlay . . . . .	2-19
Figure 2-22	Uploading Xpert Check Data Screen . . . . .	2-20
Figure 2-23	Upload incomplete Error Screen . . . . .	2-21
Figure 2-24	Downloading Xpert Check code Screen . . . . .	2-21
Figure 2-25	Applying Xpert Check code Screen . . . . .	2-22
Figure 2-26	Xpert Check complete Screen . . . . .	2-22
Figure 2-27	Exit the Program . . . . .	2-23
Figure 2-28	Confirm exit Screen . . . . .	2-23
Figure 2-29	Data collection in progress Screen . . . . .	2-24
Figure 2-30	Test Completion screen - Successful . . . . .	2-24
Figure 2-31	Write Xpert Check data to CD Screen - Step 1 . . . . .	2-25
Figure 2-32	Write Xpert Check data to CD Screen - Step 2 . . . . .	2-26
Figure 2-33	CD Writing Program - Opening Screen . . . . .	2-27

## List of Figures

---

Figure 2-34	CD Writing program - Insert a writable disk to continue Screen - Example . . . . .	2-28
Figure 2-35	CD Writing/Burning Progress Screen. . . . .	2-29
Figure 2-36	CD Writing Completion Screen . . . . .	2-30
Figure 2-37	Final Screen from Windows 7 After CD Writing has Completed . . . . .	2-30
Figure 2-38	Send the CD to your Authorized Service Provider Screen - Step 3. . . . .	2-31
Figure 2-39	Error writing Xpert Check data CD Screens - Two Examples . . . . .	2-31
Figure 2-40	Exiting Screen . . . . .	2-32
Figure 2-41	Data CD Labeling Example . . . . .	2-32
Figure 2-42	Home Screen, showing Enter Xpert Check Code Button. . . . .	2-34
Figure 2-43	Enter Xpert Check code Screen. . . . .	2-34
Figure 2-44	Xpert Check Code File - Example . . . . .	2-36
Figure 2-45	Applying Xpert Check code Screen Example . . . . .	2-37
Figure 2-46	Xpert Check complete Screen . . . . .	2-37
Figure 2-47	Xpert Check status Screen . . . . .	2-38
Figure 2-48	Xpert Check Summary Report Example - Page 1 . . . . .	2-39
Figure 2-49	Information key Screen . . . . .	2-41

# 1 Introduction

---

**Important**

**Read and understand this entire document before performing the data collection procedure.**

---

## 1.1 Proprietary Name

Xpert Check™

## 1.2 Common or Usual Name

Xpert Check

## 1.3 Intended Use

The Xpert Check kit is part of a check, verification, and hardware test system for 6-Color GeneXpert modules. The Xpert Check kit is used to check the optical system, verify the thermal system and perform a series of system-level tests to ensure full system functionality within Cepheid's instrument servicing specifications. One Xpert Check cartridge is usually used to check a single module in conjunction with the Xpert Check software. In certain cases where a retest is required, multiple cartridges may be necessary to test a module.

## 1.4 Summary and Explanation

The GeneXpert (GX) module is the basis for all GeneXpert instrument systems worldwide. Cepheid recommends that the system be checked for proper operation on an annual basis. Based upon the usage and care of each system, checks may be recommended more frequently. The system is designed to detect module issues with the internal assay controls.

The Xpert Check kit includes reagents for the optical checking and performance verification of the module. Probe Check Controls (PCCs) verify reagent rehydration, PCR tube filling in the cartridge, probe integrity, and reagent stability. Thermal performance is verified via proprietary thermal probe chemistries, and module hardware performance is tested and verified by a suite of subsystem-specific tests which exercise all critical elements of the GX module.

The Xpert Check process consists of two phases. The first phase is the execution of module testing using the cartridges contained within this kit. The second phase consists of a Cepheid Quality Assurance Review, followed by the issuance of an Xpert Check code to complete the Xpert Check process. **The Xpert Check process is not complete until this code is applied to the system.**

## 1.5 Reagents and Instruments

### 1.5.1 Materials Provided

The Xpert Check kit contains the following:

**Table 1-1. Kit Contents**

Description	Quantity
<b>Xpert Check cartridges with integrated reaction tubes</b>	<b>5 per kit</b>
<b>Each cartridge contains the following materials:</b>	
• Bead 1	<b>1 per cartridge</b>
• Reagent 1	<b>1.0 mL per cartridge</b>
<b>I-CORE Lens Cleaning Brush</b>	<b>4 per kit</b>
<b>PI/Software (P/N 950-0413)</b>	<b>1 per kit</b>
<b>Data CD</b>	<b>1 per kit</b>

**Note** Safety Data Sheets (SDS) are available at [www.cepheidinternational.com](http://www.cepheidinternational.com) under the **SUPPORT** tab.

### 1.5.2 Storage and Handling

- Store the Xpert Check cartridges at 2-28 °C. Wait at least 10 minutes after removal from cold storage before using, allow a cartridge to reach ambient temperature.
- Use the cartridge within 48 hours of opening the foil pouch.
- Discard cartridges that have been removed from their foil-wrapped pouches outside of the approved usage interval.
- Do not use cartridges that have passed the expiration date.
- Do not open a cartridge lid until you are ready to perform testing.
- Cartridge lid must be opened (vented) prior to use of a cartridge; however, no sample is required for testing.
- Discard all used and unused materials, including cleaning brushes and cartridges once the Xpert Check session is completed.

**Note** Contents of cartridges are non-hazardous.

### 1.5.3 Materials Required but Not Provided

- GeneXpert Dx or Infinity System with Cepheid-supplied computer and barcode scanner.
- 6-color GeneXpert instrument.

### 1.6 Limitations

- For use only with 6-Color GeneXpert modules (running GeneXpert Dx software version 4.0 and above). 4-Color GeneXpert modules (including Bio-threat modules) and Dual-Cal modules cannot run Xpert Check and must be tested by Cepheid Service.
- Use of the Xpert Check kit does not guarantee that the GeneXpert instrument will be free of hardware failures, nor does it take the place of a Cepheid Service Agreement.

### 1.7 Warnings and Precautions

---

**Support for Windows XP ended April 8, 2014. Microsoft no longer provides security updates or technical support for the Windows XP operating system. It is critical that you upgrade now to a newer operating system, such as Windows 10.**

**Important**

Please contact <https://www.microsoft.com/en-us/microsoft-365/windows/end-of-windows-xp-support> for Windows XP support information.

In addition, please contact your local Cepheid Technical Support if you have questions about using Windows XP.

---

**After January 14, 2020, Microsoft will no longer provide security updates or support for PCs running Windows 7. It is recommended that you upgrade to Windows 10.**

**Important**

Please contact <https://www.microsoft.com/en-us/microsoft-365/windows/end-of-windows-7-support> for Windows 7 support information.

In addition, please contact your local Cepheid Technical Support if you have questions about using Windows 7.

- 
- Follow your institution's safety procedures for working with chemicals.
  - Do not add sample or other reagents to the Xpert Check cartridges.
  - Do not use a cartridge that has a damaged reaction tube.
  - Do not use cartridges from visibly damaged or compromised foil pouches.
  - Contact your local Cepheid Technical Support office for replacement of damaged kit contents.
  - Do not use a cartridge if it is dropped.

②

- Each single-use Xpert Check cartridge is used to process one test. Do not reuse spent cartridges.
- Each cleaning brush is intended for use in a single module. Do not re-use brushes in multiple modules.
- Do not open a cartridge package or break the lid seal until you are ready to perform testing.
- Allow the Xpert Check cartridge to come to ambient temperature prior to use if it has been placed in cold storage. Wait at least 10 minutes after removal from cold storage before using.
- Do not store single cartridges. Cartridges left over from an Xpert Check session, including pouched/unopened cartridges should be discarded along with spent cartridges.
- Do not use cartridges whose shelf life has expired. The system will detect expired cartridges and abort the test.
- Once a cartridge barcode has been scanned, do not substitute another cartridge in place of the scanned cartridge.
- If using an internet-enabled Xpert Check, up-to-date anti-virus software must be installed on the desktop or laptop computer with updated virus definition files, prior to executing Xpert Check.
- Prior to running Xpert Check, ensure that the environmental operating temperature is within the correct limits (15 °C–30 °C). Xpert Check will render a system's modules unavailable if the internal temperature is above 40 °C. The internal temperature can be verified in the Maintenance section of the GeneXpert DX software. Do not proceed under these conditions.
- Xpert Check expects the same computer to be used throughout the entire process. The computer installed with the GeneXpert system should be used, and not another computer from a different GeneXpert system.
- The Xpert Check code will expire if not applied within 45 days of completion of running Xpert Check.

## 1.8 Chemical Hazards

According to Regulation (EC) No. 1272/2008 (CLP), this material is not considered hazardous.

## 1.9 Assistance and Contact Information

For a complete listing of Cepheid technical support, service support, sales support, and headquarters contacts, please see [Technical Assistance](#), in the Preface of this document.

## 1.10 Software Buttons, Icons and Symbols

Table 1-2. Software Buttons, Icons and Symbols

Symbol	Definition
	<b>Information.</b> Click on this icon to obtain additional information. Displays the Information Key workspace screen which has an explanation of the various module icon displays.
	<b>Continue.</b> This icon is located at the bottom of most screens. Click on this icon to advance the display to the next screen.
	<b>Continue to End.</b> Clicking on this icon moves the user to the last screen.
	<b>Exit.</b> Exits the Xpert Check application.
	<b>About.</b> Brings up the About screen which brings up the name of the software, the software version number, copyright notice, etc.
	<b>Home.</b> Go to the Home screen.
	<b>Repeat/Retry.</b> Retry loading an Xpert check cartridge to attempt to check a module that has had an unsuccessful test of a minor nature or if the cartridge has not been vented by the user. Used on the 'Check Test' screen.
	<b>Back.</b> Clicking on this icon takes the user to the previous screen.
	<b>Cancel.</b> Cancel the current operation. In most cases this will mean going back to the previous screen. In some cases, it may mean going back to the screen before the one that started the current operation.

**Table 1-2. Software Buttons, Icons and Symbols**

Symbol	Definition
	<p><b>Select none of the modules for check.</b> Deselects all modules for checking. If you only want to check a few modules, you may deselect ALL of them, and then reselect only the ones you wish to check.</p>
	<p><b>Select all of the modules for checking.</b> The default setting for the system.</p>
	<p><b>Connectivity Status.</b> Indicates the system is able to reach the Xpert Connectivity Center.</p>
	<p><b>Connectivity Status.</b> Indicates the system is not able to reach the Xpert Connectivity Center.</p>
	<p><b>Module unsupported for Xpert Check.</b> Skip the current module and do NOT attempt to check the current module.</p>
	<p><b>Module selected for Xpert Check.</b> Module will be included when Xpert Check is run.</p>
	<p><b>Skip Current Module.</b> Skip the current module and do not attempt to Xpert Check the current module. Used on the 'Load Xpert Check Cartridges' screen.</p>
	<p><b>Skip Remaining Modules.</b> Skip all the remaining modules and do NOT attempt to Xpert Check them. Used on the 'Load Xpert Check Cartridges' screen.</p>

Table 1-2. Software Buttons, Icons and Symbols

Symbol	Definition
	<b>Module not selected for Xpert Check.</b> Module will not be included when Xpert Check is run.
	<b>Module unavailable for Xpert Check.</b> Module will not be included when Xpert Check is run.
	Indicates a module with data collection in progress.
	Indicates data collection complete.
	<b>Retest required.</b> Indicates an incomplete Xpert Check data collection. A message will notify the user that the test must be rerun. A further message will indicate if the existing cartridge can be reused for the test or if a new cartridge must be used.
	<b>Service required.</b> Contact the Cepheid Authorized Service Provider (ASP) or your local Cepheid Technical Support office.

**Table 1-2. Software Buttons, Icons and Symbols**

Symbol	Definition
	<p><b>Lost communication.</b> Contact the Cepheid Authorized Service Provider (ASP) or your local Cepheid Technical Support office.</p>
	<p><b>Burn.</b> Burn a CD containing the collected Xpert check information (for users without an active internet connection).</p>
	<p><b>Collect Xpert Check Data.</b> Leads the user through the Xpert Check data collection process.</p>
	<p><b>Enter Xpert Check code.</b> Go to the 'Enter Xpert Check Code' screen.</p>
	<p><b>Xpert Check Status.</b> Go to the Xpert Check Status screen to review Xpert Check status.</p>
	<p><b>Upload Xpert Check Data File.</b> Go to the 'Upload Xpert Check Code Data File' screen.</p>
	<p><b>Upload Xpert Check Data CD.</b> Go to the 'Upload Xpert Check Code Data CD' screen.</p>
	<p><b>Write Xpert Check Code.</b> Write an Xpert Check Code to a file.</p>
	<p><b>Read Xpert Check Code.</b> Open a file to read the Xpert Check code.</p>
	<p><b>Scan.</b> Turn the barcode scanner on, and accept the next scanned input.</p>
	<p><b>View and Print.</b> Launch the Adobe Reader so you can view and then print a PDF file.</p>

## 2 Procedure

---

### 2.1 System Preparation

**Note**

---

Prepare the system for Xpert Check by following one of the three procedures listed in this section for the GeneXpert Dx, the Infinity-48, Infinity 48s or the Infinity-80.

---

**Important**

---

**Authorized Service Providers (ASPs) who perform Xpert Check but won't be on-site when the Xpert Check code numbers come back (non-internet connection sites), should leave the user name and password for the users to log in later to enter the codes (see [section 2.3.2](#)).**

---

#### 2.1.1 GeneXpert Dx Preparation

1. Create an Administrator or Basic level User Name and Password in the GeneXpert software if one does not exist. Xpert Check requires this logon credential to be established prior to starting.
2. Have your Authorized Service Provider (ASP) code available before continuing to the next steps.
3. Exit the GeneXpert Dx software.
4. Go to [Section 2.2, GeneXpert and Infinity Optics Cleaning](#).

#### 2.1.2 Infinity-48 Preparation

1. Create an Administrator or Basic level User Name and Password in the GeneXpert software if one does not exist. Xpert Check requires this logon credential to be established prior to starting.
2. Have your Authorized Service Provider (ASP) code available before continuing to the next steps.
3. Restart the Xpertise software and switch the system from Automatic mode to Manual mode. Follow the instructions in the GeneXpert Infinity System Operator Manual for the Infinity-48.
4. Go to [Section 2.2, GeneXpert and Infinity Optics Cleaning](#).

### 2.1.3 Infinity-48s or Infinity-80 Preparation

1. Create an Administrator or Basic level User Name and Password in the GeneXpert software if one does not exist. Xpert Check requires this logon credential to be established prior to starting.
2. Have your Authorized Service Provider (ASP) code available before continuing to the next steps.
3. Exit the Infinity System software.
4. Open the glass doors following the instructions in the *Infinity Operator Manual*.
5. Go to [Section 2.2, GeneXpert and Infinity Optics Cleaning](#).

## 2.2 GeneXpert and Infinity Optics Cleaning

This procedure describes the method for removing dust and tube debris from the surface of rod lenses of the excite and detect blocks for GeneXpert Dx and Infinity modules prior to performing the Xpert Check procedure.

---

**Note** This procedure applies only to GeneXpert 6-color modules.

---

### Materials Required or Recommended for Cleaning

- 300-8330 –Applicator brush (Quantity of four Included in the Xpert Check kit)
- Disposable gloves

**Estimated Cleaning Time: 30 Seconds per module.**

### 2.2.1 Lens Cleaning Procedure

1. Select the module to be checked and manually open the door of the module.
2. If necessary, remove the cartridge from the module.

**Biological Risks**


---

**Remove the cartridge from the GeneXpert modules prior to cleaning. Failure to remove a cartridge could result in personnel being exposed to biological hazards and/or liquid biological materials spilling into the instrument and causing damage to the instrument.**

---

3. Locate the brush provided in the Xpert Check kit (see [Figure 2-1](#)).

**Figure 2-1. Lens Cleaning Brush (300-8330)**

**Note**

The brush is designed so that it will easily insert into the I-CORE slit and make contact with the rod lenses of the excite and detect blocks.

**Biological Risks**

**Make sure you wear disposable gloves for the cleaning process. Wearing gloves prevents you from being exposed to biologically hazardous materials.**

4. Wearing disposable gloves, insert the brush into the I-CORE slit in a tilted manner up to the shank insertion shoulder, as shown in [Figure 2-2](#).

**Note**

Make sure that all the bristles are fully inserted (up to the shoulder of the plastic shank of the brush) so that it does not cause unnecessary damage to the brush.

**Caution**

**Do not insert any objects into the I-CORE slit except the provided brush. Inserting any other object may damage the I-CORE.**


**Caution**

**Do not apply any solution (such as ethanol or bleach) onto the brush bristles. The brush must be completely dry when inserting it into the I-CORE slit.**

**Important**

**The brush is intended for single-use and should not be used on more than one module. Use a new brush for each module to be cleaned.**

**Figure 2-2. Inserting the Cleaning Brush into the I-CORE Slit**

## Procedure

---

5. Insert the brush into the I-CORE slit completely up to the plastic shank (shoulder) of the brush. Hold the brush firmly in the I-CORE slit, and perform cleaning of the rod lenses as described below. The entire cleaning process should take approximately 30 seconds per module.

### Note

---

Cleaning is done by moving the brush in an up and down direction within the I-CORE slit. Brush rotation, even if it has to be done, is not the main action that results in optics cleaning.

---

- A. Begin by brushing from the top of the I-CORE slit to the bottom, making sure to apply a uniform pressure when brushing from the top to the bottom of the I-CORE slit. This will ensure that most of the tube debris and dust is brushed off from the surface of the lenses.
  - B. Rotate the brush from left to right and back again, approximately 180°.
  - C. Brush once more from the top of the I-CORE slit to the bottom.
  - D. Rotate the brush again from left to right and back again, approximately 180°.
  - E. Finally, brush again from the top of the I-CORE slit to the bottom.
6. When lens cleaning is complete, remove and discard the used brush and gloves as hazardous waste.

### Important

---

**Dispose of gloves and brushes according to your institution's safety policies and procedures for hazardous waste.**

---

7. Proceed to [Section 2.3, Data Collection Procedure: GeneXpert Dx and Infinity](#).

---

---

## 2.3 Data Collection Procedure: GeneXpert Dx and Infinity

**Important**

---

Before collecting data, be sure to prepare the system for checking as described in [Section 2.1, System Preparation](#). Internet-connected users should verify their system's connectivity status prior to beginning the Xpert Check process.

---

**Note**

---

Use care in inserting CD1 into the DVD drive. Be sure the CD is fully seated in the tray before closing the drive door.

---

1. Place Software CD1 in the computer connected to the GeneXpert Dx. For the Infinity, connect the DVD drive following the instructions in the Infinity Operator Manual and insert the CD into the DVD drive.
2. This step varies with the operating system installed on your computer:
  - **Windows XP:** On the computer desktop, right-click the My Computer icon and a drop-down menu will appear. Click **Explore**, then right-click on the applicable drive letter for your DVD drive. Select **Explore** from the drop-down menu, and the files located on the CD will then be displayed. Find and right-click the XpertCheck.exe application, and when the drop-down menu appears, click **Run** to install as Administrator. When the software has been installed, a “wrench” icon will appear on the desktop.
  - **Windows 7, Windows 10:** On the computer desktop, right-click the Computer icon and a drop-down menu will appear. Click **Open**, then right-click on the applicable drive letter for your DVD drive. Select **Open** from the drop-down menu, and the files located on the CD will then be displayed. Find and right-click the XpertCheck.exe application, and when the drop-down menu appears, click **Run** to install as Administrator. When the software has been installed, a “wrench” icon will appear on the desktop.

**Note**

---

The software may take some time to load from the CD.

---

3. Double-click the “wrench” icon to launch the Xpert Check program.
4. The Terms of Service screen appears first. Use the scroll bar to read through the entire document. You will be asked to click the check box (bottom of the screen) to verify that you have read and agree to the Terms of Service before continuing. See [Figure 2-3](#).

**Note**

---

The Xpert Check software runs on Windows XP, Windows 7 or Windows 10. The screens shown in this manual are from Xpert Check software running on Windows 7. Screens for Xpert Check software running on Windows XP or Windows 10 will be similar.

---

Support for Windows XP ended April 8, 2014. Microsoft no longer provides security updates or technical support for the Windows XP operating system. It is critical that you upgrade now to a newer operating system, such as Windows 10.


**Important** Please contact <https://www.microsoft.com/en-us/microsoft-365/windows/end-of-windows-xp-support> for Windows XP support information.

In addition, please contact your local Cepheid Technical Support if you have questions about using Windows XP.

After January 14, 2020, Microsoft will no longer provide security updates or support for PCs running Windows 7. It is recommended that you upgrade to Windows 10.

**Important** Please contact <https://www.microsoft.com/en-us/microsoft-365/windows/end-of-windows-7-support> for Windows 7 support information.

In addition, please contact your local Cepheid Technical Support if you have questions about using Windows 7.


Scroll down to read the entire document.

Note: A copy of these Terms of Service is located on CD1.

Figure 2-3. Terms of Service Screen


5. After agreeing to the Terms of Service, the Login screen will appear. Log in with your GeneXpert Dx or Infinity designated Administrator level USER NAME and PASSWORD (previously assigned to you by your system administrator). After entering your login information, click the forward arrow button at the bottom of the screen to advance to the Xpert Check Home screen. See Figure 2-4.

**Note** The user name and password are the same ones you used for the GeneXpert Dx or Xpertise software.

**Figure 2-4. Xpert Check Login Screen**

In case of a login error, the following screen will appear. See [Figure 2-5](#).

**Figure 2-5. Login Error Screen**

6. If a login error occurs, recheck the **USER NAME** and **PASSWORD** entries for errors. If necessary, reenter the information and retry. After entering your login information, click the forward arrow button at the bottom of the screen to advance to the Xpert Check Home screen.
7. Obtain a sufficient number of cartridges for the number of modules to be tested.


**Important**

**Do not open cartridge packages until you are ready to scan the cartridge barcodes (in Step 17).**

**Note**

When determining the number of cartridges that will be needed for this test, the user should be aware of the number of modules that they will be checking.

8. Click the Collect Xpert Check Data icon on the Home screen (See Figure 2-6). After a few seconds, the first Contact Information screen (Figure 2-7) will appear.

**Figure 2-6. Home Screen**

9. When the first of two Contact Information screens appear (see Figure 2-7 and Figure 2-8), fill out the fields in the two screens. Use the large navigation arrows at the bottom of the screens to move between the two screens. Note that fields marked with “\*” (at the right of the entry area) are mandatory fields.



Figure 2-7. Contact Information Screen - Page 1

In this mandatory field, enter your Authorized Service Provider's ID code. The ID code is four characters, as described in

Clicking on this icon will take the user back to the first Contact Information

Clicking on this icon will take the user to the Open Module Doors screen

Figure 2-8. Contact Information Screen - Page 2

**Note**

The ASP-provided ID code for the Service Provider on the Contact Information screen consists of four characters. (As examples: US01, 1203, etc.)

10. When all information has been entered, click the forward arrow button at the bottom of page 2 of the Contact Information screen. The Open Module Doors screen will appear. See [Figure 2-9](#). Manually open all module doors to enable cartridge loading.


Figure 2-9. Open Module Doors Screen

11. After opening all the module doors, click the forward arrow button at the bottom of the screen. The Module preparation screen may appear, showing the message *Wait while modules are being prepared.* (See [Figure 2-10.](#))

**Important**

**Note that the Module preparation screen will appear only if the firmware in the modules is not 3.0.3. The screen indicates that the software is upgrading/downgrading the firmware to the modules. The next screen you see will be the screen shown in [Figure 2-11](#), the Select Modules screen.**



Figure 2-10. Module preparation Screen

12. Follow the on-screen software instructions in [Figure 2-11](#). By default, all detected modules will be marked as selected for checking.  
 On this screen, the user can click on individual modules to exclude them from being checked, if required. The modules will disappear as they are excluded.

**Note**

For excluded modules (not selected for checking), the door position (open or closed) does not matter.


**Figure 2-11. Select modules for data collection Screen**

**Important**

**GX-XVI and Infinity systems only:** When selecting modules on screen, make a note of which module lights are blinking on the system as you select each bank, to ensure cartridges are placed in the correct modules for testing.

13. After confirming the module selection shown in [Figure 2-11](#), click the white arrow at the bottom of the screen overlay, to begin scanning cartridges. If the module selection shown is incorrect, click the red X at the bottom left corner of the screen to return to the Select Modules screen and change your selection. See [Figure 2-12](#).

Module selection is shown here. In this example, one module is selected for checking.

Figure 2-12. Confirm module selection Screen

- 14. In case of an error in the preceding step, in which either no modules have been selected, or all modules have been excluded, one of the following screens will appear (Figure 2-13). Follow the on-screen instructions to select a module, or start over by returning to the Home screen or exiting the program.


Figure 2-13. Error Screen Examples

15. After confirming your module selection, you will advance to the Scan cartridge screen, where you will be prompted to scan the barcode on the Xpert Check cartridge.

**Note** Verify you have enough cartridges on hand to perform the check procedure for the desired number of modules.

16. Remove the test kit cartridge from the package for the module you've previously selected, opening only one cartridge at a time.

**Important** Allow the cartridge to reach ambient temperature before proceeding. Do not remove a cartridge from refrigerated storage and immediately use the cartridge to run this test.


17. Scan the cartridge barcode. [Figure 2-14](#) shows a cartridge barcode being scanned. Do not substitute a cartridge with another after it's been scanned.

**Note** If the barcode cannot be scanned, skip the cartridge and contact your ASP or local Cepheid Technical Support office for a replacement cartridge, if necessary. If the barcode scanner is damaged, missing or incorrectly configured, contact your ASP or local Cepheid Technical Support office for guidance.

Scanning the cartridge barcode

**Note:** To avoid inserting a wrong cartridge into a module, do not set the cartridge aside after it has been scanned. In one operation, scan the cartridge barcode, vent the cartridge and insert it into the next available (lighted) module.

**SKIP button.** Click this icon to skip the cartridge just scanned.

**Figure 2-14. Scanning the Cartridge Barcode Screen**


- A. After scanning the barcode of the cartridge, ensure you open (vent) the cartridge lid and then close it for each cartridge as directed by the software in [Step B](#) through [Step E](#) below.

**Important** Do Not add a sample or reagent to the cartridge. Use ONLY the cartridges in the Xpert Check kit provided.

**Note** After a cartridge barcode is scanned a green light will blink on the system above the module door where the cartridge is to be loaded. (See [Figure 2-17.](#))

**Note**

If, for some reason, you want to skip the cartridge just scanned, click the **SKIP** button at the bottom of the screen. An overlay, shown in [Figure 2-15](#), will appear, asking for confirmation on skipping the cartridge. To **SKIP** the cartridge, click the forward arrow at the bottom of the confirmation screen. To proceed without skipping the cartridge, click the "X" icon at the left bottom corner of the screen. You are urged to rescan a cartridge (or substitute a new cartridge if necessary) to ensure a module is not skipped.


Click the "X" to "undo" the SKIP selection and proceed with opening, venting and loading the cartridge into a module.

Click the arrow to confirm you want to skip the cartridge just scanned.

**Figure 2-15. Confirm skip Screen**

- B. Venting the cartridge (shown in [Figure 2-16](#)), for two seconds is sufficient. This screen is animated, showing the cartridge lid being opened and closed. After venting, click the forward arrow at the bottom of the screen to continue.


**Important: If a cartridge is not vented before testing, the test will be unsuccessful.**

**Figure 2-16. Venting the Cartridge by Opening and Closing the Cartridge Lid - Animated Screen**

- C. Close the cartridge lid and ensure the module door is fully opened to receive the cartridge.
- D. Load the cartridge into the module (with the cartridge reaction tube (tab) facing away from you), as directed by the animated software screens. See [Figure 2-17](#).

**Note**

Be sure to load scanned cartridges in sequence in the next available module. This will avoid loading cartridges in the wrong location or leaving modules empty.

If a module door will not close and latch, click the SKIP button at the bottom to abort testing the module and move on to the next module to continue testing.

**Important:** Be sure to remove the cartridge from a module whose door will not latch before proceeding. Notify your ASP or local Cepheid Technical Support office of the malfunctioning module door latch so it can be serviced.

**Note** that the door on a module without a cartridge loaded will not latch. That door can be left in the open position when running the other loaded modules.

**Figure 2-17. Loading the Cartridge into the Module**

- E. If you are checking additional modules, continue by scanning the next cartridge. Place each individually scanned cartridge into the next selected open module, pressing the module door securely closed until it latches. **As each module door is closed and latched, data collection will automatically start on that specific module.** The blinking green light above the module will then become steady green, indicating that checking has started.

If a module door is not closed completely (until it latches) after loading a cartridge, the screen will continue to display a message to insert a cartridge, and the check will not run.


**Important**

If you are unable to close and latch a module door after several tries, press the SKIP button at the bottom of the screen to skip the module with the faulty door and move ahead. Notify your ASP or local Cepheid Technical Support office so the module can be serviced.

18. Checking will take approximately 20 minutes to complete after the final module has been loaded for testing. When checking begins, the Data collection in progress screen appears, as shown in [Figure 2-18](#).

**Important**

If necessary, a retest of a previously run module may be started without waiting for the present module to complete its test, as described in steps 19b through 19d.

Countdown clock, showing estimated time to test completion (21 minutes).

Note: In this example, one module is being tested.

Do not exit the software program while data collection is in progress!

**Figure 2-18. Data Collection Screen**


**Important**

If you do not have an internet connection, skip to section 2.2.1 for the remainder of this procedure. If you have an internet connection, continue with step 19.

19. After test completion, the module door will open and the light above the module door will turn off. Screens similar to those shown in [Figure 2-19](#) or [Figure 2-20](#) will appear. Press the right arrow to continue.
  - A. [Figure 2-19](#) shows the completion of a successful Xpert Check data collection. When the test is complete, click the forward button at the bottom of the screen to begin uploading Xpert Check test results to the Xpert Connectivity Center.

**Important**

When uploading test results, especially multiple files, verify the selected folder destination is correct.

Click to continue to the next screen.

Note: In this example, one module is being tested.

Figure 2-19. Test Completion Screen - Successful

- B. If the test was unsuccessful, the screen shown in Figure 2-20 will appear, showing module status. A test retry must be performed. Click the Retry icon in the lower left-hand corner of the screen.

Module test unsuccessful. Test must be rerun using either the same or a new cartridge, as instructed by the on-screen display.

Continue button. Click to proceed without retesting. See Figure 2-21.

Retry Icon. Click to retest the module.

Note: In this example, one module is being tested.



Figure 2-20. Test Completion Screen - Unsuccessful Module Checking Example

- C. If the Continue arrow at the bottom of the screen is pressed when there is an unsuccessful module test displayed (as shown in Figure 2-20), the Confirm continue screen will appear. See Figure 2-21.

To return to the Press retry or continue Screen to retest the flagged module, click the “X” icon.

To continue without retesting the flagged module, click on the right arrow at the bottom of the Confirm continue screen.

**Figure 2-21. Confirm Continue Screen Overlay**

You have the option of continuing by pressing the right arrow on the Confirm continue screen overlay. Choosing this option will result in the flagged module not being retested, and you will begin uploading check data as described in [Step 20](#).

Another option is to return to the Press retry or continue Screen to Retry (retest) the flagged module by clicking the red “X” icon at the bottom left of the Confirm continue screen. The Retry procedure is described in [Step D](#) which follows.

- D. If the Retry icon (shown above in [Figure 2-21](#) at the bottom of the screen) appears, click the Retry icon and you will return to the Scan Barcode screen ([Figure 2-14](#)) to complete the retest on the affected module(s).

Note that the retest can be of two possible types:

- 1) Retry with the same cartridge: For example, a message may appear telling you to vent the cartridge, rescan it, and put it back in the module.
- 2) Retry with a new cartridge: If the cartridge was defective, or had already been used, you will be asked to replace it by scanning the barcode on a new cartridge, venting it, and loading it into the module.

**Note**

During the course of running retests, modules may need to be skipped if the user runs out of Xpert Check cartridges. Please contact your ASP or local Cepheid Technical Support office for additional Xpert Check cartridges. Rerun Xpert Check on any modules that were skipped.

**Note**


At the completion of the Xpert Check data collection process, modules determined to require service will be flagged with an orange module icon (See [Figure 2-20](#)). Please contact your local ASP or local Cepheid Technical Support office for further assistance in servicing or replacing modules.

20. After successful test completion and Xpert Check data collection, click the forward arrow to display the screen shown in [Figure 2-22](#), if you have an active internet connection. However, if you have never been internet connected, or have lost your functioning internet connection sometime during the Xpert Check test, a Download Xpert Check code error screen or an Upload incomplete error screen ([Figure 2-23](#)) may appear instead, instructing you to write Xpert Check data to a data CD to send to your ASP or local Cepheid Technical Support office. In this case, continue to the instructions beginning at [Step 4](#) (under [Section 2.3.1](#)) of this procedure to continue as a user without an internet connection.

**Note**


With a functioning internet connection, the system should proceed normally (with [Step 21](#)), and the Xpert Check code should begin downloading, as shown in [Figure 2-24](#).

---

**Figure 2-22. Uploading Xpert Check Data Screen**


Click this icon to move to the Write Xpert Check Data to CD screen. See Figure 2-31.

**Figure 2-23. Upload incomplete Error Screen**

21. When the Xpert Check data has finished uploading, a Quality Assurance check will be performed on the data. If the check is acceptable, the Xpert Check code will automatically download. See Figure 2-24.

If the test is not acceptable, the affected module(s) will require service or replacement and will be flagged with an orange icon. Please contact Cepheid or your local ASP or the local Cepheid Technical Support office for further assistance.

**Figure 2-24. Downloading Xpert Check code Screen**

22. After the Xpert Check test results have downloaded, the Xpert Check code will be applied to each successfully tested module, and those modules will then be identified with a + symbol. See Figure 2-25. As shown here, one module is being checked.


Figure 2-25. Applying Xpert Check code Screen

**Note**

In the screen shown in Figure 2-25, some modules may display the service required icon or may be grayed out if they were skipped.

- 23. After all the Xpert Check codes have been applied to the successfully-tested modules (those green modules which appear with the plus symbols applied), the Xpert Check complete screen will appear. See Figure 2-26. This screen shows the location of the Xpert Check Data report, which is available for review, if desired.




Figure 2-26. Xpert Check complete Screen

24. Remove and discard all Xpert Check cartridges. Do not save partial kits (all unused cartridges must be discarded). When complete, click the Exit icon at the top or bottom of the screen to exit the program. See [Figure 2-27](#).

**Figure 2-27. Exit the Program**

25. The screen shown in [Figure 2-28](#) appears only if you click the exit arrow in the upper right of the screen.

**Figure 2-28. Confirm exit Screen**

This completes the Xpert Check Test for an internet-connected user.


**Note** Contact your Authorized Service Provider or the local Cepheid Technical Support office concerning modules requiring service.

**Note** To view your Xpert Check results, see [Step 22](#) and [Figure 2-26](#), which shows the file path and location of Xpert Check results and the Xpert Check Summary report.

### 2.3.1 Xpert Check Completion For Non-internet Connected Users


For Non-internet connected users, you should have completed [Step 1](#) through [Step 17](#) of section 2.2 to collect data before starting this section.

1. This section begins with the Data collection in progress screen, which is similar to step 17 in section 2.2, and the screen shown in [Figure 2-18](#) for internet-connected users.

**Figure 2-29. Data collection in progress Screen**

2. After test completion, the module door will open and the light above the module will turn off. A screen similar to that shown in [Figure 2-30](#) will appear. Press the right arrow at the bottom of the screen to advance to the next screen.

**Figure 2-30. Test Completion screen - Successful**

- When the Write Xpert Check data to CD screen appears (Figure 2-31), you will be prompted to press the Eject button to remove the existing Xpert Check Software CD so you can insert the blank data CD.


**Note**

In the following step, use care in inserting the blank CD into the DVD drive. Be sure the CD is fully seated in the tray before closing the drive door.

**Important**

If you have been running this test as an internet-connected user and then lost your internet connection and received an error screen (Figure 2-23), resume your procedure beginning with the following Step 4, continuing through Step 12.

- Insert the blank CD into the DVD drive of the computer and close the DVD drive tray fully to ensure the CD will be recognized.


**Figure 2-31. Write Xpert Check data to CD Screen - Step 1**

- After inserting the blank CD, the screen will change briefly, indicating the CD has been recognized. See Figure 2-32. This screen will remain displayed until the CD writing process is complete.

**Note**

It is not necessary for the user to locate the file to write because that process is automatic.


**Figure 2-32. Write Xpert Check data to CD screen - Step 2**

6. The CD Writing Wizard or Burn to Disc screen (Figure 2-34) will then appear as an overlay of the screen shown above, in Figure 2-32.

The next screens (Figure 2-33 through Figure 2-37) show the CD writing program screens as you progress through the writing process.

- Windows® XP users: Follow the screens on the left side of the figure.
  - Windows® 7 users: Follow the screens on the right side of the figure.
  - Windows® 10 users: Follow the screens at the bottom of the figure.
- A. On the first screen, after successful recognition of the blank CD, you will be asked to provide a name for the CD that you will be writing. DO NOT simply click the Next button (as suggested on the screen) to continue the writing process with the default name that appears. Instead, type in your facility's name, such as "XYZ Hospital," in the space provided and click Next. See Figure 2-33.


**Figure 2-33. CD Writing Program - Opening Screen**

- B. If the CD is not recognized, the screen shown in [Figure 2-34](#) may appear, instead of the screen in [Figure 2-35](#), asking you to insert a writable disc to continue. Writable discs, in this case, are CDs on which you can store files. Writable discs can only be written to once, meaning that once any files are copied to the disc, they are there permanently.


A disc that has data on it is not considered to be a writable disc and will result in an error screen, as shown in [Figure 2-39](#).

**Note**


If you are unsuccessful with any part of the CD writing process, you may contact your ASP or local Cepheid Technical Support office for assistance. It is safe for you to close the Xpert Check software now because the Xpert Check files have been saved to the hard drive and you will not lose data.

Insert a Disc Screen - Windows XP


Insert a Disc Screen - Windows 7


Insert a Disc Screen - Windows 10

Figure 2-34. CD Writing program - Insert a writable disk to continue Screen - Example


- C. After successful recognition and naming of the CD, click the Next button to continue. The writing process will begin automatically.
- D. During the writing/burning process, a progress bar will appear on the screen. See Figure 2-35.

**File Writing Progress Screen - Windows XP**

**File Burning Progress Screens - Windows 7**


**File Burning Progress Screen - Windows 10**

**Figure 2-35. CD Writing/Burning Progress Screen**


- E. When the writing of the CD is complete, the screen shown in [Figure 2-36](#) will appear. Click the Finish button to exit the CD writing program.

Completion Screen - Windows XP

Completion Screen - Windows 7

Completion Screen - Windows 10

Figure 2-36. CD Writing Completion Screen

- F. On a Windows 7 computer, you may see the screen displayed in [Figure 2-37](#) after a successful CD write. Click the OK button as many times as necessary for it to disappear, before continuing.



Figure 2-37. Final Screen from Windows 7 After CD Writing has Completed


- 7. After clicking the Finish button on the CD writing screen, the Send the CD to your Authorized Service Provider Screen will appear (see [Figure 2-38](#)). Remove the completed Xpert Check data CD from the disk drive and prepare the label, as described in [Step 10](#).

**Figure 2-38. Send the CD to your Authorized Service Provider Screen - Step 3**


8. If a problem has occurred anytime during the CD writing process, an error code screen may appear. (See [Figure 2-39](#)).
  - If a CD you have inserted already contains data as shown in the error screen below at the left, remove the CD and insert a blank CD, and then press the Retry button.
  - In the case of a read or write error, the screen shown at the right may appear and you must exit the program. Contact your ASP or the local Cepheid Technical Support office for assistance, if necessary.

**Disc is not Writable (Already Contains Data)**

Retry Button


**General Write Failure**

Exit the Program

**Figure 2-39. Error writing Xpert Check data CD Screens - Two Examples**


- After test completion, the Exiting screen will appear with the message Remove and discard all Xpert Check cartridges (see Figure 2-40).

Click to exit the program.

Figure 2-40. Exiting Screen

- Use a felt-tip pen to write on the label of the Xpert Check data CD you have just created by writing the date, instrument identification and facility/location of the test performed. See a label example in Figure 2-41.

**Note:** It is the responsibility of the user to properly enter the date, instrument identification, facility name and location on the Data CD label.

Failure to properly label your Data CD may result in lost data and a delay in obtaining your module Xpert Check code test.

Figure 2-41. Data CD Labeling Example

11. You have the option to copy the calibration\_info.gxc data file (located on the Xpert Check data CD just written) and E-mail the data file directly to your ASP or the local Cepheid Technical Support office instead of mailing the CD. If E-mail is not an option, place the Xpert Check CD2 into the provided CD shipping envelope and mail it to your local Authorized Service Provider (ASP) or the local Cepheid Technical Support office for data quality assurance checking and the issuing of your Xpert Check code.
12. Your ASP or the local Cepheid Technical Support office will perform the quality assurance review and, if successful, send back your Xpert Check code either by E-mail or regular mail, depending on what method you have previously set up with them.

**Note**

---

Discard all remaining materials from the kit. DO NOT save unopened kit pouches for later use. DO NOT discard your Software CD. For users who E-mailed their file and have not shipped their data CD: DO NOT discard your Data CD.

---

13. Restart your GeneXpert Dx or Infinity system and computer.

**Note**

---

You can continue to use your system and GeneXpert DX software while awaiting your Xpert Check code.

---

### 2.3.2 Obtaining the Xpert Check Code for Non-Internet Connected Users

**Note**

---

Ensure the system is in the same configuration as when Xpert Check was run (i.e., no software updates or changes have been made and no new GX systems have been moved to or from this computer). In the case of any module servicing and/or replacement that may occur between data collection and application of the Xpert Check Code, new or modified modules will be ignored for the purposes of the Xpert Check testing process.

---

**Note**

---

In the following step, use care in inserting the CD into the DVD drive. Be sure the CD is fully seated in the tray before closing the drive door.

---

1. Exit the GX Dx software.
2. To finish the Xpert Check process, place the Software CD in the DVD drive of the computer connected to the GeneXpert Dx instrument or in the kiosk computer for the Infinity.
3. Click on My Computer, then double-click on the applicable drive letter for your DVD drive. The files located on the CD will then be displayed. Find and double-click the XpertCheck.exe application/shortcut to launch the software.
4. Log in with your GeneXpert Dx or Infinity designated **USER NAME** and **PASSWORD** (see the IMPORTANT note in section 2.2). Also see [Figure 2-4](#) for the Login screen. After entering your login information, click the forward arrow button at the bottom of the screen to advance to the next screen (the Xpert Check Home screen).

**Note**

The user name and password are the same ones you used for the GeneXpert Dx or Xpertise software. If an ASP (FSE) previously performed Xpert Check and is not now on site, the user name and password should have been provided for this step to enter the code. If the user name or password are not now available, contact your ASP or your local cepheid Technical Support office.

5. Click on the Enter Xpert Check Code button. See Figure 2-42. The Enter Xpert Check code screen will appear. See Figure 2-43.


Figure 2-42. Home Screen, showing Enter Xpert Check Code Button


Figure 2-43. Enter Xpert Check code Screen

6. Enter your Xpert Check code as described below.

In this step, there are various ways to enter the Xpert Check code. Your four options are listed below, beginning with the simplest (recommended) method.

- A. Option 1: Use your scanner to input the barcode as follows: First, click on the icon located in the bottom center of the screen. The Xpert Check Code File (Figure 2-44) will appear on your screen. Position your scanner to scan the barcode on the Code form, using care to avoid any reflection on the monitor that may interfere with your scanner. See Figure 2-44 for an example of an Xpert Check Code File.
- B. Option 2: Print a copy of the Xpert Check Code File and place the copy on a flat surface, facing up. Use your scanner to scan the barcode on the printed page. See Figure 2-44 for an example of a Xpert Check Code File.
- C. Option 3: Copy and paste the code string into the Enter code screen from the screen's display. The code string is visible on Figure 2-44.
- D. Option 4: Type in the code string manually using the information on your screen or printed page.

When you have successfully entered the code, click the forward arrow at the bottom of the screen to continue. The Applying Xpert Check code screen will appear. See Figure 2-45.

## Xpert Check Code File

Here is the Xpert Check code for the recent data collection of your modules for the system identified below.

Xpert Check data collection performed on 30 October 2018 15:02:31 PST


GX Instrument Name:	My GeneXpert
Cepheid System ID:	123456
Software Version:	Xh1.4
Data Collected By:	admin1
Institution Name:	Institution1
Laboratory Name:	Lab1
Street Address:	100 Main Street, Suite 202
City:	New York
State/Province:	NY
Postal Code:	10001
Country:	USA
Email:	user@institution.com
Facility Phone Number:	408 400-0000
Extension:	
Mobile:	
ASP Code:	US01

Scan or enter the Xpert Check code to complete the Xpert Check process.


1. Cepheid recommends that system performance should be evaluated annually using Xpert Check.
2. Cepheid declares that the I-CORE modules in the GeneXpert® Instrument were checked using an Xpert Check product. NIST traceable qualification standards are used to control the parameters for the fluorescence standards of concentration, brightness, and spectrum. Cepheid products are manufactured, quantified and controlled under a Quality System compliant with ISO 13485 and QSR requirements.

**Figure 2-44. Xpert Check Code File - Example**

**Figure 2-45. Applying Xpert Check code Screen Example**

- E. After the Xpert Check code has been applied, the Xpert Check Complete screen will appear with the location of the Xpert Check Report displayed in the Xpert Check Data Directory area. Write down the file path and location of the Xpert Check Report file, as shown. See [Figure 2-46](#).

**Figure 2-46. Xpert Check complete Screen**

- F. Click the Review Xpert Check Status button (see [Figure 2-46](#)).
- G. The Xpert Check status screen will appear. See [Figure 2-47](#). In the Xpert Check status screen, the successfully checked modules are indicated by a + symbol on a green module.


Figure 2-47. Xpert Check status Screen

**Note**

If the Xpert Check report on the computer has been deleted, contact your ASP or the local Cepheid Technical Support office for assistance.

7. Identify the generated Xpert Check Report file in the folder C:\GeneXpert\XpertCall\Reports.
8. Identify the generated Xpert Check Summary Report file in the folder C:\GeneXpert\XpertCall\Reports.
  - A. See [Figure 2-48](#) for an example of a Xpert Check Summary Report.

The Xpert Check Summary Report lists the modules that had an unsuccessful test and require retesting or service.

The modules requiring retesting or service are listed by serial number in Table 1 on the form in [Figure 2-48](#). When requesting service, provide these listed serial numbers to your ASP or the local Cepheid Technical Support office.

Gateway information is provided in Table 2 of the form.

## Xpert Check Summary Report

Please Note:

Xpert Check data collection performed on 30 October 2018 09:37:04 PDT  
 All modules that DID NOT pass Xpert Check are listed in Table 1: Modules Requiring Service.  
 Gateway Informations are provided in Table 2.  
 Complete test results for each module are listed in Table 3: Detailed Test Results by Module.

GX Instrument Name:	My 6Color
Cepheid System ID:	12345
Data Collected By:	User1
Institution Name:	Institution1
Laboratory Name:	Laboratory Sunnyvale
Street Address:	123 Main Street
City:	Sunnyvale
State/Province:	CA
Postal Code:	90001
Country:	USA
Email:	User1@Institution1.com
Facility Phone Number:	408-400-XXXX
Extension:	
Mobile:	
ASP Code:	US03

**Table 1: Modules Requiring Service**

Module Serial Number / Location	Module Status
639563/A2	Skipped and Retest required
639565/A1	Requiring Service

**Table 2: Gateway Information**

Gateway Serial Number	MAC Address
804471	00:21:38:00:2E:1B
804470	00:21:38:00:2E:1A

## 2.4 Return System to Normal Operation

**Note**

Return the system to normal operation by following one of the three procedures listed in this section for the GeneXpert Dx, the Infinity-48, the Infinity 48s, or the Infinity-80.

---

### 2.4.1 GeneXpert Dx

Ensure all Xpert Check cartridges and CDs have been removed from the GeneXpert Dx.

1. Restart your GeneXpert system and computer. Follow the instructions in the GeneXpert Dx System Operator Manual.
2. The system will be ready for full operation.

### 2.4.2 Infinity-48

Ensure all Xpert Check cartridges and CDs have been removed from the Infinity-48.

1. Restart the Xpertise software and switch the system from Manual mode back to Automation mode. Follow the instructions in the GeneXpert Infinity System Operator Manual for the Infinity-48.
2. The system will be ready for full operation.

### 2.4.3 Infinity-48s or Infinity-80

Ensure all Xpert Check cartridges and CDs have been removed from the Infinity-48s or Infinity-80.

1. Restart the Xpertise software. Follow the instructions in the GeneXpert Infinity System Operator Manual.
2. The system will be in Automation mode, ready for full operation.

## 2.5 Information Key Screen


Figure 2-49. Information key Screen

### 2.5.1 Reasons to Repeat Xpert Check with a New Cartridge

If the onscreen instructions direct you to retest, repeat the test according to the instructions in [Step B.](#) on [page 2-18.](#)

### 2.5.2 Reasons to Repeat Xpert Check with the Same Cartridge

If software reports that the cartridge film seal was not broken, remove the original cartridge, rescan the cartridge barcode, open the lid, close the lid, and reinsert the cartridge. Restart the Xpert Check procedure for the affected module.

### 2.5.3 Application of Xpert Check Code

Xpert Check is not complete until the Cepheid-supplied Xpert Check code is applied to the system being tested. Upon receipt of the Quality Assurance Xpert Check Code from Cepheid, apply the code to your system using the Xpert Check Software to complete the Xpert Check process.

# Xpert<sup>®</sup> MTB/XDR

**REF GXMTB/XDR-10**

## **Trademark, Patents and Copyright Statements**

Cepheid<sup>®</sup>, the Cepheid logo, GeneXpert<sup>®</sup>, and Xpert<sup>®</sup> are trademarks of Cepheid.  
Windows<sup>®</sup> is a trademark of Microsoft Corporation.

THE PURCHASE OF THIS PRODUCT CONVEYS TO THE BUYER THE NON-TRANSFERABLE RIGHT TO USE IT IN ACCORDANCE WITH THIS PACKAGE INSERT. NO OTHER RIGHTS ARE CONVEYED EXPRESSLY, BY IMPLICATION OR BY ESTOPPEL. FURTHERMORE, NO RIGHTS FOR RESALE ARE CONFERRED WITH THE PURCHASE OF THIS PRODUCT.

**Copyright © Cepheid 2020. All rights reserved.**

Cepheid AB  
Röntgenvägen 5  
SE-171 54 Solna,  
Sweden

# Xpert<sup>®</sup> MTB/XDR

---

For *In Vitro* Diagnostic Use

## 1 Proprietary Name

Xpert<sup>®</sup> MTB/XDR

## 2 Common or Usual Name

Xpert MTB/XDR Assay

## 3 Intended Use

The Xpert MTB/XDR Assay, performed on the GeneXpert Instrument Systems, is a nested real-time polymerase chain reaction (PCR) *in vitro* diagnostic test for the detection of extensively drug resistant (XDR) *Mycobacterium tuberculosis* (MTB) complex DNA in unprocessed sputum samples or concentrated sediments prepared from sputum. In specimens where MTB is detected, the Xpert MTB/XDR Assay can also detect isoniazid (INH) resistance associated mutations in the *katG* and *fabG1* genes, *oxyR-ahpC* intergenic region and *inhA* promoter; ethionamide (ETH) resistance associated with *inhA* promoter mutations only; fluoroquinolone (FLQ) resistance associated mutations in the *gyrA* and *gyrB* quinolone resistance determining regions (QRDR); and second line injectable drug (SLID) associated mutations in the *rrs* gene and the *eis* promoter region.

The Xpert MTB/XDR Assay is intended for use as a reflex test for a specimen (unprocessed sputum or concentrated sputum sediments) that is determined to be MTB positive. This test is intended as an aid in the diagnosis of XDR tuberculosis (TB) when used in conjunction with clinical and other laboratory findings..

## 4 Summary and Explanation

Tuberculosis (TB), a disease caused by *Mycobacterium tuberculosis*, remains one of the deadliest diseases in the world. In 2018, there were an estimated 10 million new cases of TB and about half a million new cases of rifampicin-resistant TB, of which 78% had multidrug-resistant TB (MDR-TB)<sup>1</sup>. MDR-TB, defined as resistance to isoniazid and rifampicin (two of the most effective first lines drugs), continues to be a public health threat and new treatment guidelines calling for rapid drug susceptibility testing are released by the World Health Organization (WHO)<sup>2,3</sup>. Nevertheless, in 2018, the global number of MDR/RR-TB cases notified was still only 39% of the estimated incident cases and the number of people enrolled in treatment was equivalent to 32%<sup>1</sup>. Likewise, there is also a rising concern of undiagnosed and untreated isoniazid-resistant, rifampicin-susceptible TB. Without easy access to INH-resistance testing, countries struggle to identify patients and implement the 2018 WHO treatment recommendations for Hr-TB<sup>4</sup>. The most worrisome cases of TB are caused by MDR MTB strains which acquired additional resistances to fluoroquinolones and any one of the second line injectable drugs, amikacin (AMK), kanamycin (KAN), or capreomycin (CAP). These highly resistant strains are termed extensively drug resistant TB (XDR-TB). XDR-TB is very difficult to treat and can lead to high rates of mortality, especially when an XDR-TB diagnosis is missed and appropriate treatment is delayed<sup>5</sup>.

Culture and phenotypic drug susceptibility testing of MTB are time consuming, and labor-intensive and present a serious biohazard to laboratory workers, resulting in fewer accredited facilities in countries where MTB is endemic<sup>2</sup>. Even when available, culture-based susceptibility testing can take from weeks to months to complete. MTB may also be tested for drug resistance using fast, sensitive, and safer genotypic assays, which detect resistance by identifying mutations known to confer resistance to the first- and second-line drugs in a majority of clinical strains<sup>2</sup>. Genotypic testing approaches that can be reduced to a few manual steps are more amenable for near patient care, which can dramatically expand their availability to medically underserved populations in low and high endemic settings<sup>5</sup>.

## 5 Principle of the Procedure

The Xpert MTB/XDR Assay is an automated *in vitro* diagnostic test for detection of XDR MTB complex DNA and resistance associated mutations. The assay is performed on Cepheid GeneXpert Instrument Systems equipped with GeneXpert 10 color modules.

The GeneXpert Instrument System integrate and automate sample processing, nucleic acid amplification, and detection of the target sequences in samples using nested real-time PCR and melt peak detection. The GeneXpert Instrument Systems consists of an instrument, personal computer, barcode scanner, and preloaded software for running tests on collected samples and viewing the results. The system requires the use of single-use disposable Xpert cartridges that contain target specific polymerase chain reaction (PCR) reagents and hosts the PCR process and melt peak detection. Because the Xpert cartridges are self-contained, risk of cross-contamination between samples is minimized. For a full description of the system, see the *GeneXpert Dx System Operator Manual*.

The Xpert MTB/XDR Assay cartridge includes reagents for the detection of XDR MTB profile and sample processing control (SPC) to control for adequate processing of the target bacteria and to monitor the presence of inhibitor(s) in the PCR reaction. The Probe Check Control (PCC) verifies reagent rehydration, PCR tube filling in the cartridge, probe integrity, and dye stability.

The Xpert MTB/XDR Assay cartridge has all reagents on board, except sample reagent (SR) which requires the user to add the SR to the specimen prior loading the treated specimen into the cartridge. The test is intended to be run as a reflex test for MTB positive samples.

The results are interpreted by the GeneXpert software from measured fluorescent signals and embedded calculation algorithms and are shown in the View Results window in tabular and graphic formats. It also reports if the test is invalid, has encountered an error or produces no result. The Xpert MTB/XDR detects XDR MTB with resistance to INH, ETH, FLQs, and SLIDs directly from unprocessed sputum or from concentrated sediment from sputum in less than 90 minutes.

## 6 Reagents and Instruments

### 6.1 Material Provided

The Xpert MTB/XDR kit contains sufficient reagents to process 10 patient or quality-control specimens. The kit contains the following items:

<b>Xpert MTB/XDR Cartridges with Integrated Reaction Tubes</b>	<b>10 per kit</b>
• Bead 1, Bead 2, Bead 3, Bead 4, and Bead 5 (freeze-dried)	1 of each per cartridge
• Sample Processing Control Bead (freeze-dried)	1 of each per cartridge
• Reagent 1	4.0 mL per cartridge
• Reagent 2	4.0 mL per cartridge
<b>Disposable transfer pipettes</b>	<b>1 bag of 12 per kit</b>
<b>Sample Reagent</b>	<b>10 x 8 mL per bottle</b>
<b>CD</b>	<b>1 per kit</b>
• Assay Definition Files (ADF)	
• Instructions to import ADF into the GeneXpert software	
• Instructions for Use (Package Insert)	

---

**Note** Sample Reagent (SR) can be colorless to yellow to amber. Color may intensify with time, but color has no effect on performance.

---

**Note** Safety Data Sheets (SDS) are available at [www.cepheid.com](http://www.cepheid.com) or [www.cepheidinternational.com](http://www.cepheidinternational.com) under the SUPPORT tab.

---

**Note** The bovine serum albumin (BSA) in the beads within this product was produced and manufactured exclusively from bovine plasma sourced in the United States. No ruminant protein or other animal protein was fed to the animals; the animals passed ante- and post-mortem testing. During processing, there was no mixing of the material with other animal materials.

---

**Note** The transfer pipettes have a single mark representing the minimum volume of treated sample necessary to transfer to the cartridge. Use only for this purpose. All other pipettes must be provided by the laboratory.

---

## 7 Storage and Handling

- Store the Xpert MTB/XDR kit contents at 2–28 °C until expiration data provided on the label.
- Do not open a cartridge lid until you are ready to perform testing.

- Start the test within 2.5 hours of adding SR to the specimen or within 4 hours if stored at 2–8°C
- Do not use reagents or cartridges that have passed the expiration date.
- Do not use a cartridge that has leaked.

## 8 Materials Required but Not Provided

- GeneXpert Dx system: GeneXpert instrument equipped with GeneXpert 10 color modules., computer, barcode scanner, and operator manual
  - For GeneXpert Dx system: Software version 6.2 or higher
  - Printer: If a printer is required, contact Cepheid Sales Representative to arrange for the purchase of a recommended printer.
- Sterile screw-capped sample container
- Disposable gloves
- Labels and/or indelible labeling marker
- Sterile pipettes for sample processing

## 9 Warnings and Precautions

### 9.1 General

- For *In Vitro* Diagnostic Use

- Treat all biological specimens, including used cartridges, as if capable of transmitting infectious agents. Because it is often impossible to know which might be infectious, all biological specimens should be treated with standard precautions.
- Guidelines for specimen handling are available from the U.S. Centers for Disease Control and Prevention<sup>3</sup> and the Clinical and Laboratory Standards Institute.<sup>6,7,8</sup>
- Follow your institution's safety procedures for working with chemicals and handling biological samples.
- Wear protective disposable gloves, laboratory coats and eye protection when handling specimens and reagents. Wash hands thoroughly after handling specimens and test reagents.

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures. If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines<sup>9</sup>.
- Sample Reagent contains sodium hydroxide (pH > 12.5) and isopropanol. Harmful if swallowed (H302), causes severe skin burns and eye damage (H314). Flammable liquid and vapor (H226).
- Performance characteristics of this test have been established with the specimen types listed in the Intended Use Section only. The performance of this assay with other specimen types or samples has not been evaluated.
- Follow your institution's safety procedures for working with chemicals and handling biological samples.

### 9.2 Specimen

- Specimen collection and handling procedures require specific training and guidance.
- Maintain proper storage conditions during specimen transport to ensure the integrity of the specimen (see Section 12, Procedure). Specimen stability under shipping conditions other than those recommended has not been evaluated.
- Reject specimens with obvious food particles or other solid particulates.
- Proper sample collection, storage, and transport are essential for correct results

### 9.3 Assay/Reagent

- Do not substitute Xpert MTB/XDR Assay reagents with other reagents.
- Do not open the Xpert MTB/XDR Assay cartridge lid except when adding sample.
- Do not use a cartridge that has been dropped after removing from the kit or shaken after the cartridge lid has been opened. Shaking or dropping the cartridge after opening the lid may yield false or non-determinate results.
- Do not place the sample ID label on the cartridge lid or on the barcode label.
- Do not use a cartridge that has a damaged reaction tube.
- ② • Each single-use Xpert MTB/XDR Assay cartridge is used to process one test. Do not reuse spent cartridges.
- ② • A single-use disposable pipette is used to transfer one specimen. Do not reuse spent disposable pipettes.
- Do not use a cartridge if it appears wet or if the lid seal appears to have been broken.
- Good laboratory practices, including changing gloves between handling patient specimens, are recommended to avoid contamination of specimens or reagents.
- In the event of a spill of specimens or controls, wear gloves and absorb the spill with paper towels. Then, thoroughly clean the contaminated area with a 1:10 dilution of freshly prepared household chlorine bleach. Final active chlorine concentration should be 0.5% regardless of the household bleach concentration in your country. Allow a minimum of two minutes of contact time. Ensure the work area is dry before using 70% denatured ethanol to remove bleach residue. Allow surface to dry completely before proceeding. Or, follow your institution's standard procedures for a contamination or spill event. For equipment, follow the manufacturer's recommendations for decontamination of equipment.
- The assay has been validated using Cepheid GeneXpert Dx software version 6.2 or higher.

## 10 Chemical Hazards<sup>9,10</sup>

### Sample Reagent:

- Contains Isopropyl Alcohol
- Contains Sodium Hydroxide
- Signal Word: DANGER
- UN GHS Hazard Pictograms: 
- **UN GHS Hazard Statements**
  - Flammable liquid and vapor.
  - Causes severe skin burns and eye damage.
  - Causes severe eye damage.
  - Suspected of causing genetic defects.
  - Suspected of damaging fertility or the unborn child.
  - May cause damage to organs through prolonged or repeated exposure.
- **UN GHS Precautionary Statements**
- **Prevention**
  - Obtain special instructions before use.
  - Do not handle until all safety precautions have been read and understood.
  - Keep away from heat, sparks, open flames and/or hot surfaces. - No smoking.
  - Keep container tightly closed.
  - Do not breath mists, vapours, and/or spray.
  - Wash thoroughly after handling.
  - Wear protective gloves, protective clothing, eye protection, face protection.
  - Use personal protective equipment as required.
- **Response**
  - In case of fire: Use appropriate media for extinction.

- IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- Immediately call a POISON CENTER or doctor/physician.
- IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- Wash contaminated clothing before reuse.
- Specific treatment, see supplemental first aid information.
- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- IF SWALLOWED: Rinse mouth. DO NOT induce vomiting.
- IF exposed or concerned: Get medical advice/attention.
- Get medical advice/attention if you feel unwell.
- **Storage/Disposal**
  - Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

## 11 Specimen Collection, Transport and Storage

Specimens can be collected following the user institution's standard procedures.

Proper specimen collection, storage, and transport are critical to the performance of this test. Specimen stability under shipping and storage conditions other than those listed below have not been evaluated with the Xpert MTB/XDR Assay.

### 11.1 Transport

Specimens are recommended to be transported at 2–8 °C whenever possible.

Unprocessed sputum specimens can be kept at 2–35°C.

### 11.2 Untreated Specimens

Unprocessed sputum specimens can be stored at 2–35°C for 7 days (including shipping time)

Decontaminated/concentrated and resuspended sputum sediment can be stored at 2–8 °C for up to 7 days until testing is performed on the GeneXpert.

When testing unprocessed sputum or decontaminated/concentrated sputum sediment refer to Table 1 below to determine adequate specimen volume.

**Table 1. Required Specimen Volume**

Specimen Type	Minimum Volume for One Test	Maximum sample volume	Specimen to Sample Reagent (SR) Ratio
Sputum sediment	0.5 mL	2.5 mL	1:3 <sup>a</sup>
Unprocessed sputum	1.0 mL	4.0 mL	1:2

a. 1:2 sample to SR ratio should be used with sample volume of 0.7 mL or greater for one test.

### 11.3 Leftover Specimens Treated with SR

Xpert MTB/XDR Assay can be used to test left over SR treated specimen from Xpert MTB/RIF or Xpert MTB/RIF Ultra Assay. However, in such cases, the volume of the leftover SR treated specimen must be ≥2mL and the mix should be stored 2–8 °C for no longer than 4 hours or up to 35 °C for no longer than 2.5 hours.

## 12 Procedure

### 12.1 Procedure for Unprocessed Sputum

**Important** Start the test within 2.5 hours of adding SR to the specimen or within 4 hours if stored at 2–8 °C.

**Note** Reject specimens with obvious food particles or other solid particles.


**Volume Requirements:** ≥1mL of unprocessed sputum is required.

- Carefully open the lid of the leak-proof sputum collection container. See Figure 1.

**Figure 1. Opening the sputum collection container**

- Pour approximately 2 times the volume of the SR to the sputum (2:1 dilution, SR:sputum). See Figure 2.

**Figure 2. Examples of 2:1 dilutions**

**Note** Discard the leftover SR and bottle in an appropriate waste container according to your institution's standard practices.

- Secure the lid on the sample container.
- Shake vigorously 10 to 20 times or vortex for at least 10 seconds.

**Note** One back-and-forth movement is a single shake.

- Incubate for 10 minutes at room temperature, and then shake the specimen vigorously 10 to 20 times or vortex for at least 10 seconds.
- Incubate the sample at room temperature for an additional 5 minutes.

## 12.2 Procedure for Decontaminated Concentrated Sputum Sediments

**Important** Start the test within 2.5 hours of adding SR to the specimen or within 4 hours if stored at 2–8 °C.

**Note** Reject specimens with obvious food particles or other solid particles.

**Volume Requirements:** Sputum sediments prepared according to the NALC-NaOH (N-Acetyl-L-Cysteine–Sodium Hydroxide) method described in Kent and Kubica<sup>11</sup> and re-suspended in 67 mM Phosphate/H<sub>2</sub>O buffer can be tested using the Xpert MTB/XDR Assay. After resuspension, keep at least 0.5 mL of the resuspended sediment for the Xpert MTB/XDR Assay. For all volumes less than 0.7 mL perform steps 1 through 5 to prepare samples. These steps require 3 parts SR to 1 part sediment in order to generate adequate volume for the optimum performance of the assay. If the sample volume is equal to or greater than 0.7 mL, adequate test volume can be produced by adding 2 parts SR to 1 part sediment. In this example 1.4 mL of SR would be added to 0.7 mL sediment. These volumes scale at a ratio of 2 parts SR to 1 part sediment.

1. Transfer 0.5 mL of the total resuspended pellet to a conical, screw-capped tube labeled with the sample and/or patient ID using a transfer pipette.

**Note** Store re-suspended sediments at 2–8 °C if they are not immediately processed. Do not store for more than 7 days.

2. Add 1.5 mL of Sample Reagent (SR) to 0.5 mL of resuspended sediment.
3. Shake vigorously 10 to 20 times or vortex for at least 10 seconds.

**Note** One back-and-forth movement is a single shake.

4. Incubate for 10 minutes at room temperature, and then shake the specimen vigorously 10 to 20 times or vortex for at least 10 seconds.
5. Incubate the sample at room temperature for an additional 5 minutes.

## 12.3 Preparing the Cartridge

**Important** Ensure a module is ready to accept cartridge. Start the test as soon as possible and within 2.5 hours of adding the Sample Reagent-treated sample to the cartridge or within 4 hours if stored at 2–8 °C.

Obtain the following items: Xpert cartridge, transfer pipette (provided), and an appropriately collected and labeled test sample.

1. Remove a cartridge from the package.
2. Inspect the cartridge for damage. If damaged, do not use it.
3. Bring the cartridge to room temperature. Label each Xpert MTB/XDR cartridge with the Sample ID. See Figure 3.

**Figure 3. Write on Side of Cartridge.**

**Note** Write on the side of the cartridge or affix an ID label. Do not put the label on the lid of the cartridge or over the existing 2D barcode on the cartridge.

4. Open the cartridge lid, and then open the sample container.

- Using the provided transfer pipette, aspirate the liquefied sample to the line on the pipette. Do not process the sample further if there is insufficient volume. See Figure 4.

**Figure 4. Aspirating to the line on the pipette**

- Dispense the sample slowly to minimize the risk of aerosol formation. See Figure 5.

**Figure 5. Xpert MTB/XDR Assay Cartridge**

- Close the cartridge lid.

#### 12.4 Starting the Test

---

**Important** Before starting the test, make sure that the Xpert MTB/XDR Assay definition file is imported into the software. This section lists the basic steps of running the test. For detailed instructions, see the *GeneXpert Dx System Operator Manual*.

---

**Note** The steps you follow can be different if the system administrator changed the default workflow of the system.

---

- Turn on the GeneXpert instrument:
  - If using the GeneXpert Dx instrument, first turn on the instrument and then turn on the computer. The GeneXpert Dx software will launch automatically or may require double-clicking the GeneXpert Dx shortcut icon on the Windows® desktop.
- Log on to the GeneXpert Instrument System software using your user name and password.
- In the GeneXpert Dx System window, click **Create Test**. The Create Test window appears.

4. Scan in Patient or Sample ID or type in the Patient or Sample ID. If typing the Sample ID, make sure the Sample ID is typed correctly. The Sample ID is shown on the left side of the View Results window and is associated with the test results.
5. Scan the barcode on the Xpert MTB/XDR Assay cartridge. Using the barcode information, the software automatically fills the boxes for the following fields: Reagent Lot ID, Cartridge S/N, and Expiration Date. See Figure 6.

**Note** If the barcode on the Xpert MTB/XDR cartridge does not scan, then repeat the test with a new cartridge.

**Figure 6. GX Dx Create Test Window**

6. Click **Start Test**. Type your password in the dialog box that appears.
7. For the GeneXpert Dx Instrument:
  - A. Open the instrument module door with the blinking green light and load the cartridge.
  - B. Close the door. The test starts and the green light stops blinking. When the test is finished, the light turns off.
  - C. Wait until the system releases the door lock before opening the module door and removing the cartridge.
8. Dispose of used cartridges in the appropriate specimen waste container according to your institution's standard practices.

### 13 Viewing and Printing Results

This section lists the basic steps for viewing and printing results. For more detailed instructions on how to view and print the results, see the *GeneXpert Dx System Operator Manual*.

- Click the **View Results** icon to view results.
- Upon completion of the test, click the **Report** button of the View Results window to view and/or generate a PDF report file.

## 14 Quality Control

### 14.1 Built-in Quality Controls

#### CONTROL

Each test includes a Sample Processing Control (SPC) and Probe Check Control (PCC).

- **Sample Processing Control (SPC)**— The SPC verifies that sample processing is adequate. Additionally, this control detects sample-associated inhibition of the real-time PCR assay, ensures the PCR reaction conditions (temperature and time) are appropriate for the amplification reaction, and that the PCR reagents are functional. The SPC should be positive in a negative sample and can be negative or positive in a positive sample. The SPC passes if it meets the assigned acceptance criteria.
- **Probe Check Control (PCC)**—Before the start of the PCR reaction, the GeneXpert System measures the fluorescence signal from the probes to monitor bead rehydration, reaction tube filling, probe integrity, and dye stability. The PCC passes if it meets the assigned acceptance criteria.
- **Sample Volume Adequacy (SVA) control**—Before sample processing, the GeneXpert system measures if adequate volume of specimen is present in the sample chamber. If SVA check fails, it implies that adequate volume of sample required for testing has not been added to the sample chamber.

## 15 Interpretation of Results

The GeneXpert Instrument System generates the results from a combination of measured fluorescent signals and melting temperature ( $T_m$ ) values. Mutations and wild type sequences are detected by the GeneXpert System using  $T_m$  values. Susceptibility or resistance determination depends on where the  $T_m$  values fall within the wild type or mutant window respectively for a particular analyte. Positive results for the Xpert MTB/XDR Assay can be MTB DETECTED and all resistance targets are NOT DETECTED or MTB DETECTED and one or more of the resistance targets is DETECTED or MTB DETECTED and/or one or more of the following resistance targets is INDETERMINATE. See Table 2 for a list of possible results for each target.

**Table 2. Possible Test Results for Each Target in the Xpert MTB/XDR Assay**

Drug Class	Result Call
N/A	INVALID/ERROR/NO RESULT
	MTB DETECTED
	MTB NOT DETECTED
Isoniazid	Low INH Resistance DETECTED
	INH Resistance DETECTED
	INH Resistance NOT DETECTED
	INH Resistance INDETERMINATE
Fluoroquinolone	Low FLQ Resistance DETECTED
	FLQ Resistance DETECTED
	FLQ Resistance NOT DETECTED
	FLQ Resistance INDETERMINATE
Amikacin	AMK Resistance DETECTED
	AMK Resistance NOT DETECTED
	AMK Resistance INDETERMINATE
Kanamycin	KAN Resistance DETECTED
	KAN Resistance NOT DETECTED
	KAN Resistance INDETERMINATE

**Table 2. Possible Test Results for Each Target in the Xpert MTB/XDR Assay (Continued)**

Drug Class	Result Call
Capreomycin	CAP Resistance DETECTED
	CAP Resistance NOT DETECTED
	CAP Resistance INDETERMINATE
Ethionamide <sup>a</sup>	ETH Resistance DETECTED
	ETH Resistance NOT DETECTED

a. Ethionamide will not provide an indeterminate by assay design.

Table 3 summarizes the genes targeted by the Xpert MTB/XDR Assay and codon region and nucleotides covered for each of genes interrogated to identify or infer drug resistance.

**Table 3. Drug resistance determining regions targeted interrogated**

Drug	Gene Target	Codon Regions	Nucleotide
Isoniazid	<i>inhA</i> promoter	NA	-1 to -32 intergenic
	<i>katG</i>	311-319	939-957
	<i>fabG1</i>	199-210	597-630
	<i>oxyR- ahpC</i> intergenic region	NA	-5 to -50 intergenic (or -47 to -92) <sup>12,13</sup>
Ethionamide	<i>inhA</i> promoter <sup>a</sup>	NA	-1 to -32 intergenic
Fluoroquinolones	<i>gyrA</i>	87-95	261-285
	<i>gyrB</i>	531-544 (or 493-505) <sup>12,14</sup>	1596-1632
Amikacin, Kanamycin, Capreomycin	<i>rrs</i>	NA	1396-1417
	<i>eis</i> promoter	NA	-6 to -42 intergenic

a. The absence of mutations in the *inhA* promoter region does not exclude ETH resistance. Mutations conferring ETH resistance are reported to be present in genomic regions not targeted by the Xpert MTB/XDR assay.<sup>15</sup>

See Table 4 for examples of possible results and corresponding interpretation. Figure 7 through Figure 15 are examples of possible Xpert MTB/XDR Assay results.

**Table 4. Examples of Xpert MTB/XDR Assay Results and Interpretation**

Result	Interpretation
MTB DETECTED; INH Resistance NOT DETECTED FLQ Resistance NOT DETECTED AMK Resistance NOT DETECTED KAN Resistance NOT DETECTED CAP Resistance NOT DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>Mutations leading to INH, FLQs, AMK, KAN, CAP, or ETH resistance are not detected.</li> <li>SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>Probe Check: PASS. All probe check results pass.</li> </ul>

Table 4. Examples of Xpert MTB/XDR Assay Results and Interpretation (Continued)

Result	Interpretation
MTB DETECTED; INH Resistance DETECTED FLQ Resistance DETECTED AMK Resistance DETECTED KAN Resistance DETECTED CAP Resistance DETECTED ETH Resistance DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations contributing to INH resistance have been detected in one or more of the following genes: <i>katG</i>, <i>fabG1</i>, <i>oxyR-ahpC</i> intergenic region and <i>inhA</i> promoter</li> <li>• Mutations contributing to FLQ resistance have been detected in one or more of the following genes: <i>gyrA</i> and <i>gyrB</i> quinolone resistance determining regions (QRDR)</li> <li>• Mutations contributing to AMK resistance have been detected in one or more of the following genes: <i>rrs</i> gene and <i>eis</i> promoter</li> <li>• Mutations contributing to KAN resistance have been detected in one or more of the following genes: <i>rrs</i> gene and <i>eis</i> promoter</li> <li>• Mutations contributing to CAP resistance have been detected in the following gene: <i>rrs</i> gene</li> <li>• Mutations contributing to ETH resistance have been detected in the following gene: <i>inhA</i> promoter</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB DETECTED; INH Resistance DETECTED FLQ Resistance NOT DETECTED AMK Resistance NOT DETECTED KAN Resistance NOT DETECTED CAP Resistance NOT DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations leading to FLQs, AMK, KAN, CAP, and ETH resistance are not detected.</li> <li>• Mutations contributing to INH resistance have been detected in one or more of the following genes: <i>inhA</i> promoter, <i>katG</i>, <i>fabG1</i> and <i>oxyR-ahpC</i> intergenic region</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB DETECTED; INH Resistance DETECTED FLQ Resistance INDETERMINATE AMK Resistance NOT DETECTED KAN Resistance NOT DETECTED CAP Resistance NOT DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations leading to AMK, KAN, CAP, and ETH resistance are not detected.</li> <li>• Mutations contributing to INH resistance have been detected in one or more of the following genes: <i>inhA</i> promoter, <i>katG</i>, <i>fabG1</i> and <i>oxyR-ahpC</i> intergenic region</li> <li>• Mutations contributing to FLQ resistance could not be determined due to the detection of only WT Tm from one or more probes and missing Tms from one or more probes targeting one or more of the following genes: <i>gyrA</i> or <i>gyrB</i>. "OR" no Tm from any of the probes targeting <i>gyrA</i> and <i>gyrB</i> genes.</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB DETECTED; Low INH Resistance DETECTED FLQ Resistance NOT DETECTED AMK Resistance NOT DETECTED KAN Resistance NOT DETECTED CAP Resistance NOT DETECTED ETH Resistance DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations leading to FLQ, AMK, KAN, and CAP resistance are not detected.</li> <li>• Mutations contributing to low INH resistance have been detected in <i>inhA</i> promoter region</li> <li>• Mutations contributing to ETH resistance have been detected in the <i>inhA</i> promoter region</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>

Table 4. Examples of Xpert MTB/XDR Assay Results and Interpretation (Continued)

Result	Interpretation
MTB DETECTED; INH Resistance NOT DETECTED Low FLQ Resistance DETECTED AMK Resistance NOT DETECTED KAN Resistance NOT DETECTED CAP Resistance NOT DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample; low level FLQ, resistance is detected: <ul style="list-style-type: none"> <li>• Mutations leading to INH, AMK, KAN, CAP and ETH resistance are not detected.</li> <li>• Mutations contributing to low FLQ resistance have been detected in the following genes: <i>gyrA</i></li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB DETECTED; INH Resistance DETECTED FLQ Resistance NOT DETECTED AMK Resistance DETECTED KAN Resistance DETECTED CAP Resistance DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations leading to FLQ and ETH resistance are not detected.</li> <li>• Mutations contributing to INH resistance have been detected in one or more of the following genes: <i>katG</i>, <i>fabG1</i>, <i>oxyR-aphC</i></li> <li>• Mutations contributing to AMK resistance have been detected in one or more of the following genes: <i>rrs</i> gene; <i>eis</i> promoter</li> <li>• Mutations contributing to KAN resistance have been detected in one or more of the following genes: <i>rrs</i> gene; <i>eis</i> promoter</li> <li>• Mutations contributing to CAP resistance have been detected in the following gene: <i>rrs</i> gene</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB DETECTED; INH Resistance DETECTED Low FLQ Resistance DETECTED AMK Resistance NOT DETECTED KAN Resistance DETECTED CAP Resistance NOT DETECTED ETH Resistance NOT DETECTED	The MTB target is present within the sample: <ul style="list-style-type: none"> <li>• Mutations leading to AMK, CAP, and ETH resistance are not detected.</li> <li>• Mutations contributing to INH resistance have been detected in one or more of the following genes: <i>katG</i>, <i>fabG1</i>, <i>oxyR-ahpC</i> intergenic region and <i>inhA</i> promoter</li> <li>• Mutations contributing to Low FLQ resistance have been detected in the following gene: <i>gyrA</i></li> <li>• Mutations contributing to KAN resistance have been detected in the <i>eis</i> promoter region</li> <li>• SPC: NA (not applicable). An SPC signal is not required because MTB amplification can compete with this control.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
MTB NOT DETECTED	The MTB target is not detected within the sample: <ul style="list-style-type: none"> <li>• SPC: PASS. The SPC met the acceptance criteria.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>
INVALID	The presence or absence of MTB cannot be determined. The SPC does not meet the acceptance criteria, the sample was not properly processed, or PCR was inhibited. Repeat the test. See the Retest Procedure section of this document. <ul style="list-style-type: none"> <li>• MTB: INVALID. The presence or absence of MTB DNA cannot be determined.</li> <li>• SPC: FAIL. The MTB target result is negative, and the SPC Cycle Threshold (Ct) is not within valid range.</li> <li>• Probe Check: PASS. All probe check results pass.</li> </ul>

**Table 4. Examples of Xpert MTB/XDR Assay Results and Interpretation (Continued)**

Result	Interpretation
ERROR	<p>The presence or absence of MTB cannot be determined. Repeat the test. See the Retest Procedure section of this document.</p> <ul style="list-style-type: none"> <li>• MTB: NO RESULT</li> <li>• SPC: NO RESULT</li> <li>• Probe Check: FAIL. All or one of the probe check results failed.</li> </ul> <p><b>Note:</b> If the probe check passed, the error may be caused by a system component failure, operator error or cartridge integrity issue.</p>
NO RESULT	<p>The presence or absence of MTB cannot be determined. Repeat the test. See the Retest Procedure section of this document. A NO RESULT indicates that insufficient data was collected. For example, the operator stopped a test that was in progress.</p> <ul style="list-style-type: none"> <li>• MTB: NO RESULT</li> <li>• SPC: NO RESULT</li> <li>• Probe Check: NA (not applicable)</li> </ul>

**Note** The following figures provide representative results including melt peak tab that can be expected with the Xpert MTB/XDR Assay in the GeneXpert Dx Detailed User View. Not all possible combinations of results are shown.

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR IUO		Version	3		
Test Result	<b>MTB DETECTED;</b> INH Resistance NOT DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED					
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.3	292.5				
katG-melt	73.8	107.0				
fabG1-melt	71.5	242.0				
ahpC-melt	68.7	41.3				
gyrA1-melt	76.2	73.9				
gyrA2-melt	70.4	75.8				
gyrA3-melt	71.0	129.8				
gyrB2-melt	69.5	77.8				
rrs-melt	75.0	188.7				
eis-melt	68.5	145.3				
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Figure 7. MTB DETECTED; INH, FLQ, AMK, KAN, CAP, and ETH Resistance NOT DETECTED


Figure 8. MTB DETECTED; INH, FLQ, AMK, KAN, CAP, and ETH Resistance DETECTED

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR IUO Version 3						
Test Result	<div style="background-color: red; color: white; padding: 2px;">MTB DETECTED;</div> <div style="background-color: red; color: white; padding: 2px;">INH Resistance DETECTED;</div> <div style="background-color: green; color: black; padding: 2px;">FLQ Resistance NOT DETECTED;</div> <div style="background-color: green; color: black; padding: 2px;">AMK Resistance NOT DETECTED;</div> <div style="background-color: green; color: black; padding: 2px;">KAN Resistance NOT DETECTED;</div> <div style="background-color: green; color: black; padding: 2px;">CAP Resistance NOT DETECTED;</div> <div style="background-color: green; color: black; padding: 2px;">ETH Resistance NOT DETECTED</div>					
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
	Analyte Name		Melt Peak Temperature		Melt Peak Height	
	inhA-melt		76.6		284.9	
	katG-melt		74.0		105.2	
	fabG1-melt					
	ahpC-melt		69.0		35.4	
	gyrA1-melt		76.6		65.2	
	gyrA2-melt		70.4		64.9	
	gyrA3-melt		71.4		92.2	
	gyrB2-melt		69.7		84.7	
	rrs-melt		75.3		146.8	
	eis-melt		68.7		124.2	
	inhA-mut melt					
	katG-mut melt					
	fabG1-mut melt		75.9		178.0	
	ahpC-mut melt					
	gyrA1-mutA melt					
	gyrA1-mutB melt					
	gyrA1-mutC melt					
	gyrA2-mutA melt					
	gyrA2-mutB melt					
	gyrA3-mutA melt					
	gyrA3-mutB melt					
	gyrA3-mutC melt					
	gyrB2-mut melt					
	rrs-mut melt					
	eis-mutA melt					
	eis-mutB melt					

Figure 9. MTB DETECTED; INH Resistance DETECTED


Figure 10. MTB DETECTED; INH and KAN Resistance DETECTED; AMK and CAP INDETERMINATE

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR IUO	Version	3			
Test Result	<p>MTB DETECTED;            INH Resistance DETECTED;            Low FLQ Resistance DETECTED;            AMK Resistance NOT DETECTED;            KAN Resistance NOT DETECTED;            CAP Resistance NOT DETECTED;            ETH Resistance NOT DETECTED</p>					
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.5	313.1				
katG-melt						
fabG1-melt	71.7	211.5				
ahpC-melt	69.0	47.2				
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt	69.6	81.1				
rrs-melt	75.2	248.1				
eis-melt	68.8	158.2				
inhA-mut melt						
katG-mut melt	68.4	184.6				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt	72.3	125.0				
gyrA1-mutC melt						
gyrA2-mutA melt	76.0	207.9				
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt	76.5	128.0				
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Figure 11. MTB DETECTED; INH and Low FLQ Resistance DETECTED


Figure 12. MTB DETECTED; INH, FLQ, AMK, and KAN Resistance DETECTED

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR		Version 3			
Test Result	MTB NOT DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt						
katG-melt						
fabG1-melt						
ahpC-melt						
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt						
rrs-melt						
eis-melt						
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Figure 13. MTB NOT DETECTED

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR IUO		Version	3		
Test Result	INVALID					
For Investigational Use Only.						

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.8	102.1				
katG-melt						
fabG1-melt	71.7	53.1				
ahpC-melt	69.1	34.9				
gyrA1-melt	76.6	71.4				
gyrA2-melt						
gyrA3-melt	71.5	40.7				
gyrB2-melt	70.2	38.9				
rrs-melt						
eis-melt	68.6	109.4				
inhA-mut melt						
katG-mut melt	68.5	49.4				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Figure 14. INVALID

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR IUO		Version 3			
Test Result	<b>ERROR</b>					
For Investigational Use Only.						

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt						
katG-melt						
fabG1-melt						
ahpC-melt						
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt						
rrs-melt						
eis-melt						
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Figure 15. ERROR

## 16 Retests

### 16.1 Reasons to Repeat the Test

If any test results mentioned below occur, repeat the test according to the instructions in Section 16.2, Retest Procedure.

- An **INVALID** result indicates that the SPC failed. The sample was not properly processed, or PCR is inhibited or the sample was not properly collected.
- An **ERROR** result could be due to, but not limited to, Probe Check Control failed or the maximum pressure limits were exceeded.
- A **NO RESULT** indicates that insufficient data were collected. For example, the operator stopped a test that was in progress or a power failure occurred.
- An **INDETERMINATE** result indicates that resistance to a given drug could not definitively be concluded based on the assay algorithm (see Limitations for further explanations). Retesting with a different sample may or may not lead to a different result.

## 16.2 Retest Procedure

For retest, use a new cartridge (do not re-use a cartridge). If you have leftover sputum (should be  $\geq 1.0$  ml) or reconstituted sediment (should be  $\geq 0.5$  ml), always use new SR to decontaminate and liquefy the sputum before running the assay. Follow sample processing instructions according to Section 12.1, Procedure for Unprocessed Sputum or Section 12.2, Procedure for Decontaminated Concentrated Sputum Sediments.

If sufficient leftover SR-treated sample is available that has been stored for no longer than 2.5 hours up to 35 °C or has been stored no longer than 4 hours at 2–8 °C of the initial addition of SR to the sample, the leftover SR treated sample can be processed using a new cartridge. When retesting, always use a new cartridge and start the test within 30 minutes of adding processed sample to cartridge. See Section 12.3, Preparing the Cartridge.

## 17 Limitations

- The performance of the Xpert MTB/XDR Assay was validated using the procedures provided in this package insert. Modifications to XDR test procedure should be interpreted in conjunction with other laboratory and clinical data available to the clinician.
- The performance of the Xpert MTB/XDR Assay is dependent on operator proficiency and adherence to assay procedures. Assay procedural errors may cause false positive or false negative results. All device operators should have appropriate device and assay training.
- Because the detection of MTB complex DNA is dependent on the number of organisms present in the sample, reliable assay results are dependent on proper specimen collection, handling, and storage. Erroneous test results might occur from improper specimen collection, failure to follow the recommended sample collection procedure, handling or storage, technical error, sample mix-up, or an insufficient concentration of starting material. Careful compliance to the instructions in this insert is necessary to avoid erroneous results.
- Test results might be affected by antecedent or concurrent antibiotic therapy. Therefore, therapeutic success or failure cannot be assessed using this test because DNA might persist following tuberculosis therapy.
- A positive test result does not necessarily indicate the presence of viable organisms. It is however, presumptive for the presence of MTB complex DNA including mutations associated with INH, FLQ, AMK, KAN, CAP and ETH resistance.
- Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown XDR-MTB strains resulting in a drug-sensitive result.
- The Xpert MTB/XDR Assay does not provide confirmation of susceptibility to INH, FLQ, AMK, KAN, CAP and ETH since mechanisms of resistance other than those detected by the assay may exist that may be associated with a lack of clinical response to treatment.
- Testing of blood, cerebral spinal fluid (CSF), gastric aspirate, stool, tissue, urine has not been evaluated for use in Xpert MTB/XDR Assay.
- Although induced sputum specimens were not included in the clinical performance evaluation of the Xpert MTB/XDR Assay, isotonic or hypertonic solutions, bronchodilators, and inhaled bronchodilators commonly used in the collection of induced sputum were tested and do not interfere with the assay. Saline induction may result in insufficient number the organisms recovered and could affect detection of *M. tuberculosis*.
- Concentrated sputum sediments used in the performance evaluation of the Xpert MTB/XDR Assay were prepared following the NALC-NaOH method described in Kent and Kubica<sup>11</sup>. Use of other methods of sediment preparation may alter the performance of the test.
- A negative test does not exclude the possibility of isolating MTB complex DNA from the sputum sample. The Xpert MTB/XDR Assay may be used in conjunction with mycobacterial culture to address the risk of false negative results and to recover the organism for further characterization and susceptibility testing.
- Specimens with “MTB Trace DETECTED” results when tested with the Xpert MTB/RIF Ultra Assay are expected to be below the Limit of Detection of the MTB/XDR Assay and are not recommended for testing with the Xpert MTB/XDR Assay.
- The Xpert MTB/XDR Assay by design does not differentiate between the species of the MTB-complex (i.e., *MTB*, *M. bovis*, *M. africanum*, *M. canettii*, *M. microti*, *M. caprae*, *M. pinnipedi*, *M. mungi*, and *M. orygis*). In addition, culture must also be performed to determine if an NTM strain is present in addition to MTB-complex.
- Lower sensitivity has been reported in the literature in pediatric patients due to the diffuse nature of MTB infection in the lungs of this patient group, and difficulties encountered in obtaining adequate specimens<sup>16,17</sup>.

- A trained health care professional should interpret assay results in conjunction with the patient’s medical history, clinical signs and symptoms, and the results of other diagnostic tests.
- Mixed infections with MTB and *M. marinum* may result in “INDETERMINATE” results for FLQ at >10<sup>4</sup> CFU/mL of *M. marinum* in presence of ≤408 CFU/mL of MTB.
- In rare instances, the *rrs* primers and probes may cross-react with environmental microbes or sputum microflora which may result in “INDETERMINATE” results for AMK, KAN and CAP.
- The Xpert MTB/XDR Assay determines ETH resistance associated only with mutations in the *inhA* promoter region. The absence of mutations in the *inhA* promoter region does not exclude ETH resistance. Mutations conferring ETH resistance are reported to be present in genomic regions not targeted by the Xpert MTB/XDR assay.<sup>15</sup>
- The association of mutations in the *oxyR-ahpC* and *gyrB* genes with INH and FLQ resistance respectively, has not yet been conclusively established; however published studies have reported these mutations are found in INH and FLQ resistance strains<sup>18,19</sup>.
- Presence of deletions or rare mutations in any of the target genes could lead to “INDETERMINATE” results for a particular drug.
- In case of samples with a mixed population of both susceptible and resistant strains, there is a likelihood that the Xpert MTB/XDR Assay may not detect the mutation, if the resistant population is present at undetectable levels for the assay.
- In samples with very low bacterial load or a mixture of both susceptible and resistant strains, the Xpert MTB/XDR Assay may not reliably distinguish between low and high FLQ resistance.

## 18 Clinical Performance

A blinded clinical study was conducted to evaluate the performance of the Xpert MTB/XDR test relative to microbiological and molecular reference methods. i.e. phenotypic drug susceptibility (pDST) testing and sequencing for the detection of drug resistance to INH, ETH, FLQs and SLID (AMK, KAN and CAP). In addition, the clinical performance of the Xpert MTB/XDR test was compared to the Xpert MTB/RIF or the Xpert MTB/RIF Ultra test for the detection of MTB. Two sites with known high prevalence for MDR and XDR TB provided frozen archived unprocessed sputum or concentrated sputum sediment specimens known to be positive or negative by MTB culture

Table 5 shows the sensitivity and specificity of the Xpert MTB/XDR test relative to pDST for drug resistance. The sensitivity was >90% for INH, FLQ and AMK, >85% for KAN and CAP, and >64% for ETH; the specificity was >98% for all drugs.

**Table 5. Xpert MTB/XDR Test vs. pDST for Drug Resistance**

Drugs	N	TP	FN	TN	FP	Sensitivity (%)	95%CI	Specificity (%)	95% CI
INH	478	244	23	209	2	91.4	87.4 – 94.2	99.1	96.6 – 99.7
FLQ	417	148	11	254	4	93.1	88.0 - 96.1	98.5	96.1 – 99.4
AMK	405	79	7	317	2	91.9	84.1 – 96.0	99.4	97.7 – 99.8
KAN	343	58	8	276	1	87.9	77.9 – 93.7	99.6	98.0 – 99.9
CAP	167	21	4	142	0	84.0	65.3 – 93.6	100.0	97.4 – 100.0
ETH	230	75	41	112	2	64.7	55.6 – 72.8	98.3	93.8 – 99.5

Table 6 shows the sensitivity and specificity of the Xpert MTB/XDR test relative to sequencing for drug resistance. The sensitivity was >93% for FLQ and greater than 96% for INH, AMK, KAN, CAP and ETH; the specificity was 100.0% for all drugs listed in the table except INH which was 98.7%.

**Table 6. Xpert MTB/XDR Test vs. Sequencing for Drug Resistance**

Drug	N	TP	FN	TN	FP	Sensitivity (%)	95%CI	Specificity (%)	95%CI
INH	471	241	3	224	3	98.8	96.5 - 99.6	98.7	96.2 - 99.5
FLQ	469	152	11	306	0	93.3	88.3 - 96.2	100.0	98.8 – 100.0
AMK	463	81	3	379	0	96.4	90.0 - 98.8	100.0	99.0 – 100.0
KAN	463	88	3	372	0	96.7	90.8 - 98.9	100.0	99.0 – 100.0
CAP	463	78	3	382	0	96.3	89.7 - 98.7	100.0	99.0 – 100.0
ETH	473	104	3	366	0	97.2	92.1 – 99.0	100.0	99.0 – 100.0

Table 7 shows the positive percent agreement (PPA) and the negative percent agreement (NPA) of the Xpert MTB/XDR test relative to the Xpert MTB/RIF test for MTB detection to be 98.9% and 93.8%, respectively.

**Table 7. Xpert MTB/XDR Test vs. Xpert MTB/RIF Test for MTB Detection**

		Xpert MTB/RIF Test		
		MTB Detected	MTB Not Detected	Total
Xpert MTB/XDR Test	MTB Detected	273	2a	275
	MTB Not Detected	3b	30	33
	Total	276	32	308
		PPA	98.9% (95%CI: 96.9-99.6)	
		NPA	93.8% (95%CI: 79.9-98.3)	

Table 8 shows the PPA and NPA of the Xpert MTB/XDR test relative to the Xpert MTB/RIF Ultra test for MTB detection to be 99.5% and 100.0%, respectively.

**Table 8. Xpert MTB/XDR Test vs. Xpert MTB/RIF Ultra Test for MTB Detection**

		Xpert MTB/RIF Ultra Test		
		MTB Detected	MTB Not Detected	Total
Xpert MTB/XDR Test	MTB Detected	207	0	207
	MTB Not Detected	1 <sup>a</sup>	14	15
	Total	208	14	222
		PPA	99.5% (95%CI: 97.3-99.9)	
		NPA	100.0% (95%CI: 78.5-100.0)	

a. The Xpert MTB/RIF Ultra result was MTB Detected Trace.

Of the 531 Xpert MTB/XDR runs performed in conjunction with this study, 15 gave non-determinate (“Error”, “Invalid”, or “No Result”) results on the first attempt. Upon retest of these 14 specimens one result remained non-determinate. The initial non-determinate rate was 2.8% (15/531) and the overall non-determinate rate was 0.2% (1/531).

## 19 Analytical Performance

### 19.1 Analytical Sensitivity (Limit of Detection)

Studies were performed to determine the analytical limit of detection (LoD) of the Xpert MTB/XDR Assay with two lots of reagents across three testing days. An MTB positive result is based on the detection of the single copy *inhA* target. The higher LoD observed per strain and per lot as determined by probit analysis was selected for verification. Verification of the estimated LoD claim was performed on one reagent lot across a minimum of three testing days. The LoD was established using a representative MTBC member, *Mycobacterium bovis* BCG (*Bacille Calmette-Guerin*) spiked into a MTB negative, unprocessed sputum and into a MTB negative, decontaminated/concentrated sputum sediment.

The LoD is the lowest concentration reported in CFU/mL that can be reproducibly distinguished from negative samples with  $\geq$  95% confidence. Replicates of 20 were evaluated at five to eight concentrations with two different reagent lots over the 3 days and the LoD was determined using Probit analysis.

The higher LoD observed for each specimen type and lot as determined by probit analysis was selected for verification. Verification of the estimated LoD claim was performed on one reagent lot across a minimum of three testing days with a claim based on a minimum of 19 of 20 positive replicates. The LoD point estimates in CFU/mL are provided in Table 9.

**Table 9. Analytical Sensitivity (Limit of Detection)**

Specimen Type	LoD Point Estimate, CFU/mL
Unprocessed Sputum	136
Sediment	86

## 19.2 Analytical Specificity (Exclusivity)

The analytical specificity of the Xpert MTB/XDR Assay was evaluated by testing a panel of panel of 57 organisms consisting of 21 bacteria, 1 fungus, 7 viruses and 28 Non- tuberculous mycobacteria (NTM) representing common respiratory pathogens or those potentially encountered in the respiratory tract and/or oropharyngeal flora. Three replicates of each bacterial and yeast strain were tested at concentrations of  $\geq 1 \times 10^6$  CFU/mL. All viruses were tested at  $\geq 1 \times 10^5$  (Tissue Culture Infectious Dose) TCID<sub>50</sub>/ml. DNA or RNA were tested for 2 bacterial and 1 fungal strain at concentrations of  $\geq 10^6$  copies/ml, as whole organisms were not available or could not be accessed due to biosafety restrictions. Three replicates of each virus were tested at concentrations of  $\geq 1 \times 10^5$  TCID<sub>50</sub>/mL. The analytical specificity was 100%. The organisms tested are listed in Table 10, Table 11, and Table 12. None of the organisms tested resulted in cross-reactivity with the MTB detection probe generating “MTB NOT DETECTED” result for all the organisms and for all the replicates tested. The table below lists the organisms tested for the analytical specificity assay. *Aspergillus fumigatus* was analytically tested and showed no interference or cross reactivity. Cross reactivity with any other fungal species is not evident by *in silico* analysis.

**Table 10. Analytical Specificity of the Xpert MTB/XDR Assay (Bacteria/Fungi)**

Organism
Acinetobacter baumannii
Chlamydophila pneumoniae <sup>a</sup>
Citrobacter freundii
Corynebacterium xerosis
Enterobacter cloacae
Escherichia coli
Haemophilus influenzae
Klebsiella pneumoniae
Moraxella catarrhalis
Neisseria meningitidis <sup>a</sup>
Neisseria mucosa
Nocardia asteroides
Pseudomonas aeruginosa
Staphylococcus aureus
Staphylococcus epidermidis
Stenotrophomonas maltophilia
Streptococcus agalactiae
Streptococcus mitis
Streptococcus mutans
Streptococcus pneumoniae
Streptococcus pyogenes
Aspergillus fumigatus <sup>a</sup>

a. Genomic DNA

**Table 11. Analytical Specificity of the Xpert MTB/XDR Assay (Viruses)**

Organism
Coronavirus 229E
Human metapneumovirus (hMPV) 16 Type A1
Parainfluenza Virus Type 1
Parainfluenza Virus Type 2
Parainfluenza Virus Type 3
Respiratory Syncytial Virus
Rhinovirus 1A

**Table 12. Analytical Specificity of the Xpert MTB/XDR Assay (NTM)**

Organism
Mycobacterium asiaticum
Mycobacterium avium NJH
Mycobacterium celatum
Mycobacterium chelonae
Mycobacterium flavescens
Mycobacterium fortuitum subsp. Fortuitum
Mycobacterium gastrii
Mycobacterium gordonae
Mycobacterium gordonae
Mycobacterium gordonae
Mycobacterium genavense
Mycobacterium haemophilum
Mycobacterium malmoense
Mycobacterium marinum
Mycobacterium phlei
Mycobacterium scrofulaceum
Mycobacterium simiae
Mycobacterium szulgai
Mycobacterium terrae
Mycobacterium thermoresistibile
Mycobacterium triviale
Mycobacterium vaccae
Mycobacterium xenopi
Mycobacterium avium
Mycobacterium intracellulare
Mycobacterium abscessus
Mycobacterium marinum
Mycobacterium kansasii

### 19.3 Analytical Reactivity (Inclusivity)

The analytical reactivity (inclusivity) of the Xpert MTB/XDR Assay was evaluated using a phylogenetically diverse panel consisting of susceptible and drug resistant MTB strains to evaluate the accuracy of the drug susceptibility results of the assay. The panel of twenty-two (22) MTB-complex (MTBC)-strains included eight (8) drug susceptible strains with wild-type target genes (Table 13) and fourteen (14) well characterized drug resistant strains (Table 14). All strains were tested in triplicate at concentrations at or near 3 X LoD of the *inhA* promoter target. The copy number tested for genomic DNA lysates was based on a fluorescent dye binding assay specific for double-stranded DNA (dsDNA).

The drug susceptible strains were tested and include five strains of MTB (AR2, GD139, AH1, HR36, H37Rv) and three MTB-complex mycobacterial species (*M. bovis*, *M. canettii* and *M. microti*). The MTB strains were selected to broadly represent the range of genetic diversity and include one representative from each of the major phylogenetic lineages based on SNP-cluster groups (SCGs)<sup>20</sup>.

The 14 drug resistant MTB strains were tested using genomic DNA lysates from well characterized specimens which contain 16 clinically significant canonical mutations with at least one of each of the eight regions targeted by the assay. These mutations are commonly present in multi-drug resistant or extensively drug resistant strains of MTB worldwide with the exception of a mutation in the *gyrB* gene.

Table 13 summarizes the results with drug susceptible strains showing number of correct results for each of the individual analytes in the assay. All panel members generated “MTB DETECTED; RESISTANCE NOT DETECTED” The Xpert MTB/XDR Assay correctly identified all replicates of the strains tested near the limit of detection with wild type results for all probes except *oxyR-ahpC*. Since the *oxyR-ahpC* target has a higher LoD than the other targets in the assay, some replicates tested did not yield Tm results.

The results in Table 14 shows the assay also correctly identified expected resistance mutations in all 14 strains resistant to Isoniazid with mutations in *inhA* promoter, *katG* and *oxyR-ahpC* intergenic region; SLIDs resistance with mutations *rrs* and *eis* promoter region; and FLQ resistance with mutations in *gyrA*.

**Table 13. Analytical Reactivity (Inclusivity) for Drug Susceptible Strains**

Sample	Strain Lineage	<i>inhA</i>	<i>katG</i>	<i>fabG1</i>	<i>ahpC</i> <sup>a</sup>	<i>gyrA1</i>	<i>gyrA2</i>	<i>gyrA3</i>	<i>gyrB2</i>	<i>rrs</i>	<i>eis</i>
( <i>M.bovis</i> BCG)	1	PASS	PASS	PASS	FAIL	PASS	PASS	PASS	PASS	PASS	PASS
<i>M.bovis</i>	1	PASS	PASS	PASS	FAIL	PASS	PASS	PASS	PASS	PASS	PASS
<i>MTB</i> (AR2)	2	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS
<i>MTB</i> (GD139)	3	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS
<i>MTB</i> (AH1)	4	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS
<i>MTB</i> (HR36)	5	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS
<i>MTB</i> (HR37Rv)	4	PASS	PASS	PASS	FAIL	PASS	PASS	PASS	PASS	PASS	PASS
<i>M.canetti</i>	Not assigned	PASS	PASS	PASS	FAIL	PASS	PASS	PASS	PASS	PASS	PASS
<i>M.microti</i>	Not assigned	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

a. The LoD for *oxyR-ahpC* is higher than that of *inhA* used for determination of MTB positivity. "PASS" indicates all the replicates tested generated the expected wild type Tm; "FAIL" indicates at least one or more replicates generated no Tm values.

**Table 14. Analytical Reactivity (Inclusivity) for Drug Resistant Strains (# positive results / total tested)**

Strain ID	Gene	Expected Mutation	MTB Detected	Mutant Probe Tm Detected (# positive/tested)	Correct RESISTANCE DETECTED Calls (# positive/tested)
Clinical	<i>gyrA</i>	GAC 94 TAC	3 / 3	<i>gyrA1</i> -MutB (3/3); <i>gyrA3</i> -MutC (3/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>fabG1</i>	G609A		<i>fabG1</i> Mut (3/3)	INH [3/3]
Clinical	<i>gyrA</i>	GGC 88 GCC, GCG 90 GTG, TCG 91 CCG	3 / 3	<i>gyrA1</i> -MutB (2/3), <sup>a</sup> <i>gyrA1</i> -MutC (2/3), <i>gyrA2</i> -MutA (3/3), <i>gyrA3</i> -MutB (1/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>rrs</i>	A1410G		<i>rrs</i> -Mut (3/3)	AMK, CAP, KAN [3/3]
Clinical	<i>gyrA</i>	GAC 94 GGC	3 / 3	<i>gyrA3</i> -MutB (3/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>rrs</i>	A1410G		<i>rrs</i> -Mut (3/3)	AMK, CAP, KAN [3/3]
14-14194	<i>gyrA</i>	GAC 94 GCC	3 / 3	<i>gyrA1</i> -MutA, <i>gyrA2</i> -MutA	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>inhA</i>	C -15 T		<i>inhA</i> -Mut (3/3)	INH, ETH [3/3]
15-14175	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>eis</i>	-10G/A		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]
15-14191	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>eis</i>	-10G/A		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]
16-05612	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>inhA</i>	C -15 T		<i>inhA</i> -Mut (3/3)	INH, ETH [3/3]
	<i>eis</i>	-12C/T		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]

Table 14. Analytical Reactivity (Inclusivity) for Drug Resistant Strains (# positive results / total tested) (Continued)

Strain ID	Gene	Expected Mutation	MTB Detected	Mutant Probe Tm Detected (# positive/tested)	Correct RESISTANCE DETECTED Calls (# positive/tested)
16-05613	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
	eis	-12C/T		eis-Mut (3/3)	AMK, KAN [3/3]
14-13764	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	ahpC	-48G/A		ahpC-Mut (3/3)	INH [3/3]
14-13806	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	ahpC	-48G/A		ahpC-Mut (3/3)	INH [3/3]
Clinical	gyrA	GCG 90 GTG, GAC 94 GGC	3 / 3	gyrA3-MutB (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH [3/3]
	ahpC	G-6A		ahpC (2/3) <sup>b</sup>	INH [3/3]
Clinical	katG	AGC 315 ACC	3 / 3	katG Mut (3/3)	INH [3/3]
Clinical	gyrB2	ACC 539 AAC	3 / 3	gyrB2 WT <sup>c</sup>	*No resistance detected [0/3]
	rrs	A1410G		rrs-Mut (3/3)	AMK, CAP, KAN [3/3]
	gyrA	GCG 90 GTG		gyrA1 MuB (3/3), gyrA2 MutA (3/3), gyrA3 MutB (3/3)	FLQ [3/3]
	ahpC	g -6 a		ahpC Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA Mut (3/3)	INH, ETH [3/3]
Clinical	gyrA	TCG 91 CCG	3 / 3	gyrA1-MutB (3/3), gyrA2-MutA (3/3), gyrA3-MutC (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA Mut (3/3)	INH, ETH [3/3]

- This sample containing three different mutations in the *gyrA* gene did not generate mutant Tms for all the three *gyrA* probes all the time. However, for the correct resistance call to be made, at least one probe needs to generate a mutant Tm, call was correctly made for all the replicates, since at least one *gyrA* probe always generated at least one mutant Tm when tested.
- This sample is a *katG* / *ahpC* double mutant. The replicate with a missed *ahpC* mutant Tm was called INH-R due to the presence of the *katG* mutation, which was detected by the assay.
- This specific mutation is not detected by the assay. However, there is limited clinical evidence that this mutation may actually contribute to FLQ resistance (Low confidence mutation for FLQ-Resistance).

#### 19.4 Interfering Substances Study

Performance of the Xpert MTB/XDR Assay was evaluated in the presence of 36 potentially interfering substances that may be present in the sputum. Potentially interfering substance classes include endogenous substances that are may be present in the specimen and exogenous substances that might be introduced into the specimen. Isotonic or hypertonic solutions, bronchodilators, and inhaled bronchodilators commonly used in the collection of induced sputum were tested and do not interfere with the assay. Saline induction may result in insufficient number the organisms recovered and could affect detection of *M. tuberculosis*.

The substances tested are listed in Table 15 with active ingredients and concentrations tested shown. Negative samples (n = 8) were tested per each substance to determine the effect on the performance of the sample processing control (SPC). Positive samples (n = 8) *Mycobacterium bovis*, *Bacille Calmette-Guerin (BCG)* spiked at 3x the analytical Limit of Detection for TB positivity were tested per substance. All substances were tested in MTB-negative pooled human sputum background included in this study. All positive and negative replicates were identified correctly using the Xpert MTB/XDR Assay, except for Zicam gel (50% w/v; resulted in "MTB NOT DETECTED" in 11.1% of the replicates tested).

Table 15. Potentially Interfering Substances in the Xpert MTB/XDR Assay

Substance/Class	Description / Active Ingredient	Concentration Tested
Blood (human)	Blood 5% (v/v)	5% (v/v)
Human DNA/Cells	HELA 229 cell line	10 <sup>6</sup> cells/mL
White Blood Cells (human)	WBC/Pus matrix (30% buffy coat; 30% plasma; 40% PBS)^	100% (v/v)
Antimycotic; Antibiotic	Nystatin 500KU (100%)	20% (v/v)
Germicidal Mouthwash	Chlorhexidine gluconate (0.12%) oral rinse, USP	20% (v/v)
Specimen Processing Reagents	Cetylpyridinium chloride, 1% in 2% NaCl	0.5% (v/v) in 1% NaCl
Specimen Processing Reagents	Cetylpyridinium chloride, 1% in 2% NALC	0.5% (v/v) in 1% NALC
Specimen Processing Reagents	Cetylpyridinium chloride, 1% in 2% NALC plus 25 mM Citrate	0.5% (v/v) in 1% NALC plus 12.5mM Citrate
Gastric Acid	pH 3 to 4 solution in water, neutralized with sodium bicarbonate	100% (v/v)
Anesthetics (endotracheal intubation)	Lidocaine HCl 4%	4% (v/v)
Nebulizing solutions	NaCl 5% (w/v)	5% (w/v)
Mucin	Mucin 5% (w/v)	5% (w/v)
Antibacterial, systemic	Levofloxacin 25 mg/mL	5 mg/mL
Nasal corticosteroids	Fluticasone 500 mcg/spray	5 µg/mL;
Inhaled bronchodilators	Albuterol Sulfate (2 mg/5mL)	100 µg/mL
Oral anesthetics	Orajel (20% Benzocaine)	5% (w/v)
Anti-viral drugs	Acyclovir	50 µg/mL
Antibiotic, nasal ointment	Neosporin (400U Bactracin, 3.5mg Neomycin, 5000U Polymyxin B)	5% (w/v)
Tobacco	Nicogel 40% tobacco extract	0.5%
Anti-tuberculosis drugs	Streptomycin 1mg/mL	25µg/mL
Anti-tuberculosis drugs	Ethambutol 1mg/mL	50 µg/mL
Anti-tuberculosis drugs	Isoniazid 50mg/5ml	50 µg/mL
Oral expectorants	Guaifenesin (400mg/tablet)	5 mg/mL
Anti-tuberculosis drugs	Pyrazinamide (500mg/tablet)	100 µg/mL
Nasal gel (Homeopathic)	Zicam gel	50% (w/v)
		20% (w/v)
Nasal spray	Phenylephrine 1%	0.5% (v/v)
Anti-tuberculosis drugs	Rifampicin (300mg/tablet)	25 µg/mL
Allergy relief medicine (Homeopathic)	100% Pure Tea tree oil (<5% Cineole, >35% Terpinen-4-01)	0.5% (v/v)
Nebulizing solutions	Pentamidine isethionate	300ng/ml
Anti-tuberculosis drugs	Amoxicillin	25 µg/mL
Bronchodilator	Epinephrine	1mg/mL
Anti-tuberculosis drugs	Amikacin	70ug/ml
Anti-tuberculosis drugs	Capreomycin	50ug/ml
Anti-tuberculosis drugs	Kanamycin	50ug/ml
Anti-tuberculosis drugs	Ethionamide	50ug/ml
Flu Mist Qual Nasal	Influenza Virus Vaccine Live-nasal	5%

### 19.5 Carry-over Contamination Study

A study was conducted to demonstrate that carry-over, cross contamination does not occur when using the single-use, self-contained Xpert MTB/XDR cartridges. The study consisted of processing a negative sample immediately following processing a high concentration of *Mycobacterium bovis-Bacille Calmette-Guerin* (BCG) at  $1 \times 10^{+6}$  CFU/mL in human sputum the same Gene Xpert module. This testing scheme was repeated at least 20 times in two GeneXpert modules producing a total of 41 runs resulting in 20 positives and 21 negatives per module.

All 20 positive samples were correctly reported as MTB DETECTED; INH Resistance NOT DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED. All 21 negative samples were correctly reported as MTB NOT DETECTED. Under the conditions of this study, there was no evidence of any carry-over contamination when testing with very high positive BCG sample at the concentration of  $1.0 \times 10^{+6}$  CFU/mL.

### 19.6 Competitive Interference Study

Competitive interference of the assay caused by the presence of high concentrations of high concentrations of non-tuberculous Mycobacteria (NTM) on the detection of low levels of MTB in the Xpert MTB/XDR Assay was evaluated by testing the representative member of the MTBC. BCG at  $\sim 3 \times$  LoD (411 CFU/mL) in the presence of different NTM strains at  $1 \times 10E+06$  CFU/mL concentration in a background of negative control buffer. MTB positivity is based on detection of *inhA* promoter valid melt peak height and melt peak temperature. Resistance detection is based on valid mut melt peak height and mut melt peak temperature for individual analytes (*inhA*, *katG*, *gyrA1*, *gyrA2*, *gyrA3*, *gyrB2* and *eis*). *oxyR-ahpC* and *fabG1* analytes were excluded due to lower sensitivity and *rrs* was excluded due to known interference with microflora. All samples containing BCG should have results as MTB DETECTED; INH Resistance NOT DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED.

Four replicates of each NTM/ BCG competitive mixture test condition along with a positive control condition with only BCG at  $\sim 3 \times$  LoD were tested. None of the NTM strains tested interfered with the detection of 411 CFU/mL of BCG and generated the correct result as mentioned above. However, under the conditions of this study, competitive inhibitory effects were observed in the presence of only one of the two strains of *M.marinum* (ATCC 0927) tested. Interference with *gyrA2* probes was observed only at challenge concentrations  $>10^4$  CFU/mL resulting in FLQ resistance INDETERMINATE calls at these high challenge concentrations. Refer to Section 17, Limitations for further information..

**Table 16. Competitive Interference by NTM on MTB detection and drug susceptibility detection**

Test Condition / NTM Strain ID	NTM CFU/mL	MTB Detected	INH	FLQ	AMK	KAN	CAP	ETH
MTB + <i>M. avium</i> / (NJH)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.gastir</i> / (ATCC 15754)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.gordonae</i> / (NJH)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.gordonae</i> / (ATCC 14470)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.gordonae</i> / (ATCC 35760)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.marinum</i> / (NJH)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.marinum</i> / (ATCC 0927)	10E+06	PASS	PASS	FAIL	PASS	PASS	PASS	PASS
	10E+05	PASS	PASS	FAIL	PASS	PASS	PASS	PASS
	10E+04	PASS	PASS	PASS	PASS	PASS	PASS	PASS
	10E+03	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.xenopi</i> / (ATCC 700084)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.avian</i> / (ATCC 15769)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.intracellulare</i> / (ATCC 35771)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS
MTB + <i>M.abscessus</i> / (ATCC 19977)	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS

**Table 16. Competitive Interference by NTM on MTB detection and drug susceptibility detection (Continued)**

Test Condition / NTM Strain ID	NTM CFU/mL	MTB Detected	INH	FLQ	AMK	KAN	CAP	ETH
<i>MTB + M.kansasii / (ATCC 12478)</i>	10E+06	PASS	PASS	PASS	PASS	PASS	PASS	PASS

"PASS" indicates all the replicates tested generated the expected "RESISTANCE NOT DETECTED" result for the relevant drugs;  
"FAIL" indicates at least one or more replicates generated "RESISTANCE INDETERMINATE" result for the particular drug.

### 19.7 Inactivation of Mycobacteria in Sputum Samples

The disinfection capability of the Xpert MTB/XDR Sample Reagent was determined using a standardized quantitative tuberculocidal culture method.<sup>21</sup> Samples of sputum were spiked with a high concentration of viable *M. bovis*, mixed with sample reagent at a ratio of 2:1 and incubated for 15 minutes. Following incubation, the sample reagent/sputum mixture was neutralized by dilution and filtration and then cultured. The viability of the *M. bovis* organisms from the treated sputum was reduced by at least 6 logs relative to the un-treated control.

Each laboratory must determine the effectiveness of the sample reagent disinfection properties using their own standardized methods and must adhere to recommended biosafety regulations.

## 20 Precision and Reproducibility

The precision and reproducibility of the Xpert MTB/XDR test was established in a multicenter (three sites), blinded study utilizing a multi-factor nested design. The study consisted of a five-member sample panel and each panel member was prepared by spiking an MTB wild type (WT) strain and an MTB mutant (MUT) strain into artificial sputum matrix. The WT and MUT strains were made from plasmids carrying either MTB XDR wild type or mutant sequences for the genes targeted by the assay, encapsulated in killed, chemically fixed *E. coli*.

The panel members were prepared at ~1xLoD and ~3xLoD using the melt temperatures ( $T_m$ ) of the *inhA* promoter target in the Xpert MTB/XDR test, which generates the MTB DETECTED/NOT DETECTED result depending on the presence or absence of the wildtype or mutant *inhA* promoter specific  $T_m$ . Testing was conducted for six days with three lots of Xpert MTB/XDR cartridges. Each site had two operators (Op1 and Op2) who performed two runs each with two replicates/run each day. A replicate was a single cartridge test. The percent agreement for each panel member is presented in Table 17.

**Table 17. Percent Agreement of Xpert MTB/XDR Test for MTB and *inhA* Detection**

Sample	Site 1			Site 2			Site 3			Total Agreement by Sample
	OP 1	OP 2	Subtotal	OP 1	OP 2	Subtotal	OP 1	OP 2	Subtotal	
MTB MUT 1xLoD	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	95.8% (23/24)	97.9% (47/48)	91.7% (22/24)	91.7% (22/24)	91.7% (44/48)	96.5% (139/144)
MTB MUT 3xLoD	95.8% (23/24)	100% (24/24)	97.92% (47/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	99.3% (143/144)
MTB WT 1xLoD	100% (24/24)	91.67% (22/24)	95.8% (46/48)	91.7% (22/24)	91.7% (22/24)	91.7% (44/48)	91.7% (22/24)	100% (24/24)	95.8% (46/48)	94.4% (136/144)
MTB WT 3xLoD	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
NEG	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	95.8% (23/24)	97.9% (47/48)	99.3% (143/144)

The performance of the Xpert MTB/XDR test in MTB WT and MUT strains at low (~1x) and moderate (~3x) LoD panel samples for each gene target where MTB was detected is presented in Table 18.

**Table 18. Percent Agreement of Xpert MTB/XDR Test in MTB MUT and WT Types Specimens**

Drug	Percent Concordance			
	MTB MUT 1x LoD (95% CI) [n agree/total n]	MTB MUT 3xLoD (95% CI) [n agree/total n]	MTB WT 1x LoD (95% CI) [n agree/total n]	MTB WT 3x LoD (95% CI) [n agree/total n]
INH	100.00% (97.3-100) [139/139]	100.00% (97.4-100.0) [143/143]	89.1% (82.6-93.4) [115/129]	99.3% (96.2-99.9) [143/144]
FLQ	87.80% (81.3-92.2) [122/139]	100.00% (97.4-100.0) [143/143]	81.4% (73.8-87.2) [105/129]	95.8% (91.2-98.1) [138/144]
ETH	100.00% (97.3-100) [139/139]	100.00% (97.4-100.0) [143/143]	99.2% (95.7-99.9) [128/129]	100.0% (97.4-100.0) [144/144]
AMK	100.00% (97.3-100) [139/139]	100.00% (97.4-100.0) [143/143]	91.5% (85.4-95.2) [118/129]	98.6% (95.1-99.6) [142/144]
CAP	99.30% (96.3-99.0) [138/139]	100.00% (97.4-100.0) [143/143]	98.4% (94.5-99.6) [127/129]	99.3% (96.2-99.9) [143/144]
KAN	100.00% (97.3-100) [139/139]	100.00% (97.4-100.0) [143/143]	91.5% (85.4-95.2) [118/129]	98.6% (95.1-99.6) [142/144]

## 21 References

1. WHO. Global Tuberculosis Report 2019. [www.who.int/tb/publications/global\\_report](http://www.who.int/tb/publications/global_report)
2. WHO. 2018, Rapid Communication: Key changes to treatment of multidrug- and rifampicin-resistant tuberculosis.
3. WHO, 2019, Rapid Communication: Key changes to the treatment of drug-resistant tuberculosis
4. Sulis G, Pai M (2020) Isoniazid-resistant tuberculosis: A problem we can no longer ignore. *PLoS Med* 17(1): e1003023
5. Soumitesh Chakravorty, Sandy S. Roh, et al Department of Medicine, New Jersey Medical School, Newark, New Jersey, USA Detection of Isoniazid-, Fluoroquinolone-, Amikacin-, and Kanamycin-Resistant Tuberculosis in an Automated, Multiplexed 10-Color Assay Suitable for Point-of-Care Use *J Clin Microbiol* Jan 2017 Volume 55 Issue 1.
6. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Chosewood, LC and Wilson, DE (eds) (2009). HHS Publication number (CDC) 21-1112.
7. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (refer to latest edition).
8. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Laboratory Detection and Identification of Mycobacteria; Document M48A(refer to latest edition).
9. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC (amending Regulation (EC) No 1907/2007).
10. Occupational Safety and Health Standards, Hazard Communication, Toxic and Hazard Substances (March 26, 2012) (29 C.F.R., pt. 1910, subpt. Z.
11. Kent PT, Kubica GP 1985. Public Health Mycobacteriology—A Guide for Level III Laboratory, Centers of Disease Control, Atlanta, Publication no. PB 86-216546
12. Camus JC, et al. 2002. Re-annotation of the genome sequence of *Mycobacterium tuberculosis* H37Rv. *Microbiology* 148:2967–73.)
13. Zhang et al. 1996. Molecular basis for the exquisite sensitivity of *Mycobacterium tuberculosis* to isoniazid. *Proc Natl Acad Sci U S A*. 1996;93(23):13212–13216
14. Maruri et al. 2012. A systematic review of gyrase mutations associated with fluoroquinolone-resistant *Mycobacterium tuberculosis* and a proposed gyrase numbering system. *J Antimicrob Chemother*.
15. Vadwai V, Ajbani K, Jose M, Vineeth VP, Nikam C, Deshmukh M, Shetty A, Soman R, Rodrigues C. Can inhA mutation predict ethionamide resistance? *Int J Tuberc Lung Dis*. 2013 Jan;17(1):129-30. doi: 10.5588/ijtld.12.0511. Epub 2012 Nov 8. PubMed PMID: 23146620
16. Dunn JJ, Starke JR, Revell PA. 2016. Laboratory diagnosis of *Mycobacterium tuberculosis* infection and disease in children. *J Clin Microbiol* 54:1434–1441. doi:10.1128/JCM.03043
17. Diagnostic Standards and Classification of Tuberculosis in Adults and Children, *Am J Respir Crit Care Med* Vol 161. pp 1376–1395, 2000.
18. Seifert M, Catanzaro D, Catanzaro A, Rodwell TC (2015) Genetic Mutations Associated with Isoniazid Resistance in *Mycobacterium tuberculosis*: A Systematic Review. *PLoS ONE* 10(3): e0119628. doi:10.1371/journal.pone.0119628
19. Farhat MR, Jacobson KR, Franke MF, Kaur D, Sloutsky A, Mitnick CD, Murray M. 2016. Gyrase mutations are associated with variable levels of fluoroquinolone resistance in *Mycobacterium tuberculosis*. *J Clin Microbiol* 54:727–733. doi:10.1128/JCM.02775-15.
20. Gagneux S, Small PM. Global phylogeography of *Mycobacterium tuberculosis* and implications for tuberculosis product development. *Lancet Infect Dis*. 2007 May;7(5):328-37.
21. Banada, P. et.al. Containment of Bioaerosol Infection Risk by the Xpert MTB/RIF Assay and Its Applicability to Point of Care Settings. *Journal of Clinical Microbiology*. 2010.48:10. 3551-3557.

## 22 Cepheid Headquarters Locations

### Corporate Headquarters

Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089  
USA  
Telephone: +1 408.541.4191  
Fax: +1 408.541.4192  
www.cepheid.com

### European Headquarters

Cepheid Europe SAS  
Vira Solelh  
81470 Maurens-Scopont  
France  
Telephone: +33 563 825 300  
Fax: +33 563 825 301  
www.cepheidinternational.com

## 23 Technical Assistance

Before contacting Cepheid Technical Support, collect the following information:

- Product name
- Lot number
- Serial number of the instrument
- Error messages (if any)
- Software version and, if applicable, Computer Service Tag number

### Contact Information

United States  
Telephone: + 1 888 838 3222  
Email: techsupport@cepheid.com

France  
Telephone: + 33 563 825 319  
Email: support@cepheideurope.com

Contact information for all Cepheid Technical Support offices is available on our website:  
[www.cepheid.com/en/CustomerSupport](http://www.cepheid.com/en/CustomerSupport).

## 24 Table of Symbols

Symbol	Meaning
	Catalog number
	<i>In vitro</i> diagnostic medical device
	CE marking – European Conformity
	Do not reuse
	Batch code
	Consult instructions for use
	Manufacturer
	Contains sufficient for <n> tests
	Control
	Expiration date
	Temperature limitation
	Biological risks
	Caution
	Flammable Liquids
	Skin Corrosion
	Severe Health Hazards
	Country of manufacture

Cepheid AB  
 Röntgenvägen 5  
 SE-171 54 Solna,  
 Sweden

GeneXpert.  
Powered By CEPHEID INNOVATION

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Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR IUO Version 3						
<b>Test Result</b> <b>MTB DETECTED;</b> INH Resistance NOT DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED						
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.3	292.5				
katG-melt	73.8	107.0				
fabG1-melt	71.5	242.0				
ahpC-melt	68.7	41.3				
gyrA1-melt	76.2	73.9				
gyrA2-melt	70.4	75.8				
gyrA3-melt	71.0	129.8				
gyrB2-melt	69.5	77.8				
rrs-melt	75.0	188.7				
eis-melt	68.5	145.3				
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

7 ilustracija. MTB APTIKTA; atsparumas INH, FLQ, AMK, KAN, CAP ir ETH NEAPTIKTAS

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name		MTB-XDR IUO		Version		3
Test Result		MTB DETECTED; INH Resistance DETECTED; FLQ Resistance DETECTED; AMK Resistance DETECTED; KAN Resistance DETECTED; CAP Resistance DETECTED; ETH Resistance DETECTED				
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
	Analyte Name		Melt Peak Temperature			Melt Peak Height
	inhA-melt					
	katG-melt					
	fabG1-melt					
	ahpC-melt					
	gyrA1-melt		76.1			90.0
	gyrA2-melt		69.6			39.7
	gyrA3-melt					
	gyrB2-melt					
	rrs-melt					
	eis-melt					
	inhA-mut melt		70.9			259.6
	katG-mut melt		68.4			214.0
	fabG1-mut melt		75.9			181.1
	ahpC-mut melt		66.2			68.2
	gyrA1-mutA melt					
	gyrA1-mutB melt					
	gyrA1-mutC melt					
	gyrA2-mutA melt					
	gyrA2-mutB melt					
	gyrA3-mutA melt					
	gyrA3-mutB melt		76.0			125.0
	gyrA3-mutC melt					
	gyrB2-mut melt		66.0			103.2
	rrs-mut melt		71.0			125.7
	eis-mutA melt		71.4			163.9
	eis-mutB melt					

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Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR IUO		Version	3		
Test Result	<b>MTB DETECTED;</b> <b>INH Resistance DETECTED;</b> FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED					
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
	Analyte Name		Melt Peak Temperature			Melt Peak Height
	inhA-melt		76.6			284.9
	katG-melt		74.0			105.2
	fabG1-melt					
	ahpC-melt		69.0			35.4
	gyrA1-melt		76.6			65.2
	gyrA2-melt		70.4			64.9
	gyrA3-melt		71.4			92.2
	gyrB2-melt		69.7			84.7
	rrs-melt		75.3			146.8
	eis-melt		68.7			124.2
	inhA-mut melt					
	katG-mut melt					
	fabG1-mut melt		75.9			178.0
	ahpC-mut melt					
	gyrA1-mutA melt					
	gyrA1-mutB melt					
	gyrA1-mutC melt					
	gyrA2-mutA melt					
	gyrA2-mutB melt					
	gyrA3-mutA melt					
	gyrA3-mutB melt					
	gyrA3-mutC melt					
	gyrB2-mut melt					
	rrs-mut melt					
	eis-mutA melt					
	eis-mutB melt					

9 ilustracija. MTB APTIKTA; atsparumas INH APTIKTAS

Test Result Analyte Result Detail Melt Peaks Errors History Support

Assay Name MTB-XDR RU0 Version 4

Test Result

MTB DETECTED;  
 INH Resistance DETECTED;  
 H.Q. Resistance NOT DETECTED;  
 AMK Resistance INDETERMINATE;  
 KAN Resistance DETECTED;  
 CAP Resistance INDETERMINATE;  
 ETH Resistance NOT DETECTED

For Research Use Only

Test Result Analyte Result Detail Melt Peaks Errors History Support

Analyte Name	Melt Peak Temperature	Melt Peak Height
inhA-melt	76.3	265.9
katG-melt		
fabG1-melt	71.5	183.2
ahpC-melt	68.9	50.6
gyrA1-melt	76.3	83.8
gyrA2-melt	70.3	63.8
gyrA3-melt	71.2	84.2
gyrB2-melt	69.5	100.1
rrs-melt		
eis-melt		
inhA-mut melt		
katG-mut melt	68.3	168.7
fabG1-mut melt		
ahpC-mut melt		
gyrA1-mutA melt		
gyrA1-mutB melt		
gyrA1-mutC melt		
gyrA2-mutA melt		
gyrA2-mutB melt		
gyrA3-mutA melt		
gyrA3-mutB melt		
gyrA3-mutC melt		
gyrB2-mut melt		
rrs-mut melt		
eis-mutA melt		
eis-mutB melt	63.7	57.0

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Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR IUO Version 3						
<b>Test Result</b> MTB DETECTED; INH Resistance DETECTED; Low FLQ Resistance DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED						
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name			Melt Peak Temperature			Melt Peak Height
inhA-melt			76.5			313.1
katG-melt						
fabG1-melt			71.7			211.5
ahpC-melt			69.0			47.2
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt			69.6			81.1
rrs-melt			75.2			248.1
eis-melt			68.8			158.2
inhA-mut melt						
katG-mut melt			68.4			184.6
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt			72.3			125.0
gyrA1-mutC melt						
gyrA2-mutA melt			76.0			207.9
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt			76.5			128.0
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

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Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name		MTB-XDR IUO		Version 3		
Test Result		<p>MTB DETECTED;                      INH Resistance DETECTED;                      FLQ Resistance DETECTED;                      AMK Resistance DETECTED;                      KAN Resistance DETECTED;                      CAP Resistance NOT DETECTED;                      ETH Resistance NOT DETECTED</p>				
For Investigational Use Only.						
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.6	278.9				
katG-melt						
fabG1-melt	71.7	226.6				
ahpC-melt	69.0	42.9				
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt	69.8	68.7				
rrs-melt	75.3	198.7				
eis-melt						
inhA-mut melt						
katG-mut melt	68.5	204.1				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt	72.9	88.0				
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt	69.1	113.4				
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt	71.6	183.4				
eis-mutB melt						

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14 iliustracija. NEGALIOJANTIS

The screenshot displays the Xpert MTB/XDR software interface. At the top, there are tabs for 'Test Result', 'Analyte Result', 'Detail', 'Melt Peaks', 'Errors', 'History', and 'Support'. Below the tabs, the 'Assay Name' is 'MTB-XDR IUO' and the 'Version' is '3'. The 'Test Result' field shows 'ERROR' in a yellow box. Below this, there is a section labeled 'For Investigational Use Only'. Below the error message, there is a table with columns 'Analyte Name', 'Melt Peak Temperature', and 'Melt Peak Height'. The table lists various analytes, including 'inhA-melt', 'katG-melt', 'fabG1-melt', 'ahpC-melt', 'gyrA1-melt', 'gyrA2-melt', 'gyrA3-melt', 'gyrB2-melt', 'rrs-melt', 'eis-melt', and their mutant counterparts.

Analyte Name	Melt Peak Temperature	Melt Peak Height
inhA-melt		
katG-melt		
fabG1-melt		
ahpC-melt		
gyrA1-melt		
gyrA2-melt		
gyrA3-melt		
gyrB2-melt		
rrs-melt		
eis-melt		
inhA-mut melt		
katG-mut melt		
fabG1-mut melt		
ahpC-mut melt		
gyrA1-mutA melt		
gyrA1-mutB melt		
gyrA1-mutC melt		
gyrA2-mutA melt		
gyrA2-mutB melt		
gyrA3-mutA melt		
gyrA3-mutB melt		
gyrA3-mutC melt		
gyrB2-mut melt		
rrs-mut melt		
eis-mutA melt		
eis-mutB melt		

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Padermės ID	Genas	Tikėtina mutacija	MTB aptikta	Zonde aptikta mutacijos Tm (# teigiami/tirta)	Teisingas ATSPARUMAS APTIKTA (# teigiami/tirta)
Klinikinis	<i>gyrA</i>	GAC 94 TAC	3/3	<i>gyrA1</i> -MutB (3/3); <i>gyrA3</i> -MutC (3/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>fabG1</i>	G609A		<i>fabG1</i> Mut (3/3)	INH [3/3]
Klinikinis	<i>gyrA</i>	GGC 88 GCC, GCG 90 GTG, TCG 91 CCG	3/3	<i>gyrA1</i> -MutB (2/3), <sup>a</sup> <i>gyrA1</i> -MutC (2/3), <i>gyrA2</i> -MutA (3/3), <i>gyrA3</i> -MutB (1/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>rrs</i>	A1410G		<i>rrs</i> -Mut (3/3)	AMK, CAP, KAN [3/3]
Klinikinis	<i>gyrA</i>	GAC 94 GGC	3 / 3	<i>gyrA3</i> -MutB (3/3)	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> Mut (3/3)	INH [3/3]
	<i>rrs</i>	A1410G		<i>rrs</i> -Mut (3/3)	AMK, CAP, KAN [3/3]
14-14194	<i>gyrA</i>	GAC 94 GCC	3 / 3	<i>gyrA1</i> -MutA, <i>gyrA2</i> -MutA	FLQ [3/3]
	<i>katG</i>	AGC 315 ACC		<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>inhA</i>	C -15 T		<i>inhA</i> -Mut (3/3)	INH, ETH [3/3]
15-14175	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>eis</i>	-10G/A		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]
15-14191	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>eis</i>	-10G/A		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]
16-05612	<i>katG</i>	AGC 315 ACC	3 / 3	<i>katG</i> -Mut (3/3)	INH [3/3]
	<i>inhA</i>	C -15 T		<i>inhA</i> -Mut (3/3)	INH, ETH [3/3]
	<i>eis</i>	-12C/T		<i>eis</i> -Mut (3/3)	AMK, KAN [3/3]

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Padermės ID	Genas	Tikėtina mutacija	MTB aptikta	Zonde aptikta mutacijos Tm (# teigiami/tirta)	Teisingas ATSPARUMAS APTIKTA (# teigiami/tirta)
16-05613	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
	eis	-12C/T		eis-Mut (3/3)	AMK, KAN [3/3]
14-13764	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	ahpC	-48G/A		ahpc-Mut (3/3)	INH [3/3]
14-13806	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	ahpC	-48G/A		ahpc-Mut (3/3)	INH [3/3]
Klinikinis	gyrA	GCG 90 GTG, GAC 94 GGC	3 / 3	gyrA3-MutB (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH [3/3]
	ahpC	G-6A		ahpC (2/3) <sup>b</sup>	INH [3/3]
Klinikinis	katG	AGC 315 ACC	3 / 3	katG Mut (3/3)	INH [3/3]
Klinikinis	gyrB2	ACC 539 AAC	3 / 3	gyrB2 WT <sup>c</sup>	*Atsparumas neaptiktas [0/3]
	rrs	A1410G		rrs-Mut (3/3)	AMK, CAP, KAN [3/3]
	gyrA	GCG 90 GTG		gyrA1 MuB (3/3), gyrA2 MutA (3/3), gyrA3 MutB (3/3)	FLQ [3/3]
	ahpC	g -6 a		ahpC Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA Mut (3/3)	INH, ETH [3/3]
Klinikinis	gyrA	TCG 91 CCG	3 / 3	gyrA1-MutB (3/3), gyrA2-MutA (3/3), gyrA3-MutC (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA Mut (3/3)	INH, ETH [3/3]

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1. WHO. Global Tuberculosis Report 2019. [www.who.int/tb/publications/global\\_report](http://www.who.int/tb/publications/global_report)
2. WHO. 2018, Rapid Communication: Key changes to treatment of multidrug- and rifampicin-resistant tuberculosis.
3. WHO, 2019, Rapid Communication: Key changes to the treatment of drug-resistant tuberculosis
4. Sulis G, Pai M (2020) Isoniazid-resistant tuberculosis: A problem we can no longer ignore. *PLoS Med* 17(1): e1003023
5. Soumitesh Chakravorty, Sandy S. Roh, et al Department of Medicine, New Jersey Medical School, Newark, New Jersey, USA Detection of Isoniazid-, Fluoroquinolone-, Amikacin-, and Kanamycin-Resistant Tuberculosis in an Automated, Multiplexed 10-Color Assay Suitable for Point-of-Care Use *J Clin Microbiol* Jan 2017 Volume 55 Issue 1.
6. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Chosewood, LC and Wilson, DE (eds) (2009). HHS Publication number (CDC) 21-1112.
7. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (refer to latest edition).
8. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Laboratory Detection and Identification of Mycobacteria; Document M48A(refer to latest edition).
9. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC (amending Regulation (EC) No 1907/2007).
10. Occupational Safety and Health Standards, Hazard Communication, Toxic and Hazard Substances (March 26, 2012) (29 C.F.R., pt. 1910, subpt. Z).
11. Kent PT, Kubica GP 1985. Public Health Mycobacteriology—A Guide for Level III Laboratory, Centers of Disease Control, Atlanta, Publication no. PB 86-216546
12. Camus JC, et al. 2002. Re-annotation of the genome sequence of *Mycobacterium tuberculosis* H37Rv. *Microbiology* 148:2967–73.)
13. Zhang et al. 1996. Molecular basis for the exquisite sensitivity of *Mycobacterium tuberculosis* to isoniazid. *Proc Natl Acad Sci U S A*. 1996;93(23):13212–13216
14. Maruri et al. 2012. A systematic review of gyrase mutations associated with fluoroquinolone-resistant *Mycobacterium tuberculosis* and a proposed gyrase numbering system. *J Antimicrob Chemother*.
15. Vadwai V, Ajbani K, Jose M, Vineeth VP, Nikam C, Deshmukh M, Shetty A, Soman R, Rodrigues C. Can inhA mutation predict ethionamide resistance? *Int J Tuberc Lung Dis*. 2013 Jan;17(1):129-30. doi: 10.5588/ijtld.12.0511. Epub 2012 Nov 8. PubMed PMID: 23146620
16. Dunn JJ, Starke JR, Revell PA. 2016. Laboratory diagnosis of *Mycobacterium tuberculosis* infection and disease in children. *J Clin Microbiol* 54:1434–1441. doi:10.1128/JCM.03043
17. Diagnostic Standards and Classification of Tuberculosis in Adults and Children, *Am J Respir Crit Care Med* Vol 161. pp 1376–1395, 2000.
18. Seifert M, Catanzaro D, Catanzaro A, Rodwell TC (2015) Genetic Mutations Associated with Isoniazid Resistance in *Mycobacterium tuberculosis*: A Systematic Review. *PLoS ONE* 10(3): e0119628. doi:10.1371/journal.pone.0119628
19. Farhat MR, Jacobson KR, Franke MF, Kaur D, Sloutsky A, Mitnick CD, Murray M. 2016. Gyrase mutations are associated with variable levels of fluoroquinolone resistance in *Mycobacterium tuberculosis*. *J Clin Microbiol* 54:727–733. doi:10.1128/JCM.02775-15.
20. Gagneux S, Small PM. Global phylogeography of *Mycobacterium tuberculosis* and implications for tuberculosis product development. *Lancet Infect Dis*. 2007 May;7(5):328-37.
21. Banada, P. et.al. Containment of Bioaerosol Infection Risk by the Xpert MTB/RIF Assay and Its Applicability to Point of Care Settings. *Journal of Clinical Microbiology*. 2010.48:10. 3551-3557.

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