CERTEST C. difficile antigen GDH POSITIVE CONTROL
CERTEST C. difficile Toxin A POSITIVE CONTROL
CERTEST C. difficile Toxin B POSITIVE CONTROL

INTENDED USE

The CerTest C. difficile antigen GDH Positive Control is an external quality control for CerTest Clostridium difficile GDH+Toxin A+B combo card test. **Only for the strip A (GDH).**
 The CerTest C. difficile Toxin A Positive Control is an external quality control for CerTest Clostridium difficile GDH+Toxin A+B combo card test. **Only for the strip B (Toxin A).**
 The CerTest C. difficile Toxin B Positive Control is an external quality control for CerTest Clostridium difficile GDH+Toxin A+B combo card test. **Only for the strip C (Toxin B).**
 Use those positive control swabs to check that the extraction reagents and the test are working properly.

COMPOSITION

CerTest C. difficile antigen GDH Positive Control: Clostridium difficile Glutamate Dehydrogenase (GDH) antigen dried in swab containing a preservative.
 CerTest C. difficile Toxin A Positive Control: inactivated Toxin A dried in swab containing a preservative.
 CerTest C. difficile Toxin B Positive Control: inactivated Toxin B dried in swab containing a preservative.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The reagent is stable through the expiration date printed on the sealed pouch. The reagent must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- The swab must be carried out within 2 hours of opening the sealed bag.
- Do not use the swab if the aluminium pouch has been damaged during storage.
- The swab should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

MATERIALS PROVIDED (REF.: GX872001VC)

- CerTest C. difficile antigen GDH Positive Control swab
- CerTest C. difficile Toxin A Positive Control swab
- CerTest C. difficile Toxin B Positive Control swab
- CerTest Clostridium difficile GDH+Toxin A+B combo card tests
- Stool collection tubes with diluent (REF.: GX120V, blue cap)
- Positive control reagent tubes (REF.: GX12, red cap)
- Testing tubes
- Disposable pipettes
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves
- Timer


MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves
- Timer

QUALITY CONTROL TESTING PROCEDURE

Allow tests and positive control swabs to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


1. Remove the CerTest Clostridium difficile GDH+Toxin A+B combo card test from its sealed bag just before using (1).
2. Take the positive control reagent tube, cut the end of the cap (2). Use a separate positive control reagent tube, pipette and testing tube for each control.

(1) Remove the combo card test

CerTest C. difficile antigen GDH Positive Control swab procedure


3. Add 15 drops (3) of positive control reagent into a testing tube and immediately introduce the C. difficile antigen GDH Positive Control swab into the testing tube.
4. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (4). Discard the swab.
5. Dispense 3 drops from the testing tube in the circular window marked with the letter A (5).
6. Read the result at 10 minutes. Do not read the test result later than 10 minutes.

- (3) Add 15 drops of the positive control reagent
- (4) Introduce the CerTest C. difficile antigen GDH Positive Control swab into the testing tube, rotating 1 minute and extract the liquid
- (5) Add 3 drops in the circular window GDH (strip A)

CerTest C. difficile Toxin A Positive Control swab procedure


7. Take a different positive control reagent tube, and add 15 drops (6) of positive control reagent into a different testing tube and immediately introduce the C. difficile Toxin A Positive Control swab into another testing tube.
8. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (7). Discard the swab.
9. Dispense 3 drops from the testing tube in the circular window marked with the letter B (8).
10. Read the result at 10 minutes. Do not read the test result later than 10 minutes.

- (6) Add 15 drops of the positive control reagent
- (7) Introduce the CerTest Toxin A Positive Control swab into the testing tube, rotating 1 minute and extract the liquid
- (8) Add 3 drops in the circular window Toxin A (strip B)

CerTest C. difficile Toxin B Positive Control swab procedure

11. Take a different positive control reagent tube, and add 15 drops (9) of positive control reagent into a different testing tube and immediately introduce the C. difficile Toxin B Positive Control swab into another testing tube.
12. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (10). Discard the swab.
13. Dispense 3 drops from the testing tube in the circular window marked with the letter C (11).
14. Read the result at 10 minutes. Do not read the test result later than 10 minutes.

INTERPRETATION OF RESULTS (please refer to the illustration below)

VALID RESULT

Strip A: Clostridium difficile antigen GDH POSITIVE: A GREEN line (control line (C)) and a RED line (GDH test line (T)) in the results window appear at the same time. That means the test is working properly.

Strip B: C. difficile Toxin A POSITIVE: A GREEN line (control line (C)) and a RED line (Toxin A test line (T)) in the results window appear at the same time. That means the test is working properly.

Strip C: C. difficile Toxin B POSITIVE: A GREEN line (control line (C)) and a RED line (Toxin B test line (T)) in the results window appear at the same time. That means the test is working properly.


CerTest Clostridium difficile antigen GDH Positive Control is a control for strip A in CerTest Clostridium difficile GDH+Toxin A+B combo card test.

CerTest C. difficile Toxin A Positive Control is a control for strip B in CerTest Clostridium difficile GDH+Toxin A+B combo card test.

CerTest C. difficile Toxin B Positive Control is a control for strip C in CerTest Clostridium difficile GDH+Toxin A+B combo card test.

INVALID RESULTS


Any other result is considered to be invalid. Total absence of the control coloured lines (GREEN) regardless the appearance or not of the test lines (RED), or only the control (C) colored lines (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.

A: Green-Red→Positive GDH
 B: Green-Red→Positive Toxin A
 C: Green-Red→Positive Toxin B

	A GDH	B Toxin A	C Toxin B	Interpretation of the results
1	+	+	+	GDH (A) POSITIVE , Toxin A (B) POSITIVE and Toxin B (C) POSITIVE : A GREEN line (control line (C)) and a RED line (test line (T)) in the three strips, in the results window, appear at the same time due to presence of Clostridium difficile in the three strips (A, B, C).
2,3 4,5, 6,7, 8,9	GREEN- RED	GREEN- RED	GREEN- RED	Invalid result, we recommend repeating the assay with another test and positive control swabs.

SIMBOLS FOR IVD PRODUCTS AND REAGENTS

CERTEST C. difficile antigen GDH POSITIVE CONTROL
CERTEST C. difficile Toxin A POSITIVE CONTROL
CERTEST C. difficile Toxin B POSITIVE CONTROL

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