CERTEST *Rotavirus* POSITIVE CONTROL
CERTEST *Adenovirus* POSITIVE CONTROL
CERTEST *Norovirus* GII POSITIVE CONTROL

INTENDED USE

The CerTest *Rotavirus* Positive Control is an external quality control for CerTest ROTA+ADENO+NORO combo card test. **Only for the strip A (Rota).**
 The CerTest *Adenovirus* Positive Control is an external quality control for CerTest ROTA+ADENO+NORO combo card test. **Only for the strip B (Adeno).**
 The CerTest *Norovirus* GII Positive Control is an external quality control for CerTest ROTA+ADENO+NORO combo card test. **Only for the strip C (Noro).**
 Use those positive control swabs to check that the extraction reagents and the test are working properly.

COMPOSITION

CerTest *Rotavirus* Positive Control: inactivated *Rotavirus* antigens dried in swab containing a preservative.
 CerTest *Adenovirus* Positive Control: inactivated *Adenovirus* antigens dried in swab containing a preservative.
 CerTest *Norovirus* GII Positive Control: inactivated *Norovirus* GII antigens dried in swab containing a preservative.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The reagent is stable through the expiration date printed on the sealed pouch. The reagent must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- The swab must be carried out within 2 hours of opening the sealed bag.
- Do not use the swab if the aluminium pouch has been damaged during storage.
- The swab should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

MATERIALS PROVIDED (REF.: RA872001VC/RA871001VC)

MATERIALS PROVIDED

- CerTest *Rotavirus* Positive Control swab
- CerTest *Adenovirus* Positive Control swab
- CerTest *Norovirus* GII Positive Control swab
- CerTest ROTA+ADENO+NORO combo card tests
- Stool collection tubes with diluent (REF.: RA110V, blue cap)
- Positive control reagent tubes (REF.: RA12, red cap)
- Testing tubes
- Disposable pipettes
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves
- Timer


MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves
- Timer


QUALITY CONTROL TESTING PROCEDURE

Allow tests and positive control swabs to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Remove the CerTest ROTA+ADENO+NORO combo card test from its sealed bag just before using (1).
2. Take the positive control reagent tube, cut the end of the cap (2). Use a separate positive control reagent tube, pipette and testing tube for each control.

(1) Remove the combo card test



(2) Cut the cap

CerTest *Rotavirus* Positive Control swab procedure

3. Add 15 drops (3) of positive control reagent into a testing tube and immediately introduce the ***Rotavirus* Positive Control swab** into the testing tube.
4. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (4). Discard the swab.
5. Dispense 4 drops from the testing tube in the **circular window marked with the letter A** (5).
6. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

Wait 10 minutes

(3) Add 15 drops of the positive control reagent

(4) Introduce the **CerTest *Rotavirus* Positive Control swab** into the testing tube, rotating 1 minute and extract the liquid



(5) Add 4 drops in the circular window **Rota (strip A)**

CerTest *Adenovirus* Positive Control swab procedure

7. Take a different positive control reagent tube, and add 15 drops (6) of positive control reagent into a different testing tube and immediately introduce the ***Adenovirus* Positive Control swab** into another testing tube.
8. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (7). Discard the swab.
9. Dispense 4 drops from the testing tube in the **circular window marked with the letter B** (8).
10. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

Wait 10 minutes

(6) Add 15 drops of the positive control reagent

(7) Introduce the **CerTest *Adenovirus* Positive Control swab** into the testing tube, rotating 1 minute and extract the liquid


(8) Add 4 drops in the circular window **Adeno (strip B)**

CerTest *Norovirus* GII Positive Control swab procedure

11. Take a different positive control reagent tube, and add 15 drops (9) of positive control reagent into a different testing tube and immediately introduce the ***Norovirus* GII Positive Control swab** into another testing tube.
12. Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn (10). Discard the swab.
13. Dispense 4 drops from the testing tube in the **circular window marked with the letter C** (11).
14. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.


(9) Add 15 drops of the positive control reagent

(10) Introduce the **CerTest *Norovirus* GII Positive Control swab** into the testing tube, rotating 1 minute and extract the liquid

Wait 10 minutes

(11) Add 4 drops in the circular window **Noro (strip C)**

INTERPRETATION OF RESULTS (please refer to the illustration below)

VALID RESULT

Strip A: *Rotavirus* POSITIVE: A GREEN line (control line (C)) and a RED line (Rota test line (T)) in the results window appear at the same time. That means the test is working properly.

Strip B: *Adenovirus* POSITIVE: A GREEN line (control line (C)) and a RED line (Adeno test line (T)) in the results window appear at the same time. That means the test is working properly.

Strip C: *Norovirus* GII POSITIVE: A GREEN line (control line (C)) and a RED line (Noro test line (T)) in the results window appear at the same time. That means the test is working properly.

CerTest *Rotavirus* Positive Control is a control for strip A in CerTest ROTA+ADENO+NORO combo card test.

CerTest *Adenovirus* Positive Control is a control for strip B in CerTest ROTA+ADENO+NORO combo card test.

CerTest *Norovirus* GII Positive Control is a control for strip C in CerTest ROTA+ADENO+NORO combo card test.

INVALID RESULTS


Any other result is considered to be invalid. Total absence of the control coloured lines (GREEN) regardless the appearance or not of the test lines (RED), or only the control (C) colored lines (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.

1. POSITIVE

A: Green-Red→Positive Rota
B: Green-Red→Positive Adeno
C: Green-Red→Positive Noro


2. INVALID


3. INVALID


4. INVALID


5. INVALID


6. INVALID

7. INVALID

8. INVALID

	A Rota	B Adeno	C Noro	Interpretation of the results
1	+	+	+	<i>Rotavirus</i> (A) POSITIVE , <i>Adenovirus</i> (B) POSITIVE and <i>Norovirus</i> GII (C) POSITIVE : A GREEN line (control line (C)) and a RED line (test line (T)) in the three strips, in the results window, appear at the same time due to presence of antigens in the <i>Rotavirus</i> Positive Control swab in strip A, presence of antigens in <i>Adenovirus</i> Positive Control swab in strip B and presence of antigens in <i>Norovirus</i> GII Positive Control swab in strip C.
2,3 4,5, 6,7, 8,9	Any other result			Invalid result, we recommend repeating the assay with another test and positive control swabs.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

In vitro diagnostic device

Consult instructions for use

Keep dry

Temperature limitation

Use by

Positive Control

Batch code

Catalogue number

Contains sufficient for <n> tests

Manufacturer

DIL

Diluent

CerTest
BIOTEC S.L.

Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840, San Mateo de Gállego, Zaragoza
(SPAIN) www.certest.es

October 2014. Revision: 00

IU-RA287 rev. 00

CERTEST Rotavirus POSITIVE CONTROL
CERTEST Adenovirus POSITIVE CONTROL
CERTEST Norovirus GII POSITIVE CONTROL

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