Rotavirus (ROTA)

Description

Rotavirus is highly infectious even with a small number of fewer than 100 copies, causing severe dehydration symptoms due to frequent watery diarrhea accompanied by fever, vomiting, and abdominal pain, especially in infants and young children, and in severe cases can lead to death. ichroma™ ROTA test diagnoses infected people earlier prevent outbreaks and allow quick and appropriate treatment.

Indication

Acute gastroenteritis, Rotavirus gastroenteritis.

Diagnostics Value

ROTA tests do not require ELISA or PCR. This makes it an ideal choice for neonatal and children's clinics for detecting ROTA infection promptly before they spread to the general population.

Selection guide

Item / Platform	ichroma™ II	ichroma™ III
ichroma™ ROTA	●	●

Specification

Sample type	Feces
Sample vol.	Spot
Assay time	12 min
Working range	Qualitative
Storage	4-30°C
Shelf life	20 months

Performance

Accuracy	Positive: 99.02% Negative: 95.23% Overall: 97.10%
----------	---

Ordering information

Cat #	CFPC-75
Test/Kit	25T/Kit

Rotavirus and Adenovirus (ROTA/Adeno)

Description

Rotavirus causes severe diarrhea in infants and infants under the age of 2 worldwide, and intestinal adenovirus causes gastroenteritis mainly in children under 3 years old, immunosuppressed patients and bone marrow transplant patients. Differential diagnosis with ichroma™ Rota/Adeno reduces antibiotic abuse and misdiagnosis and helps with proper treatment.

Indication

Acute gastroenteritis, Rotavirus gastroenteritis, Adenovirus gastroenteritis.

Diagnostics Value

ichroma™ Rota/Adeno test kits identify the virus type of the rotavirus and the adenovirus sharing identical symptoms. These POC tests can help reduce the likelihood of misuse of antibiotics in neonatal clinics and children's hospitals.

Selection guide

Item / Platform	ichroma™ II	ichroma™ III
ichroma™ Rota/Adeno	●	●

Specification

Sample type	Feces
Sample vol.	Spot
Assay time	12 min
Working range	Qualitative
Storage	4-30°C
Shelf life	20 months

Performance

Accuracy	Sensitivity	Rota 98.2%, Adeno 95%
	Specificity	Rota 95.7%, Adeno 96.2%

Ordering information

Cat #	CFPC-79
Test/Kit	25T/Kit

Mycoplasma

Description

Mycoplasma causes acute pneumonia with high fever and complications through upper respiratory infection, and it must be treated with specific antibiotics with severe side effects. Therefore, it is necessary to distinguish mycoplasma from common bacterial and viral infections. This test detects mycoplasma specifically and helps to shorten the patient's recovery time by prescribing appropriate antibiotics.

Indication

Mycoplasma infection, Mycoplasma pneumonia.

Diagnostics Value

Mycoplasma is highly contagious and can be dangerous for the elderly and children. One needs specific antibiotics for the treatment. Mycoplasma tests will help clinicians with the diagnosis and the treatment planning against Mycoplasma by providing an early diagnosis.

Selection guide

Item / Platform	ichroma™ M2	ichroma™ II	ichroma™ III	ichroma™-50
ichroma™ Mycoplasma	●	●	●	●

Specification

Sample type	Throat swab
Assay time	10 min
Working range	Qualitative
Storage	1-30°C
Shelf life	18 months

Performance

Accuracy	Clinical sensitivity: 81.9% Clinical specificity: 98.8%
----------	---

Ordering information

Cat #	CFPC-94
Test/Kit	25T/Kit

IGRA-TB (IFN- γ release assay-tuberculosis)

Description

Patients with latent tuberculosis or reduced immunity can be converted into active tuberculosis patients with lethal transmission, so proactive treatment and response are required. ichroma™ IGRA-TB is a unique product that combines a highly specific cell-mediated immune response (CMI) against *Mycobacterium tuberculosis* and a highly sensitive lateral flow detection system. It can help fight tuberculosis by detecting latent tuberculosis in the simplest, fastest, and most accurate way.

Indication

Latent tuberculosis, LTB.

Diagnostics Value

IGRA-TB provides easy-to-use, rapid, and cost-effective diagnostics for latent TBs and can replace the time-consuming and costly ELISA based tests. These features would facilitate the utility at the level of public health centers or equivalent and third-party laboratories.

Selection guide

Item / Platform	ichroma™ II	ichroma™ III	ichroma™-50
ichroma™ IGRA-TB	●	●	●

Specification

Sample type	Whole blood
Sample vol.	1 mL × 3 tubes
Assay time	15 min
Working range	Qualitative
Storage	4-30°C
Shelf life	20 months

Performance

Accuracy (Compared to Qunatiferon-Plus)	Positive percent agreement: 81.6% Negative percent agreement: 97.5%
--	--

Ordering information

Cat #	CFPC-86	CFPC-86-1
Test/Kit	100T/Kit (300 cartridges)	50T/Kit (150 cartridges)

Autoimmune

Rheumatoid Arthritis (RF IgM)	84
Total IgE	85
Anti-CCP Plus	86

Rheumatoid Arthritis (RF IgM)

Description

It helps to diagnose rheumatoid arthritis, distinguish it from other types of arthritis and other inflammatory diseases, and ensure proper treatment by rapidly and accurately quantifying RF IgM in blood.

Indication

Rheumatoid arthritis, autoimmune disease.

Diagnostics Value

It helps diagnose rheumatoid factor (RA) in the blood with high sensitivity and distinguishes it from other types of arthritis and other inflammatory diseases, thereby helping active treatment to minimize complications and tissue damage.

Selection guide

Item / Platform	ichroma™	ichroma™ II	ichroma™ III
ichroma™ RF IgM	●	●	●

Specification

Sample type	Whole blood, Serum, Plasma
Sample vol.	10/5/5 µL
Assay time	5 min
Working range	8-200 IU/mL
Storage	2-8°C for DB/4-30°C for Cartridge
Shelf life	20 months

Performance

Accuracy	$y=1.082x - 12.038$ $R=0.9035$
----------	--------------------------------

Ordering information

Cat #	CFPC-39
Test/Kit	25T/Kit

Total IgE

Description

It quantifies IgE in blood to support diagnosing of acute allergic diseases.

Indication

Allergic disease (Type 1 hypersensitivity reaction), autoimmune disease, Job's syndrome.

Diagnostics Value

It can be usefully used in small and medium hospitals including emergency rooms, as it can specifically detect immunoglobulin E in blood to diagnose excessive acute allergies and provide appropriate treatment to patients.

Selection guide

Item / Platform	ichroma™ II	ichroma™ III
ichroma™ Total IgE	●	●

Specification

Sample type	Whole blood, Serum, Plasma
Sample vol.	100/50/50 µL
Assay time	12 min
Working range	1-1,000 IU/mL
Storage	2-8°C for DB/4-30°C for Cartridge
Shelf life	20 months

Performance

Accuracy	$y=0.9954x + 1.0943$ $R=0.9982$
----------	---------------------------------

Ordering information

Cat #	CFPC-91
Test/Kit	25T/Kit

Anti-CCP Plus

Description

It helps to diagnose and provide appropriate treatment for chronic systemic autoimmune diseases and rheumatism by quantifying Anti-CCP in the blood.

Indication

Rheumatoid arthritis, autoimmune disease.

Diagnostics Value

It supports to diagnose rheumatoid arthritis at an earlier stage by quantifying anti-citrulline peptide antibody (Anti-CCP) in the blood to minimize complications and tissue damage by distinguishing it from other types of arthritis and other inflammatory diseases. Furthermore, it is useful for monitoring the effectiveness of treatment and treatment as appropriate.

Selection guide

Item / Platform	ichroma™ II	ichroma™ III
ichroma™ Anti-CCP Plus	●	●

Specification

Sample type	Whole blood, Serum, Plasma
Sample vol.	5 µL
Assay time	12 min
Working range	3.5-300 U/mL
Storage	2-8°C for DB/4-30°C for Cartridge
Shelf life	20 months

Performance

Accuracy	Positive: 93.9% Negative: 96.0% Overall: 94.9%
----------	--

Ordering information

Cat #	CFPC-97
Test/Kit	25T/Kit

Other test items

Calprotectin	88
Ferritin	89
<i>Helicobacter pylori</i> Stool Antigen (H. Pylori SA)	90
Vitamin D	91

Calprotectin

Description

The elevated level of fecal calprotectin may indicate a case of intestinal inflammation including IBD. This test can tell apart the IBD from the non-inflammatory bowel disease: the level stays low with the latter. It measures the level of fecal calprotectin to help differentiate those seemingly similar diseases.

Indication

Inflammatory bowel disease (IBD): Ulcerative Colitis (UC), Crohn's disease (CD)
Irritable bowel syndrome (IBS)

Diagnostics Value

The test value can render the status of the inflammatory bowel disease (IBD) and the response to the therapy. The calprotectin test would provide the essential information for the diagnosis and monitoring of the therapeutic progress.

Selection guide

Item / Platform	ichroma™	ichroma™ II	ichroma™ III	ichroma™-50
ichroma™ Calprotectin	●	●	●	●

Specification

Sample type	Feces
Sample vol.	10 mg
Assay time	10 min
Working range	10-1,000 mg/kg feces
Storage	4-30°C
Shelf life	20 months

Performance

Accuracy	$y=1.0602x - 2.9932$ $R=0.9866$
----------	---------------------------------

Ordering information

Cat #	CFPC-83
Test/Kit	25T/Kit

Ferritin

Description

It helps in the diagnosis of diseases related to iron metabolisms, such as iron deficiency and excess, as well as inflammation and tumors.

Indication

Iron deficiency anemia, anemia of chronic disease, Hemolytic anemia, Hereditary hemochromatosis, acute-phase reactants (APR), tumor marker.

Diagnostics Value

The Ferritin test measures the level of ferritin in blood to provide answers on iron metabolic disorders such as iron deficiency anemia and iron excess, as well as inflammation and tumor development, helping to provide appropriate treatment to patients.

Selection guide

Item / Platform	ichroma™	ichroma™ II	ichroma™ III
ichroma™ Ferritin	●	●	●

Specification

Sample type	Serum, Plasma
Sample vol.	30 µL
Assay time	10 min
Working range	10-1,000 ng/mL
Storage	2-8°C for DB/4-30°C for Cartridge
Shelf life	20 months

Performance

Accuracy	$y=0.99198x + 0.56317$ R=0.9897
----------	---------------------------------

Ordering information

Cat #	CFPC-32
Test/Kit	25T/Kit

H. pylori SA

Description

H. pylori SA is an antigen produced by the immune system to cope with the *H. Pylori* infection. The level of the antigen examined through the stool serves as a more reliable indicator of the infection than by the one via a blood test, which does not tell if the infection is still active.

Indication

Chronic gastritis, gastric ulcers conditions, duodenal ulcers, stomach cancer etc.

Diagnostics Value

It is designed to be POC compatible, enabling a physician to make the diagnostics on-site, without the lengthy delay with the conventional test. With the immediate result on hand, the clinician can deliver diagnostics and treatment plans in the shortest time possible to maximize the patient's convenience.

Selection guide

Item / Platform	ichroma™	ichroma™ II	ichroma™ III
ichroma™ H. pylori SA	●	●	●

Specification

Sample type	Feces
Sample vol.	10 mg
Assay time	12 min
Working range	Qualitative
Storage	4-30°C
Shelf life	20 months

Performance

		Reference ELISA result		
		Positive	Negative	Total
ichroma™ H. pylori SA	Positive	102	1	103
	Negative	3	54	57
	Total	105	55	160

- Percent positive agreement = 97.1%
- Percent negative agreement = 98.2%
- Overall percent agreement = 97.5%

Ordering information

Cat #	CFPC-81
Test/Kit	25T/Kit

Vitamin D

Description

Vitamin D is closely related to bone formation, as well as autoimmune regulation, metabolic function, and cancer prevention. This product is useful to check for association with vitamin D when abnormal calcium, phosphorus, and parathyroid hormone (PTH) levels are observed in a patient. In addition, when starting drug treatment for bone disease, bone weakness, or osteoporosis, it is useful for diagnosis and monitoring of vitamin D prescription treatment through periodic quantification of vitamin D. The quality of this product is monitored by participating in CDC's Vitamin D Standardization Certification Program (VDSCP) every quarter.

Indication

Deficiency of Vitamin D, Excess of Vitamin D.

Diagnostics Value

Vitamin D plays an essential role in making strong bones. The vitamin D test kits come in dry format, eliminating the need for cold chain logistics. With the room temperature storage and its operation on a POC platform, one can integrate the vitamin D test into a routine on-site diagnostic service.

Selection guide

Item / Platform	ichroma™	ichroma™ II	ichroma™ III
ichroma™ Vitamin D	●	●	●

Specification

Sample type	Serum, Plasma
Sample vol.	50 µL
Assay time	28 min
Working range	8-70 ng/mL
Storage	4-30°C
Shelf life	20 months

Performance

Accuracy	$y=1.1629x - 0.2829$ $R=0.9951$
----------	---------------------------------

Ordering information

Cat #	CFPC-47
Test/Kit	25T/Kit

Molecular diagnostics

NuActor	94
ExAmplar	95
NuActor Viral RNA extraction kit	96
NuActor Viral DNA extraction kit	97
NuActor sputum DNA extraction kit	98
ExAmplex COVID-19 PCR 3-gene Lyo kit	99
ExAmplex COVID-19 PCR 3-gene kit	100
ExAmplex Influenza A/B SARS-CoV-2 PCR kit	101
ExAmplar COVID-19 real-time PCR kit (L)	102
ExAmplar Quanti-HBV v1.0 60T	103
ExAmplar Quanti-HCV v1.0 60T	104
ExAmplar Quanti-HIV v1.0 60T	105
ExAmplar MTB kit v2.0 72T	106

Fully Automated Extraction Solution

NuActor

Fully automatic nucleic acid (RNA/DNA) extraction system that can quickly and easily extract high-purity nucleic acids from small amount of samples with a ready-to-use cartridge and Tabletop style extractor.

- ▶ 8 extraction/run in < 13 min: quick result
- ▶ Reagent-ready Cartridge; Ready-to-use cartridge
- ▶ Room Temperature Storage: stable performance
- ▶ 200 μ L sample to 100 μ L elution: concentrated nucleic acid
- ▶ Viral RNA, DNA and sputum DNA extraction kits are available (Exclusive NuActor cartridges)

Specification

Item	NuActor Specification
Dimensions	240 mm (L) x 240 mm (W) x 260 mm (H)
Weight	6.98 Kg
Power supply	100-240V AC(50/60 Hz), 60W
Throughput	8 samples
Operation temperature	15-40°C
Memory	N/A
Interface	Touch-screen LCD
Cat#	FPRR032

Fast and Simple Viral-Load Analyzer

ExAmplar

A cost-effective and compact analyzer based on PCR that detects virus/bacteria quantitatively or qualitatively in a short time using ready-to-use cartridges.

- ▶ 8 Reaction/run in 30 min (DNA) / 45-55 min (RNA): quick result
- ▶ Compact, light and low power consumption (~2.88 kg)
- ▶ Small footprint
- ▶ Stand-alone

Specification

Item	ExAmplar Specification
Dimensions	200 mm (L) x 240 mm (W) x 180 mm (H)
Weight	2.88 Kg
Power supply	100-240V AC(50/60 Hz), 110W
Data output	SD card
Operation temperature	15-40°C
Memory	+5,000 test results
Interface	Touch-screen LCD
Cat#	FPRR030

NuActor Viral RNA extraction kit

Description

This product is provided as ready-to-use cartridges for extracting viral RNA. Samples are placed in a cartridge, and the cartridge is mounted on a NuActor instrument, and it is automatically processed to extract highly purified viral RNA quickly and simply.

Indication

Applicable for Serum / Plasma / VTM / UTM samples.

Diagnostics Value

Stable results in molecular biological diagnostic tests depend on the purity and stability of the extracted RNA used. This product automatically extracts high-purity virus RNA from human blood (serum/plasma) and respiratory swab samples (VTM, UTM) by introducing a cartridge-type all-in-one extraction kit.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
ExAmplex COVID-19 PCR 3-gene Lyo kit		●
ExAmplex Influenza A/B SARS-CoV-2 PCR kit		●
ExAmplar COVID-19 real-time PCR kit (L)		●
Quanti-HCV v1.0 60T		●
Quanti-HIV v1.0 60T		●

Specification

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene Lyo kit	ExAmplex Influenza A/B SARS-CoV-2 PCR kit	ExAmplar COVID-19 real-time PCR kit (L)	Quanti-HCV v1.0 60T	Quanti-HIV v1.0 60T
Sample type	Serum/Plasma VTM/UTM	Swab, sputum	Swab, sputum	Swab, sputum	Serum/Plasma	Serum/Plasma
Sample vol.	200 µL	5 µL	5 µL	5 µL	5 µL	5 µL
Assay time	<13 min	<55 min (RNA)	<55 min (RNA)	<55 min (RNA)	<45 min (RNA)	<45 min (RNA)
Storage	15-35°C	2-30°C	2-8°C	2-8°C	-20°C	-20°C
Shelf life	12 months	12 months	12 months	12 months	12 months	12 months

Ordering information

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene Lyo kit	ExAmplex Influenza A/B SARS-CoV-2 PCR kit	ExAmplar COVID-19 real-time PCR kit (L)	Quanti-HCV v1.0 60T	Quanti-HIV v1.0 60T
Cat #	UFPK-1	UFPK-7	UFPK-6	UFPK-4	UFPK-11	UFPK-12
Test/Kit	96	96	96	48	60 (+ 36 controls)	60 (+ 36 controls)

NuActor Viral DNA extraction kit

Description

This product is provided as ready-to-use cartridges for extracting viral DNA. Samples are placed in a cartridge, and the cartridge is mounted on a NuActor instrument, and it is automatically processed to extract highly purified viral DNA quickly and simply.

Indication

Viral DNA purification from human serum/plasma.

Diagnostics Value

Stable results in molecular biological diagnostic tests depend on the purity and stability of the extracted DNA used. This product introduced a cartridge-type all-in-one extraction kit that automatically extracts high-purity viral DNA from human blood (serum/plasma) and provides stable results.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral DNA Extraction kit	●	
Quanti-HBV v1.0 60T		●

Specification

	Viral DNA Extraction kit	Quanti-HBV v1.0 60T
Sample type	Serum/Plasma	Serum/Plasma
Sample vol.	200 µL	5 µL
Assay time	<13 min	<30 min (DNA)
Storage	15-35°C	-20°C
Shelf life	12 months	12 months

Ordering information

	Viral DNA Extraction kit	Quanti-HBV v1.0 60T
Cat #	UFPK-2	UFPK-2
Test/Kit	96	60 (+36 controls)

NuActor sputum DNA extraction kit

Description

Tuberculosis is the most deadly infectious disease in human history. About 10 million patients infected by tuberculosis occur every year. This kit is a DNA extraction kit required for Tuberculosis PCR testing. It is used with NuActor, an automatic nucleic acid extraction device, to easily extract DNA of Mycobacterium tuberculosis from sputum samples.

Indication

Extraction of bacterial DNA from human Sputum sample.

Diagnostics Value

Sputum samples for tuberculosis testing are difficult to handle due to the high risk of infection. This kit helps to extract bacterial DNA quickly and easily from sputum samples.

Selection guide

Item / Platform	NuActor	ExAmplar
Sputum DNA Extraction Kit	●	
MTB v2.0 72T		●

Specification

	Sputum DNA extraction Kit	MTB v2.0 72T
Sample type	Sputum	Sputum
Sample vol.	250 µL	2.5 µL
Assay time	<13 min	<30 min (DNA)
Storage	15-35°C	2-8°C
Shelf life	12 months	12 months

Ordering information

	Sputum DNA extraction Kit	MTB v2.0 72T
Cat #	UFPK-3	UFPK-5
Test/Kit	96	72 (+12 controls)

ExAmplex COVID-19 PCR 3-gene Lyo kit

(for open system)

Description

This product detects two COVID-19 confirmation genes, and one Pan-corona gene with high accuracy. This product is provided as a single-test unit with ready-to-use form that can be stored at room temperature. It can be used immediately without the process of reaction mixture preparation by the user. After dissolving the dried reaction bead with a diluent, user can simply perform PCR for COVID-19 testing by just adding RNA extracted beforehand.

Indication

SARS-CoV-2 infection (ORF1ab gene, N gene, Sarbeco E gene detection), Coronavirus disease-19 (COVID-19).

Diagnostics Value

The SARS-CoV-2 virus shares common symptoms with other respiratory diseases. It is contagious even in the asymptomatic stage, which makes it imperative to diagnose the virus early to provide timely treatment and to prevent further transmissions. Since this product is designed to increase diagnostics accuracy by detecting 2 confirmation gene of SARS-CoV-2 as well as 1 Pan-corona gene, it helps to prevent spread of disease by early detection of COVID-19. In addition, this kit is stable at room temperature, so it can be easily used even in environments where frozen delivery conditions are not available or freezer storage facilities are not available.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
ExAmplex COVID-19 PCR 3-gene Lyo kit		●

Specification

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene Lyo kit
Sample type	Swab (VTM/UTM)	Swab, Sputum
Sample vol.	200 µL	5 µL
Assay time	<13 min	<55 min (RNA)
Storage	15-35°C	2-30°C (Lyophilized)
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene Lyo kit
Cat #	UFPK-1	UFPK-7
Test/Kit	96	96

ExAmplex COVID-19 PCR 3-gene kit (Liquid)

(for open system)

Description

This product detects two COVID-19 confirmation genes, and one Pan-corona gene with high accuracy. Users make reaction mix by mixing primer mix and enzyme mix solution and can perform PCR for COVID-19 testing by adding RNA extracted beforehand to the reaction mix.

Indication

SARS-CoV-2 infection (ORF1ab gene, N gene, Sarbeco E gene detection), Coronavirus disease-19 (COVID-19).

Diagnostics Value

The SARS-CoV-2 virus shares common symptoms with other respiratory diseases. It is contagious even in the asymptomatic stage, which makes it imperative to diagnose the virus early to provide timely treatment and to prevent further transmissions. Since this product is designed to increase diagnostics accuracy by detecting 2 confirmation genes of SARS-CoV-2 as well as 1 Pan-corona gene, it helps to prevent spread of disease by early detection of COVID-19. This kit is provided as liquid form which needs to be stored in -20°C .

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
ExAmplex COVID-19 PCR 3-gene kit		●

Specification

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene kit
Sample type	Swab (VTM/UTM)	Swab, Sputum
Sample vol.	200 μL	5 μL
Assay time	<13 min	<55 min (RNA)
Storage	15-35 $^{\circ}\text{C}$	-20 $^{\circ}\text{C}$
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	ExAmplex COVID-19 PCR 3-gene kit
Cat #	UFPK-1	UFPK-9
Test/Kit	96	96

ExAmplex Influenza A/B SARS-CoV-2 PCR kit

(for open system)

Description

A highly stable, ready-to-use PCR kit that can diagnose COVID-19 and Influenzas (Flu A, Flu B) simultaneously from a single sample collection.

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Diagnostics Value

This is a product that maximizes user convenience and storage so that RT-PCR products can be used even in areas where frozen delivery is difficult or in an environment that does not have a freezer storage facility. It helps to enable accurate antiviral treatment by simultaneously diagnosing single or multiple infections of COVID-19 and Flu (A & B) from one sample easily.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
ExAmplex Influenza A/B SARS-CoV-2 PCR kit		●

Specification

	Viral RNA extraction kit	ExAmplex Influenza A/B SARS-CoV-2 PCR kit
Sample type	Swab (VTM/UTM)	Swab, Sputum
Sample vol.	200 µL	5 µL
Assay time	<13 min	<55 min (RNA)
Storage	15-35°C	2-8°C
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	ExAmplex Influenza A/B SARS-CoV-2 PCR kit
Cat #	UFPK-1	UFPK-6
Test/Kit	96	96

ExAmplar COVID-19 real-time PCR kit (L)

Description

Highly stable and ready-to-use PCR kit for COVID-19 diagnosis.

Indication

SARS-CoV-2 infection
(RdRp, E gene detection),
Coronavirus disease-19 (COVID-19).

Diagnostics Value

This is the user-friendly RT-PCR kit for COVID-19 detection which is highly stable during storage as well as delivery. It provides a very easy operation protocol and is highly acceptable to decentralized laboratories where freezing storage or cold-chain logistics is limited. This can help to diagnose SARS-COV-2 infection and to isolate the infected person so that it can be treated quickly and prevent transmission.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
ExAmplar COVID-19 real-time PCR kit (L)		●

Specification

	Viral RNA extraction kit	ExAmplar COVID-19 real-time PCR kit (L)
Sample type	Serum/Plasma/VTM/UTM	Swab, Sputum
Sample vol.	200 µL	5 µL
Assay time	<15 min	<55 min (RNA)
Storage	15-35°C	2-8°C (Lyophilized)
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	ExAmplar COVID-19 real-time PCR kit (L)
Cat #	UFPK-1	UFPK-4
Test/Kit	96	up to 48

ExAmpliar Quanti-HBV

Description

The viral loads of hepatitis B virus (HBV) in the blood are an important indicator used to determine the patient's prognosis, whether antiviral therapy is administered, and to evaluate the response after treatment. This product can perform PCR to quantify the viral load of HBV immediately by simply dissolving the dry reagent with a diluent and adding DNA extracted beforehand.

Indication

Acute and chronic HBV infection, antiviral therapy monitoring.

Diagnostics Value

Since hepatitis B is easily transmitted through asymptomatic chronic carriers and vertical infection through pregnant women is possible, it is important to prevent transmission and chronicity through appropriate antiviral treatment. Since this product maximizes user convenience and stability, it is possible to diagnose the hepatitis B viral infection quickly and accurately by quantitatively measuring the viral load of HBV even in environments where frozen shipping and frozen storage facilities are not available.

Selection guide

Item / Platform	NuActor	ExAmpliar
Viral DNA extraction kit	●	
Quanti-HBV v1.0 60T		●

Specification

	Viral DNA extraction kit	Quanti-HBV v1.0 60T
Sample type	Serum/Plasma	Serum/Plasma
Sample vol.	200 µL	5 µL
Assay time	<13 min	<30 min (DNA)
Storage	15-35°C	-20°C
Shelf life	12 months	12 months

Ordering information

	Viral DNA extraction kit	Quanti-HBV v1.0 60T
Cat #	UFPK-2	RAB-003
Test/Kit	96	60 (+36 controls)

ExAmpliar Quanti-HCV

Description

The titer of hepatitis C virus (HCV) in the blood is an important index used to evaluate the therapeutic effect of direct-acting antiviral therapy (DAA) for confirming the prognosis of patients and for treatment. This product can perform PCR to quantify the viral load of HCV immediately by simply dissolving the dry reagent with a diluent and adding DNA extracted beforehand.

Indication

Acute and chronic HCV infection, DAA (Direct Acting Antivirals) therapy monitoring.

Diagnostics Value

DAA treatment, which has recently emerged as a treatment therapy for hepatitis C, is very expensive and has side effects such as drug interaction and drug resistance. By maximizing user convenience and storage, this product helps to identify hepatitis C virus infection and maximize treatment effects by quantitatively measuring the viral load of HCV easily and reliably even in environments where frozen shipping and frozen storage facilities are not available.

Selection guide

Item / Platform	NuActor	ExAmpliar
Viral RNA extraction kit	●	
Quanti-HCV v1.0 60T		●

Specification

	Viral RNA extraction kit	Quanti-HCV v1.0 60T
Sample type	Serum/Plasma	Serum/Plasma
Sample vol.	250 µL	2.5 µL
Assay time	<13 min	<45 min (RNA)
Storage	15-35°C	-20°C
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	Quanti-HCV v1.0 60T
Cat #	UFPK-1	RAC-004
Test/Kit	96	60 (+36 controls)

ExAmplar Quanti-HIV

Description

Highly stable, ready-to-use, and sensitive PCR kit for diagnosis of HIV infection and monitoring of patients under ART (Anti-Retroviral Treatment) by HIV viral load quantification.

Indication

HIV infection and/or monitoring patients under ART (Anti-Retroviral Treatment)
> HIV infection, antiretroviral treatment (ART) monitoring.

Diagnostics Value

Since this product is an RT-PCR kit for HIV detection that maximizes user convenience and storage, it is possible to quantitatively measure the viral load of HIV easily and reliably even in areas where frozen delivery is difficult or in an environment that does not have a freezer storage facility.

Selection guide

Item / Platform	NuActor	ExAmplar
Viral RNA extraction kit	●	
Quanti-HIV v1.0 60T		●

Specification

	Viral RNA extraction kit	Quanti-HIV v1.0 60T
Sample type	Serum/Plasma	Serum/Plasma
Sample vol.	200 µL	5 µL
Assay time	<13 min	<45 min (RNA)
Storage	15-35°C	-20°C
Shelf life	12 months	12 months

Ordering information

	Viral RNA extraction kit	Quanti-HIV v1.0 60T
Cat #	UFPK-1	RAC-013
Test/Kit	96	60 (+36 controls)

ExAmpliar MTB kit v2.0

Description

Tuberculosis is a common and fatal respiratory infectious disease caused primarily by infection with *Mycobacterium tuberculosis* (MTB). MTB complexes such as *M. bovis*, *M. africanum*, *M. microti*, *M. caprae* and *M. canettii*, including *M. tuberculosis*, have 99.9% similar genomes to each other. Therefore, rapid and accurate detection of infection causative bacteria is very important for the prevention of transmission of tuberculosis and selection of an appropriate antituberculosis agent for treatment. This product can perform PCR to diagnose the MTB complex immediately by simply dissolving the dry reagent with a diluent and adding DNA extracted beforehand.

Indication

Mycobacterium tuberculosis complex
(*M. tuberculosis*, *M. bovis*, *M. africanum*,
M. microti, *M. caprae* and *M. canettii*)

Diagnostics Value

This product is a PCR kit for MTB complex detection that maximizes user convenience and storage. Users can stably store the product even in an environment that does not have a refrigeration facility for delivery and storage. It is possible to accurately identify *M. tuberculosis* by conducting an examination.

Selection guide

Item / Platform	NuActor	ExAmpliar
Sputum DNA Extraction Kit	●	
MTB v2.0 72T kit		●

Specification

	Sputum DNA Extraction Kit	MTB v2.0 72T kit
Sample type	Sputum	Sputum
Sample vol.	250 µL	2.5 µL
Assay time	<13 min	<30 min (DNA)
Storage	15-35°C	2-8°C (Lyophilized)
Shelf life	12 months	12 months

Ordering information

	Sputum DNA Extraction Kit	MTB v2.0 72T kit
Cat #	UFPK-3	UFPK-5
Test/Kit	96	72 (+24 controls)

Catalogue number index

Test items/Molecular diagnostics/ Instruments/Accessories	108
Controls/Calibrators	109

Test items

Test item	Catalog No.
	ichroma™
Adeno	CFPC-96
AFP	i-CHROMA AFP-25
AFP Plus	CFPC-73
AMH	CFPC-89
Anti-CCP Plus	CFPC-97
Anti-HBs	CFPC-52
Anti-HCV	CFPC-31
ASO	CFPC-46
BNP	CFPC-121
Calprotectin	CFPC-83
Cardiac Triple	CFPC-78
CEA	13013
CEA Plus	CFPC-72
CK-MB	CFPC-33
Cortisol	CFPC-24
COVID-19 Ab	CFPC-114
COVID-19 Ag	CFPC-115
COVID-19/Flu Ag Combo	CFPC-117
COVID-19 nAb	CFPC-120
CRP	i-CHROMA CRP-25
Cystatin C	CFPC-43
D-Dimer	CFPC-25
Dengue IgG/IgM	CFPC-60
Dengue NS1 Ag	CFPC-62
Ferritin	CFPC-32
FSH	CFPC-35
H. pylori SA	CFPC-81
HbA1c	CFPC-38
HBsAg	CFPC-29
hsCRP	CFPC-6
iFOB Neo	CFPC-15-1
IGRA-TB	CFPC-86/CFPC-86-1
IL-6	CFPC-116
Influenza A+B	CFPC-61
Influenza A+B/RSV	CFPC-80
LH	13010
Microalbumin	i-CHROMA MAU-25
Mycoplasma	CFPC-94
Myoglobin	CFPC-37
NORO	CFPC-76
NT-proBNP	CFPC-77
PCT	CFPC-23-1
PCT Plus	CFPC-64
PRL	CFPC-27
Progesterone	CFPC-21
PSA	i-CHROMA PSA-25
PSA Plus	CFPC-71
RF IgM	CFPC-39
ROTA	CFPC-75
Rota/Adeno	CFPC-79
RSV	CFPC-88
ST2	CFPC-100
Strep A	CFPC-74

Test item	Catalog No.
	ichroma™
T3	CFPC-44
T4	CFPC-26
Testosterone	13012
Tn-I	13011
Tn-I Plus	CFPC-65
Total IgE	CFPC-91
Toxo IgG/IgM	CFPC-112
Troponin T	CFPC-122
TSH	CFPC-22
TSH Plus	CFPC-45
Vitamin D	CFPC-47
β-HCG	CFPC-36
β-HCG Plus	CFPC-66

Instruments

Instrument	Catalog No.
ichroma™ II	FPRR021
ichroma™ III	FPRR037
ichroma™-50	FPRR022
ichroma™ M2	FPRR031
ichroma™ M3	FPRR035
hemochroma PLUS	FPRR016
i-chamber	FPRR009

Molecular Diagnostics

Molecular Diagnostics	Catalog No.
NuActor	FPRR032
ExAmplar	FPRR030
NuActor Viral RNA extraction kit	UFPK-1
NuActor Viral DNA extraction kit	UFPK-2
NuActor sputum DNA extraction kit	UFPK-3
ExAmplex COVID-19 PCR 3-gene Lyo kit	UFPK-7
ExAmplex COVID-19 PCR 3-gene kit	UFPK-9
ExAmplex Influenza A/B SARS-CoV-2 PCR kit	UFPK-6
ExAmplar COVID-19 real-time PCR kit (L)	UFPK-4
ExAmplar Quanti-HBV v1.0 60T	RAB-003
ExAmplar Quanti-HCV v1.0 60T	RAC-004
ExAmplar Quanti-HIV v1.0 60T	RAC-013
ExAmplar MTB kit v2.0 72T	UFPK-5

Controls

Control item	Catalog No.
Boditech AFP Control	CFPO-248
Boditech AFP Plus Control	CFPO-249
Boditech AMH Control	CFPO-214
Boditech Anti-CCP Plus Control	CFPO-288
Boditech Anti-HBs Control	CFPO-144
Boditech Anti-HCV Control	CFPO-143
Boditech BNP Control	CFPO-304
Boditech Calprotectin Control	CFPO-221
Boditech Cardiac Control	CFPO-98
Boditech Cardiac Triple Control	CFPO-204
Boditech CEA Control	CFPO-246
Boditech CEA Plus Control	CFPO-247
Boditech CK-MB Control	CFPO-243
Boditech Cortisol Control	CFPO-236
Boditech COVID-19 Ab Control	CFPO-292
Boditech COVID-19 Ag Control	CFPO-293
Boditech COVID-19 nAb Control	CFPO-303
Boditech COVID-19/Flu Ag Control	CFPO-298
Boditech CRP Control	CFPO-100
Boditech D-Dimer Control	CFPO-101
Boditech Dengue IgG/IgM Control	CFPO-280
Boditech Dengue NS1 Ag Control	CFPO-282
Boditech Ferritin Control	CFPO-99
Boditech FSH Control	CFPO-230
Boditech H.pylori Ag Control	CFPO-222
Boditech HbA1c Control	CFPO-96
Boditech HBsAg Control	CFPO-142
Boditech hCG Control	CFPO-232
Boditech hCG Plus Control	CFPO-233
Boditech Hormone Control	CFPO-95
Boditech IGRA-TB Control	CFPO-294
Boditech iFOB Neo Control	CFPO-14
Boditech IL-6 Control	CFPO-296
Boditech LH Control	CFPO-234
Boditech MAU Control	CFPO-4
Boditech Myoglobin Control	CFPO-244
Boditech NORO Control	CFPO-165
Boditech NT-proBNP Control	CFPO-245
Boditech PCT Control	CFPO-97
Boditech PCT Plus Control	CFPO-225
Boditech PRL Control	CFPO-226
Boditech Progesterone Control	CFPO-238
Boditech PSA Control	CFPO-250
Boditech PSA Plus Control	CFPO-251
Boditech RF IgM Control	CFPO-103
Boditech Rota Control	CFPO-170
Boditech Rota/Adeno Control	CFPO-164
Boditech ST2 Control	CFPO-289
Boditech T3 Control	CFPO-240
Boditech T4 Control	CFPO-237
Boditech Testosterone Control	CFPO-239
Boditech Tn-I Control	CFPO-241
Boditech Tn-I Plus Control	CFPO-212
Boditech Total IgE Control	CFPO-219

Control item	Catalog No.
Boditech Troponin T Control	CFPO-306
Boditech TSH Control	CFPO-228
Boditech TSH Plus Control	CFPO-229
Boditech Tumor Marker Control	CFPO-94
Boditech Vitamin D Control	CFPO-102

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
LITERATŪROS ŠARŠAS

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. Mol Immunol 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. Circulation 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C-reactive protein and Lipid Screening. Clin. Chem. 2001; 47:28-30.
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. Clin. Chem. 2001; 47(3): 403-411.
5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. Circulation 1999; 99:855-860.
6. Taubes G. Does inflammation cut to the heart of the matter? Science 2002; 296:242-245.
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. N Engl J Med 2000;342(12): 836-843.
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. Clin Chem 1999; 45:1676-1678.
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of-care testing. Clin Chim Acta 2005; 356:172-177.
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J. Lab. Clin Med 1976;87:120-128.

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ichroma™ CRP

INTENDED USE

ichroma™ CRP is a fluorescence Immunoassay (FIA) for the quantitative determination of CRP in human whole blood/serum/ plasma. It is useful as an aid in management and monitoring of autoimmune diseases and infectious processes, such as rheumatoid arthritis.

For *in vitro* diagnostic use only.

INTRODUCTION

The C-Reactive Protein (CRP) is synthesized by the liver in response to interleukin-6 and well known as one of the classical acute-phase reactants and as a marker of inflammation. CRP is the first acute-phase protein to be described and is an exquisitely sensitive systemic marker of inflammation and tissue damage. The acute-phase response comprises the nonspecific physiological and biochemical responses of endothermic animals to most forms of tissue damage, infection, inflammation, and malignant neoplasia. The serum CRP level may rise from a normal level of <5 mg/L to 500 mg/L during the body's general, non-specific response to infectious and other acute inflammatory events. For some time, the measurement of CRP concentration has been used as a clinical tool for monitoring autoimmune diseases and infectious processes, such as rheumatoid arthritis.^{1,2}

PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes, and migrate onto nitrocellulose matrix to be captured by the other immobilized-antibodies on test strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for ichroma™ tests to show CRP concentration in the sample.

COMPONENTS

ichroma™ CRP consists of 'cartridges', 'detection buffers' and 'sample collectors'.

- The cartridge contains the membrane called a test strip which has anti human CRP at the test line, and streptavidin at the control line. All cartridges are individually sealed in an aluminum foil pouch containing a desiccant in a box.
- The detection buffer contains anti-human CRP fluorescence conjugate, BSA-biotin-fluorescence conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide in phosphate buffered saline (PBS) as a preservative.
- The detection buffer is pre-dispensed in a tube. Detection buffers are packaged in detection buffer box and further

packed in a Styrofoam box with ice-pack for the shipment.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instruction for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge, detection buffer and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges or detection buffer tubes. A cartridge should be used for testing one sample only. A detection buffer tube should be used for processing of one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use the cartridge, if pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- Allow the cartridge, detection buffer and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for ichroma™ tests may generate slight vibration during use.
- Used cartridges, detection buffers and sample collectors should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- **ichroma™ CRP** will provide accurate and reliable results subject to the below conditions.
 - Use **ichroma™ CRP** should be used only in conjunction with instrument for ichroma™ tests.
 - Have to use recommended anticoagulant sample.

Recommended anticoagulant

K₂ EDTA, K₃ EDTA, Sodium heparin

STORAGE AND STABILITY

Component	Storage condition	
	Storage Temperature	Shelf life
Cartridge	2 - 30 °C	20 months
Detection buffer	2 - 8 °C	20 months

- After the cartridge pouch is opened, the test should be performed immediately.

LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is the most common if the epitope is masked by some unknown

components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.

- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

MATERIALS SUPPLIED

[REF] i-CHROMA CRP-25

Components of **ichroma™ CRP**

- Cartridge Box:
 - Cartridge 25
 - ID Chip 1
 - Instruction for Use 1
 - Sample Collector 25
- Box containing Detection Buffer
 - Detection buffer 25

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **ichroma™ CRP**.

Please contact our sales division for more information.

- Instrument for **ichroma™** tests
 - **ichroma™ Reader** **[REF]** FR203
 - **ichroma™ II** **[REF]** FPRR021
 - **ichroma™-50** **[REF]** FPRR022
 - **ichroma™ III** **[REF]** FPRR037
- Printer **[REF]** FPRR007
- **Boditech CRP Control** **[REF]** CFPO-100

SAMPLE COLLECTION AND PROCESSING

The sample type for **ichroma™ CRP** is human whole blood/serum /plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.
- Samples may be stored for a week at 2-8°C prior to being tested. If testing will be delayed more than a week, samples should be frozen at -20°C.
- Samples stored frozen -20°C for 6 months showed no performance difference.
- However, the whole blood sample should not be kept in a freezer in any case.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.

TEST SETUP

- Check the contents of **ichroma™ CRP**: Sealed Cartridges, Detection buffers, Sample collectors, ID Chip and Instruction for use.
- Ensure that the lot number of the cartridges matches that of the detection buffers as well as an ID chip.
- If the sealed cartridge and the detection buffer have been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for **ichroma™** tests.
 (Please refer to the 'Instrument for **ichroma™** tests Operation Manual' for complete information and operating instructions.)

TEST PROCEDURE

► **ichroma™ Reader/ ichroma™ II/ ichroma™ III**

<Multi Mode>

- 1) Make a puncture on the top of a detection buffer by inserting an empty sample collector.
- 2) Draw 10 µL (Human whole blood/serum/ plasma/control) of sample with a sample collector.
- 3) Assemble the sample collector and the detection buffer into one.
- 4) Shake 10 times or more until the sample out of the sample collector by inversion. The mixture of the buffer and the sample has to be used within 30 seconds.
- 5) Remove the cap off the top of assembled tube. Discard two drops of reagent onto the paper towel before applying to a cartridge.
- 6) Load only two drops of the mixture onto the sample well of the cartridge.
- 7) Leave the cartridge at room temperature for 3 min before inserting the device into the holder.
 ⚠ Scan the sample-loaded cartridge immediately when the incubation time is over. If not, it will cause inaccurate test result.
- 8) To scan the sample-loaded cartridge, insert it into the cartridge holder of the instrument for **ichroma™** tests. Ensure proper orientation of the cartridge before pushing it all the way inside the cartridge holder. An arrow is marked on the cartridge especially for this purpose.
- 9) Press the 'Select' or Tap the 'START' button on the instrument for **ichroma™** tests to start the scanning process.
- 10) The instrument for **ichroma™** tests will start scanning the sample-loaded cartridge immediately.
- 11) Read the test result on the display screen of the instrument for **ichroma™** tests.

<Single Mode >

- 1) The test procedure is same with "Multi test 1) – 4)".
- 2) Remove the cap off the top of assembled tube. Discard two drops of reagent onto the paper

3 t.s.

towel before applying to a cartridge.

- 2) Load only two drops of the mixture onto the sample well of the cartridge.
- 3) Inserting the device into the holder of the instrument for ichroma™ tests. Ensure proper orientation of the cartridge before pushing it all the way inside the cartridge holder. An arrow is marked on the cartridge especially for this purpose.
- 4) Press the 'Select' or Tap the 'START' button on the instrument for ichroma™ tests.
- 5) The cartridge goes inside the Instrument for ichroma™ tests and will automatically start scanning the sample-loaded cartridge after 3 min.
- 6) Read the test result on the display screen of the instrument for ichroma™ tests.

► **ichroma™-50**

- 1) Insert the tip array in the tip station.
- 2) Insert the detection buffer array in the reagent station and cover the reagent station.
- 3) Open the detection buffer and insert it in the detection buffer station.
- 4) Open the cover of the magazine station and pull and lift the cartridge magazine.
- 5) Insert the cartridges in the cartridge magazine individually.
- 6) Insert the cartridge loaded cartridge magazine into the magazine station and close the cover of the magazine station.
- 7) Insert the sample tube into the blood collection tube rack and load the blood collection tube rack into the sampling station (loading part).
- 8) Tap the button located in the upper side of the No. of test cartridge region to select ID chip what you want to use.
- 9) When the selected cartridge slot is activated, set the number of test cartridge by tapping.
- 10) Tap the button located in the upper side of the No. of reagent region to select ID chip what you want to use.
- 11) When the selected slot is activated, set the number of detection buffer by tapping.
- 12) Set the number of pipette tips by tapping.
- 13) Tap the 'START' button on the left upper of the main screen to start test.
 (Please refer to the ichroma™-50 operation manual for complete information and operation instructions.)

INTERPRETATION OF TEST RESULT

- Instrument for ichroma™ tests calculates the test result automatically and displays CRP concentration of the test sample in terms of mg/L.
- **The cut-off (reference value): 10 mg/L**
- Working range: 2.5-300 mg/L
- Effect of Hematocrit
 The CRP Whole Blood of ichroma™ Reader is calibrated to read the CRP serum concentration of a blood sample with a

hematocrit of 40%. If the actual hematocrit value deviates from 40%, the result should be corrected by multiplying with the respective factor in the table: deviates from 40%, the result should be corrected by multiplying with the respective factor in the table:

Hct (%)	Factor	Hct (%)	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6
43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

Reference range, HCT:

- Woman: 35-44%
- Men: 39-48%

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **ichroma™ CRP**. For more information regarding obtaining the control materials, contact [Boditech Med Inc.'s Sales Division for assistance](#).
 (Please refer to the instruction for use of control material.)

PERFORMANCE CHARACTERISTICS

- **Analytical sensitivity**
 Analytical sensitivity of **ichroma™ CRP** was determined by testing 10 times with three lots of reagents. Analytical sensitivity of **ichroma™ CRP** system was 0.5 mg/L.
- **Specificity**
 Biomolecules such as hemoglobin, CEA, AFP, ALT, Troponin I, CK-MB, Albumin and serum amyloid P were added to the test sample(s) at concentration much higher than their normal physiological levels in the blood. ichroma CRP test results did not show any significant cross-reactivity with these biomolecules.
- **Precision**

- Intra-assay
 The intra-assay precision was calculated by one evaluator, who tested different concentration of control standard twenty times each with three different lots of **ichroma™ CRP**.
- Inter-assay
 The inter-assay precision was confirmed by 3 different evaluators with 3 different lots, testing ten times each different concentration for ten days.


CRP (mg/L)	Intra-assay			Inter-assay		
	Mean	SD	CV (%)	Mean	SD	CV (%)
5.0	5.06	0.3	6.7	5.11	0.4	7.1
40.0	39.95	3.0	7.6	40.18	3.1	7.6
150.0	150.46	5.9	3.9	150.15	4.5	3.0

■ **Linearity**

The high pool (300 mg/L) was diluted with the low pool (2.5 mg/L) to the following final percentages; 100%, 75%, 50%, 25%, 10%, 5% and 0%. Sample was assayed in triplicate in one analytical run at each CRP level. The coefficient of linear regression was R=0.997.

■ **Comparability**

CRP concentrations of 166 samples were quantified independently with **ichroma™ CRP (ichroma™ Reader)** and Comparator A as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and coefficient of correlation (R). Linear regression and coefficient of correlation between the two tests were $Y=0.99X + 1.0081$ and $R = 0.9933$ respectively.

Note: Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance; please contact:
Boditech Med Inc.'s Technical Services
 Tel: +82 33 243-1400
 E-mail: sales@boditech.co.kr

REFERENCES

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. *J Clin. Invest* 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. *Mol Immunol* 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. *Circulation* 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C-reactive protein and Lipid Screening. *Clin. Chem.* 2001; 47:28-30.
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. *Clin. Chem.* 2001; 47(3): 403-411.
5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. *Circulation* 1999; 99:855-860.
6. Taubes G. Does inflammation cut to the heart of the matter? *Science* 2002; 296:242-245.
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. *N Engl J Med* 2000;342(12): 836-843.
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. *Clin Chem* 1999; 45:1676-1678.
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of-care testing. *Clin Chim Acta* 2005; 356:172-177.
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. *J. Lab. Clin Med* 1976;87:120-128.

Boditech Med Incorporated
 43, Geodudanji 1-gil, Dongnae-myeon,
 Chuncheon-si, Gang-won-do, 24398
 Republic of Korea
 Tel: +(82)-33-243-1400
 Fax: +(82)-33-243-9373
 www.boditech.co.kr

Obelis s.a
 Bd. Général Wahis 53,
 1030 Brussels, BELGIUM
 Tel: +(32)-2-732-59-54
 Fax: +(32)-2-732-60-03
 E-Mail: mail@obelis.net

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Boditech CRP Control

INTENDED USE

Boditech CRP Control is intended for *in vitro* diagnostic use in the quality control of CRP Assay Kit. **For *in vitro* diagnostic use only.**

INTRODUCTION

The use of Boditech CRP Control may be considered as an objective assessment of the precision of CRP Assay Kits and is an integral part of Good Laboratory Practices. Boditech CRP Control is provided in liquid form.

COMPONENTS

Boditech CRP Control consists of 'Boditech CRP Control level 1', 'Boditech CRP Control level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The control contains CRP antigen stock and Horse serum
- Each control vial packed in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech CRP Control should not be used past the expiration date.
- Boditech CRP Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and CRP Assay Kits.
- Human source materials from which Boditech CRP Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

STORAGE AND STABILITY

- Opened: Store refrigerated (+2 to +8 °C). Once reconstituted, stable for 28 days at (+2 to +8 °C) if kept capped in original container and free from contamination. After use, any residual product should not be returned to the original tube.
- Unopened: Store refrigerated (+2 to +8°C).
- Boditech CRP Control is stable to expiration date on the label.
- Bacterial contamination of reconstituted Boditech CRP Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the tube should be discarded and a fresh tube needs to be reconstituted.

INSTRUCTIONS FOR USE

Boditech CRP Control is supplied in liquid form.

1. Thoroughly mix the contents of the tube before each use by gently inverting for several times.
2. Follow the procedure according to the instruction for use provided with the kit.

Please refer to package inserts of the test cartridges for detailed test procedure.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

In the event of damage to the package, contact the

Boditech Med Inc.'s Technical Services.

MATERIALS SUPPLIED

REF CFPO-100

Boditech CRP Control Box (2 tubes)	
Boditech CRP Control level 1 (0.5 mL)	1
Boditech CRP Control level 2 (0.5 mL)	1
Instruction For Use	1
Control value & Barcode Sheet	1

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

For Technical Assistance

Boditech Med Inc.'s Technical Services at

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Boditech Med Inc. techninė pagalba