Improving the quality of life

The Matrix Patch™ equine pericardial patch

Our mission:


Leading in cardiovascular tissue engineering

We are a biotechnology company that focuses on the development and manufacture of biological implants. Our cell-free implants minimize the need for long-term medication, which helps to reduce possible negative side effects. This can significantly improve the quality of life, especially for younger patients.


Our technology:

Optimizing immunoreactivity by reducing DNA content


We use a proprietary glutaraldehyde-free process to decellularize and sterilize our tissue implants. This technique reduces residual DNA content to a minimum. After the process, RNA is no longer detectable. Subsequently, the transmission of pathogenic microorganisms is impossible. The complete removal of all xenogenic cells ensures significantly reduced immunoreactivity against our implants.

Native tissue consists of cells, cellular components and an extracellular matrix (ECM). The ECM consists of structural proteins which determine the mechanical properties of the tissue. The cells and cellular components are responsible for the rejection of xenogenic tissue implants.

To prevent rejection and degeneration of the implants, autotissue uses the naturally occurring detergent deoxycholic acid (DOA) during processing to remove all cells and cellular components from the tissue.

Decellularized tissue is characterized by the absence of cell nuclei and a very low residual DNA and RNA content. Due to the very low RNA and DNA content, autotissue implants reveal a high degree of viral safety.

Case Report of Prof. Dr. Boulos Asfour: Deutsches Kinderherzzentrum

Diagnosis

Critical congenital aortic stenosis

• May 2005

Catheter intervention directly after birth:

Balloon valvuloplasty.

• October 2006

1st surgery at age of 1 year:

Aortic valve repair.

Replacement of right coronary cusp with autologous pericardium.

• September 2012

2nd surgery at age of 7 years:


Redo aortic valve repair due to failure of the autologous pericardial patch. Replacement

of the right coronary cusp with a Matrix Patch equine.

• September 2017

3rd surgery at age of 12 years:

Aortic valve replacement with a mechanical valve (On-X 21 mm) due to retraction of the two native cusps. The Matrix Patch is still pliable and macroscopically free of calcifications.

Intraoperative view of the opened aorta demonstrating the pliable and nicely grown-in Matrix Patch™ equine in situ after 5 years.

Macroscopic view of half of the explanted right-coronary patch (Matrix Patch™ equine).

The von-Kossa staining confirms complete absence of any calcification.

The von-Willebrand staining shows a thin endothelial layer.

Vergnat M, Asfour B, Arenz C et al. Aortic stenosis of the neonate: A single-center experience. *J Thorac Cardiovasc Surg.* 2019;157:318-326.e1.

Vergnat M, Asfour B, Arenz C et al. Contemporary results of aortic valve repair for congenital disease: lessons for management and staged strategy. *Eur J Cardiothorac Surg.* 2017;52:581-587.

Hraška V, Sinzobahamvya N, Haun C et al. The long-term outcome of open valvotomy for critical aortic stenosis in neonates. *Ann Thorac Surg.* 2012;94:1519-1526.

Histology: Courtesy of Prof. Dr. Matthias Sigler
Department of Cardiac Surgery, University of Göttingen.

The Matrix Patch™ is available in the following sizes

Nr. 37	• 10 x 5 cm	Article No. 42 505 306 041 55
	• 10 x 2 cm	Article No. 42 505 306 041 24
Nr. 36	• 8 x 8 cm	Article No. 42 505 306 048 89
Nr. 35	• 8 x 4 cm	Article No. 42 505 306 048 41
Nr. 34	• 5 x 5 cm	Article No. 42 505 306 040 56
	• 3 x 3 cm	Article No. 42 505 306 043 39

Visit our new website
for videos and presentations:

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
Matrix Patch^{CE}₀₁₉₇

equine pericardial patch

autotissue 
 leading in cardiovascular tissue engineering


The equine Matrix Patch™ is a cell-free pericardial patch for use in cardiac surgery. Compared to native or conventional pericardial patches treated with glutaraldehyde, the patch shows similar tear strength but is considerably more flexible and features lower tendency for bleeding from stitch channels.

Aortic valve reconstruction in a patient with congenital aortic valve stenosis.
 Courtesy of MU Dr Matej Nosál'
 BC Children's Heart Centre, Bratislava, Slovakia

- Nr. 34
- Nr. 35
- Nr. 36
- Nr. 37
- Nr. 38


The equine Matrix Patch™ is manufactured using a unique technology without the use of glutaraldehyde, in which the accumulation of calcium and fat in tissue is reliably prevented.

Plastic enlargement of the ascending aorta and pulmonary artery surgery in a patient with Williams-Beuren Syndrome.
 Courtesy of Prof. Dr Mikhail Abramyan
 Research Centre for Obstetrics, Gynecology and Perinatology, Moscow, Russia

Pericardium has been used in heart surgery for more than 40 years. The primary advantages of decellularised equine pericardium include outstanding mechanical properties, biocompatibility and the absence of transferable diseases. Moreover, no calcification or shrinkage has been observed over the years.

To date, the equine Matrix Patch™ has been successfully used in more than 5500 surgeries in Europe, Turkey and Russia. It is delivered sterile in antibiotic solution and can be stored refrigerated at 4-8°C for up to 18 months. Thus, it is ready for use at any time without lengthy preparation and can be used immediately after one minute rinsing in physiological saline to remove the antibiotic solution.

Tricuspid valve reconstruction.
 Courtesy of Prof. Dr Robert Cesjnevár
 Paediatric Cardiac Surgery at the Erlangen University Hospital

Nr. 38

Every equine Matrix Patch™ is inspected by trained specialists in the quality assurance department. Validated testing monitors and ensures biocompatibility, absence of cells and sterility of the product.

The equine Matrix Patch™ is 0.17 to 0.37 mm thick and comes in the following sizes:

10 x 10 cm, 10 x 5 cm, 10 x 2 cm, 8 x 8 cm, 8 x 4 cm, 8 x 2 cm, 8 x 1 cm, 6 x 1 cm, 5 x 5 cm and 3 x 3 cm.

TÜV Rheinland LGA Products GmbH certifies the approval of the equine Matrix Patch™ in all available sizes. Auto Tissue Berlin GmbH is also a certified manufacturer in accordance with EN ISO 13485.

Auto Tissue Berlin GmbH
 Goerzallee 305d
 14167 Berlin
 Germany
 Email: info@autotissue.de

Tel: +49 (0)30 3302 925 80
 Fax: +49 (0)30 3302 925 89
 Email: info@autotissue.de
 Internet: http://www.autotissue.de

Management System
 EN ISO 13485:2016
 www.tuv.com
 ID 0000038078

Version 2020-10-27

Matrix Patch™

Equine pericardial patch

Instructions for Use

English (en)

Auto Tissue Berlin GmbH

Goerzallee 305 d

14167 Berlin

Germany

Phone: +49 30 330 292 580

Fax: +49 30 330 292 589

	Caution! Read this instruction for use carefully!
	Caution! Store refrigerated at 4 to 8 °C!
	Caution! Matrix Patch™ equine pericardial patch must not be used if packaging has been damaged or seal is broken!
	Caution! Matrix Patch™ equine pericardial patch must not be re-sterilised!
	Caution! Matrix Patch™ equine pericardial patch for single use only!
STERILE*	Chemically sterilised and validated according to DIN EN ISO 14160

Device description

The Matrix Patch™ equine pericardial patch is a decellularized structure made from equine pericardial tissue, which is manufactured in accordance with a patented procedure of Auto Tissue Berlin GmbH using controlled time, pH and temperature conditions. The treatment kills and eliminates all living cells (connective tissue cells and endothelial cells, as well as bacteria, viruses, fungi and mycoplasma). This produces a matrix that is made predominantly of collagen and elastin. The decellularization process reduces the DNA and RNA content of the tissue so that the possibility of transmitting endogenous retroviruses to the recipient is excluded.

Note: Any slight changes in the colour of the tissue have no influence on the quality and performance of the device.

Magnetic resonance imaging

The Matrix Patch™ equine pericardial patch contains no metallic or ferromagnetic materials and therefore poses no additional risk for patients undergoing MRI examinations.

Available models

The Matrix Patch™ equine pericardial patch is available in the following sizes:

10 cm x 10 cm, 8 cm x 8 cm, 10 cm x 5 cm, 8 cm x 4 cm, 5 cm x 5 cm, 10 cm x 2 cm, 8 cm x 2 cm, 4 cm x 4 cm, 3 cm x 3 cm, 8 cm x 1 cm, 6 cm x 1 cm and 2 cm x 2 cm.

Special features of this prosthesis are given in the section "device description".

The serial number of each Matrix Patch™ equine pericardial patch is given on the container and on separate identification labels provided within the package.

Expiry date

The expiry date of the prosthesis is 18 months from the date of manufacture as stated on the package label.

Delivery

The Matrix Patch™ equine pericardial patch is delivered sterile in a sealed container. The container together with identification labels, the implant card, Instructions for Use and Freeze Watch™ temperature indicator are packed in a cardboard box. Every cardboard box is transported in a Styrofoam outer package with a second temperature indicator.

The Matrix Patch™ equine pericardial patch is shipped and stored in an antibiotic preservation solution.

Temperature indicator

Each shipping unit contains two different kinds of temperature indicator. First one in the outer package for high and low temperatures to track the temperature to which the device was exposed during transportation.

The temperature indicator in the Styrofoam outer package is intended only for transport of the patch and is not suitable for monitoring temperatures during storage. Immediately after delivery the customer should remove the Styrofoam and check the temperature probe. In case the sensor has been activated, the patch must not be used. In this case, please contact Auto Tissue Berlin GmbH or the local distributor for return and replacement of the device.

Note: In cases in which more than 48 hours have elapsed since delivery, and activation of the temperature sensor is noticed, the patch can be replaced only at the expense of the customer as storage in this case is in the control of the customer.

In every cardboard box, a second temperature indicator (Freeze Watch™) monitors the temperature during storage. Further information is given in the section "Preparation of the patches for implantation".

Storage

To minimise the risk of contamination and to protect the patch please keep the prosthesis and the instructions for use until use in the cardboard box. Store dry and protect against contamination in a refrigerator at 4 to 8 °C. The customer is advised to check the stock on a regular basis to make sure that prostheses are used before the expiry date given on the cardboard box.

Caution! Do not freeze! Store this device always in a dry refrigerator. Any prosthesis that has been frozen or may have been frozen must not be implanted in patients.

Sterility

The Matrix Patch™ equine pericardial patch is sterile as long as the container is not damaged, and the storage solution is clear. The external surface of the container is not sterile and must not come in contact with sterile instruments.

Caution! Patches removed from their container, but not implanted, are no longer sterile and must not be used at a later date. **Patches must not be re-sterilised!** Resterilisation can change the properties of the device, whereby a safe application of the device can no longer be guaranteed. Do not use the patch if container has been opened or damaged.

Intended use

The Matrix Patch™ equine pericardial patch is intended for the extension, augmentation and repair of tissue defects during cardiovascular surgery.

Indications for use

The Matrix Patch™ equine pericardial patch is indicated for congenital and acquired disease for the extension, augmentation and repair of defects during cardiac, arterial and venous surgical operations such as augmentation and repair of the pulmonary artery, use as a transannular patch, reconstruction of heart valves, augmentation and repair of veins or arteries and closure of intracardiac septal defects.

Contraindications

Since the choice of which graft material to use and the use in each patient is of utmost importance, a physician should analyse the risk / benefit ratio for the application of a pericardial patch based on the condition of each patient and his / her symptoms.

In patients with allergies to Penicillin, Streptomycin or Amphotericin, the Matrix Patch™ equine pericardial patch should be used only after careful consideration of the advantages of the Matrix Patch™ equine pericardial patch. In these patients, it is strongly recommended that the saline rinsing procedure, which is a standard procedure for all tissue implants before implantation, be intensified.

General considerations

Preparation of the patches for implantation

The Freeze Watch™ inside the package must be checked. Dark staining of the Freeze Watch™ indicates that the patch was exposed to temperatures below the acceptable storage temperature. The patch must not be used if the Freeze Watch™ is activated (darkly stained).

The Matrix Patch™ equine pericardial patch is delivered sterile in a container. Before opening this container, it should be carefully inspected for leakage or damage in the form of cracks. Patches from damaged or leaking containers or containers without sufficient solution to completely immerse the Matrix Patch™ equine pericardial patch should not be used.

Unscrew the container to take the Matrix Patch™ equine pericardial patch out of the container. To avoid contamination, the contents of the container must be handled aseptically. Once the container has been opened, the patch must be used or disposed of. Re-use of the patch is not permitted. The contents of the container are sterile, and it is recommended that the Patch be handled with sterile atraumatic forceps.

Caution! Do not touch the Matrix Patch™ equine pericardial patch itself with surgical forceps. The patches must be handled with atraumatic forceps. During the handling of the tissue, extreme caution is necessary.

Thorough rinsing of the bioprosthesis in physiologic saline before the implantation is necessary to clear the tissue of antibiotics in the storage solution. For rinsing, the patch should be placed for 1 minute in 500 cc of sterile normal saline. The patch should be immersed completely in the rinsing solution, and the rinsing solution should be discarded after rinsing the patch is complete.

Caution! No chemical solutions, drugs, chemical compounds, antibiotics etc. should be added to the storage and rinsing solution as changes to the tissue may result, which may have an adverse effect on long-term performance on the patch.

Caution! Always keep the patch moist to prevent tissue dehydration.

Implantation

All implants require special care and precise handling.

If the Matrix Patch™ equine pericardial patch is not sterile, is damaged or is handled incorrectly, it may not be used for human implantation.

The implantation should be carried out using a suitable suture technique. For vascular surgery, the patch must be oriented such that the smooth side is aligned to the bloodstream. Because of the complexity and diversity of surgical procedures in which the Matrix Patch™ equine pericardial patch may be used, it is up to the surgeon to determine which surgical technique to choose after weighing the warnings, precautions and relative risks and benefits of particular implantation techniques.

Standard intra-operative anticoagulation appropriate for the standard of care of cardiac surgery at each site should be used. No additional anti-thrombotic therapy, beyond standard anticoagulation for cardiac surgery, is necessary when using the Matrix Patch™ equine pericardial patch. No anti-thrombotic therapy need be stopped after using the patch.

There are no reports of negative interactions between drugs and the Matrix Patch™ equine pericardial patch.

Follow-up examination

The Matrix Patch™ equine pericardial patch does not require a device specific follow-up examination. The necessity and frequency of follow-up examinations should be determined by the underlying pathology of the patient in consultation with the surgeon who implanted the patch.

Explantation

Due to the diverse applications of the Matrix Patch™ equine pericardial patch, Auto Tissue Berlin GmbH cannot make any general recommendations regarding the procedure of explanting the patch.

Disposal

The Matrix Patch™ equine pericardial patch may be disposed of as regular biological waste.

Complications

All cardiac interventions bear the risk of severe complications that may lead to a patient's death.

Possible complications during or following use of cardiovascular implants include, but are not limited to, the following:

- Stitch channel bleeding
- Suture pull out
- Rupture
- Thrombosis
- Inflammatory reaction
- Thickening
- Shrinking
- Calcification
- Aneurysms formation
- Infection

The cardiovascular surgeon should provide each patient or his legal representative with appropriate information.

Documentation

For documentation at the customers every Matrix Patch™ equine pericardial patch is delivered with sic (6) labels.

Please fill out the customer questionnaire and send or fax it to Auto Tissue Berlin GmbH.

Counselling: For any other technical or surgical questions, please contact:

info@autotissue.de

Product specialists from Auto Tissue Berlin GmbH will contact you as soon as possible.

Warranty

Auto Tissue Berlin GmbH guarantees that the Matrix Patch™ equine pericardial patch is designed, manufactured and packaged in accordance with its EN ISO 13485-certified quality management system with all possible care, using the most suitable state-of-the-art processes and applying the principle of integrating safety into its design and manufacture. This system is designed to guarantee the safe use of device when it is employed under the conditions and for the purposes for which it is intended, in accordance with precautions described in the previous paragraphs, and to minimise the risks associated with the use of the device as far as possible, although they cannot be totally eliminated.

The device must only be used under the responsibility of specialist medical personnel, taking account of the unavoidable risks and the possible side effects and complications of the treatment for which, it is designed, including those referred to in other sections of these Instructions for Use.

In view of the technical complexity and delicacy of the related medical decisions and the methods of applying the device, Auto Tissue Berlin GmbH is unable to give any express or implied guarantee that the use of the device will produce good results or effectively cure a medical disorder. The results, in terms of clinical profile and the functionality and duration of the device, depend on numerous factors outside the manufacturer's control, including the patient's condition, the surgical implantation or procedures used, and the way in which the device is handled after the packaging is opened.

In view of these factors, the liability of Auto Tissue Berlin GmbH is strictly limited to replacing the device if it should prove defective at the time of delivery. Under such circumstances, the customer shall deliver the device to Auto Tissue Berlin GmbH, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty solely covers replacement of a device found to be defective with another of same or equivalent type manufactured by Auto Tissue Berlin GmbH.

The warranty shall only apply if the device is returned to Auto Tissue Berlin GmbH, correctly packed, together with a detailed written report on the defects complained of and, in the case of an implanted device, the reasons for its removal from the patient.

If the device is found to be faulty and replaced, Auto Tissue Berlin GmbH will refund the costs directly incurred by the customer. Auto Tissue Berlin GmbH declines all liability in the event of failure to comply with the directions for use or precautions specified in this instruction booklet, or in the case of use after the expiry date marked on the packaging.

Auto Tissue Berlin GmbH further declines all liability for the consequences of medical decisions and methods of use or application of the device, and shall therefore not be liable for any material, biological or non-pecuniary loss or injury of any kind resulting from the application of the device.

Our current business terms and conditions apply.

Versija 2020-10-27

Matric Patch™

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Naudojimo instrukcija

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Auto Tissue Berlin GmbH

Goerzallee 305 d

14167 Berlynas

Vokietija

Telefonas: +49 30 330 292 580

Faksas: +49 30 330 292 589

autotissue 
Matrix Patch 
0197
equine pericardial patch

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