

ECO.next & ECO.LUX

**ISTRUZIONI PER L'USO
OPERATING INSTRUCTIONS
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES**

La casa costruttrice declina ogni responsabilità per le possibili inesattezze contenute nel presente manuale, imputabili ad errori di trascrizione e si riserva il diritto di apportare ai propri prodotti quelle modifiche che ritenesse necessarie o utili, senza pregiudicarne le caratteristiche essenziali.

Da conservare per futuri riferimenti.

The manufacturer declines all and any responsibility for possible inexactitudes due to typing errors contained in this manual and reserves the right to modify his products without changing their main characteristics, if he considers such modifications necessary or useful.

To be preserved for future consultation.

Le fabricant décline toute responsabilité pour les possibles inexactitudes contenues dans ce manuel, imputables à des erreurs de transcription ou traduction et se réserve le droit de modifier ses produits, sans modifier les caractéristiques essentielles, s'il le considère nécessaire ou utile.

À conserver pour l'éventuelle consultation future.

El fabricante declina toda responsabilidad por las inexactitudes que pudiera contener el presente manual, imputables a errores de transcripción y se reserva el derecho de aportar a sus productos las modificaciones que considere oportunas o necesarias sin perjudicar las características esenciales de los mismos.

Conserven este manual para futuras consultas.

Tutti i componenti in poliuretano prodotti da Tecnodent sono:

All polyurethane parts produced by Tecnodent are:

Tous les composants en polyuréthane produits par Tecnodent sont :

Todos los componentes de poliuretano producidos por Tecnodent son:

HCFC / Free

POLTRONA TIPO
CHAIR TYPE
FAUTEUIL TYPE
SILLON MODELO

ECO.next ECO.LUX

ISTRUZIONI PER L'USO
OPERATING INSTRUCTIONS
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES

TECN**DENT**
ERGONOMIC INNOVATION

INDEX

1.	Introduction	20
1.1.	Definitions.	20
1.2.	General warning.	20
1.3.	Documentation purpose.	20
1.4.	Document keeping.	21
1.5.	Limitation of responsibility.	21
1.6.	Transport and storage.	21
1.7.	Packing disposal.	21
1.8.	End-of-life equipment disposal information.	21
2.	Equipment description	22
2.1.	Intended use.	22
2.2.	General features.	22
2.3.	Warranty.	22
2.4.	Identification and labelling.	23
2.5.	Technical specifications.	23
3.	Installation	26
3.1.	Unpacking.	26
3.2.	Supply status.	26
3.3.	Handling.	26
3.4.	Installation.	27
4.	Using the product	28
4.1.	Typical use.	28
4.2.	Controls.	28
4.3.	Chair operations.	29
5.	Maintenance	30
5.1.	General information.	30
5.2.	Cleaning.	30
5.3.	Main fuses replacement.	30
5.4.	Upholsteries replacement.	31
6.	Troubleshooting	32
6.1.	Procedures.	32
6.2.	Acoustic signals.	33
7.	Accessories	34
7.1.	Headrests.	34
7.2.	Armrests.	35
7.3.	Backrests.	35
7.4.	Rotation.	35
7.5.	Programmes.	36
7.6.	Trendelenburg position.	36
7.7.	Remote control.	36
7.8.	Mobile foot controls.	36

1. INTRODUCTION

1.1. DEFINITIONS.

This manual adopts the following graphic and linguistic conventions.

NOTE.



Contains important information, which should be highlighted with respect to the text.

CAUTION.



Before conducting certain procedures, this message may be displayed. If not observed, it could cause damage to the equipment.

WARNING.



Before conducting certain procedures, this message may be displayed. If not observed, it could cause damage to the operator and the equipment.

Double function control(s).

Foot operated switches which control the multiple chair functions, depending on how they are activated. For example, continuous activation: manual chair movement; short activation (tap): automatic movement of a pre-programmed working position.

Synchronized movement.

Simultaneous and interdependent positioning of the chair parts, which may be activated in an automatic or manual setting (e.g., raising the chair and simultaneously leaning the backrest).

Compensated movement.

Position variation of two chair parts acted by a motion system which takes into consideration ergonomic needs (e.g., backrest leaning movement).

PCB.

Electronic card (printed circuit board).

CN.

Electric connector.

PWR.

Power supply.

LED.

Light-emitting diode.

1.2. GENERAL WARNING.

Upon user request, Tecnodent agrees to provide circuit diagrams, list of component parts, calibration instructions or other information available, which will assist the user's technical personnel to conduct repairs of those parts of the chair which are considered repairable by the manufacturer.

WARNING.



The manufacturer will accept responsibility for the safety, reliability and performance of the equipment on the following conditions:

- The installation operations, modifications or repairs have been performed by authorised personnel;
- The electrical system of the installation site conforms to the legal regulations and standards;
- The equipment is used in accordance with the operating instructions.

NOTE.



Please note that in accordance with Art. 14 of EEC Directive 85/374 "Responsibility for damages arising from defective products" implemented in Italy by Decreto del Presidente della Repubblica on May 24, 1988, n. 224: "Right to compensation shall cease at the expiration of ten years from the day the manufacturer or the importer within the EU Nations has first marketed the product, which has caused the claimed damage".

1.3. DOCUMENTATION PURPOSE.

This manual shows the ECO.next and ECO.LUX chairs and contains information regarding their correct use, functions, performances, maintenance, troubleshooting and corresponding solutions.

This manual is intended for the final user, the professional that uses the equipment to perform his job.

1.4. DOCUMENT KEEPING.

This manual is part of the equipment and must be conserved until its disposal. It should be stored in a protected location, but near the device, for prompt usage if necessary; if this documentation is lost, the user should request a replacement copy.

1.5. LIMITATION OF RESPONSIBILITY.

The Manufacturer must be relieved of any responsibility in the following cases:

- Improper use or usage by untrained personnel for professional use;
- Use contrary to the applicable standards;
- Improper installation, if it is not specifically included in the supply contract;
- Power supply failure or defects;
- Lack of maintenance or servicing;
- Unauthorised modification or interventions;
- Use of non-original spare parts, or not model-specific;
- Failure to comply with the instructions (even partially);
- Exceptional events.

1.6. TRANSPORT AND STORAGE.

The equipment, if contained in its original packaging, may be exposed for a period of no more than 5 weeks to environmental conditions within the following ranges:

	Min	Max
Temperature	-40° C	+70° C
Relative humidity (condensation included)	10%	95%
Atmospheric pressure	500 hPa (mbar)	1060 hPa (mbar)

1.7. PACKING DISPOSAL.

All materials used for packaging are recyclable and respect the environment:

- Wooden pallet with fumigation treatment;
- Cardboard box;
- Polythene film or bubblewrap.

Please deliver the packing to an authorized landfill, for the recovery of these materials.

1.8. END-OF-LIFE EQUIPMENT DISPOSAL INFORMATION.

According to the article 13 of D.lgs. no. 151, dated July 25th, 2005, "Implementation of the Directives 2002/95/EEC, 2002/96/EEC and 2003/108/EEC, concerning the reduction of hazardous substances in electric and electronic equipment, as well as waste disposal", and subsequent modifications such as D.lgs. no. 49, dated March 14, 2014, "Implementation of the Directive 2012/19/EEC concerning waste of electric and electronic equipment", Tecnodent states as follows.



- This symbol, placed in the identification label, indicates that the product at the end of its useful life must be disposed separately from other waste.
- The separate collection of end-of-life equipment is organized and managed by the manufacturer. The owner who wishes to dispose this equipment should contact the manufacturer or one of its representatives to allow separate collection of this end-of-life equipment.
- The correct separate collection for recycling, processing and environmentally friendly disposal helps avoiding potential adverse effects on human health, and promotes reuse and recycling of the equipment materials.
- Illegal disposal of the product by the owner will lead to penalties, as determined by the laws in force.
- If the product is placed out of use at the end of its working life, it is necessary to definitely put it out of service, by disconnecting the plug from the socket and cutting the power supply cable.

2. EQUIPMENT DESCRIPTION

2.1. INTENDED USE.

The product described by this manual is a dental patient chair.

The equipment is designed for lifting and holding the patient in a proper working position during dental diagnosis or therapy.

The chair can be also suitable for more generic use, in relation to the basic positioning of the patient. In case of any doubt, contact the manufacturer for further details.

It can be used by specialists in private clinics, hospital structures, specialized wards.

The apparatus requires the use of a compliant electrical installation, as specified in the section "installation". It is not advised for usage in close proximity with devices which monitor patient vital functions.

Essential performance of the device is holding the patient in any normal use situation, in fault condition too.

WARNING.



The equipment was developed and manufactured only for dental use. It can be also utilised for different generic fields. For more information, ask the manufacturer. It is furthermore prohibited to perform any modification of the equipment or its parts, without explicit written authorization of the manufacturer.

2.2. GENERAL FEATURES.

2.2.1. Basic version.

A patient chair model ECO.next or ECO.LUX presents the following functions:

- Completely electromechanical seat lifting/lowering and backrest tilting movements;
- Compensated backrest tilting, up to the horizontal position;
- Anatomical narrow backrest;
- ELLE 2 MOBILE type headrest;
- Complete symmetry, for right-handed and/or left-handed professionals;
- Safety STOPS on the backrest, on the base and on the lifting arm;
- Return to patient exit ("zero") position automatic programme;
- Rinse position automatic program;
- Low voltage controls (5 V).

2.2.2. Accessories.

The chairs can also be equipped with different accessories, depending on the working needs of the professional, as well as the comfort needs of the patient to be treated. In detail, the following options are available:

- "UNI tappezzato" type headrest;
- Anatomical large backrest;
- Magnetic pillow model C95;
- Magnetic pillow model C2002;
- Right or left armrest(s);
- Device for rotation around the vertical axis, installed under the seat;
- Device for programming three different working positions;
- Trendelenburg movement;
- Foot control, on the base;
- Mobile foot control;
- Infra-red remote control.

The operation and the performance of these accessories are described in chapter 7.

2.3. WARRANTY.

NOTE.



Upholstery damage caused by contact with clothes or clothing accessories (ex: belts) having metal, sharp, pointed or similar parts will not be replaced under warranty.

CAUTION.



The use of disinfectants or detergents with a strong alcoholic content on the apparatus damages the surface of the product, especially the upholsteries.

2.4. IDENTIFICATION AND LABELLING.

The identification label is placed in the base of the chair under the seat. It indicates the chair model, serial number and other information. Figure 1 highlights the label position and a facsimile.




FIG. 1

	Manufacturer.		Manufacturing date.
	Model. It identifies unambiguously the device of the product family.		Serial number of the device.
	Medical device.		UDI code. UDI coding may not be shown on the label.
	Read the operating instruction carefully.		Contains type B applied parts.
	CE marking.		Dispose of the product in compliance with the procedures provided for electric and electronic equipment.

2.5. TECHNICAL SPECIFICATIONS.

2.5.1. Sector-specific standard compliance.

The patient chairs, for dental or generic use, model ECO.next and ECO.LUX comply with the Regulation (EU) 2017/745 (MDR), concerning medical devices.

Specifically, design and manufacture of the equipment comply with the following international standards:

EN 60601-1:2006 / A11 A1

EN 60601-1:2015

EN 80601-2-60:2015

2.5.2. Technical data.

ELECTRICAL SPECS.

Voltage	V	230 ± 10%	110 ± 10%
Frequency	Hz	50 – 60	
Input power	VA	1000	1000
Main fuses		2 x T 6,3A; 250V	2 x T 10A; 110V
Operating mode		Intermittent loading 25s ON 300s OFF	Intermittent loading 20s ON 300s OFF
Insulation class		I	
Applied part type		B	
Controls voltage	V	5	
Medical device risk class		I	

PERFORMANCES.

Chair type		Without rotation	With rotation
Minimum space required for installation	m	3 x 2	3 x 3
Max. height	mm	780	795
Min. height	mm	380	395
Rotation angle	deg	-	±30°
Trendelenburg (if present)	deg	15° (±2°)	
Width	mm	660	
Length (with completely elongated headrest)	mm	2050	
Packing dimensions	mm	1430x700x1105 (h)	
Gross weight	Kg	170	
Net weight	Kg	140	
Maximum lifting capability (patient + unit)	Kg	135 + 75	
External controls interface		Yes	




FIG. 2 Chair without rotation

Chair with rotation

2.5.3. Electromagnetic compatibility.

ECO.next and ECO.LUX chairs do not cause influences of an electromagnetic nature or other type with other equipment, nor are them susceptible to influences from other devices, and are compliant with the requirements of the harmonized standard EN 60601-1-2.

The following tables provide details on the electromagnetic environment in which the product can be used.

Emission aspects.

Emission test	Conformity	Electromagnetic environment – guide
RF Emissions Cispr 11	Group 1	The product uses RF energy for its internal function only. Therefore, its RF emissions are extremely low and not likely to cause interference with nearby electronic devices.
RF Emissions Cispr 11	Class B	It is possible to use the equipment in all buildings, including domestic use buildings, and those directly connected to a public power supply network with a low voltage supplying buildings for domestic use, taking installation steps, and guaranteeing an increased installation distance for potentially sensitive devices.
Harmonic emissions IEC 61000-3-2	Class A Conformity	It is possible to use the equipment in all buildings, including domestic buildings, and those directly connected to a public power supply network with a low voltage supplying buildings for domestic use.
Emission fluctuations voltage/flicker IEC61000-3-3	Conformity	

Immunity aspects.

ECO.next and ECO.LUX chairs are intended for use in an electromagnetic environment under specification. The customer or user should ensure that it is used correctly in such an environment.

Immunity test	Test level EN 60601-1-2	Conformity level	Electromagnetic environment – guide
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Burst / Fast transient EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (>95% dip UT) for 0,5 cycles 40% UT (60% dip UT) for 5 cycles 70% UT (30% dip UT) for 25 cycles < 5% UT (>95% dip UT) for 5 seconds	< 5% UT (>95% dip UT) for 0,5 cycles 40% UT (60% dip UT) for 5 cycles 70% UT (30% dip UT) for 25 cycles < 5% UT (>95% dip UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be that of a typical commercial or hospital environment.

Immunity aspects at RF.

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an electromagnetic environment.

Immunity test	Test level EN 60601-1-2	Compliance level	Electromagnetic environment – guide
RF conducted EN 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff from 150kHz to 80MHz	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: d = 1,2 · √P from 150kHz to 80MHz d = 1,2 · √P from 80 MHz to 800 MHz d = 2,3 · √P from 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
RF radiated EN 61000-4-3	3 Veff from 80MHz to 2,5GHz	3 Veff from 80MHz to 2,5GHz	
Field strengths from fixed RF transmitters; as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:			

Recommended separation distances between portable and mobile RF communication equipment and the device.

ECO.next and ECO.LUX chairs are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150kHz to 80MHz d = 1,2 · √P	From 80MHz to 800MHz d = 1,2 · √P	From 800MHz to 2GHz d = 2,3 · √P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

Notes:
 (1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
 (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

3. INSTALLATION

3.1. UNPACKING.

To remove the cardboard box, it is necessary to cut the straps and to unscrew the screws located in the lower part of the box.




FIG. 3

3.2. SUPPLY STATUS.

The following elements must be contained inside the packaging:

- An ECO.next or ECO.LUX chair;
- An operating instruction manual;
- A bottle of fabric cleaner;
- Two transportation handles.

In case of mistake or missing parts, contact the dealer.

3.3. HANDLING.

To remove the chair from the wooden pallet, unscrew the two screws that fix the chair base to the pallet (see figure 3).

To move the chair from the wooden pallet, employ two people.

To move the chair, hold it by the handles.

After positioning the chair, remove the handles as follows (see figure 4):

- Unscrew the hexagonal screws 1 and 2;
- Remove the "U" shaped handle, near the lower body;
- Unscrew the hexagonal screws 3 and 4;
- Remove the "L" shaped handle, near the upper body.




FIG. 4

3.4. INSTALLATION.

3.4.1. Typical installation.

Proceed as follows.

1. Be sure that the installation site conforms to local regulatory technical requirements for medical studios or outpatient clinics.
2. Be sure that the electrical system of the cabinet has a safe earthing and a residual-current device system.
3. Be sure that the main supply of the electrical wiring of the cabinet corresponds to the one of the identification label.
4. Check that the fuses, located at the front side of the chair near to the base, are correctly fixed, because during transport they could have become loose.
5. Insert the plug into the socket.
6. Activate the main switch to turn on the chair. Now the chair is ready for use. In case of non-operation of the chair, refer to section 6 about troubleshooting.
7. If needed, drive the chair to a comfortable height and assemble, if provided, the unit bracket and the unit itself, following the instructions of the manufacturer.

NOTE.




Equipment must not be used in potentially explosive environment, or in rooms not protected against weather conditions. Not suitable for use in altered atmospheres.

3.4.2. Position choice and stability.

To ensure a good chair installation, please verify the stability and flatness of the floor in the location where the product will be installed: if necessary, a damaged or unstable floor should be repaired.

1.1.2. The chair is designed to remain stable without floor fixing. However, according to the unit characteristics (dimensions, space required, weight) it could be necessary to fix the chair to the floor, using the appropriate holes in the base (figure 5, locations A and B).

For further information, ask the installer.



WARNING.



Tecnodent declines all and any responsibility for damage caused by the non-compliance with the above instructions.

4. USING THE PRODUCT

The advanced technology applied to the production of this apparatus has allowed the design of a patient chair with multiple functions and performances, all of them easily usable: all regulations and functions of the chair are coordinated, processed and stored in a "data bank", by a microprocessor inside the chair. The patient chairs ECO.next and ECO.LUX are designed to be easily controlled by foot controls, allowing the professional's hands to be completely free to operate in the patient's mouth.

WARNING.



Before using the information contained in this section, make sure that the chair is correctly installed, by the service technician.

CAUTION.



Use in a clean place (specialized surgery, clinic) and avoid accidental spilling or falling of liquid substances on the product.

CAUTION.



In case of a product malfunction that is not described in this documentation, put equipment out of service (shut it down, unplug main supply, mark the equipment as "out of order" or similar) and call the technical service.

4.1. TYPICAL USE.

To use the product, the operator must be near the chair, typically on the back of the chair or sideways, while acting on the controls. Care should be taken throughout the operation of the product, as the presence of automatic movements could cause hazardous situations, even if moderate level ones.

The chair is provided with several emergency stop systems, e.g., in case of presence of obstacles during the movements, or if the operator checks any unsuitable situation.

WARNING.



Do not overload the leg rest, e.g., by sitting on the front part of the upholstery.

4.2. CONTROLS.

The controls are shown in figure 6.

NOTE.



Depending on the configuration chosen, the chair may not be provided with foot controls as described below, because the chair is also designed to be controlled by other devices (e.g., a dental unit).


Please refer to section 7 for more information on optional controls

CHAIR WITHOUT PROGRAMMES.

- P Main switch (Power on/off).
- L Foot control, placed on the base. It controls raising, lowering and backrest tilt manual movements of the chair and the automatic return to patient exit ("zero") position.
- S Control the automatic movement of rinse position, in order to allow the patient's mouth wash and the return to the last working position.

CHAIR WITH PROGRAMMES.

- P Main switch (Power on/off).
- L Double function foot control, placed on the base. It controls the raising, lowering and backrest tilt manual movements of the chair, as well as the 3 programmable working positions and the automatic reset to patient exit ("zero") position.
- S Control the automatic movement of rinse position, in order to allow the patient's mouth wash and the return to the last working position.
- M Stores the programmed working positions.



4.3. CHAIR OPERATIONS.

4.3.1. Chair without programmes.

Connect the plug in the socket and switch on the chair with the main switch. The chair is now ready to be used.

"L" FOOT CONTROL OPERATION

The rear base foot controls are provided with a joystick and two buttons (see figure 7).

To control the main movements, move the joystick with the foot in the allowed directions.

- While keeping the joystick downward, the chair lowers.
- While keeping the joystick upward, the chair lifts.
- While keeping the joystick to the right of the chair, the backrest moves to the vertical position.
- While keeping the joystick to the left of the chair, the backrest moves to the supine position.
- While tapping the joystick downward, the chair moves automatically to the patient exit position.

"S" BUTTON OPERATION

While pushing "S" the rinse position function is activated; the chair backrest moves up to vertical position to allow the mouthwash. While pushing "S" again the chair backrest moves back to its last position.

NOTE.



Any switch, if pushed during an automatic movement or during the automatic reset to "zero" position movement, works like a safety STOP.




FIG. 7

4.3.2. Chair with programmes.

Connect the plug in the socket and switch on the chair with the main switch. The chair will emit a double sound signal ("BEEP – BEEP"), which indicates that the chair is ready for use.

"L" FOOT CONTROL OPERATION

- While keeping the joystick downward, the chair lowers.
- While keeping the joystick upward, the chair lifts.
- While keeping the joystick to the right of the chair, the backrest moves to the vertical position.
- While keeping the joystick to the left of the chair, the backrest moves to the supine position.
- While tapping the joystick upward, left or right, the chair moves automatically to the programmed working positions.
- While tapping the joystick downward, the chair moves automatically to the patient exit position.

NOTE.



To distinguish between "keep pressed" and "tap", the chair controller needs a delay of approximately 0.5 seconds.

"S" BUTTON OPERATION

While pushing "S" the rinse position function is activated; the chair backrest moves up to vertical position to allow the mouthwash. While pushing "S" again the chair backrest moves back to its last position.

NOTE.



Any switch, if pushed during an automatic movement or during the automatic reset to "zero" position movement, works like a safety STOP.

PROGRAMMING PROCEDURE

This function is only available on chairs equipped with such feature.

1. Set the chair to the patient exit position by tapping the joystick downward.
2. Configure the chair to the desired position by activating the corresponding foot controls (see paragraph 4.3).
3. While keeping the M button pressed, tap the joystick upward ("P1"), right ("P2") or left ("P3"), according to the selected program number.

An acoustic signal ("BEEP") indicates that the work position has been memorised.

Every time that the joystick is tapped upward, right or left, the chair will automatically configure to the selected pre-programmed working position.

NOTE.



To avoid losing the programmed position, it is recommended to regularly set the chair to its "zero" position, by using the "0" function at any patient change.

5. MAINTENANCE

5.1. GENERAL INFORMATION.

ECO.next and ECO.LUX chairs were conceived to not require any servicing during their working life. Consequently, it is not required to perform any adjustment of the equipment. If the chair will not be used for some time, it is recommended to cover it with a cloth, to turn off the main switch, and to disconnect the plug from the socket.

If you use the ECO.next and ECO.LUX chairs according to the technical specifications and operating instructions in this manual, their useful life has been calculated to be 10 years, based on laboratory simulations.

After that time, it is necessary to have the chair serviced by qualified personnel, especially the electronic parts and those subject to wear.

WARNING.



Any technical adjustment, as well as any repair of the equipment not specifically listed in this chapter, must be carried out by qualified technicians.

5.2. CLEANING.

For a better and longer life of the equipment, it is necessary to perform an accurate, methodical and periodical general cleaning of the chair. It is recommended to proceed as follows.

CAUTION.



The use of disinfectants or detergents with a strong alcoholic content on the apparatus damages the surface of the product, especially the upholsteries.

It is recommended to avoid the use of any dissolvent, strong detergent or abrasive agent, even to remove "difficult" stains.

5.2.1. Upholsteries.

The fabric that covers the chair upholstery has to be cleaned using the liquid contained in the bottle kit supplied with the chair, using a soft cloth to avoid any surface scratching and to guarantee a better elasticity and smoothness of the surface.

NOTE.



Once the liquid detergent has been consumed, please contact the nearest assistance centre for new product.

5.2.2. Plastic or metal parts.

The polyurethane parts, as well as the metal parts (painted and not) have to be cleaned with a soft cloth dipped in water and soap.

5.3. MAIN FUSES REPLACEMENT.

ECO.next and ECO.LUX chairs are protected against electrical power overrating by two fuses, placed as shown in figure 8. To replace them, proceed as follows:

1. Turn off the chair's main switch;
2. Disconnect the plug from the socket;
3. Unscrew the fuse holder cap, by using a crosshead screwdriver of medium size;
4. Replace the damaged fuses;
5. To re-install the fuses, follow the inverse sequence of the above steps.




FIG. 8

FUSES

Tecnodent opera una politica di continuo sviluppo dei prodotti.
In questa ottica, Tecnodent si riserva il diritto di variare in qualsiasi momento e senza
previa comunicazione le caratteristiche tecniche ed estetiche dei prodotti, qui riportate,
che non risultano pertanto vincolanti o impegnative.

Tecnodent pursues a policy of continuous product development.
Therefore, Tecnodent reserves all rights to change the technical and
aesthetical characteristics of the product at any time.

Tecnodent suit une philosophie de recherche continue sur les produits.
Tecnodent se réserve donc le droit de changer à tout moment et sans préavis
les caractéristiques techniques et esthétiques des produits présents dans ce manuel
qui n'engagent cependant en rien ou ne sont obligatoires.

Tecnodent sigue una política de constante desarrollo de sus productos.
Tecnodent se reserva el derecho de modificar en cualquier momento y sin
necesidad de aviso previo las características técnicas y estéticas de los productos
que aparecen en este manual por lo que no deben considerarse vinculantes.



UNI EN ISO 9001:2015
UNI EN ISO 13485:2016

TECNO⁺DENT
ERGONOMIC INNOVATION


Tecnodent S.r.l.

Via 63° Brigata Bolero – 24 – 40033 – Casalecchio di Reno (BO) – Italy

Tel. +39 051 613 11 43 – Fax +39 051 575 402 – P.IVA 02807791203

<http://www.tecnodent.com> – e-mail: info@tecnodent.com

Pat L. Costas



ECO.next chair is the evolution of a great classic that mixes high performances to the absolute pleasure of refined and timeless Italian design.

ECO.next has been designed to offer maximum anatomical comfort to the patient, complete ergonomics and perfect hygiene to the medical team.

As a result of the combination of different materials and technologies, the **ECO.next** chair matches comfort and high performances, giving a unique experience to the professional and to the patients.

The chair's structure is almost completely covered in integral polyurethane to assure the maximum strength, the compact base improves the functionality and the two lateral wings make the chair perfectly stable, the joystick pedal controls directly in a simple and efficient way the movements and allows to reach the programmed positions.

Control	Function	Unit connection	Programs	Armrest	Backrest	Pillows	Headrests
• Standard	• Small integrated utility box	• Knee-break legrest	• Knee program	• Return to "zero" program	• 17° shaped right armrest	• E/LE/FSSO	• E/LE/FSSO
• 2009new	• Large integrated utility box	• Security	• Outside side drilled upper body	• Left armrest	• 17° shaped left armrest	• E/LE/EL2 MOBILE	• E/LE/EL2 MOBILE
• ECO19	• Joystick	• Handshaking	• Fixed central arm	• Rotation	• 17° shaped right armrest	• UNI and Light	• UNI and Light
• ECO.next	• Foot control on the base	• Rotation with central arm	• Fixed central arm	• Handshaking	• 17° shaped left armrest	• C2002	• C2002
• Linda3	• Mobile foot control	• Handshaking	• Fixed central arm	• Handshaking	• 17° shaped right armrest	• Narrow	• Narrow
• Sting	• Infrared remote control	• Handshaking	• Fixed central arm	• Handshaking	• 17° shaped left armrest	• Large	• Large

HEADRESTS

Tecnodent manual Headrests are available in 4 different models: E/LE/FSSO, E/LE/EL2 MOBILE, UNI and Light. Every headrest was made of ensuring every operator's needs.



E/LE/EL2

E/LE/EL2, E/LE/EL2 FSSO has only the movement along its vertical axis. **E/LE/EL2 MOBILE** is an anatomical headrest, characterized by movement of inclination of the upper part (double articulation) for a correct cervical support.



C95



C2002

Magnetic pillows. E/LE/EL2 type, both fix and mobile version, can be combined with magnetic pillow C95 or C2002. Either pillow can be applied directly on the backrest, for ensure a better comfort, independency of the patient's height.



UNI




Light

UNI and Light. The UNI and Light headrests are multi-adjusted and allows positioning of the headrest along 3 axes.

TECNODENT
ERGONOMIC INNOVATION


Tecnodent S.p.A. Via 637 Brigata Bolero, 24 - 40073 Castelcchio di Reno (BO), Italy
Tel. +39 051 6131143 - Fax +39 051 573402 - info@tecnodent.com - www.tecnodent.com



ECO.next

TECNODENT
ERGONOMIC INNOVATION

— / —



ECO.net

- Electromechanical lifting, lowering and tilting movements;
- Compressed backrest tilting, up to the horizontal position;
- Narrow and slim anatomical backrest for the maximum liberty of action of the professional;
- Multi-articulated, anatomical, headrest, ELE2 MOBILE for a patient's comfortable posture and for a perfect visibility during the treatment;
- Wrapping seamless upholstery in high quality leatherette fabric;
- Possible customization with different colors of polyurethane and leatherette fabric;
- Integrated utility box;
- Complete symmetry for right-handed or left-handed professionals;

- Joystick control on the base for fast and efficient movements of the chair without using the hands, assuring a perfect hygiene;
- Security STOPS on the backrest, on the base and on the pantograph arm;
- Return to "zero" position automatic program;
- Rise position automatic program;
- Possibility to program 3 different working position (OPTION);
- Possibility to equip the chair with one/two rotating armrests to improve the accessibility of the chair and to assure a natural support for the arms, allowing a relaxed position during the treatment (OPTION).

Narrow backrest

This backrest was studied to allow the doctor the closest approach to the patient during the treatment, also guaranteeing an excellent comfort to the patient.

Tendelhubung Movement

At 15° with automatic compression of the backrest for a correct blood flow during the treatment.

Mobile foot control "Easy Touch"

Easy and functional foot control which provides different precise positions for the dentist and doctor's handle foot control.

Infrared remote control

Natural control of the movements and set to each programmed position.

Central arm (OPTION)

Unit fixing and rotating arm for ambidextrous functionality.



Large anatomical backrest

(OPTION)
Guarantees the maximum comfort, wrapping the patient's back in a pleasant relax.

Rotation (OPTION)

Rotating device which allows the chair to rotate around the vertical axis by a lever under the seat.

Right and/or left armrest

(OPTION)
Rotating armrests help the patient to get out of the chair.

Upright attachment

(OPTION)
Dedicated body for fastening dental unit lamp.

Polyurethane colors

Possibility to choose between 5 different colors.

TECHNICAL FEATURES

Maximum height	mm	Without rotation	780	With rotation	795
Minimum height	mm	Without rotation	380	With rotation	395
Backrest tilting angle	°		-15/+45		
Rotation angle	°		±30		
Width	mm		660		
Length (without headrest)	mm		1570		
Length (with headrest)	mm		2050		
Cross width	mm		170		
Net weight	kg		141		
Maximum lifting capability (patient + unit)	kg		135+75		
Classification					B

DICHIARAZIONE CE DI CONFORMITÀ

EU declaration of conformity

Il fabbricante

The manufacturer

TECNODENT S.r.l.

Via 63a Brigata Bolero, 24

40033 Casalecchio di Reno (BO) – ITALY

SRN: IT-MF-000027920

dichiara sotto la propria totale responsabilità che il seguente prodotto

declares under its own exclusive responsibility that the following product

Modello *Product Name*

ECO.NEXT

Poltrona dentale. Poltrona ad uso generico per applicazioni correlate al semplice posizionamento del paziente

Destinazione d'uso

Intended use

Dental patient chair; General purpose chair for simple patient positioning - related applications

UDI di Base (*basic UDI*): **8055728120000PDENTALN4**

È conforme ai requisiti essenziali di sicurezza come da allegati I ed alle prescrizioni di MDR 2017/745

Complies to the essential safety requirements and prescriptions of Medical Devices Regulation 2017/745



Classe del dispositivo medico *Medical device class*

I (ref. MDR2017/745, All. VIII regola 13)
(ref. MDR2017/745 Annex VIII rule 13)

Allegato II + Allegato III

Procedura di certificazione *Certification procedure*

Annex II & Annex III

L'azienda dichiara inoltre che i prodotti indicati nel presente dichiarazione di conformità sono:

- Conformi alla Direttiva 2011/65/UE (Restriction of Hazardous Substances), sulla restrizione dell'uso di sostanze pericolose nelle apparecchiature elettriche ed elettroniche, comprese le disposizioni della Direttiva Delegata (UE) 2015/863 emanata dalla Commissione Europea.

TECNODENT moreover hereby states, that the products mentioned in this Declaration of Conformity are compliant with Directive 2011/65/EU (Restriction of Hazardous Substances) concerning the restrictions on the use of hazardous substances in electrical and electronic equipment, including the provisions of Commission Delegated Directive (EU) 2015/863 issued by the European Parliament and the Council.

Casalecchio di Reno (BO)

14/09/2022

Il Consigliere Delegato

Massimiliano Tabacchi
(Member of the Board)

TO WHOM IT MAY CONCERN

We, company TECNODENT Srl, located in Casalecchio di Reno (Bologna) Italy, Via 63a Brigata Bolero 24, as manufacturer of dental/ medical chairs and stools we hereby authorize our partner LIMETA Jsc. Graiciuno str. 4, LT-02241 Vilnius, Lithuania (Company code: 221906050 - VAT code: LT219060515) to supply and install our production range, provide after sale, warranty services, promote and train his employees on Tecnodent products in its territory of competence.

Best regards.

Casalecchio di Reno, March 21st 2022


TECNODENT Srl
Elisa Guidubaldi
General Sales Manager

TECNODENT

ERGONOMIC INNOVATION

TO WHOM IT MAY CONCERN

We, company TECNODENT Srl, located in Casalecchio di Reno (Bologna) Italy, Via 63° Brigata Bolero 24, as manufacturer of dental / medical chairs and stools we hereby state that:

ECO NEXT patient chair footrest section is equipped with a transparent protective plastic cover that can be easily removed for cleaning purposes.

Best regards.

Casalecchio di Reno, January 9th 2023



TECNODENT SRL
Elisa Guidubaldi
General Sales Manager



TECNODENT S.r.l. Società Unipersonale

Via 63a Brigata Bolero, 24
40033 CASALECCHIO DI RENO - BOLOGNA (ITALIA)
Tel. 39 - 051 61 31 143 - Telefax 39 - 051 57 54 02
Internet: www.tecnodent.com - E-mail: info@tecnodent.com

Capitale Sociale € 100.000 int. vers.
P. I.V.A. IT02807791203
Cod. Fisc. e Reg. Imp. di Bo N. 02807791203



STOMADENT®

CERTIFICATE

for authorized service technician:

A. Blokhin

According to attended training and successfully carried out examinations in STOMADENT training center we hereby certify the authorization to carry out service operations of under-mentioned products:

Dental Units:

**STOMADENT IMPULS
STOMADENT GLANC
STOMADENT NEO**

Dental Chair:

**STOMADENT DC-70
STOMADENT NEO**

For this purpose we are pleased to hand over this certificate



Valid to January 18st, 2026

STOMADENT SK s.r.o.
HUMENNÉ, January 18st, 2021



STOMADENT®

CERTIFICATE

for authorized service technician:

Vidas Sileikis

According to attended training and successfully carried out examinations in STOMADENT training center we hereby certify the authorization to carry out service operations of under-mentioned products:

Dental Units:

**STOMADENT IMPULS
STOMADENT GLANC
STOMADENT NEO**

Dental Chair:

**STOMADENT DC-70
STOMADENT NEO**

For this purpose we are pleased to hand over this certificate



Valid to January 18st, 2026

STOMADENT SK s.r.o.
HUMENNÉ, January 18st, 2021



3EC International a.s., Hraničná 18, 821 05 Bratislava, Slovakia
Notified Body No. 2265

EC CERTIFICATE

No. 2019-MDD/QS-054

issued in compliance with the Council Directive 93/42/EEC as amended by 2007/47/EC,
which is implemented by the Slovak Government Decree No. 582/2008 Coll. as amended by 215/2013
Coll. certifies that the medical device of Class IIa.

Dental Unit
STOMADENT IMPULS
STOMADENT NEO
STOMADENT HARMONY

manufactured by company

STOMADENT SK s.r.o.
300, Lackovce 066 01, Slovak Republic

is manufactured under conditions fulfilling the quality system requirements of Annex VI, of the
Directive 93/42/EEC as amended by 2007/47/EC.

The Notified Body No. 2265 has performed an audit of the above device quality system. The product
quality assurance has been assessed and found that it meets the requirements above. The quality system
is subject to continuous surveillance according to Annex VI, Sections 3.3. and 4, of the Directive
93/42/EEC as amended 2007/47/EC. The detailed description of the system, requirements and measures
applied by the manufacturer are presented in the Audit Report No. 310431, and the Final protocol No.
310431/2019.

This certificate is issued under the following conditions:

It applies only to the quality system maintained in the manufacture of the above referenced model of
medical device and it does not substitute the design or type-examination procedures, if requested. The
certificate remains valid until the manufacturing conditions or the quality system are changed but until May
26th, 2024 at the latest. The certificate validity is conditional upon positive results of regular surveillance
audits and fulfilment of relevant legal and other requirements by manufacturer.



Dr. Katarína Tomin Srdošová
Responsible to act on behalf of NB 2265

In Bratislava, on September 15th, 2019

/Logotipas: 3EC International/

3EC International a.s., Hranična 18, 821 05 Bratislava, Slovakijos Respublika
Notifikuotosios įstaigos Nr. 2265

EC SERTIFIKATAS

Nr. 2019-MDD/QS-054

išleistas pagal Tarybos direktyvą 93/42/EEB, kaip atnaujinta su 2007/47/EB, kuri yra pritaikyta Slovakijos Vyriausybės dekretu Nr. 582/2008 Coll, kaip atnaujinta su 215/2013 Coll., patvirtina, kad medicinos prietaisai yra IIa klasės

Odontologinis prietaisas
STOMADENT IMPULS
STOMADENT NEO
STOMADENT HARMONY

pagamino įmonė

STOMADENT SK s.r.o.
Lackovce 188, 066 01 Humenne, Slovakijos Respublika

gamina pagal sąlygas, atitinkančias kokybės sistemos reikalavimus, pateiktus Direktyvos 93/42/EEB, kaip atnaujinta su 2007/47/EB VI Priede.

Notifikuotoji įstaiga Nr. 2265 atliko aukščiau paminėto prietaiso kokybės sistemos auditą. Buvo įvertinta pilna kokybės užtikrinimo sistema ir rasta, kad atitinka paminėtus reikalavimus. Kokybės sistema gali būti nuolat tikrinama pagal Direktyvą 93/42/EEB, kaip atnaujinta su 2007/47/EB. Detalus sistemos, reikalavimų ir būtų, kuriuos taiko gamintojas, aprašymas yra pateiktas Audito ataskaitoje Nr. 310431 ir Galutiniame protokole Nr. 310431/2019.

Šis sertifikatas išleistas pagal šias sąlygas:

Taikomas tik kokybės sistemai, kurią taiko gamintojas aukščiau paminėtiems medicininių prietaisų modeliams ir nepakeičia dizaino ar tipo tyrimo procedūrų, jei reikalaujama. Sertifikatas galioja kol nepakeičiamos kokybės sistemos gaminimo sąlygos, bet ilgiausiai iki 2024 m. gegužės 26 d. Sertifikato galiojimas priklauso nuo reguliarių priežiūros auditų teigiamų rezultatų ir atitikimo atitinkamiems įstatymams ir kitiems gamintojo reikalavimams.

/Antspaudas/

/Parašas/

Dr. Katerina Tomin Srdošova
Atsakinga už dokumentą ir NB 2265 vardu

Bratislava, 2019 m. rugsėjo 15 d.

Vertėjas Siantas Bantashevicius

1. ISTRUZIONI D'USO PER L'ODONTOIATRIA/PARAMEDICO

- **1.1 GENERALITÀ**
La "SIRINGA MINIMATE" è un dispositivo ad uso esclusivamente odontoiatrico, nato con lo scopo di insuflare aria ed acqua (da sole o combinate) per mantenere costantemente pulito ed asciutto il campo operatorio.
- **1.2 CARATTERISTICHE GENERALI**
La SIRINGA MINIMATE è stata progettata con i più moderni concetti ergonomici per un facile uso ed una immediata pulizia e sterilizzazione. Sia il puntale che l'impugnatura sono facilmente estraibili per una perfetta disinfezione e sterilizzazione in autoclave a 135°C. Sono disponibili pertanto più impugnature di forma diversa: il dentista può scegliere la forma che preferisce in base alla necessità: a gomito o a stilo.
- **1.3 MARCATURA CE**
Ogni prodotto è marcato CE.

2. DATI IDENTIFICATIVI E GARANZIA

- **2.1 LOTTO DI COSTRUZIONE**
Ogni prodotto è caratterizzato da un numero riportato all'interno della parte terminale, numero che identifica univocamente il lotto di produzione. Mediante questo numero è sempre possibile rintracciare il periodo di costruzione e le relative schede di controllo.
- **2.2 GARANZIA**
Il prodotto è garantito dalla nostra società per 12 mesi dalla data del documento di consegna. Qualsiasi manipolazione o modifica non autorizzata, oltre a far decadere immediatamente la garanzia, solleva la nostra società da ogni responsabilità per danni a persone, animali o cose che potrebbero verificarsi in conseguenza di tale manutenzione. Per ogni controversia il Foro competente è quello di Milano (Italy).

3. CARATTERISTICHE TECNICHE

Pressione acqua di alimentazione max	Bar	2,5
Pressione di alimentazione aria max	Bar	4,5
Portata aria	l/min	10
Portata acqua	Cc/min	110

SCHEDA SEGNALAZIONE ANOMALIE SIRINGA MINIMATE

Prodotto _____
 Tipo _____
 Tipo di segnalazione: Anomalia Suggerimento

Descrizione

Spedire a: LUZZANI DENTAL SRL - VIA TORINO, 3
 SENAGO (MI) ITALY
 oppure spedire al seguente fax: +39 0299010379
 e-mail: info@luzzani.com

4. INSTALLAZIONE E COLLEGAMENTI

- **4.1** Data la natura dello strumento l'installazione e riparazione deve essere eseguita SOLO dal fabbricante di riuniti o da un suo tecnico autorizzato (il fabbricante di riuniti è dotato di un apposito manuale di installazione).

5. USO NORMALE

- **5.1** Per insuflare acqua fredda nel campo operatorio basta premere il pulsante sinistro dell'impugnatura.
- **5.2** Per insuflare aria fredda nel campo operatorio basta premere il pulsante destro dell'impugnatura.
- **5.3** Per insuflare acqua ed aria fredde in modo combinato (spray), basta premere contemporaneamente i due pulsanti dell'impugnatura.

6. PULIZIA E STERILIZZAZIONE

Dopo ogni intervento su un paziente si deve, al fine di garantire la massima igiene, pulire e sterilizzare la siringa. A tal fine occorre eseguire le fasi seguenti:

- staccare il puntale (svitando il terminale ferma punta) e/o l'impugnatura completa (premiendo il pulsante posto nella parte inferiore della stessa e tirando verso l'alto).
- pulirli con un panno pulito, rimuovendo eventuali macchie.
- porli in autoclave a vapore d'acqua a 135°C 4 minuti minimo.

7. MANUTENZIONE

Non è prevista né necessaria alcuna manutenzione specifica dell'apparecchio, se non la normale pulizia e sterilizzazione descritte al punto precedente. Ogni operazione di lubrificazione è da escludere in quanto può causare danni irreparabili alla siringa.

8. SMALTIMENTO E ROTTAMAZIONE

Il prodotto non contiene componenti pericolosi o tossicologici, né viene in contatto con prodotti di questo tipo durante il ciclo di lavoro.

9. SEGNALAZIONE ANOMALIE

Per il soddisfacimento dei requisiti di conformità alla dir. 93/42 CEE, l'azienda ha istituito una procedura di sorveglianza post-vendita su eventuali problemi generati dall'utilizzo dei nostri prodotti. Vi chiediamo gentilmente di segnalarci eventuali anomalie inviandoci il tagliando sottostante.

SCHEDA SEGNALAZIONE ANOMALIE SIRINGA MINIMATE

Data _____
 Società _____
 Indirizzo _____

Tipo di segnalazione: Anomalia Suggerimento

Descrizione

Spedire a: LUZZANI DENTAL SRL - VIA TORINO, 3
 SENAGO (MI) ITALY
 oppure spedire al seguente fax: +39 0299010379
 e-mail: info@luzzani.com

1. DIRECTIONS FOR USE BY DENTIST/PARAMEDIC

- **1.1 INTRODUCTION**
The "MINIMATE SYRINGE" is a device for dental use only, designed to insufflate air and water (separately or combined) to keep the operating field constantly clean and dry.

2. GENERAL CHARACTERISTICS

- **2.1 IDENTIFICATION DATA AND GUARANTEE**
A number marked on the inside of the terminal section of each product identifies the production batch. This number can be used to trace the period of manufacture and the associated control cards.
- **2.2 GUARANTEE**
The product is guaranteed by our company for 12 months from the date of the delivery document. Any unauthorised manipulation or modification will immediately terminate the guarantee and exonerate our company from all liability for injury or damage to persons, animals or chattels that may be caused by such interference. All disputes shall be subject to the jurisdiction of the courts of Milan (Italy).

3. TECHNICAL SPECIFICATIONS

Max. water input pressure	Bars	2,5
Max. air input pressure	Bars	4,5
Air flow rate	l/min	10
Water flow rate	Cc/min	110

4. INSTALLATION AND CONNECTIONS

- **4.1** In view of the nature of the instrument, installation and repairs may ONLY be performed by the dental unit manufacturer or its authorised technicians (the dental unit manufacturer holds the specific installation manual).

5. NORMAL USE

- **5.1** To insufflate cold water into the operating field, press the left-hand button on the handpiece.
- **5.2** To insufflate cold air into the operating field, press the right-hand button on the handpiece.
- **5.3** To insufflate a combination of cold air and water (spray), press the two buttons on the handpiece simultaneously.

6. CLEANING AND STERILISATION

After each use on a patient, the syringe must be cleaned and sterilised to guarantee the maximum hygiene. For this purpose, proceed as follows:

- detach the tip (by unscrewing the tip-retaining terminal) and/or the complete handpiece (by pressing the button in the lower part of the handpiece and pulling upwards).
- wipe with a clean cloth and remove any stains.
- place in a steam autoclave at 135°C for 4 minutes.

7. MAINTENANCE

No specific maintenance of the instrument is required apart from the normal cleaning and sterilisation described in the previous paragraph. Lubrication must not be performed, as it can cause irreparable damage to the syringe.

8. DISPOSAL AND SCRAPPING

The product does not contain any dangerous or toxic/hazardous components, nor does it come into contact with such products during the operating cycle.

9. MALFUNCTION REPORTS

To meet the requirements of Directive 93/42/ECC, the company has instituted a post-marketing surveillance procedure to monitor any problems generated by the use of our products. We should therefore be grateful if you would use the form below to inform us of any malfunctions.


MINIMATE SYRINGE MALFUNCTION REPORT

Date _____
 Signature _____
 Address _____

Type of report: Malfunction Suggestion


Description

Please send by post to: LUZZANI DENTAL SRL - VIA TORINO, 3
 SENAGO (MI) ITALY
 or fax to: +39 0299010379
 e-mail: info@luzzani.com




Shawmut 250
Pacheco's -
- 250 Street Co


Stamant Isibatis - Plakats - wotshar ka




FORMEX




CLINE K




CLINE



DENTAL



DENTAL PL-K



CLINE FL

ST 1123	PLU2
ST 5061	VI02
ST 1017	ORA4
ST 5141	ORA3
ST 3096	BLU5
ST 5154	STL1
ST 3104	LB01
ST 1210	GRN4
ST 5153	TYR1
ST 6156	
ST 6100	YEL3
ST 6021	LGV5
ST 6127	GRV1
NEFRYT	BLU6
ST 5161	GRN5
ST 7000	RED3

Please note that this colour chart is for illustration purposes only.

Catalogo

CLINICAL MICROMOTORS

83

CLINICAL
MICROMOTORS

82

NLXplus

Optic Internal Water Spray

Electric Micromotor Integrated System **BRUSHLESS**
 Speed : 100 – 40,000 min⁻¹, Torque : 4.0 Ncm, with Endo Function



Features

- Torque control & Auto reverse feature
- 3 working modes: Auto-reverse on, Auto-stop or Auto-reverse off
- Autoclavable titanium body
- LED illumination
- Installable built-in type conforms to most dental equipments by employing original NLX BF power module

Integrated System Individual Components



- MODEL NLX BF** **ORDER CODE Y1002122**
- Power module with lamp and air function for integrated system
 - Dimensions with case : W 115 x D 88.5 x H 40 (mm)



- MODEL NLS-A ADP** **ORDER CODE U1020050**
- 24 Vdc - 32 Vdc converter for integrated system
 - Dimensions : W 60.2 x D 81.4 x H 45.2 (mm)



- MODEL NLS-A SEL** **ORDER CODE Y1010403**
- Selector module for integrated system
 - Dimensions : W 68.5 x D 56.5 x H 38.5 (mm)



- MODEL NLX plus** **ORDER CODE E1040051**
- NLX plus micromotor only (Cord not included)
 - Solid Titanium Body
 - Autoclavable up to 135°C



- MODEL plus CD** **ORDER CODE E1040062**
- NLX plus micromotor cord for integrated system (2.2 m)



- MODEL Multi Pad** **ORDER CODE Y1001845**
MODEL Multi Pad Auto-Select **ORDER CODE Y1001846**
- Multi-Control Panel for built-in system
 - Dimensions : W 95 x D 138 x H 31 (mm)

NBX

Optic Internal Water Spray

Electric Micromotor Integrated System **BRUSH**
 Speed : 60 – 40,000 min⁻¹



Integrated System Individual Components



- MODEL iMD BS** **ORDER CODE Y1002845**
- Power module for integrated system
 - Dimensions : W 96 x D 57 x H 33 (mm)



- MODEL iMD SEL** **ORDER CODE Y1002848**
- Selector module for integrated system
 - Dimensions : W 88.5 x D 56.5 x H 38.5 (mm)



- MODEL NBX** **ORDER CODE E1059051**
- Micromotor only (Cord not included)
 - Solid Titanium Body



- MODEL NBX CDB** **ORDER CODE E1059061**
- NBX micromotor cord for integrated system



- MODEL Multi Pad** **ORDER CODE Y1001845**
MODEL Multi Pad Auto-Select **ORDER CODE Y1001846**
- Multi-Control Panel for built-in system
 - Dimensions : W 95 x D 138 x H 31 (mm)

EC Certificate
Directive 93/42/EEC Annex II, excluding Section 4
Full Quality Assurance System
Medical Devices

Registration No.: HD 60148508 0001

Report No.: 12031343 009

Manufacturer: Nakanishi Inc.
700 Shimohinata, Kanuma,
Tochigi,
322-8666 Japan

Products: Instruments and Equipment in Dental and Surgical Fields
(see attachments for products included)
Replaces Approval, Registration No.: HD 60144337 0001

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

Effective Date: 2020-10-03

Date: 2020-10-03



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

**Attachment to
Certificate**

Registration No.: HD 60148508 0001
Report No.: 12031343 009

Manufacturer: **Nakanishi Inc.**
700 Shimohinata, Kanuma,
Tochigi,
322-8666 Japan

Products included:

- Pneumatic dental drilling system handpiece
- Diamond dental bur, reusable
- Dental drilling system attachment
- Line-powered dental drilling system handpiece
- Pneumatic dental scaling system handpiece
- Pneumatic dental scaling system handpiece tip
- Air Motor
- Remotely-driven dental drilling system motor
- Endodontic enlarger
- Dental-professional prophylaxis motor
- Endodontic apex locator
- Dental abrasive air jet system handpiece
- Dental abrasive powder
- Sawing power tool attachment, reciprocating
- Sawing power tool attachment, oscillating



Notified Body

Date: 2020-10-03

M. Aihara
M.Sc. M. Aihara

TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

**Attachment to
Certificate**

Registration No.: HD 60148508 0001

Report No.: 12031343 009

Manufacturer: Nakanishi Inc.
700 Shimohinata, Kanuma,
Tochigi,
322-8666 Japan

Products included:

- Sawing power tool attachment, sagittal
- Reciprocating surgical saw blade, reusable
- Oscillating surgical saw blade, reusable
- Sagittal surgical saw blade, reusable
- Dental implant system
- Surgical screwdriver, reusable
- Ultrasonic dental scaling system
- Ultrasonic dental scaling system handpiece tip
- Dental ultrasonic surgical system handpiece tip
- Dental ultrasonic surgical system
- Total Surgical System
- Sterile Products:
 - Sterile Drills
 - Sterile Burs
 - Sterile Blades
 - Sterile Rasps
 - Sterile Files
- Dental suction system, non-surgical



Notified Body

M. Aihara

Date: 2020-10-03

M.Sc. M. Aihara

TÜV Rheinland
LGA Products GmbH
Tillystraße 2, 90431 Nürnberg

**Attachment to
Certificate**

Registration No.: HD 60148508 0001
Report No.: 12031343 009

Manufacturer: Nakanishi Inc.
700 Shimohinata, Kanuma,
Tochigi,
322-8666 Japan

Products included:

- Pneumatic Surgical Drilling Machine
- Electric drill system for dental surgery, dental implant surgery and foot- and ankle- surgery
- External Coolant Straight Handpiece for foot surgery and oral surgery

Site included:

A1 Factory
990 Fukahodo, Kanuma, Tochigi, 322-0302, Japan




Notified Body

M. Aihara

Date: 2020-10-03

M.Sc. M. Aihara



嘱託人株式会社ナカニシ宮崎友美子の提出した別添証書の謄本は、その原本と対照し、原本と符合することを認めた。

よって、これを認証する。

令和3年 2 月 26 日、本公証人役場において
東京都台東区東上野1丁目7番2号

東京法務局所属

公証人
Notary

畑野隆二



HATANO Ryuji


証 明

上記署名は、東京法務局所属公証人の署名に相違ないものであり、かつ、その押印は、真実のものであることを証明する。

令和3年 2 月 26 日

東京法務局長

山西宏紀



APOSTILLE

(Convention de La Haye du 5 octobre 1961)

1. Country: JAPAN
This public document
 2. has been signed by HATANO Ryuji
 3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
 4. bears the seal/stamp of HATANO Ryuji ,Notary
- Certified
5. at Tokyo
 6. FEB. 26. 2021
 7. by the Ministry of Foreign Affairs
 8. 21- No 012132
 9. Seal/stamp:
 10. Signature



Tanaka Toshie

TANAKA Toshie

For the Minister for Foreign Affairs

Registered No. 59 of 2021

NOTARIAL CERTIFICATE

This is to certify that the attached copy of the documents
produced before me by MIYAZAKI Yumiko , NAKANISHI
INC. exactly corresponds with the seen original.

畑野 隆二



Notary, HATANO Ryuji

Ueno Notary's Office of Tokyo Legal Affairs

Bureau , 1-7-2 Higashiueno , Taito-ku ,

Tokyo, Japan

BROSIVERA

Bien Air⁺
Dental



 BIEN-AIR
SWISS
MADE

CONTRA-ANGLES AND HANDPIECES
UNFAILING PRECISION



5



ACCU-SPRAY

Research conducted by Bien-Air engineers has clearly demonstrated the superior cooling power of six separate air and water nozzles. With the ACCU-SPRAY system, the air disperses the water to a highly precise point before contact with the bur. The ultra-precise nebulisation cools the bur and dental material as efficiently as possible. It also ensures perfect visibility in the working area. The separate water and air ducts also prevent blockages.



COOL TOUCH

Bien-Air is the only manufacturer in the world to offer a push-button bur change mechanism with an anti-heating safety system. This reduces the risk of burning your patients with the head of your instrument. A tungsten carbide ball is placed between the pushbutton and the top of the bur locking mechanism. It absorbs heat in the event of contact when rotating. COOL TOUCH. An exclusive invention for your comfort and your patients' safety.



TECHNICAL DATA

micro
SERIES



CA 1:5 L Micro-Series

Type	Contra-angle with internal spray, anti-retraction valve
Ratio	1:5 (speed increasing factor of 5)
Coupling	Compatible with Micro-Series
Locking system	Push-button bur locking with anti-heating system
Bur dimension	Bur shank Ø 1.6 mm
Light	Multi-strand optical glass conductor light
Cleaning/sterilisation	Thermal washer/steam autoclave up to 135°C (275°F)
Ref.	1600690-001

micro
SERIES



CA 1:1 L Micro-Series

Type	Contra-angle with internal spray, anti-retraction valve
Ratio	1:1 (direct)
Coupling	Compatible with Micro-Series
Locking system	Push-button bur locking
Bur dimension	Bur shank Ø 2.35 mm
Light	Multi-strand optical glass conductor light
Cleaning/sterilisation	Thermal washer/steam autoclave up to 135°C (275°F)
Ref.	1600691-001

micro
SERIES



CA 20:1 L Micro-Series / CA 20:1 L KM Micro-Series

Type	Contra-angle with internal spray
Ratio	20:1 (speed reducing factor of 20)
Coupling	Compatible with Micro-Series
Locking system	Push-button bur locking
Bur dimension	Bur shank Ø 2.35 mm
Light	Multi-strand optical glass conductor light
Cleaning/sterilisation	Thermal washer/steam autoclave up to 134°C (273°F)
CA 20:1 L Micro-Series ref.	1600692-001
CA 20:1 L KM Micro-Series ref.	1600786-001

micro
SERIES




PM 1:1 Micro-Series

Type	Straight handpiece with internal spray
Ratio	1:1 (direct)
Coupling	Compatible with Micro-Series
Locking system	Locking ring
Bur dimension	Bur shank Ø 2.35 mm
Cleaning/sterilisation	Thermal washer/steam autoclave up to 135°C (275°F)
Ref.	1600693-001

CA 1:5 L / CA 1:5

Type	Contra-angle with internal spray, anti-retraction valve
Ratio	1:5 (speed increasing factor of 5)
Coupling	E-type as per ISO 3964
Locking system	Push-button bur locking with anti-heating system
Bur dimension	Bur shank Ø 1.6 mm
Cleaning/sterilisation	Thermal washer/steam autoclave up to 135°C (275°F)
Ref. CA 1:5 L	1600386-001
Light	Multi-strand optical glass conductor light
Ref. CA 1:5	1600325-001



CA 1:1 L / CA 1:1

Type	Contra-angle with internal spray, anti-retraction valve
Ratio	1:1 (direct)
Coupling	E-type as per ISO 3964
Locking system	Push-button bur locking
Bur dimension	Bur shank Ø 2.35 mm
Cleaning/sterilisation	Thermal washer/steam autoclave up to 134°C (273°F)
Ref. CA 1:1 L	1600384-001
Light	Multi-strand optical glass conductor light
Ref. CA 1:1	1600424-001



CA 10:1 L / CA 10:1

Type	Contra-angle with internal spray
Ratio	10:1 (speed reducing factor of 10)
Coupling	E-type as per ISO 3964
Locking system	Push-button bur locking
Bur dimension	Bur shank Ø 2.35 mm
Cleaning/sterilisation	Thermal washer/steam autoclave up to 134°C (273°F)
Ref. CA 10:1 L	1600385-001
Light	Multi-strand optical glass conductor light
Ref. CA 10:1	1600425-001



CAP 15:1 Prophy

Type	Contra-angle without light or spray, screw or clip cup interchangeability
Ratio	15:1 (speed reducing factor of 15)
Cleaning/sterilisation	Thermal washer/steam autoclave up to 134°C (273°F)
Ref.	1600290-001



PMP 10:1 Prophy

Type	Handpiece without light or spray, for standard disposable heads
Ratio	10:1 (speed reducing factor of 10)
Cleaning/sterilisation	Thermal washer/steam autoclave up to 134°C (273°F)
Ref.	1600289-001



Home / Products / Contra-Angles & Straight Handpieces / Restorative Classic / CA 1:1

CA 1:1 HERITAGE OF RELIABILITY.



There are some devices still earning respect and admiration forty years in, and Bien-Air's classic CA 1:1 and CA 1:5 electric handpieces are two of them. Forever evolving, their dominating productivity and remarkable durability make them the most popular models in the lineup. The series' fully redesigned bur-locking and rotation drive mechanism ensures micron-precise concentricity and stability, while a 60 to 200,000 rpm speed range lets you complete the most intricate restorative work with incredible accuracy – and on your own terms. CA 1:1 and CA 1:5 epitomize the reliability and quality Bien-Air has built its reputation on.

Reference:

CA 1:1 L Standard 1600384-001

CA 1:1 L Micro-Series 1600691-001

CA 1:1 Standard 1600424-001



TECHNICAL DATA

	CA 1:1 L Standard	CA 1:1 L Micro-Series	CA 1:1 Standard
Transmission ratio	1:1	1:1	1:1
Maximum speed (rpm)	40'000	40'000	40'000
Head	-	-	-
Head diameter (mm)	9.8	9.8	9.8
Head size (mm)	15.3	15.3	15.3
Weight (g)	86	79	88
Noise level (dBa)	57	57	57
Spray system	2x3 separate air/water sprays	2x3 separate air/water sprays	2x3 separate air/water sprays
Lighting system	Multi-strand optical glass technology	Multi-strand optical glass technology	-
Anti-heating system	-	-	-
Coupling system	ISO 3964 & Bien-Air Micro-Series	ISO 3964 & Bien-Air Micro-Series	ISO 3964 & Bien-Air Micro-Series
Warranty	2 to 3 years ¹	2 to 3 years ¹	2 to 3 years ¹
Reference	CA 1:1 L Standard 1600384-001	CA 1:1 L Micro-Series 1600691-001	CA 1:1 Standard 1600424-001

CA 1:1 L Standard

Transmission ratio	1:1
Maximum speed (rpm)	40'000
Head	-
Head diameter (mm)	9.8
Head size (mm)	15.3
Weight (g)	86
Noise level (dBa)	57
Spray system	2x3 separate air/water sprays
Lighting system	Multi-strand optical glass technology
Anti-heating system	-
Coupling system	ISO 3964 & Bien-Air Micro-Series
Warranty	2 to 3 years ¹
Reference	CA 1:1 L Standard 1600384-001

CA 1:1 L Micro-Series

Transmission ratio	1:1
Maximum speed (rpm)	40'000
Head	-

CA 1:1 L Micro-Series

Head diameter (mm)	9.8
Head size (mm)	15.3
Weight (g)	79
Noise level (dBa)	57
Spray system	2x3 separate air/water sprays
Lighting system	Multi-strand optical glass technology
Anti-heating system	-
Coupling system	ISO 3964 & Bien-Air Micro-Series
Warranty	2 to 3 years ¹
Reference	CA 1:1 L Micro-Series 1600691-001

CA 1:1 Standard

Transmission ratio	1:1
Maximum speed (rpm)	40'000
Head	-
Head diameter (mm)	9.8
Head size (mm)	15.3
Weight (g)	88
Noise level (dBa)	57
Spray system	2x3 separate air/water sprays
Lighting system	-
Anti-heating system	-
Coupling system	ISO 3964 & Bien-Air Micro-Series
Warranty	2 to 3 years ¹
Reference	CA 1:1 Standard 1600424-001

¹ 2-year standard warranty and 1-year optional warranty available through [Bien-Air's PlanCare extended warranty program](#).

Bien-Air Medical Technologies develops, manufactures, and markets instrumentation devices and systems across multiple dental and surgery specialities. The company operates eight subsidiaries from its Bienne, Switzerland headquarters and is backed by a solid network

Information

Where to buy?

Where to repair?

Maintenance of the products

Privacy

General Terms and Conditions

Certificates of Quality

Newsletter

Keep up on our always evolving product features and technology. Enter your e-mail and subscribe to our newsletter.

Email Address

Contact Us

Bien-Air Dental SA

Länggasse 60

Case Postale

CH-2500 Bienne / Suisse

Phone: +41(0) 32 344 64 64

Fax: +41(0) 32 344 64 91

Email: office@bienair.com

created by Ide

ORIGINAL



LETTER OF AUTHORIZATION

To whom it may concern

Bienne, February 25th 2016

We, Bien-Air Dental S.A. located in Länggasse 60, case Postale , CH-2500 Bienne 6, Switzerland hereby confirm that the company:

Medical Trade, s.r.o
Kreslicka 1
101 00 Prague 10
Czech Republic

is an authorized sales and service partner for all the Bien-Air Dental products in the Czech Republic.

With best regards.

Michel NOEL
Sales Manager

A handwritten signature in black ink, appearing to read "Michel Noel", written over a faint circular stamp.

Bien-Air Dental SA
Länggasse 60
Case postale
CH-2500 Bienne 6

A member of the Bien-Air Group

Prague, April 19th 2021

To whom it may concern

We, Medical Trade, s.r.o., located in Křeslická 301/1, 101 00 Prague 10, Czech Republic,
as an authorized **sale and service partner for all Bien-Air Dental S. A. products in
the Czech Republic,**

hereby confirm

that **CA 1:5 (L) and CA 1:1 (L) contra-angles**, produced by Bien-Air Dental S. A., are
equipped with **internal cooling system and 2 x 3 separate air/water sprays.**

With best regards

Ing. Vladimír Valenta
Executive Director



Medical Trade, s. r. o.
Křeslická 1
101 00 Prague 10
Czech Republic
Tel.: +420 602 350 661

 **MEDICAL** TRADE[®] s.r.o.
Křeslická 1, 101 00 Praha 10
Tel.: 272 765 043
IČO: 25790579, DIČ: CZ25790579



Prague, April 8th 2021

To whom it may concern

We, Medical Trade located in Křeslická 301/1, 101 00 Prague 10, Czech Republic,
who are **sale and service partner for all Bien-Air Dental S. A. products in the Czech Republic,**

hereby confirm that following **Bien-Air Dental S. A. turbines and handpieces are sterilizable in Class B autoclaves up to 135°C:**

All versions of turbines: **Bora, Prestige, Tornado**

All versions of contra-angles: **CA 1:5 (EVO), CA 1:1 (EVO), CA 10:1, CA 20:1**


All versions of straight handpieces: **PM 1:1**

With best regards

Bc. Martina Pokrývková Valentová
Export Sales Manager

Medical Trade, s. r. o.
Křeslická 1
101 00 Prague 10
Czech Republic
Tel.: +420 724 380 730

MEDICALTRADE s.r.o.
Křeslická 1, 101 00 Praha 10
Tel: 272 765 043
ICO: 25790579, DIČ: CZ25790579



MK-dent

COLLECTION

2019 / 2020

TURBINES CLASSIC LINE

Reference Number	HC20KL	HC21KL	HC20K	HC21K	HC22KL	HC22K
Light	V	V	-	-	V	-
Head Size	Standard head	Standard head	Standard head	Standard head	Small head	Small head
Head Dimensions (mm)	Ø12.5 x H 14.3	Ø12.5 x H 14.3	Ø12.5 x H 14.3	Ø12.5 x H 14.3	Ø10.8 x H 12.5	Ø10.8 x H 12.5
Coupling Type	Kavo	Kavo	Kavo	Kavo	Kavo	Kavo
KEY FEATURES						
Bur Guide	Carbide insert	Carbide insert	Carbide insert	Carbide insert	Carbide insert	Carbide insert
Back Cap	Push button	Push button	Push button	Push button	Push button	Push button
Bearing	Ceramic	Ceramic	Ceramic	Ceramic	Ceramic	Ceramic
Spray Outlet	4 asymmetrical mixed	4 asymmetrical mixed	4 asymmetrical mixed	4 asymmetrical mixed	3 asymmetrical mixed	3 asymmetrical mixed
TECH SPECS						
Power (W)	28	23	28	23	18	18
Non-load Speed (rpm)	375 000	375 000	375 000	375 000	420 000	420 000
Sound Level (db)	Max. 57	Max. 57	Max. 57	Max. 57	Max. 60	Max. 60
Weight (g)	63	63	63	63	63	63
Tubing	6-pin	6-pin	4-pin	4-pin	6-pin	4-pin
Coupling	QC6016KW (LED) QC6016K (LED) QC5016KW (Xenon) QC5016K (Xenon)	QC6016KW (LED) QC6016K (LED) QC5016KW (Xenon) QC5016K (Xenon)	QC4014K QC4012K	QC4014K QC4012K	QC6016KW (LED) QC6016K (LED) QC5016KW (Xenon) QC5016K (Xenon)	QC4014K QC4012K
Warranty (months)	18	18	18	18	18	18

Nr. 1.



BARGTEHEIDE, 18.08.2021

MK-dent

Marie-Curie-Str. 2, 22941 Bargteheide Germany

To whom it may concern,

We, the MK-dent GmbH, located in Marie-Curie-Strasse 2 in 22941 Bargteheide GERMANY, confirm that our turbines fulfill the criteria of internal water-air cooling across the entire product range.

Please feel free to contact us for further information

With kind regards

FABIAN KLEENE
VICE PRESIDENT - MARKETING & SALES
MK-DENT GMBH

+49 4532 400 49 42
FA.KLEENE@MK-DENT.COM

PHONE +49 4532 400 49-0 FAX +49 4532 400 49-10 INFO@MK-DENT.COM MK-DENT.COM

BANK: HAMBURGER SPARKASSE BIC: HASPDE3333 EUR IBAN: DE53 2005 0550 1354 1204 36 USD IBAN: DE21 2005 0550 1640 1326 90
AG LÜBECK HRB 4272 AH MICHAEL KLEENE (MANAGING DIRECTOR) VAT NO.: DE207950833 TAX NO.: 30 293 2502 9

MK-DENT GMBH, MARIE-CURIE-STR. 2, 22941 BARGTEHEIDE, GERMANY



MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide, Germany


Bargteheide, 20.03.2020

Product confirmation

To whom it may concern, We, MK-dent GmbH located at Marie-Curie-Strasse 2 in 22941 Bargteheide in Germany, confirm that our manufactured dental turbines and contra angles are:


- 1) compressed air driven,
- 2) connected to sleeve according to ISO standard Ritter / Midwest,
- 3) can be autoclaved in Type B autoclaves.

MK-dent[®] GmbH
 Marie-Curie-Str. 2
 22941 Bargteheide
 Tel: +49 321 40049-0 Fax: +49 321 40049-10



Martin Kryniecki
Partner Development Manager
MK-dent GmbH

MK-dent GmbH	Phone: +49 (0)4532 40049-0	Bank: Hamburger Sparkasse	Amtsgericht Ahrensburg 90 HRB 4272
Marie-Curie-Str. 2	Fax: +49 (0)4532 40049-10	Swift Code: HASP DE HHXXX	Michael Kleene
22941 Bargteheide	E-mail: info@MK-dent.com	Bank Code: 200 505 50	VAT-No.: DE207950833
Germany	Internet: www.MK-dent.com	Account No.: 1354/120436	TAX-No.: 30 293 2502 9



1 Explanation of symbols

	See chapter 2. Safety
	Important information for the user
	Product can be sterilized in a steam sterilizer (Autoclave)
	Can be used in thermal disinfectors
	CE-mark – Confirms that this product fulfills all requirements for medical products
REF	Reference number
SN	Serial number

2 Safety

2.1 Description of danger levels

	CAUTION
	CAUTION is used if a deficiency in adequate care may result in an endangering of the patient, the user or others.
	WARNING
	WARNING is used if improper use may result in serious injuries.

2.2 Safety notes

	WARNING
	Risk of injury for the patient and the operator Concerning damages like irregular running noise, irregular vibrations, an unspecific rise in temperature and/or other defects. → Stop operating immediately and contact our service.
	WARNING
	Risk of infection In case of an injury of the tissue in the oral area, do not proceed to work with air-operated instruments because of a high risk of an infection.



Product Service

EC Certificate

Production Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class IIa, IIb or III)

No. G2 18 05 51317 005

Manufacturer: **MK-dent GmbH**

Marie-Curie-Straße 2
22941 Bargteheide
GERMANY



Facility(ies):

MK-dent GmbH
Marie-Curie-Straße 2, 22941 Bargteheide, GERMANY

**Product
Category(ies):**

**Turbines, high and low speed handpieces,
prophylaxis handpieces, motors and couplings
including spare parts and accessories for
the dental area**

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex V. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class IIb and III devices an additional Annex III certificate is mandatory. See also notes overleaf.

Report No.: 713132582

Valid from: 2018-08-05

Valid until: 2023-08-04



Date, 2018-06-08

Stefan Preiß

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

Page 1 of 1



MK-dent GmbH, Marie-Curie-Str.2, 22941 Bargteheide, Germany

Bargteheide ,01.06.2020

Letter of Authorization

This is to certify that
MK-dent GmbH
Marie-Curie-Strasse 2
22941 - Bargteheide
Germany

Appoints
Limeta Jsc
Graiciuno Str. 4
02241 Vilnius
Lithuania

, as an authorized dealer for importation, distribution and repairs of **MK-dent** products in:

Lithuania

This appointment is effective since 01.06.2020. This contract is valid for 1 year from the date of issue and it is automatically renewed year by year unless cancelled by one of both parties at least 30 days before each natural expiry.

MK-dent GmbH

(authorized signature)

Fabian Kleene
Vice President
Sales & Marketing

MK-dent GmbH	Phone:	+49 (0)4532 40049-0	Bank:	Hamburger Sparkasse	Amtsgericht Ahrensburg 90 HRB 4272
Marie-Curie-Str.2	Fax:	+49 (0)4532 40049-10	Swift Code	HASP DE HHXX	Michael Kleene
22941 Bargteheide	E-mail:	info@mk-dent.com	Bank Code:	200 505 50	VAT-No.: DE207950833
Germany	Internet:	www.mk-dent.com	Account No.:	1354/120436	TAX-No.: 30 293 2502 9

mectron

 PRODUCT CATALOGUE

2021

→ TECHNICAL DATA

operating frequency:
24.000 – 36.000 Hz
power types:
adjustable via the dental chair
power supply:
24 VAC, 50-60 Hz / 32 VDC
water supply:
via the dental unit

warranty:
2 years, handpiece 1 year

→ COMPACT PIEZO LED

Equipment set
1 electronic module
1 LED ultrasonic handpiece
1 basic kit equipped with: scaling inserts S1
and S6, perio inserts P1 and P10

Ref. no. 05090008

→ COMPACT PIEZO P2K

Equipment set
1 electronic module
1 ultrasonic handpiece
1 basic kit equipped with: scaling inserts S1, S2,
S3 and 1 high-efficiency scaling insert S6

Ref. no. 05090003

→ COMPACT PIEZO LED
→ COMPACT PIEZO

Compact efficiency or directly built-in.



Including feedback-system for constant power.

SPARE PARTS

High quality spare parts for your ultrasound device.



DSC adaptor for ultrasound inserts



scaler insert adaptor

ADAPTORS FOR THERMODISINFECTION

Product	Ref. no.
DSC Adaptor (for scaler handpiece irrigation line)*	04610007
Spare filter for DSC Adaptor*	04590006
Scaler insert Adaptor*	04610009

COMMON SPARE PARTS

dynamometric torque wrench K10	02900137-001
dynamometric torque wrench K6	02900074
dynamometric wrench K7 (for tipholder DB1)	02900081
water filter	00420004
insert tray	02900072
universal cart	03540011
universal cart with water & air connections	03540012

* not disposable



dynamometric torque wrench K10



dynamometric torque wrenches K6 and K7



water filter

ULTRASOUND



→ every insert
is supplied with its
own torque wrench

→ SCALING INSERTS

Product	Ref. no.
insert S1 universal curette with semicircular diameter for considerable tartar removal	02960001-001
insert S1-S slim universal curette with triangled slightly curved surface for removal of supragingival tartar, subgingival concretions and biofilm	02960009-001
insert S2 universal curette with triangled, slightly curved surface, efficient in the interdental spaces and posterior surfaces	02960002-001
insert S3 flat, with rounded edges for considerable supragingival tartar removal	02960003-001
insert S4 universal curette with 45° angled, triangled and slightly curved surface for tartar removal on mesial and distal interdental surfaces in the posterior area	02960004-001
insert S5 similar shape to S1 insert but longer and thinner for gentle supra- and subgingival tartar removal and for gingivitis	02960005-001

→ INSERTS

For a wide range of applications

→ SCALING



→ SCALING HIGH EFFICIENCY



→ every insert is supplied with its own torque wrench

→ SCALING HIGH EFFICIENCY INSERTS

Product	Ref. no.
insert S6 contra-angled universal curette with semicircular diameter, powerful insert (twice the power of S1 insert) for considerable tartar removal	02960006-001
insert S7 contra-angled universal curette with triangled, slightly curved surface, powerful insert (twice the power of S2 insert) efficient in the interdental spaces and posterior surfaces	02960007-001
insert S8 contra-angled, flat working surface with round edges, powerful insert (twice the power of S3 insert) for considerable supragingival tartar removal	02960008-001

DICHIARAZIONE DI CONFORMITA'**Nome dispositivo:** COMPACT PIEZO LED**Numero di serie / serial number:** 108003863**Tipologia / Typology:** Ablatore ad Ultrasuoni/Ultrasonic Scaler

La Ditta Mectron S.p.A., fabbricante delle tipologie di dispositivi medici:

"Dispositivi per chirurgia ossea, Abiatori ad Ultrasuoni, Apparecchi combinati ablatore e rimozione della placca, Apparecchi per la rimozione della placca, Lampade Polimerizzatrici ed Inserti per chirurgia ossea", dichiara sotto la propria responsabilità che tali apparecchi soddisfano le disposizioni applicabili delle seguenti direttive e dei corrispondenti decreti di attuazione:

- a) Direttiva 93/42/CEE - recepita in Italia con D.Lgs n. 46 del 24-2-1997, modificata dalla Direttiva 2007/47/CEE - recepita in Italia con D.Lgs n. 37 del 25-1-2010.

*Mectron S.p.A., manufacturer for the following Medical Devices:**"Bone surgery units, Ultrasonic Scalers, Ultrasonic Scalers and Jet-Cleaner Polishers, Jet-Cleaner Polishers, Light Curing Units and Inserts for bone surgery", declares on its responsibility that those apparatus comply with the applicable requirements of the following directives:*

- a) *Directive 93/42/EEC - D.Lgs. 24 February 1997 n. 46, amended by Directive 2007/47/EEC - D.Lgs. 25 January 2010 n. 37.*


A tale scopo Mectron S.p.A. dichiara e garantisce che:

- 1) Il dispositivo in oggetto soddisfa i requisiti essenziali richiesti dall'Allegato I delle Direttive menzionate al punto a);
- 2) La progettazione del dispositivo in oggetto soddisfa le disposizioni applicabili dell'Allegato I delle Direttive menzionate al punto a);
- 3) La fabbricazione soddisfa le disposizioni applicabili dell'Allegato II (Sistema Completo di Garanzia di Qualità) escluso il punto 4 delle Direttive menzionate al punto a);
- 4) E' garantito un sistema di Assicurazione della Qualità in conformità alla norma UNI EN ISO 9001-2015 ed UNI EN ISO 13485-2016;
- 5) Il dispositivo in oggetto è da considerarsi appartenente alla classe IIa.

For this purpose Mectron S.p.A. declares and guarantees that:

- 1) *This device complies with the main requisites required by Annex I Directives listed in a);*
- 2) *The planning of this device comply with the applicable requirements of Annex I Directives listed in a);*
- 3) *The manufacture complies with the applicable requirements of Annex II (Full Quality Assurance System) without point 4 of the Directives listed in a);*
- 4) *A system of Quality Assurance is guaranteed in conformity with the UNI EN ISO 9001-2015 and UNI EN ISO 13485-2016;*
- 5) *This device is to be considered as belonging to the class IIa.*

Si dichiara pertanto che i dispositivi in oggetto sono conformi a quanto prescritto dalle Direttive menzionate al punto a), con approvazione dell'Ente Notificato e sono stati immessi in commercio con la marcatura CE, secondo quanto disposto dall'Articolo 17 delle suddette Direttive e Articolo 16 del D.Lgs n. 37 del 25-01-2010.

*Therefore we declare that this device complies with all that is prescribed by the Directives listed in a), with approval of Notified Body and they have been put on the market with the EEC mark according to the Article 17 of the above-mentioned Directives and Article 16 of D.Lgs. 25 January 2010 n. 37.***Ente:** KIWA CERMET ITALIA 0476
Via Cadriano 23**Certificato:** MED 21014A**Data / Date:** 02/07/2021**Il Responsabile**
Riccardo Valle

Anti-Forgery statement

Dear Users:

Thank you for trusting and supporting WOODPECKER products!

Recently we found FAKE woodpecker products on the market, due to the quality is not insured, the fake product can easily break your device and will bring potential risk for patients' safety. In order to protect our users' benefits, we particularly inform that:

1. Finding those who counterfeit and sell imitations, please email to woodpecker@mailgl.cn to report. **We offer 60% of the award judged** by the court to the resultful reporter, and your information will be kept strictly confidential.
2. Authenticity verification of WOODPECKER products can be checked from www.glwoodpecker.com and click "Anti-Forgery Query" for checking.
3. From 15th September 2012, all detachable handpiece as well as all scaling tips that sold separately are pasted with Anti-Forgery labels. Please login into our official website to check.


Scan and Login website
for more information





CURING LIGHT BUILT-IN C USER'S MANUAL

(Please read this manual before operating)



- Certified Management System
- EN ISO 9001
- EN ISO 13485

GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

www.glwoodpecker.com

Content

1. Installation.....	1
2. Principle and usage.....	1
3. Structure and components.....	1
4. Technical specifications.....	1
5. Assembly and disassembly.....	2
6. Operation.....	4
7. Cautions.....	5
8. Contraindications.....	5
9. Maintenance.....	5
10. After service.....	6
11. Trouble shooting.....	6
12. Storage and transportation.....	6
13. Environmental protection.....	7
14. Packing list.....	7
15. Manufacturer's right.....	7
16. For technical data, please contact.....	8
17. Symbol instruction.....	8
18. Statement.....	9
19. Declaration of conformity.....	10

1. Installation

Guilin Woodpecker Medical Instrument Co., Ltd. is a high-tech enterprise in researching, developing, and producing dental equipment, and has a perfect quality assurance system, main products including ultrasonic pizeo scaler, curing light, micro motor, apex locator and ultrasurgery etc.

2. Principle and usage

2.1 BUILT-IN C adopts the principle of ray radiation to solidify the light-sensitive resin by shooting at it in a short time.

2.2 This product is used to restore teeth and solidify material for whitening teeth.

3. Structure and components

BUILT-IN C is composed mainly of high power LED, optical fiber, main unit.

4. Technical specifications

4.1 Power supply: 24V~ 50/60Hz

4.2 Rated current: 0.3A

4.3 Applied part: optical fiber

4.4 Lightsource:

4.4.1 Blue light

4.4.2 Wave length: 420nm-480nm

4.4.3 Light intensity: 1000mW/cm²~1200mW/cm²

4.5 Working condition:

4.5.1 Environment temperature: 5°C to 40°C

4.5.2 Relative humidity: ≤80%

4.6 Dimensions: Φ23mm×250mm

4.7 Net weight: 160g

4.8 Consumption power: ≤8W

4.9 Protection type against electrical shock: class II

4.10 Protection against electrical shock: type B

4.11 Protection against harmful ingress of water or particular matter: Ordinary equipment (IPX0)

4.12 Safety in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide: not suitable under this condition.


4.13 Applied part: optical fiber

Notice: the equipment needs to be installed in the medical power supply, and to ensure that the final product meets the requirements of IEC 60601-1.

5. Assembly and disassembly

5.1 Connect the power line which marked with 24V~ to the dental unit 24V~ output, lace the nylon thread to the dental unit, it's ready to be operated.

Caution: when assembling, the power of the dental unit should be cut off. The power line in the cable should be longer than the nylon thread, in case break the power line when operating.



picture 1

5.2 Take off the red cap from the optical fiber and then insert the metal part into the front of BUILT-IN C revolving, make sure to screw the fiber to the end.

5.3 To install the light hood on as show in picture 1.

5.4 Uninstall the LED, just reverse the procedure above.

6. Operation

6.1 Press the time button to choose the solidification time.4 working time modes are available: 5,10, 15, 20 seconds.

6.2 Lightly press the mode key. Following three modes are available.

6.2.1 Full power mode: blue light shine in full power.

6.2.2 Ramping mode: The power of the blue light turns from weak to stronger, and reaches the highest power in 5 seconds.

6.2.3 Pulse mode: blue light work in the mode of pulse.

6.3 During the operation, aim blue light at the position needing solidification. The machine will gives off a sound "di", when it is on. It is strong blue light after 2 seconds with week blue light. Then it counts down to "0" second to end the solidification.

6.4 After the operation, please clean the fiber with calico in order not to affect the light intensity.

6.5 The depth of solidification of composite is no less than 4mm

per 10 seconds.

6.6 The optical fiber can be spinned off by 360° and autoclaved to 135°C and 0.22MPa.

7. Cautions

7.1 When in clinical operation, make sure the light source be aimed at the resin directly improperty position will affect the solidification.

7.2 Avoid aiming at eyes directly.

① WARNING: If the curing light works for 40s continuously, the temperature of the top of optical fiber may reach 56°C.

② WARNING: Do not modify this equipment without authorization of the manufacturer.

8. Contraindications

Heart disease patients, pregnant women and children should be cautious to use this equipment.

9. Maintenance

9.1 Only the optical fiber can be autoclaved under high temperature and pressure, other parts should be waped by clean water or disinfectants if necessary, never immersed.

9.2 After operation each time, clean the optical fiber.

10. After service

Two year warranty according to the warranty card.

11. Trouble shooting

Fault	Cause	Solution
Non-indication Non-act	1. The main unit doesn't connect well with the unit or the power doesn't on.	1. Check with the connection of the main unit and the unit. 2. Make sure the power is on.
Light intensity insufficient	1. The optical fiber is not inserted till the botton. 2. The optical fiber has cracked. 3. There is resin remained on the surface of optical fiber.	1. Insert the optical fiber again correctly. 2. Change the optical fiber. 3. Wipe off the resin.

If any malfunction case was found, please contact with the dealer the unit was purchased or our company.

12. Storage and transportation

12.1 This equipment should be handled carefully, kept away from shaking point, installed or stored at shadowy, dry, cool and

ventilated places.

12.2 Don't store it together with articles that are combustible, poisonous, caustic and explosive.

12.3 This equipment should be stored in the environment where the relative humidity is $\leq 80\%$, the atmosphere pressure is 70kPa to 106kPa and the temperature is -10°C to $+55^{\circ}\text{C}$.

12.4 Excess impact or shake should be avoided during transportation.

12.5 Don't mix it with dangerous articles during transportation.

12.6 Keep it away from sun or snow or rain during transportation.

13. Environmental protection

There is not any harmful element in this equipment. It can be disposed of according to the local laws.

14. Packing list

The components of the equipment are listed in the packing list.

15. Manufacturer's right


We reserve the rights to change the design of the equipment, the technique, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some

differences between blueprint and real equipment, take the real equipment as the norm.

16. For technical data, please contact


EC REP Wellkang Ltd (www.CE-Marking.eu)
29 Harley St., LONDON, W1G 9QR, UK

17. Symbol instruction

	Trademark		CE marked product
	Type B applied part		FDA marked product
IPX0	Ordinary equipment		Class II equipment
	Date of manufacture		Alternating current
	Manufacturer		Recovery
	Used indoor only		Keep dry
	Screw inside/ outside		Handle with care

 Temperature limitation

 Humidity limitation

 Atmospheric pressure for storage

 Appliance compliance WEEE directive

 Consult the accompanying documents

 Authorised Representative in the EUROPEAN
COMMUNITY



- Certified Management System
- EN ISO 9001
- EN ISO 13485

Got the quality management system certification and CE certification issued by TüV Rheinland

18. Statement

All rights of modifying the product are reserved to the manufacturer without further notice. The pictures are only for reference. The final interpretation rights belong to GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD. The industrial design, inner structure, etc, have claimed for several patents by WOODPECKER, any copy or fake product must take legal responsibilities.

19. Declaration of conformity

19.1 Product conforms to the following standards:

EN 60601-1:2006	EN 1041:2008
EN 60601-1-2:2007	EN ISO 14971:2009
EN 61000-3-2:2006	EN ISO 7405:2008
EN 61000-3-3:2008	EN ISO 17664:2004
EN 60601-1-4:1996	EN ISO 17665-1:2006
EN 60825-1:2007	EN ISO 10993-1:2009
EN 980:2008	EN ISO 10993-5:2009
ISO 9687:1993	EN ISO 10993-10:2010

19.2 EMC - Declaration of conformity

Guidance and manufacturer's declaration - electromagnetic emissions		
The model BUILT-IN C is intended for use in the electromagnetic environment specified below. The customer or the user of the model BUILT-IN C should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model BUILT-IN C uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The model BUILT-IN C is suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance & Declaration — electromagnetic immunity			
The model BUILT-IN C is intended for use in the electromagnetic environment specified below. The customer or the user of the model BUILT-IN C should assure that it is used in such an environment.			
Immunity test	IEC 60501 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±8 kV air	±8 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model BUILT-IN C requires continued operation during power mains interruptions, it is recommended that the model BUILT-IN C be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance & Declaration - Electromagnetic Immunity

The model BUILT-IN C is intended for use in the electromagnetic environment specified below. The customer or the user of the model BUILT-IN C should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3V 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model BUILT-IN C, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1.2 \times P^{1/2}$ $d=1.2 \times P^{1/2} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d=2.3 \times P^{1/2} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model BUILT-IN C is used exceeds the applicable RF compliance level above, the model BUILT-IN C should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model BUILT-IN C.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and the model BUILT-IN C**

The model BUILT-IN C is intended for use in electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the model BUILT-IN C can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model BUILT-IN C as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2 \times P^{1/2}$	80MHz to 800MHz $d=1.2 \times P^{1/2}$	800MHz to 2,5GHz $d=2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

Scan and Login website
for more information



Guilin Woodpecker Medical Instrument Co., Ltd.
Information Industrial Park, National High-Tech
Zone, Guilin, Guangxi, 541004 P. R. China

Tel:

Europe Sales Dept.: +86-773-5873196, +86-773-2125222

North America, South America &

Oceania Sales Dept.: +86-773-5873198, +86-773-2125123

Asia & Africa Sales Dept.: +86-773-5855350, +86-773-2125896


Fax: +86-773-5822450

E-mail: woodpecker@mailgl.cn, sales@glwoodpecker.cn

Website: <http://www.glwoodpecker.com>



Wellkang Ltd (www.CE-Marking.eu)
29 Harley St., LONDON, W1G 9QR, UK



200 Patents; 19 invention patents; Exported to 134 countries

WOODPECKER

Since **1989**

CATALOGUE

GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

Spare parts for curing light



Beveled Optical Fiber

- Easy to touch any part of teeth.
- Can radiate vertically at any angles.
- Solve the problem of weak light intensity when radiate vertically with old optical fiber.
- The light intensity is 40% higher than old optical fiber on certain spot when radiate vertically.

Optical fiber
 Black and white optional
 Can be autoclaved under high temperature and high pressure
 1#Dimensions: $\Phi 8\text{mm}$ - $\Phi 8\text{mm}$ / $\Phi 5\text{mm}$
 2#Dimensions: $\Phi 8\text{mm}$
 1#Applicable for:
 LED.B, LED.C, LED.D, LED.F
 LED.G, LED.P
 2#Applicable for: LED.H

Power supply
 Power Input:
 AC100V-240V 50/60Hz
 Power Output / Current:
 DC5.0V / 1.0A
 Applicable for:
 LED.B, LED.C, LED.D
 LED.F, LED.H

Battery

①: ICR18490A Battery Voltage: 3.7V Applicable for: LED.C	②: ICR18490B Battery Voltage: 3.7V Applicable for: LED.D	③: ICR18650A Battery Voltage: 3.7V Applicable for: LED.F
④: ICR18650B Battery Voltage: 3.7V Applicable for: LED.H	⑤: ICR18650C Battery Voltage: 3.7V Applicable for: LED.B	

1# 2# 3#
Light hood

Lightmeter

Bathymeter

Bleaching fiber for LED.F

LED lamp
 3W, 5W

Glasses



WOODPECKER

GUILIN WOODPECKER MEDICAL INSTRUMENT CO.,LTD.

We GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

With legal address at INFORMATION INDUSTRIAL PARK, GUILIN
NATIONAL HIGH-TECH ZONE, GUILIN, GUANGXI, P.R.CHINA

authorize STOMADENT SK s.r.o. with the legal address at: Lackovce, 300,
06601 Lackovce, Slovakia

to be our official dealer to promote, sell and provide service, install, train
employees and perform warranty services the woodpecker series
products manufactured by us as the part of the dental units
manufactured by STOMADENT SK s.r.o.

Jason Lau



Add: Information Industrial Park, Guilin National High-Tech Zone,
Guilin, Guangxi, 541004 P. R. China

Tel: +86-773-2350599,2350527

Fax: +86-773-5822450

E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com

Website: <http://www.glwoodpecker.com>



EC Declaration of Conformity

Manufacturer:

whose single Authorized Representative:

Guilin Woodpecker Medical Instrument Co., Ltd.

MedNet EC-Rep GmbH • Borkstrasse 10 • 48163
Muenster • Germany

We, the manufacturer, herewith declare that the products
Curing Light UMDNS-Code: 16386

MODEL: LED.B、LED.C、LED.D、LED.E、LED.F、BUILT-IN C、BUILT IN-H、LED.H、LUX V、
LED.G、LUX E、LUX VI、i Led、X-Cure、O-Light、B-Cure

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class I according to Annex IX of the Directive 93/42/EEC. It bears the mark



The product concerned has been manufactured under a quality management system according to Annex VII of Directive 93/42/EEC.

following the procedure relating to the EC Declaration of Conformity set out in Annex V and Annex VII of Directive 93/42/EEC.

This Declaration of Conformity covers all medical devices as specified in the product list belonging to this declaration and is only valid in connection with a batch specific Certificate of Compliance for all products concerned bearing the CE mark

The above mentioned declaration of conformity is exclusively under the responsibility of

Company: Guilin Woodpecker Medical Instrument Co., Ltd.

Address: Information Industrial Park, Guilin National High-Tech Zone, Guilin, GuangXi,
541004, P.R.China

杨芸凤 2020.3.20


Preparation, date

王波 2020.3.20

Review, date



Legally binding signature, Function



MANUALE ISTRUZIONI
OPERATOR'S HANDBOOK
MANUEL D'UTILISATION
GEBRAUCHSANWEISUNG
MANUAL DE INSTRUCCIONES

UNI-JET 75 CARENATO

UNI-JET 75 WITH BOX
UNI-JET 75 CAPOTE
UNI-JET 75 MIT GEHÄUSE
UNI-JET 75 CARENADO



- **General running data
Dental Aspirator**

Model	Uni-Jet 75 with box and control panel
Rated voltage	230 V \sim \pm 5%
Rated frequency	50 Hz
Rated current	3,1 A
Insulation class	Class I
Type of appliance	B
Operating conditions	continuous operation
Protection against ingress of liquids	ordinary
Degree of protection against electric shock	type B \uparrow
The motor is protected by a thermal switch	
Output power	0,4 kW
Maximum flow	1250 l/min
Maximum head for continuous operation	1300 mm H ₂ O
Sound pressure level	58 dB(A)*
Other voltages available on request: 240 V \sim 50 Hz 2,95 A - 220 V \sim 60 Hz 3,5 A - 120 V \sim 60 Hz 6,0 A - 110 V \sim 60 Hz 7,0 A - 100 V \sim 60 Hz 7,0 A	
This appliance cannot work in the presence of an anaesthetic mixture flammable with air, oxygen or nitrous oxide	
A bacteriological filter is also available to filtrate the exhausted air	

\sim	Alternating current	IEC 417-5032
\oplus	Earthing	IEC 417-5019
\uparrow	Appliance Type B	IEC 878-02-02
\bigcirc	OPEN (disconnected from the main electrical supply)	IEC 417-5008
I	CLOSED (connected to the main electrical supply)	IEC 417-5007
\downarrow	Dangerous voltage	IEC 878-03-01

*Sound pressure level with canalized air tested according to the standard ISO 3746-1979 (E).
Parameters: r=1 - Background noise \leq 51 dB (A) - Instrument: Brüel & Kjær Type 2232.

Manufactured by CATTANI S.P.A. - PARMA - ITALY



Integrating your business

FROM PRODUCTS TO DISTRIBUTION

Air and vacuum for dentistry

Compressed air and vacuum for a wide range of medical devices located in dental surgeries and clinics, generally in class 2.4.1 in accordance with ISO 8573-1 (5 microns solid particles filtration, dry air at pressure dew point of +3°C / -20°C of atmospheric dew point, without content of oil vapors).



COMPRESSORS COMPRESSORS WITH SUCTION SUCTION ACCESSORIES

Advantages of Ekom compressors

- Oil-free air
- A wide range of air treatment accessories
- Externely quiet operation
- Long service life
- Easy maintenance
- High durability, efficient heat removal system
- Attractive design, easy to fit into treatment rooms

00-41

Areas of use

- Dental chairs/units
- CAD/CAM systems
- Dental laboratories
- Disinfection and sterilization
- Air distribution for small and medium-sized clinics



DENTISTRY ▾

Compressors ▾

DK50 SIMPLE

[DK50-10](#)

DK50 PLUS

[DK50 2V](#)

[DK50 2V/50](#)

[DK50 4VR/50](#)

[DK50 2X2V/110](#)

[DK50 2X4VR/110](#)

[DK50 4x2VT/M](#)

[DK50 3X4VR/M](#)

[DK50 6x2VT/M](#)

[DK50 4X4VRT/M](#)

[DK50 6X4VRT/M](#)

[DK50 9X4VRT/M](#)

Compressors with suction ▾

Suction ▾

Accessories ▾

Air for dentistry

HEALTHCARE ▾

LABORATORY AND DIAGNOSTICS ▾

INDUSTRY ▾

SHEET METAL PROCESSING ▾

DK50-10

Small-size, but powerful compressor with 10 l air tank. Performance up to 75 l·min⁻¹ and stable pressure up to 8 bar make it ideal for one* dental unit and other air-driven devices in a dental clinic.



PDF Leaflet

DK50-10 Z
without noise reduction cabinet

DK50-10 Z/M
without noise reduction cabinet, with dryer

DK50-10 S
with noise reduction cabinet

DK50-10 S/M
with noise reduction cabinet, with dryer

TECHNICAL PARAMETERS

[DIAGRAMS](#)

[DOWNLOADS](#)

[PHOTOGALLERY](#)

WITHOUT DRYER

WITH MEMBRANE DRYER

WITH ADSORPTION DRYER

	DK50-10 Z	DK50-10 S	DK50-10 S/M
Dryer		X	X
Noise reduction cabinet		X	✓
Number of dental units		1	1
Nominal voltage / frequency		230 V / 50 (60) Hz 115 V / 60 Hz	
Motor performance		0.55 kW	
Air tank capacity		10 l	
Compressor efficiency (working pressure 5 - 7 bar)		75 l·min ⁻¹ / 5 bar (50 Hz) 85 l·min ⁻¹ / 5 bar (60 Hz)	



EKOM spol. s r.o.
Priemyselná 5031/18
921 01 Piešťany
Slovak Republic

**EC DECLARATION OF CONFORMITY
(GB)**

Document No: DK50/2021/EU
Article: 112000526-000

DÉCLARATION DE CONFORMITÉ CE (F)

N° de document: DK50/2021/EU
Article: 112000526-000

Manufacturer: Ekom spol. s r.o.
Priemyselná 5031/18
921 01 Piešťany
Slovak Republic

Fabricant: Ekom spol. s r.o.
Priemyselná 5031/18
921 01 Piešťany
Slovak Republic

SRN: P58533

Product name: COMPRESSOR

Model: DK50B, DK50BS, DK50-10Z, DK50-10Z/M, DK50-10S, DK50-10S/M, DK50 PLUS, DK50 PLUS/M, DK50 PLUS S, DK50 PLUS S/M, DK50, DK50 PLUS MOBILE, DK50 PLUS/M MOBILE, DK50 2V, DK50 2V/M, DK50 2VS, DK50 2VS/M, DK50 2V/50, DK50 2V/50/M, DK50 2V/50S, DK50 2V/50S/M, DK50 2x2V/110, DK50 2x2V/110/M, DK50 2x2V/110S, DK50 2x2V/110S /M, DK50 2V MOBILE, DK50 2V/M MOBILE, DK50 4VR/50, DK50 4VR/50/M, DK50 4VR/50S, DK50 4VR/50S/M, DK50 2x4VR/110, DK50 2x4VR/110 /M, DK50 2x4VR/110S, DK50 2x4VR/110S /M

Le numéro d'enregistrement unique : P58533

Nom du produit : COMPRESSEUR

Modèle: DK50B, DK50BS, DK50-10Z, DK50-10Z/M, DK50-10S, DK50-10S/M, DK50 PLUS, DK50 PLUS/M, DK50 PLUS S, DK50 PLUS S/M, DK50, DK50 PLUS MOBILE, DK50 PLUS/M MOBILE, DK50 2V, DK50 2V/M, DK50 2VS, DK50 2VS/M, DK50 2V/50, DK50 2V/50/M, DK50 2V/50S, DK50 2V/50S/M, DK50 2x2V/110, DK50 2x2V/110/M, DK50 2x2V/110S, DK50 2x2V/110S /M, DK50 2V MOBILE, DK50 2V/M MOBILE, DK50 4VR/50, DK50 4VR/50/M, DK50 4VR/50S, DK50 4VR/50S/M, DK50 2x4VR/110, DK50 2x4VR/110 /M, DK50 2x4VR/110S, DK50 2x4VR/110S /M

Risk class: class I

Intended purpose: The compressor is used as a source of clean, oil-free compressed air to power active medical devices where the parameters and properties of the compressed air are suitable for the specific application.

Registration code: P68338

Classe de risque : classe I

Usage prévu : Le compresseur est utilisé comme source d'air comprimé propre et sans huile pour alimenter les dispositifs médicaux actifs où les paramètres et les propriétés de l'air comprimé conviennent à l'application spécifique.

Code d'enregistrement : P68338

Noted product is in conformity with technical requirements and applicable regulations:

Regulation: MDR 2017/745

Quality Assurance Standards: EN ISO 13485:2016

Procedural Standards: EN 60601-1:2006 /A1:2013, EN 60601-1-2:2015, EN ISO 14971:2019

Le produit noté est conforme aux exigences techniques et aux réglementations en vigueur :

Règlement : MDR 2017/745

Normes d'assurance-qualité : EN ISO 13485:2016

Normes procédurales : EN 60601-1:2006 /A1:2013, EN 60601-1-2:2015, EN ISO 14971:2019

This declaration of conformity is issued the sole responsibility of Ekom spol. s r.o. We hereby declare that the medical device(s) specified above meet the provision of the Regulation (EU) MDR 2017/745 for medical device. This declaration is supported by the Quality System approval to ISO 13485 issued by DNV GL PRESAFE AS. Certificate No. 279404-2018-AQ CZS-NA-PS Rev. 0.0 valid until: January 30, 2022.

Cette déclaration de conformité est la responsabilité exclusive de Ekom spol. s r.o. Nous déclarons par la présente que le dispositif médical détaillé ci-dessus répond aux normes du Règlement (EU) MDR 2017/745 portant sur les dispositifs médicaux.

Cette déclaration est soutenue par l'approbation de système-qualité de l'ISO 13485 émise par DNV GL PRESAFE AS. N° de certificat 279404-2018-AQ-CZS-NA-PS Rév. 0.0 valide jusqu'au : 30 janvier 2022.



Piešťany, 26.5.2021

.....
Ing. Jozef Kováč
Regulatory Affairs Specialist



Piešťany, 26.5.2021

.....
Ing. Jozef Kováč
Spécialiste des affaires réglementaires



V.A. Graičiūno g. 4, LT-02241 Vilnius, tel. (8~5) 2649696, faks. (8~5) 2602055, el.paštas vilnius@limeta.lt
Kodas 221906050, PVM mokėtojo kodas LT219060515, Lietuvos Respublikos Juridinių asmenų registras.

VŠĮ Jonavos pirminės sveikatos priežiūros centras
Siunčiama per CVP IS

2023-01-10, Nr. 2-014

PATVIRTINIMAS

Nr. 644338

- Mes UAB „LIMETA“, vadovaudamiesi įdiegta kokybės vadybos sistema, užtikriname,
kad pirkimo dokumentuose vertinimo kriterijų sąlygos:
- užtikrinime papildomą garantiją 24 mėn. pratęsimą (po privalomų 24 mėn.).

Viešųjų pirkimų specialistė

Aušra Silickienė