## Menisko siuvimo sistemos

### FAST-FIX

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


2 Deploy preloaded implant 2.


THE TECHNIQUE

The **ULTRA FAST-FIX<sup>®</sup> Meniscal Repair System** is a unique all-inside implant system, offering the fixation strength of an open-vertical mattress stitch without the invasive surgical procedure normally required for suture-based repairs.

When the original FAST-FIX<sup>®</sup> Meniscal Repair System was introduced in 2001, it set the benchmark for minimally-invasive, all-inside repairs. Thanks to its preloaded implants, pre-tied sliding knot, and innovative pusher/cutter device, this innovative system lets you deploy two implants vertically or horizontally on either side of the meniscus, lighten the suture and trim the excess.

Building on its proven clinical success, the all-new **ULTRA FAST-FIX** system adds easier knot sliding and stronger suture. The result is a faster, more secure meniscal repair system that will help maximize the chances of a successful meniscus tear recovery.



1.5

**3** Slide the pre-tied knot, tighten the suture, and trim the excess.

*↑ is antisto peruvito, sistentis mego.*

**4** Use the reverse-curved needle to repair tears on the inferior surface (tibial side).


**5** Completed repair with vertical and horizontal mattress stitches.

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**Fast and easy**

Unlike conventional suture-based repair systems, the ULTRA FAST-FIX system is an implant system with a pre-tied, self-sliding knot that eliminates the need for intra-articular knot tying.

**Strong and clinically proven**

Provides a strong, reproducible and reliable meniscal repair with biomechanical properties equal to that of the open, vertical mattress suture technique.

**Contains no hard device heads**

Minimizes trauma to articular cartilage

**Easy knot sliding with ULTRABRAID® suture**

Offers advantages over traditional polyester suture, including higher knot-breaking strength, increased lubricity and a stronger resistance to fraying.

**PEEK-OPTIMA® high strength non-absorbable, polymer implants**

PEEK-OPTIMA polymer lends confidence that the implant will resist breakage when used with high-strength ULTRABRAID suture

**PLLA absorbable implants**

Bioabsorbable version of ULTRA FAST-FIX.


**Curved and reverse-curved needles**

Curved needles promote safer and easier access to a multitude of tear sites. The reverse-curved needle is designed for repairing tears on the inferior surface. Because the needle's point is on the opposite side of the curve, it can safely enter the inferior area without skiving the meniscus or the tibial plateau.

*1.6*  
*↑ lenktos ir reverseine' adatos*

*1.3 Chemini' sudetis "PeeK-Optima"*


ULTIMATE STRENGTH, N

ELONGATION

STIFFNESS, N/mm

After 1000 cycles between 5N and 20N. All values are means +/- standard errors. Bovine meniscus model.

An evaluation of the mechanical properties of the ULTRA FAST-FIX Meniscal Repair System showed a significant increase in ultimate strength and stiffness, and a decrease in elongation compared to ETHIBOND®, RAPIDLOC® and the conventional, clinically successful FAST-FIX device\*. This in vivo mechanical evaluation of meniscal

1 Zontan T, Lyggers A, Masahl V, Weimann A, Petersen W. Cyclic testing of flexible all-inside meniscus suture anchors. AM J Sports Med. 2005; 33:1-7  
2 Data on file at Smith & Nephew

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Setup

*A.L.*

Each ULTRA FAST-FIX® device contains two 5 mm polymer integrated anchors (resorbable or bio-inert), with a pre-tied, self-sliding knot comprised of #2 non-absorbable, UHMW polyethylene ULTRABRAID® co-braid suture. The entire system is packaged in an easy-to-insert, integrated delivery needle. The anchors are placed into the meniscus sequentially, seated safely beyond the capsule, and are then tightened in a simple manner without the need for arthroscopic knot tying (Photos 2 and 3).

*1.7* The dark blue sheath comes preset to a depth of 25 mm from the tip of the needle and 17 mm from the back of the implant, which has been shown to avoid neurovascular injury while allowing predictable meniscocapsular placement<sup>2</sup>.

Peripheral, popliteal, hiatal, and mid-1/3 medial meniscus tears may require penetration less than the 17 mm allowed by the dark blue sheath. Use of the meniscal depth probe, in conjunction with the trimmable depth penetration limiter (white plastic sheath), allows controlled penetration (Figures 1a and 1b).

If the trimmable depth penetration limiter is used with the split cannula, then the split cannula should be completely split before inserting it over the white depth penetration limiter, to allow for easier removal.


Photo 2. Completed repair - horizontal

Photo 3. Completed repair - vertical x 2

*Arthroscopic -  
Ultra Fast Fix  
anchors suture  
has is only  
5 mm polymer  
integrated (bio-  
resorbable or  
inert), so it is  
permitted system.  
The meniscus, suture  
is #2 non-  
absorbable  
UHMW polyete-  
lene. Ultra fast  
is built into  
sheath. This sys-  
tem is superior  
to the previous  
systems, integrated  
to provide  
anchors.*

*Intercan you post  
cannula with  
meniscus most distal  
suture is top  
suture, it has  
intercan pop-  
liteal back, be  
polyester  
meniscus  
suture*


Figure 1a. Trimmable depth penetration limiter and meniscal depth probe


Figure 1b. Trimmable depth penetration limiter

*5*  
[Circular stamp]  
[Signature]

**Caution:** The opening of the split cannula must be at 90° to the curved ULTRA FAST-FIX® needle to prevent the needle from slipping out as it is introduced into the knee (Figure 3). If using the slotted cannula, make sure the curve of the needle faces down (Figure 3a).

To minimize needle bending, grasp the cannula on the shaft and hold it like a dart or pencil when passing through the fat pad (Photo 4).

**Note:** The pretied, self-sliding knot, included in the ULTRA FAST-FIX device, slides from the first implant (T1) to the second implant (T2). Therefore, placing T1 further away than T2 will facilitate sliding of the knot.

**Note:** Maintaining the needle insertion tip within the arthroscopic view at all times avoids suture tangling.

2. For a horizontal repair, place the first implant (T1) farthest away and advance the needle into the outer meniscal fragment until the implant pops through the meniscus.

For a vertical repair, place the superior implant first and advance the needle into the outer meniscal fragment (bisecting the fragment) until the implant pops through the meniscus (Photo 6).

Using the curved ULTRA FAST-FIX device may facilitate initial penetration.

Using the slotted cannula minimizes needle skiving when accessing more anterior tears. Leave the cannula in to help steer the needle tip.

3. Oscillate the needle approximately 5° and pull the needle out of the meniscus, releasing T1 behind the meniscus (Figure 4).

To reduce the amount of suture in the field of view, slowly pull back on the needle after deploying implant 1. Use a forefinger for control upon withdrawal. Piercing the meniscus by 2-3 mm prior to advancing T2 can also help with suture management.


Figure 3. Blue split cannula positioned 90° to the curve of the delivery needle


Figure 3a. Needle curve facing down


Photo 6. Implant 1 placed superior to tear


Figure 4. Placement of implant 1

1.3

1.5

1.2

1.4

1.2

is untested previous  
 for by insertion  
 needles, create  
 Ultra Fast Fix  
 procedure, suture  
 now placed T  
 implants prior  
 anterior T implants  
 to T2 will posi  
 tioned T1  
 anterior now T2  
 but posterior  
 untested suture  
 suture manage.

positioned T1 as anterior

7  
 8


Photo 7. Implant 2 fully advanced to the ready position

4. Using the tip (rather than the volar pad) of the thumb, slide the gold trigger forward to advance the second implant into the ready position (Photo 7 and Figure 5).

**Note:** It is normal to encounter resistance prior to achieving the ready position. A snap or click is heard when the trigger is fully advanced, ensuring that the implant is fully seated at the end of the needle (Figures 6 and 7).


Figure 5. Implant 2 fully advanced to the ready position

*"U" James J. Esposito  
M.D.*


Figure 6. Proper positioning of implant 2


Figure 7. Improper positioning of implant 2

© 2003 by the American Academy of Orthotronics  
 10000 W. 16th Ave., Suite 1000  
 Odessa, FL 33464  
 888-888-8888  
 www.aao.org

*[Signature]*

9

"U" former fibrocyte

5. Insert the delivery needle to release implant 2.

- For a horizontal repair, insert the needle between the entry point and the first implant, approximately 4–5 mm inferior from implant 1 (Figure 8, Photo 8).
- For a vertical repair, insert the needle approximately 4–5 mm inferior from implant 1 (Figure 8, Photo 9).

A.4



Figure 8. Implant 2 ready for release

Photo 8. Horizontal mattress suture

Photo 9. Vertical mattress suture

yr vadybinke  
Olga Raklevičienė

Handwritten signature


Figure 9. Prior to tightening suture construct

6 Remove the delivery needle from the knee, leaving the free end of the suture. Pull the free end of the suture to advance the sliding knot and reduce the meniscal tear (Figure 9, Photos 10 and 11). It is normal to encounter considerable resistance as the knot is snugged down. It is important to pull the free end of the suture in a line directly perpendicular to the tear site.

Avoid suture breakage by wrapping the suture around several fingers and using the tibia as a fulcrum to provide a tactile feel. Apply slow, increasing tension. In most cases, this steady pulling of the suture will cinch the knot down. As the knot is tightened, it may strangle the free leg of suture, creating a loop of suture. If controlled tightening does not eliminate the loop, place a probe under the tight leg of suture and use it as a pulley.

7. To further snug down the suture construct, thread the free end of the suture through the ULTRA FAST-FIX\* Knot Pusher/Suture Cutter. Both curved and straight knot pushers/suture cutters are available. This threading can be facilitated with the use of the suture funnel.


Photo 10. Hand-tightened suture construct - vertical mattress


Photo 11. Hand-tightened suture construct - horizontal mattress

Handwritten signature and stamp area. The signature is written in cursive and appears to be 'Marta Kalkoviciene'. There is a circular stamp to the left of the signature.

Handwritten initials 'MA' at the bottom right corner.

U" James Fitzpatrick 1.4

8. While holding the suture taut, gently slide the knot pusher/suture cutter to the meniscus to achieve the desired tension (Figure 10). The knot pusher should engage the suture in a direct line and perpendicular to the repair. A manual suture "pull"/"push" maneuver is suggested.
9. Rest the tip against the knot to allow for a 2-3 mm suture tail. Cut the suture by sliding the gold trigger forward (Photos 12 and 13, Figure 11). Alternatively, trim the suture with arthroscopic scissors.

To reduce puckering that may result from the femoral surface repair, the implants can subsequently be placed on the tibial side of the meniscus to help pull down the meniscal flap. The reverse curve ULTRA FAST-FIX® device is recommended for tibial side tears.


Figure 10. Suture construct tensioning


Figure 11. Suture cutting

### Postoperative Care

Reestablish full extension and quadriceps activation early, along with joint kinematics and proprioception involving the entire kinetic chain progression to full weight bearing; limit flexion to 90° for three weeks and to torsion for six weeks. Running is indicated at 8-12 weeks. Cutting activities are indicated at 10-12 weeks. Return to full activity is indicated at 3-6 months. Individualization is based on the stability of the tear, repair construct security, and associated pathology.


Photo 12. Suture cutting - horizontal mattress


Photo 13. Suture cutting - vertical mattress

James Fitzpatrick  
 Orthopedic Surgeon  
 Ocala, Florida  
 [Signature]

## Additional Instruction

Prior to performing this technique, consult the Instruction for Use documentation provided with individual components – including indications, contraindications, warnings, cautions, and instructions.

## References

1. Caborn, Borden, Nyland, Pienkowski: Biomechanical Comparison of the FAST-FIX Meniscal Repair Suture System with Vertical Mattress Sutures and Meniscus Arrows. *The American Journal of Sports Medicine*, Vol 31, No. 3, 2003.
2. Coen, Caborn, Urban, et al: *Arthroscopy* 1998.

Courtesy of Smith & Nephew, Inc.,  
Endoscopy Division

Caution: U.S. Federal law restricts this device to sale  
by or on the order of a physician.

1.2 →

## Pearls

- Prepare site and assess geometry of reduction.
- Reduce tear center-to-center to avoid "dog ears"/gapping/ruffles.
- Approach tear from contralateral portal/view from ipsilateral portal.
- Use curve of needle to increase safety zone of vector and improve insertion positioning.
- Hold delivery needle like a dart to stabilize.
- Insert curved delivery needle within split cannula with convexity up.
- Ease insertion using metal slotted cannula.
- Vertical mattress suture: T1 goes posterior and superior; T2 goes anterior and inferior.
- Insert T2 (implant 2) 4 mm to 5 mm from T1 (implant 1).
- Advance gold slide trigger and implant 2 completely to tip until a click is heard (requires force). Use the tip of the thumb on slide trigger rather than volar pad of thumb.
- If implant 1 does not deploy, it is most likely NOT inserted through the entire meniscal tissue: advance deeper.
- If implant 2 does not deploy, it is most likely NOT advanced to the deployment position at the tip of the delivery needle.
- Thread suture onto knot pusher/suture cutter with suture funnel.
- If the knot does not cinch smoothly, it usually requires a more forceful steady pull which is facilitated by wrapping the suture around several fingers like a pulley and applying traction.
- Cinch knot to obtain compression of the suture across the tear but avoid over-cinching or puckering the tissue.
- Alternate divergent femoral side and tibial (tensile) side suture placement.
- Consider reverse curved devices for tibial side fixation.

\*Trademarks of Smith & Nephew, registered U.S. Patent & Trademark Office.

Endoscopy  
Smith & Nephew, Inc  
Andover, MA 01810  
USA

www.smith-nephew.com  
+1 978 749 1000  
+1 978 749 1108 Fax  
+1 800 343 5717 U.S. Customer Service

©2002, 2008 Smith & Nephew, Inc.  
All rights reserved  
02/2008 1061031 Rev. C

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## ACL/PCL Drill Guide System Components and Accessories

### Cruciate Pins and Wires

#### Drill-Tip Passing Pins

With eyelet for graft passage

Reference #	Description
72202009	2.7 mm x 15" (38 cm) Graduated Drill-Tip Passing Pin, 2.7 mm diameter head, 2.4 mm diameter shaft, sterile, single use (6 per package)
014508*	2.7 mm x 15" (38 cm) Drill-Tip Passing Pin, 2.7 mm diameter head, 2.4 mm diameter shaft, sterile, single use (6 per package)
014485	2.7 mm x 15" (38 cm) Drill-Tip Passing Pin, 2.7 mm diameter head, 2.4 mm diameter shaft, non-sterile, single use (6 per package)
72202817	2.4 mm x 15" (38 cm) Drill-Tip Passing Pin, sterile, single use (6 per package)
7208678	2.4 mm x 15" (38 cm) Graduated Drill-Tip Passing Pin, sterile, single use
1107089	2.4 mm x 12" (30 cm) Beath Pin, two eyelets, sterile, single use (3 per package)

#### Spade-Tip Passing Pin

With eyelet for graft passage

Reference #	Description
013110*	3 mm Spade-Tip Passing Pin, 2.4 mm diameter shaft x 11" (28 cm) long, non-sterile, single use (6 per package)
013111	4 mm Spade-Tip Passing Pin, 2.4 mm diameter shaft x 11" (28 cm) long, non-sterile, single use (6 per package)

#### Trocar-Tip Passing Pin

With eyelet for graft passage

Reference #	Description
014395	2.4 mm x 11" (28 cm) Trocar-Tip Passing Pin, sterile, single use (1 per box)
7207220*	2.4 mm x 17" (43 cm) Trocar-Tip Passing Pin, sterile, single use (1 per box)
72202816	2.4 mm x 11" (28 cm) Trocar-Tip Passing Pin, sterile, single use (6 per package)
7207702	2.4 mm x 17" (43 cm) Trocar-Tip Passing Pin, non-sterile, single use (6 per package)

#### Cruciate Guide Wires

For use with ACUFEX® cruciate guides and cannulated drill bits

Reference #	Description
72202872*	2.4 mm Drill-Tip Guide Wire, 10" (25 cm) long, no eyelet, sterile, single use (6 per package)
014396	2.4 mm Drill-Tip Guide Wire, 10" (25 cm) long, no eyelet, sterile, single use (1 per box)
7207221	2.4 mm Drill-Tip Guide Wire, 12" (30 cm) long, no eyelet, sterile, single use (1 per box)
7207284	PCL Safety Guide Wire, 2.4 mm x 10" (25 cm) long, no eyelet, sterile, single use, calibrated for use with the ACUFEX DIRECTOR PCL Safety Stop (1 per box)
72202389	2.4 mm Drill-Tip Guide Wire, 12" (30 cm) long, no eyelet, sterile, single use (6 per box)

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Kopija tisto

## ACL/PCL Drill Guide System Components and Accessories

### Cruciate Pins and Wires

#### Drill-Tip Passing Pins

With eyelet for graft passage

Reference #	Description
72202009	2.7 mm x 15" (38 cm) Graduated Drill-Tip Passing Pin, 2.7 mm diameter head 2.4 mm diameter shaft, sterile, single use (6 per package)
014508*	2.7 mm x 15" (38 cm) Drill-Tip Passing Pin, 2.7 mm diameter head, 2.4 mm diameter shaft, sterile, single use (6 per package)
014485	2.7 mm x 15" (38 cm) Drill-Tip Passing Pin, 2.7 mm diameter head, 2.4 mm diameter shaft, non-sterile, single use (6 per package)
72202817	2.4 mm x 15" (38 cm) Drill-Tip Passing Pin, sterile, single use (6 per package)
7208678	2.4 mm x 15" (38 cm) Graduated Drill-Tip Passing Pin, sterile, single use
1107089	2.4 mm x 12" (30 cm) Beath Pin, two eyelets, sterile, single use (3 per package)

#### Spade-Tip Passing Pin

With eyelet for graft passage

Reference #	Description
013110*	3 mm Spade-Tip Passing Pin, 2.4 mm diameter shaft x 11" (28 cm) long, non-sterile, single use (6 per package)
013111	4 mm Spade-Tip Passing Pin, 2.4 mm diameter shaft x 11" (28 cm) long, non-sterile, single use (6 per package)

3

#### Trocar-Tip Passing Pin

With eyelet for graft passage

Reference #	Description
014395	2.4 mm x 11" (28 cm) Trocar-Tip Passing Pin, sterile, single use (1 per box)
7207220*	2.4 mm x 17" (43 cm) Trocar-Tip Passing Pin, sterile, single use (1 per box)
72202816	2.4 mm x 11" (28 cm) Trocar-Tip Passing Pin, sterile, single use (6 per package)
7207702	2.4 mm x 17" (43 cm) Trocar-Tip Passing Pin, non-sterile, single use (6 per package)

#### Cruciate Guide Wires

For use with ACUFEX® cruciate guides and cannulated drill bits

Reference #	Description
72202872*	2.4 mm Drill-Tip Guide Wire, 10" (25 cm) long, no eyelet, sterile, single use (6 per package)
014396	2.4 mm Drill-Tip Guide Wire, 10" (25 cm) long, no eyelet, sterile, single use (1 per box)
7207221	2.4 mm Drill-Tip Guide Wire, 12" (30 cm) long, no eyelet, sterile, single use (1 per box)
7207284	PCL Safety Guide Wire, 2.4 mm x 10" (25 cm) long, no eyelet, sterile, single use, calibrated for use with the ACUFEX DIRECTOR PCL Safety Stop (1 per box)
72202389	2.4 mm Drill-Tip Guide Wire, 12" (30 cm) long, no eyelet, sterile, single use (6 per box)

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7207222 pravedimo viela 38 cm

17

# ACL/PCL Drill Guide System Components and Accessories


## Drill Bits, Routers and Bone Graft Components

### Cannulated Drill Bits

Used with ACL and PCL Drill Guide Systems, 2.4 mm I.D., 7" (17.75 cm) length.


Reference #	Description	Reference #	Description
013498	5 mm	013546	10 mm
7207491	5.5 mm	7207487	10.5 mm
013542	6 mm	013547	11 mm
7207483	6.5 mm	7207488	11.5 mm
013543	7 mm	013548*	12 mm
7207484	7.5 mm	7207489	12.5 mm
013544	8 mm	013549	13 mm
7207485	8.5 mm	7207490	13.5 mm
013545	9 mm	013551	14 mm
7207486	9.5 mm		

### Endoscopic Cannulated Drill Bits

Used with Endoscopic ACL Drill Guide Systems for Femoral Tunnel Drilling, 2.4 mm I.D., 7" (17.75 cm) length, calibrated.

Reference #	Description
7207315	4.5 mm (use with ENDOBUTTON®, ENDOBUTTON CL or RCI screws), sterile, disposable, single use
013184*	4.5 mm (use with ENDOBUTTON, ENDOBUTTON CL or RCI screws), non-sterile, disposable, single use
72202971	4.5 mm CLANCY Flexible Drill

Reference #	Description	Reference #	Description
7208686	5 mm	7207495	9.5 mm
7207482	5.5 mm	013663*	10 mm
013499	6 mm	7207496	10.5 mm
7207492	6.5 mm	013664	11 mm
013660	7 mm	7207497	11.5 mm
7207493	7.5 mm	013665	12 mm
013661	8 mm	7207498	12.5 mm
7207494	8.5 mm	013666	13 mm
013662	9 mm		

UAB "Ostecca" (UAB "Ostecca")  
 UAB "Ostecca" (UAB "Ostecca")  
 UAB "Ostecca" (UAB "Ostecca")

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








G.I.

**Forefoot Set**

A new modular and scalable forefoot instruments set:

- Bold 2.5 – 3.0 – 3.7 Compression screw.
- QWIX 3.0 Fixation screw.
- Spin 2.0 – 2.7 Spin break off screws.
- Solustaple Standard staple.
- Uni-clip Compression Staple.

**Frame**

- 229996 Frame - 2 floors
- 229990 Frame - 3 floors
- 229970 Forefoot II- lid- generic

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20

**Hand Fractures and Deformity Repair Set**

6.1

It is a set to be used for fractures.

Only one instrumentation module in order to set up 4 screws:

- Build Compression Screw:
  - > 2.5 mm
  - > 3.0 mm
  - > 3.7 mm
- Qwix Fixation Screw:
  - > 3.0 mm


} 6.1

And two associated optional modules :

- Implant module: for non sterile screws.
- Generic module: for other personal instrumentation and support for 2 or 3 floors for the chosen modules.

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
*[Handwritten signature]*

**Bold® Compression Screw Dia. 2.5 mm / 3.0 mm / 3.7 mm** (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Totally intra-osseous implants. *6.3.*
- Large range of sizes:
  - > 10 to 30 mm for Bold® 2.5 screws by 2 mm increment.
  - > 10 to 34 mm for Bold® 3.0 screws by 2 mm increment.
  - > 14 to 34 mm for Bold® 3.7 screws by 2 mm increment.
- K-Wire guided drilling and insertion.
- Self-tapping.
- True & controlled compression.
- Dual threads. *6.3.*
- Material: Titanium alloy TA6V. *6.2.*

**Bold® dia 2.5 mm**

Reference	Length
200010(S)	L. 10 mm
200012(S)	L. 12 mm
200014(S)	L. 14 mm
200016(S)	L. 16 mm
200018(S)	L. 18 mm
200020(S)	L. 20 mm
200022(S)	L. 22 mm
200024(S)	L. 24 mm
200026(S)	L. 26 mm
200028(S)	L. 28 mm
200030(S)	L. 30 mm

**Bold® dia 3.0 mm**

Reference	Length
111010(S)	L. 10 mm
111012(S)	L. 12 mm
111014(S)	L. 14 mm
111016(S)	L. 16 mm
111018(S)	L. 18 mm
111020(S)	L. 20 mm
111022(S)	L. 22 mm
111024(S)	L. 24 mm
111026(S)	L. 26 mm
111028(S)	L. 28 mm
111030(S)	L. 30 mm
111032(S)	L. 32 mm
111034(S)	L. 34 mm

**Bold® dia 3.7 mm**

Reference	Length
131014(S)	L. 14 mm
131016(S)	L. 16 mm
131018(S)	L. 18 mm
131020(S)	L. 20 mm
131022(S)	L. 22 mm
131024(S)	L. 24 mm
131026(S)	L. 26 mm
131028(S)	L. 28 mm
131030(S)	L. 30 mm
131032(S)	L. 32 mm
131034(S)	L. 34 mm

**QWIX® Fixation Screw, Dia. 3.0 mm** (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Self-drilling and self-tapping.
- Totally intra-osseous fixation.
- K-Wire guided drilling and insertion.
- Material: Titanium alloy TA6V.


**QWIX® Dia. 3.0 mm**

Reference	Length
111312(S)	L. 12 mm
111314(S)	L. 14 mm
111316(S)	L. 16 mm
111318(S)	L. 18 mm
111320(S)	L. 20 mm
111322(S)	L. 22 mm
111324(S)	L. 24 mm
111326(S)	L. 26 mm
111328(S)	L. 28 mm
111330(S)	L. 30 mm
111332(S)	L. 32 mm
111334(S)	L. 34 mm

(S) means available in sterile and non-sterile

## Implant Details


Projekto vadovas:  
Evaldas Ročys

## Indications

- For fixation of bone fractures or for bone reconstruction.
- Examples include:
  - Fixation of osteotomies for Hallux Valgus treatment (such as Scarf, Chevron, etc.)
  - Mono or Bi-cortical osteotomies in the foot or hand.
  - Distal or proximal metatarsal or metatarsal osteotomies.
  - Arthrodesis in hand or foot surgery.
  - Fixation of small bone fragments, in long bones or small bones fractures.

6.4.  
6.5.  
6.6.

## Instruments Details

Screws and K-wire

As straight handle

Instrument kit for fixation of small bone fragments and osteotomies

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
Surgical technique of fixation of small bone fragments and osteotomies  
Products for Sale: Europe, Middle East and Africa only

24

Reference of the manufacturer of the device: please refer to the product and technical manual for the correct usage of the device for use as a specific, stable fixation. The patient's doctor should always be consulted for advice on the correct use of the device. An example of the correct use of the device is shown in the illustration.

## Surgical Technique

- K-wire Insertion**  
Secure the bone fragments with a K-wire (fig. 1). The position should be checked under fluoroscopy.

This K-wire guides the screw. Due to the traumatic aspects of the K-wire, it can be withdrawn and replaced until the desired position is achieved.

1. K-wire insertion into the bone fragment. Diagrams of the device



Screw	Reference	Diameter	Length
Bold: 2.5 mm	1800045	2.5 mm	20 mm
Bold: 3.0 mm	1800046	3.0 mm	20 mm
Bold: 3.7 mm	1800047	3.7 mm	20 mm
Quick: 3.0 mm	1800048	3.0 mm	20 mm

1. K-wire insertion into the bone fragment

## Measurement

Insert the appropriate cannulated screwdriver on the K-wire (fig. 2-4) and read the screw length to be used directly on the scale (fig. 4-2).

**Note**  
If the surgeon wishes a monocortical fixation, subtract 1.5 mm to determine the appropriate screw length.

2. K-wire insertion into the bone fragment. Diagrams of the device

2. K-wire insertion into the bone fragment. Diagrams of the device

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