## VIDAS<sup>®</sup> Anti-HCV (HCV)

VIDAS Anti-HCV is an automated qualitative test for use on the instruments of the VIDAS family, for the detection of IgG antibodies to hepatitis C virus (anti-HCV) in human serum or plasma (heparin) using the ELFA technique (Enzyme Linked Fluorescent Assay). The detection of these specific antibodies, in conjunction with other clinical information, aids in the diagnosis of infection in persons with symptoms of hepatitis and in persons at risk for hepatitis C infection.

### SUMMARY AND EXPLANATION

The Hepatitis C virus (HCV) discovered in 1989 using advanced molecular biology techniques, was rapidly found to account for the majority of those patients with non-A non-B hepatitis. HCV represents a major worldwide public health problem requiring global action for the diagnosis, treatment and prevention of this infection (1).

HCV is primarily parenterally transmitted through direct blood-to-blood contact between two people: use of unsterilized injection devices and transfusion of unscreened blood or blood products (2). The disease frequently progresses to chronic hepatitis C (80%), exposing patients to a greater risk of hepatic complications such as cirrhosis or hepatocellular carcinoma. (3).

The current standard of treatment for HCV is a combination of two drugs: pegylated interferon and ribavirin, but due to the high genetic variability of HCV (4), it is still only partially effective: viral eradication in less than 50% of patients infected with genotype 1 hepatitis C virus against approximately 80% of patients infected with genotype 2 or 3. New therapeutic options are under study to offer more effective and safer personalized treatments (5,6).

Diagnosis of patients infected with HCV can be performed using two categories of virological tests: indirect tests, and direct tests (7). Indirect serological tests are third-generation enzyme immunoassays that detect antibodies to HCV. The antigens used in the tests to detect antibodies are from the structural and non-structural regions of the HCV (8) (capsid, protein, cofactors, polymerase, etc.). The presence of anti-HCV antibodies indicates that an individual may have been infected with HCV in the past or may have an ongoing HCV infection. To confirm the presence of active HCV infection, a positive serological test can be completed using direct tests (e.g.: molecular assays that detect RNA genomes). The results will be used to guide patient management and determine the optimal duration of treatment.

The VIDAS Anti-HCV assay is a third-generation test using antigens corresponding to the HCV core, NS3 and NS4 proteins for the qualitative detection of anti-HCV antibodies.

### PRINCIPLE

The assay principle combines a two-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR<sup>®</sup>) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and are pre-dispensed in the sealed reagent strips.

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During the first step, the sample is diluted, and then cycled in and out of the SPR several times. The anti-HCV antibodies present in the sample will bind to the antigens representing the HCV core, NS3 and NS4 proteins coated on the interior of the SPR. Unbound sample components are washed away.

During the second step, mouse monoclonal anti-human IgG antibodies in Fab form, conjugated to recombinant alkaline phosphatase (yeast) are cycled in and out of the SPR several times and will bind to the human Ig bound to the molecules on the solid phase. Further wash steps remove unbound components.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibody present in the sample. At the end of the assay, the results are automatically calculated by the instrument in relation to the Standard S1 stored in memory, and then printed out.

**CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS:**

60 Strips HCV	STR	Ready-to-use.
60 SPRs HCV (2 x 30)	SPR®	Ready-to-use. Interior of SPRs coated with antigens representing the HCV core, NS3 and NS4 proteins.
HCV Positive control 1 x 1.9 ml (liquid)	C 1	Pooled human serum or plasma* containing anti-HCV IgG in a phosphate buffer + BSA + preservatives Index: The confidence interval is indicated on the MLE card after the following mention: "Control C1 (+) Test Value Range".
HCV Negative control 1 x 1.9 ml (liquid)	C 2	Phosphate buffer + BSA + preservatives Index: The confidence interval is indicated on the MLE card after the following mention: "Control C2 (-) Test Value Range".
Standard 1 x 1.9 ml (liquid)	S 1	Pooled human serum or plasma* containing anti-HCV IgG in a phosphate buffer + BSA + preservatives The confidence interval in "Relative Fluorescence Value (RFV)" is indicated on the MLE card after the following mention: "Standard (S1) RFV Range".
1 MLE card (Master Lot Entry)		Specifications for the factory master data required to calibrate the test: to read the MLE data, please refer to the User's Manual
1 Package Insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

\* This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1 and HIV2. The product has been inactivated. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

**The SPR**

The interior of the SPR is coated during production with the antigens representing the HCV core, NS3 and NS4 proteins. Each SPR is identified by the code "HCV". Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening.**

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**Description of the HCV strip**

Wells	Reagents
1	Sample well.
2	Sample diluent: TRIS buffered saline + Tween 20 + BSA + preservatives (600 µl)
3 – 4 – 5 – 7 - 8	Wash buffer: TRIS buffered saline + Tween 20 + preservatives (600 µl)
6	Conjugate: mouse monoclonal anti-human IgG antibodies conjugated to recombinant ALP in Phosphate buffered saline + protein stabilizer + preservatives (400 µl)
9	Reactive diluent: Phosphate buffered saline + preservative (400 µl)
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

**\* IRRITANT reagent:**

- **R 36:** Irritating to eyes.
- **S 26:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

For further information, refer to the Safety Data Sheet available on request.

## MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS family.

## WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPR<sup>®</sup>s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the risk phrase "R" and the precautions "S" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

## STORAGE CONDITIONS

- Store the VIDAS Anti-HCV kit at 2-8°C.
- Do not freeze reagents.
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the desiccant inside and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

## SPECIMENS

### Specimen type and collection

Human serum or plasma

#### **Types of tubes validated:**

- Plain tube,
  - Tube with lithium heparin,
  - Tube with sodium heparin,
  - Tube with lithium heparin and separation gel,
  - Plastic tube with clot activator,
  - Plastic tube with clot activator and separation gel.
- The use of heat-inactivated sera has not been validate.

**Note:** Depending on the manufacturer, blood collection tubes may contain materials and additives that could generate different test results.

It is the responsibility of each laboratory to validate use of these tubes in accordance with the manufacturer's recommendations for use.

### Specimen preparation

Plain tubes: wait for samples to coagulate and **centrifuge** according to the tube manufacturer's recommendations to eliminate fibrin.

Other tubes: follow the tube manufacturer's recommendations for use.

Frozen-stored samples: after thawing, these samples must be homogenized before analysis.

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- hemolysis (after spiking samples with hemoglobin: 0 to 300 µmol/l (monomer),
- lipemia (after spiking samples with lipids: 0 to 30 mmol/l equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 220 mg/ml or 376 µmol/l).

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### **Do not inactivate samples.**

### Specimen stability

Serum and plasma samples separated from the clot, can be stored at 2-8 °C in stoppered tubes for 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C.

Do not exceed 3 freeze/thaw cycles.

A study performed on samples frozen for 12 months, showed that the quality of results is not affected.

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### VIDAS PTC protocol data entry

When using the assay for the first time, **and before reading the MLE data**, scan the bar code(s) (at the end of the package insert) using the instrument's external bar code reader. This reading will allow VIDAS PTC protocol data to be transferred to the instrument software for its update. These data should only be read the first time the assay is used.

### Master lot data entry

**Note: When using the assay for the first time, enter the VIDAS PTC protocol (bar codes at the end of the package insert) before reading the MLE data. If the MLE data have been read before the VIDAS PTC protocol, read the MLE data again.**

Before each new lot of reagents is used, specifications (or factory master calibration data) must be entered into the instrument using the MLE data. If this operation is not performed **before initiating the tests**, the instrument will not be able to print results. The master lot data need only be entered once for each lot.

It is possible to enter MLE data manually or automatically depending on the instrument (refer to the User's Manual).

### Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 28 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by "S1", must be tested **in duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

### Procedure

1. **Only remove the required reagents from the refrigerator. They can be used immediately.**
  2. Use one "HCV" strip and one "HCV" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
  3. The test is identified by the "HCV" code on the instrument. The standard must be identified by "S1", and tested in duplicate. If the positive control is to be tested, it should be identified by "C1". If the negative control is to be tested, it should be identified by "C2".
  4. If necessary, clarify samples by centrifugation.
  5. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
- |  |
|--|
| 6. <b>For this test, the standard, control, and sample test portion is 100 µl.</b> |
|--|
7. Insert the "HCV" SPRs and the "HCV" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
  8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.

9. Restopper the vials and return them to 2-8°C after pipetting.
10. The assay will be completed within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument.

The patient RFV is interpreted as follows:

Test value = (patient RFV / standard RFV).

This test value and the interpretation are also included on the result sheet. The interpretation depending on the test value is as follows:

Test value (TV)	Interpretation
<1.00	negative
≥1.00	positive

All positive patient results must be verified in duplicate. If at least one of the repeat values is positive, the patient result is considered as positive. In that case, additional tests should be performed (another immunoassay or HCV marker) on the same sample or on a second one.

**Note: In all cases, refer to current national guidelines concerning HCV diagnosis.**

**Interpretation of VIDAS Anti-HCV test results should be made taking into consideration the patient history and the results of any other tests or hepatitis C markers.**

### QUALITY CONTROL

One positive control and one negative control are included in each VIDAS Anti-HCV kit.

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## LIMITATIONS OF THE METHOD

For the diagnosis of HCV infection, the serological results should be used and interpreted taking into account the patient history, the clinical record, and further tests.

A negative test result does not exclude the possibility of exposure to HCV or infection with HCV. Anti-HCV antibodies may be undetectable in some stages of the infection (acute phase of hepatitis or presence of a serological scar) and in some clinical conditions (immunosuppression) (7,9).

Interference may be encountered with certain sera containing antibodies against reagent components.

This test has not been validated for use with any specimen matrices other than human serum or plasma.

## RANGE OF EXPECTED VALUES(1)

Hepatitis C has a worldwide prevalence of 2-3% that varies according to country:

Region of the world	Anti-HCV prevalence (%)
Europe	2.3
Africa	3.2
Americas	1.5
Australia & Oceania	1.2
Asia	2.1
Middle East	4.7
Total	2.4

## PERFORMANCE

The following study results demonstrate the conformity of VIDAS Anti-HCV to the Common Technical Specifications of 98/79/CE Directive:

### 1. Specificity for blood donor population:

5104 blood donor samples (including 2904 fresh samples with a negative status collected ≤ 24 hours previously) obtained from 2 blood transfusion centers, were tested using the VIDAS Anti-HCV assay.

VIDAS Anti-HCV	Status	
	Positive	Negative
Positive	0	20
Negative	0	5084

Diagnostic specificity of the VIDAS Anti-HCV assay on this population: 99.61%

(95% confidence interval: 99.40% - 99.76%)

### 2. Clinical specificity for hospitalized patients

200 samples with a negative status were tested using the VIDAS Anti-HCV assay.

Diagnostic specificity of the VIDAS Anti-HCV assay on this population: 99.50%.

(95% confidence interval: 97.25% - 99.99%)

### 3. Diagnostic sensitivity

439 samples with a positive status, including 102 fresh samples (collected ≤ 24 hours previously), were tested using the VIDAS Anti-HCV assay.

Genotypes 1 to 6 were tested:

Genotype	Number tested
1	21
2	21
3	23
4	22
(including non-a sub-types)	
5	6
6	2

Results on tested populations:

Population	Positive VIDAS Anti-HCV /total tested	Diagnostic sensitivity observed (95% confidence interval)
HCV/HIV-negative patient	254/254	100% [98.56% - 100%]
HCV/HIV-positive patient	60/61*	98.36% [91.20% - 99.96%]
Patient with unknown HCV/HIV status	124/124	100% [97.07% - 100%]
Total HCV population	438/439*	99.77% [98.74% - 99.99%]

\* The patient who was not detected using VIDAS Anti-HCV either had a low antibody level or was not detected using equivalent methods.

### 4. Sensitivity for seroconversion panels

Testing of 30 seroconversion panels demonstrated the precocity of detection of the VIDAS Anti-HCV assay. The results are comparable to those obtained using the most sensitive methods.

## 5. Precision

The repeatability and reproducibility were determined at two sites and calculated according to the recommendations of the CLSI® documents EP5-A2 / EP12-A2.

Four human samples were tested in duplicate using two lots of reagents. Testing was performed twice a day for 10 days on three instruments at one experimental site (N=120). Each reagent lot used a single calibration curve throughout the study. Data from this study are summarized in the following table:

Sample ID / Target value	Sample 1		Sample 2		Sample 3		Sample 4	
	0.26		0.93		1.08		1.19	
	Standard deviation	CV (%)						
Repeatability	0.01	5.6	0.04	4.3	0.05	4.8	0.06	4.7
Reproducibility	0.07	26.9	0.05	5.9	0.07	6.3	0.07	5.8

## 6. Cross-reactivity

273 samples from patients with a physiological status that can potentially interfere with the detection of hepatitis C antibodies, were tested using VIDAS Anti-HCV. All of the samples were found to be negative with another EIA method (except one CMV IgG+ sample). No disease-related interference was observed for VIDAS Anti-HCV.

	VIDAS Anti-HCV
HSV +	0/10
VZV +	0/10
EBV +	0/10
HIV +	0/10
CMV IgG +	1/11*
LYME Ig+	0/10
HAV IgG +	0/10
HVB (HBcT +)	0/8
HVB (Ag HBs +)	0/10
Syphilis	0/10
Rubella IgG +	0/10
Toxoplasmosis IgG +	0/10
Rheumatoid factor	0/10
Anti-Nuclear Antibody	0/10
Anti-E. coli antibody	0/10
Anti-Pichia Antibody	0/10
Pregnant women**	0/114

\* The reference EIA method also showed one false positive sample, but on a different sample.

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**REFERENCES BIBLIOGRAPHIQUES**

1. LAVANCHY D. Evolving epidemiology of hepatitis C virus - Clin Microbiol Infect 2011; 17: 107-115.
2. ALTER MJ. Epidemiology of hepatitis C virus infection, World J Gastroenterol 2007; 13(17): 2436-2441.
3. GURTSEVITCH VE. Human Oncogenic Viruses: Hepatitis B and Hepatitis C Viruses and Their Role in Hepatocarcinogenesis, Biochemistry (Moscow), 2008, 73(5): 504-513.
4. SIMMONDS P. et al. Consensus Proposals for a Unified System of Nomenclature of Hepatitis C Virus Genotypes, HEPATOLOGY 2005;42: 962-973.
5. VERMEHREN J, SARRAZIN C. New HCV therapies on the horizon, Clin Microbiol Infect 2011; 17: 122-134.
6. WEBSTER D.P, KLENERMAN P., COLLIER J., JEFFERY K. JM. Development of novel treatments for hepatitis C, Lancet Infect Dis 2009; 9: 108-17.
7. CHEVALIEZ S. Virological tools to diagnose and monitor hepatitis C virus infection, Clin Microbiol Infect 2011; 17: 116-121.
8. PENIN F, DUBUISSON, REY F-A, MORADPOUR D, ET PAWLOTSKY JM. Structural Biology of Hepatitis C Virus, HEPATOLOGY 2004; 39: 5-19.
9. THIO CL, NOLT KR, ASTEMBORSKI J, VLAHOV D, NELSON KE et THOMAS DL. Screening for Hepatitis C Virus in Human Immunodeficiency Virus-Infected Individuals, Journal of clinical microbiology, Feb. 2000, 38(2): 575-577.

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests

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BIOMERIEUX, the blue logo, VIDAS and SPR are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

CLSI is a trademark belonging to Clinical and Laboratory Standards Institute Inc.

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1. LAVANCHY D. Evolving epidemiology of hepatitis C virus - Clin Microbiol Infect 2011; 17: 107-115.
2. ALTER MJ. Epidemiology of hepatitis C virus infection, World J Gastroenterol 2007; 13(17): 2436-2441.
3. GURTSEVITCH VE. Human Oncogenic Viruses: Hepatitis B and Hepatitis C Viruses and Their Role in Hepatocarcinogenesis, Biochemistry (Moscow), 2008, 73(5): 504-513.
4. SIMMONDS P. et al. Consensus Proposals for a Unified System of Nomenclature of Hepatitis C Virus Genotypes, HEPATOLOGY 2005;42: 962-973.
5. VERMEHREN J, SARRAZIN C. New HCV therapies on the horizon, Clin Microbiol Infect 2011; 17: 122-134.
6. WEBSTER D.P, KLENERMAN P., COLLIER J., JEFFERY K. JM. Development of novel treatments for hepatitis C, Lancet Infect Dis 2009; 9: 108-17.
7. CHEVALIEZ S. Virological tools to diagnose and monitor hepatitis C virus infection, Clin Microbiol Infect 2011; 17: 116-121.
8. PENIN F, DUBUISSON, REY F-A, MORADPOUR D, ET PAWLITSKY JM. Structural Biology of Hepatitis C Virus, HEPATOLOGY 2004; 39: 5-19.
9. THIO CL, NOLT KR, ASTEMBORSKI J, VLAHOV D, NELSON KE et THOMAS DL. Screening for Hepatitis C Virus in Human Immunodeficiency Virus-Infected Individuals, Journal of clinical microbiology, Feb. 2000, 38(2): 575-577.

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**VIDAS<sup>®</sup> HBs Ag Ultra (HBS)**

VIDAS HBs Ag Ultra (HBS) is an automated qualitative test for use on the VIDAS family instruments for the detection of hepatitis B surface antigen (HBs Ag) in human serum or plasma, using the ELFA technique (Enzyme Linked Fluorescent Assay).

**SUMMARY AND EXPLANATION**

The hepatitis B virus is responsible for acute and chronic hepatitis infections. Acute hepatitis can be asymptomatic or present symptoms of varying severity which may progress to fulminant hepatitis in 0.1 to 0.5% of cases. Chronicity occurs in 5 to 10% of cases in adults, but up to 90% of cases in infants following perinatal transmission. Currently, approximately 350 million people worldwide are chronic carriers of the virus (1). Chronic hepatitis may be asymptomatic and lead to liver lesions of varying severity, possibly evolving to cirrhosis, with an evolution in 5% of cases to hepatocellular carcinoma (2). The hepatitis B virus can be transmitted by parenteral or perinatal pathways or through sexual contact. Persons most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, as well as close friends and family of an infected subject, and newborns of an infected mother (2).

The discovery of the Australia antigen in 1965 – later known as HBs antigen - combined with viral hepatitis (3, 4) was a major breakthrough in the diagnosis of hepatitis B.

HBs antigen appears several days to several weeks after contact with the virus and can persist for several months. Persistence of HBs antigen for more than 6 months serologically defines chronic HBV infection. Disappearance of the HBs antigen is normally followed by the appearance of anti-HBs antibodies, which is a sign of recovery. The anti-HBs antibody assay is performed to confirm efficacy of the HBV vaccine. Anti-HBc antibodies are normally detected at the onset of the disease (2). Anti-HBc IgM titers are elevated (> 100 UPEI/ml) during acute hepatitis, then their titer decreases or disappears. However, during chronic hepatitis, the appearance of low IgM titers, reflecting hepatic cytolysis, confirms the active phase of the disease (5). Total anti-HBc antibodies, mainly IgG, are detected during acute and chronic hepatitis and persist after recovery.

HBe antigen is a circulating protein with a sequence very similar to that of HBc antigen, but with distinct antigenicity. During acute or chronic hepatitis, the presence of HBe antigen is generally associated with intense viral replication. Its disappearance coincides with the appearance of anti-HBe antibodies. Anti-HBe seroconversion is generally an indicator of a favorable prognosis of recovery (6). Nevertheless, the possible selection of a mutant strain incapable of synthesizing Hbe antigen may also lead to anti-HBe seroconversion, thereby limiting the prognostic value of HBe and anti-HBe markers (7). In this case, only the viral DNA assay is capable of directly detecting viral replication.

**PRINCIPLE**

The VIDAS HBs Ag test is an enzyme-linked fluorescent immunoassay (ELFA) that is performed in the automated VIDAS system (see User's Manual). This assay can be performed according to 2 protocols: HBL long protocol (90 minutes), HBS short protocol (60 minutes).

The Solid Phase Receptacle (SPR<sup>®</sup>) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready to use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

After a preliminary washing step, the antigen present in the sample will bind simultaneously to the monoclonal antibody (8) coating the interior of the SPR and to the antibody conjugated with biotin. Unbound sample components are washed away. The antigen bound to the solid phase and to the biotinylated antibody is in contact with streptavidine conjugated with alkaline phosphatase, which will bind with biotin. Another wash step follows and removes unbound components.

During the final detection step, the substrate (4-Methylumbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methylumbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample.

At the end of the assay, results are analyzed automatically by the instrument and are expressed as an index calculated using a standard.

**CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS:**

60 HBS strips	STR	Ready-to-use.
60 HBS SPRs 2 x 30	SPR <sup>®</sup>	Ready-to-use. Interior of SPRs coated with two monoclonal anti-HBs Ag antibodies (mouse) (5).
HBS standard 3 x 1 ml (lyophilized)	S1	Serum base* supplemented with inactivated human plasma HBs antigen + preservatives. The standard must be reconstituted with 1 ml of sterile distilled water (measured exactly). Allow to dissolve for at least 20 mins, then mix using a vortex. After reconstitution, the standard must be stored in aliquots at $-25 \pm 6^{\circ}\text{C}$ for up to 6 months, by freezing within an hour of reconstitution. Avoid successive freezing and thawing.
HBS positive control 1 x 1.5 ml (liquid)	C1	Ready-to-use. Serum base* supplemented with inactivated human plasma HBs antigen + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1(+) Test Value Range)
Negative control 1 x 1.9 ml (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives.
Specifications for the factory master data required to calibrate the test:		
<ul style="list-style-type: none"> <li>• MLE data (Master Lot Entry) provided in the kit,</li> </ul> or <ul style="list-style-type: none"> <li>• MLE bar code printed on the box label.</li> </ul>		
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

\* This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

**The SPR**

The interior of the SPR is coated during production with monoclonal anti-HBs Ag antibody (mouse). Each SPR is identified by the HBS code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

**The Reagent Strip**

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

**Description of the HBS strip**

Wells	Reagents
1	Sample well.
2	Conjugate: buffer containing goat serum + biotin-labeled polyclonal anti-HBs Ag antibodies (goat) + 1 g/l sodium azide (300 µl).
3 - 4 - 7 - 8 - 9	Wash buffer: buffer containing diethanolamine (DEA*: 0.85 mol/l or 9.1%) + Tween 20 + 1 g/l sodium azide (600 µl).
5	Tracer: alkaline phosphatase-labeled streptavidine + 0.9 g/l sodium azide (400 µl).
6	Prewash solution: TRIS buffer (50 mmol/l) (pH = 7.4) + protein and chemical stabilizers + 1 g/l sodium azide (600 µl).
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + DEA** (0.62 mol/l or 6.6%) pH 9.2 + 1 g/l sodium azide (300 µl).

\* Signal Word: **DANGER**

Hazard statement

H318 : Causes serious eye damage.

H373 : May cause damage to organs through prolonged or repeated exposure.

H315 : Causes skin irritation.

H302 : Harmful if swallowed.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

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H318 : Causes serious eye damage.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

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For further information, refer to the Material Safety Data Sheet.

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- Pipette with disposable tip to dispense 1 ml and 150 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- **Sample to sample contamination through contact with gloves: since high concentrations of HBs Ag may be encountered, it is strongly recommended to keep tips on a stand prior to use. It is also recommended to keep one tube of sample aside for testing HBs antigen.**
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The wash buffer (wells 3, 4, 7, 8, 9) contains a harmful agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- The optical cuvette with substrate (well 10) contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

**STORAGE CONDITIONS**

- Store the VIDAS HBs Ag Ultra kit at 2-8°C.
- **Do not freeze reagents, with the exception of S1 standard after reconstitution.**
- Store all unused reagents at 2-8°C, with the exception of S1 standard which must be kept frozen after reconstitution.
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

## SPECIMENS

### Specimen type and collection

Use sera (plain tube, tube with separator gel, plain tube with beads) or plasmas collected in lithium heparin. Serum and plasma should be stored separated from the pellet. Samples containing impurities must be clarified by centrifugation.

The results obtained were not found to be influenced by: icteric samples (bilirubin concentrations up to 500 µmol/l), hemolyzed samples (hemoglobin concentrations up to 270 µmol/l of monomer) and lipemic samples (up to 30 mg/ml).

Do not inactivate specimens.

### Specimen stability

Samples can be stored for 5 days in stoppered tubes at 2-8°C. If longer storage is required, freeze the sera or plasma at -25 ± 6°C. A study performed on samples frozen for 2 months, showed that the quality of results is not affected.

## INSTRUCTIONS FOR USE

For complete instructions, see the User's Manual.

### Reading Master lot data

Before each new lot of reagents is used, enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

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**Note: the master lot data need only be entered once for each lot.**

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User's Manual).

### Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested **in duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

### Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HBS" strip and one "HBS" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBS" code on the instrument, or by the "HBL" code for the long protocol. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.

4. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

**5. For this test, the calibrator, control, and sample test portion is 150 µl.**

6. Insert the "HBS" SPRs and "HBS" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the required temperature after pipetting.
9. The assay will be completed within approximately 60 to 90 minutes, depending on the protocol selected. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet. The patient RFV is interpreted by the VIDAS system as follows:

$$i = \text{test value} = \text{patient RFV} / \text{standard RFV}$$

The test value and interpretation are also indicated on the result sheet. Interpretation of results according to test values is as follows:

Test value		Interpretation
Short protocol	Long protocol	
$i < 0.13$	$i < 0.10$	Negative
$i \geq 0.13$	$i \geq 0.10$	Positive

**If a positive result is obtained for a patient with no previous history, the assay must be repeated and confirmed using a neutralization test (VIDAS HBs Ag Ultra Confirmation ref. 30 317) or other tests.**

A positive sample should be confirmed using the same assay protocol (short or long) as was initially used to test the sample.

Before retesting, samples should be centrifuged again so as to eliminate any interference caused by fibrin fragments or cell elements.

Interpretation of test results should be made taking into consideration the patient history, and the results of any other tests performed.

VIDAS HBs Ag Ultra is calibrated against the panel of the Société Française de Transfusion Sanguine (French Blood Transfusion Society) (adw2/ayw3 mixture expressed in ng/ml).

The analytical sensitivity of VIDAS HBs Ag Ultra, which was determined using this panel, is less than 0.20 ng/ml with the short protocol (HBS), and less than 0.15 ng/ml with the long protocol (HBL).

#### QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HBs Ag Ultra kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

Results cannot be validated if the control values deviate from the expected values.

#### Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

#### LIMITATIONS OF THE METHOD

- As interference (i.e. anti-idiotypic antibody) may be encountered with certain sera, the test should only be declared positive after taking into account the patient's history and the results of the other hepatitis B markers.
- A negative Ag HBs result does not allow infection by the hepatitis B virus to be excluded. The Ag HBs serum concentration may be below the analytical sensitivity of the reagent. The presence of a modified HBs antigen (variant) cannot be excluded; the antigen may, in this case, have been incorrectly recognized or not recognized by the antibodies in the reagent. The results of this test must be interpreted taking into consideration the patient's history and the results of any other tests performed (confirmation by neutralization, HBV DNA etc.).
- In rare cases it is possible to detect both HBs antigen and anti-HBs antibodies.
- This assay has been validated for serum and plasma and should not be used for other biological fluids such as saliva, CSF or urine.
- This assay should not be used with specimens collected post-mortem.
- Do not use serum pools.

#### RANGE OF EXPECTED VALUES (9)

The world wide prevalence of hepatitis B varies depending on the countries:

Region of the world	Prevalence of HBs Ag (%)
Europe, North America, Australia	0.2 - 0.5
Eastern Europe, Mediterranean basin, Russia, southwest Asia, South America	2 - 7
Southeast Asia, tropical Africa	8 - 20

#### PERFORMANCE

The results of the studies, which demonstrated the conformity of VIDAS HBs Ag Ultra to the Common Technical Specifications of Directive 98/79/EC, are as follows:

#### 1. Specificity for blood donor population

5059 blood donor samples from 2 blood transfusion centers were tested using the short protocol and the long protocol. No discrepancies were observed between the two protocols.

VIDAS HBs Ag Ultra (HBS/HBL)	Final interpretation of other EIA techniques	
	Positive	Negative
Positive	0	0
Negative	0	5059

Relative specificity with VIDAS HBs Ag Ultra for this population: 100.00%  
(95% confidence interval: 99.87% - 100.00%).

#### 2. Clinical specificity

##### a) for hospitalized patients:

200 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.  
(95% confidence interval: 98.04% - 100.00%).

##### b) for outpatients in a screening center:

100 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.  
(95% confidence interval: 96.38% - 100.00%).

#### 3. Analytical sensitivity

An external study performed using the panel of the French Blood Transfusion Society, showed a sensitivity of 0.12 ng/ml with the short protocol (HBS) and 0.08 ng/ml with the long protocol (HBL).

The sensitivity of VIDAS HBs Ag Ultra, determined using the international standard NIBSC 00/588, was estimated at 0.05 IU/ml with the long protocol (HBL) and 0.075 IU/ml with the short protocol (HBS).

#### 4. Diagnostic sensitivity

30 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the short protocol (HBS). 30 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the short protocol (HBS).

31 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the long protocol (HBL). 31 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the long protocol (HBL).

## 5. Detectability

A study was performed using 506 samples characterized as positive, including 28 subtyped or genotyped samples and 46 samples from patients with acute hepatitis. Detectability with the long protocol (threshold of 0.08 ng/ml) and with the short protocol (threshold of 0.12 ng/ml), is 100.00 %.

(95% confidence interval: 99.22% - 100.00%).

## 6. Sensitivity of seroconversion panels

During a study, 32 seroconversion panels were tested using the two VIDAS HBs Ag Ultra protocols.

With the short protocol, HBs Ag was detected earlier by the VIDAS HBs Ag Ultra than with the comparison method in 13 out of 32 seroconversion panels (12 panels one specimen collection earlier, 1 panel two specimen collections earlier). HBs Ag was detected in one panel, one specimen collection later.

With the long protocol, 14 of the 32 panels tested were detected earlier than with the comparison method (9 panels one specimen collection earlier, 5 panels two specimen collections earlier).

## 7. Sensitivity of HBs Ag mutants

A panel of 27 recombinant proteins imitating the most significant mutations in HBs Ag amino acid sequences and 4 native samples harboring HBs Ag mutants were tested successfully (10).

## 8. Precision

Intra-assay and inter-assay reproducibility were determined at two sites and calculated according to the recommendations of the NCCLS EP5-T2 document, volume 12-4.

### Intra-assay reproducibility

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	3.41	3.07	2.68
Weak positive	72	0.21	4.98	0.15	3.94
Negative	80	0.01	0.00 *	0.02	0.00 *

\* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

### Inter-assay reproducibility

The total precision takes into account all the sources of variability.

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	11.58	3.07	6.17
Weak positive	72	0.21	12.41	0.15	6.75
Negative	80	0.01	0.01 *	0.02	0.01 *

\* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

## 9. Cross-reactivity

367 samples from patients whose physiological status is likely to affect the detection of HBs antigen, were tested using both the short and the long protocols. All the samples were found to be negative with another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

	VIDAS HBs Ag Ultra positive samples
HCV +	0/20
EBV +	0/10
HIV +	0/10
CMV IgG +	0/10
HAV IgG +	0/10
HSV +	0/10
Syphilis	0/10
Rubella IgG +	0/10
Toxoplasmosis IgG +	0/10
Rheumatoid factor	0/9
Anti-nuclear antibodies	0/9
Non-viral hepatitis	0/5
Dialysis patients	0/10
Children less than 15 years old	0/10
Vaccinated: Anti-HBs +	0/10
Pregnant women*	0/214**

\* including 22 multipara

\*\* including 2 false positive results which were not repeated with the long protocol.

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It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardfulness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

**LITERATURE REFERENCES**

1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. *Journal of Clinical Virology* 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in *Fields Virology*, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. *JAMA*, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, *Proc Natl Acad Sci USA*, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIANTONI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, *Journal of Hepatology*, 1993, 19, 431-436.
6. ALDERSHVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, *Journal of Infectious Diseases*, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in *The Molecular Medicine of Viral Hepatitis*, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - *Journal of Medical Virology* - 2001, vol. 65, p.241-249.
9. ZUCKERMAN AJ., Hepatitis Viruses. In: Baron S, eds. *Medical Microbiology*, 4th ed. The University of Texas M Branch at Galveston, 1996: 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, *Journal of Virological Methods* 135, 2006, 109-117.

**WARRANTY**

*bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.*

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Release date	Part Number	Change Type	Change Summary
2015/01	11728I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE
2016/05	11728J	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. *Journal of Clinical Virology* 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in *Fields Virology*, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. *JAMA*, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, *Proc Natl Acad Sci USA*, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIANTONI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, *Journal of Hepatology*, 1993, 19, 431-436.
6. ALDERSHVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, *Journal of Infectious Diseases*, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in *The Molecular Medicine of Viral Hepatitis*, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - *Journal of Medical Virology* - 2001, vol. 65, p.241-249.
9. ZUCKERMAN AJ., Hepatitis Viruses. In: Baron S, eds. *Medical Microbiology*, 4th ed. The University of Texas M Branch at Galveston, 1996 : 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, *Journal of Virological Methods* 135, 2006, 109-117.

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**VIDAS<sup>®</sup> B·R·A·H·M·S PCT<sup>™</sup> (PCT)**

VIDAS<sup>®</sup> B·R·A·H·M·S PCT<sup>™</sup> is an automated test for use on the VIDAS<sup>®</sup> family of instruments for the determination of human procalcitonin in human serum or plasma (lithium heparin) using the ELFA (Enzyme-Linked Fluorescent Assay) technique.

Used in conjunction with other laboratory findings and clinical assessments, VIDAS<sup>®</sup> B·R·A·H·M·S PCT<sup>™</sup> aids in the risk assessment of critically ill patients on their first day of ICU admission, for progression to severe sepsis and septic shock.

Used in conjunction with other laboratory findings and clinical assessments, VIDAS<sup>®</sup> B·R·A·H·M·S PCT<sup>™</sup> also aids in decision making on antibiotic therapy for patients with lower respiratory tract infections (LRTI) (including community acquired pneumonia, exacerbation of chronic obstructive pulmonary disease, acute bronchitis) seen during medical consultations, including at the Emergency Department.

**SUMMARY AND EXPLANATION**

Procalcitonin (PCT) is the prohormone of calcitonin. Whereas calcitonin is only produced in the C cells of the thyroid gland as a result of hormonal stimulus, PCT is secreted by different types of cells from numerous organs in response to proinflammatory stimulation, particularly bacterial stimulation (1).

Depending on the clinical background, a PCT concentration above 0.1 ng/mL can indicate clinically relevant bacterial infection, requiring antibiotic treatment (2). At a PCT concentration > 0.5 ng/mL, a patient should be considered at risk of developing severe sepsis or septic shock (3, 4).

Sepsis is an excessive reaction of the immune system and coagulation system to an infection (5). The diagnosis and monitoring of infected patients are major problems for physicians. It has been proven that PCT levels increase precociously, specifically in patients with a bacterial infection. For laboratory diagnosis, PCT is therefore an important marker enabling specific differentiation between a bacterial infection and other causes of inflammatory reactions (2). Moreover, the resorption of the septic infection is accompanied by a decrease in the PCT concentration which returns to normal with a half-life of 24 hours (6, 7).

Several randomized clinical trials have shown that the use of procalcitonin to guide the initiation and also the duration of antibiotic treatment in patients with LRTI, significantly reduced antibiotic consumption across different LRTI diagnoses (8, 9). The randomized trials were conducted in different settings, including primary care (10), the Emergency Department (11) or the ICU (12). In the ED or the primary care setting, a cut-off concentration at 0.25 ng/ml was used to either withhold or stop antibiotics in LRTI patients (13, 9, 14).

The reduction of antibiotic consumption was clinically safe in this framework, as no higher rates of mortality or treatment failure were associated with a PCT-guided antibiotic therapy (9).

The safe reduction of antibiotic use through PCT-guided therapy was confirmed in an observational quality surveillance study that considered consecutive LRTI patients, recruited without exclusion criteria, who were seen at an emergency department or a physician's office (14).

In certain situations (newborns, polytrauma, burns, major surgery, prolonged or severe cardiogenic shock, etc.) PCT elevation may be independent of any infectious aggression. The return to normal values is usually rapid. Viral infections, allergies, autoimmune diseases and graft rejection do not lead to a significant increase in PCT (15). A localized bacterial infection can lead to a moderate increase in PCT levels (2, 16).

The evaluation of VIDAS<sup>®</sup> B·R·A·H·M·S PCT<sup>™</sup> assay results must always be performed taking into consideration the patient's history and the results of any other tests performed.

If discrepancies are found between the laboratory findings and the clinical signs, additional tests should be performed.

**PRINCIPLE**

The assay principle combines a one-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR<sup>®</sup>) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument. The sample is transferred into the wells containing anti-procalcitonin antibodies labeled with alkaline phosphatase (conjugate). The sample/conjugate mixture is cycled in and out of the SPR several times. This operation enables the antigen to bind with the immunoglobulins fixed to the interior wall of the SPR and the conjugate to form a sandwich. Unbound components are eliminated during washing steps.

Two detection steps are then performed successively. During each step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR<sup>®</sup>. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample. At the end of the assay, the results are automatically calculated by the instrument in relation to two calibration curves stored in memory corresponding to the two detection steps. A fluorescence threshold value determines the calibration curve to be used for each sample. The results are then printed out.

**CONTENT OF THE KIT - RECONSTITUTION OF REAGENTS (60 TESTS):**

60 PCT Strips	STR	Ready-to-use.
60 PCT SPR®s 2 x 30	SPR®	Ready-to-use. Interior of SPR®s coated with mouse monoclonal anti-human procalcitonin immunoglobulins.
PCT Controls C1 Control 2 x 2 mL (lyophilized) C2 Control 2 x 2 mL (lyophilized)	C1  C2	Reconstitute with 2 mL of distilled water. Wait for 5 to 10 minutes and then mix. Stable after reconstitution for 8 hours at 2-8°C, or until the expiration date on the kit at - 25 ± 6°C. 5 freeze/thaw cycles are possible. TRIS NaCl buffer (pH 7.3) + recombinant human PCT + preservatives. MLE data indicate the acceptable range in ng/mL ("Control C1 Dose Value Range" or "Control C2 Dose Value Range").
PCT Calibrators S1 Calibrator 2 x 2 mL (lyophilized) S2 Calibrator 2 x 2 mL (lyophilized)	S1  S2	Reconstitute with 2 mL of distilled water. Wait for 5 to 10 minutes and then mix. Stable after reconstitution for 8 hours at 2-8°C, or until the expiration date on the kit at - 25 ± 6°C. 5 freeze/thaw cycles are possible. TRIS NaCl buffer (pH 7.3) + recombinant human PCT + preservatives. MLE data indicate the concentration in ng/mL ("Calibrator (S1) Dose Value" or "Calibrator (S2) Dose Value") and the acceptable range in "Relative Fluorescence Value" (Calibrator (S1) RFV Range or Calibrator (S2) RFV Range).
Specifications for the factory master data required to calibrate the test:		
<ul style="list-style-type: none"> <li>MLE data (Master Lot Entry) provided in the kit. or</li> <li>MLE bar codes printed on the box label.</li> </ul>		
1 package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a> .		

**The SPR**

The interior of the SPR® is coated during production with mouse monoclonal anti-human procalcitonin immunoglobulins. Each SPR is identified by the PCT code. Only remove the required number of SPR®s from the pouch and **carefully reseal the pouch after opening**.

**The Reagent Strip**

The strip consists of 10 wells covered with a labeled foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

**Description of the PTC strip**

Well	Reagents
1	Sample well.
2 - 3 - 4	Empty wells.
5	Conjugate: alkaline phosphatase-labeled mouse monoclonal anti-human procalcitonin immunoglobulins + preservative (400 µL).
6 - 7 - 8	TRIS NaCl Tween (pH 7.3) + preservative (600 µL).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/L) + diethanolamine* (DEA*) (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

\* Signal Word: **DANGER**

**Hazard statement**

H318: Causes serious eye damage.

**Precautionary statement**

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information, please refer to the Material Safety Data Sheet.

**MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED**

- Pipette with disposable tip to dispense 2 mL and 200 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instrument of the VIDAS® family.

**ADDITIONAL REAGENT**

- Serum free (ref. 66 581)

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- The kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPR®s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the kit label.
- Do not mix reagents (or disposables) from different lots.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- The substrate in well 10 contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" indicated above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

**STORAGE CONDITIONS**

- Store the VIDAS® B·R·A·H·M·S PCT™ kit at 2-8°C.
- **Do not freeze reagents, with the exception of calibrators and controls after reconstitution.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR® pouch is correctly sealed and undamaged. If not, do not use the SPR®s.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPR®s and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

**SPECIMENS****Specimen type and collection**

Human serum or plasma (lithium heparin).

Since EDTA causes a decrease in the values measured, **plasma collected on EDTA should not be used.**

**For a given patient, the PCT assays must be performed on the same type of sample tube.**

**Specimen Type**

- **Dry tubes:** wait for samples to coagulate and **centrifuge** according to the tube manufacturer's recommendations to eliminate fibrin.
- **Other tubes:** follow the tube manufacturer's recommendations for use.
- **Frozen-stored samples:** after thawing, these samples must be clarified by centrifuging before testing.

**Note:** Blood sampling tube results may vary from one manufacturer to another depending on the materials and additives used.

It is the responsibility of each laboratory to validate the type of sample tube used and to follow the manufacturer's recommendations for use.

**Sample preparation**

Follow the tube manufacturer's recommendations for use.

The pre-analytical step, including the preparation of blood samples, is an essential first step when performing medical analyses. In accordance with Good Laboratory Practice, this step is performed under the responsibility of the laboratory manager.

Insufficient clot time can result in the formation of fibrin with micro-clots that are invisible to the naked eye. The presence of fibrin, red blood cells, or suspended particles can lead to erroneous results.

Samples containing suspended fibrin particles or erythrocyte stroma should be centrifuged before testing.

**For serum specimens, ensure that complete clot formation has taken place prior to centrifugation.**

Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting times.

**Sample stability**

Samples separated from the clot can be stored at 2-8°C in stoppered tubes for up to 48 hours; if longer storage is required, freeze the sera or plasma at  $-25 \pm 6^\circ\text{C}$ . Six-month storage of frozen samples does not affect the quality of results. Three freeze/thaw cycles have been validated.

**Test sample volumes less than 200 µL**

Test sample volumes between 50 µL and 200 µL can be tested after performing a manual dilution up to 1/4 (1 volume of test sample + 3 volumes of Serum Free reagent (ref. 66 581)) and no more than two hours after dilution.

**Sample-related interferences**

None of the following factors have been found to significantly influence this assay:

- hemolysis (after spiking samples with hemoglobin: up to 347 µmol/L (monomer)),
- lipemia (after spiking samples with lipids, up to 30 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: up to 574 µmol/L).

However, it is recommended not to use samples that are clearly hemolyzed, lipemic or icteric and, if possible, to collect a new sample.

## INSTRUCTIONS FOR USE

For complete instructions, see the Instrument User Manual.

Reading VIDAS® Protocol Test Change (PTC) data and MLE data

### When using the assay for the first time:

With the external instrument barcode reader,

1. Scan the PTC barcode(s) at the end of the package insert or downloadable from [www.biomerieux.com/techlib](http://www.biomerieux.com/techlib). This reading allows VIDAS® PTC protocol data to be transferred to the instrument software for its update.
2. Scan the MLE data on the box label.

**Note: If the MLE data have been read before the VIDAS® PTC protocol, read the MLE data again.**

### When opening a new lot of reagents:

Enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results.

**Note: the master lot data need only to be entered once for each lot.**

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the VIDAS® User Manual).

### Calibration

Calibration, using the **two calibrators** provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered, and then **every 28 days**. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

**The calibrators identified by S1 and S2**, must be tested **in duplicate** (see User Manual) in the same run. The calibration values must be within the set RFV (Relative Fluorescence Value). If this is not the case, **recalibrate using S1 and S2**.

### Procedure

1. **Remove the required reagents from the refrigerator.**
2. Use one "PCT" strip and one "PCT" SPR® for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPR®s have been removed.**
3. The test is identified by the "PCT" code on the instrument. The calibrators must be identified by "S1" and "S2", and tested **in duplicate**. If the controls need to be tested, they should be identified by "C1" and "C2" and tested singly.

4. Mix the calibrators and controls using a vortex-type mixer.
5. To obtain optimum results, refer to all the paragraphs in the **SPECIMENS** section.
6. Before pipetting, ensure that samples, calibrators and controls are free of bubbles.

**7. For this test, the test portion of the calibrators, controls and samples is 200 µL.**

8. Insert the SPR®s and the strips into the instrument. Check to make sure the color labels with the assay code on the SPR®s and the Reagent Strips match.
9. **Initiate the assay immediately.** All the assay steps are performed automatically by the instrument.
10. Reclose the vials and return them to the required temperature after pipetting.
11. The assay will be completed within approximately **20 minutes**. After the assay is completed, remove the SPR®s and strips from the instrument.
12. Dispose of the used SPR®s and strips into an appropriate recipient.

## RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer using two calibration curves which are stored by the instrument; the concentrations are expressed in ng/mL.

As no international standard is available, VIDAS® B·R·A·H·M·S PCT™ is calibrated against an internal panel of human sera with known procalcitonin concentrations. In case of patient follow-up, it is recommended to use the same PCT assay technique.

Samples with procalcitonin concentrations > 200 ng/mL should be retested after dilution by 1/10 (1 volume of sample + 9 volumes of PCT negative sample or Serum Free reagent (ref. 66 581)). The final result should take into account the original dilution factor.

### **Special case of a sample volume < 200 µL:**

If the result obtained after dilution is below the measuring range for the test (0.05 ng/mL) then the VIDAS® instrument will report an "INVALID" result with the mention "over-diluted". The final result cannot be calculated and should be reported as **less than (0.05 ng/mL x dilution factor)**. In this case, a sample which is tested with a dilution at 1/4 should be reported as **less than 0.2 ng/mL** ( $0.2 = 0.05 \times 4$ ).

If the dilution factor has not been entered when the Work List was created (see User Manual), multiply the result by the dilution factor to obtain the sample concentration.

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed.

## QUALITY CONTROL

Two controls are included in each VIDAS® B·R·A·H·M·S PCT™ kit. These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2. Results cannot be validated if the control values deviate from the expected values.

### Note

It is the responsibility of the user to perform Quality Control in accordance with any applicable local regulations.

## LIMITATIONS OF THE METHOD

Interference may be encountered with certain samples containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

## RANGE OF EXPECTED VALUES

### - Risk assessment for progression to severe sepsis and septic shock

In agreement with the literature (3, 4), the results obtained with VIDAS® B·R·A·H·M·S PCT™ during a study performed on patients admitted to intensive care units (refer to the "Clinical performance" section) are as follows:

- a concentration < 0.5 ng/mL represents a low risk of severe sepsis and/or septic shock.
- a concentration > 2 ng/mL represents a high risk of severe sepsis and/or septic shock.

Nevertheless, concentrations < 0.5 ng/mL do not exclude an infection, on account of localized infections (without systemic signs) which can be associated with such low concentrations, or a systemic infection in its initial stages (< 6 hours). Furthermore, increased procalcitonin can occur without infection. PCT concentrations between 0.5 and 2.0 ng/mL should be interpreted taking into account the patient's history. It is recommended to retest PCT within 6-24 hours if any concentrations < 2 ng/mL are obtained.

### - Decision making on antibiotic therapy for patients with lower respiratory tract infections

In agreement with the literature (8, 9, 10, 11, 12, 13), and validated using VIDAS® B·R·A·H·M·S PCT™ (14), recommended cut-off values are as follows:

PCT concentration	Analysis / Recommendation	Note
< 0.10 ng/mL	Indicates absence of bacterial infection. Antibiotic therapy strongly discouraged.	<ul style="list-style-type: none"> <li>• If antibiotics are withheld, repeat PCT measurement within 6-24 hours (also in outpatients if symptoms persist/worsen).</li> </ul>
0.10-0.25 ng/mL	Bacterial infection unlikely. Antibiotic therapy discouraged.	<ul style="list-style-type: none"> <li>• Antibiotic therapy should be considered for:               <ul style="list-style-type: none"> <li>○ Respiratory or hemodynamic instability, severe comorbidities, ICU admission</li> <li>○ <b>PCT &lt;0.1 ng/mL:</b> CAP with PSI V or CURB-65&gt;3, COPD with GOLD IV</li> <li>○ <b>PCT 0.1-0.25 ng/mL:</b> CAP with PSI IV &amp; V or CURB-65&gt;2, COPD with GOLD III &amp; IV</li> </ul> </li> </ul>
0.26-0.50 ng/mL	Bacterial infection is possible. Antibiotic therapy recommended.	<ul style="list-style-type: none"> <li>• Follow up samples can be tested at regular intervals and antibiotic therapy may be discontinued using the same cut-off values per this table.</li> <li>• If PCT remains high, consider treatment failure</li> </ul>
> 0.50 ng/mL	Suggestive of presence of bacterial infection. Antibiotic therapy strongly recommended.	

- CAP: **C**ommunity-**A**cquired **P**neumonia
- PSI: **P**neumonia **S**everity **I**ndex
- CURB-65: **C**onfusion – **U**rea – **R**espiratory Rate – **B**lood pressure – age > **65** years
- COPD: **C**hronic **O**bstuctive **P**ulmonary **D**isease
- GOLD: **G**lobal initiative for chronic **O**bstuctive **L**ung **D**isease

**PERFORMANCE**

Studies performed using VIDAS® B·R·A·H·M·S PCT™ gave the following results:

**Measurement range**

The VIDAS® B·R·A·H·M·S PCT™ assay measurement range is 0.05-200 ng/mL.

**Detection limits**

The Limit of Blank (LoB), the Limit of Detection (LoD) and the Limit of Quantitation (LoQ) were determined according to the CLSI® EP17-A2 recommendations:

Limit of Blank (LoB)	0.01 ng /mL
Limit of Detection (LoD)	0.03 ng /mL
Limit of Quantitation (LoQ)	0.05 ng /mL

The Limit of Quantitation (LoQ) is the lowest concentration of PCT measured with a within-site precision of 20% CV.

**Hook effect**

No hook effect was found up to procalcitonin concentrations of 2 600 ng/mL.

**Normal values**

These results are given as a guide; it is recommended that each laboratory establish its own reference values from a rigorously selected population.

A study was performed using the VIDAS® B·R·A·H·M·S PCT™ test on serum samples from apparently healthy male (N=98) and female (N=102) individuals. The normal values corresponding to the 95<sup>th</sup> and 99<sup>th</sup> percentiles were respectively found at < 0.05 ng/mL and at 0.09 ng/mL.

**Precision**

A precision study was performed according to the recommendations of CLSI® document EP5-A3.

Nine samples were tested in duplicate in 2 runs per day, over 20 days using 3 VIDAS® instruments installed in 3 laboratories (N=240 values for each sample).

Two reagent lots were used: 10 days of tests and 2 calibrations were performed for each lot (5 test days per calibration).

The repeatability, within-laboratory precision and reproducibility/total precision (between-laboratory precision) were estimated for each sample and are reported in the following table:

Sample	N	Mean Concentration (ng/mL)	Repeatability		Within-laboratory Precision		Reproducibility/Total precision	
			Standard deviation (ng/mL)	CV (%)	Standard deviation (ng/mL)	CV (%)	Standard deviation (ng/mL)	CV (%)
Sample 1	240	0.12	0.011	9.4%	0.017	14.7%	0.017	14.7%
Sample 2	240	0.15	0.010	6.4%	0.019	12.5%	0.019	12.5%
Sample 3	240	0.20	0.010	5.1%	0.017	8.3%	0.017	8.3%
Sample 4	240	0.53	0.013	2.4%	0.023	4.3%	0.023	4.3%
Sample 5	240	2.12	0.027	1.3%	0.079	3.7%	0.079	3.7%
Sample 6	240	23.09	0.501	2.2%	1.067	4.6%	1.067	4.6%
Sample 7	240	92.17	3.099	3.4%	7.252	7.9%	7.252	7.9%
Sample 8	240	128.37	5.155	4.0%	12.578	9.8%	12.578	9.8%
Sample 9	240	162.79	7.233	4.4%	18.934	11.6%	18.934	11.6%

**Specificity**

The following compounds, tested at the concentrations indicated in the table, do not affect the VIDAS® B·R·A·H·M·S PCT™ assay.

Tested compound	Tested concentration
Protein (albumin)	4 g/dL
Human Calcitonin	60 ng/mL
Human Katalcalcin	10 ng/mL
Human a-CGRP*	10 µg/mL
Human b-CGRP*	10 µg/mL

\*Calcitonin Gene Related Peptide

**Drug interference**

The following drugs, at the concentrations indicated in the table, do not affect the VIDAS® B•R•A•H•M•S PCT™ assay:

<b>Tested drug</b>	<b>Tested concentration</b>
Acetaminophen (paracetamol)	1 324 µmol/L
Acetylsalicylic Acid	3.62 mmol/L
Alcohol	86.8 mmol/L
Amoxicillin	206 µmol/L
Ampicillin	152 µmol/L
Azithromycin	15.3 µmol/L
Beclometasone dipropionate	1.00 µg/mL
Caffeine	308 µmol/L
Cefotaxime	673 µmol/L
Ceftriaxone	1 416 µmol/L
Celecoxib	240 µg/mL
Cetirizine HCl	7.71 µmol/L
Cromolyn	24 µg/mL
Dextromethorphan	3.70 µmol/L
Dopamine	5.87 µmol/L
Dobutamine	1 500 ng/mL
Epinephrine (adrenaline)	1.8 µg/mL
Fluticasone	0.30 µg/mL
Formoterol	29 ng/mL
Furosemide	181 µmol/L
Heparin	3 000 IU/mL
Ibuprofen	2 425 µmol/L
Imipenem	180 µg/mL
Levofloxacin	48.6 µmol/L
Linezolid	480 µg/mL
Loratadine	0.78 µmol/L
Naproxen	2 170 µmol/L
Nicotine	6.2 µmol/L
Noradrenaline	2.1 ng/mL
Oxymetazoline HCl	90 ng/mL
Phenylephrine	180 ng/mL
Prednisolone	8.31 µmol/L
Salmeterol	60 ng/mL
Theophylline	222 µmol/L
Tiotropium	22 ng/mL
Vancomycin	69 µmol/L

## Accuracy

The test linearity was studied according to a procedure taken from the CLSI® EP6-A guideline. The test is linear over the complete measurement range.

Three samples were diluted in a PCT-negative serum pool and tested in triplicate. The ratio of the mean concentration measured over the expected mean concentration is expressed as a mean recovery percentage.

Samples	Dilution factor	Expected mean concentration (ng/mL)	Measured mean concentration (ng/mL)	Mean recovery percentage (%)
1	1/1	137.07	137.07	100.0
	1/2	68.54	71.54	104.4
	1/3	45.69	49.56	108.5
	1/4	34.27	37.26	108.7
	1/8	17.13	19.50	113.8
	1/16	8.57	8.75	102.2
	1/20	6.85	7.73	112.8
2	1/1	38.67	38.67	100.0
	1/2	19.34	19.75	102.1
	1/3	12.89	13.90	107.8
	1/4	9.67	9.79	101.3
	1/8	4.83	4.96	102.5
	1/16	2.42	2.26	93.3
	1/20	1.93	1.85	95.7
3	1/1	7.58	7.58	100.0
	1/2	3.79	4.17	110.1
	1/3	2.53	2.70	107.0
	1/4	1.90	1.98	104.7
	1/8	0.95	0.94	99.2
	1/16	0.47	0.51	108.4
	1/20	0.38	0.37	98.5

### Concordance with the B•R•A•H•M•S PCT LIA method

A concordance study between VIDAS® B•R•A•H•M•S PCT™ and B•R•A•H•M•S PCT LIA was performed using 204 samples with cut-off values at 0.5 ng/mL and 2 ng/mL.

VIDAS® B•R•A•H•M•S PCT™	B•R•A•H•M•S PCT LIA		
	≤ 0.5 ng/mL	> 0.5 ng/mL	Total
≤ 0.5 ng/mL	74	1	75
> 0.5 ng/mL	5	124	129
<b>Total</b>	<b>79</b>	<b>125</b>	<b>204</b>

VIDAS® B•R•A•H•M•S PCT™	B•R•A•H•M•S PCT LIA		
	≤ 2 ng/mL	> 2 ng/mL	Total
≤ 2 ng/mL	109	4	113
> 2 ng/mL	8	83	91
<b>Total</b>	<b>117</b>	<b>87</b>	<b>204</b>

The percentages of concordance between the 2 techniques for the cut-off values at 0.5 and 2 ng/mL are respectively 97.1% and 94.1%.

## Clinical Performance

### • Risk assessment for progression to severe sepsis and septic shock

A study performed at four (4) sites (2 in France and 2 in the USA) determined the clinical performance of the VIDAS® B•R•A•H•M•S PCT™ assay. This study included 229 patients (141 males and 88 females), who were consecutively admitted to the medical intensive care unit (MICU). The data represents first day admission testing. Patients admitted for trauma, surgery, burns, or prolonged or severe cardiogenic shock were excluded from the study.

Based on criteria from the consensus conference of the American College of Chest Physicians/Society of Critical Care Medicine (5), patients were classified into 5 categories: no infection, SIRS (Systemic Inflammatory Response Syndrome), sepsis, severe sepsis and septic shock. The classification was reviewed by an independent expert.

The number, range and mean age in each category were as follows:

- no infection: 27 patients aged between 22 and 92 years (mean 64.4 years)
- SIRS: 62 patients aged between 18 and 87 years (mean 59.0 years)
- sepsis: 42 patients aged between 21 and 92 years (mean 64.2 years)
- severe sepsis: 48 patients aged between 19 and 89 years (mean 66.3 years)
- septic shock: 50 patients aged between 33 and 88 years (mean 68.2 years)

PCT values for the groups of patients with no infection or SIRS or sepsis versus severe sepsis or septic shock with cut-off values at 0.5 ng/mL and 2.0 ng/mL are shown in the tables below:

- Results obtained with a cut-off value at 0.5 ng/mL:

	No infection/SIRS/sepsis	Severe sepsis/septic shock	Total
PCT ≤ 0.5 ng/mL	88	3	91
PCT > 0.5 ng/mL	43	95	138
Total	131	98	229

- Results obtained with a cut-off value at 2 ng/mL:

	No infection/SIRS/sepsis	Severe sepsis/septic shock	Total
PCT ≤ 2 ng/mL	115	19	134
PCT > 2 ng/mL	16	79	95
Total	131	98	229

#### • Decision making on antibiotic therapy for patients with lower respiratory tract infections

A clinical study (14) established clinical performance data for the VIDAS® B.R.A.H.M.S PCT™ assay in terms of decision making on antibiotic therapy for patients with lower respiratory tract infections.

This study shows that the **duration of antibiotic therapy is significantly reduced** (from 7.4 to 5.9 days) compared to medical care that does not include the evaluation of PCT concentrations (difference of -1.51 days; 95% confidence interval [-2.04 ; -0.98]; p < 0.01). The study population was composed of patients with community-acquired pneumonia, exacerbation of chronic obstructive pulmonary disease or acute bronchitis.

Furthermore, **no increase in adverse medical outcomes** (relapse, hospitalization, side effects of antibiotic therapy, and/or mortality) is associated with withholding antibiotic therapy for patients with low PCT values (≤ 0.25 ng/mL) on hospital admission, and with stopping antibiotic therapy.

#### WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

**LITERATURE REFERENCES**

1. DANDONA P , NIX D, et al. Procalcitonin increase after endotoxin injection in normal subjects, *J Clin Endocrinol Metab.* 1994; 79: 1605-1608.
2. CHRIST-CRAIN M, JACCARD-STOLZ D, et al. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections : cluster-randomised single-blinded intervention trial. *Lancet* 2004; 363: 600-607.
3. MULLER B, BECKER KL, et al. Calcitonin precursors are reliable markers of sepsis in medical intensive care unit. *Crit. Care Med.* 2000;28: 977-983.
4. HARBARTH S, HOLECKOVA K, et al. Diagnostic value of procalcitonin, interleukin-6 and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med.* 2001;164: 396-402.
5. American College of Chest Physicians/Society of Critical Care Medicine: Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. *Crit Care Med* 1992; 20: 864-874.
6. LUYT CE, GUERIN V, et al. Procalcitonin kinetics as a prognostic marker of ventilator-associated pneumonia. *Am J Respir Crit Care Med* 2005;171: 48-53.
7. BRUNKHORST FM, HEINZ U, et al. Kinetics of procalcitonin in iatrogenic sepsis. *Intensive Care Med.* 1998; 24:888-892.
8. SCHUETZ P, CHRIST-CRAIN M, et al. Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: the ProHOSP randomized controlled trial. *JAMA* 2009;302:1059-1066.
9. SCHUETZ P, MULLER B, et al. Procalcitonin to initiate or discontinue antibiotics in acute respiratory tract infections. *Cochrane Database Syst Rev.* 2012;9:CD007498.
10. BURKHARDT O, EWIG S, et al. Procalcitonin guidance and reduction of antibiotic use in acute respiratory tract infection. *Eur Respir J.* 2010;36:601-7.
11. CHRIST-CRAIN M, JACCARD-STOLZ D, et al. Effect of Procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. *Lancet* 2004;363:600-7.
12. BOUADMA L, LUYT CE, et al. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): a multicentre randomised controlled trial. *Lancet* 2010; 375:463-474.
13. STOLZ D, CHRIST-CRAIN M, et al. Diagnostic value of signs, symptoms and laboratory values in lower respiratory tract infection. *SwissMedWkly.* 2006;136:434-440.
14. ALBRICH WC, DUSEMUND F, et al. Effectiveness and Safety of Procalcitonin-Guided Antibiotic Therapy in Lower Respiratory Tract Infections in "Real Life". *Arch Intern Med.* 2012;172:715-722.
15. MEISNER M. Procalcitonin (PCT) – A new, innovative infection parameter. Biochemical and clinical aspects. Thieme Stuttgart, New York 2000, ISBN: 3-13-105503-0
16. CHRIST-CRAIN M, MULLER B. Procalcitonin in bacterial infections--hype, hope, more or less? *Swiss Med Wkly* 2005;135:451-460.
17. SCHUETZ P, ALBRICH W, et al. Procalcitonin for guidance of antibiotic therapy. *Expert Rev Anti Infect Ther* 2010;8:575-587.

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalog number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

**LIMITED WARRANTY**

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the "System") other than as set forth in the IFU.

**REVISION HISTORY**

Change type categories:

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
<b>Note:</b>	<i>Minor typographical, grammar, and formatting changes are not included in the revision history</i>

Release date	Part Number	Change Type	Change Summary
2015/01	13207F	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical change	CONTENT OF THE KIT - RECONSTITUTION OF REAGENTS (60 TESTS); WARNINGS AND PRECAUTIONS
2015/10	13207G	Technical change	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2016/11	13207H	Technical change	SUMMARY AND EXPLANATION CONTENT OF THE KIT - RECONSTITUTION OF REAGENTS (60 TESTS); MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED SPECIMENS INSTRUCTIONS FOR USE RESULTS AND INTERPRETATION QUALITY CONTROL RANGE OF EXPECTED VALUES PERFORMANCE LITERATURE REFERENCES LIMITED WARRANTY

*Vidas* B•R•A•H•M•S PCT


Reagent developed in collaboration with B•R•A•H•M•S

BIOMERIEUX, the blue logo, VIDAS and SPR are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

B•R•A•H•M•S PCT™ is the property of Thermo Fisher Scientific Inc and its subsidiaries.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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1. DANDONA P, NIX D, et al. Procalcitonin increase after endotoxin injection in normal subjects, *J Clin Endocrinol Metab.* 1994; 79: 1605-1608.
2. CHRIST-CRAIN M, JACCARD-STOLZ D, et al. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections : cluster-randomised single-blinded intervention trial. *Lancet* 2004; 363: 600-607.
3. MULLER B, BECKER KL, et al. Calcitonin precursors are reliable markers of sepsis in medical intensive care unit. *Krit. Care Med.* 2000;28: 977-983.
4. HARBARTH S, HOLECKOVA K, et al. Diagnostic value of procalcitonin, interleukin-6 and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med.* 2001;164: 396-402.
5. American College of Chest Physicians/Society of Critical Care Medicine: Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. *Crit Care Med* 1992; 20: 864-874.
6. LUYT CE, GUERIN V, et al. Procalcitonin kinetics as a prognostic marker of ventilator-associated pneumonia. *Am J Respir Crit Care Med* 2005;171: 48-53.
7. BRUNKHORST FM, HEINZ U, et al. Kinetics of procalcitonin in iatrogenic sepsis. *Intensive Care Med.* 1998; 24:888-892.
8. SCHUETZ P, CHRIST-CRAIN M, et al. Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: the ProHOSP randomized controlled trial. *JAMA* 2009;302:1059-1066.
9. SCHUETZ P, MULLER B, et al. Procalcitonin to initiate or discontinue antibiotics in acute respiratory tract infections. *Cochrane Database Syst Rev.* 2012;9:CD007498.
10. BURKHARDT O, EWIG S, et al. Procalcitonin guidance and reduction of antibiotic use in acute respiratory tract infection. *Eur Respir J.* 2010;36:601-7.
11. CHRIST-CRAIN M, JACCARD-STOLZ D, et al. Effect of Procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. *Lancet* 2004;363:600-7.
12. BOUADMA L, LUYT CE, et al. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): a multicentre randomised controlled trial. *Lancet* 2010; 375:463-474.
13. STOLZ D, CHRIST-CRAIN M, et al. Diagnostic value of signs, symptoms and laboratory values in lower respiratory tract infection. *SwissMedWkly.* 2006;136:434-440.
14. ALBRICH WC, DUSEMUND F, et al. Effectiveness and Safety of Procalcitonin-Guided Antibiotic Therapy in Lower Respiratory Tract Infections in "Real Life". *Arch Intern Med.* 2012;172:715-722.
15. MEISNER M. Procalcitonin (PCT) – A new, innovative infection parameter. Biochemical and clinical aspects. Thieme Stuttgart, New York 2000, ISBN: 3-13-105503-0
16. CHRIST-CRAIN M, MULLER B. Procalcitonin in bacterial infections--hype, hope, more or less? *Swiss Med Wkly* 2005;135:451-460.
17. SCHUETZ P, ALBRICH W, et al. Procalcitonin for guidance of antibiotic therapy. *Expert Rev Anti Infect Ther* 2010;8:575-587.

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
**Quality Control VIDAS<sup>®</sup> (QCV)**

Automated test for use on the VIDAS<sup>®</sup> system to detect abnormal operation of the VIDAS and mini VIDAS instrument pipette mechanisms and optical systems.

**SUMMARY AND EXPLANATION**

Quality Control VIDAS (QCV) is used to detect abnormal operation of pipette mechanisms which may affect the results of biological tests. It is also intended for checking that the optical system is capable of measuring high fluorescence levels.

**The QCV test must be run at least once a month in each position** of the VIDAS and mini VIDAS instruments, or any time that one of the above-mentioned instrument problems is suspected.

**PRINCIPLE**

The test corresponds to successive aspirations/dilutions of fluorescent substrate (4-Methyl-umbelliferone) solutions which are standardized and have varying levels of concentration. The aspirations are performed at different speeds to check pump aspiration capability.

The results obtained at the end of the test correspond to:

- a fluorescence ratio "Test Value 1" (TV1 for mini VIDAS) for checking that the pipette mechanism is operating correctly,
- a fluorescence reading (Reading 3 for VIDAS and R3 for mini VIDAS) for checking that the optical system is intact.

The values obtained must be within the acceptable range defined later in this package insert.

**CONTENT OF THE KIT (60 TESTS):**

60 QCV Strips	STR	Ready-to-use.
60 QCV SPRs (2 x 30)	SPR	Ready-to-use.
1 Package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com/techlib">www.biomerieux.com/techlib</a>		

**DESCRIPTION****The SPR<sup>®</sup>:**

The SPR is used as a pipetting device and is identified by the code QCV. **Only remove the required number of SPRs from the pouch.**

**The strip:**

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The last well is a cuvette in which the fluorometric reading is performed. The wells in the center of the strip contain the various reagents required for the assay.

**Description of the QCV strip:**

Wells	Reagents
1	<b>Empty well: no sample is needed for this test.</b>
2	400 µL of 4-Methyl-umbelliferone solution (100 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
3	Empty well.
4	400 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
5 - 6	600 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
7	Empty well.
8 - 9	400 µL of 4-Methyl-umbelliferone solution (26 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
10	Optical cuvette with substrate: 300 µL of 4-Methyl-umbelliferone (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.

**WARNINGS AND PRECAUTIONS**

1. **Caution: the strips are photosensitive and must be kept in their box to be protected from light.**
2. Do not use reagents after the expiration date indicated on the label.
3. Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.

4. It is recommended to wear gloves as if handling biological tests.

**STORAGE**

- Store the VIDAS QCV kit at 2 - 8°C. Do not freeze reagents. **Store all unused reagents in their box at 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Do not use the reagents after the expiration date.

## INSTRUCTIONS FOR USE

Your instrument should be checked **at least once a month** using the VIDAS QCV test, and each time abnormal operation of the pipette mechanism or optical system is suspected.

**It is essential to perform the VIDAS QCV test on all the positions to ensure that the whole instrument has been checked.**

It is essential to archive patient work lists and their results on a monthly basis in order to track results obtained between two QCV tests.

1. Remove necessary components from the kit and return all unused components to storage at 2-8°C.

2. Create the work list: type QCV and indicate the number of tests to be run.
3. Load the reagent strips and SPRs into the instrument (**Note:** no desiccant is needed in the SPR pouch). Be sure to close all SPR® doors and tray covers.
4. Run the test (see User's Manual). All steps will be performed automatically by the instrument. The test will be completed in approximately 20 minutes.
5. **After the test is completed, print the results and analyze them as instructed later in this package insert. Once interpretation of results has been validated, dispose of the used SPRs and strips into an appropriate recipient.**

## RESULTS AND INTERPRETATION

The first reading is performed at the beginning of the assay to check that the reagent is intact:

reading 1 for VIDAS  
R1 for mini VIDAS

The second reading is taken after successive aspirations which are performed at different speeds; they are used to check the ability of the pump to pipette properly:

reading 2 for VIDAS  
R2 for mini VIDAS

The third reading is taken after a small quantity of highly concentrated substrate is pipetted. This operation is performed to check that the optical system is capable of measuring high fluorescence levels:

reading 3 for VIDAS  
R3 for mini VIDAS


When the test is completed, the computer automatically calculates Test Value 1 and 2 (**TV1** and **TV2** for mini VIDAS) using the 3 fluorescence readings. The calculations appear on a report.

$$\text{VIDAS: Test Value (TV1)} = \frac{\text{reading 2}}{\text{reading 1}} \quad \text{Test Value (TV2)} = \frac{\text{reading 3}}{\text{reading 1}}$$

$$\text{mini VIDAS: TV1} = \frac{R2}{R1} \quad \text{TV2} = \frac{R3}{R1}$$

The values for reading 3 (R3 for mini VIDAS) and Test Value 1 (TV1 for mini VIDAS) must be within the acceptable ranges indicated below (see diagram):

- **Acceptable range for Test Value 1 (TV1 for mini VIDAS):** the Test Value 1 for each position must be  $\geq$  the value indicated on the kit label. If the result of a particular position is outside the range, two new VIDAS QCV tests must be run successively in all the positions in the section concerned (for VIDAS and mini VIDAS). If at least one other non-compliant result is produced in the same section, independent of the position, take the section offline and contact bioMérieux.
- **Acceptable range for reading 3 (R3 for mini VIDAS):** reading 3 for each position must be  $\geq$  4100 RFU. If the result of a particular position is outside the range, take the instrument offline and contact bioMérieux.

**LIMITATIONS OF THE METHOD**

- VIDAS QCV may produce an invalid result when the reading 1 or R1 value is < 200 RFU. In this case, perform the test again using a new QCV strip from a new lot (if possible).
- VIDAS QCV should only be used to diagnose problems with the instrument pipette mechanism or optical system. It is not a calibrator.
- VIDAS QCV does not diagnose temperature problems or leakages.
- VIDAS QCV should only be used to test the instrument.

**INDEX OF SYMBOLS**

Symbol	Meaning
	Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Protect from light

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BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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**VIDAS® PTH (1-84) (PTH)****Intended Use**

**VIDAS® PTH (1-84) (PTH)** is a third generation automated quantitative assay for use on the VIDAS® family of instruments, for the quantitative measurement of the biologically active parathyroid hormone PTH (1-84) in human serum or plasma using the ELFA (Enzyme Linked Fluorescent Assay) technique.

Used in conjunction with other laboratory findings and clinical assessments, this assay is intended for use as follows:

- As an aid in the diagnosis of hyper or hypoparathyroidism.
- As an aid for the monitoring of calcium homeostasis in patients with chronic kidney disease.

**Summary and Explanation**

Parathyroid hormone (PTH) is the main hormone involved in calcium and phosphorus homeostasis. PTH is formed in the parathyroid glands (essentially four, located on the back of the thyroid gland) and secreted as an 84 amino-acid peptide [1-84 according to the universal numbering convention beginning with the N-terminus] when the extracellular calcium concentration decreases.

The main function of the PTH is to maintain the ionized calcium level within a narrow range by a powerful hypercalcemic effect, by stimulating calcium release from bone and distal tubular reabsorption of calcium by the kidney, and by increasing renal production of 1,25-dihydroxyvitamin D-3 (calcitriol) which stimulates intestinal calcium absorption.

PTH increases urine phosphate levels and decreases serum phosphate levels. PTH stimulates bone remodeling by mobilizing osteoblasts, thus justifying its use in osteoporosis treatment to reduce the risk of vertebral fractures.

The secretion of PTH is promoted by low calcium concentrations and high phosphate concentrations, and is inhibited by high calcium concentrations and low phosphate concentrations. It is also inhibited by calcitriol, PTH (7-84) fragment, and magnesium deficiency.

In the bloodstream, the half-life of PTH (1-84) does not exceed 4 minutes. PTH is rapidly degraded by the liver into non-(1-84) C-terminal fragments, essentially (7-84) and (53-84) with long half-lives, and which are eliminated by the kidneys. These fragments accumulate in the blood when the kidneys fail to function.

Parathyroid gland disorders lead to either hyperparathyroidism (primary, secondary or tertiary) or hypoparathyroidism. Apart from bone fragility<sup>(1)</sup> or nephrolithiasis, the associated symptoms may not be very specific, and asymptomatic parathyroid disorder may be unexpectedly detected during a routine blood test<sup>(2)</sup>.

This biological assay contributes in determining a patient's status regarding the function of the parathyroid glands as normal or hyper- or hypoparathyroidism, when interpreted with other biological parameters (such as calcium, phosphate, vitamin D) and clinical signs.

PTH measurement is indicated when performing a differential diagnosis for hypercalcemia or hypocalcemia. It enables the diagnosis of primary or secondary hyperparathyroidism<sup>(3)</sup>, and the diagnosis of hypoparathyroidism<sup>(4)</sup> mainly after parathyroid surgery. It is also prescribed for patients with chronic kidney disease, in particular in the follow-up of dialysis patients<sup>(5)</sup>.

In case of kidney failure, medical monitoring is guided by the KDOQI (Kidney Disease Outcomes Quality Initiative) and KDIGO (Kidney Disease Improving Global Outcomes) recommendations.

In non-kidney disorders, PTH levels in blood are one of the results of routine biological tests investigated by a healthcare specialist (for example: Endocrinology, Rheumatology, Geriatrics, Metabolic Diseases).

**Principle**

The assay principle combines a 2-step sandwich enzyme immunoassay method with a final fluorescence detection (ELFA). VIDAS® PTH (1-84) is a third generation assay traced to the WHO International Standard for Parathyroid Hormone 1-84, recombinant, coded 95/646.

The single-use Solid Phase Receptacle (SPR) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and pre-dispensed in the sealed single-use reagent strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR device several times.

The sample is collected and then diluted in a suitable diluent. The mixture is cycled in and out of the SPR device several times. This operation enables PTH to bind with the capture antibodies fixed to the interior wall of the SPR device.

Unbound components are eliminated during washing steps.

The conjugate mixture is cycled in and out of the SPR device several times. This operation enables the conjugate to bind with the immune complexes that are fixed to the interior wall of the SPR device and to form a sandwich. Unbound components are eliminated during washing steps.

During the final detection step, the substrate (4-Methylumbelliferyl phosphate) is cycled in and out of the SPR device. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of PTH in the sample.

At the end of the assay, the results are automatically calculated and reported in pg/mL by the instrument according to the calibration curve stored in memory. The results can then be printed out.

### Content of the Kit (30 tests)

30 PTH Strips <sup>(a)</sup>	STR	Ready-to-use.
30 PTH Solid Phase Receptacles 1 x 30	SPR	Ready-to-use. Interior of SPR device coated with mouse monoclonal anti-PTH immunoglobulins (specific to the N-terminal (1-6) fragment).
PTH Calibrator <sup>(b)</sup> 1 x 4.4 mL (liquid)	S1	Ready-to-use. Buffer containing recombinant PTH (1-84) + Stabilizer of animal origin + preservative. MLE (Master Lot Entry) data indicate the acceptable range in "Relative Fluorescence Value" ("Calibrator (S1) RFV Range").
PTH Control <sup>(b)</sup> 1 x 2.3 mL (liquid)	C1	Ready-to-use. Buffer containing recombinant PTH (1-84) + Stabilizer of animal origin + preservative. MLE data indicate the acceptable range in pg/mL ("Control C1 Dose Value Range").
Specifications for the factory master data required to calibrate the assay: MLE barcode printed on the box label.		
1 package insert provided in the kit or downloadable from <a href="http://www.biomerieux.com">www.biomerieux.com</a> .		

(a) DANGER  
P305 + P351 + P338

WARNING

EUH208 / H317 / H318 / P261 / P280 / P302 + P352 /

(b) WARNING

EUH208 / H317 / P261 / P280 / P302 + P352

#### Hazard statements

EUH208: Contains 2-methyl-2H-isothiazolin-3-one. May produce an allergic reaction.

H317: May cause an allergic skin reaction.

H318: Causes serious eye damage.

#### Precautionary statements

P261: Avoid breathing dust/fume/gas/mist/vapours/spray.

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P302 + P352: IF ON SKIN: Wash with plenty of soap and water.

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For further information, consult the Safety Data Sheet.

### The SPR device

The interior of the SPR device is coated during production with mouse monoclonal anti-PTH N-ter (1-6) fragment immunoglobulins. Each SPR device is identified by the "PTH" code.

Only remove the required number of SPR devices from the pouch and carefully reseal the pouch after opening.

### Reagent Strip

The strip consists of 10 wells covered with a labeled foil seal. The label comprises a barcode which mainly indicates the assay code, kit lot number, and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

### Description of the VIDAS® PTH (1-84) (PTH) strip

Well	Reagents
1	Sample well: dispense 300 µL of calibrator, control, or sample.
2	Sample diluent: Buffer + Protein stabilizer of animal origin + preservative (300 µL).
3	Empty well.
4 - 5	Wash buffers: Buffer + preservative (600 µL).
6	Conjugate: Buffer containing anti-(7-84) PTH antibodies + Stabilizer of animal origin + preservative (400 µL).
7 - 8	Wash buffers: Buffer + preservative (600 µL).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine (DEA) (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

### Materials and Disposables Required but Not Provided

- Single-use pipette and/or micropipettes to dispense the appropriate volumes.
- Powderless disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instruments of the VIDAS® family.

### Warnings and Precautions

- **For *in vitro* diagnostic use only.**
- **For professional use only, by qualified laboratory personnel, in clinical laboratories.**
- This kit does not contain products of human origin.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest; do not inhale).
- Handle all patient specimens as potentially infectious material and follow the usual safety precautions.
- Although the PTH (1-84) recombinant proteins present in the S1 Calibrator and C1 Control are inactive, handle the tubes as potentially infectious products.
- Do not use the SPR devices if the pouch is pierced or if the dot sealing a SPR device has come unstuck.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.

- Do not mix reagents (or disposables) from different lots.
- VIDAS® PTH (1-84) assay reagents are only for use with the instruments of the VIDAS® family.
- Use powderless gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- Refer to the hazard statements “H” and precautionary statements “P” indicated above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (refer to the User Manual for user and preventive maintenance operations).

## Storage Conditions

---

- Store the VIDAS® PTH (1-84) kit at +2°C/+8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at +2°C/+8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPR devices.
- **After use, carefully reseal the pouch with the desiccant inside to maintain stability of the SPR devices, and return the complete kit to +2°C/+8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the box label.

## Samples

---

### Sample type and collection

Human serum or plasma.

### Types of tubes validated:

- Plastic tube with clot activator
- Plastic tube with clot activator and separation gel
- Plastic tube with lithium heparin
- Plastic EDTA tube

It is recommended to validate collection tubes before use as some contain substances which interfere with test results.

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### Sample Preparation

The current WHO/DIL/LAB/99.1 document provides recommendations for sample preparation. <sup>(6)</sup>

For use of sample tubes, refer to the tube manufacturer's recommendations for use.

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For serum samples, ensure that complete clot formation has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting times.

### Sample stability

Samples (serum and plasma) can be stored in closed primary tubes at +18°C/+25°C for up to 8 hours and at +2°C/+8°C for up to 48 hours. If longer storage is required, freeze the serum or plasma at -19°C/-31°C.

A study performed on samples frozen for 6 months at -19°C/-31°C showed that the quality of results is not affected. Three freeze/thaw cycles were validated.

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- hemolysis (after spiking samples with hemoglobin, up to 5 g/L (monomer)),
- lipemia (after spiking samples with lipids, up to 30 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin, up to 0.3 g/L).

It is recommended not to use hemolyzed, lipemic, icteric samples, and, if possible, to collect a new sample.

Refer to the section **PERFORMANCE – Study of drugs and other potentially interfering substances** for the compounds tested.

## Instructions for Use

---

**For complete instructions, see the Instrument User Manual.**

### Reading VIDAS® PTC (Protocol Test Change) data and MLE data

#### When using the assay for the first time:

With the external instrument barcode reader, **scan the barcodes (PTC and MLE) in the following order:**

1. According to the instrument used, scan the PTC barcode(s) downloadable from [www.biomerieux.com](http://www.biomerieux.com). This reading allows VIDAS® PTC protocol data to be transferred to the instrument software for its update.
2. Scan the MLE data on the box label.

#### When opening a new lot of reagents:

With the external instrument barcode reader, scan the MLE data on the box label before performing the test.

**Note: The master lot data need only be entered once for each lot.**

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User Manual).

#### Calibration

The calibrator is traceable to the WHO International Standard 95/646 for PTH (1-84).<sup>(7)</sup>

Calibration, using the Calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the MLE data have been entered, and then every **84 days**. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf life of the kit.

The Calibrator identified by S1 must be tested in **duplicate** (refer to the User Manual). The Calibrator value must be within the set RFV (Relative Fluorescence Value) range. If this is not the case, recalibrate using S1.

**Caution:** The assay Calibrator (S1) and Control (C1) are ready-to-use reagents made liquid by the addition of a tension-active stabilizer. When pipetting the 300 µL calibration and calibration control portion, make sure that the entire 300 µL volume is dispensed out of the pipet tip.

#### Kit controls

One Control is included in each VIDAS® PTH (1-84) (PTH) kit. This Control must be performed immediately after opening a new kit to ensure that reagent performance has not been altered.

Each calibration must also be checked using this Control. The instrument will only be able to check the control value if it is identified by C1.

Results cannot be validated if the control value deviates from the expected values.

**Note:** The aim of the VIDAS® PTH (1-84) Quality Control (C1) is to validate calibrations. The use of the VIDAS® PTH (1-84) Quality Control (C1) as an Internal Quality Control is under the customer's responsibility.

#### Procedure

1. **Remove the kit from storage at +2°C/+8°C and take out the required reagents. Carefully reseal the SPR pouch and return the kit to +2°C/+8°C. The reagents can be used immediately.**

2. Use one "PTH" strip and one "PTH" SPR device for each sample, control, or calibrator to be tested. Make sure the SPR pouch has been carefully resealed after the required Solid Phase Receptacles have been removed.
3. The test is identified by the "PTH" code on the instrument. The calibrator must be identified by "S1" and tested in duplicate. The control identified by "C1" must be tested singly.
4. Mix the calibrator, control, and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
5. For optimal results, refer to all the paragraphs in the **SAMPLES** section.
6. Before pipetting, ensure that the samples, calibrator, control, and diluent are free of bubbles.
7. **For this test, the calibrator, control, and sample test portion is 300 µL.**
8. Insert the "PTH" Solid Phase Receptacles and "PTH" strips into the instrument. Check to make sure that the color labels with the assay code on the Solid Phase Receptacles and the Reagent Strips match.
9. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
10. Close the vials and return them to the required temperature after pipetting.
11. The assay will be completed within approximately **24 minutes**. After the assay is completed, remove the Solid Phase Receptacles (SPR) and strips from the instrument.
12. Dispose of the used Solid Phase Receptacles (SPR) and strips into an appropriate recipient.

## Quality Control

---

Additional quality controls can be performed in accordance with local regulations or requirements related to accreditation, as well as requirements defined in the laboratory's quality control procedure.

**Note:** It is the responsibility of the user to perform Quality Control in accordance with any applicable local regulations.

## Results and Interpretation

---

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR device is introduced into the substrate.

The second reading is taken after incubating the substrate with the enzyme bound to the interior of the SPR device. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (predefined mathematical model) and are expressed as pg/mL.

Interpretation of test results should be made taking into consideration the patient's clinical history, and the results of any other tests performed.

The instrument displays the VIDAS® PTH (1-84) (PTH) assay results from 4.0 to 1500.0 pg/mL.

### Sample dilution

Samples with VIDAS® PTH (1-84) concentrations greater than 1500.0 pg/mL must be retested after 1/2 dilution (dilution factor = 2) in saline solution (1 volume of sample + 1 volume of saline solution) in order to obtain a result within the measurement range.

**The dilution must be performed manually on the VIDAS®, MINI VIDAS®, and VIDAS® 3 instruments.**

If the dilution factor (= 2) was entered when the Work List was created, the result is calculated automatically. If the dilution factor was not entered, multiply the result by the dilution factor to obtain the sample concentration.

## Limitations of the Method

---

- Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's clinical history and the results of any other tests performed.
- Refer to the section **PERFORMANCE – Study of drugs and other potentially interfering substances** for the compounds tested.
- Any results that do not correspond to the patient's clinical history may be due to inadequate instrument maintenance (see the Instrument User Manual).

## Reference Values

These results are given as a guide. It is recommended that each laboratory establish its own reference values from a rigorously selected population.

Reference values were determined in a clinical study using 491 samples from apparently healthy individuals. The population was divided into 3 groups, depending on the vitamin D status. For each group, a 95% reference value range (or reference interval, according to CLSI EP28-A3c Guideline) was determined, as well as for the total of the apparently healthy population as described below:

	N	Median (in pg/mL)	PTH (1-84) reference interval (in pg/mL)	90% CI*of the Lower reference limit (in pg/mL)	90% CI of the Upper reference limit (in pg/mL)
Healthy population	491	20.8	9.2 - 44.6	[8.7 ; 9.7]	[42.5 ; 46.7]

Healthy population by Vitamin D level	N	Median (in pg/mL)	PTH (1-84) reference interval (in pg/mL)	90% CI of the Lower reference limit (in pg/mL)	90% CI of the Upper reference limit (in pg/mL)
Deficient ≤ 20 ng/mL	176	21.5	8.9 - 45.3	[8.0 ; 9.9]	[42.1 ; 48.8]
Insufficient > 20 and < 30 ng/mL	213	20.7	8.8 - 47.8	[7.6 ; 9.9]	[41.6 ; 57.3]
Normal ≥ 30 ng/mL	102	21.3	9.6 - 47.3	[8.7 ; 10.6]	[41.6 ; 54.1]

\* CI: Confidence Interval.

According to the KDIGO guidelines, the PTH concentration of the patients treated by dialysis should be maintained within two and nine times the upper normal limit of the assay.

## Performance

Studies performed using the VIDAS® PTH (1-84) (PTH) assay gave the following results.

### Measurement range

The measurement range (or analytical measuring interval, according to CLSI EP36-Ed1 Guideline) is the range of values corresponding to the acceptable performance limits (precision and linearity).

The measurement range of the VIDAS® PTH (1-84) (PTH) assay is 4.0 to 1500.0 pg/mL.

### Linearity

Linearity was evaluated according to CLSI EP06-A recommendations. The VIDAS® PTH (1-84) assay is linear between 4.0 and 1500.0 pg/mL.

### Detection limits

Limit of Blank (LoB)	0.8 pg/mL
Limit of Detection (LoD)	2.2 pg/mL
Limit of Quantitation (LoQ)	4.0 pg/mL

LoB, LoD and LoQ values were determined according to CLSI EP17-A2 recommendations.

The Limit of Blank (LoB) is the concentration below which 95% of analyte-free samples are found.

The Limit of Detection (LoD) is the lowest concentration of analyte in a sample that can be distinguished from the analyte-free sample with a probability of 95% (observed result greater than the LoB with a 95% probability).

The Limit of Quantitation (LoQ) is the lowest concentration of analyte that can be detected and measured with an acceptable level of precision. All results under 4.0 pg/mL are printed  $\leq$  4.0 pg/mL. For the VIDAS® PTH (1-84) assay, the acceptable level of precision corresponds to within-lot precision fixed at 20% CV.

#### Hook effect

The assay uses a 2-step immunoassay sandwich method. Hook effect can therefore be excluded by design.

#### Precision

A precision study was performed according to CLSI EP05-A3 recommendations. A panel of human samples representing 6 concentration levels in the measurement interval was analyzed on the VIDAS® instrument to include the following main sources of variability: repeatability, run, day, calibration, and lot.

Repeatability (within-run precision), intermediate precision (total within-lot) and within-laboratory precision (between-lot within-instrument) were estimated for each sample. The values obtained are reported in the following table, as a guide.

Sample	Concentration level (pg/mL)	Repeatability		Intermediate precision		Within-laboratory	
		Standard deviation (pg/mL)	CV (%)	Standard deviation (pg/mL)	CV (%)	Standard deviation (pg/mL)	CV (%)
Sample 1	14.7	0.9	6.0	1.1	7.3	1.2	7.9
Sample 2	46.3	2.8	6.1	3.0	6.5	5.7	12.4
Sample 3	106.7	4.0	3.8	5.1	4.8	5.6	5.2
Sample 4	299.0	10.3	3.4	13.1	4.4	16.2	5.4
Sample 5	590.0	15.3	2.6	25.4	4.3	30.7	5.2
Sample 6	1130.1	42.5	3.8	51.1	4.5	62.4	5.5

These results are an indicator showing the ability to monitor chronic kidney disease patients in the routine mode.

#### Analytical specificity

The analytical specificity of the VIDAS® PTH (1-84) (PTH) assay was established by testing cross-reactive compounds according to CLSI document EP7-A2 recommendations. Cross-reactivity was evaluated by spiking serum samples with cross-reactive compounds (about 50 pg/mL and about 400 pg/mL of PTH).

The results of this study are reported in the following table.

Tested substance	Tested concentration (pg/mL)	Cross-reactivity %
PTH (7-84)	100 000	0.11%
PTH (1-34)	100 000	0.04%
PTH (13-34)	100 000	0.0004%
PTH (39-68)	100 000	0.0004%
PTH (44-68)	100 000	0.009%
PTH (39-84)	100 000	0.002%
PTH (53-84)	100 000	0.006%
Calcitonin	100 000	0.005%
Osteocalcin	100 000	0.014%
C-telopeptid	100 000	0.023%

#### Study of drugs and other potentially interfering substances

Potential interference by commonly used drugs and other substances was studied according to CLSI EP7-A2 recommendations. No significant interference was detected up to the maximum concentrations tested.

Tested drug	Maximum concentration	Tested drug	Maximum concentration
Acetaminophen	35 mg/dL	Etidronate	105 mg/dL
Acetylsalicylic acid	65 mg/dL	Ibuprofen	50 mg/dL (242.5 $\mu$ mol/dL)
Alendronate	8 mg/dL	Lanthanum chloride	40 mg/dL
Alfacalcidol	2.5 $\mu$ g/mL	Magnesium chloride	40 mg/dL
Aluminum sulfate	40 mg/dL	Pamidronate	18 mg/dL

Tested drug	Maximum concentration	Tested drug	Maximum concentration
Biotin	2 µg/mL	Risedronate	6 mg/dL
Calcitriol	1 ng/mL	Salicylic acid	60 mg/dL
Calcium acetate	40 mg/dL	Vitamin D2	1800 ng/mL
Calcium citrate	40 mg/dL	Vitamin D3	240 ng/mL

Tested substance	Maximum concentration
Bilirubin (unconjugated)	0.3 g/L (510 µmol/L)
Cholesterol	4.21 g/L
Hemoglobin	5 g/L (310 µmol/L)
Human Albumin	60 g/L
Human Anti Mouse Antibodies (HAMA)	2 µg/mL
Lipids (intralipids)	30 g/L
Rheumatoid factor	800 IU/mL

### Method comparison

The comparison studies were performed according to CLSI EP09-A3 recommendations.

All the intended populations covering the measuring interval from 4.0 to 1500.0 pg/mL have been tested with VIDAS® PTH (1-84) (PTH) assay. Among the 215 patients enrolled, the following populations were represented:

- Adenoma
- Bone disorder
- Dialysis
- Non parathyroidian hypercalcemia
- Parathyroidectomy
- Renal failure stage 3
- Renal failure stage 4
- Renal failure stage 5
- Renal transplanted
- Apparently healthy population with  $20 < \text{Vit D} < 30 \text{ ng/mL}$
- Apparently healthy population with  $\text{Vit D} \geq 30 \text{ ng/mL}$

The comparison of the VIDAS® PTH (1-84) assay (Y) with another commercially available PTH (1-84) immunoassay (X) gave the following results:


Number of samples tested: 215

Samples tested between 4.1 pg/mL and 1179.4 pg/mL for the VIDAS® PTH (1-84) (PTH) assay.

Samples tested between 6.68 pg/mL and 1200.00 pg/mL for the other commercially available immunoassay.

Equation for Passing-Bablok regression:  $Y = 1.0460 X - 2.2489$

Correlation Coefficient: 0.9889

### Monitoring performance

A total number of 21 patients were enrolled with 3 different time points. In conclusion a total of 63 samples were collected and tested with VIDAS® PTH (1-84) (PTH) assay (Y) and another commercially available PTH (1-84) immunoassay (X).

According to the KDIGO guidelines, the PTH concentration of the patients treated by dialysis should be maintained within two and nine times the upper normal limit of the assay. For information, the highest reference value established in a clinical study was 44.6 pg/mL, meaning an interval of PTH values between 89.2 and 401.4 pg/mL. Correspondence was defined by a concordant interpretation between both methods.

The global (all subjects and time points mixed) frequency and percent of concordant interpretation between the 2 assays regarding the KDIGO criterion was calculated.

Criteria	% Concordance
Concordant interpretation between the 2 assays regarding the KDIGO criterion	88.89 %

### Waste Disposal

Dispose of used or unused reagents, as well as any other contaminated disposable materials, following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced, according to their nature and degree of hazardousness, and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

### Literature References

1. Camacho PM *et al.* – American Association of Clinical Endocrinologists and American College of Endocrinology Clinical Practice Guidelines for the Diagnosis and Treatment of Postmenopausal Osteoporosis - 2016. – *Endocrine Practice* – 2016 September; 22 (Suppl 4): 1-42.
2. Khan AA *et al.* – Primary hyperparathyroidism: review and recommendations on evaluation, diagnosis, and management. A Canadian and international consensus. – *Osteoporos Int.* – 2017; 28: 1-19.
3. Wilhelm SM *et al.* – The American Association of Endocrine Surgeons Guidelines for Definitive Management of Primary Hyperparathyroidism. – *JAMA Surg.* – 2016 Oct; 151 (10): 959-68.
4. Bollerslev J *et al.* – European Society of Endocrinology Clinical Guideline: Treatment of chronic hypoparathyroidism in adults. – *European Journal of Endocrinology* – 2015; 173, G1-20.

5. KDIGO 2017 Clinical Practice Guideline Update for the Diagnosis, Evaluation, Prevention, and Treatment of Chronic Kidney Disease-Mineral and Bone Disorder (CKD-MBD). – *Kidney International Supplements* – 2017 July; 7 (1): 1-59.
6. WORLD HEALTH ORGANIZATION, USE OF ANTICOAGULANTS IN DIAGNOSTIC, LABORATORY INVESTIGATIONS, 2002, WHO/DIL/LAB/99.1 Rev.2
7. WORLD HEALTH ORGANIZATION, Parathyroid Hormone 1-84, human, recombinant (1st International Standard) - WHO\_BS\_09.2115.

## Index of Symbols

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

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## Revision History

Change type categories

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

**Note:** *Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release Date	Part Number	Change Type	Change Summary
2018-10	050049-01	N/A	Not applicable (First publication)

---

Release Date	Part Number	Change Type	Change Summary
2018-12	050049-02	Correction	PERFORMANCE - Method comparison
2019-11	050049-03	Technical	Content of the Kit (30 tests) / Warnings and Precautions
2022-02	050049-04	Technical change	Materials and Disposables Required but Not Provided

---

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries or one of its companies.

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CLSI yra „Clinical Laboratory and Standards Institute, Inc.“ priklausantis prekės ženklas.