

A surgeon in a blue operating room, wearing a blue surgical cap, clear safety glasses, and a blue surgical mask. The surgeon is wearing blue scrubs and white gloves, and is holding a handpiece in their right hand. The background is dark, and the overall lighting is blue.

**MIDAS REX™ MR8™**  
HIGH SPEED  
DRILL SYSTEM

Advancing  
surgical possibilities

**Medtronic**  
Further, Together

## MR8™ electric and electric plus motors

### Multifunction FCU

**Note:** Conduct a system check by pressing the foot pedal to briefly run the motor and confirm proper function prior to any procedure. If motor overheats, discontinue use of the motor.


**Important:** By default, the FCU is activated by pressing the corresponding button for at least 100 milliseconds (mS). Use the IPC™ system touch screen settings screen to change the default value.

To use the multifunction FCU (Figure 9) to control the motor, do the following:

- To select FWD or REV mode, press the mode button (Figure 9, Number 1). 
  - To change the motor, press the motor button (Figure 9, Number 2). 
  - To toggle between the start/stop mode and variable speed mode, press the control button (Figure 9, Number 3). 
- Note:** Functionality of the control button may be changed in the motor Defaults menu to pause irrigation. Refer to **Change System Settings**, in the **Pre-Operating** instructions section of the **IPC™ User's Guide**.
- To start or adjust the speed of a motor in variable mode, press the foot pedal (Figure 9, Number 5).

### Cleaning the FCU

For cleaning instructions refer to the section, **Clean the FCU**, in the **Cleaning and sterilization** section of the **Integrated Power Console (IPC™) System User's Guide**.



### Technical specifications and performance characteristics

**Table 1: MR8™ electric motor technical specifications and performance characteristics**

MR8™ electric motor (EM800)			
Size	Weight	Speed	Duty cycle for applied part
Length: 9.73 cm (3.83 in) x Diameter: 1.55 cm (0.61 in) Length of motor cable: 460 cm (181 in)	87 grams	FWD/REV, variable speed of 200–75,000 rpm	For operating room temperatures up to 40°C (104°F), the electric motor is rated for 3 minutes (03:00) of continuous cutting time at 75,000 rpm, followed by 25 minutes (25:00) of rest.  For normal operating room temperatures (typically 20°C/68°F), the electric motor is rated for continuous cutting at 75,000 rpm.

*Handwritten notes: 1.2, 1.1, 1.4, 1.5*

**Table 2: MR8™ electric plus motor technical specifications and performance characteristics**

MR8™ electric plus motor (EM850)			
Size	Weight	Speed	Duty cycle for applied part
Length: 9.73 cm (3.83 in) x Diameter: 2.10 cm (0.83 in) Length of motor cable: 460 cm (181 in)	171 grams	FWD/REV, variable speed of 200–75,000 rpm	For operating room temperatures up to 40°C (104°F), the electric plus motor is rated for 3 minutes (03:00) of continuous cutting time at 75,000 rpm, followed by 25 minutes (25:00) of rest.  For normal operating room temperatures (typically 20°C/68°F), the electric plus motor is rated for continuous cutting at 75,000 rpm.



ACCESS /  
CLOSURE

# THE POWER TO DO MORE

## THE INTEGRATED POWER CONSOLE (IPC®) PLATFORM

Inspired by listening to your needs, the IPC® platform delivers maximum torque and ease of use, it's the first multispecialty powered surgery console from Medtronic for cranial, spine, and ENT. Our innovative IPC system drives more handpieces with more power than ever before.

IntelliFlow™  
irrigation  
remote control

Legend EHS Stylus®  
motor delivers  
>50% more torque

Color coded  
connectors  
for ease of use



Helpful  
graphics for  
each handpiece

Easy-to-use  
touchscreen

Connect up to  
four handpieces  
at once

1.6

**The Midas Rex™ MR8™ Electric High-Speed Drill System**




**MR8™ electric motors**

The electric motor (Figure 6) is a smaller and more compact version of the electric plus motor (Figure 7). The electric touch motor (Figure 8) is a small compact motor that includes a rotating and removable finger lever.

The electric motors are high-speed, high-torque, reversible electric motors used to dissect bone and biomaterial at selectable speeds from 200 to 75,000 rpm. These motors are designed to interface with a series of interchangeable attachments and dissecting tools utilizing a standard locking mechanism. The electric motors provide the power to operate disposable surgical dissecting tools, intended for use in various surgical procedures.

The electric motors can be controlled either by using the FCU or through the finger lever. The finger-controlled electric motors operate in the same manner as the foot-controlled electric motors. Pressing the finger lever increases the speed and lifting off of the finger lever reduces the speed to a stop. The electric motors are reusable devices that are supplied non-sterile and require cleaning and sterilization prior to each surgical use.

**Note:** All motor cables are integrated and cannot be removed from the motor.



1.8

**Use**

The following instructions for the electric motors are in addition to the general assembly instructions found in the **Integrated Power Console (IPC™) System User's Guide**. To obtain a copy of the **IPC™ System User's Guide**, please contact Medtronic or your local distributor.


Use an IPC™ system that has software version V2.7.4.0 or later, or update the IPC™ system to the latest software by contacting your Medtronic Neurosurgery Group sales representative. If outside the USA, contact your Medtronic regional distributor or Medtronic Neurosurgery Group sales representative.

**Connection to IPC™ system**

Locate the motor and FCU ports on the IPC™ system connector panel (Figure 1) and insert the multi-pin connector into either of the two motor ports provided.

**Note:** Prior to use, check for bent or missing pins in the cable connectors.

**Note:** To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.



1. Electric, electric plus, and electric touch motors
2. Intelliflow irrigation remote control
3. FCU

**IPC™ system pump detection**

The IPC™ system incorporates the electric, electric plus, and electric touch motor irrigation tubing sets at Pump 1. If you do not use irrigation for the motors, manually change Pump 1 to **None**. Refer to the **Integrated Power Console (IPC™) System User's Guide, Set up and Prime Pumps** for more information.

When the IPC™ system detects the electric or electric plus motor, Pump 2 defaults to **None**.




When the IPC™ system detects both the electric and the electric plus motor, the system defaults Pump 2 to the Shared configuration. Use the pumps screen to override the Shared default by selecting the electric motor for Pump 1. Refer to the **Integrated Power Console (IPC™) System User's Guide, Set up and Prime Pumps** for more information.

Control the operation of the electric, electric plus and electric touch motors with the IPC™ system touchscreen and the multifunction FCU.

**IPC™ system touchscreen controls**

1.8

To set or adjust the electric (Figure 2-1), the electric plus (Figure 2-2), or electric touch (Figure 2-3) motor controls, on the IPC™ system touchscreen:



In the touchscreen control box, adjust or set motor speed, irrigation flow rate, motor acceleration/deceleration, motor rotation, and SAFE mode option.

5





# SURGICAL TECHNIQUE

## OTHER TIPS

The drill can cut or tear the dura mater when it unlocks. Check some conditions before performing the procedure, such as the existence of adhering dura mater, high intracranial pressure, or other anomalies adjacent to the drilling site.

- Do not use the drill in a hole or an area that has already been drilled, or the surrounding areas. The composition of reconstructed bone differs from that of healthy bone, increasing the risk of unwanted breakage. Leave a minimum spacing of 15 mm from one hole to another.

- It is essential to keep the drill perpendicular (90°) at the predetermined point of the skull to be drilled, as an excessive deviation from perpendicularity may cause the product to fail and lead to serious patient injury. The perforator can nick or tear the dura or brain if not held at a 90° angle to the skull.
- Select a drill bit suitable for the bone thickness according to the table below. This prevents the drill from tearing the bone or brain tissue (similar effect when the bit is not positioned at 90°).
- Take extra care with the pressure applied when drilling the skull of infants, children, and the elderly or pathological bone.

## ORDERING INFORMATION

Item #	Colour	Outer Diam.	Inner Diam.	Shelf	Minimum Bone Thickness
DM0008FAA		9.0 mm	6.0 mm	3.0 mm	3.25 mm
DM0011FAA		11.0 mm	7.0 mm	3.0 mm	3.25 mm
DM0210FAA		14.0 mm	11.0 mm	1.5 mm	1.75 mm
DM0010FAA		14.0 mm	11.0 mm	3.0 mm	3.25 mm
Item #	Description				
*AD03	Midas Rex® Perforator Driver Attachment Manufactured by Medtronic				

## Medtronic

**Europe**  
Medtronic International  
Trading Sàrl  
Route du Mallet, 31  
Case postale  
CH-1137 Tolochenaz  
Tel: +41 (0)21 302 70 00  
Fax: +41 (0)21 302 73 00

**United Kingdom/Ireland**  
Medtronic UK Ltd  
Building 9  
Crawley Green Business Park  
Widford  
Hertfordshire WD13 3WW  
UK  
www.medtronic.co.uk  
Tel: +44 (0)1923 212313  
Fax: +44 (0)1923 241004

Manufactured by:



Distributed by Medtronic

UC201701133EE ©2016 Medtronic  
All rights reserved. Printed in Europe

Rx only. Refer to product instruction manual/  
package insert for instructions, warnings,  
precautions and contraindications.

For more information, contact your Medtronic  
Neurosurgery sales representative or refer to  
www.MedtronicNeurosurgery.com.

For a listing of indications, contraindications,  
precautions, warnings, and potential adverse  
events, please refer to the Instructions for Use.

### Jacobs chuck attachments

The Jacobs chuck attachment is a non-cannulated, 5/32" chuck with key for drilling.


**Note:** A drill bit may be installed in the attachment before the Jacobs chuck attachment is installed on the motor.

#### Assembly


1. Slide the attachment over the motor collet aligning the triangular markers (Figure 26-1). An audible click, heard and perceptible by touch, confirms that the attachment is fully seated.
2. Rotate the attachment to the locked position (Figure 26-2).
3. To install a drill bit, rotate the Jacobs chuck key counterclockwise to open the jaw (Figure 26-3).
4. Insert the drill bit and rotate the Jacobs key clockwise to tighten the jaws until the drill bit is secure.

#### Disassembly

To open the jaws, turn the Jacobs chuck key counterclockwise. Remove and discard the drill bit. Rotate the Jacobs chuck attachment to the unlocked position and lift the attachment off of the motor.




**Small Bore Attachments**  
ACCOMMODATE HRS™  
2.4 mm DIAMETER TOOLS



**Large Bore Attachments**  
ACCOMMODATE HRS™  
3.3 mm DIAMETER TOOLS



**Telescoping Attachments**




**Footed Attachments**



**Small Bore Variable Exposure Attachments**

Variable exposure attachments allow you to vary the setting and exposure. Small bore variable attachments can increase exposure up to 10mm in 10 use only. (Check for safety & waste exposure requirements)



**Large Bore Variable Exposure Attachments**

Large bore variable exposure attachments allow you to vary the setting and exposure up to 12mm in 12 use only. (Check for safety & waste exposure requirements)



**Perforator Driver and Jacobs™ Chuck**

The Perforator Driver has a Platinum coated bit so it is easy to use. The Perforator Driver has a 1/2" diameter bit. (Check for safety & waste exposure requirements)



**J-Latch**

The J-Latch is a new feature that allows you to use the J-Latch to hold the bit in place. (Check for safety & waste exposure requirements)



**Metal Cutting Attachments**

Metal cutting attachments allow you to cut through metal. (Check for safety & waste exposure requirements)







AF01 AF02 AF03 AF01R AF02R

Rotating

## FOOTED ATTACHMENTS

Also referred to as craniotome attachments, footed attachments are used in cranial and spinal procedures.

The footed attachments allow the surgeon to protect soft tissue when cutting through bone.



AC16

## CONTRA ANGLED ATTACHMENT

The contra angled attachment is designed for limited access right-angled applications.



ASMC

## METAL CUTTING ATTACHMENT

This attachment can accommodate tungsten carbide or diamond tools for the trans-section of metal.



ASJL

## J-LATCH ATTACHMENT

ASJL universal attachment that accepts all manufacturer tools with a J shaped end; compatible with all Midas Rex motors.

7

AD02



AD01 (800 RPM) AD03 (1000RPM)

## PERFORATOR DRIVER AND JACOBS® CHUCK

The Perforator Driver is used to drive perforator bits with a Hudson shank, while the non-cannulated Jacobs Chuck drives pins and wires.

11

## Metal cutting attachments

### Warnings:

- Do not use metal cutter dissecting tools on bone.
- For metal transection, observe the following safety guidelines:
  - a. Wear eye protection.
  - b. Irrigate copiously to cool the cutting surfaces
  - c. Protect the wound site from metal debris.
  - d. Use a clamp or grasping device to control loose fragments during transection of any metal component.
- Use only dissecting tools specifically designed for use with this drill system. Match the nomenclature and color code on the dissecting tool packaging to the same nomenclature and color band on the attachment.

### Notes:

- The metal cutting attachments use the tungsten carbide or diamond wheel dissecting tools.
- A dissecting tool may be installed and locked in the attachment before the metal cutting attachment is installed on the motor.
- The Legend system metal cutter dissecting tool can be used with the MR8-ASMC9 attachment. Match the color code on the metal cutter Legend dissecting tool packaging to the same color band on the attachment.

All metal cutter dissecting tools have an "MC" attachment prefix in their nomenclature (for example MR8-9MC30). Metal cutter dissecting tools should not be installed into other attachments.

The metal cutting attachments come in a double locking (DK) configuration.

### Assembly:

1. Slide the metal cutting attachment over the motor collet aligning triangular marks on the attachment and the motor case. An audible click, heard and perceptible by touch, confirms that the attachment is fully seated.
2. Rotate the attachment to the locked position on the motor case.
3. To insert the tool, make sure that the tool locking ring is in the unlocked position, and insert the dissecting tool into the top of the tube (Figure 23).
4. Rotate the dissecting tool until it drops into position and is fully seated (Figure 23).
5. Rotate the tool locking ring to the locked position. Gently pull on the shaft of the dissecting tool to verify proper installation.

### Disassembly:

1. Rotate the tool locking ring to the unlocked position. Pull the dissecting tool from the attachment.
2. Discard the tool.
3. Push the attachment distally, then rotate the attachment to the unlocked position on the motor case and lift the attachment off the motor.

23

MR8™ metal cutting attachment



9

**Telescoping attachments and tubes**

Telescoping attachments provide support to the rotating dissecting tool. Telescoping tubes are disposable following multiple uses and should be discarded when heat or excessive vibration is noticed or insertion of tools becomes difficult. Telescoping attachment MR8-AT10ADK is provided in a double locking (DK) configuration. The straight telescoping attachment (MR8-AT10) is available in both a DK and non-DK configuration.

**Warnings:** Use only dissecting tools specifically designed for use with this drill system. Match the nomenclature and color code on the dissecting tool packaging to the same nomenclature and color band on the attachment.

- Notes:**
- To assemble DK attachments, follow the instructions for the respective attachment. To disassemble DK attachments, push the attachment distally before rotating.
  - The Legend system 12 cm, 14 cm, 15 cm, and 18 cm telescoping tubes and dissecting tools can be used with the MR8™ system. Match the color code on the telescoping Legend dissecting tool packaging to the same color band on the telescoping tube.

**Assembly**


1. Slide the attachment over the motor collet aligning the triangular markers. A tactile click confirms the attachment is fully seated.
 

**Note:** A dissecting tool and telescoping tube may be installed and locked in the attachment before the telescoping attachment is installed on the motor.
2. Rotate the attachment to the locked position.
3. Rotate the tube locking ring toward the unlocked icon (Figure 24-1).
4. Insert the base of the telescoping tube into the attachment (Figure 24-2).
5. To lock the tube in place, rotate the tube locking ring towards the locked icon. Verify that the tube is secure by gently pulling on the tube. Do not over tighten the tube locking ring.
6. To insert the tool, make sure that the tool locking ring is in the unlocked position, and insert the dissecting tool into the top of the tube (Figure 24-3). A tactile click confirms that the tool is fully seated.
7. Rotate the tool locking ring to the locked position. Gently pull on the shaft of the dissecting tool to verify proper installation (Figure 24-4).
8. If the tube position needs to be changed, rotate the tube locking ring toward the unlocked icon and reposition the tube. Then, rotate the tube locking ring toward the locked icon.

**Disassembly**

1. To remove the attachment, rotate the tool locking ring and the tube locking ring to the unlocked position, and pull the telescoping tube and tool out of the attachment.
2. Rotate the attachment to the unlocked position and lift the attachment off the motor.
 


**Note:** Telescoping tubes should be discarded when heat or excessive vibration is noticed, or when insertion of the tool becomes difficult.



13

# Telescoping Attachments

The MR8™ system telescoping attachments offer a range of straight, angled, curved and hooded tubes.





# EASYDRILL™ CRANIAL PERFORATORS

The EasyDrill™ Cranial Perforators combine efficient bone cutting performance with the protection of AutoStop technology.

## PROTECTION

Employs a clutch mechanism to automatically disengage once it ceases to find resistance to bone.

## PERFORMANCE

Drill tip design helps prevent slipping, while sharp edges cut quickly through bone.

## VARIETY

Four color-coded sizes to meet a wide range of cranium types.

## ORDERING INFORMATION

23.1 in 24.1

Item #	Colour	Outer Diam.	Inner Diam.	Shelf	Minimum Bone Thickness
DM0008FAA	Orange	9.0 mm	6.0 mm	3.0 mm	3.25 mm
DM0011FAA	Blue	11.0 mm	7.0 mm	3.0 mm	3.25 mm
DM0210FAA	Green	14.0 mm	11.0 mm	1.5 mm	1.75 mm
DM0010FAA	White	14.0 mm	11.0 mm	3.0 mm	3.25 mm

23.2

Item #	Description
*AD03	Midas Rex® Perforator Driver Attachment Manufactured by Medtronic

## Medtronic

Europe  
Medtronic International  
Trading Sdn Bhd  
Route 06, Malimau 57  
Cusack Centre  
C-4-1131, Tropicana 2  
Tel: +41 (0)21 802 75 00  
Fax: +41 (0)21 802 73 00

United Kingdom/Ireland:  
Medtronic UK Ltd  
Building 9  
Crowley Green Business Park  
Weston  
Hertfordshire WD18 8WW  
UK  
www.medtronic.co.uk  
Tel: +44 (0)1923 213213  
Fax: +44 (0)1923 241004

Manufactured by:



Distributed by Medtronic

©2014/05-2015 Medtronic  
All rights reserved.  
Prices in Europe

www.medtronic.eu



Rx only. Refer to product instruction manual/  
package insert for instructions, warnings,  
precautions and contraindications.

For more information, contact your Medtronic  
Neurosurgery sales representative or refer to  
[www.MedtronicNeurosurgery.com](http://www.MedtronicNeurosurgery.com).

For a listing of indications, contraindications,  
precautions, warnings, and potential adverse  
events, please refer to the Instructions for Use.

# SURGICAL TECHNIQUE

## THE EASYDRILL™ CRANIAL PERFORATORS

The EasyDrill Cranial Perforator device is a bone cutting and drilling instrument used in conjunction with a surgical motor and a Hudson Chuck - speed reducer attachment, to drill access holes through a patient's skull.

When properly used, the EasyDrill Cranial Perforator employs a clutch mechanism to automatically disengage once perforation is accomplished and as the drill ceases to find resistance to bone.

### EQUIPMENT NEEDED



1. EasyDrill Cranial Perforator™



2. Attachment Driver (AD03)



3. Drill

### EASYDRILL CRANIAL PERFORATOR FUNCTIONAL TEST

Perform the following function test before each trephination. Use sterile disposable gloves to perform the test.

**CAUTION:** If the sterile protective cover is not available, use sterile gauze (as a protective barrier) when handling the cutting end of the EasyDrill Cranial Perforator.

**NOTE:** Do not use the product and contact customer service if either of the two function tests below fail to present the indicated results:

### RECOMMENDED SPEED RANGE

We recommend a speed range of 800 to 1000 RPM for the drilling with EasyDrill Cranial Perforators.



1. Firmly hold the EasyDrill Cranial Perforator's plastic drive shank (Hudson shank) and the protective cover, turning the cutting end of the EasyDrill Cranial Perforator device in a clockwise direction. The cutting end of the device should turn freely.



2. Next, grasp the EasyDrill Cranial Perforator by the plastic drive shank (Hudson shank) and press the protected cutting end of the device into the opposite hand, turn the drill clockwise until it locks; to unlock, cease to apply pressure.

## CRANIAL BONE DRILLING

1. After the functional tests, attach the EasyDrill Cranial Perforator to the Hudson-type chuck of the motorized driver. Verify that the EasyDrill Cranial Perforator is firmly fixed to the Hudson-type chuck before proceeding.
2. Remove the protective cover from the cutting end of the EasyDrill Cranial Perforator device before drilling.
3. Position and press the cutting end of the EasyDrill Cranial Perforator gently against the skull to be perforated.
4. Manually turn the EasyDrill Cranial Perforator until you feel the drill engagement.
5. Position the EasyDrill Cranial Perforator perpendicular to the skull (90°) on the point of perforation. Supply power to the drilling system and begin trephination.
6. Maintain a gentle and continuous pressure while drilling.
7. Irrigate the point of perforation with physiological serum while drilling to prevent bone necrosis.
8. At the drilling depth where the EasyDrill Cranial Perforator device ceases to find resistance to bone, it automatically disengages the cutting end of the device, which stops rotating. Note: At this stage, the plastic shank of the device may appear to continue to rotate.
9. Carefully remove the EasyDrill Cranial Perforator from the skull.
10. There may still be a thin bone layer in the hole that must be carefully removed with appropriate surgical instrumentation. If this layer is adhered to the brain material, take appropriate care not to injure the dura.
11. Before performing another perforation, wipe the cutting end of the EasyDrill Cranial Perforator removing the residues from the previous drilling. Perform testing procedures again.



MR8™ Twist Drill

Head Diameter (mm)	Shank Diameter (mm)	Shank Length (mm)	Head Diameter (mm)	Head Length (mm)	Large Bore
0.8	4.0	4.0	1.5	6.0	MRB-7TD156
1.0	4.0	4.0	1.5	7.0	MRB-7TD17
1.0	6.0	6.0	1.5	8.0	MRB-7TD18
1.1	4.0	4.0	1.5	12.0	MRB-7TD114
1.1	6.0	6.0	1.5	20.0	MRB-7TD123B
1.2	4.0	4.0	1.7	4.0	MRB-7TD174
1.2	6.0	6.0	1.7	6.0	MRB-7TD176
1.2	8.0	8.0	1.7	8.0	MRB-7TD178
1.5	4.0	4.0	1.0	4.0	MRB-7TD264
1.1	6.0	6.0	2.0	6.0	MRB-7TD206
1.5	4.0	4.0	2.2	8.0	MRB-7TD188

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Hole Saw

Head Diameter (mm)	Head Length (mm)	Large Bore
6.4	48.3	MRB-9HS54
12.8	33.0	MRB-9HS108
12.8	33.0	MRB-9HS115

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Hole Maker

Head Diameter (mm)	Head Length (mm)	Large Bore
6.5	10.4	MRB-3HM6
9.5	10.4	MRB-3HM6
10.0	10.4	MRB-3HM150
12.6	10.4	MRB-3HM128

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Metal Cutter

Head Diameter (mm)	Head Length (mm)	Head Diameter (mm)	Head Length (mm)
3.0	6.4	MRB-5MC16	MRB-5MC16
3.0	8.0	MRB-5MC24	MRB-5MC24
5.2	17.9	MRB-5MC30	MRB-5MC30
20.4	13.2	MRB-5MC32	MRB-5MC32

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Telescoping - Ball

Head Diameter (mm)	Head Length (mm)	Head Diameter (mm)	Head Length (mm)
2.0	MRB-T12BA2D	MRB-T12BA2D	MRB-T12BA2D
3.0	MRB-T12BA4D	MRB-T12BA4D	MRB-T12BA4D
4.0	MRB-T12BA6D	MRB-T12BA6D	MRB-T12BA6D

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Telescoping - Ball Diamond

Head Diameter (mm)	Head Length (mm)	Head Diameter (mm)	Head Length (mm)
2.0	MRB-T12BA2D	MRB-T12BA2D	MRB-T12BA2D
3.0	MRB-T12BA4D	MRB-T12BA4D	MRB-T12BA4D
4.0	MRB-T12BA6D	MRB-T12BA6D	MRB-T12BA6D

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Telescoping - Reverse Taper

Head Diameter (mm)	Head Length (mm)	Head Diameter (mm)	Head Length (mm)
4.0	MRB-T12RT4D	MRB-T12RT4D	MRB-T12RT4D
5.4	MRB-T12RT6D	MRB-T12RT6D	MRB-T12RT6D
6.8	MRB-T12RT8D	MRB-T12RT8D	MRB-T12RT8D

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

MR8™ Telescoping - Tapered

Head Diameter (mm)	Head Length (mm)	Head Diameter (mm)	Head Length (mm)
1.1	MRB-T12TA11	MRB-T12TA11	MRB-T12TA11
1.5	MRB-T12TA15	MRB-T12TA15	MRB-T12TA15
2.0	MRB-T12TA20	MRB-T12TA20	MRB-T12TA20

ATTACHMENT

MRB-A127  
MRB-A128  
MRB-A129  
MRB-A130  
MRB-A131  
MRB-A132

29  
30  
31

35  
36  
38  
40

34  
37  
39

\*MRB-7TD114 has a coating, full length of 7.0mm and no length stop.

## ACCESSORIES

### Sterilization Trays and Cases

- CA800 MR8 Instrument Tray 1/2 DIN
- CA800S\* MR8 Sterilization Case 1/2 DIN
- CA850 MR8 Instrument Tray 3/4 DIN
- CA850S\* MR8 Sterilization Case 3/4 DIN
- CA875 MR8 Instrument Tray Navigated 3/4 DIN

### Electric Motor Accessories

- EA815 MR8 Electric Control Switch

### On-Drill Irrigation Tubing

- IRB800 MR8 On-Drill Irrigation Tubing, Single
- IRB820 MR8 On-Drill Irrigation Tubing, 5-Pack

### ClearView Irrigation Tubing

- IRB870 MR8 ClearView Irrigation Tubing, Single
- IRB875 MR8 ClearView Irrigation Tubing, 5-Pack



CA800S



EA815



CA850



EA815

19

- Pneumatic Motor Accessories**
- PA700 Pneumatic Lubricant Diffuser Cartridge
- PC700 Pneumatic Foot Control Unit
- PC720 Pneumatic Regulator Hose
- PA200 Pneumatic Triton Adapter
- PC100 N2 DISS to Male Schrader Adapter
- PC110 N2 DISS to Female Schrader Adapter
- PC120 N2 DISS to Air DISS Adapter
- PC130 N2 DISS to Female Schrader Adapter
- PC140 N2 DISS to Male 3/8 Adapter (AUS)
- PC150 N2 DISS to Male 3/8 Adapter (AUS)

### General Accessories

- PA115 Tool Storage Rack
- PA130 Motor Wrench
- PA305 Small 1.2mm Cleaning Brush 5-Pack
- PA325 Medium 1.6mm Cleaning Brush 5-Pack
- PA325 Large 1.2mm Cleaning Brush 5-Pack



PC700



PC110



PC111



PC120



PC140



PC150



PA700



PA200

# Midas Rex ClearView® Irrigation Tubing Set 33

**EC REP**

Medtronic B. V.  
Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands



Medtronic Xomed, Inc.,  
6743 Southpoint Drive North  
Jacksonville, FL 32216-0980 USA 33

**R**Only **2**

**LOT**

**0223544736**


**STERILE EO**



**2023-06-30**



**DEHP**




M0000030B050 B 0097

20

# Medtronic

REF IRD800 X1

Midas Rex™ MR8™  
Irrigation Tubing Set



MR8

LOT 0228019536

STERILE EO

2028-01-02



DEHP



CE2797



ROnly

M026477C001 A 0106

32

LOT 0228019536



2024-01-03

M026477C001 A 0106

68G1156 B



21





MR8™ Clearview™ LP Ball Diamond

TOOL TYPE	Head Diameter (mm)	5 cm	7 cm
	0.5	MR8-LP08BA0D	
	0.8	MR8-LP08BA0D	
	1.1	MR8-LP08BA0D ✓	MR8-LP07BA0D ✓
	1.2	MR8-LP08BA0D	
	1.5	MR8-LP08BA0D	
	2.0	MR8-LP08BA0D	MR8-LP07BA0D
	3.0	MR8-LP08BA0D	MR8-LP07BA0D
	4.0	MR8-LP08BA0D	MR8-LP07BA0D

48  
50  
52



MR8™ Clearview™ LP Ball

TOOL TYPE	Head Diameter (mm)	5 cm	7 cm
	1.0	MR8-LP08BA0	MR8-LP07BA0
	1.5	MR8-LP08BA0	MR8-LP07BA0
	2.0	MR8-LP08BA0	MR8-LP07BA0
	3.0	MR8-LP08BA0	MR8-LP07BA0
	4.0	MR8-LP08BA0	MR8-LP07BA0

47  
49  
51  
44  
46

MR8™ Clearview™ LP Ball Diamond Coarse

TOOL TYPE	Head Diameter (mm)	5 cm	7 cm
	2.0	MR8-LP08BA0DC	MR8-LP07BA0DC
	3.0	MR8-LP08BA0DC	MR8-LP07BA0DC
	4.0	MR8-LP08BA0DC	MR8-LP07BA0DC



MR8™ Clearview™ SD/SP Matchhead

TOOL TYPE	Head Diameter (mm)	12 cm	14 cm
	2.5	MR8-SD12MH25	
	2.5	MR8-SD12MH25	
	3.0	MR8-SD12MH25	
	3.0	MR8-SD12MH25	
	2.5	MR8-SD12MH25	MR8-SP14MH25
	2.5	MR8-SD12MH25	MR8-SP14MH25
	3.0	MR8-SD12MH25	MR8-SP14MH25
	3.0	MR8-SD12MH25	MR8-SP14MH25 ✓

MR8™ Clearview™ SD/SP Ball

TOOL TYPE	Head Diameter (mm)	12 cm	14 cm
	1.0	MR8-SD12BA10	
	4.0	MR8-SD12BA10	
	1.0	MR8-SD12BA10	MR8-SP14BA10
	4.0	MR8-SD12BA10	MR8-SP14BA10

MR8™ Clearview™ SD/SP Matchhead Diamond

TOOL TYPE	Head Diameter (mm)	12 cm	14 cm
	2.5	MR8-SD12MH25D	
	3.0	MR8-SD12MH25D	
	2.5	MR8-SD12MH25D	MR8-SP14MH25D
	3.0	MR8-SD12MH25D	MR8-SP14MH25D

23

MR8™ Match Head - Diamond Coarse



TOOL TYPE	Head Dia. (mm)	Head Length (mm)	5 cm Large Base	10 cm Small Base	14 cm Large Base	15 cm Small Base
	3.0	3.8	MRS-SPH00C	MRS-1PH00C MRS-1PH00D MRS-1PH00E MRS-1PH00F MRS-1PH00G MRS-1PH00H	MRS-1PH00C MRS-1PH00D MRS-1PH00E MRS-1PH00F MRS-1PH00G MRS-1PH00H	MRS-1PH00C MRS-1PH00D MRS-1PH00E MRS-1PH00F MRS-1PH00G MRS-1PH00H

ATTACHMENT

MR8™ Match Head - Symmetri



TOOL TYPE	Head Dia. (mm)	Head Length (mm)	14 cm Large Base
	3.0	3.8	MRS-1PH00T MRS-1PH00U MRS-1PH00V MRS-1PH00W

ATTACHMENT

MR8™ Match Head - Fluted



Navigation Compatible

TOOL TYPE	Head Dia. (mm)	Head Length (mm)	7 cm Small Base	9 cm Large Base	10 cm Small Base	12 cm Small Base	14 cm Large Base	15 cm Small Base	18 cm Small Base	21 cm Large Base	26 cm Large Base
	1.7	2.2	MRS-7PH07		MRS-10PH17			MRS-15PH07			
	2.2	3.8	MRS-7PH02		MRS-10PH02	MRS-12PH02	MRS-14PH02	MRS-15PH02	MRS-18PH02		
	3.0	3.8	MRS-7PH00	MRS-9PH00	MRS-10PH00	MRS-12PH00	MRS-14PH00	MRS-15PH00	MRS-18PH00	MRS-21PH00	MRS-26PH00
	4.0	5.0	MRS-7PH07	MRS-9PH00	MRS-10PH00	MRS-12PH00	MRS-14PH00	MRS-15PH00	MRS-18PH00	MRS-21PH00	MRS-26PH00

ATTACHMENT

MR8™ Match Head - Diamond



TOOL TYPE	Head Dia. (mm)	Head Length (mm)	7 cm Small Base	9 cm Large Base	10 cm Small Base	12 cm Small Base	14 cm Large Base	15 cm Small Base	18 cm Small Base	21 cm Large Base	26 cm Large Base
	1.7	2.2	MRS-7PH07D		MRS-10PH07D			MRS-15PH07D			
	2.2	3.8	MRS-7PH02D		MRS-10PH02D		MRS-14PH02D	MRS-15PH02D	MRS-18PH02D		
	3.0	3.8	MRS-7PH00D	MRS-9PH00D	MRS-10PH00D	MRS-12PH00D	MRS-14PH00D	MRS-15PH00D	MRS-18PH00D	MRS-21PH00D	MRS-26PH00D
	4.0	5.0	MRS-7PH07D	MRS-9PH00D	MRS-10PH00D	MRS-12PH00D	MRS-14PH00D	MRS-15PH00D	MRS-18PH00D	MRS-21PH00D	MRS-26PH00D

ATTACHMENT



MR8™ Ball - Carbide



TOOL TYPE	Head Diameter (mm)	7 cm Small Base	10 cm Small Base	14 cm Large Base
	1.0	MR8-7BA5C		
	2.0	MR8-7BA6C		
	3.0	MR8-7BA7C		
	4.0	MR8-7BA8C	MR8-10BA8C	
	5.0	MR8-7BA9C		MR8-14BA9C
	6.0	MR8-7BA0C		
	7.0	MR8-7BA1C		
	8.0	MR8-7BA2C		
<b>ATTACHMENT</b>				
		MR8-7BA3C MR8-7BA4C MR8-7BA5C MR8-7BA6C MR8-7BA7C MR8-7BA8C	MR8-10BA3C MR8-10BA4C MR8-10BA5C MR8-10BA6C MR8-10BA7C MR8-10BA8C	MR8-14BA3C MR8-14BA4C MR8-14BA5C MR8-14BA6C MR8-14BA7C MR8-14BA8C

MR8™ Ball - Symmetri



TOOL TYPE	Head Diameter (mm)	7 cm Small Base	9 cm Large Base	10 cm Small Base	12 cm Small Base	14 cm Large Base	15 cm Small Base	18 cm Small Base
	3.0	MR8-3BA3T	MR8-9BA3T	MR8-12BA3T		MR8-14BA3T	MR8-15BA3T	
	4.0	MR8-3BA4T	MR8-9BA4T	MR8-12BA4T		MR8-14BA4T	MR8-15BA4T	MR8-18BA4T
	5.0	MR8-3BA5T	MR8-9BA5T			MR8-14BA5T		MR8-18BA5T
	6.0	MR8-3BA6T	MR8-9BA6T			MR8-14BA6T		
	7.0	MR8-3BA7T	MR8-9BA7T			MR8-14BA7T		
<b>ATTACHMENT</b>								
		MR8-3BA8T MR8-3BA9T MR8-3BA0T MR8-3BA1T MR8-3BA2T	MR8-9BA8T MR8-9BA9T MR8-9BA0T MR8-9BA1T MR8-9BA2T	MR8-12BA8T MR8-12BA9T MR8-12BA0T MR8-12BA1T MR8-12BA2T	MR8-14BA8T MR8-14BA9T MR8-14BA0T MR8-14BA1T MR8-14BA2T	MR8-15BA8T MR8-15BA9T MR8-15BA0T MR8-15BA1T MR8-15BA2T	MR8-18BA8T MR8-18BA9T MR8-18BA0T MR8-18BA1T MR8-18BA2T	MR8-21BA8T MR8-21BA9T MR8-21BA0T MR8-21BA1T MR8-21BA2T

MR8™ Acorn - Fluted



TOOL TYPE	Head Diameter (mm)	7 cm Small Base	9 cm Large Base	10 cm Small Base	14 cm Large Base	15 cm Small Base	21 cm Large Base	25 cm Large Base
	5.0	MR8-7AC50	MR8-9AC50	MR8-10AC50	MR8-14AC50	MR8-15AC50		
	6.0	MR8-7AC60	MR8-9AC60	MR8-10AC60	MR8-14AC60	MR8-15AC60	MR8-21AC60	MR8-25AC60
	7.5	MR8-7AC75	MR8-9AC75		MR8-14AC75		MR8-21AC75	
	9.0	MR8-7AC90	MR8-9AC90		MR8-14AC90		MR8-21AC90	MR8-25AC90
<b>ATTACHMENT</b>								
		MR8-7AC00 MR8-7AC10 MR8-7AC20 MR8-7AC30 MR8-7AC40 MR8-7AC50	MR8-9AC00 MR8-9AC10 MR8-9AC20 MR8-9AC30 MR8-9AC40 MR8-9AC50	MR8-10AC00 MR8-10AC10 MR8-10AC20 MR8-10AC30 MR8-10AC40 MR8-10AC50	MR8-14AC00 MR8-14AC10 MR8-14AC20 MR8-14AC30 MR8-14AC40 MR8-14AC50	MR8-15AC00 MR8-15AC10 MR8-15AC20 MR8-15AC30 MR8-15AC40 MR8-15AC50	MR8-21AC00 MR8-21AC10 MR8-21AC20 MR8-21AC30 MR8-21AC40 MR8-21AC50	MR8-25AC00 MR8-25AC10 MR8-25AC20 MR8-25AC30 MR8-25AC40 MR8-25AC50

56

26



MR8™ Ball - Diamond Extra Coarse

TOOL TYPE	Ball Diameter (mm)	7 cm Small Ball	9 cm Light Ball	10 cm Small Ball	12 cm Light Ball	14 cm Light Ball	15 cm Small Ball	21 cm Large Ball	28 cm Large Ball
	1.0	MRB-7BA00DEL							
	1.5	MRB-7BA1SDVEL							
	2.0	MRB-7BA00DEL							
	3.0	MRB-7BA00DEL							
	4.0	MRB-7BA00DX	MRB-9BA0DX	MRB-10BA0DX	MRB-12BA0DX	MRB-14BA0DX	MRB-15BA0DX	MRB-21BA0DX	MRB-28BA0DX
	5.0	MRB-7BA0DX	MRB-9BA0DX	MRB-10BA0DX	MRB-12BA0DX	MRB-14BA0DX	MRB-15BA0DX	MRB-21BA0DX	MRB-28BA0DX
	6.0	MRB-7BA0DX	MRB-9BA0DX	MRB-10BA0DX	MRB-12BA0DX	MRB-14BA0DX	MRB-15BA0DX	MRB-21BA0DX	MRB-28BA0DX
	7.0			MRB-10BA0DX				MRB-21BA0DX	MRB-28BA0DX
	8.0								
ATTACHMENT									
		MRB-A01E MRB-A02M MRB-A027 MRB-A030M MRB-A031M	MRB-A02E MRB-A02M MRB-A027 MRB-A030M MRB-A031M	MRB-A030 MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D



MR8™ Ball - Diamond Coarse

TOOL TYPE	Ball Diameter (mm)	7 cm Small Ball	9 cm Light Ball	10 cm Small Ball	14 cm Light Ball	15 cm Small Ball
	1.0	MRB-7BA00CL				
	1.5	MRB-7BA1SDCL				
	2.0	MRB-7BA00CL				
	2.5	MRB-7BA00CL				
	3.0	MRB-7BA00CL	MRB-9BA0CL	MRB-10BA0CL	MRB-14BA0CL	MRB-15BA0CL
	3.5	MRB-7BA00CL	MRB-9BA0CL	MRB-10BA0CL	MRB-14BA0CL	MRB-15BA0CL
	4.0	MRB-7BA00CL	MRB-9BA0CL	MRB-10BA0CL	MRB-14BA0CL	MRB-15BA0CL
	5.0	MRB-7BA00CL	MRB-9BA0CL	MRB-10BA0CL	MRB-14BA0CL	MRB-15BA0CL
	6.0	MRB-7BA00CL	MRB-9BA0CL	MRB-10BA0CL	MRB-14BA0CL	MRB-15BA0CL
	7.0					
ATTACHMENT						
		MRB-A01E MRB-A02M MRB-A027 MRB-A030M MRB-A031M	MRB-A02E MRB-A02M MRB-A027 MRB-A030M MRB-A031M	MRB-A030 MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D	MRB-A031E MRB-A031M MRB-A031S MRB-A031T MRB-A031D

## Midas Rex® MR7 Motor Kits

Product	Description	Quantity
PK700	Midas Rex® MR7 Foot Control Motor Kit	

*Kit includes the following:*

PM700	Midas Rex® MR7 Pneumatic Motor	1
PC700	Midas Rex® MR7 Pneumatic Foot Control	1
PC710	Midas Rex® MR7 Regulator Hose	1
PA700	Midas Rex® MR7 Lubricant/Diffuser Cartridge	4-pk
PA600	Midas Rex® MR7 Instrument Case 3/4 DIN	1
175031EN	Midas Rex® Instructions for Use Manual, English	1
175032EN	Midas Rex® MR7 Preventative Maintenance Manual, English	1
PA300	Cleaning Brush, 1.2 mm	1
PA310	Cleaning Brush, 2.4 mm	1
PA320	Cleaning Brush, 3.2 mm	1

Product	Description	Quantity
PK710	Midas Rex® MR7 Finger Control Motor Kit	

*Kit includes the following:*

PM710	Midas Rex® MR7 Touch Pneumatic Motor	1
PC700	Midas Rex® MR7 Pneumatic Foot Control	1
PC710	Midas Rex® MR7 Regulator Hose	1
PA700	Midas Rex® MR7 Lubricant/Diffuser Cartridge	4-pk
PA600	Midas Rex® MR7 Instrument Case 3/4 DIN	1
175031EN	Midas Rex® Instructions for Use Manual, English	1
175032EN	Midas Rex® MR7 Preventative Maintenance Manual, English	1
PA300	Cleaning Brush, 1.2 mm	1
PA310	Cleaning Brush, 2.4 mm	1
PA320	Cleaning Brush, 3.2 mm	1

58.2

## MR7 accessories

### MR7 pneumatic control unit

The pneumatic control unit (Figure 8-1) provides variable speed motor control through a foot pedal. It also allows the user to switch between finger and foot control of the motor (if applicable).

### Regulator

The Legend Regulator (Figure 8-2) controls the delivery pressure of compressed gas to the Legend pneumatic control unit. The pressure gauges indicate cylinder pressure (right gauge) and delivery pressure (left gauge).

**Note:** Outlet pressure gauge accurate to +/- 12 psi.

### Regulator hose

Connects the gas source to the pneumatic control unit to deliver compressed gas.

### MR7 instrument case

The instrument case (Figure 8-3) is used to organize equipment and to secure instruments during sterilization and transport.

### Motor wrench

Used to align arrows on motor collet flats prior to installation of a Legend attachment (Figure 8-4).

### Lubricant/diffuser cartridge (PA700)

The lubricant/diffuser cartridge (Figure 8-5) provides lubrication to the motor and filters oil from exhausted air. Can be used before sterilization of the instrument, for internal protection of the instrument.

**Caution:** The Legend lubricant/diffuser cartridge and the MR7 lubricant/diffuser cartridge are not interchangeable.

58.1

