# VICRYL™ Mesh

## STERILE

VICRYL™ Mesh STERILE	<b>en</b>	Rete VICRYL™ STERILE	<b>it</b>
Χιρουργικό платно VICRYL™ СТЕРИЛНО	<b>bg</b>	VICRYL™ tinklelis STERILUS	<b>lt</b>
Sítka VICRYL™ STERILNÍ	<b>cs</b>	VICRYL™ tīkliņš STERILS	<b>lv</b>
VICRYL™ mesh STERIL	<b>da</b>	VICRYL™-mesh STERIEL	<b>nl</b>
VICRYL™ Netzimplantat STERIL	<b>de</b>	VICRYL™ nett STERILT	<b>no</b>
Πλέγμα VICRYL™ ΑΠΟΣΤΕΙΡΩΜΕΝΟ	<b>el</b>	Siatka VICRYL™ JAŁOWA	<b>pl</b>
Malla VICRYL™ ESTÉRIL	<b>es</b>	Rede VICRYL™ ESTÉRIL	<b>pt</b>
Võrk VICRYL™ STERIILNE	<b>et</b>	Plasă VICRYL™ STERILĂ	<b>ro</b>
VICRYL™-verkko STERIILI	<b>fi</b>	Sietka VICRYL™ STERILNÁ	<b>sk</b>
Treillis VICRYL™ STÉRILE	<b>fr</b>	Mrežica VICRYL™ STERILNO	<b>sl</b>
Mrežica VICRYL™ STERILNA	<b>hr</b>	VICRYL™ nät STERILT	<b>sv</b>
VICRYL™ háló STERIL	<b>hu</b>	VICRYL™ Mesh STERİL	<b>tr</b>

## Instructions for Use

# VICRYL™ Mesh

## STERILE

### Introduction

This package insert is designed to provide instructions for use of VICRYL™ Mesh. It is not a comprehensive reference to surgical technique.

The selection of mesh for any given patient is a function of numerous factors including, but not limited to, the patient's past medical and surgical history, current medical condition (i.e., comorbidities), surgical technique, and size and location of the wound or organ support. The healthcare professional is advised to consult the medical literature regarding techniques, complications, and adverse reactions before selecting a mesh.

### Description

VICRYL™ Mesh is a synthetic absorbable sterile mesh made of 90% glycolide and 10% L-lactide (Polyglactin 910). The empirical formula of the copolymer is  $(C_2H_2O_2)_m(C_3H_4O_2)_n$ . The copolymer is identical in composition to that used in VICRYL™ synthetic absorbable suture.

VICRYL™ Mesh has a maximum weight per area of 800 mg per 100 cm<sup>2</sup>.

VICRYL™ Mesh is intended for use only by healthcare professionals who are trained in the surgical procedures and techniques required for the implantation of absorbable synthetic meshes for wound or organ support.

VICRYL™ Mesh is intended for use as a buttress to provide temporary support during the healing process.

The clinical benefit to be expected is secondary to treatment of temporary wound or organ support where it is required. A summary of safety and clinical performance can be found at the following link: <https://ec.europa.eu/tools/eudamed>.

### Indications / Intended Use

VICRYL™ Mesh may be used wherever temporary wound or organ support is required, particularly in instances in which compliant and

stretchable support is desired. VICRYL™ Mesh may be cut to the shape or size desired for each specific application.

### Contraindications

Because VICRYL™ Mesh is absorbable, it should not be used where extended wound or organ support is required. VICRYL™ Mesh should not be used in conjunction with breast implants and/or breast expanders.

### Patient Target Group(s)

The selection of mesh for any given patient is a function of numerous factors including, but not limited to, the patient's past medical and surgical history, current medical condition (i.e., comorbidities), surgical technique, and size and location where temporary wound or organ support is required.

### Warnings

- Failure to properly follow instructions may result in improper functioning of the device and could lead to injury. Please read all information carefully.
- This device is indicated for temporary wound or organ support and not for gynecologic procedures. Gynecologic procedures should be performed using devices indicated for gynecologic repairs.
- VICRYL™ Mesh is indicated for temporary wound or organ support. VICRYL™ Mesh is not recommended to be used unaided for clinical applications that require long term support. If long term support is needed a device that is intended for long term support should be used.
- As with any implant, an acute foreign body response will occur. In some patients, this response can result in one or more of the adverse reactions listed below.
- Insufficient or improper fixation may increase the risk of postoperative complications. Consult the Application/Instructions for Use section.

- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and anyone coming in contact with the device.
- Inspect the mesh carefully before implantation. Do not use the device if it is damaged.
- Discard unintentionally opened devices, partially used devices, partially used packages, and used devices and packaging.
- In the event of a product malfunction before use such as product damage, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at healthcare professional's discretion. In case the product is to be disposed, disposal of the device and packaging should be according to the user's facility policies and standard operating procedures concerning biohazardous materials and waste.

### Precautions

- The safety and effectiveness of pretreating VICRYL™ Mesh with solutions (e.g., saline, medications) prior to implantation have not been studied.
- The safety and effectiveness of VICRYL™ Mesh in neural tissue and in cardiovascular tissue has not been established.

### Adverse Reactions / Undesirable Side Effects

Potential adverse reactions are those typically associated with surgically implantable materials, including infection, inflammation (foreign body sensation/ reaction, local reaction, granuloma, seroma), acute or chronic pain, hematoma, impaired healing (delayed wound healing, adhesion, postprocedural fistula, wound dehiscence), nerve injury, soft tissue injury, extrusion/erosion, excessive contraction or shrinkage of the tissue surrounding the mesh, mesh failure, and treatment failure.

Method of mesh fixation may carry the risk of additional adverse events.

One or more revision surgeries may be necessary to treat the above-mentioned adverse reactions. Revision surgery may not resolve the adverse reactions and may pose a risk of additional adverse reactions. The healthcare professional should convey adverse reactions, undesirable side effects, and risks associated with the product and the procedure to the patient prior to undertaking mesh implantation and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

The healthcare professional is advised to inform the patient about any residual risks, contraindications, and any foreseeable undesirable side

effects prior to undertaking implantation of surgical meshes. If the patient experiences any of the Adverse Events, contact the healthcare professional immediately.

### Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

VICRYL™ Mesh is Magnetic Resonance (MR) Safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350), or reproductive toxicant (H360) based on human evidence and animal studies.

### Application / Instructions for Use

It is recommended that absorbable sutures be placed at least 1 cm from the edge of the mesh. Some surgeons prefer to suture into position a mesh larger than the wound. The edges are then sutured to assure proper closure under correct tension. When all margin sutures have been placed, the excess mesh is trimmed away, leaving at least 6 mm of mesh extending beyond the suture line. The use of thermal cutting devices is not recommended, because it has not been tested.

#### Disposal

Discard unintentionally opened, partially used, and used devices and packaging according to your facility's policies and procedures concerning biohazardous materials and waste.

### Performance / Actions

Two important characteristics describe the *in vivo* function and behavior of VICRYL™ Mesh: reinforced wound strength and the rate of absorption (loss of mass). *In vivo* studies indicate that the dehiscence force of healing abdominal wounds closed with size 4-0 absorbable sutures and reinforced with VICRYL™ Mesh is significantly greater than the sutured incisional wound alone. Explanted VICRYL™ Mesh was found to have 80% of its original burst strength remaining after fourteen days *in vivo*.

VICRYL™ Mesh is essentially absorbed within 90 days.

### Sterility

VICRYL™ Mesh is sterilized using ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged.

### Storage

No special storage conditions required. Do not use after expiry date.

### How Supplied

VICRYL™ Mesh is available in single-use, sterile packets in a variety of sizes.

## Traceability

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted. The following specific information can be found on the device packaging label: Catalogue number, batch code, expiry and manufacturing date, manufacturer name, address, and the website, and a Unique Device Identification bar code with the Global Trade Item Number information.

## Information to be Conveyed to the Patient

### Device Description

VICRYL™ Mesh is made of a manmade material. It is the same material that is used in VICRYL™ sutures. It will help support your tissues while they heal.

### Device Materials

VICRYL™ Mesh is a sterile mesh device made of manmade material called VICRYL™. It is a mixture of two polymers. One is called glycolide. The other is called L-lactide. It is the same material that is used in VICRYL™ Sutures. VICRYL™ Mesh is absorbed by the body.

VICRYL™ Mesh has a maximum weight per area of 800 mg per 100 cm<sup>2</sup>.

### Information for Safe Use

There is always a risk of adverse reactions with any type of surgery. Some reactions could be due to your overall health condition. Some reactions could be related to just having the surgery. Some reactions could be related to the materials being placed into the body. It is difficult to know the source of some adverse reactions. Contact your doctor about any signs of infection or adverse reactions. Be sure to keep up with your doctor visits after surgery.

This device will not cause any issues with a medical test called an MRI. MRI means magnetic resonance imaging.

### Expected Device Lifetime and Follow-up

VICRYL™ Mesh is compatible with your tissues while it is being absorbed. It still has about 80% of its strength after 14 days. It will be essentially absorbed within 90 days.

### Patient Information Portal

[ic.jnjmedicaldevices.com](http://ic.jnjmedicaldevices.com)

You can find any updates about the device on this website. Your implant card will contain information about your specific device. This can also be found in your medical record.

Caution

Do not re-use

Do not re-sterilize

Do not use if package is damaged and consult instructions for use

Manufacturer

Packaging unit

Catalogue number

Batch code

Use-by date

Sterilized using ethylene oxide

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Date of manufacture

MR Safe


Medical Device / for implant card: Device Name

TEAR HERE

Single sterile barrier system with protective packaging inside

Absorbable Mesh Device

Consult instructions for use or consult electronic instructions for use.  
 EU: Call paper on demand help desk to get paper copies free of charge within 7 days.

Unique Device Identifier

Patient information website / for implant card:  
Information website for patients

Patient Identification / for implant card: Patient Name or  
Patient ID

Date / for implant card: Date of Implantation

Health care centre or doctor / for implant card: Name and  
address of the implanting healthcare institution / provider

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