
5 Pirliimo dālis

smith&nephew

Q-FIX^o

All-Suture Implants

The next generation in all suture implant technology


Koofia tikra


siilinis bambulas,
kuris ir nry purty
dolydai remiasi is
kaulq.

Q-FIX[®] all-suture implants


The Q-FIX family of implants is the next generation in all suture implant technology, providing the benefits of an all suture implant with performance characteristics equivalent or better than traditional anchor designs.¹

Biomet JuggerKnot

Smith & Nephew Q-FIX*


Stryker Iconix

sūlīnīs bumbulės

Specifications of the 1.8 and 2.8mm implants

5.1.


Kopija tikra

Ancillaries

Consistent deployment

Q-FIX[®] provides consistent deployment due to the radially expanding implant design and delivery system that acts like a knot pusher.


Inserter is 2mm below bone surface

Radial expanding deployment


Superior fixation¹

The Q-FIX delivery system with automated suture pull to a target of 140 newtons ensures full deployment of the implant within the bone hole, giving the surgeon reproducible bone fixation performance they can count on.


Load (N) to failure

Anchor Displacement

Load (N) @ 2mm Displacement

References


- 1 P/N 54231-01 Rev. A; P/N 49193-01 Rev. A; P/N 51963-01 Rev. A
 - * Reference P/N 49193-01 Rev. A
 - ** Reference P/N 51800
 - *** Reference P/N 51963-01 Rev. A
 - **** Reference P/N 51801
 - ***** Reference P/N 54231-01 Rev. A
- Tested in 20pcfl foam with synthetic cortical shell, including ultimate tensile strength Q-FIX N=10, Juggerknot N=2, 4.5mm Corkscrew FT N=3

Streamlined technique

Implant deployment with a turn of the dial through the drill guide streamlines implant placement.

Kopija tikra

mentartinio naudojimo
gudimo instrumentas

S.1.



Intero ismetovimas

Q-FIX® Applications

5.1.




SLAP/Bankart repair

✓ The 1.8mm Q-FIX® implant and disposable kit is a great option for SLAP and Bankart repairs with implant performance superior to commonly used all-suture anchors and traditional anchors. The small implant that marries with the drill guide ensures accurate implant delivery with limited bone removal.


Rotator cuff repair

✓ The 2.8mm Q-FIX implant's fixation performance and small footprint makes it a great option for rotator cuff repair. With performance similar to traditional corkscrew anchors more than twice the size, the 2.8mm Q-FIX is ideal for medial row fixation where performance is needed and space is limited.



Q-FIX 2.8mm versus Arthrex 4.5mm Corkscrew FT

Acetabular labral repair

With a 1.8mm footprint, Q-FIX is optimal for regions of the acetabular rim with limited bone stock. Pre-loaded with a #2 suture allows the surgeon to maximize points of fixation without adding additional implants. By reducing instrument passes and focusing on consistent deployment, Q-FIX helps streamline acetabular labral fixation in a highly technical procedure.


Q-FIX 1.8mm versus
Arthrex 2.4mm
Biosuturetak


Kopija tikra


Disposable kits

Disposable kit for 1.8mm Q-FIX® implant

Disposable kit for 2.8mm Q-FIX implant

Disposable Kit, Hip, for 1.8mm Q-FIX implant

Ordering Information

Q-FIX® All-Suture Implants

Reference #	Description
25-1800	1.8MM Q-FIX ALL-SUTURE ANCHOR ✓
25-2800	2.8MM Q-FIX ALL-SUTURE ANCHOR ✓
25-1810	DISPOSABLE KIT FOR 1.8MM Q-FIX IMPLANT, INCLUDES DRILL, DRILL GUIDE AND OBTURATOR
25-1811	DISPOSABLE KIT, HIP, FOR 1.8MM Q-FIX IMPLANT, INCLUDES DRILL, DRILL GUIDE AND OBTURATOR
25-2810	DISPOSABLE KIT FOR 2.8MM Q-FIX IMPLANT, INCLUDES DRILL, DRILL GUIDE AND OBTURATOR
25-2820	REUSABLE BONE PUNCH FOR 2.8MM Q-FIX IMPLANT
25-1812	REUSABLE 1.8MM PATHFINDER® OBTURATOR
25-2812	REUSABLE 2.8MM PATHFINDER OBTURATOR

ArthroCare Corporation
 7000 West William Cannon Drive
 Austin, TX 78735
 USA

www.smith-nephew.com

Order Entry: 1-800-343-5717
 Order Entry Fax: 1-888-994-2782

© 2015 Smith & Nephew, Inc. Trademark of Smith & Nephew. Reg. US Pat. & TM Office.

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Rx Only

Rx Only

Q-FIX® Soft Suture Anchor Instructions for Use

DESCRIPTION

The Q-FIX® Soft Suture Anchor System consists of an implant (anchor) and associated instruments which facilitate the attachment of tissue to bone. The implant, which is available in 1.8mm (single loaded) and 2.8mm (double loaded) sizes, is manufactured from braided polyester with ultra high molecular weight polyethylene suture and is pre-loaded into a disposable tool designed to facilitate direct insertion into a pre-drilled/pre-punched bone hole.

INDICATIONS FOR USE

The Q-FIX Soft Suture Anchor System is intended to be used for soft tissue to bone fixation for:

- Shoulder: Bankart lesion repair; SLAP lesion repair; acromio-clavicular repair; capsular shift/capsulolabral reconstruction; deltoid repair; rotator cuff tear repair; biceps tenodesis.
Foot & Ankle: Medial/Lateral repair and reconstruction; midfoot and forefoot repair; Hallux valgus reconstruction; Metatarsal ligament/tendon repair or reconstruction; Achilles tendon repair.
Elbow: Ulnar or radial collateral ligament reconstruction; lateral epicondylitis repair; biceps tendon reattachment.
Knee: Extra-capsular repair: medial collateral ligament (MCL), lateral collateral ligament (LCL) and posterior oblique ligament; iliotibial band tenodesis (IBT); patellar tendon repair; vastus medialis obliquus advancement (VMO); joint capsule closure.
Hip: Acetabular labral repair.

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS FOR USE section.
Pathologic bone conditions such as cystic changes or severe osteopenia, which would impair secure fixation of the implant.
Pathological changes in the affected soft tissues, which would prevent secure fixation.
Insufficient quality or quantity of bone, inadequate bone stock or comminuted bone surface(s) which would compromise secure fixation of the implant.
The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this device and placement of implants must not bridge, disturb or disrupt the growth plate.
Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitations, active or previous infections, etc.
Foreign body sensitivity to implant materials. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
Conditions such as senility, mental illness or alcoholism, which may impact the patient's ability or willingness to restrict activities or comply with physician instructions during the healing period.
The implant is not designed for and should never be used to attach artificial ligaments or other implants.
Use of suture other than ArthroCare® USP #2 MAGNUMWIRE® suture.

WARNINGS AND PRECAUTIONS

- Do not use if package is opened or damaged.
Detailed instructions on the use and limitations of the device should be given to the patient.
Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases. The implant is a single use device to be used only once for a single patient.
Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete insertion or poor bone quality may result in implant pullout or suture breakage.
Do not bend, apply excessive torque, or

Caution: Do not use the drill or punch without the drill guide as the bone surface can be damaged.

Caution: Care must be taken to not create too deep or too shallow a bone hole. To accomplish this, be sure to drill the hole until the shoulder on the drill bottoms in the drill guide. If using the punch, insure the handle lines up with the proximal edge of the drill guide. Care must be taken to not advance the Drill Guide into the bone as the implant retention force may be compromised.

Caution: Confirm that bone hole is clear of material by advancing and retracting the drill until the flutes of the drill are visually clear of bone fragments upon withdrawal. Failure to clear material from the bone hole may result in difficulty inserting the implant which could compromise retention force.

- 5) Remove the drill or punch. A PATHFINDER® may be used to confirm bone hole location and depth. Caution: Do not attempt to insert the implant into the bone without first creating a bone hole.
6) Insert the Q-FIX implant through the drill guide into the bone hole (see figures, Item 5). Be sure that the implant inserter is fully engaged with the drill guide throughout the deployment of the anchor (see figures, Item 6). This can be accomplished by holding the drill guide and inserter together by hand or by applying downward pressure to maintain the engagement during deployment. Also, be sure to maintain the alignment of the drill guide to the bone hole during insertion.
Caution: Do not use excessive force or pound the inserter into the bone hole as damage to the inserter/implant may occur or result in injury to the patient. If resistance is met on insertion, turn the inserter handle clockwise while keeping a loose hold on the drill guide. If continued resistance is met, remove the inserter and re-drill the existing bone hole.
Caution: Use care to properly align the implant and inserter handle with the bone hole during insertion. A PATHFINDER may be used to confirm bone hole location and depth. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy bent or damaged inserter.
Caution: Be sure to advance the inserter fully into the bone hole. Incomplete insertion or poor bone quality may result in implant pullout.
7) Rotate the activation knob clockwise on the proximal end of the implant inserter to deploy the implant (see figures, Item 7). Rotate the knob until a hard stop is reached and the suture cleat is exposed at the proximal end of the inserter (see figures, Item 8).
8) Release the suture from the cleat (see figures, Item 9) and remove the inserter and drill guide. Pull on the suture tails to ensure that the implant is set in the bone hole.
9) If placing more than one implant, assure that the anchors are located at least 7 mm apart and repeat steps 1 through 9.

MATERIALS

- Anchor implant: Braided Polyester.
Suture: Ultra High Molecular Weight Polyethylene.
Inserter Handle: Surgical grade stainless steel and medical grade plastic.
The product is not made with natural rubber latex.
The product does not contain DEHP.

HOW SUPPLIED

- The Device is provided STERILE for SINGLE USE ONLY. The device is only sterile if the packaging is not opened, damaged or broken.
Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.
The device is sterilized using ethylene oxide.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require

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

REF 25-1800

LOT 2125405

1.8mm Q-FIX[®] Suture Anchor
with MAGNUMWIRE[®] Suture
Blue Co-Braid¹

SIZE 2
Metric 5

5.1.
Sterilize
package

2026-07-05

QTY: (1)

Rx Only

STERILE EO

CE
0123

Consult instructions for use

Keep away from sunlight

Do not re-use

Keep dry

Do not resterilize

Use-by date

Do not use if package is damaged

Manufacturer

EC REP

Authorized Representative in the European Community

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Rx Only

CE 0123

Q-FIX[®] Soft Suture Anchor Instructions for Use

DESCRIPTION

The Q-FIX[®] Soft Suture Anchor System consists of an implant (anchor) and associated instruments which facilitate the attachment of tissue to bone. The implant, which is available in 1.8mm (single loaded) and 2.8mm (double loaded) sizes, is manufactured from braided polyester with ultra high molecular weight polyethylene suture and is pre-loaded into a disposable tool designed to facilitate direct insertion into a pre-drilled/pre-punched bone hole.

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Shoulder: Bankart lesion repair; SLAP lesion repair; acromio-clavicular repair; capsular shift/capsulolabral reconstruction; deltoid repair; rotator cuff tear repair; biceps tenodesis.

Foot & Ankle: Medial/Lateral repair and reconstruction; midfoot and forefoot repair; Hallux valgus reconstruction; Metatarsal ligament/tendon repair or reconstruction; Achilles tendon repair.

Elbow: Ulnar or radial collateral ligament reconstruction; lateral epicondylitis repair; biceps tendon reattachment.

Knee: Extra-capsular repair: medial collateral ligament (MCL), lateral collateral ligament (LCL) and posterior oblique ligament; Iliotibial band tenodesis (IBT); patellar tendon repair; vastus medialis obliquus advancement (VMO); joint capsule closure.

Hand & Wrist: Collateral ligament repair; Scapholunate ligament reconstruction; Tendon transfers in phalanx; Volar plate reconstruction.

Hip: Acetabular labral repair.

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS FOR USE section.

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- 5) Remove the drill or punch. A PATHFINDER[®] may be used to confirm bone hole location and depth. **Caution:** Do not attempt to insert the implant into the bone without first creating a bone hole.
- 6) Insert the Q-FIX implant through the drill guide into the bone hole (see figures, Item 5). Be sure that the implant inserter is fully engaged with the drill guide throughout the deployment of the anchor (see figures, Item 6). This can be accomplished by holding the drill guide and inserter together by hand or by applying downward pressure to maintain the engagement during deployment. Also, be sure to maintain the alignment of the drill guide to the bone hole during insertion.

Caution: Do not use excessive force or pound the inserter into the bone hole as damage to the inserter/implant may occur or result in injury to the patient. If resistance is met on insertion, turn the inserter handle clockwise while keeping a loose hold on the drill guide. If continued resistance is met, remove the inserter and re-drill the existing bone hole.

Caution: Use care to properly align the implant and inserter handle with the h...

68