



Overview

- ◆ **Simultaneously 12 leads Acquisition** 2.
- ◆ Applicable people: Adult and pediatric
- ◆ **Built-in protection against defibrillator**
- ◆ 10.4"HD LED Full Touch Screen with standard grid background
- ◆ Sampling Modes:
 - 9.2 **Real-time/Pre-sampling/Trigger sampling/ Period sampling**
 - ◆ **Manual/Auto/Rhythm Modes** 6.
 - ◆ 3x4, 3x4+1R, 3x4+3R, 6x2, 6x2+1R, 6x2+3R and **12x1** display and simultaneously printing. 9.1
 - ◆ Support Multiple Languages
 - ◆ Electrode Disconnection Alarm
 - ◆ **Automatic measurement & analysis** 13.1
 - ◆ **Support PACE detection** 16.
 - ◆ Max 60 s waveform freezing and review function.
 - ◆ 1 min – 5 minutes Rhythm Analysis.
 - ◆ ZQECG/PDF/PNG/XML/HL7/DICOM-WAVEFORM/DICOM-IMAGE/DICOM-PDF/GDT format report transfer.
 - ◆ Support to connect with scanner, USB mouse, USB keyboard, card reader and so on peripheral products.
- ◆ **Support WiFi** and LAN cable connecting to the third party ECG network system, like: HIS, EMR and **PACS** system without ECG workstation. 10.5
- ◆ Transfer method: FTP/HTTP/Samba/**DICOM** etc.
- ◆ Full alphanumeric capability QWERTY keyboard.

Technique Specification

- ◆ **Sampling rate:** 32000Hz (pcs)
- ◆ **Pace Sample Rate:** 32000Hz
- ◆ **Sample frequency of signal analysis and storage:** 1000 samples/second/channel
- ◆ Time Constant: $\geq 3.2s$
- ◆ **Frequency Response:** 0.01 ~ 350Hz 19.

- ◆ **CMRR:** $\geq 130dB$
- ◆ **Patient Leakage Current:** $<10\mu A$
- ◆ **Input Circuit:** $<50\mu A$
- ◆ **RTI, peak-to-valley value:** $\leq 15\mu V$
- ◆ **Sensitivity valve:** 20 μV
- ◆ **DC bias voltage range:** $\pm 950mV$
- ◆ **Calibration Voltage:** 1mV
- ◆ **A/D conversion:** 24 bit
- ◆ **Input Impedance:** $\geq 60 M\Omega(10 Hz)$ 21.
- ◆ **Sensitivity:** 2.5, 5, 10, 20, 40, Auto mm/mV 4.
- ◆ **Paper Speed:** 5, 6.25, 10, 12.5, 25, 50mm/s 3.
- ◆ **Recording Mode:** Thermal array printer 8 dot/ mm (Vertical), 40 dot/ mm (Horizontal)
- ◆ **Paper Type:** Z-Fold Paper 9.3
- ◆ **Paper Size:** 210mm/216mm \times 140mm--20m 210mm/216mm \times 150mm (A4)--20m (Option)
- ◆ **Filter Setting:** AC: off/50/60Hz.
- ◆ **EMG:** off/25/35/45/75/100/150/250/350Hz 5.
- ◆ **Baseline Drift:** off/0.01/0.05/0.3/0.6Hz
- ◆ **Heart Rate Range:** 30-350 bpm
- ◆ **WiFi Security Standards:**
 - ◆ WEP
 - ◆ WPA/WPA2 PSK
 - ◆ WPA2-Enterprise / ieee802.1x
- ◆ **Compliant to IEEE 802.11b/g/n** 10.1

Physical Specification

- ◆ **Dimension & weight:**
 - ◆ **Equipment:** 285mm*360mm*94mm 3.7kg 24.
 - ◆ **Package:** 345mm *195mm *430mm 6.2kg
- ◆ **Monitor:**
 - ◆ 10.4 inch(L*W=21.13cm*15.85cm, Diagonal: 26.42cm) HD Color LED. 7.
 - ◆ Full Touch Screen
 - ◆ 1024*768 pixel resolution
- ◆ **Electrical Power:**
 - ◆ AC: 100-240V, 50Hz/60Hz, 120VA 22.
 - ◆ DC: Built-in rechargeable battery
- ◆ **Battery Capacity:** 14.52V/5200mAh
- ◆ **Continuous Operating:** DC $\geq 10.5h$ 22.2

Recharge time: The battery is charged to 90% for no more than 3 hours, and charged to 100% for no more than 3.5 hours

Data Storage:

- ◆ **18.5** > 1500 reports(sampling time 10s, 6*2)
- ◆ > 2000 reports(sapling time 2.5s, 6*2)

Environmental Conditions:

- ◆ Transportation Temperature: -20°C~+55°C
- ◆ Relative Humidity: 10%~95%
- ◆ Atmospheric Pressure: 500hPa~1060hPa
- ◆ Operation Temperature: +5°C~+40°C
- ◆ Relative Humidity: 25%~85%
- ◆ Atmospheric Pressure: 570hPa~1060hPa

System Inputs & Outputs

- ◆ 2*USB2.0/ LAN/ SD Interface

Standard Package

- ◆ iMAC 120 ECG Machine (with battery)
- ◆ **Patient Cable** 25.1
- ◆ **Chest Electrodes** 25.2
- ◆ **Limb Electrodes** 25.3
- ◆ **Thermal Printing Paper** 25.5
- ◆ Power Cord
- ◆ Grounding Wire
- ◆ User Manual

Optional Accessories

- ◆ Series Electrodes
- ◆ Thermal Printing Paper
- ◆ SD Memory Card
- ◆ Lithium Battery
- ◆ Workstation Software
- ◆ **Trolley** 25.6
- ◆ **Cable Arm for trolley**
- ◆ Soft bag
- ◆ Laser Printer
- ◆ **Scanner** 11
- ◆ Card Reader
- ◆ **ECG Gel** 25.4

Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Add: 380, High-tech 2nd road, Eastlake High-tech Development Zone, Wuhan, Hubei (430206), China.
Tel: 86-27-87770203 Email: info@zoncare.com
Fax: 86-2787770581 Website: www.zoncare.com/en/

Advanced Function

- | | | |
|----|--|--|
| 11 | <ul style="list-style-type: none">● External laser printer to output the A4 report● Scanner, card reader and other patient information input devices● Wireless WiFi function● DICOM function (Worklist, Storage, MPPS, Storage Commit)● QTc formula: Hodge, Bazett, Fridericia, Framingham● Measurement Matrix● Average Template● R-R histogram, R-R chart● Bradycardia and Tachycardia range adjustable● Lead System: Standard, Posterior wall, Right chest, Right chest posterior wall, Prev. intercostal space, Next intercostal space, CABRERA, Custom● GLASGOW ECG Algorithm can make precise diagnosis according to patient age, gender, racial, taking medicine history, etc., especially for AMI, Ventricular Fibrillation, Ventricular Tachycardia, quarterly Bradycardia, malignant arrhythmia, etc. | <p>Comply with
93/42/EEC As amended by
2007/47/EC</p> <ul style="list-style-type: none">EN ISO 13485:2016EN ISO 14971:2012EN ISO 15223-1:2016EN 1041:2008ISO 10993-1:2009ISO 10993-5:2009ISO10993-10:2010IEC60601-1:2005/A1:2012EN 60601-1-2:2015IEC 60601-2-25:2011EN 62304:2006/A1:2015EN 62366-1:2015EN 60601-1-6:2010/A1:2015 |
|----|--|--|

13.2

Certificates

Quality Standards

ISO9001
ISO13485
CE

NOTE:

- The specifications of this system may change without any prior notification.
- This document contains intellectual property information that is proprietary to Zoncare and is protected by law. Neither the document nor the information contained therein should be used or reproduced in whole or partially, without prior written agreement consent of Zoncare.

Dec. 2021

Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Add: 380, High-tech 2nd road, Eastlake High-tech Development Zone, Wuhan, Hubei (430206), China.
Tel: 86-27-87770203 Email: info@zoncare.com
Fax: 86-2787770581 Website: www.zoncare.com/en/

Product Information:

Name: Digital Multi-Channel Electrocardiograph Machine

Type: iMAC 120

Manufacturer: Wuhan Zoncare Bio-medical Electronics Co., Ltd

Registered Address: #380, High-tech 2nd road, Eastlake high-tech district,
Wuhan, Hubei, P. R. China

Manufacture Address: #380, High-tech 2nd road, Eastlake high-tech district,
Wuhan, Hubei, P. R. China

Production Date: See the label

Operating Life: 10 years

Tel: +86(27)87770581

Tel /Fax: +86(27)87770203

Post Code: 430206

After-sale Service by: Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Website: <http://www.zoncare.com>

Authorized Representative:

Company Name: Well Kang Limited

Company Address: The Black Church, St. Mary's Place, Dublin 7, Ireland

Tel: +353(1)4433560

Fax: +353(1)6864856

Statement

- Wuhan Zoncare Bio-medical Electronics Co., Ltd (hereinafter referred to as Zoncare) reserves the copyright of this unofficially-published user manual.
- No part of this manual may be reproduced, transmitted in any form or by any means, electronic or mechanic, including photocopying, recording, foreign language translating or by any information storage or retrieval system, without prior written permission from Zoncare.
- Zoncare is not responsible for any accidental or consequential damage because of wrongful installation or misoperation. For any consequences by infringing copyright law or violating the third party rights, Zoncare shall hold no responsibility.
- All contents of this manual are regarded as correct but not guaranteed with any responsibility. Parts of the contents in this manual may be changed without prior notice.
- All illustrations in this manual are provided as examples only. Depending on system configuration, screens in the manual may differ from the screens on your system.

Responsibility of the Manufacturer

Zoncare is responsible for the safety, reliability and performance of hardware supplied by Zoncare only if the following conditions are met:

- Assembly operations, expansions, readjustments, modifications or repairs are performed by persons authorized by Zoncare.
- The electrical safety of the room where the unit is installed complies with the requirements of appropriate local, state, and other government regulations.
- The unit is used in accordance with operation requirements.

Safety Conventions in this manual

A Hazard is a source of potential injury to a person, property, or the device.

This manual uses the term WARNING, CAUTION, and NOTICE to point out hazards and to designate a degree or level of seriousness. Familiarize yourself with the following definitions and their significance.

Safety Convention	Definition
 WARNING	◆ It indicates a potential hazard or unsafe practice, which, if not avoided, could result in death or serious injury.
 CAUTION	◆ It indicates a potential hazard or unsafe practice, which, if not avoid, could result in moderate or minor injury.
 NOTICE	◆ It indicates a potential hazard or unsafe practice, which, if not avoided, could result in the loss or destruction of property or data. .
	◆ Mandatory to refer to this manual for equipment use.

Contents

Chapter1 Safety Information and Introduction.....	5
1.1 Safety Information	5
1.1.1 Hazards	5
1.1.2 Warning	5
1.1.3 Cautions	7
1.2 Device Symbols.....	9
1.3 Introduction.....	10
1.3.1 Precautions	10
1.3.2 Product Performance	10
1.3.3 Product Composition.....	10
1.3.4 Intended Application.....	11
1.3.5 Contraindication	11
1.3.6 Theory of Operation	11
1.3.7 Main Features	12
1.3.8 Intended Places and Users.....	13
1.3.9 Safety Standards and Requirements.....	13
Chapter 2 Structure Identification	15
2.1 Identification of Front Panel, Buttons, Symbols and Expansion Slot.....	15
2.1.1 Top View	15
2.1.2 Side View	17
2.1.3 Back View.....	18
2.1.4 Bottom View	18
2.2 Modes.....	19
2.2.1 Standard Mode.....	19
2.2.2 Standby Mode	19
2.2.3 Demo Mode.....	20
2.3 Interface Display.....	21
2.3.1 Waveform Acquisition Interface	21
2.3.2 “Standard ECG”, “Record Analyze”, “Lead System” Interface	24
2.3.3 Patient Info Interface.....	29
2.3.4 Paper Speed Interface	32
2.3.5 Gain Amplitude Interface.....	33
2.3.6 Lead Layout Interface	34
2.3.7 Report Management Interface	35
2.3.8 Application Interface.....	36
Chapter 3 Installation	37

3.1 Preparation for Installation	37
3.1.1 Open-package Inspection	37
3.1.2 Environment Requirements.....	38
3.2 Power Selection.....	40
3.2.1 Connecting to the AC Power.....	40
3.2.2 Battery Power.....	40
3.3 Installing the Recording Paper	41
3.3.1 Installing the Paper Block.....	41
3.3.2 Installation of recording paper.....	42
3.4 Connecting the Patient Cable.....	44
3.5 Power on.....	45
3.5.1 Checks before Power on.....	45
3.5.2 Power on	45
3.6 ECG Machine Setting.....	46
3.7 Power off.....	46
Chapter 4 System Applications.....	47
4.1 Entering Main Menu	47
4.2 Application Interface.....	48
4.2.1 General–Machine Setting.....	48
4.2.2 General-Report Setting	50
4.2.3 General-Recorder Setting	52
4.2.4 General-Filter Setting	55
4.2.5 General-Heartbeat Warning.....	57
4.2.6 General-Institution Information.....	58
4.2.7 General-About the Machine	59
4.2.8 General-Transfer Setting.....	60
4.2.9 General-DICOM Setting.....	64
4.2.10 E-Mail setting.....	66
4.2.11 LAN setting.....	67
4.2.12 WIFI setting	68
4.2.13 Standard ECG setting	68
4.3 Machine Setting.....	73
4.3.1 First time setup.....	73
4.3.2 Setup before Use	73
Chapter 5 Connecting ECG Cable	74
5.1 Environmental Requirements.....	74
5.2 Preparation	75

5.3 Electrode Selection and Usage	75
5.3.1 Examinee's Skin Preparation	75
5.3.2 Electrode Selection	76
5.4 Electrode Placement	80
5.4.1 ECG Patient Cable	80
5.4.2 Limb Lead Placement	81
5.4.3 Chest Lead Placement	81
5.4.4 Pediatric Electrode Placement	82
5.5 Electrode Connection	82
5.6 Lead-off Handling	82
5.7 Entering Patient Information	82
Chapter 6 ECG Acquisition and Recording	83
6.1 Acquisition Preparation	83
6.2 Acquisition and Recording	84
6.2.1 Recording Setup	84
6.2.2 Recording Report	85
6.2.3 ECG Waveform Description	87
6.3 Freeze	88
6.4 Report Storage	90
6.4.1 Report Storage Location	90
6.4.2 Report Storage Mode	90
Chapter 7 Report Management	92
7.1 Report Storage	92
7.2 Report Management	92
7.2.1 Selecting Reports	92
7.2.2 Searching Reports	93
7.2.3 Editing Reports	94
7.2.4 Printing Reports	95
7.2.5 Deleting Reports	95
7.3 Report Transmitting	96
7.4 Report Refresh	98
Chapter 8 Troubleshooting	99
8.1 Interference Problem	99
8.1.1 AC Interference	99
8.1.2 EMG Interference	100
8.1.3 Baseline Drift	101
8.2 Recorder Failure	102

Chapter 9 Maintenance	81
9.1 Cleaning and Disinfection.....	81
9.1.1 Cleaning	82
9.1.2 Disinfection	84
9.1.3 Sterilization	84
9.2 Routine Inspection and Test.....	84
9.2.1 Daily Inspection	84
9.2.2 Regular Inspection	85
9.3 Battery Usage & Maintenance.....	86
9.3.1 Overview	86
9.3.2 Battery Charging	87
9.3.3 Battery Replacement.....	87
9.3.4 Battery Guidance.....	87
9.3.5 Battery Maintenance	88
9.3.6 Battery Recycling	90
9.4 Usage and Maintenance of Recording Paper	90
9.5 Maintenance of Electrodes and Lead Wires	90
Chapter 10 After-sale Service information.....	92
Chapter 11 Accessories.....	94
Appendix I	95
I.1 Performance Index	95
I.2 Safety Index.....	99
I.3 Power Specifications	99
I.4 Appearance Parameters.....	100
I.5 Environmental Conditions	100
I.6 Adherence to Standards.....	100
Appendix II Electromagnetic (EMC)	102
Appendix III Environmental Statement	111

Chapter1 Safety Information and Introduction

1.1 Safety Information

1.1.1 Hazards

This product does not involve any information of hazard level.

1.1.2 Warning



WARNING

- ◆ *This ECG machine can only be operated on one single patient at one time.*
 - ◆ *The unit, its accompanying cables and accessories shall be checked prior to use in order to guarantee that they work normally and safely.*
 - ◆ *Explosion hazard. Do not use the ECG machine in the presence of flammable anesthetics, gases or chemicals. Otherwise it will incur explosion or fire.*
 - ◆ *The unit can only be connected to AC power outlet with grounding protection. If proper grounding cannot be guaranteed, you shall operate the unit on built-in rechargeable battery instead of AC power. The unit shall also be well grounded to avoid the risk of electric shock. Please place the unit where is easy for proper grounding.*
 - ◆ *The ECG machine is to work in an environment free from interferences caused by high-voltage cable, X-ray machine, Ultrasound scanner, and Electrotherapeutic equipment. Do not use the unit in an environment with high static electricity. Otherwise the unit might be affected by electromagnetic interference.*
 - ◆ *Do not open the unit cover, or there might be risk of electrical shock. Only service personnel trained and authorized by the manufacturer can repair or upgrade the unit.*
 - ◆ *Keep the ECG machine away from water. Do not install or store the unit in the*
-



WARNING

place where chemicals are stored or ventilation is poor. Keep the unit away from excessive humidity, temperature, dust, salt, and sulphate.

- ◆ *The ECG machine should be placed gently on a stable platform and be protected from tilting, excessive vibration, and/or shocking in the process of transportation.*
- ◆ *Since excessive leakage current in total will harm the patient, only IEC 60601-1 class I equipment is allowed to be connected to the ECG machine. Therefore, the manufacturer of connected equipment shall hold relevant responsibility for leak current monitoring. When the unit is being used together with other instruments, attention shall be paid to good connections so as to avoid incorrect diagnosis. If necessary, you shall consult a professional technician for advice.*
- ◆ *Electrodes and connectors can only be in touch with the patient, but shall not be in touch with other conductor parts, including the earth.*
- ◆ *The operators shall not leave the exam room when the ECG machine is in operation. They shall keep careful observation on the patient and, if necessary, turn off the power or disconnect the electrodes to assure patient safety. If an accident happens during operation, please power off the unit immediately and check it*
- ◆ *Chemicals from a broken LCD display panel are toxic when ingested. Be cautious when handling an ECG machine with a broken display panel.'*
- ◆ *According to Standard IEC 60601-1, this ECG machine belongs to type CF defibrillation-proof equipment, therefore its applied part can be connected conductively to human heart.*
- ◆ *Used together with a defibrillator, the unit's defibrillator protection is only guaranteed with the manufacturer recommended defibrillator protected electrodes and cables (for specifications, see chapter 11-Accessories). If the defibrillation takes more than 5 seconds or the unit is used with the high*



WARNING

frequency equipment, please use standard disposable electrodes so as to prevent metal electrodes from burning the patient's skin. While used together with other electrical stimulator, the unit shall be operated under instructions from professionals at present.

- ◆ *Do not touch the patient during defibrillation. Otherwise it will incur serious injury or death.*
 - ◆ *ECG signal acquisition may be affected by special environment, incorrect operation of ECG machine and patient condition. For safety information, please refer to the corresponding chapter in this manual.*
 - ◆ *Using unspecified patient cable, limb clamp, and suction bulb may degrade anti-interference performance of ECG machine. Connectivity of patient cable should be checked periodically, at least once a month.*
 - ◆ *It is advised to use the unit with paper suggested by the manufacturer, as this is the only way the printer head life time and clear ECG waveform can be guaranteed.*
 - ◆ *The ECG machine is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms. Physiological waveform and parameters displayed in this ECG machine are for the doctors' reference only, and cannot be used as the basis for clinical treatment.*
-

1.1.3 Cautions



CAUTION

- ◆ *Please used the accessories specified in this manual.*
 - ◆ *When the unit and accompanying will exceed their life time, dispose of them in accordance with relevant local laws and regulations or the regime of local hospitals.*
-

-
- ◆ *Electromagnetic field will affect the performance of this unit. Therefore other equipment used in the vicinity of this unit must comply with relevant EMC requirements.*
 - ◆ *Before connecting the unit to an AC power outlet, please verify its power voltage and frequency comply with the label on the unit or the requirements specified in this manual.*
 - ◆ *Please properly install and carry this unit to prevent it from dropping, collision, strong oscillation or damage by other external mechanical forces.*
 - ◆ *Please install the unit in the place available to observe, operate and maintain.*
 - ◆ *Place this manual near the unit so that it is available when necessary.*
 - ◆ *In order to describe and record electrocardiographs more accurately, the ECG machine should be placed and used in a quiet and comfortable environment.*
-

1.2 Device Symbols

The following table describes symbols or icons that may be on your device and its packaging.

Symbol	Description	Symbol	Description
	Attention! Consult accompanying documentation		On/Off
	DC power supply		AC power supply
	Charging battery		SD card
	Type CF equipment equipped with protector against defibrillation.		Equipotential terminal
	Ethernet port		USB port
	Manufacturing date		Serial number
	China Metrology Accreditation		EU authorized representative
	Indoor use		CE Certification
	Refer to this manual for equipment use		Environment-friendly use period
	It indicates this device contains electronic or electrical components that must not be disposed of as unsorted municipal waste but collected separately. Contact an authorized representative of the manufacturer for information concerning the decommissioning of your device.		



CAUTION

- ◆ *Do not damage any label on the unit.*
 - ◆ *These labels provide important information for the safety and operation of the unit. Damaging or moving the labels may lead to disoperation.*
-

1.3 Introduction

1.3.1 Precautions

You are required to read through the operating instructions before use the ECG machine so as to ensure proper operation of the unit.

1.3.2 Product Performance

Power voltage: AC100V-240V ;

Power frequency: 50/60Hz \pm 1Hz , 75VA;

Continuous operation time: more than 10.5 h;

Front-end acquisition mode: A/D sampling. A/D digits are not less than 24 digits.

The effective sampling is not less than 32000 samples per second (or 32000 Hz/channel); Gain: the equipment is available in 40mm/mV, 20mm/mV, 10mm/mV, 5mm/mV and 2.5mm/mV and automatic six gears.; Gain accuracy is \pm 3%;

Paper speed: the equipment is available in 5mm/s, 6.25mm/s, 10mm/s, 12.5mm/s, 25mm/s, and 50mm/s. Accuracy is \pm 3%;

General safety for the product should conform to IEC 60601-1-2:2012 standard.

Particular safety for the product should conform to IEC 60601-2-25:2011 standard.

1.3.3 Product Composition

The ECG machine is mainly composed of the host, patient cable, limb electrodes, and chest electrodes.

1.3.4 Intended Application

The unit is intended to extract ECG wave group from human body for its morphological and rhythm analysis.

1.3.5 Applicable people

Adult and Children

1.3.6 Contraindication

None

1.3.7 Side effects

No side effects.

1.3.8 Theory of Operation

Schematics:

Only qualified field service engineer will be provided schematics and spare parts list for the ECG machine.

Theory of Operation:

The ECG machine acquires a microvolt level signal from human body surface via patient cable and electrodes. It amplifies the signal by the Amplification Module before Analog to Digital (A/D) converting. Following A/D conversion, the signal is processed by the CPU of the Keyboard Control Module. The CPU outputs the signal to the thermal printer. Precision control program from the Keyboard Control Module is used to drive the stepping motor so as to make the recording paper run at a constant speed. By controlling temperature of the above mentioned thermal emitting components, relevant ECG trace and character could be printed out on the thermal recording paper. In addition, the Keyboard Control Module also processes the keyboard signal, and controls the trace display, the real time clock, and etc. The Power Supply Module provides other modules of the ECG machine with power sources, of which AC power is of priority to the rechargeable battery. When the

ECG machine is powered by AC power, battery is charged provided there is no operation of the unit.

1.3.9 Main Features

- It adopts high-resolution thermal array printout system, uses thermal Zip fold paper of 216mm or 210mm in width, and records clear accurate ECG waveform and information about lead marks, gain(sensitivity), time reference (or paper speed), filter status etc.
- It adopts a unique high-precision digital filter to prevent baseline drift and other interference without causing waveform distortion. It enhances the ability of anti-baseline drift, which is convenient to interpret the waveform.
- It has man-machine interaction interface, independent full alphanumeric keyboard and full touchscreen. Also it can be connected with mouse and keyboard externally. It operates flexibly and easily, reduces workload, and improves working efficiency.
- Support English input with more complete information;
- Support various file output formats, and meet the needs of clinical information.
- Color LCD display, which can record 12-channel ECG waveform and information.
- Real-time waveform freezing;
- Age grading function;
- External USB disk and SD card are supported to store as many patient reports as possible;
- Standby mode, reduces power consumption and extends LCD life ;
- It is designed according to the safety standard of IEC class I type CF. ECG amplifier is fully floated and has good safety performance.
- It has both AC and DC power supply modes. Rechargeable environmentally-friendly lithium-ion battery is installed in the machine and is easy to replace. And it has a dedicated battery charging circuit and a perfect battery management and protection system.

1.3.10 Intended Places and Users

- It can be used as measuring instrument in relevant hospital sections or patient rooms.
- It can be used for mass exam.
- ECG signal measurement is intended for users of different ages, but for the children under 2 years old, the interpretations may have deviation.
- When use the unit on a pacemaker patient, please activate the PACE Detection with reference to *Section 2.3.2.2*.



WARNING

- ◆ *Wherever the ECG machine is used, you must verify that it is connected with reliable dedicated grounding wire.*
-

1.3.11 Safety Standards and Requirements

- Strictly in accordance with IEC 60601-1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance; IEC60601-2-25 Medical electrical equipment-Part1-2:General requirements for basic safety and essential performance - Collateralstandard: Electromagnetic compatibility-Requirements and tests. The security type is type I CF.
- Power supply conditions for the room where the ECG machine is to work shall be suitable for standard three-plug socket with the grounding plug properly grounded. Otherwise, you are required to ground the unit using the accompanying grounding cable, with one end connected to the grounding post and the other end to the ground.



WARNING

- ◆ *Grounding shall be performed in accordance with relevant standards or under the guidance of experienced electricians.*
-

- When the unit is used in combination with other medical devices, you are required to connect the grounding cable together with that of the others so as to protect the patient from possible shock.
- Connect one end of the grounding cable to the equipotential post of the ECG machine and other end to the ground. Avoid using water piping or other piping as a grounding conductor, otherwise the precautions for grounding protection of the unit shall lose effectiveness and the patient may have the risk of electrical shock.
- The ECG machine is an instrument of continuous operation and ordinary equipment. Avoid ingress of liquids into the unit. Explosion hazard. Do not use the unit in the presence of flammable anesthetics or gases.

Classification :

Type of protection against electrical shock	Class I internally powered equipment
Degree of protection against electrical shock	Type CF defibrillation-proof applied part
Degree of protection against harmful ingress of liquids	Ordinary equipment (enclosed equipment without protection against ingress of liquids).
Degree of safety of application in the presence of flammable gas	Equipment is not appropriate for use in the presence of flammable gas.
Signal input and output:	With input and output parts
Mode of operation	Continuous operation

Chapter 2 Structure Identification

2.1 Identification of Front Panel, Buttons, Symbols and Expansion Slot

2.1.1 Top View




Fig. 2.1 Top View

1. Touch Screen: Displaying ECG parameters and waveforms, patient and system information; for more interface information, please see *Section 2.3*. —
2. Indicator Light: Indicating power supply status and battery charge status; for details, see the table below.
3. Operation Panel: Operation keys; for details, see the table below.

Identification of Front Panel Keys	
Key	Function
	“ON/OFF”: Provided the unit is powered on, key “ON/OFF” could be used to change its condition between “ON” and “OFF”.

Identification of Front Panel Keys	
Key	Function
F1~F7	Real-time function keys, which correspond to the hotkeys displayed at the bottom of the screen
	“Mode”: For switch over the unit’s operation mode between Automatic and Manual.
	Switching leads in manual mode
	Switching leads in manual mode
	Quickly save the ECG
	“Patient ID”: Provided on Waveform Acquisition interface, key “Patient ID” could be used to quick enter into Patient Info interface.
	“ECG” could be used to print the ECG trace.
Keyboard	Composite keyboard for information input and some function keys

Identification of Front Panel Buttons	
Indicator Light	Function
	AC power indicator light On: AC power connected Off: AC power disconnected
	DC power indicator light On: power supplied by battery Off: power not supplied by battery, or no battery installed
	Charging status indicator light On: Battery being charged Off: No battery installed or battery fully charged

2.1.2 Side View




Fig. 2.2 Right Side View

1. Patient Cable Connector: For connection with the patient cable.

2. SD Card Slot: External SD storage.

3. USB Port: To connect external USB devices.

18.4

4. USB Port: To connect external USB devices.

5. LAN Connector: For Network connection.




Fig. 2.3 Left Side View

1. Printer Door Switch: Toggle the printer door upwards to open it.

2. Printer Door: To close the printer.

2.1.3 Back View




Fig. 2.4 Back View

18.2

1 Equipotential terminal: When other equipment is used together with the electrocardiograph, the other equipment and the equipotential terminal of the electrocardiograph should be connected by wires to eliminate the potential difference between different devices to ensure safety.

2 AC power connector: Connect the power cable to provide AC power to the ECG.

2.1.4 Bottom View




Fig.2.5 Bottom View

Battery Compartment: To contain the battery.

2.2 Operation Mode

2.2.1 Standard Mode

Automatically enters the standard ECG mode after power on. In this mode you can perform ECG measurements, record waveform and measured value results. You can also perform system setup, export data and perform configuration management.

2.2.2 Standby Mode

If the user does not have any operation within the set time, the ECG machine automatically enters the standby state.

Follow the steps below to set the time to automatically enter standby mode:

1. Directly click “” on Waveform Acquisition interface to enter into Application interface;



2. Click “” to enter the general interface;
3. Select “machine setting”;
4. The display turns off in standby mode, which reduces power consumption and extends display life. Click the “” button in the upper left corner of the screen to return to the system application interface, and continue to click the “” button to return to the normal mode.

2.2.3 Demo Mode

In this mode, the ECG can demonstrate features of the machine without connecting the patient cable and accessories. Click the "☰" button in the upper left corner to enter [Record Analyze] → [Demo Mode] → [ON], click the "✕" button to return to enter the demo mode.




Fig. 2.6 Demo mode

In the demo mode, there are demonstration waveforms for demonstration, which can be printed. Click the "☰" button in the upper left corner to enter [Record Analyze] → Demo Mode → Off. Click the "✕" button to exit the demo mode.



WARNING

- ◆ *Demo function is mainly used to display the machine performance and provide training for users. When the machine is connected to the patient in clinical practice, it's prohibited to use the Demo function in case that the medical staff shall mistake the demo waveform as the patient's, thus affecting the patient's measurement and delaying his or her treatment. Before use, the user must check the unit, its cables and accessories in order to make sure that all of them*

will work safely and properly.



CAUTION

- ◆ Once entering into the Demo Mode, the system cannot exit automatically. Even the ECG machine is Reboot after shutdown, it is still in demo mode. You need to click the "" button in the upper left corner to enter [Record Analyze] to close the demo mode.
-

2.3 Interface Display

2.3.1 Waveform Acquisition Interface

There is locking function in this interface. After pressing the keyboard, press the key combination "CTRL+ALT+L". A lock dialog box will pop up, which can lock the paper feed speed, gain, layout and report printing functions. After locking, the option button is gray; if you need to change the option settings, press the key combination "CTRL+ALT+L" again to unlock the option.




Fig. 2.7 Waveform Acquisition Interface

1 Patient and system information zone :

- a  Mode setting, press it to quickly enter the mode setting dialog window;
- b  Patient ID;
- c  Heart Rate Indication (tachycardia/bradycardia); Click it to quick enter into Heart Rate Indication menu;
- d  External printer indication. The icon being dark indicates that it is not connected;
- e  Lead indication: light-on icon indicates that the patient cable is wellconnected. Click it to quick enter into Lead Indication menu;
- f  Storage information; Click it to quick check the storage device;
- g  LAN connection state; Dark icon indicates network disconnection, click it to set LAN connection;
- h  WIFI connection state; Dark icon indicates network disconnection, click it to set WIFI connection;
- i  Battery Information. The white bar indicates the battery level; Click it to enter Power Management menu;
- j  Time display; Click it and quick set the date and time;

18.3

2 Waveform zone: Displaying ECG waveforms. Red lead indicates the lead falls off, while white indicates the lead is well connected;

3 Hot keys zone: Displaying real-time function keys, which respectively correspond to F1~F7 keys. Each key's function is listed in the table below:

Name	Key	Function Illustration
F1		Enter into Patient Info interface
F2		Change the paper speed
F3		Change the gain amplitude
F4		Change the lead layout
F5		Freeze to generate the report (for details, please see <i>Chapter 6</i>)
F6		Enter into Report Management interface
F7		Enter into Application interface

2.3.2 “Standard ECG”, “Record Analyze”, “Lead System” Interface

On Waveform Acquisition interface, click “” at upper left of screen to quick enter into “Standard ECG”, “Record Analyze”, or “Lead System” interface.

2.3.2.1 Standard ECG




Fig. 2.8 Standard ECG Interface

Standard ECG			
Menu	Drop-down Menu	Defaults	Illustration
Sampling mode	Real-time Pre-sampling Trigger Sampling Period Sampling	Real-time	Real-time: real-time acquisition and measurement of ECG signals; Pre-sampling: Start acquisition several seconds before you decide the start point of recording. If the time of recording signals before the start point is insufficient, the signals can be supplemented by those acquired after the start point. Trigger sampling: Once triggered by preset

Standard ECG			
Menu	Drop-down Menu	Defaults	Illustration
			<p>conditions, start acquisition 2 seconds before the start point.</p> <p>Period sampling: sampling in periodic intervals</p>
Real-time sampling	Auto-simultaneous Auto-sequential Manual Rhythm	Auto-simultaneous	<p>Auto: When recording the ECG waveform, the system automatically records each lead according to the sampling time and automatically switches the record;</p> <p>Simultaneous: record the ECG waveform of the 12-lead at the same time;</p> <p>Sequential: The 12-lead is divided into 4 averaging periods according to the record format layout (such as 3×4), and recorded in the lead sequence. The lead recording time of each column is set by the "printing duration" item;</p> <p>Manual: Manually change the recording when recording ECG waveform. Manual recording is only available in real time sampling;</p> <p>External printer is not supported in manual mode.;</p> <p>Only sequential recording is available in manual mode;</p> <p>Rhythm: Record a single rhythm lead waveform.</p>
Real-time rhythm analysis	1-5 min	1min	Set the single rhythm lead sampling time.
Periodic sampling	3-100min	12min	Set the total sampling time.

Standard ECG			
Menu	Drop-down Menu	Defaults	Illustration
total time			
Periodic sampling interval	3-100min	3min	Set the sampling time for each cycle
Printing duration	2.5 sec、 5 sec 、 7 sec、 10 sec、 Two pages、 Three pages、 Four pages	Two pages	The time of each lead being recorded when the recording mode is auto-simultaneous, auto-sequential
Freezing time	20 sec、 30 sec、 60 sec	20 sec	The ECG signal waveform is frozen 20 seconds before the button is pressed to observe the measurement result. If the time is not enough, only take the effective time.



CAUTION

- ◆ *When the sampling mode is pre-sampling, triggering sampling or periodic sampling, the recording mode defaults to auto-simultaneous with no auto-sequential option.*
-

2.3.2.2 Record Analyze




Fig. 2.9 Record Analyze Interface

Record Analyze			
Menu Item	Drop-down Menu	Defaults	Illustration
Demo Mode	Off / On	Off	Turn on/off the Demo mode.
Auto Diagnosis	Off / On	On	Turn on/off the Auto Diagnosis After turning off, the ECG report does not show the diagnosis conclusion
Auto Measurement	Off / On	On	Turn on/off the Auto Measurement. After turning off, the ECG report does not show the diagnosis parameters.
Minnesota Code	Off / On	On	Turn on/off the Minnesota Code. If it is turned off or can not acquire the code, the ECG report does not show the Minnesota Code
PACE Detection	Off / On	On	Turn Off / On the PACE Detection

2.3.2.3 Lead System

Lead system can freely change among Standard, Posterior Wall, Right Chest, Right Chest Posterior Wall, Previous Intercostal Space, Next Intercostal Space, CABRERA, and Custom.




Fig. 2.10 Lead System Interface

Lead System	
Menu Item	Lead Label
Standard	I , II , III , aVR , aVL , aVF , V1 , V2 , V3 , V4 , V5 , V6
Posterior Wall	I , II , III , aVR , aVL , aVF , V1 , V2 , V3 , V7 , V8 , V9
Right Chest	I , II , III , aVR , aVL , aVF , V1 , V2 , V3 , V3R , V4R , V5R
Right Chest Posterior Wall	I , II , III , aVR , aVL , aVF , V3R , V4R , V5R , V7 , V8 , V9
Previous Intercostal Space	I , II , III , aVR , aVL , aVF , V`1 , V`2 , V`3 , V`4 , V`5 , V`6
Next Intercostal Space	I , II , III , aVR , aVL , aVF , V.1 , V.2 , V.3 , V.4 , V.5 , V.6
CABRERA	aVL , I , -aVR , II , aVF , III , V1 , V2 , V3 , V4 , V5 , V6
Custom	I , II , III , aVR , aVL , aVF , V1 , V2 , V3 , V4 , V5 , V6

2.3.3 Patient Info Interface

In the waveform acquisition interface, you can click **F1** key on panel or the “ **Patient Info** ” or “  ” button at the bottom of the touch screen to enter the patient information interface and then enter the information by keyboard. After the input is completed, click the “ **Confirm** ” button to save and return. Click the “ **Previous** ” button to automatically generate the last patient information.

8.5

8.5




Fig. 2.11 Patient Info Interface

Click the “  ” button in the Patient information interface, enter the server's IP address, port and path (for details, please consult your network administrator). You can synchronize the patient information data in the server, and directly import the patient information into the list. Click on one of the patients to automatically generate patient information. Before importing data, please test whether it can be connected. After the test is successful, click the “Save” button, and then click the “Sync” button.

If the automatic list synchronization function is enabled, the system will automatically synchronize the server data each time you enter the list interface;

If the function of the ID card module is enabled, and the resident ID card reading machine can be connected to the electrocardiograph through the USB interface, the patient ID is automatically generated by brushing the patient ID card.




Fig 2.12 Patient Information Interface - Settings




Fig 2.13 Patient Information Interface - List

2.3.4 Paper Speed Interface

On Waveform Acquisition interface, touch the key **25mm/s** or click the F2 key on the panel to change the paper speed. Paper speed can be selected among 5mm/s , 6.25mm/s , 10 mm/s , 12.5mm/s , 25mm/s , 50mm/s.



Fig. 2.14 Paper Speed Interface

2.3.5 Gain Amplitude Interface

On Waveform Acquisition interface, touch the key **10mm/mv** or click the F3 key on the panel to change the gain amplitude. Gain amplitude can be selected among 2.5mm/mv, 5mm/mv, 10mm/mv, 20 mm/mv, 40 mm/mv and Auto.




Fig. 2.15 Gain Amplitude Interface

2.3.6 Lead Layout Interface

In the waveform acquisition interface, you can directly click the “6x2” button on the touch screen or click the F4 key on the panel to switch the lead layout. Lead layout options are: 3x4 , 3x4+1 , 3x4+3 , 6x2 , 6x2+1 , 12x1 , 6x2+3. Take 3x4+1 for example, When printing, the 12-lead ECG waveform is arranged in a 3-row, 4-column format with a rhythm lead below when printing. (The lead layout will be saved after shutdown.)




Fig. 2.16 Lead Layout Interface

2.3.7 Report Management Interface

In the waveform acquisition interface directly clicks the "Report" button on the touch screen or click the F6 key on the panel to enter the report management interface. See Chapter 7, Report Management for details. This interface has locking function. After the external keyboard is pressed, press the key combination "CTRL+ALT+L", a lock dialog box will pop up, and the report printing, report transmission and report deletion functions can be locked. After locking, the option button is grayed out; if you need to change the option settings, press the key combination "CTRL+ALT+L" again to unlock the option.




Fig. 2.17 Report Management Interface

2.3.8 Application Interface

On Waveform Acquisition interface, directly touch the key  or click the F7 key on the panel to enter the Application interface.




Fig. 2.18 Application Interface

Chapter 3 Installation

3.1 Preparation for Installation



CAUTION

- ◆ *This device shall be installed by the personnel authorized by the manufacturer.*
 - ◆ *Don't open the device cover. Otherwise there might be risk of electrical shock. Only the service personnel authorized and trained by the manufacturer can maintain or update the device.*
 - ◆ *This device includes software protected by copyright and international laws, all right reserved to the manufacturer, No part of the device may be modified, reproduced or transmitted in any form or by any means, without prior written permission from the manufacturer.*
 - ◆ *All the analog digital devices connected to this device must be approved by the designated standard (such as IEC60601-1 Safety of Medical Electrical Equipment). And all equipment must be connected in accordance with the valid version of the system standard IEC60601-1-1. The person who is in charge of connecting the additional devices to the input and output signal ports shall be responsible for whether the system conforms to the standard IEC60601-1-1 or not. For any questions, please contact the manufacturer.*
 - ◆ *When this device and another electrical device are connected into a conjunction with certain function, if it cannot be determined whether this conjunction is hazardous (for example, electrical shock caused by leakage current crowding) or not in terms of each device's specification, please contact the manufacturer or related experts in hospitals in order to guarantee that safety of all devices will not be breached.*
-

3.1.1 Open-package Inspection

Before opening the package box, carefully inspect it. If any damage is found, please contact the shipping company immediately.

Please open the package in a correct way. Carefully remove the device and other components from the box and check them one by one per the packing list. Check whether the device is damaged mechanically, or the goods are complete. For any questions, please contact the manufacturer immediately.



WARNING

- ◆ *Keep packaging materials out of reach of children. When dispose of the packaging materials, please follow local laws and regulations or regime of medical waste disposal in local hospitals.*
 - ◆ *The device may be contaminated by microbes during storage, transportation and usage. Please verify that the package is intact before use, especially disposable accessories. If any damage is found, please stop using.*
-



CAUTION

- ◆ *We'll keep the packaging boxes and materials for future shipment or storage.*
-

3.1.2 Environment Requirements

This device must be used in an environment complying with environment specifications in this manual.



- ◆ *Reasonably avoid using the device in the presence of noise, vibration, dust, corrosive or flammable and explosive substances. If the device is installed in a box, please make sure that the front and back space is enough for operation, maintenance and service. In order to allow unimpeded air circulation for a good cooling effect, at least 5cm space should be saved around the device*
 - ◆ *During the process of device movement from one environment to another, it might cause condensation due to differences in temperature or humidity. At this moment, you can use it until the condensation disappears.*
-



WARNING

- ◆ *Please ensure that the device works under required specified environment. Otherwise it will not comply with the technical specifications alleged in this manual, which may lead to unpredictable consequences (damaging the device, etc).*
 - ◆ *Do not use the device in oxygen-rich environment, or the presence of flammable or explosive substances (anesthetics, etc) in case of fire or explosion.*
 - ◆ *Electromagnetic fields can affect the performance of this unit. Therefore other devices used in the vicinity of the unit must comply with EMC requirements. Mobile phone, X-ray or MRI equipment are potential interference sources, because they could emit electromagnetic radiation of high intensity.*
 - ◆ *Power plug is used to separate the ECG circuit from power mains. Do not put the ECG machine in a place where is difficult to handle the plug.*
 - ◆ *Be sure to connect the AC power cord to a hospital-grade three-core socket with a ground wire to guarantee reliable grounding*
 - ◆ *Before the unit is connected to the AC power, please verify that the power's voltage and frequency comply with its label or the requirements specified in this manual.*
-



CAUTION

- ◆ *If strong electromagnetic radiation exists in surroundings, it will produce different levels of interference to the ECG machine. Please make sure there are no high-voltage lines and heavy-load power cables passing-by near the unit and patient bed.*
 - ◆ *While examining the patient, prevent irrelevant individuals from contacting the machine or the patient in case that interference affects the interpretations.*
-

3.2 Power Selection

This ECG machine can work by 220V AC power or lithium-ion rechargeable battery power.

3.2.1 Connecting to the AC Power

Plug one end of accompanying three-core power cord into the power jack at the back of the machine, and the other end into the three-core socket with a grounding cable. Turn on the power switch at the back of the machine, then the AC input light is on, which indicates that AC power has been connected.



WARNING

- ◆ *Use the dedicated adapter power cord provided by the manufacturer only. If the power cord is damaged, please contact the manufacturer to buy a new one for replacement.*
 - ◆ *If proper grounding cannot be guaranteed, you shall operate the device with the built-in rechargeable battery. Otherwise it may incur electrical shock to the patient and operator.*
-

3.2.2 Battery Power

The ECG machine has a built-in rechargeable lithium battery pack that can be used to power the unit during transport or when AC power is not available. For battery use and maintenance, please refer to the related contents in *Section 9.3* of this manual.



WARNING

- ◆ *Make sure that the ECG machine is powered by dedicated rechargeable battery. Before use, please refer to the contents in *Section 9.3* of this manual. Safe and proper use of the battery shall be guaranteed to prevent current leakage, heat or explosion.*
 - ◆ *Battery replacement shall be carried out by the manufacturer authorized service engineer. For battery replacement, please contact the service engineer*
-

authorized by the manufacturer.



CAUTION

- ◆ *To prevent data loss caused by accidental AC power interruption, a battery must always be installed in ECG machine.*
 - ◆ *Whenever the unit is connected to AC power and AC power is on, the battery is being charged. Therefore, it is recommended that the unit remain connected to AC power when not in use. This will ensure a fully-charged battery whenever it is needed.*
-


3.3 Installing the Recording Paper

This ECG machine is designed to use Zip fold thermal recording paper of 216 mm or 210 mm in width.

3.3.1 Installing the Paper Block

Prior to install the 216mm-wide paper, please install the accompanying paper block first.

A: Push the printer door switch, and pull out the door, then put in the paper block as shown.






Fig. 3.1 Sketch of Installing the Recording Paper

3.3.2 Installation of recording paper

The method of installing the 216mm and 210mm wide Z-fold thermal recording paper is the same as the following figure:




Fig 3.2 Schematic of Installing the Recording Paper



- ◆ The black-marked end of recording paper should stay close to the black-marked detecting end of the recorder.
 - ◆ The grid side of recording paper should face up to the printer head.
 - ◆ The end page of recording paper is marked red. Make sure it is placed on the bottom.
 - ◆ Recording paper should be placed above the plastic pull tab.
 - ◆ Place the recording paper naturally while installing. Keep the lower edge of recording paper close to that of the printer door. Don't center it deliberately.
 - ◆ When close the printer door, parallel the vertical grid line of recording paper to the tearing margin as possible in case of the paper jam caused by huge placement deviation while printing.
 - ◆ Please use the paper recommended by the manufacturer.
-

3.4 Connecting the Patient Cable

Properly connect the patient cable to the lead connector at the right side of the machine.

Tighten the screws on patient cable connector and attach the cable to the ECG machine.

The other end of the cable shall be connected to the patient via electrodes. For more details, please refer to *Chapter 5*.



CAUTION

- ◆ *Please use the dedicated patient cable configured by the manufacturer. If the cable is damaged, please contact the manufacturer in time to replace or purchase a new one.*
-

3.5 Power on

3.5.1 Checks before Power on

In order to ensure safe exam and stable printout of ECG, you are required to make the above-mentioned checks before operating of the ECG machine.

- **Checks on Operation Environment**

Verify that the ground wire is well connected, the ground bolt is tight, and ground wire connector and grounding wire are properly connected.

Operation environment of the machine should be free of X-ray equipment, short-wave devices, and the like, which may impose interference on the ECG machine.

The machine shall be operated in warm indoors (room temperature should be no less than 18°C) in order to avoid myoelectric interference caused by coldness.

Verify that the power cord is properly connected, and disentangled with other cables.

- **Checks on Power Supply**

When the machine is to operate on AC power, please check whether the power voltage is the same as local voltage used and whether the power cord is firmly connected with the machine. Please use a properly grounded AC outlet. If a battery is used, please inspect whether it is fully charged.

- **Checks on Patient Cable**

Check whether the patient cable is firmly connected with the ECG machine.

Verify that the patient cable plugs are correctly and reliable connected with the related electrodes.

- **Checks on Recording Paper**

Verify that recording paper is sufficient and properly installed.

3.5.2 Power on

After installation and checks, connect the power cable and turn on the ECG machine, then start to acquire and record the patient's ECG.



CAUTION

- ◆ *If the device is damaged or does not work, it cannot be used to acquire and record the patient ECG. Please contact the service personnel or the manufacturer immediately.*
-

3.6 ECG Machine Setting

There are some settings need to be done for the first time use, such as time, report storage ways etc. During ECG waveform acquisition and measurement, you need to set such parameters as recording modes, sampling time, gain, paper speed, filter etc. For more details, please refer to *Section 2.3* and *Chapter 4*.

3.7 Power off

Please follow the steps below to power off the ECG machine:

- 1) Verify that the patient ECG acquisition and recording can be ended.
- 2) Disconnecting all the ECG electrodes with the patient.
- 3) Press the power switch, a dialog appears with three options, then select "Power Off".



CAUTION

- ◆ *If you cannot shut down the machine normally, or special circumstances arise, please long press the power switch for 10 seconds to force a shutdown. Forced shutdown may incur ECG data loss, thus usually not recommended except in special circumstances.*
-

Chapter 4 System Applications

4.1 Entering Main Menu

In the waveform acquisition interface, click the “  ” button on the touch screen or **F7** key on panel to enter into the system application interface. The system application interface layout is described in Section 2.3.8. The touch screen allows you to directly click on the item you want to view.

Key operation mode:

- Tab key to make the focus switch to the next control, Shift+Tab to make the focus switch to the next control;
- In the control options, you can press the “↑” and “↓” direction keys to switch vertically. The “←” “→” direction keys switch horizontally, and the space bar confirms the selection;
- You can touch the screen to select; or connect an external USB mouse and left click to select.

The settings made in the system settings interface will be saved as default user settings and will remain active the next time you turn the phone on.

The general setting interface has a locking function. Pressing the key combination “CTRL+ALT+L” will pop up the lock dialog box, which can lock the machine settings, report settings, recorder settings, filter settings and heartbeat warning settings. After locking, the option button is grayed out; if you need to change the option settings, press the key combination "CTRL+ALT+L" again to unlock the option.

4.2 Application Interface

4.2.1 General–Machine Setting




Fig. 4.1 Machine Setting Interface

Machine Setting			
Menu Item	Option	Defaults	Illustration
Brightness	/	/	Adjust the backlight intensity of the display.

Standby	Off, 5 min, 10 min, 15 min, 20 min, 25 min, 30 min	Off	Set the time when the ECG machine automatically enters standby mode. If there is no user operation within the set standby time, the ECG machine enters the standby state. When entering standby mode, the screen is turned off. The standby time you set cannot be longer than the automatic shutdown time.
Shut down	Off, 5 min, 10 min, 15 min, 20 min, 25 min, 30 min	Off	Set the auto power off time. If the user does not perform any operation during the set auto power off time, the ECG will automatically shut down.
Mouse	Show/Hide	Hide	When "Display" is turned on, the screen displays the current position of the mouse.
Language	Simplified Chinese, English,etc.	Chinese	Select the interface language.
Date and time	/	/	Set the date and time of the machine.
Factory settings	Restore factory settings	/	To restore factory settings
Burn verification	Burn verification (MD5 file) / system upgrade	/	Burn verification is used to verify that the current version is successfully burned.



CAUTION

◆ *The factory settings will be restored and the machine will restart automatically.*

4.2.2 General-Report Setting




Fig. 4.2 General-Report setting

Report setting			
Menu Item	Option	Defaults	Illustration
Storage location	Local/ SD / USB	Local	See section 6.4.1 for details.
Storage method	Auto/Manual	Auto	See section 6.4.2 for details.
Report template	Waveform report / Measurement matrix / Average template	Waveform report	<p>Check the waveform and display the ECG waveform on the printed report.</p> <p>Check the measurement matrix and display the ECG measurement matrix on the printed report (when printing, multiple periodic waveform of each lead will be combined into a single cycle waveform.)</p> <p>Check the average template and display the ECG average template on the printed report (when</p>

			<p>printing, the horizontal display shows 12 leads; the vertical display shows the parameters of each lead, such as P, QRST wave start and end points, P, QRS, T wave group period, etc.)</p> <p>Built-in print report in real-time sampling mode, only supports waveform report; other sampling mode built-in printing or ECG preview interface printing regular ECG report can support all report templates, you can select multiple or single.</p> <p>The report template for regular ECG external printing supports single or multiple selection by printer model. (See 4.2.3 General Settings - Recorder Settings)</p>
Local/USB/SD Card	/	/	Display storage capacity

4.2.3 General-Recorder Setting




Fig. 4.3 General-Recorder setting

Recorder Setting			
Menu Item	Options	Default Value	Illustration
Printer	Built-in/external/paperless	Built-in	In paperless mode, click  key, the system will acquire and save ECG waveforms. No matter there is recording paper in the printer or not, it will not print the waveforms.
Graded gain	Off/On	Off	When set to On, the gain amplitude of the first six leads is twice the gain of the last lead. For example, when the gain amplitude is set to 20 mm/mv, the first six-lead gain amplitude is 20 mm/mv, and the last lead amplitude is 10 mm/mv.

Recorder Setting			
Menu Item	Options	Default Value	Illustration
			If the waveform acquisition interface selects automatic gain, the gain amplitude of each lead is automatically assigned.
Printing Grid	Off/On	On	Set the grid when export PDF reports. When it is set "On", the waveform zone and footer zone will display grid background; When set "Off", no grid will appear on the report. Printing grid is only available for an External printer.
Line width	1/2/3	2	The width of the waveform line displayed on the thermal paper.
Printing	Current page / all / two pages / three pages / four pages	Two pages	Preview the contents of the print when the it is in built-in printer.

Recorder Setting			
Menu Item	Options	Default Value	Illustration
Supported external printer models	/	/	<p>Supported printer models:</p> <p>HP Dest_Jet P1102 (only report template radio is supported)</p> <p>HP Dest_Jet P1108 (only report template radio is supported)</p> <p>HP Dest_Jet P1008 (only report template radio is supported)</p> <p>LaserJet Pro M202n (supports report template multiple selection and single selection)</p> <p>If the report template is multi-selected, the connection only supports external printers that report single-selection, and the waveform report is printed by default.</p>

4.2.4 General-Filter Setting




Fig. 4.4 General-filter setting

Filter setting			
Menu Item	Options	Default Value	Illustration
AC Filter (Hz)	OFF、 50、 60	50	Select whether to turn on the AC filter to filter out AC power interference. AC Filter should always be turned on, unless it is necessary to turn it off.
Baseline drift (Hz)	OFF、 0.01、 0.05、 0.3、 0.6	0.6	Select the baseline drift filter (BDR) frequency. Baseline drift filter rejects most of the baseline drift interference and makes the ST segment undistorted. Baseline drift filtering is not performed when set to [Off]. (This filter is high-pass filtering)

EMG (Hz)	OFF、 25、 35、 45、 75、 100、 150、 250、 350	35	Set the frequency of EMG. Turning on the EMG filter can filter out EMG interference, but it may reduce the bandwidth and cause the ECG waveform shape to change. The EMG filter is low-pass filtering. signals with a frequency higher than the set value will be filtered out. When set to 35 Hz, the system only displays signals below 35 Hz and below 35 Hz, and signals beyond 35 Hz will be filtered out.
------------	---	----	--



CAUTION

- ◆ *The settings of the baseline drift filter and the EMG filter affect the frequency response range of the electrocardiograph. If the baseline drift filter setting value is too high and the EMG filter setting value is too low, it will affect the accuracy of the ECG machine for waveform restoration. For example, when the baseline drift filter frequency is set to 0.6 Hz and the EMG filter frequency is set to 25 Hz, the system frequency response range is 0.6 Hz to 25 Hz.*
-

4.2.5 General-Heartbeat Warning




Fig. 4.5 Heartbeat Warning Interface

Heartbeat Warning			
Menu Item	Options	Default Value	Illustration
Bradycardia	40-80	60	Set the bradycardia threshold.
Tachycardia	81-140	100	Set the tachycardia threshold.

4.2.6 General-Institution Information

The screenshot displays the 'Settings' application interface. At the top, there is a navigation bar with a '< Back' button, the title 'Settings', and system status icons (Bluetooth, USB, Wi-Fi, battery) along with the date '2019-08-01' and time '09:00:09'. On the left side, a vertical menu lists various settings: 'Machine setting', 'Report setting', 'Recorder setting', 'Filter setting', 'Heart beat warning', 'Institution information' (highlighted in blue), 'About the machine', 'Transfer setting', and 'iHE DICOM settings'. The main content area on the right is titled 'Institution information' and contains a form with the following fields:

SN code	8231AB6002869
Device model	iMAC 120pro
Device code	<input type="text"/>
Institution name	<input type="text"/>
Institution code	<input type="text"/>
Dept. name	<input type="text"/>
Dept. code	<input type="text"/>

Fig. 4.6 Institution Information Interface

4.2.7 General-About the Machine




Fig. 4.7 About the Machine Interface

4.2.8 General-Transfer Setting

Enter the transfer setting interface. Enter the unlock password according to the prompt. The password defaults to 4000400499.




Fig. 4.8 Transfer Setting Setting Interface

Transfer Setting			
Menu Item	Options	Default Value	Illustration
Report transmission format	ZQECG、PDF、PNG、XML、HL7、DICOM-WAVEFORM、DICOM-IMAGE、DICOM-PDF GDT	ZQECG	The transmission format of the ECG file.
Naming rule	SN/PID/StudyID/Time + XXX + SN/PID/StudyID/Time + XXX + SN/PID/StudyID	SN+Time. ZQECG	The naming rules can be a combination of SN, PID, StudyID, Time, or can be named separately. The file name suffix is automatically generated based

	/Time .		on the selected report transfer format.
Auto transfer	OFF/FTP/SAMBA /HTTP/DICOM	OFF	<p>When it is OFF, you can manually select the ECG report transmission on the [Report Management] interface.</p> <p>When it is on, the report can be transmitted in real time in the FTP/SAMBA/HTTP/DICOM transmission mode when the ECG report is printed and frozen. (You should ensure that the network is properly connected and the FTP/SAMBA/HTTP/DICOM transmission mode is successfully enabled)</p> <p>Note: The corresponding button will only be highlighted if the FTP/SAMBA/HTTP/DICOM transfer mode is enabled, otherwise it will be grayed out.</p>
FTP setting	IP address, port number, username, password, directory	/	<p>1) The FTP username and password that you set must be the username and password that can be logged into the FTP server.</p> <p>2) The FTP path that you set must be a subdirectory that already exists in the root directory of the FTP server.</p> <p>Note: For more information about the FTP server, please consult your network administrator.</p> <p>3) After setting, you can click the “” button to check if</p>

			<p>it can be connected.</p> <p>4) After the test is passed, check the “” button to enable it.</p> <p>(Note: It can only be used in the transmission of the [Report Management] interface after it is enabled.)</p>
HTTP setting	IP address, port number, username, password, directory	/	<p>1) The HTTP username and password that you set must be the username and password that can be logged into the FTP server.</p> <p>2) The HTTP path that you set must be a subdirectory that already exists in the root directory of the HTTP server.</p> <p>Note: For more information about the HTTP server, please consult your network administrator.</p> <p>3) After setting, you can click the “” button to check if it can be connected.</p> <p>4) After the test is passed, check the “” button to enable it.</p> <p>(Note: It can only be used in the transmission of the [Report Management] interface after it is enabled.)</p>
Remote diagnosis settings	IP address, port number, username, password, directory	/	<p>1) The remote diagnosis username and password that you set must be the username and password that can be logged into the FTP server.</p>

			<p>2) The remote diagnosis path that you set must be a subdirectory that already exists in the root directory of the remote diagnosis server.</p> <p>Note: For more information about the remote diagnosis server, please consult your network administrator.</p> <p>3) After setting, you can click the  button to check if it can be connected.</p> <p>4) After the test is passed, check the  button to enable it.</p> <p>(Note: It can only be used in the transmission of the [Naming Rule] interface after it is enabled.)</p>
SAMBA setting	IP address, port number, username, password, directory	/	<p>1) The SAMBA username and password that you set must be the username and password that can be logged into the FTP server.</p> <p>2) The SAMBA path that you set must be a subdirectory that already exists in the root directory of the SAMBA server.</p> <p>Note: For more information about the SAMBA server, please consult your network administrator.</p> <p>3) After setting, you can click the  button to check if it can be connected.</p> <p>4) After the test is passed, check</p>

			<p>the “” button to enable it. (Note: It can only be used in the transmission of the [Report Management] interface after it is enabled.)</p>
--	--	--	--

4.2.9 General-DICOM Setting

Enter the transfer setting interface. Enter the unlock password according to the prompt.

The password defaults to 4000400499.

10.2, 10.3





Fig. 4.9 DICOM Setting Interface

DICOM SETTING			
Menu Item	Options	Default Value	Illustration
Dicom Storage	IP address, remote port number, remote AE	/	The ECG file is uploaded to this server in the "dicom" transfer

(upload configuration)	directory, local AE directory		format. 1) After setting, you can click the “  ” button to check if it can be connected. 2) After the test is passed, check the “  ” button to enable it. (Note: It can only be used in the transmission of the [Report Management] interface after it is enabled.)
Worklist (Server settings)	IP address, remote port number, remote AE directory, local AE directory	/	The ECG machine imports patient information data from this server. 1) After setting, you can click the “  ” button to check if it can be connected. 2) After the test is passed, check the “  ” button to enable it.
MPPS (Server settings)	IP address, remote port number, remote AE directory, local AE directory	/	The ECG machine will send three states of acquisition, start to collect, complete the acquisition, and cancel the acquisition, to this server. 1) After setting, you can click the “  ” button to check if it can be connected. 2) After the test is passed, check the “  ” button to enable it.
Storage Commit (Storage server)	IP address, remote port number, remote AE directory, local AE directory	/	A server that confirms whether the ECG file has been uploaded. 1) After setting, you can click the “  ” button to check if it can be connected.

configuration)			2) After the test is passed, check the "  " button to enable it.
Character set	ASCII/GB1 8030/ISO-8859-1~ ISO-8859-9/ UTF-8	ASCII	Select the appropriate character set to avoid garbled characters in the transmitted ECG file.

4.2.10 E-Mail setting



1. Click the "  " button to open the mail dialog. Note: For more information on configuring your network and interfaces, please consult your network administrator. Use wired Ethernet or WIFI transmission:

- 1) Manually add the sender, recipient's account number, server address and port number.
- 2) Select and confirm the appropriate sender, recipient.

3) Click the "  " button to check the files you want to upload. The system will automatically fill in the email subject.

4) Click the "  " button to send the mail.

2. The recipient can browse the png format ECG file through the mailbox.




Fig 4.10 Email Interface

4.2.11 LAN setting



Click the "LAN" button to open the Wired Settings dialog box to make the appropriate network settings. (Note: For details on configuring the network, please consult your network administrator.) After setting, click the "Save" button to save and take effect. You can automatically obtain the local IP address or manually enter the local IP address.




Fig. 4.11 LAN Interface

4.2.12 WIFI setting



Click the "WIFI" button to open the Wireless Settings dialog box, you can choose to turn WiFi on or off. Click the "Add internet" button to manually connect hidden WiFi and support WEP encryption type hotspots.




Fig. 4.12 WIFI Interface

4.2.13 Standard ECG setting



Click the "Standard ECG" button to enter the waveform acquisition interface.

4.2.13.1 Auto mode

The automatic mode refers to the automatic acquisition and printing of waveforms by the electrocardiograph. It is a common mode of the electrocardiograph and is used for routine electrocardiogram examination. The user simply presses the "ECG" button to automatically capture and print the waveform. (Does not include pre-sampling mode).



Specific operation method:

1. In the waveform acquisition interface, set the layout, gain, and paper speed according to actual needs;
2. Click the  button in the upper left corner of the waveform acquisition interface to enter the [Standard ECG] setting interface, select the auto-simultaneous or auto-sequential recording mode, set the sampling mode, acquisition duration, and printing duration;
3. Enter the [12-lead] interface and select the lead system and the single or three-rhythm lead label;
4. According to your own needs, set other parameters, and exit the [Standard ECG] interface after the setting is completed.
5. Click the  button on the panel to print the automatic mode ECG report.

 **CAUTION**

- ◆ *When the ECG signal is just connected or the device receives the overload noise, the waveform will be disordered, the baseline drift is serious, and the signal waveform amplitude may exceed the maximum display width. At this time, wait for the device connection and the patient as well as the waveform to stabilize. Then start measuring and recording.*
 - ◆ *When the ECG is overloaded or any part of the amplifier is saturated, the ECG machine is in an abnormal working state. At this time, the interface only displays the baseline. In order to obtain accurate measurement results, wait for the device connection and the patient as well as the waveform to stabilize. Then start measuring and recording.*
 - ◆ *If the waveform becomes cluttered or unstable during the patient's ECG signal acquisition, see Chapter 8.*
-

4.2.13.2 Manual mode

Manual mode means that the user can manually control the duration of the ECG machine to collect and print the ECG. It is generally used by the user to collect and print ECG waveforms of any length according to clinical needs.

Specific operation method:

1. Click the “” button in the upper left corner of the waveform acquisition interface to enter the [Standard ECG] setting interface, select the real-time sampling mode and manual recording mode;
2. Enter the [12-lead] interface and select the lead system and the single the single or three-rhythm lead label;
3. In the waveform acquisition interface, set the layout, gain, and paper speed according to actual needs;
4. According to your needs, set other parameters, and exit the [Standard ECG] setting interface after the setting is completed.
5. Click the “” button on the panel to print the manual mode ECG report.
6. The built-in /paperless can be switched through the floating button in the lower left corner of the interface. Note: External printer is not supported in manual mode.




Fig. 4.13 Manual mode waveform acquisition Interface

4.2.13.3 Rhythm mode

The rhythm mode is for user to collect and print a single lead for a long time to observe and capture sporadic or frequent arrhythmia.

Specific operation method:

1. Click the  button in the upper left corner of the waveform acquisition interface to enter the [Standard ECG] setting interface, and select the real-time sampling mode, rhythm recording mode and rhythm analysis time;

2. In the waveform acquisition interface directly touche the screen to select the rhythm lead label.

3. After setting, click the  button on the screen, the sampling time countdown will appear. A complete rhythm waveform will be recorded after the countdown ends.

During the recording process, press the  button to stop recording if necessary.

4. Click the  button on the panel, and the sampling time countdown will appear. A complete rhythm waveform will be printed after the countdown ends. During the

recording process, press the  button to stop recording or stop printing the rhythm analysis report if necessary.




Fig. 4.14 Rhythm Analysis waveform acquisition Interface




Fig. 4.15 2 minute rhythm analysis report Interface

4.3 Machine Setting

4.3.1 First time setup

When using the ECG for the first time or use it after repair or update, be sure to check and improve the following settings:

Fill in the relevant device information in [Institution Information] for the convenience of centralized management and maintenance.

Set the date and time in [Machine Settings] to ensure accurate date and time for the convenience of patient number generation and report management.

Set the backlight brightness in [Machine Setting].

Set the standby time in [Machine Setting]

Set the language type in [Machine Setting].

Set the report save location and save mode in [Report Setting].



- ◆ *When the equipment is repaired or upgraded, the equipment settings will usually be restored to the factory values and need to be reset at that time.*
-

4.3.2 Setup before Use

When using the ECG for the first time or use it after repair or update, be sure to check and improve the following settings:

In the [12-lead] interface , select the lead system and the single or three rhythm lead label;

In the [General-Filter Setting] setting interface, set the filter frequency;

In the waveform acquisition interface, set the layout, gain, and paper speed according to actual needs;

In the [Standard ECG] setting interface, select sampling mode, printing time, freezing time and other parameters according to actual needs.



- ◆ *In the measurement, the paper speed, gain, layout, etc. can be quickly set by the real-time shortcut key. For details, see section 2.3.*
-

Chapter 5 Connecting ECG Cable



- ◆ *During defibrillation, don't touch the patient, electrodes, patient cable and lead terminals. Otherwise it may result in serious injury or death.*
 - ◆ *For a pacemaker patient, ECG machine may interpret and record the pacemaker pulse as the QRS complex wave. Please carefully inspect the ECG waveform recorded by ECG machine.*
 - ◆ *For a pacemaker patient, pacemaker detection should be activated when set the machine. Please refer to Section 2.3.2.2 for details.*
 - ◆ *Please verify that all electrodes are connected to the correct sites on patient body. Avoid electrodes (including neutral electrodes) and the patient contacting the ground or any other electric conductors.*
 - ◆ *As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.*
 - ◆ *Suction ball of chest electrodes contains natural rubber, which may cause allergy. Please pay close attention to the skin placed with electrodes, if allergy occurs, please change other types of electrodes.*
 - ◆ *Automatic measurement and diagnosis are for the doctors' reference only, and cannot be directly used as the basis for clinical treatment.*
-

5.1 Environmental Requirements

- It is required to keep warm indoors in winter, and the room temperature is no lower than 18 °C so as to avoid EMG interference caused by coldness. In summer, open the air conditioner to control room temperature and to prevent insecure placement of electrodes cause by the patient's sweating.
- AC power outlet must be connected with dedicated and reliable ground wire. .
- When place the ECG machine, try to keep its power cord away from the patient's bed and patient cable. Do not place other appliances or power cord near the bed.

- The patient's bed should be of suitable size to ensure that the patient can naturally lie down, with hands and feet stretching naturally.

5.2 Preparation

To acquire accurate ECG signals, the patient should be explained about the following information.

- Before ECG examination, the patient should not do strenuous exercise, drink alcohol, smoke, have full meal, have tea, have uncooked and cold food.
- Inform the patient of preparation for ECG examination, and make a good explanation to him or her so as to eliminate his or her mental nervousness.
- Before examination, the patient should have a good rest and stay calm.
- During examination, the patient should lie down naturally, relax, and breathe calmly and evenly.
- During examination, the patient should not move or turn over casually, with hands and feet no touching metal objects, such as metals at the bed edges.
- Open the ECG machine to warm it up, and then set its parameters. For startup and system setup, please refer to *Chapter 3* and *Chapter 4*. Power off the ECG machine when connecting the patient cable and electrodes.

5.3 Electrode Selection and Usage

5.3.1 Examinee's Skin Preparation

When environmental factors are correct and the examinee is ready for the exam, place electrodes on examinee body. To obtain accurate ECG signals, properly prepare the examinee's skin to be placed with electrodes.

Remove body hair on intended sites.

Wipe the intended sites with alcohol to degrease and remove dead skin cells.

Dry each site with a dry cotton ball.

5.3.2 Electrode Selection

The ECG machine is equipped with reusable chest bulbs and limb clamps. For your convenience, you can also purchase disposable electrode slices from the manufacturer. Use chest bulbs or disposable electrode slices to connect the chest leads (C1 ~ C6), while use limb clamps or disposable electrode slices to connect the limb leads (R, F, L, RF).

5.3.2.1 Chest Bulb

The chest bulb is composed of suction bulb and metal electrode. There is a connection hole on metal electrode, which is used to connect the lead wire of $\Phi 4.0$ mm connector.




Figure 5.1 Chest Bulb

1. Suction bulb;
2. Connection hole ($\Phi 4.0$ mm);
3. Metal electrode;

Placing Chest Electrodes

- Check the chest electrodes, and verify that they are clean;
- Respectively connect 6 chest electrodes with its corresponding lead wire. Straighten out the lead wires, and avoid twisting and twining;
- Unsnap the patient's chest button to expose electrode sites;
- Prepare the skin;
- Apply a thin layer of conductive cream on the surface of electrode sites.
- Evenly apply a thin layer of conductive cream on the edges of chest bulbs,;

- Well place the electrodes: squeeze the chest suction bulb to make the lower edge of the metal electrode tightly connected with the skin. Loosen the bulb to make the electrode adsorbed on the patient's skin;
- Ensure that the ECG machine and its cable, electrodes and lead wires are firmly connected.

5.3.2.2 Limb Clamp

The limb clamp is composed of the clamp and metal electrode. There is a connection hole on the metal electrode, which is used to connect the lead wire of $\Phi 4.0$ mm connector.




Fig 5.2 Limb Clamp

1. Connection hole ($\Phi 4.0$ mm);
2. Metal electrode;
3. Clamp

Placing Limb Clamps

Limb electrodes should be placed above the forearm wrist joint and the calf ankle joint medial to make them closely contact with skin, bypassing the bones.

- Check the electrodes to verify that they are clean,
- The four limb electrodes are respectively connected with the corresponding lead wire. Straighten out the lead wires, and avoid twisting,
- Roll up sleeves and trousers to expose electrode sites.
- Prepare the skin,
- Apply a thin layer of conductive cream on the surface of electrode sites;
- Apply a thin layer of conductive cream on metal electrodes,
- Well place the electrodes: clamp the prepared sites by limb clamps;

- Ensure the ECG machine and its patient cable, the electrodes and the lead wires are firmly connected.

5.3.2.3 Disposable Electrode Slice

Disposable electrode slices are not accompanied with this ECG machine. Please contact the manufacturer or authorized representatives if needed.

Steps of placing disposable electrode slices:

- Roll up the patient's sleeves and trousers, and loosen his or her chest button to expose electrode sites,
- Prepare the skin,
- Attach the electrodes in correct position. Limb electrodes should be placed above the forearm wrist joint and the calf ankle joint medial to make them closely contact the skin, bypassing the bones.
- Straighten out the lead wires, and avoid twisting and twining. Connect the lead wires with electrode slices;
- Ensure the ECG machine and its patient cable, the electrodes and lead wires are firmly connected.



- ◆ *Please don't mix up electrodes of different types and brands. Otherwise it may cause large baseline drift or longer baseline recovery after defibrillation.*
 - ◆ *Disposable electrodes can be used only once, repeated use may result in performance degradation or cross infection.*
 - ◆ *Reusable electrodes must be cleaned and disinfected before use.*
 - ◆ *Suction bulb of chest electrode contains natural rubber, which may cause allergy. Please pay close attention to the electrode sites, if allergy occurs, please change another type of electrodes.*
 - ◆ *If reusable electrodes are damaged after long-term use, please contact our after-sale service in time for purchase and replacement. We receive order of*
-

electrodes only by set. Old and new electrodes cannot be mixed up.



CAUTION

- ◆ *In order to achieve satisfactory ECG recording results, metal electrodes must be closely contacted with the skin evenly.*
 - ◆ *Metal electrodes must be clean. And the prepared sites contacting metal electrodes must be clean, free of grease and perspiration.*
 - ◆ *When placing chest electrodes, don't let metal electrodes contact with each other, or the conductive cream zone of adjacent electrodes overlap.*
 - ◆ *When placing the four limb electrodes, don't injury the patient' hands and feet. After placement, check whether the electrodes are too loose or too tight.*
 - ◆ *Frequent plugging and unplugging may cause the metal plates of limb electrodes moving or loosening, which need adjustment during use.*
 - ◆ *Reusable electrodes should be immediately cleaned after each use.*
-

5.4 Electrode Placement

5.4.1 ECG Patient Cable

Contrast Table of Electrodes and Lead Wires:

Placement	Symbol	Color Code
Right Arm	R	Red
Left Arm	L	Yellow
Left Foot	F	Green
Right Foot	RF	Black
Chest	C1	Red
Chest	C2	Yellow
Chest	C3	Green
Chest	C4	Brown
Chest	C5	Black
Chest	C6	Purple

5.4.2 Limb Lead Placement

R —Right Arm

L —Left Arm

RF —Right Foot

F —Left Foot




Fig 5.3 Chest & Limb Lead Electrode Placement

5.4.3 Chest Lead Placement

C1: Fourth intercostal space at the right sternal border.

C2: Fourth intercostal space at the left sternal border

C3: Midway between location C2 and C4.

C4: Left mid-clavicular line in the fifth intercostal space.

C5: Left anterior axillary line on the same horizontal level as C4.

C6: Left mid-axillary line on the same horizontal level as C4.

5.4.4 Pediatric Electrode Placement

When acquiring pediatric ECG, C3 lead should be placed in C4R position rather than the site where standard C3 lead is placed, as shown in right figure.




Fig 5.4 Pediatric Chest Lead Electrode Placement

5.5 Electrode Connection

Plug chest lead wires and limb lead wires respectively into connection holes of the chest electrode bulbs and limb electrode clamps. Adjust the contact sites to guarantee compact connection. Pay attention to the electrode placement.

5.6 Lead-off Handling

When the electrode is loose, or the electrode is not properly connected to the lead wire, or the ECG cable is loose from the machine side, the lead name on the ECG display will be

red. The "" icon is illuminated (white) after all leads are connected.

5.7 Entering Patient Information

Some of the patient's information directly affects ECG measurements. Patient information must be checked before starting ECG acquisition. The user can enter into Patient

Information page with **F1** key on panel or hotkey  at bottom right of numeric keyboard. Enter patient information by keyboard; select the desired blank by direction key. For details, please refer to *Chapter 2*. After filling out the information, please confirm and exit.

Chapter 6 ECG Acquisition and Recording

6.1 Acquisition Preparation

To secure the patient and to record stable accurate electrocardiogram (ECG), please carefully check the following items before you power on the machine for measurement.

- 1) Check whether the exam room is appropriate;
 - ◆ Whether there are apparatuses such as X-ray machine, ultrasound system and so forth in the patient room, as they may interfere with each other;
 - ◆ Whether the ground is well connected;
 - ◆ Whether room temperature and humidity are suitable; temperature would better be in the range of 20 ~ 25°C, ambient humidity would better be in the range of 30% ~ 60%.
- 2) Whether the power is well connected.
 - Whether the power plug is loose;
 - Whether the power cord is intertwined;
 - Whether the power is sufficient if supplied by battery.
- 3) Whether the lead wires are well connected.
 - Whether the plug is loose;
 - Whether the patient cable stays too close to the power cord;
 - Whether the lead wires and electrodes are well connected;
 - Whether electrodes are installed loosely, or adjacent electrodes have contact with each other.
- 4) How is the patient?
 - Whether the patient is too nervous; whether he or she moves or talks
 - Whether the patient's hands or feet contact metal objects such as bed edges.
- 5) Whether instruments are in good condition.
 - Whether the ECG machine is damaged;
 - Whether ECG machine is inspected and maintained regularly;
 - Whether recording paper is sufficient.

6.2 Acquisition and Recording

6.2.1 Recording Setup

After the checks are all right, power on the ECG machine, it enters into Waveform Acquisition interface, and then you can observe the waveforms after they become stable. Set the sampling mode, layout, gain, and paper speed on actual needs, and then press




key to print ECG waveforms. On Waveform Acquisition interface, you can set the

ECG machine with **F1~F7** key or by clicks on touchscreen. Press the key



to switch between auto mode and manual mode. In manual mode, you can change leads with



“” key. If you select a layout with rhythm lead, set it in **12-Lead**. For details, please refer to *Chapter 2*.



CAUTION

- ◆ *When ECG signal is just connected or the machine receives overload noise, waveforms will be chaotic, baseline drift will be severe, and waveform amplitude may exceed the maximum width. At this moment, please wait until the machine is connected, the patient calms down, and waveforms displayed on the interface become stable, then start measuring and recording.*
 - ◆ *When the ECG machine is overloaded or any part of the amplifier is saturated, or ECG machine works abnormally, at this moment, only baseline is displayed on the interface. In order to achieve accurate results, please wait until the machine is connected, the patient calms down, and waveforms displayed on the interface become stable, then start measuring and recording.*
 - ◆ *During the process of ECG signal acquisition, if waveforms become cluttered or unstable, please refer to Chapter 8.*
-

6.2.2 Recording Report

The following are the reports printed in real-time auto-simultaneous mode:




Fig. 6.1 ECG Recording Report(a)




Fig. 6.1 ECG Recording Report(b)

The ECG report above consists of two parts (a) and (b) with a print layout of 6×2.

This ECG report contains the following information: The ECG waveform 6×2 layout, patient information, date of inspection, and measurement information.

Measurement information:

HR: heart rate

P duration: average value of the P-wave duration of the average heart beat of each lead

PR interval: average value of the PR interval of the average heart beat of each lead

QRS duration: average value of the QRS duration of the average heart beat of each lead

QT/QTc interval: average/normalized value of the QT interval of the average heart beat per lead of each lead

P/QRS/T electric axis: the dominant direction of the average integrated ECG vector on the frontal plane.

RV5/SV1 amplitude: the maximum amplitude of the R and R' waves in the average heart beat on the lead V5 / The maximum value of the absolute amplitude of the S and S' waves in the average heart beat on the lead V1

RV5+SV1 amplitude: RV5 and SV1 sum

Minnesota code (optional) ECG code

Average template (optional) : multiple periodic waveform for each lead will be combined into an average of a single periodic waveform

Measurement matrix (optional) : the horizontal display shows 12 leads; the vertical display shows the parameters of each lead, such as PQRST wave start and end points, P, QRS, T wave group interval, etc.

The dotted line on the ECG waveform is the position marker, which marks the start and end points of the P wave, the start and end points of the QRS wave, and the T wave end point.

Diagnostic results: the diagnostic results show the results of the automated diagnosis.

Top information; name of medical institution

Bottom information: 0.6~35Hz (0.6Hz baseline drift filter, 35Hz low pass filter)

AC50 (50Hz AC filter)

25mm/s (paper feed speed)

10mm/mV (gain)

iMAC 120pro (machine model)

1.0.25 (software version)

Simultaneous (recording order)

Inspection date

Confirm report and sign



WARNING

- ◆ *This machine's auto-measurement precision conforms to the measurement standard of ECG machine (please refer to Appendixes). Diagnostic structure is made of related parameters in each wave segment. The interpretations are for doctors' reference only, which cannot be used as the basis of clinical treatment.*
-

6.2.3 ECG Waveform Description

Standard ECG waveform is shown as follows:




Fig 6.2 Standard ECG Waveform

Meaning and description of each segment of ECG waveform:

- ◆ P wave is of blunt round shape, and its amplitude is lower than T wave's.
- ◆ PR interval: time from the start point of P wave to the start point of QRS wave group. It represents the depolarization time from atrial to ventricle. The older the patient is or the slower the heart rate is, the longer PR interval will be. Abnormally-prolonged PR interval indicates atrioventricular conduction disturbance.
- ◆ QRS waveform group: it indicates the changes of potential and time during ventricular

muscle depolarization. R wave is high and narrow without notch, and its segments are completely above or below the base point.

- ◆ S-T segment: a horizontal line from the end point of QRS wave group (J point) to the start point of T wave is called S-T segment. Down deviation of any normal lead's S-T segment should not exceed 0.05mV. When S-T segment's down deviation exceeds the standard range, it is common in myocardial ischemia or strain. Normally, if S-T segment deviates upwards, limb lead and precordia leads of V4-6 should not exceed 0.1mV and precordia leads of V1-3 should not exceed 0.3mV. When S-T segment's up deviation exceeds the standard range, it is common in acute myocardial infarction or pericarditis.
- ◆ T wave is of blunt round shape. Its amplitude is lower than 1/3 of R wave's. But it takes longer time. The direction of T wave is usually the same as that of the main wave of QRS wave group. I, II, V4-6 leads are upright while aVR is upside down. Other leads can be upright, bidirectional or upside down. If V1 is upright, V3 could not be upside down. In the leads of QRS wave group whose main wave is up, when T wave is low and smooth or upside down, it is common in myocardial ischemia or hypokalemia.



CAUTION

- ◆ *Here is a simple description of ECG waveforms. For details, please see relevant references.*
-

6.3 Freeze

Before ECG signal acquisition, please set the freeze time on Standard ECG interface. For details, please refer to *Section 2.3.2.1*. On Freeze interface, use the real-time hotkeys to observe ECG waveforms as needed.




Fig 6.3 Freeze Interface

6.4 Report Storage

6.4.1 Report Storage Location

Configure the report storage location in Application→ Convention→ General with options such as Local/SD card/U disk.

Local: When sampling time is 10s, and layout is 6×2, you can store no less than 1500 reports.

When sampling time is 2.5s, and layout is 6×2, you can store no less than 2000 reports.

U disk (SD card): you can select U disk (SD card) only when it is inserted. Reports are stored in .ECG format and can be viewed from this ECG Machine or dedicated software workstation.



CAUTION

- ◆ *Local storage has limited capacity, thus U disk (SD card) is recommended to store the reports.*
-

6.4.2 Report Storage Mode

In the [Report Management] interface, you can set the report save mode to automatic / manual.

Auto Storage

In auto recording mode, after each measurement, the system will automatically store the current report to designated location.

When freezing the ECG waveforms, the system will automatically generate the report and store it in designated location.



CAUTION

- ◆ *In auto recording mode, the system will store the report only when the measurement is finished. If it stops recording abnormally, the report will not be saved.*
-

-
- ◆ *In manual recording mode, reports are automatically saved every 10s.*
 - ◆ *When a lead is off the patient, freezing loses its effect, still it can print, but the system cannot store the recording results automatically.*
-

Manual Storage

When the save mode is set to manual, when the recording ends, the system will pop up the "Save the ECG file?" dialog box .

After freezing the ECG waveform, click the Back button on the Freeze interface. The system will pop up the "Save the ECG file?" dialog box .

Chapter 7 Report Management

7.1 Report Storage

The report is default to be saved locally, but the local storage capacity is limited. It is recommended to use USB storage. The USB has large storage capacity, which is convenient for category management and storage.

To avoid damage or loss of external storage devices, you can back up the reports in the storage device to a dedicated computer server.

7.2 Report Management

In the waveform acquisition interface directly click the “” button at the bottom of the screen or click **F6** key enter into the **[Report Management]** interface. For the layout of the interface, see section 2.3.7.

The general setting interface has a locking function. Pressing the key combination “CTRL+ALT+L” will pop up a lock dialog box, which can lock the functions of report printing, report transmission, report deletion and so on. After locking, the option button is grayed out; if you need to change the option settings, press the key combination "CTRL+ALT+L" again to unlock the option.

7.2.1 Selecting Reports

Ways of selecting reports:

- Select a storage device to view reports;
- Touch select the desired report or control it by an external USB mouse;
- Select all reports on the current page;
- Check the Select All button to select all reports.

7.2.2 Searching Reports

Enter relevant information (such as ID number, name) in the “” query dialog box to query related reports. After the query, the list area only displays the reports that meet the requirements, which is convenient for finding and operating. Figure 7.1.




Fig 7.1 Searching Reports

After the query is completed, the query information is deleted, and the content of the report list before the query can be restored.

You can also filter the report in the following ways:

(1) Click the up arrow next to the “ID” button, all reports will be sorted in ascending order by ID number; click the down arrow, all reports will be sorted by ID number in descending order; click the “ID” button, all reports will be sorted in descending order of report modification time.

(2) Click the “Type” button to filter the report by report type. They are: all, standard 12 leads, rhythm, etc.

(3) Click the “Status” button to filter the report by report upload status. They are: all, not uploaded, uploaded, failed to upload, diagnosed, etc.

7.2.3 Editing Reports

After selecting the report, click the “” button to enter the ECG report preview interface, you can edit the patient information in the report, edit or re-analyze the diagnostic parameters and diagnosis results in the report. Click the “” button, the “Case dictionary” page will pop up, and the diagnosis result can be quickly edited; the medical record can also be customized and added to the “Case dictionary”. Click the “” button, the “Previous diagnosis” and “Current diagnosis” comparison pages will pop up. Click the “” button to give the diagnosis parameters and diagnosis results of the ECG waveform displayed on the current ECG preview interface. After editing, click the “” button to save and return to the report management interface.




Fig 7.2 Editing Reports




Fig. 7.3 Reanalysis report

7.2.4 Printing Reports

When preview a report, you can print it with  key.



CAUTION

◆ When preview and edit the report, you can select only one report.

7.2.5 Deleting Reports

Select a report and press the key  to delete it.

- (1) Check to select all the reports in **[Report Management]**
- (2) Check to select all the reports in the current page.




Fig 7.4 Delete Reports

7.3 Report Transmitting

Here are steps of transmitting reports:

1. Click the key  at bottom of **Report Management** interface.
2. Select a time option (today, three days, one week, all)
3. Select the report to be transmitted
4. Click the "" button on the report transmission interface to select the transfer mode (FTP/HTTP/SAMBA/DICOM/USB/SD card/local); (see 4.2.8 [General Setting – Transfer Settings] section and 4.2.9 for details). **[General Settings – DICOM Settings]** section. Note: The transfer mode FTP/HTTP/SAMBA/DICOM can only be used in the [Report Management] interface transmission after the settings interface is enabled.)
5. Continue to click the "" button to complete the report transmission. (Support multi-page PDF transmission. For example, if the report template is set to select all, upload the report PDF document, the first page displays the waveform, the second page

displays the measurement matrix, and the third page displays the average template. The PNG format only supports a single report template options.)

6. Status bar :
-  means not uploaded ;
 -  means upload successful ;
 -  means upload failed ;
 -  means report reload successful;

8.1



The screenshot shows the 'Report management' interface. At the top, there are filters for 'Today', 'Three days', 'One week', and 'All' (selected). A search bar is labeled 'Query ID/name'. Below the filters, there is a table with columns: ID, Name, Gender, Type, Status, and a checkbox. The table contains three rows:

ID	Name	Gender	Type	Status	Checkbox
20190801085922			ECG-12		<input checked="" type="checkbox"/>
20190801085813	AA	M	ECG-12		<input type="checkbox"/>
20190801084820			ECG-12		<input type="checkbox"/>

Below the table, there is a confirmation dialog box that says 'Transfer 1 reports?' with 'Yes' and 'No' buttons. On the left side of the interface, there are several sections:

- Diagnostic parameters:** HR: 80 bmp, P/QRS/T: 21/16/21°, PR: 166 ms, QRS: 82 ms, QT/QTc: 354/389 ms, RV5/SV1: 1.109/0.487 mv, RV5+SV1: 1.596 mv.
- Diagnostic result:** Normal ECG, Sinus rhythm, --- Interpretation made without knowing patient's gender/age ---
- Report information:** Inspection time: 2019-08-01 08:59:22, Reporting time: 20 s, Lead type: 12Lead.
- Status description:** Not uploaded (with a yellow upload icon and a red 'X' over it).

At the bottom of the interface, there are navigation buttons: Preview, Last page, Next page, Delete, Local, Transfer, and Refresh.

Fig 7.5 Report Transmission

7.4 Report Refresh

Report brush operation steps:

- 1) Successfully enable HTTP and remote diagnosis on the [**General Settings - Transfer Settings**] screen. (For details, see section 4.2.8 [**General Settings - Transfer Settings**]);
- 2) Select the report that needs to be transmitted;
- 3) Click the "" button on the report management interface to select the HTTP transmission mode;
- 4) Continue to click the "" button to complete the report transmission.
- 5) After the remote diagnosis is completed, click the "" button to return the ECG report after the diagnosis has been updated.

Chapter 8 Troubleshooting

To record a stable and accurate ECG, when a failure occurs, please find out its cause, and solve it with effective solutions.



WARNING

- ◆ *ECG machine cover should be opened only by qualified service personnel.*
There are no user-serviceable parts inside the ECG machine.
-

8.1 Interference Problem

During use, ECG machine will inevitably be disturbed by the environment, itself, human static electricity etc. The ECG machine is desired with functions of myoelectric filter, baseline drift filter, and frequency filter. As the filter band is limited, interference signals cannot be filtered out completely. Therefore, please avoid the interferences caused by the environment or non-standard operation during use.

8.1.1 AC Interference




Figure 8.1 ECG with AC interference

1) Environment Cause:

- Both ECG machine and metal bed are properly grounded.
- Avoid electrical devices of large power working in the vicinity of the ECG machine, such as X-Ray machine or ultrasound instrument etc.

2) Patient Cause:

Inform the patient of no touching the wall or metal bed edges. Don't let other people contact the patient.

3) Electrode Cause:

- Check whether the electrodes or lead wires are connected correctly,
Electrodes and skin are well applied with conductive cream, clean the patient's electrode sites with medical alcohol, apply conductive cream on the sites evenly, conductive cream on each electrode can't be cross-linked
- Check whether the patient cable is too close to or intertwined with the power cord.
- Check whether the metal part on the connection of lead wire and electrode is rusty or dirty, if it is, please clean it.
- Check whether the patient cable has poor contact, please replace a new cable and try again.

If the interference can't be cleared out by the solutions above, please make sure whether frequency filter is activated.

8.1.2 EMG Interference




Fig 8.2 EMG Interference

1) Environment Cause:

- Check whether the exam room is comfortable,
- Check whether the indoor temperature is too low,
- Check whether the bed is small and narrow.

2) Patient Cause:

- Explain to the examinee that ECG examination is very simple, which will not injury his or her body, or have sequelae;
- Make the patient relax physically and mentally, and breathe gently.
- Do not let the patient move or talk.

3) Electrode Cause:

- Check whether limb electrodes are installed too tightly, which makes the patient

feel uncomfortable,

- Check the metal part on the connection of the lead wire and electrode is rusty or dirty, if it is, please clean it.

8.1.3 Baseline Drift



Fig 8.3 Baseline Drift Waveform Graph

1) Environment Cause:

- Check whether the exam room is comfortable,
- Check whether the indoor temperature is too low,
- Check whether the bed is small and narrow.

2) Patient Cause:

- Explain to the patient that ECG examination is very simple, which will not injury his or her body, or have sequelae;
- Make the patient relax physically and mentally, and breathe gently.
- Let the patient not move or talk.

3) Electrode Cause:

- Check whether limb electrodes are installed too tightly, which makes the patient feel uncomfortable,
- Check whether the electrode is loose or poorly connected.
- Check whether the metal part on the connection of the lead wire and electrode is rusty or dirty.
- Make sure that all electrodes are of the same specification; mixed use of new and old batteries will also cause interference.

If the interference can't be cleared out by the solutions above, please make sure whether frequency filter is activated.

8.2 Recorder Failure

Failure	Possible Cause	Solutions
The paper feeds slowly and unevenly.	As the paper-feeding device has been used for a long time, its transmission ability is degraded by worn gear or loose connector.	Tighten the transmission unit, and apply some lubricating oil on the gear and both ends of paper shaft.
	As the paper-feeding device has been used for a long time, its transmission resistance increases.	Contact our service department for maintenance or replacement.
	The recorder is deformed by external force collision, thus affecting the paper speed.	Contact our service department for maintenance or replacement.
	The paper is out of specification, thus the resistance becomes over-large.	Select and use the specified paper.
	The paper has been installed for a long time, it gets heated or moistened, which makes local viscosity increase, thus affecting the paper speed.	Replace the paper
	ECG machine is not well cleaned and maintained. The recorder's transmission unit is dusty, thus degrading the transmission ability.	Inspect and clean the ECG machine to remove moisture and dust.
The paper doesn't feed, while paper is detected.	The motor is damaged.	Contact our service department for replacement.
	Main control board failure.	Return to depot maintenance.
The printer works with noises but the paper doesn't feed.	Transmission gear is stuck by some hard object.	Clear out the hard object
	Transmission gear teeth are damaged.	Contact our service department for replacement.
It's detected lack of paper.	Recording paper is not well placed or the recorder's printer door is not well closed.	Place the paper again and well close the printer door.
	Paper detector transducer is dusty.	Clear the transducer with Anhydrous ethanol.

Failure	Possible Cause	Solutions
It prints unclear or with breakpoints	Recording paper is out of specification.	Replace it with our paper or with better paper of the same specification.
	Paper shaft is dusty.	Clean the paper shaft.
	Print head is dusty.	Clean the print head.
After pressing "Stop", the recorder still works, but prints nothing.	Recording paper is installed backwards. Black label direction is wrong	Reinstall the recording paper.
	Recording paper is out of specification.	Select the recording paper with black label.
	Black label detection sensor head is dusty.	Clean the transducer head with a cotton swab dipped in medical alcohol.
It prints empty	Recording paper is installed backwards	Properly install the recording paper, with grid side right facing to the print head.

The solutions above can solve common printing failures. If there are still some issues unsolved, please contact our service department, or return the ECG machine to us for maintenance or replacement.

Chapter 9 Maintenance

9.1 Cleaning and Disinfection 23.

Please keep the ECG machine and its accessories clean. And in order to avoid damaging the ECG machine, please follow the regulations below:

- Please dilute the cleaner and disinfectant according to the manufacturer's instructions, or use the cleaner and disinfectant whose concentration is as low as possible;
- Do NOT immerse the device into liquid;
- Do NOT dump any liquid onto the device or its accessories;
- Do NOT let any liquid enter into the device ;
- Do NOT use abrasive materials such as steel wool or silver polisher, or any strong solvents such as acetone or acetone detergent.



- ◆ *You must turn off the power, and disconnect the power cord and the outlet before cleaning and disinfecting the machine;*



WARNING

- ◆ *The ECG machine can be cleaned or disinfected only by the materials and methods listed in this chapter. We will not provide warranty for any damage or accident caused by using other materials or methods;*
 - ◆ *We will not assume any responsibility for the effectiveness of using the listed chemicals or methods as infection control ways. For ways of infection control, please consult the infection prevention department in hospitals or epidemiologists.*
-



CAUTION

- ◆ *If you accidentally dump liquid onto the equipment or its accessories, cause damage, please contact our service department.*
-

9.1.1 Cleaning 23.

Available detergents for cleaning the host are listed as follows:

NaClO (Bleaching)

Oxydol (3%)

Ethanol (75%)

Isopropyl alcohol (70%)

It is recommended to clean the accessories with 75% ethanol.

Cleaning the host:

The ECG machine should be cleaned regularly. In those areas where the environment is seriously polluted or the sand blows heavily, it should be cleaned more frequently. Please consult or know about the hospital regulations for cleaning the ECG machine before you clean it.

While cleaning the machine:



- ◆ *Turn off the power. Disconnect the power cord, accessories and other devices connected to this ECG machine before cleaning;*
 - ◆ *Use a soft cotton ball to wipe the LCD screen with some detergents;*
 - ◆ *Use a soft cloth to clean the surface of the machine with some detergents. Avoid the ports at the sides and rear of the machine;*
 - ◆ *Wipe off the remaining detergents with a dry cloth when necessary;*
 - ◆ *Put the machine in a place with cool ventilation to dry it naturally.*
-

Cleaning ECG cables and lead wires:



- ◆ *Please remove the cables from ECG machine before cleaning them and the lead wires.*
 - ◆ *Use a soft cloth with some 75% ethanol to wipe the surface of the cables and lead wires. Avoid the metal connection parts;*
 - ◆ *Wipe off the remaining detergent with a dry cloth if necessary;*
 - ◆ *Put the cables and lead wires in a place with cool ventilation to dry them naturally.*
-

Cleaning reusable electrodes:



- *Reusable electrodes must be cleaned after each use.*
 - ◆ *Use a soft cloth with some 75% ethanol to wipe the surface of the electrodes;*
 - ◆ *Wipe off the remaining detergent with a dry cloth if necessary;*
 - ◆ *Put the electrodes in a place with cool ventilation to dry them naturally.*
-

Cleaning the recorder head:

Stains and dirt on the surface of thermosensitive recorder head will influence the record's definition. Therefore the recorder head should be cleaned regularly (at least once a month). If you find that characters on the report are light or the recorder doesn't work, it indicates that the recorder head needs cleaning.

Please follow the steps below to clean the recorder head:



- *Turn off the ECG machine;*
 - *Push the button to open the print door, and take out the paper;*
 - *Clear out the stains and dirt on the surface of thermosensitive recorder head with a cotton swab dipped with 75% alcohol;*
-

-
- *Dry the recorder head gently with a clean cotton swab;*
 - *Dry the recorder head naturally, reinstall the recording paper and close the printer door.*
-



CAUTION

- ◆ *Please don't clean the recorder head immediately after recording as the head might be extremely hot at this time.*
-

9.1.2 Disinfection 23.

Disinfection may cause some damage to the ECG machine or its accessories. It's recommended to perform disinfection only when it is necessary for the service plan in your hospital. Perform cleaning before disinfection.

Disinfectants recommended for the host contain 75% ethanol, 70% isopropyl alcohol, Shu US sterilization agent (C/D-level) with R active oxygen. It is recommended to disinfect the accessories with 75% ethanol.

9.1.3 Sterilization

It is not recommended to sterilize the ECG machine and its accessories unless otherwise required in the manual for accessories.

9.2 Routine Inspection and Test

9.2.1 Daily Inspection

Before the first use each day, the machine appearance should be inspected. Once the ECG machine is found damaged, please stop using it immediately, and contact the engineers in your hospital or our maintainers.

Inspection items include:

- No stain is on ECG machine shell; the panel and LCD screen is not broken or damaged;

- All buttons are in good condition;
- Ports, plugs and cables are not damaged or twined;
- The power cord and ECG cable are firmly and respectively connected with the machine;
- The recording paper is installed correctly, and sufficient for use;
- The battery is installed and fully charged;
- Chest bulbs are free of cracks, and limb clamps clamp well with adequate force.

9.2.2 Regular Inspection

When used continuously for 6 to 12 months, or after maintenance or upgrading, the ECG machine should be tested completely by the qualified service personnel, ensuring that the ECG machine works normally.

Inspection items are listed as follows:

- The environment and power meet the requirements;
- The ECG machine and its accessories are not mechanically damaged;
- The power cord, ECG cable and lead wires are not worn;
- The battery performance is in good condition;
- Function test: used for inspecting the inside of the ECG machine. This test needs to be performed by our professionals or by the authorized personnel under the guidance of our technicians.



CAUTION

- ◆ *For accidents or equipment damage caused by lack of necessary maintenance, we will not assume any responsibility.*
-
-

9.3 Battery Usage & Maintenance

9.3.1 Overview

ECG machine is configured with rechargeable lithium-ion battery so as to ensure that it works normally while moving in the hospital or when the power fails. When the power fails suddenly, the system will automatically enable the battery to supply the ECG machine, thus the machine won't stop working. .

The rechargeable lithium-ion battery used in the machine has over-charge protection circuit, so it will not get over charged. Its output voltage is related to its power. When its power is low, its output voltage will decrease, but it will not affect ECG machine's normal running. There is a battery detection system inside the machine. If the battery is short of power, it will not be used in order to prevent its over-discharging. The battery plug corresponds to its socket, which could prevent misconnection of battery polarity. The battery is well sealed itself, which will not leak electrolyte or dangerous gas during use.



WARNING

- *Be sure to use and maintain the battery according to the contents in this chapter.*
 - *If the battery has sign of damage or leakage, please replace it immediately. Do not install faulty battery into ECG machine.*
-



CAUTION

- ◆ *In order to prevent the machine from working interruption caused by sudden power failure, we suggest the user always install a full-charged battery in it.*
 - ◆ *When ECG machine is supplied by battery, if the battery is short of power, the machine will crash with black screen. This is a normal phenomenon, which could be eliminated by connection with AC power or charging the battery.*
-

Battery icons on the screen indicate its status:

Those two icons  indicate that the battery works normally. The white grid shows the power capacity.

This icon  indicates that the battery power is extremely low. When battery power is extremely low, ECG machine will pop out a message “Extremely Low Battery Power”, and the recorder cannot work. At this moment, please immediately connect the machine with AC power to charge the battery. Otherwise it will power off automatically.

9.3.2 Battery Charging

When the ECG machine is connected with AC power, no matter it is turned on or off, the battery will be charged. When the battery is being charged, its light will be lit. Once fully charged, the light will go out

When charging the battery with the ECG machine turned off, in an environment with temperature range of $25^{\circ}\text{C}\pm 5^{\circ}\text{C}$, the battery is charged to 90% in no more than 3 hours, and charged to 100% in no more than 3.5 hours.

9.3.3 Battery Replacement

The battery installed in this ECG machine should be replaced by authorized service engineers. Please contact our service engineers if it is demanded for battery replacement.

9.3.4 Battery Guidance

The battery's life depends on its usage frequency and time. If the lithium-ion battery is properly maintained and stored, its life will last for about 2 years. If used improperly, its life will be shortened. We recommend replacing the battery every two years.

In order to guarantee the battery's life, please pay attention to the following guidance:

- Battery performance must be inspected once half a year. Besides, you also need to inspect the battery performance before maintenance of the ECG machine or when the battery is suspected to be faulty.

- When the battery has been used or stored for three months or when its working time obviously shortens, perform an optimization on it.
- Before the ECG machine is delivered or when it will not be used for more than 3 months, please take out the battery.
- If the ECG machine has not been used for a long time with the battery installed in it, the battery's life will shorten. The battery should be charged and discharged at least once every three months.
- When the lithium battery is laid aside with 50% of its full power, it can be stored for about 6 months. After 6 months, the battery must be charged again to full power, and then use it to supply ECG machine. When its power reduces to 50% of the full power, take it out of ECG machine and lay it aside again.
- When storing the battery, please make sure that its electrodes do not touch metal objects. If the battery needs to be stored for a long time, put it in a cool environment, which can delay battery aging. Ideally, the battery should be stored in a cool environment whose temperature is 15°C. If the battery is placed in high heat for a long time, its life will obviously shorten. Do not store the battery in the environment whose temperature is not within the range of -20°C~60°C.



CAUTION

- ◆ *Place the battery in a place out of reach of children.*
 - ◆ *Only use the battery designated by the manufacturer.*
-

9.3.5 Battery Maintenance

Battery Performance Optimization

The battery should be optimized for its initial use. A complete optimization period is continuously charging the battery to full capacity. Then discharge it until the ECG machine is power off. During use, the battery should be optimized annually to sustain its life.

Please optimize the battery by following steps:

1. Disconnect ECG machine with the patient;

2. Connect ECG machine with AC power, continuously charge the battery to its full capacity, and then the indicator light goes out.
3. Disconnect AC power, and supply ECG machine with battery power until it is power off.
4. Connect ECG machine with AC power again, and continuously charge the battery to full capacity, then the light goes out.



CAUTION

- ◆ *As the time of using the battery increases, its actual power capacity will decrease. For the used battery, full-capacity icon indicates that neither its power capacity nor supply time could still meet the manufacturer's specification. When optimizing the battery, if you find that its supply time shortens obviously, please replace it.*
-

Battery Performance Inspection

Battery performance will degrade as times of using the battery increase, thus it should be inspected once a year. Besides, it also needs inspection before servicing the ECG machine or when the battery is suspected to be faulty.

Please inspect the battery according to the following steps:

1. Disconnect ECG machine with the patient
2. Connect ECG machine with AC power, and constantly charge the battery to full power, then the light goes out.
3. Disconnect AC power, supply ECG machine by battery until it is power off.

The Battery's supply time reflects its performance. After announced charging time, if its actual supply time is obviously lower than the time declared in specification, please contact the maintainer to replace the battery.



CAUTION

- ◆ *Battery supply time depends on configuration and operation of the machine..*
 - ◆ *If the battery's supply time is too short after fully charged, obviously less than*
-

stated in the specification, it might be damaged or faulty. Please contact the maintenance staff to replace the battery.

9.3.6 Battery Recycling

If the battery is obviously damaged or cannot be charged, it should be replaced and recycled properly. When disposing of the used battery, please follow relevant laws and regulations. .



WARNING

- ◆ *Do not disassemble the battery or throw it into fire or short it out. Its burning, explosion or leakage may cause personal injury.*
-

9.4 Usage and Maintenance of Recording Paper

Please follow the rules below when storing the recording paper:

- Store it in cool dry environment free from high temperature, humidity and direct sunlight.
- Do not put it in fluorescent light for a long time.
- Do not let it contact polyvinyl chloride (PVC), which will cause its color change.
- Do not overly the used paper while storing, which may cause its printout transferring with each other.
- Using the paper provided by the manufacturer or of specification dedicated by the manufacturer. Otherwise it may shorten thermosensitive recorder head's life, recorded waveforms will become fuzzy and the paper will feed poorly.

9.5 Maintenance of Electrodes and Lead Wires

Conduction of each lead wire will directly affect ECG traces. If it conducts poorly (any one lead conducts poorly), it will cause virtual image of corresponding lead wire on ECG traces. Therefore the conduction must be inspected regularly, at least once a month.

Slightly bending or entangling the lead wire will shorten its life. Please put it in as good order as possible before use.

Electrodes must be properly stored. After long-term use, their surfaces may become oxidized and discolored because of corrosion, at this moment, it's better to replace them.

Chapter 10 After-sale Service information

1. When users begin to use the ECG machine, they should fill the details in warranty card and send it back to the manufacturer by mail or email in time, the manufacturer will build the users' profiles and regularly contact them to know about the usage, which will help provide targeted first-rate services constantly.
2. During normal use per the manual and operation notes, once the machine breaks down, please contact the manufacturer's after-sale service center immediately. Users can enjoy free service within the stipulated time on warranty card since the purchase day.
3. The manufacturer after-sale service team or local support partners may fulfill its warranty promise by ways of visiting your place, telephone guidance or delivery back to the manufacturer. .
4. Even within warranty period, the following services will be charged:
 - ①Fault and damage caused by improper operation of users;
 - ②Fault or damage caused by falling down while moving the machine after purchase;
 - ③Fault and damage caused by repairing, transforming or decomposing the machine without the manufacturer's authorization.
 - ④Fault and damage caused by fire or natural disaster after purchase;
 - ⑤Fault and damage caused by using thermal paper unspecified by the manufacturer;
 - ⑥Fault and damage caused by connection with other devices;
 - ⑦Warranty seal is broken. Users privately alter and replace the serial numbers of the machine and lead wires. .
5. If the product fails within three months and it is not caused by article 4, the company will replace the main unit free of charge, but the accessories, worn parts and consumables will not be replaced.
6. The company shall not be responsible for the failure of other connected devices directly or indirectly caused by the failure of the product. This warranty system is only valid in China.

7. If warranty label is damaged, the manufacturer has rights to exempt free service within stipulated time on warranty card.
8. For chargeable services out of warranty period, it's recommended to continue "Service Contract Rules". For details, please consult the customer service center of the manufacturer.

Chapter 11 Accessories



WARNING

- ◆ *Use the accessories stipulated in this manual only. Other accessories may damage this ECG machine or cannot meet the specification declared in this manual.*
 - ◆ *Disposable accessories can be used only once; Reuse will cause performance degradation or cross infection.*
 - ◆ *If the accessories or their packages are found damaged, please do not use the accessories.*
-

Accessories:

Name	Type
Patient Cable	Pinning fibrillation-proof ECG cable
Chest Bulbs	ECG chest electrodes(φ4)
Limb Electrodes	ECG limb electrodes (φ4)
Power Cord for Adapter	Power cord of European standard

I.4 Appearance Parameters

Size	(L × W × H) 285mm×360mm×94mm
------	------------------------------

Weight	about 3.7 kg
--------	--------------

I.5 Environmental Conditions

Operation

Ambient Temperature	+ 5°C ~ + 40°C
---------------------	----------------

Ambient Humidity	20%~85%(no condensation)
------------------	--------------------------

Atmospheric Pressure	570hPa~1060hPa
----------------------	----------------

Shipment and Storage

Ambient Temperature	-20°C ~ +55°C
---------------------	---------------

Ambient Humidity	10%~95%
------------------	---------

Atmospheric Pressure	500hPa~1060hPa
----------------------	----------------

I.6 Adherence to Standards

EN ISO 13485:2016	Medical devices-quality management system-Requirements for regulatory purpose;
-------------------	--

EN ISO 14971:2012	Medical devices-Application of risk management to medical devices
-------------------	---

EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
---------------------	--

ISO10993-1:2009/AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
-------------------------	--

ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
-----------------------------	--

ISO 10993-10:2010	Biological evaluation of medical devices-Part 10:Tests for irritation and skin sensitization
------------------------------	--

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

EN 1041:2008	Information supplied by the manufacturer of medical devices
IEC 60601-1:2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-25: 2011	Medical Electrical Equipment - Part 2-25: Particular Requirements For The Basic Safety And Essential Performance of Electrocardiographs;
EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance -Collateral standard: Electromagnetic compatibility-Requirements and tests;
EN 62304:2006/AC:2008	Medical device software - Software life-cycle processes
EN 62366:2008	Medical devices - Application of usability engineering to medical devices
EN 60601-1-6:2015	Medical electrical equipment - part 1-6: general requirements for basic safety and essential performance - collateral standard: usability

Appendix II Electromagnetic (EMC)

Electromagnetic compatibility (EMC) is defined as the ability of a product, device, or system to function properly in its electromagnetic environment without posing unacceptable electromagnetic disturbances to anything in the environment.

Anti-electromagnetic interference is the ability of a product, device, or system to function properly in the presence of electromagnetic interference (EMI).

It is designed and manufactured in accordance with existing electromagnetic compatibility standards and related requirements. Use in the presence of an electromagnetic field may cause performance degradation such as output anomalies. If this happens frequently, it is recommended to check the environment in which the ECG is used to determine possible sources of disturbance. These harassments may come from other electrical equipment used in the same room or in a nearby room, or from portable and mobile RF communications equipment such as cell phones, walkie-talkies, or from nearby radios, televisions, or microwave transmission equipment. If electromagnetic interference (EMI) interferes with the ECG, it may be necessary to move the ECG to another location or take appropriate electromagnetic interference suppression measures.

This product complies with the requirements of the EMC standard IEC 60601-1-2.

 **Warning**

- ◆ *It will not be used for the lead wire and power cord of the electrocardiograph for the electrocardiograph, which may result in an increase in the emission of the electrocardiograph or a decrease in the immunity.*
- ◆ *The electrocardiograph should not be used close to or stacked with other equipment. If it must be used close to or stacked, it should be observed that it will function properly in the configuration in which it is used.*

 **Note**

- ◆ Medical equipment has special precautions for EMC and needs to be installed and used according to the EMC information provided in the
-

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

documentation provided with the ECG.

- ◆ This section includes information on electromagnetic radiation and immunity to electromagnetic systems. Ensure that the operation of the electrocardiograph meets the conditions specified in the reference information. Operating an electrocardiograph in an environment that does not meet these conditions may degrade the performance of the system.
- ◆ To ensure electromagnetic compatibility when installing and using an electrocardiograph, follow the information and warnings contained in this and other sections.

Description

- ◆ If you operate and use an electrocardiograph in the electromagnetic environment described in "Anti-Electromagnetic Interference" below, it will work safely and provide the following basic properties:
 - 1 button works normally;
 - 2 The host continuously collects signals and displays the waveform and measured value results on the display.

Electromagnetic radiation		
The electrocardiograph is expected to be used in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in this electromagnetic environment.		
Launch Test	Compliance	Electromagnetic Environment - Guide
Radio frequency emission CISPR 1	Group 1	The ECG uses RF energy only for its internal functions. Therefore, its RF emissions are low and there is little chance of interference with nearby electronic equipment.

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

Radio frequency emission CISPR 11	Class A	The electrocardiograph is suitable for use in all facilities that are not directly connected to the home and to the public low voltage power supply network of the home.
Harmonic emission IEC61000-3-2	Not applicable	
Voltage fluctuation/ flicker emission IEC61000-3-3	Not applicable	

Approved accessories that meet electromagnetic standards

Accessories for electrocardiographs may affect their amount of radiation. The accessories listed in this section have been tested in accordance with international standards when used in electrocardiographs to confirm compliance with radiation standards. Please use only the attachments listed in this section.

When connecting the accessories to the ECG machine, the user should ensure the electromagnetic compatibility of the ECG machine. Unless otherwise stated, use only EMC-compliant equipment.

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

No.	Name	Cable length (m)	Shielding or Not
1	power cord	1.6	NO
2	Patient cable	< 3	NO

Anti-electromagnetic interference

The electrocardiograph is designed for use in the electromagnetic environment specified herein. The user or user of the electrocardiograph should ensure that it is used in such an environment.

Electromagnetic Immunity Guidelines and Statements			
The electrocardiograph is expected to be used in a defined electromagnetic environment and the purchaser or user should ensure that it is used in such an electromagnetic environment.			
Anti-interference test	IEC60601Test level	Compliance level	Electromagnetic environment- guide
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV Contact discharge ±8 kV Air discharge	±6 kV Contact discharge ±8 kV Air discharge	The floor must be wood, concrete or tile. If the floor is covered with synthetic material, the relative humidity is at least 30%.
Electrical fast transient (EFT)/burst IEC 61000-4-4	±2kV power cord ±1 kV I/O cables (Length >3m).	±2kV power cord ±1 kV input/output cables (>3m)	The network power supply should have the quality used in a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	The network power supply should have the quality used in a typical commercial or

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

			hospital environment.
Voltage dip, short interruption and voltage change on the power input line IEC 61000-4-11	< 5%UT (drop > 95% UT) 0.5cycle 40%UT (drop 60%UT 5cycle 70%UT (drop 30%UT 25cycle < 5%UT (drop > 95% UT) 5sec	< 5%UT (drop > 95% UT) 0.5cycle 40%UT (drop 60%UT 5cycle 70%UT (drop 30%UT 25cycle < 5%UT (drop > 95% UT) 5sec	The network power supply should have the quality used in a typical commercial or hospital environment. If the user of the ECG machine needs to run continuously during a power outage, it is recommended that the ECG machine be powered by an uninterruptible power supply or battery.
Power frequency magnetic field (50/60Hz) EC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic field should have a power frequency magnetic field level characteristic typical of a typical commercial or hospital environment.
Note 1: UT is the A.C. mains voltage prior to application of the test level.			

Electromagnetic interference

Electromagnetic interference can occur on the electrocardiograph in a variety of ways. These disturbances depend on the immunity of the device. In the presence of interference, be careful when continuing to use the ECG machine.

Electromagnetic Immunity Guidelines and Statements

The electrocardiograph is expected to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in such an electromagnetic environment.

Immunity test	IEC60601 Test level	Compliance level	Electromagnetic environment - guide
Radio frequency conduction 61000-4-6	3 Vrms 150k ~ 80MHz	3Vrms	<p>Portable and mobile RF communications equipment should not be used closer to any part of the electrocardiograph, including cables, than the recommended isolation distance. This distance is calculated by a formula corresponding to the transmitter frequency.</p> <p>The recommended isolation distance is calculated as:</p> $d = 1.2 \sqrt{P}$

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

<p>Radio frequency radiation 61000-4-3</p>	<p>3V/m 80MHz 2.5GHz</p>	<p>~ 3 V/m</p>	<p>The recommended isolation distance is calculated as:</p> <p>80 MHz to 800 MHz: $d = 1.2 \sqrt{P}$</p> <p>800MHz to 2.5GHz: $d = 2.3 \sqrt{P}$</p> <p>P is the transmitter's rated maximum output power, in watts, provided by the transmitter manufacturer;</p> <p>d is the recommended isolation distance in meters.</p> <p>The field strength of the RF transmitter obtained by electromagnetic field measurement must be less than the compliance level in each frequency range.</p> <p>Interference may occur near devices marked</p>  <p>with the following symbols:</p>
<p>Note 1: At 80MHz and 800MHz frequencies, the formula for the higher frequency band is used.</p> <p>Note 2: The above guidelines may not be applicable in all situations. Electromagnetic wave propagation is affected by absorption and reflection from buildings, objects and the human body.</p> <p>Note 3: This equipment is intended to receive RF electromagnetic energy and is exempt from basic performance requirements in the occupied frequency band (2395.825MHz-2487.645MHz), but remains safe.</p>			

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

For fixed-type airports, such as wireless (cellular/cordless) telephones and terrestrial mobile radio base stations, amateur radio, AM (amplitude modulation) and FM (frequency modulation) radio broadcasts, and television broadcasts, the field strength is theoretically can not be accurately predicted. In order to assess the electromagnetic environment of a stationary RF transmitter, an electromagnetic field survey should be considered. If the measured field strength of the location where the ECG is located is higher than the RF compliance level of the above application, the ECG should be observed to verify that it is functioning properly. Additional measures may be necessary if abnormal performance is observed, such as reorienting the ECG machine and repositioning it.

Recommended isolation distance between portable and mobile RF communications equipment and ECG machines

The electrocardiograph is expected to be used in an electromagnetic environment where radio frequency radiation disturbance is controlled. Depending on the maximum rated output power of the communication device, the purchaser or user can prevent electromagnetic interference by maintaining the minimum distance between the portable and mobile RF communication device (device machine) and the ECG machine as recommended below.

Maximum rated output power of the transmitter (W)	Isolated distance for different frequency transmitters (m)		
	150 kHz ~ 80 MHz $d = 1.2 \sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2 \sqrt{P}$	800 MHz ~ 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30

Technical data file (confidential)
ECG equipment iMAC120, Zoncare

100	12.00	12.00	23.00
<p>For the rated maximum output power of the transmitter not listed above, the recommended isolation distance d, in meters (m), can be determined using the formula in the corresponding transmitter frequency column, where P is provided by the transmitter manufacturer. Transmitter maximum output rated power in watts (W).</p> <p>Note 1: Formulas for higher frequency ranges are used at frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by the absorption and reflection of buildings, objects and human bodies.</p>			

Appendix III Environmental Statement

Names and Contents of Hazardous and Noxious Substances:

Name	Hazardous and Noxious Substances or Elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (CrVI)	PBB	PBDE
Built-in PCB	○	○	○	○	○	○
plug-in board	○	○	○	○	○	○
Metal parts	○	○	○	○	○	○
Shell	○	○	○	○	○	○
Display part	○	○	○	○	○	○
Package	○	○	○	○	○	○
Accessories	○	○	○	○	○	○
○: It indicates that contents of hazardous substances in the all homogeneous materials of the part are below the limits specified in SJ/T 11363-2006 standard.						
 CAUTION		The product and its spare parts should be disposed of in accordance with the local laws and regulations, and shall not be discarded as useless together with household waste.				

I.4 Appearance Parameters

Size	(L × W × H) 285mm×360mm×94mm
Weight	about 3.7 kg

I.5 Environmental Conditions

Operation

Ambient Temperature	+ 5°C ~ + 40°C
Ambient Humidity	20%~85%(no condensation)
Atmospheric Pressure	570hPa~1060hPa

Shipment and Storage

Ambient Temperature	-20°C ~ +55°C
Ambient Humidity	10%~95%
Atmospheric Pressure	500hPa~1060hPa

I.6 Adherence to Standards

EN ISO 13485:2016	Medical devices-quality management system-Requirements for regulatory purpose;
EN ISO 14971:2012	Medical devices-Application of risk management to medical devices
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO10993-1:2009/AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
ISO 10993-10:2010	Biological evaluation of medical devices-Part 10:Tests for irritation and skin sensitization

EN 1041:2008	Information supplied by the manufacturer of medical devices
IEC 60601-1:2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-25: 2011	Medical Electrical Equipment - Part 2-25: Particular Requirements For The Basic Safety And Essential Performance of Electrocardiographs;
EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance -Collateral standard: Electromagnetic compatibility-Requirements and tests;
EN 62304:2006/AC:2008	Medical device software - Software life-cycle processes
EN 62366:2008	Medical devices - Application of usability engineering to medical devices
EN 60601-1-6:2015	Medical electrical equipment - part 1-6: general requirements for basic safety and essential performance - collateral standard: usability

Appendix II Electromagnetic (EMC)

Electromagnetic compatibility (EMC) is defined as the ability of a product, device, or system to function properly in its electromagnetic environment without posing unacceptable electromagnetic disturbances to anything in the environment.

Anti-electromagnetic interference is the ability of a product, device, or system to function properly in the presence of electromagnetic interference (EMI).

It is designed and manufactured in accordance with existing electromagnetic compatibility standards and related requirements. Use in the presence of an electromagnetic field may cause performance degradation such as output anomalies. If this happens frequently, it is recommended to check the environment in which the ECG is used to determine possible sources of disturbance. These harassments may come from other electrical equipment used in the same room or in a nearby room, or from portable and mobile RF communications equipment such as cell phones, walkie-talkies, or from nearby radios, televisions, or microwave transmission equipment. If electromagnetic interference (EMI) interferes with the ECG, it may be necessary to move the ECG to another location or take appropriate electromagnetic interference suppression measures.

This product complies with the requirements of the EMC standard IEC 60601-1-2.

Warning

- ◆ *It will not be used for the lead wire and power cord of the electrocardiograph for the electrocardiograph, which may result in an increase in the emission of the electrocardiograph or a decrease in the immunity.*
- ◆ *The electrocardiograph should not be used close to or stacked with other equipment. If it must be used close to or stacked, it should be observed that it will function properly in the configuration in which it is used.*

Note

- ◆ Medical equipment has special precautions for EMC and needs to be installed and used according to the EMC information provided in the

documentation provided with the ECG.

- ◆ This section includes information on electromagnetic radiation and immunity to electromagnetic systems. Ensure that the operation of the electrocardiograph meets the conditions specified in the reference information. Operating an electrocardiograph in an environment that does not meet these conditions may degrade the performance of the system.
 - ◆ To ensure electromagnetic compatibility when installing and using an electrocardiograph, follow the information and warnings contained in this and other sections.
-

Description

- ◆ If you operate and use an electrocardiograph in the electromagnetic environment described in "Anti-Electromagnetic Interference" below, it will work safely and provide the following basic properties:
 - 1 button works normally;
 - 2 The host continuously collects signals and displays the waveform and measured value results on the display.
-

Electromagnetic radiation		
The electrocardiograph is expected to be used in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in this electromagnetic environment.		
Launch Test	Compliance	Electromagnetic Environment - Guide
Radio frequency emission CISPR 1	Group 1	The ECG uses RF energy only for its internal functions. Therefore, its RF emissions are low and there is little chance of interference with nearby electronic equipment.

Radio frequency emission CISPR 11	Class A	The electrocardiograph is suitable for use in all facilities that are not directly connected to the home and to the public low voltage power supply network of the home.
Harmonic emission IEC61000-3-2	Not applicable	
Voltage fluctuation/ flicker emission IEC61000-3-3	Not applicable	

Approved accessories that meet electromagnetic standards

Accessories for electrocardiographs may affect their amount of radiation. The accessories listed in this section have been tested in accordance with international standards when used in electrocardiographs to confirm compliance with radiation standards. Please use only the attachments listed in this section.

When connecting the accessories to the ECG machine, the user should ensure the electromagnetic compatibility of the ECG machine. Unless otherwise stated, use only EMC-compliant equipment.

No.	Name	Cable length (m)	Shielding or Not
1	power cord	1.6	NO
2	Patient cable	< 3	NO

Anti-electromagnetic interference

The electrocardiograph is designed for use in the electromagnetic environment specified herein. The user or user of the electrocardiograph should ensure that it is used in such an environment.

Electromagnetic Immunity Guidelines and Statements			
The electrocardiograph is expected to be used in a defined electromagnetic environment and the purchaser or user should ensure that it is used in such an electromagnetic environment.			
Anti-interference test	IEC60601 Test level	Compliance level	Electromagnetic environment- guide
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV Contact discharge ±8 kV Air discharge	±6 kV Contact discharge ±8 kV Air discharge	The floor must be wood, concrete or tile. If the floor is covered with synthetic material, the relative humidity is at least 30%.
Electrical fast transient (EFT)/burst IEC 61000-4-4	±2kV power cord ±1 kV I/O cables (Length >3m).	±2kV power cord ±1 kV input/output cables (>3m)	The network power supply should have the quality used in a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	The network power supply should have the quality used in a typical commercial or

			hospital environment.
Voltage dip, short interruption and voltage change on the power input line IEC 61000-4-11	<p>< 5% UT (drop > 95% UT) 0.5cycle</p> <p>40%UT (drop 60%UT) 5cycle</p> <p>70%UT (drop 30%UT) 25cycle</p> <p>< 5% UT (drop > 95% UT) 5sec</p>	<p>< 5% UT (drop > 95% UT) 0.5cycle</p> <p>40%UT (drop 60%UT) 5cycle</p> <p>70%UT (drop 30%UT) 25cycle</p> <p>< 5% UT (drop > 95% UT) 5sec</p>	<p>The network power supply should have the quality used in a typical commercial or hospital environment.</p> <p>If the user of the ECG machine needs to run continuously during a power outage, it is recommended that the ECG machine be powered by an uninterruptible power supply or battery.</p>
Power frequency magnetic field (50/60Hz) EC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic field should have a power frequency magnetic field level characteristic typical of a typical commercial or hospital environment.
Note 1: UT is the A.C. mains voltage prior to application of the test level.			

Electromagnetic interference

Electromagnetic interference can occur on the electrocardiograph in a variety of ways. These disturbances depend on the immunity of the device. In the presence of interference, be careful when continuing to use the ECG machine.

Electromagnetic Immunity Guidelines and Statements

The electrocardiograph is expected to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in such an electromagnetic environment.

Immunity test	IEC60601 Test level	Compliance level	Electromagnetic environment - guide
Radio frequency conduction 61000-4-6	3 Vrms 150k ~ 80MHz	3Vrms	<p>Portable and mobile RF communications equipment should not be used closer to any part of the electrocardiograph, including cables, than the recommended isolation distance. This distance is calculated by a formula corresponding to the transmitter frequency.</p> <p>The recommended isolation distance is calculated as:</p> $d = 1.2 \sqrt{P}$

Radio frequency radiation 61000-4-3	3V/m 80MHz 2.5GHz	~ 3 V/m	<p>The recommended isolation distance is calculated as:</p> <p>80 MHz to 800 MHz: $d = 1.2\sqrt{P}$</p> <p>800MHz to 2.5GHz: $d = 2.3\sqrt{P}$</p> <p>P is the transmitter's rated maximum output rated power, in watts, provided by the transmitter manufacturer;</p> <p>d is the recommended isolation distance in meters.</p> <p>The field strength of the RF transmitter obtained by electromagnetic field measurement must be less than the compliance level in each frequency range.</p> <p>Interference may occur near devices marked</p> <div style="text-align: center;">  </div> <p>with the following symbols:</p>
--	-------------------------	---------	--

Note 1: At 80MHz and 800MHz frequencies, the formula for the higher frequency band is used.

Note 2: The above guidelines may not be applicable in all situations. Electromagnetic wave propagation is affected by absorption and reflection from buildings, objects and the human body.

Note 3: This equipment is intended to receive RF electromagnetic energy and is exempt from basic performance requirements in the occupied frequency band (2395.825MHz-2487.645MHz), but remains safe.

For fixed-type airports, such as wireless (cellular/cordless) telephones and terrestrial mobile radio base stations, amateur radio, AM (amplitude modulation) and FM (frequency modulation) radio broadcasts, and television broadcasts, the field strength is theoretically can not be accurately predicted. In order to assess the electromagnetic environment of a stationary RF transmitter, an electromagnetic field survey should be considered. If the measured field strength of the location where the ECG is located is higher than the RF compliance level of the above application, the ECG should be observed to verify that it is functioning properly. Additional measures may be necessary if abnormal performance is observed, such as reorienting the ECG machine and repositioning it.

Recommended isolation distance between portable and mobile RF communications equipment and ECG machines

The electrocardiograph is expected to be used in an electromagnetic environment where radio frequency radiation disturbance is controlled. Depending on the maximum rated output power of the communication device, the purchaser or user can prevent electromagnetic interference by maintaining the minimum distance between the portable and mobile RF communication device (device machine) and the ECG machine as recommended below.

Maximum rated output power of the transmitter (W)	Isolated distance for different frequency transmitters (m)		
	150 kHz ~ 80 MHz $d = 1.2 \sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2 \sqrt{P}$	800 MHz ~ 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30

100	12.00	12.00	23.00
<p>For the rated maximum output power of the transmitter not listed above, the recommended isolation distance d, in meters (m), can be determined using the formula in the corresponding transmitter frequency column, where P is provided by the transmitter manufacturer. Transmitter maximum output rated power in watts (W).</p> <p>Note 1: Formulas for higher frequency ranges are used at frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by the absorption and reflection of buildings, objects and human bodies.</p>			

Appendix III Environmental Statement

Names and Contents of Hazardous and Noxious Substances:

Name	Hazardous and Noxious Substances or Elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (CrVI)	PBB	PBDE
Built-in PCB	○	○	○	○	○	○
plug-in board	○	○	○	○	○	○
Metal parts	○	○	○	○	○	○
Shell	○	○	○	○	○	○
Display part	○	○	○	○	○	○
Package	○	○	○	○	○	○
Accessories	○	○	○	○	○	○
○: It indicates that contents of hazardous substances in the all homogeneous materials of the part are below the limits specified in SJ/T 11363-2006 standard.						
 CAUTION		The product and its spare parts should be disposed of in accordance with the local laws and regulations, and shall not be discarded as useless together with household waste.				


ECG report

14.

Reporting time : 2021-06-03 11:15:20
Confirm and sign:

Interpretations :
Sinus rhythm
Normal ECG

ID : 36212120090
Name : Jonas Jonaitis
Gender : M
Age : 59 Years
Dept : 92 kg
Bed No : 180 cm
HR : 80 bpm
PR : 162 ms
QRS : 84 ms
QT/QTc : 356/391 ms
P/QRS/T : 51/31/48 °
RV5/SV1 : 0.964/0.461 mv
RV5+SV1 : 1.425 mv



ECG report

ID : 36212120090
 Name : Jonas Jonaitis
 Gender : M
 Age : 59 Years
 Dept : 92 kg
 Bed No : 180 cm

Reporting time : 2021-06-03 11:15:21
 Confirm and sign:


	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
P Ons	264	264	264	264	264	264	264	264	264	264	264	264
P Dur	94	94	94	94	94	94	94	94	94	94	94	94
QRS Ons	426	426	426	426	430	426	440	442	428	426	426	426
QRS Dur	60	82	80	78	70	82	68	68	82	84	84	68
Q Dur	16	18	21	0	0	19	0	0	12	15	17	17
R Dur	43	40	34	17	70	38	33	32	39	42	44	50
S Dur	0	23	24	41	0	24	34	35	29	25	22	0
ST Dur	174	128	128	152	172	150	166	156	122	126	144	162
T Ons	660	636	634	656	672	658	674	666	632	636	654	656
T Dur	122	146	148	126	110	124	94	116	150	146	128	126
QRS IntD	36	36	37	62	34	37	22	20	34	39	39	39
Q Amp	-49	-143	-94	0	0	-118	0	0	-35	-73	-103	-80
R Amp	582	834	250	96	166	542	184	381	847	1042	964	687
S Amp	0	-173	-130	-707	0	-151	-461	-768	-571	-236	-86	0
ST Amp	7	14	6	-10	0	11	0	3	9	10	13	7
QRS Area	694	804	109	-747	295	458	-223	-280	526	1085	1103	802
P Area	173	251	77	-211	49	166	108	135	172	291	228	182
T Area	649	951	297	-794	182	630	201	408	882	1215	929	747
P Morph	1	1	1	-1	1	1	1	1	1	1	1	1
T Morph	1	1	1	-1	1	1	1	1	1	1	1	1
STM Amp	3	5	2	-4	0	4	0	1	1	2	4	3
ST60 Amp	7	9	1	-8	3	5	0	3	8	10	9	7
STTMidAmp	2	5	3	-3	0	4	0	3	2	1	4	0

Examination time: 2021-06-03 11:10:02

ECG report

ID : 36212120090
Name : Jonas Jonaitis
Gender : M
Age : 59 Years
Dept : 92 kg
Bed No : 180 cm

Reporting time : 2021-06-03 11:15:22
Confirm and sign:



25mm/s 10mm/mv

Examination time: 2021-06-03 11:10:02


ECG report

Reporting time: 2021-06-03 11:15:25

ID : 36212120090
Name : Jonas Jonaitis
Gender : M
Age : 59 Years
Dept : 92 kg
Bed No : 180 cm
HR : 113 bpm
PR : 154 ms
QRS : 80 ms
QT/QTc : 290/383 ms
P/QRS/T : 79/-37/109 °
RV5/SV1 : 1.609/1.515 mv
RV5+SV1 : 3.124 mv

Interpretations :
Sinus tachycardia with sinus arrhythmia
Possible right atrial abnormality
Possible left anterior fascicular block
Left ventricular hypertrophy
Widespread ST-T abnormality ~ may be due to the hypertrophy and/or ischemia
Abnormal ECG

Confirm and sign:



ECG report

Reporting time: 2021-06-03 11:15:27
 ID : 36212120090
 Name : Jonas Jonaitis
 Gender : M
 Age : 59 Years
 Dept : 92 kg
 Bed No : 180 cm


	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
P Ons	280	280	280	280	280	280	280	280	280	280	280	280
P Dur	106	106	106	106	106	106	106	106	106	106	106	106
QRS Ons	444	434	434	436	436	434	434	444	438	440	440	440
QRS Dur	52	78	76	62	72	76	74	66	74	74	74	46
Q Dur	0	0	0	45	11	0	0	0	0	0	0	0
R Dur	52	40	14	16	60	15	11	13	43	43	44	46
S Dur	0	37	61	0	0	15	62	52	30	30	29	0
ST Dur	150	126	138	148	140	116	140	128	118	120	128	148
T Ons	646	638	648	646	648	626	648	638	630	634	642	634
T Dur	72	82	76	74	76	98	52	70	90	86	78	86
QRS IntD	22	34	10	51	41	32	6	10	24	26	26	28
Q Amp	0	0	0	-633	-28	0	0	0	0	0	0	0
R Amp	790	481	62	181	904	64	28	34	713	1528	1609	1414
S Amp	0	-542	-1117	0	0	-160	-1515	-635	-658	-553	-367	0
ST Amp	-56	-87	-31	71	-18	-59	51	-43	-122	-145	-133	-97
QRS Area	1029	-235	-1240	-385	1128	-736	-1826	-842	158	1280	1480	1456
P Area	1	476	475	-238	-236	474	-266	-434	-347	-158	19	94
T Area	-645	99	963	262	-779	448	749	85	-489	-1159	-1219	-1042
P Morph	1	1	1	-1	-1	1	-1	1	1	1	1	1
T Morph	-1	1	1	2	-1	1	1	1	-2	-1	-1	-1
STM Amp	-46	-83	-36	64	-5	-59	73	-28	-91	-129	-126	-105
ST60 AmI	-66	-55	10	60	-38	-22	75	-27	-98	-143	-140	-118
STTMidAr	-50	-84	-32	64	-1	-62	74	-27	-91	-128	-125	-104

Confirm and sign:

Examination time: 2021-06-03 11:08:26

ECG report


Reporting time: 2021-06-03 11:15:27
ID : 36212120090
Name : Jonas Jonaitis
Gender : M
Age : 59 Years
Dept : 92 kg
Bed No : 180 cm



25mm/s 10mm/mv

Examination time: 2021-06-03 11:08:26

Confirm and sign:



ID: 81

09-06-2021
11:11:14

Real-time sampling-Auto-Simultaneous Paperless



Patient Info ZScore's 12x1 Freeze Report Application

F1

F2

F3

F4

F5

F6


F7



iMAC 300/ iMAC 120/ iMAC 12 trolley characteristics and material


Code No.: 3902800106

Full View



Height Range: 90cm - 110cm

Size




Package size: 750 * 502 * 235 (mm)

G.W.: approx. 7.1 kg

N.W.: approx. 5.6 kg

Material



Notice: Product color matching maybe different with the above one, please refer to the actual received goods.

DECLARATION OF CONFORMITY

Manufacturer: *Wuhan Zoncare Bio-Medical Electronics Co., Ltd.*

380, High-tech 2nd road, Eastlake High-tech Development Zone, Wuhan, Hubei (430206), China.

Product: *Trolley (Medical use)*

UMDNS (Name/Code): *Carts, Instrument/10641*


Classification: *Class I, Rule 1 according to Annex IX of the MDD.*

WE HEREWITH DECLARE THAT the ECG Trolley(for models iMAC 300) conform to the safety, reliability and product design criteria referenced in regulations MDD 93/42/EEC and conform with EN 60601-1, 3rd edition.

Place, Date of Issue: WUHAN, Dec 10, 2020

Signature: 蔡海学 (Paul Tsai)


Position: Sales Director



ECG Cable Specifications



1. Model : ECG-FD12
2. Code : 3902800031
3. Material and specifications :




1	Plug		DB 26P
2	PVC	20g	55P gray-white color PVC
3	Gum	10g	LD-PE natural color
4	Copper foil	1PCS	W30*L60mm
5	Screws	2PCS	Lock bolt
6	Wire	2210mm	Ø6.5
7	Printed board	1SET	PCB Board
8	Gum	50g	
9	Gum	100g	
10	Wire	2*500mm	Ø3.2
11	Wire	2*650mm	Ø3.2
12	Internal mold	5g	PP natural color
13	PVC	4-10g	55P White PVC
14	Resistance	10PCS	Anti-defibrillation resistance 25.1
15	PVC	6-10g	55P White PVC
16	Color label set	5g	85P, 10 characters
17	Pins	10PCS	Ø4, H62 nickel plating
18	Clip chain	1PCS	W190*L220*0.1mm

19	Chip	1PCS	DS2431
20	product label	1PCS	Transparent PVC
21	Identification label	1PCS	PVC cold gray 3C
22	CE label	1PCS	
23	Magic cable tie	1PCS	1.5×15 Gray

4. Packing

Package Method



Limb electrodes, 4 units 25.3



Chest electrodes, 6 units 25.2



Confirm Patient Info Previous

Input List Settings

ID: 20210330131447

Name: Age setting ✕


Gender: 3 | Days

Dept: Birth date Year Months D

Race: Unknown

Medicine: 0 Undefined

PREV DIAG: 255



Neonatal ECG cable and disposable electrodes



EC CERTIFICATE

Full Quality Assurance System

Certificate No.:
12351-2018-CE-RGC-NA-PS Rev. 2.0

Project No.:
PRJC-217612-2010-PRC-CHN

Valid Until:
27 May 2024

This is to certify that the quality system of:

Wuhan Zoncare Bio-medical Electronics Co., Ltd.

#380, High-tech 2nd Road, Eastlake high-tech district, Wuhan, Hubei, P. R. China.

For design, production and final product inspection/testing of:

DIGITAL MULTI-CHANNEL ELECTROCARDIOGRAPH

Has been assessed with respect to:

THE CONFORMITY ASSESSMENT PROCEDURE DESCRIBED IN ANNEX II EXCLUDING SECTION 4 OF COUNCIL DIRECTIVE 93/42/EEC ON MEDICAL DEVICES, AS AMENDED

and found to comply.

Further details of the product(s) and conditions for certification are given overleaf.

Place and date:
Høvik, 06 January 2021

For:
DNV GL PRESAFE AS
Notified Body No.: 2460

Sholeh Gheissar

The certificate is digitally verified by blockchain technology. For more info, see www.dnvgl.com/assurance/certificates-in-the-blockchain.html



Certificate No.:
12351-2018-CE-RGC-NA-PS Rev. 2.0

Project No.:
PRJC-217612-2010-PRC-CHN

Valid Until:
27 May 2024

Jurisdiction

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as "Forskrift om Medisinsk Utstyr" by the Norwegian Ministry of Health and Care Services.

Certificate history:

Revision	Description	Issue Date
0.0	Original Certificate	2018-04-20
1.0	Extension in scope	2020-06-17
2.0	Recertification	2021-01-06

Products covered by this Certificate:

Product Description	Product Name	Class
Digital Multi-Channel Electrocardiograph	iMAC 120, iMAC 120pro	IIa

The complete list of devices is filed with the Notified Body

Sites covered by this certificate

Site Name	Address
Wuhan Zoncare Bio-medical Electronics Co., Ltd.	#380, High-tech 2nd Road, Eastlake high-tech district, Wuhan, Hubei, P. R. China.

EU Representative

Well kang Limited , The Black Church, St. Mary's Place , Dublin 7, Ireland

Certificate No.:
12351-2018-CE-RGC-NA-PS Rev. 2.0

Project No.:
PRJC-217612-2010-PRC-CHN

Valid Until:
27 May 2024

Terms and conditions

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform Presafe of any intended updating of the quality system and Presafe will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies the quality system. Presafe reserves the right, on a spot basis or based on suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.

Conformity declaration and marking of product

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number of Presafe.

End of Certificate

1. Declaration of Conformity

Declaration of Conformity

For the following products:

Digital Multi-channel Electrocardiograph

(Product Name)

iMAC 120, iMAC 120A, iMAC 120B, iMAC 120pro

(Model Designation)

is hereinafter confirmed to comply with the requirements set out in the Council Directive on the harmonization of the Laws of the Member States concerning Medical Device Directive (93/42/EEC As amended by 2007/47/EC)

EN ISO 13485:2016, EN ISO 14971:2012, EN ISO 15223-1:2016, EN 1041:2008, ISO 10993-1:2009, ISO 10993-5:2009, ISO 10993-10:2010, IEC 60601-1:2005/A1:2012, EN 60601-1-2:2015, IEC 60601-2-25:2011, EN 62304:2006/A1:2015, EN 62366-1:2015, EN 60601-1-6:2010/A1:2015.

Conformity Assessment Route:

Annex II excluding section 4 of Medical Device Directive

Notified Body:

DNV GL Presafe AS (NB No. 2460)

Veritasveien 3, 1363 Høvik, Norway

The following representative in Europe is responsible for making this declaration:

Company Name: Well kang Limited

Company Address: The Black Church, St. Mary's Place, Dublin 7, Ireland

The following manufacturer is responsible for making this declaration:

Company Name: Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Company Address: #380, High-tech 2nd road, Eastlake high-tech district, Wuhan, Hubei, P. R. China

Qi Wang

General Manager

(Legal Signature)

(Position/title)



(Date)



Notified Body Confirmation Letter Reference: C684165

To whom it may concern,

Confirmation of the status of a formal application, written agreement, and appropriate surveillance in the framework of Regulation EU 2023/607 amending Regulations (EU) 2017/745 and (EU) 2017/746 as regards the transitional provisions for certain medical devices and in vitro diagnostic medical devices

This letter confirms that, DNV Product Assurance AS, a Notified Body (NB) designated against Regulation (EU) 2017/745 (MDR) and identified by the number NB 2460 on NANDO, has received a formal application in accordance with Section 4.3, first subparagraph of Annex VII of MDR and has signed a written agreement in accordance with Section 4.3, second subparagraph of Annex VII of MDR with the following manufacturer:

Wuhan Zoncare Bio-medical Electronics Co., Ltd.

#380, High-tech 2nd Road, Eastlake high-tech district, Wuhan, Hubei, P.R. China.

SRN Number: CN-MF-000010439

The devices covered by the formal application and the written agreement mentioned above are identified in the Tables below. Table 1 identifies the devices for which an MDR application has been received, written agreement concluded and for which the NB is also responsible for appropriate surveillance of the corresponding devices under the applicable Directive. Table 2 identifies the devices for which an MDR application has been received and a written agreement concluded, but the NB has not yet taken the responsibility for appropriate surveillance of the corresponding devices under the applicable Directive.

In the case of devices covered by certificates issued under Directive 93/42/EEC (MDD) that expired after 26 May 2021 and before 20 March 2023, without having been withdrawn, this letter also confirms that the manufacturer signed the written agreement under MDR by the date of MDD certificate expiry; or provided evidence that a competent authority of a Member State had granted a derogation or exemption from the applicable conformity assessment procedure in accordance with Article 59(1) of MDR or Article 97(1) of the MDR respectively, by the 20 Mar 2023 for the relevant devices.

Place and date:
Høvik, 2024/04/17



For the issuing office:
DNV Product Assurance AS – Notified Body 2460
Veritasveien 1, 1363 Høvik, Norway

André Fernandes
Management Representative

Lack of fulfilment of conditions as set out in the Certification Agreement may render this letter invalid.

NOTIFIED BODY 2460: DNV Product Assurance AS, Veritasveien 1, 1363 Høvik, Norway, Tel +47 67 57 88 00, www.dnv.com

The transition timelines that apply to the devices covered by this letter, subject to the manufacturer's continued compliance to the other conditions specified in Article 120.3c of MDR (as amended by (EU) 2023/607), are shown below:

- 26 May 2026 for Class III custom-made implantable devices
- 31 December 2027 for Class III devices and Class IIb implantable devices excluding Well-established technologies (WET - sutures, staples, dental fillings, dental braces, tooth crowns, screws, wedges, plates, wires, pins, clips, and connectors)
- 31 December 2028 for other Class IIb devices, Class IIa, Class I devices placed on the market in sterile condition or have a measuring function.
- 31 December 2028 for devices not requiring the involvement of a notified body under MDD but requiring it under MDR (e.g., class I devices that qualify as re-usable surgical instruments)

Table 1: Devices covered by this letter and for which the NB is also responsible for appropriate surveillance of the corresponding devices under the applicable Directive:

Device name and Basic UDI-DI (under MDR application)	MDR Device classification (as proposed by the manufacturer and verified at the quotation request review stage)	If the MDR device is a substitute device, identification of the corresponding MDD device	MDD Certificate Reference(s) of the devices under MDR application, and the NB Identification
Device name: Digital Multi-channel Electrocardiograph iMAC 12, iMAC 12pro Basic UDI-DI: 695691348214MQ	Class IIa	NA	Certificate number: 12350-2018-CE-RGC-NA-PS Rev. 2.0 NB name: DNV Product Assurance AS, NB number: 2460
Device name: Digital Multi-channel Electrocardiograph iMAC 120, iMAC 120pro Basic UDI-DI: 695691348214MQ	Class IIa	NA	Certificate number: 12351-2018-CE-RGC-NA-PS Rev. 2.0 NB name: DNV Product Assurance AS, NB number: 2460
Device name: Digital Multi-channel Electrocardiograph iMAC 300, iMAC 300pro Basic UDI-DI: 695691348214MQ	Class IIa	NA	Certificate number: 12349-2018-CE-RGC-NA-PS Rev. 2.0 NB name: DNV Product Assurance AS, NB number: 2460

Device name and Basic UDI-DI (under MDR application)	MDR Device classification (as proposed by the manufacturer and verified at the quotation request review stage)	If the MDR device is a substitute device, identification of the corresponding MDD device	MDD Certificate Reference(s) of the devices under MDR application, and the NB Identification
Device name: Digital Multi-channel Electrocardiograph iMAC 1800, iMAC 1800pro Basic UDI-DI: 695691348217MW	Class IIa	NA	Certificate number: 9861-2017-CE-RGC-NA-PS Rev. 3.0 NB name: DNV Product Assurance AS, NB number: 2460

Table 2: Devices covered by this letter and for which the NB is NOT responsible for appropriate surveillance of the corresponding devices under the applicable Directive:

Device name and Basic UDI-DI (under MDR application)	MDR Device classification (as proposed by the manufacturer and verified at the pre-application stage)	If the MDR device is a substitute device, identification of the corresponding MDD device	MDD Certificate Reference(s) of the devices under MDR application, and the NB Identification
N/A			

Confirmation Letter Revision History

Date	NB internal reference traceable to each version of the letter	Action
2024/04/17	C684165	Initial issue

Lack of fulfilment of conditions

The following may render this letter of confirmation invalid:

- Lack of compliance to the requirements of Regulation (EU) 2023/607.
- Significant changes to design or intended purpose of the devices.
- Changes in the quality system affecting production.
- Periodical audits not held within the timeframe.

12. Statement on iMAC ECG Noise Reduction

To whom it may concern,

We, Wuhan Zoncare Bio-medical Electronics Co., Ltd., manufacturer of ECG Machine, Ultrasound Machine, Patient Monitor makes the statement regarding our iMAC 120, iMAC 12, iMAC 300 ECG machines as below:

In the case that some specific tenders require the indication of ECG noise reduction function/technology, while it's not included in our technical specifications such as datasheet, user manual, etc., we hereby state that the iMAC ECGs obtained two patents for inventions relating to noise filtering/removal:

- ◆ Apparatus and method for filtering out industrial frequency noise(Patent No. ZL 2017 1 0082475.3).
- ◆ An ECG signal filtering method for removing ringing(Patent No. ZL 2018 1 1443998.7).

The original patent certificate are attached.

This parameter is for reference and Zoncare holds the final right of explanation.



Date: July. 26th 2023



Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Add: 380, High-tech 2nd road, Eastlake High-tech Development Zone, Wuhan, Hubei (430206), China.

Tel: 86-27-87770581 Email: info@zoncare.com

Fax: 86-27-87770203 Website: www.zoncare.com

Service Certificate

We, Wuhan Zoncare Bio-medical Electronics Co., Ltd., located at the following address: #380, High-tech 2nd road, Eastlake high-tech development zone, Wuhan, Hubei, P. R. China hereby certify that, Dr. Raimondas Binkys of Litfarma, Jeruzales 45., Vilnius Lithuania, has passed the training courses and is authorized for installation, training of the personnel for work with the equipment and to perform warranty and after warranty service, repair, maintenance for the following models:

Digital Electrocardiograph:

iMAC 300, iMAC 12, iMAC 120, iMAC 1200

Wuhan Zoncare Bio-medical Electronics Co., Ltd.

Jan 3rd, 2022

Signature: 严谨 Johnny Yan

Title: Area Sales Director

