
Edwards SAPIEN 3 valve with Alterra adaptive prestant

Technical Sheet

Edwards SAPIEN 3 Transcatheter Pulmonary Valve System with Alterra Adaptive Prestant.


Device description

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant consists of the Edwards 29 mm SAPIEN 3 transcatheter heart valve (THV), the Edwards SAPIEN 3 pulmonic delivery system (PDS), the Edwards Alterra adaptive prestant system and accessories.

Edwards Alterra Adaptive Prestant System consists of an Alterra adaptive prestant that is fully loaded in an Alterra delivery system and supplied together in one package.

Edwards Alterra Adaptive Prestant is used as a docking adaptor for the 29 mm Edwards SAPIEN 3 transcatheter heart valve (THV). It is comprised of a self-expanding, radiopaque, nitinol frame assembly with polyethylene terephthalate (PET) fabric covering. The prestant has designated inflow and outflow sides. The proximal inflow section is identifiable by the presence of two triangular tabs (prestant connector) that are attached to the catheter of the delivery system. The distal outflow section is distinguished by the open cells for blood flow. The PET fabric is attached by sutures to the inside surface of the frame to create sealing at the inflow section and opening for the outflow. Sutures are also used in the center to support the middle section when an Edwards SAPIEN 3 transcatheter heart valve is implanted. Three (3) fluoroscopically visible radiopaque markers at the prestant waist help with positioning.

Edwards Alterra adaptive prestant

Designed to provide a landing zone for 29 mm SAPIEN 3 valve

Prestant sizing in RVOT landing zone

Perimeter	84.9 - 119.3 mm
Perimeter derived diameter*	27 - 38 mm
Prestant size	40 x 49 mm
SAPIEN 3 valve compatibility	29 mm

* Diameter range throughout cardiac cycle

Edwards Alterra adaptive prestant dimensions

Inflow sealing OD	39 mm
Outflow sealing OD	41 mm
Height	49 mm

Patient sizing range

Sizing recommendation for the prestant in the right ventricular outflow tract/pulmonary valve (RVOT/PV) landing zone are shown in the table on top.

Edwards

Intended use

The bioprosthesis with prestant is intended for use in patients requiring pulmonary heart valve replacement. The delivery systems and accessories are intended to facilitate placement of the bioprosthesis and prestant via the transfemoral access approach.

Indications

The 29 mm Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant is indicated for use in the management of patients with pulmonary regurgitation who have a native or surgically-repaired right ventricular outflow tract and are clinically indicated for pulmonary valve replacement.

Product Description	29 mm
Edwards SAPIEN 3 Transcatheter Pulmonary Valve System	S3PAS4045CE
Edwards SAPIEN 3 Transcatheter Heart Valve	9600TFX29
Edwards SAPIEN 3 Pulmonic Delivery System	9630PL29CE
Edwards Crimper	9600CR
Edwards eSheath+ Introducer Set (16F)	916ESP
Locking Syringe	96406
Alterra Adaptive Prestent with Alterra Delivery System	29AP4045CE


Edwards Alterra Delivery System

The delivery system includes a handle which consists of a wheel that allows for deployment, two primary shafts with a flush port to flush the delivery system, and a long tapered tip at the distal end to facilitate tracking through the vasculature. A fluoroscopically visible radiopaque marker band shows the location of the tip of the outer shaft. The prestant is fully loaded in the delivery system. A stylet is included within the guidewire lumen.

Delivery system dimensions


<u>Alterra capsule OD</u>	6.7 mm / 20 Fr
<u>Alterra shaft OD</u>	5.7 mm / 17Fr
<u>Wire compatibility</u>	0.035" / 0.89mm

Alterra adaptive prestant is compatible with the 16Fr eSheath+ introducer set or equivalent

Edwards SAPIEN 3 Transcatheter Heart Valve

The Edwards SAPIEN 3 transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) fabric skirt. The leaflets are treated according to the Carpentier-Edwards ThermaFix process.

	29 mm
Tissue	Bovine pericardium, Thermafix tissue process ¹
Frame	Cobalt-chromium alloy
Crimped height ²	31 mm
Expanded height	22.8 mm
Foreshortening	8.2 mm


¹ No clinical data are available with which to evaluate the long-term impact of the Carpentier-Edwards ThermaFix tissue process in patients

Edwards SAPIEN 3 Pulmonic Delivery System (PDS)

The Edwards SAPIEN 3 pulmonic delivery system facilitates the placement of the bioprosthesis. It consists of an inline sheath, balloon catheter for deployment of the Edwards SAPIEN 3 transcatheter heart valve, and an outer shaft and valve capsule to cover the transcatheter heart valve during insertion and tracking to the intended deployment location. The delivery system includes a tapered tip to facilitate crossing of right heart structures. The valve capsule and tapered tip are hydrophilic coated. A visual balloon shaft marker is provided to assist with balloon recapture. A stylet is included within the guidewire lumen of the delivery system. The 28 F hydrophilic coated dilator (packaged with the delivery system) is used to predilate the vessel prior to insertion of the delivery system, if necessary.

Pulmonic delivery system dimensions

Catheter capsule outer diameter	28.5 Fr / 9.5 mm
In-line sheath outer diameter	25.5 Fr / 8.5 mm
Catheter outer shaft diameter	17 Fr / 5.7 mm
Guidewire compatibility	0.035" / 0.89mm

How Supplied

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system, sheath, and crimper are supplied sterilized with ethylene oxide gas.

The Edwards Alterra adaptive prestant system is supplied pouched and sterilized by e-beam sterilization.

The THV is supplied nonpyrogenic packaged in buffered glutaraldehyde, in a plastic jar to which a tamper evident seal has been applied. Each jar is shipped in a shelf box containing a temperature indicator to detect exposure of the THV to extreme temperature. The shelf box is enclosed in Styrofoam prior to shipping.

Storage

The transcatheter heart valve must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the THV to extreme temperature. The delivery system and accessories should be stored in a cool, dry place. The prestant and delivery system must be stored in a cool, dry place.

MR Safety

MR Conditional

Non-clinical testing has demonstrated that the Non-clinical testing has demonstrated that the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is MR Conditional. A patient can be scanned safely immediately after placement of this implant in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 Tesla and 3.0 Tesla
- Spatial magnetic gradient field of 3000 Gauss/cm (30 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) scanning per sequence
- Gradient system is in normal operating mode

Under the scan conditions defined above, the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is expected to produce a maximum temperature rise of 4.0 °C or less after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 15 mm for gradient echo images when scanned using a 3.0 T MRI system. The artifact obscures the device lumen in spin and gradient echo images. The delivery system has not been evaluated for MR compatibility and is considered MR unsafe.

Qualitative and Quantitative Information related to the valve

The Alterra adaptive prestant system and SAPIEN 3 transcatheter heart valve contain the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0 Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

For the qualitative and quantitative information on the materials and substances constituting the valve please refer to the IFU.

Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information. The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the Alterra platform and the SAPIEN 3 platform.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, Alterra, Carpentier-Edwards, Edwards eSheath, Edwards SAPIEN, Edwards SAPIEN 3, eSheath, SAPIEN, SAPIEN 3, and TheraFix are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2024 Edwards Lifesciences Corporation. All rights reserved. PP-EU-9104 v1.0


Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com

Edwards

Proven technology that extends SAPIEN 3 valve benefits to more patients

The Alterra adaptive prestant—accurate, secure placement for an expanded group of patients


SAPIEN 3 valve with the Alterra adaptive prestant


Proven technology that extends SAPIEN 3 valve benefits to more patients

The Allterra adaptive prestent—accurate, secure placement for an expanded group of patients

A minimally invasive solution that helps minimize reinterventions.

Deliver excellent outcomes, delay potential reintervention, and improve quality of life³ for more of your pulmonary patients.⁷

References:

1. Hinderbach et al. Medical care and prevention in ACHD. *Cardiovasc Diagn Ther*. 2018;8(6):705-715
2. Haseck et al. Outcomes of transcatheter pulmonary SAPIEN 3 valve implantations: an international registry. *Eur Heart J*. 2024;45(2):198-210
3. Dimaev et al. Transcatheter Pulmonary Valve Implantation with the Allterra Adaptive Prestent, SAPIEN 3 Transcatheter Heart Valve and Pulmonary Delivery System: Two-year Follow-up Outcomes of the ALLTERRA Trial. Presented at PICS 2024, San Diego, CA
4. Aloufiani, Jamil. Five-year Results from the COMPRESSION III Trial of Treatment of Patients with a Dysplastic Aortic PVD: Conduit or Previously Implanted Pulmonic Valve with the SAPIEN 3 Transcatheter Heart Valve. Presented at PICS 2024, San Diego, CA
5. Edwards Lifesciences, Allterra adaptive prestent EU
6. Edwards Lifesciences SAPIEN 3 Transcatheter Pulmonary Valve System with Allterra Adaptive Prestent EU
7. FDA Summary of Safety and Effectiveness Data, December 2021

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the instructions for use (consult eu.edwards.com where applicable).

© 2024 Edwards Lifesciences Corporation. All rights reserved. PPS-10-1897-A-10
Edwards Lifesciences S.A. Route de Hérisay 20, 1202 Nyon, Switzerland • edwards.com


Unique patients Excellent outcomes

The Edwards SAPIEN 3 Pulmonic Portfolio—where heart valve innovation and living life to the fullest meet

Built on a proven platform; designed for performance

With pulmonary patients living longer¹ and having higher expectations for their health and quality of life, lifetime management is increasingly important. That's why valve choice matters, and the Edwards SAPIEN 3 pulmonary portfolio delivers.^{2,3}

Excellent clinical outcomes, patient after patient

	SAPIEN 3 valve at 5 years ¹ COMPASSION 3 Trial	SAPIEN 3 valve with the Alterra adaptive present at 2 years ¹ Alterra Pivotal Trial
Mortality	0% (n=42)	1.7% (n=188)
≤ Mild pulmonary regurgitation	95% (n=40)	92% (n=45)
Reintervention	2% (n=41)	2% (n=188)

Minimizing complications that impact the need for reintervention

	SAPIEN 3 valve at 5 years ¹ COMPASSION 3 Trial	SAPIEN 3 valve with the Alterra adaptive present at 2 years ¹ Alterra Pivotal Trial
Endocarditis	3.8% (n=41)	0% (n=188)
Stent Thrombosis	4.3% (n=18)	4.2% (n=188)
Frame fracture	0%	0.8% (n=4)
Stent fracture requiring reintervention	—	0% (n=188)

¹In patients with aortic and/or mitral outflow tract (AVRT) or aortic and/or mitral outflow tract (MVRT) patients with aortic and/or mitral outflow tract (AVRT/MVRT) patients with aortic and/or mitral outflow tract (AVRT/MVRT).

Delivery systems that offer controlled, predictable, and precise positioning
Ready to adapt and adjust based on patient needs.

Alterra Delivery System
Recapture and reposition* with ergonomic control, smooth tracking, and a system developed specifically for the Alterra adaptive present.

Pulmonic Delivery System
Covered system* designed to track through complex anatomy and deliver the SAPIEN 3 valve into the Alterra adaptive present.

Edwards Commander Delivery System
Dual articulation designed for precise positioning and deployment of the Edwards SAPIEN 3 valve.

Clinical Summary

Multicenter Pivotal Study of the Alterra Adaptive Presept for the Treatment of Pulmonary Regurgitation

Dimas VV, Babaliaros V, Kim D, et al. Multicenter Pivotal Study of the Alterra Adaptive Presept for the Treatment of Pulmonary Regurgitation. *JACC Cardiovasc Interv.* 2024;17(19):2287-2297.

Study aim

The Alterra adaptive presept was designed to expand transcatheter therapy to treat patients with dilated right ventricular outflow tract (RVOT). The study aims to report **2 years outcomes** of the ALTERRA pivotal trial in dilated native RVOTs.

Methodology

Prospective, single-arm, multicenter study design to assess the safety and effectiveness of the Alterra adaptive presept with the 29 mm SAPIEN 3 valve in patients with dysfunctional RVOT requiring treatment of moderate or severe pulmonary regurgitation (PR).

Patients characteristics

- Median age: 23years
- Gender: 55.7% male
- NYHA class: 54.1% NYHA class I
- Diagnose: Tetralogy of Fallot (70.5%); Pulmonic valve stenosis (23%); Balloon valvuloplasty (6 patients) and Surgical valvotomy (8).

The primary end point was transcatheter heart valve (THV) dysfunction at 6 months, defined as a nonhierarchical composite of RVOT/pulmonary valve gradient 35 mm Hg or greater on transthoracic echocardiography. Patients enrolled, n=60

Results

The Alterra adaptive presept had shown excellent procedural outcomes achieving **94.9% of device success**, composite of:

- Single presept deployed in the desired location (100%)
- Single THV implanted in the desired location within the presept (96.7%)
- RV-PA peak-to-peak gradient less than 35mm Hg post-implant (100%)
- Less than moderate PR on discharge TTE (or earliest evaluable TTE) (98.3%)
- Free of explantation at 24 h post implantation (100%)


Edwards

Excellent clinical outcomes were reported at 2 years follow-up. Key outcomes reported below.

2-Year Outcomes of the ALTERRA Pivotal Trial

Prospective, single-Arm, Multicenter Study of Patients with Dysfunctional Native or Patched RVOT and \geq Moderate PR by TTE

Trial Devices: Alterra adaptive prestant and 29 mm SAPIEN 3 THV	Screen-Fail Rate for Patient Population: 37.1% (36/97) N=60 patients implanted	2-Year Follow-Up: 96.7% (58/60)		
	Both devices implanted in a single procedure	98.3%	Key Outcomes at 2 Years	
	Free of explant within 24 h	100%	All cause mortality	0%
	Primary Endpoint THV dysfunction at 6 months	0% (0/59)	Fracture requiring intervention	0%
	RVOT reintervention	0	Coronary artery compression	0%
	Moderate or greater PR	0	Endocarditis	0%
	Mean gradient \geq 35 mm Hg	0	Sustained VT	0%
			RVOT/PV reintervention	1.7%
		\geq Moderate PR	7.5%	

- At 2 years, 92.5% of patients treated in the Alterra Pivotal Trial had \leq mild PR
- No deaths, endocarditis, device embolization, or explants

Total pulmonary valve regurgitation (PR) and mean right ventricular outflow tract (RVOT) gradient were assessed by a transthoracic echocardiography (TTE) core laboratory | ALTERRA = ALTERRA Pivotal Trial Multicenter Study of Congenital Pulmonic Valve Dysfunction Studying the SAPIEN 3 THV with the Alterra adaptive prestant; PV = pulmonic valve; THV = transcatheter heart valve; VT = ventricular tachycardia.

Conclusion

It is extremely favorable to see **sustained excellent clinical outcomes out to 2 years**. The patient cohort will continue to be monitored through annual follow-up assessment up to 5 years.

As of today, we can state that dilated RVOT presents unique challenges for transcatheter management. However, what we are seeing is that not only the Alterra adaptive prestant allows to expand the transcatheter therapy for this dilated RVOT patients population but also, *“the Alterra adaptive prestant and 29 mm SAPIEN 3 THV resulted in restoration of PV competency and can be safely and effectively performed, with improvement in quality of life, as well as favorable remodeling of the right ventricle.”*

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, Alterra, SAPIEN and SAPIEN 3 are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2024 Edwards Lifesciences Corporation. All rights reserved. PP--EU-9442 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com

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EC REP

Edwards Lifesciences Services GmbH
Edisonstr. 6
85716 Unterschleissheim
Germany

2022-06
10051932001 A
© Copyright 2022, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands



Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU

Edwards eSheath introducer set

The Dynamic Expansion Mechanism (DEM) is designed to reduce vascular trauma by

- Allowing for transient sheath expansion during delivery system passage
- Reduces the time the access vessel is expanded


¹ The SAPIEN 3 valve is retrievable prior to THV deployment. THV and sheath are removed together with THV in the body of the sheath.

74

Edwards eSheath introducer set*

- For 23 and 26 mm THV, no color coding for 14F sheath
- Sheath size is labeled on sheath handle
- For 29 mm THV, 16F sheath is either color coded **green** (not **orange**) or no color and sheath size is labeled on sheath handle

Note

An additional dilator may or may not be included with the system for vessel dilation. If the same dilator is used for vessel dilation and sheath insertion, check to ensure dilator has not been damaged before inserting into sheath.

THV	Sheath ID (unexpanded)	Sheath OD (unexpanded)	Minimum vessel diameter
23 mm SAPIEN 3 valve	14F (4.6 mm)	6.0 mm	5.5 mm
26 mm SAPIEN 3 valve	14F (4.6 mm)	6.0 mm	5.5 mm
29 mm SAPIEN 3 valve	16F (5.3 mm)	6.7 mm	6.0 mm
23 mm SAPIEN XT valve	16F (5.3 mm)	6.7 mm	6.0 mm
26 mm SAPIEN XT valve	18F (5.9 mm)	7.2 mm	6.5 mm
29 mm SAPIEN XT valve	20F (6.6 mm)	8.0 mm	7.0 mm

*Refer to the IFU for any contraindications

15

Edwards eSheath introducer set

Tri-seal valve technology designed for hemostasis


Seal	Hemostasis
 <p>Disc valve</p>	<p>Designed to provide hemostasis with a 5F or greater catheter</p>
 <p>Cross slit valve</p>	<p>Designed to provide hemostasis for guidewires $\geq 0.035''$</p>
 <p>Duckbill valve</p>	<p>Designed to provide hemostasis when nothing is inside sheath</p>

16

Edwards eSheath introducer set

Hydrophilic coating to facilitate easy insertion

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



SAPIEN 3 Transcatheter Heart Valve

SAPIEN XT Valve

	23 mm	26 mm	29 mm	
Tissue	Bovine Pericardium, TheraFix tissue process ¹			
Frame	Cobalt-Chromium Alloy			
Native annulus size	2D TEE diameter	18 – 22 mm	21 – 25 mm	24 – 27 mm
	3D area	314 – 415 mm ²	415 – 530 mm ²	530 – 660 mm ²
	3D area derived diameter	20 – 23 mm	23 – 26 mm	26 – 29 mm
Crimped height	17 mm	20 mm	22 mm	
Expanded height	14 mm	17 mm	19 mm	
Foreshortening	3 mm	3 mm	3 mm	
Inner skirt height³	6.7 mm	8.7 mm	11.6 mm	

¹ No clinical data are available with which to evaluate the long-term impact of the Carpentier-Edwards TheraFix tissue process in patients
² Rounded to the nearest 0.5 mm
³ Measured at bottom of zig-zag

SAPIEN 3 Valve

	23 mm	26 mm	29 mm	
Tissue	Bovine Pericardium, TheraFix tissue process ¹			
Frame	Cobalt-Chromium Alloy			
Native annulus size	2D TEE diameter	18 – 22 mm	21 – 25 mm	24 – 28 mm
	3D area	338 – 430 mm ²	430 – 546 mm ²	540 – 683 mm ²
	3D area derived diameter	20.7 – 23.4 mm	23.4 – 26.4 mm	26.2 – 29.5 mm
Crimped height⁽²⁾	24.5 mm	27 mm	31 mm	
Expanded height	18 mm	20 mm	22.5 mm	
Foreshortening	6.5 mm	7 mm	8.5 mm	
Inner skirt height³	9.3 mm	10.2 mm	11.6 mm	
Outer sealing skirt height	6.6 mm	7.0 mm	8.1 mm	