# PIEZOTOME CUBE

Adopt **PIEZOTOME®** extraction to favor bone preservation with immediate implant placement

## PIEZOTOME® CUBE

with Essential kit

Delivered with a **1.4.** multifunction footswitch, a **1.7.** PIEZOTOME® CUBE LED handpiece, **1.2.** 1 handpiece holder, an **1.3.** Essential III kit (6 tips: BS1S, BS4, SL1, SL2, SL3 et LC2), a **1.** PIEZOTOME® autoclave dynamometric wrench, **1.9.** 1 bracket, **1.5.** 2 sterilizable irrigation lines, **1.6.** 5 clips and 15 perforators.

## PIEZOTOME® CUBE

with Universal kit

Delivered with a multifunction footswitch, a **F50100** PIEZOTOME® CUBE LED handpiece, **F50113** 1 handpiece holder, an Universal kit (6 tips: BS6, LC1, LC2, CE1, CE3 et PZ1), a PIEZOTOME® autoclave dynamometric wrench, **F50113** 1 bracket, 2 sterilizable irrigation lines, 5 clips and 15 perforators.

For more information, please consult our website

Website

[www.acteongroup.com/en/products/equipment](http://www.acteongroup.com/en/products/equipment)

## PIEZOTOME® CUBE ACCESSORIES

**PIEZOTOME® CUBE LED HANDPIECE**

Compatible with 2<sup>nd</sup> generation piezo tips →

Class I/II medical device (GMD) - CE 0459 - Manufacturer: SATELEC® - France

**F12816**

**MULTIFUNCTION FOOTSWITCH  
PIEZOTOME® CUBE**

**F50109**

**HANDPIECE HOLDER**

**F50133**

**PIEZOTOME® AUTOCLAVABLE  
DYNAMOMETRIC WRENCH**

**F50150**

**IRRIGATION LINE CLIPS x10**

**F50116**

**BRACKET FOR PIEZOTOME® CUBE LED**

**F50129**

**LIGHT RING  
PIEZOTOME® CUBE LED handpiece**

**F12828**

## FOR FURTHER INFORMATION

Piezotome® Cube  
& Tips  
Catalog  
Ref: D50111

Around Piezotome® Cube video

FAST, SAFE AND MINIMALLY INVASIVE

FOR AN UNMATCHED VISIBILITY

1.8.

**Intralift**

### OSTEOTOMY - PILOT DRILLING


- ▶ **TKW1**  
Conical and diamond-coated tip for pilot drilling in very dense cortical bone (Ø 1.35mm).

### OSTEOTOMY - SECONDARY DRILLING

- ▶ **TKW2**  
Cylindrical and diamond-coated tip dedicated for drilling and fracturing the sinus floor (Ø 2.1mm).
- ▶ **TKW3**  
Cylindrical tip dedicated to drill and widen the access canal to sinus membrane (Ø 2.35mm).
- ▶ **TKW4**  
Cylindrical and diamond-coated tip intended for the receptacle preparation and widening the access canal to the sinus membrane (Ø 2.8mm).

### MEMBRANE DETACHMENT (ELEVATION)

- ▶ **TKW5**  
Tip Ø 3mm to be placed in the receptacle preparation for sinus membrane elevation by crestal approach.

**SinusLift**

### VESTIBULAR BONE WINDOW CUT OSTEOTOMY

- ▶ **SL1**  
Diamond-coated tip for vestibular bone window cut and for attenuation of sharp angles.

### OSTEOPLASTY

- ▶ **SL2**  
Diamond-coated ball tip for smoothing the vestibular bone window and precise osteoplasty.

### MEMBRANE DETACHMENT

- ▶ **SL3**  
Plateau tip served for sinus membrane elevation on the window's edges.
- ▶ **SL4**  
Spatula served for sinus membrane elevation inside the sinus.
- ▶ **SL5**  
Spatula used for sinus membrane elevation inside the sinus and for disengagement of anatomical structures.

# 11 Technical specifications of the medical device

## 11.1 Identification

Manufacturer	SATELEC, a company of Acteon group
Name of the medical device	Piezotome Cube

## 11.2 Generator

Supply voltage	100 - 240 VAC	7.
Power supply frequency	50 / 60 Hz	7.
Power consumption	150 - 150 VA	
Voltage supplied to handpiece	150 VAC	
Power setting range	D1 - D4	5.
Output frequency	28 kHz - 36 kHz	3.
Type of leakage currents	BF	
Operating mode	Intermittent: 10 minutes ON / 5 minutes OFF	
Electrical rating	I	
Fuse (mains connector)	2 fuses T2AL, 250 VAC	
Width	251 mm	
Height	160 mm, 481 mm with bracket	
Depth	271 mm	
Weight	3,500 g without accessories	
Ingress protection rating	IPX0	

## 11.3 Length of cords

Handpiece cord	2 500 mm +/- 50 mm	1.2.
Control pedal cord	2 500 mm +/- 50 mm	7.2.

## 11.4 Irrigation

Bottles or irrigation bags should not weigh more than one kilogram. A heavier container will make the medical device tip over.

Maximum volume of irrigation solution bags	1,000 ml	
Maximum weight of irrigation solution bags	1,000 g	
Nominal water output flow at the end of the handpiece Cube LED	0 ml/min to 120 ml/min	4.
Maximum water output flow at purge	120 ml/min	

## 11.5 Footswitch

Width	173 mm
Height	140 mm including arch
Depth	176 mm
Weight	1,060 g
Ingress protection rating	IPX1

# ACTEON® MAKES THE DIFFERENCE

## CUBE LED HANDPIECE

Upgraded conception based on ACTEON®  
6 ceramic rings

- ▶ More power than ever

Lighter with perfectly balanced design

- ▶ Natural gesture, rotation free for less hand fatigue

## Full white color (LED ring (100,000 Lux) 2.1.

- ▶ Great tissue distinction and posterior areas visibility

## TIPS 6.

- ▶ Many exclusive designs
- ▶ Clinical versatility

Strengthened by surface treatment

- ▶ Excellent durability

## MAKE IT EASY

### ▶ PERISTALTIC PUMP

- Fast and intuitive insertion of the irrigation line in the pump
- Small size and robust
- Precise adjustment (ml/min) and regular flow rate reaching the tip's end

### ▶ EXTERNAL IRRIGATION LINE

- Extend the handpiece lifetime
- Immediate visual control throughout the procedure
- Provide complete sterile work and ease maintenance
- Autoclavable or single use delivered sterile

## DURABILITY

Piezotome® Cube accessories are designed, with robust and durable materials like stainless steel, limiting the use of plastic. A willingness of ACTEON® to offer you a clinical accompaniment over the years.

## FAST ASSEMBLY

### 6.2. TORQUE WRENCH

- Saves precious time
- Robust stainless steel
- Dynamometric wrench

## TOTAL ADAPTABILITY

Autoclavable handpiece holder with 2 positions

- Protecting from accidental fall:
- On\* table
- On\*\* table

## ABSOLUTE CONTROL

Weighted multifunctional footswitch

## STERILISATION CONSCIENTIOUSNESS 2.2.

- Easy maintenance
- Direct access to all parts
- Meeting hygiene constraints

## 7.7 Screen of the medical device


The screen of the medical device serves as the interface display.

It is tactile such that practitioners can make adjustments by pressing the active areas.

Always make adjustments with your finger. Never use stylets or instruments, as these may damage the screen.

The touch zones are capacitive and extremely sensitive. The screen must therefore be constantly clean and dry to avoid disturbing user-defined settings.

The tactile areas, those with which you can interact with the medical device, are as shown below.

Set the irrigation flow by pressing the tactile zones + and -

Select the requisite mode by pressing the tactile zones ↑ and ↓

Activate the purge by pressing the tactile zone ⇒

## 7.8 Setting the power

The ultrasound power must be adjusted in accordance with the tip used and the required treatment.

Select the requisite mode by pressing the tactile zones ↑ and ↓

Each tip must be used in accordance with the settings defined in the power settings table for intraoral surgery ultrasonic generators [J58010].

## 7.9 Setting the irrigation flow

The medical device must be set to minimum power to adjust the irrigation flow rate. Press the footswitch until a spray appears.

Because work habits, feedback and professional training differ from one professional to another, the user must ensure that the irrigation flow is compatible with the procedure to be carried out to prevent burns to the clinical site.

Adjust the irrigation flow using the irrigation flow configuration arrows. This setting depends on the tip used and on the procedure to be carried out.

Set the irrigation flow by pressing the tactile zones + and -

### 7.9.1 Activating the purge function/starting irrigation

Press and hold the purge icon for as long as necessary.

Activate the purge by pressing the tactile zone ⇒

## 7.10 Air inlets

Air inlets ensure correct ventilation of the control unit. Leave them uncovered to allow air to circulate.

## 7.11 Control pedal

The pedal can be set either in ON/OFF mode or in progressive operation. 7.

Pressing the foot-switch automatically activates the handpiece ultrasounds and the irrigation function. **7**

The control footswitch fitted with its cord must be disconnected for daily cleaning using an alcohol disinfectant wipe.

The light function remains active for approx. 9 seconds after the pedal is released.

For further information, refer to the chapter entitled *Pedal Overview page 1*.

## 7.12 Mains Connector

The mains connector with its earthing pin is used to connect the device to the electrical network via a disconnectable mains cord.

## 7.13 Switch

The mains switch is used to switch on (position I) or to stop (position O) the medical device.

## 7.14 Fuse recess

The recess holds two fuses designed to protect the medical device in the event of overvoltage or an internal fault.

Please read the instructions listed in the chapter *Replacing the fuses page 25*

## 7.15 Irrigation lines

The autoclavable irrigation lines can be reused after cleaning, disinfection and sterilisation.

After use, sterile irrigation lines should be discarded in a biomedical waste container.

Bottles or irrigation bags should not weigh more than one kilogram. A heavier container will make the medical device tip over.

The medical device is not designed to deliver medicinal substances.

The medical device may only be used with bags or bottles of saline solution or sterile water.

# SECURE & BROADEN YOUR SCOPE IN SURGERY

6.

## Extraction

**Syndesmotomy / Osteotomy**  
For maximum bone preservation  
LC1, LC2 (K2),  
LC2L, LC2R, Ninja™ tips  
Ref. R97546

*Widening the ligamentary space. LC tips separate the tooth from its attachment system for a rapid extraction, with less luxation, keeping intact the alveolar bone edge.*

R97540: EXTRACTION PACK (A: Extraction II kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## Bone Surgery

**Osteotomy / Osteoplasty**  
Clean and thin cut for maximal bone volume  
BS1, BS2L, BS2R, BS4,  
BS5, BS6 tips  
Ref. R97509

*Performing bone graft, allows cutting, existing and remodeling bone structures reducing the risk of soft tissue lesions.*

R97500: BONE SURGERY PACK (A: Bone Surgery II kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## Crest Splitting

Rapid and minimally invasive technique for controlled expansion  
CS1, CS2, CS3, CS4, CS5,  
CS6 tips  
Ref. R97567

*The increasing tips thickness allows a smooth crest expansion limiting the risk of bone fracture.*

R97560: CREST SPLITTING PACK (A: CS kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)


## Sinus Lift

**Lateral Sinus Lift**  
Unbeatable comfort: selective and hemostatic cut  
SL1, SL2, SL3, SL4,  
SL5 tips  
Ref. R97519

*A thin cut, performed without effort, limiting any risk of sinus membrane perforation or carotid artery damage.*

R97510: SINUS LIFT PACK (A: SL11 kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## IntraLift

**Crestal Sinus Lift**  
Minimally invasive surgery for smooth sinus floor fracture  
TKW1, TKW2, TKW3,  
TKW4, TKW5 tips  
Ref. R97536

*Sinus lift by the crestal approach: the IntraLift® kit makes it possible to perform non-invasive surgery in full safety.*

R97530: IntraLift® Pack (A: IntraLift® kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## Piezocision

**Accelerated Orthodontic Surgery**  
Minimal incisions, minimally invasive  
PZ1, PZ2L, PZ2R,  
PZ3 tips  
Ref. R97576

*Treatment of malocclusions by ultrasonic corticotomies for a bone regeneration surgery in full safety. Bone healing is stimulated by micro-perforations. Tooth movement, the treatment is 3 to 4 times faster\*.*

R97570: Piezocision® Pack (A: Piezocision® kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## Crown Extension

**Crown Lengthening**  
Comparable precision and accessibility  
BS6, CE1, CE2, CE3 tips  
Ref. R97554

*Completely innocuous to the gingiva guaranteeing the integrity of bone and adjacent teeth.*

R97550: CROWN EXTENSION PACK (A: CE kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

## Essential


**The 6 fundamental ACTEON® tips**  
BS1S, BS4, SL1, SL2, SL3,  
LC2 tips  
Ref. R97528  
Kit delivered with the following references: F50510, F50100

*For osteotomy, osteoplasty, extraction and lateral sinus lift.*

R97520: ESSENTIAL Pack (A: Essential kit, a sterilization box, a Piezotome® LED handpiece and its autoclavable cord, a dynamometric wrench)

14 \*Kobayashi B, Saruhashi J, Bergman CD, Dohert S. Acceleration of orthodontic tooth movement following selective alveolar decortications: biological rationale and outcome of an innovative tissue engineering technique. *International Orthodontics*. 2008;36:232-245.


### 5.5 Attaching a tip 6.2.

The tips that can be used on this medical device are referred to as "second generation". They are recognizable by the engraving II on the base. They are incompatible with first-generation Implant Center and Piezotome tips. Conversely, older-generation tips are incompatible with Piezotome Cube, Implant Center Cube, Implant Center 2 LED, Piezotome 2 and Piezotome Solo LED.

A tip vibrates correctly when it is perfectly tightened without being forced beyond its stop point. Tighten it moderately using the wrench supplied to ensure optimum ultrasonic function. Over-tightening of the tip can result in breakage of the tip or handpiece.

| To prevent self-locking of the tip, the latter must be removed and sterilised after each use. 6.3.

The wrench is a sliding torque wrench. After a few turns, the wrench will seem to slide or turn without resistance. This means the tightening torque is reached. **6.2.**

## 4 Manual instructions

### 4.1 Pre-disinfection and cleaning – Manual method

Equipment: soft brush, soft lint-free swab, lint-free cloth, syringe, pipette or water spray, alkaline cleaner, ultrasonic cleaner.

Minimum duration of step	Cleaning instructions
1 minute	Rinse the soiled device under cold running water. Use a soft-bristled brush, a swab or a lint-free cloth to remove most of the contamination. Use a syringe, a pipette or a water spray with an alkaline or enzymatic cleaning solution to rinse the tip cannulation.
10 minutes	Immerse the medical device in a freshly prepared alkaline cleaning solution in an ultrasonic cleaner for at least ten minutes. Adhere to the manufacturer's exposure time, concentration, water quality and temperature recommendations.
1 minute	Rinse the device under cold running water. Use a syringe, a pipette or a water spray to rinse the tip cannulation.
4 minutes	Wash and disinfect the medical device using an alkaline cleaning solution. Remove surface contamination using a soft brush or a swab. Wash the medical device in water to prevent contaminants from spreading into the air
1 minute 30 seconds	Rinse the medical device in deionised or purified water. Use a syringe, a pipette or a water spray to rinse the tip cannulation.
	Visually inspect the medical device. Repeat this procedure until the medical device is visibly clean. Perform a final rinse of the device using distilled or purified water. Dry the medical device using a soft lint-free cloth or medical grade clean compressed air

### 4.2 Sterilisation 6.3.

Unless otherwise specified, non-sterile products can be resterilised using validated steam sterilisation methods (ISO 17665 or national standards). SATELEC, a company of Acteon group recommends the following:

Sterilisation exposure time	Sterilisation exposure temperature	Drying time
4 minutes	132 °C	15 minutes minimum and 20 minutes
18 minutes	134 °C	15 minutes minimum and 20 minutes
4 minutes	134 °C	15 minutes minimum and 20 minutes
3 minutes	134 °C	15 minutes minimum and 20 minutes

Saturated steam sterilisation with pre-vacuum

The drying times vary from 15 to 60 minutes according to the following criteria:

- the type of packaging material, such as a sterile barrier system or rigid reusable containers;
- steam quality;
- device materials;
- total mass;
- steriliser performance;
- usual practices for the geographical area;
- varying cool-down times.

The manufacturer accepts no responsibility for sterilisation procedures performed by the end user or the customer that are not performed according to the manufacturer's recommendations.

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ANNEX II excluding section 4 Directive 93/42/EEC concerning medical devices

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For class III devices, a EC design certificate is required

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*Medical devices and electromedical devices for dental art, electrosurgery and bone and bone substitute surgery equipments, with associated sterile and non-sterile accessories*

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GMED certifies that, on the basis of the results contained in the file referenced P159650 / P601259, the quality system - for design, manufacturing, and final inspection - of medical devices listed here above complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4

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The validity of the certificate is subject to periodic or unexpected verification

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GMED\_159650-F-10-07-2018

GMED - 21427 rev. 11  
Modifie le certificat 21427-10

On behalf of the President  
**Béatrice LYS**  
Technical Director

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**GMED 0459**

  
 DocuSigned by:  
  
 On behalf of the President  
**Béatrice LYS**  
 Technical Director

**Addendum au certificat n° 21427 rev. 11** page 2/5  
*Addendum of the certificate n° 21427 rev. 11*  
 Dossiers / Files N° P159650 / P601259

Désignation du dispositif / Accessoires marqués CE <i>Device designation / CE marked accessories</i>	Réf commerciale du dispositif ou code article <i>Device commercial reference or article code</i>	Code NBOG <i>NBOG code</i>	Classe du DM <i>MD class</i>
Pièces à main ultrasoniques de chirurgie intra-orale stériles et non stériles -- <i>Sterile and non-sterile ultrasonic handpieces for intraoral surgery</i>			
PIEZOTOME LED (Handpiece) -- PIEZOTOME SOLO (Handpiece) -- CUBE LED (Handpiece)	TE10	MD1106	Ila
Pièces à main ultrasoniques destinées à la chirurgie des os et substituts osseux stériles et non stériles -- <i>Sterile and non-sterile ultrasonic handpieces for bone and bone substitute surgery</i>			
PIEZOTOME M+ LED (Handpiece) -- PIEZOELECTRIC SYSTEM (Handpiece)	TO42	MD 1104	Ilb
Pièces à main d'aéropolissage dentaire stériles et non stériles et leurs accessoires stériles et non stériles -- <i>Sterile and non-sterile Dental air polishing handpiece with associated sterile and non-sterile accessories</i>			
AIR-N-GO EASY	TC07	MD 1106	Ila
PROPHYPEN	TC22	MD 1106	Ila
Inserts dentaires stériles et non stériles -- <i>Sterile and non-sterile dental tips</i>			
1 -- 1S -- 2 -- 3 -- 10P -- 10X -- 10Z -- P14D -- P15LD -- P15RD -- S12-70D -- AS3D -- AS6D -- AS9D -- ASLD -- ASRD -- 5AE -- PM1 -- PM2 -- PM3 -- PM4 -- C20 -- EX1 -- EX2 -- EX3 -- EXL -- EXR -- PMS1 -- PMS2 -- PMS3 -- PMV1 -- PMV2 -- PMV3 -- PMV4 -- PMV5 -- PMV6 -- H1 -- H2L -- H2R -- H3 -- H4L -- H4R -- PH1 -- PH2L -- PH2R -- P2L -- P2R -- PFU -- PFL -- PFR -- TK1-1S -- TK1-1L -- TK2-1L -- TK2-1R -- IP1 -- IP2L -- IP2R -- IP3L -- IP3R -- ET18D -- ET20 -- ET20D -- ET25 -- ET25S -- ET25L -- ET40 -- ET40D -- ETBD -- ETPR -- SO4 -- CAP1 -- CAP2 -- CAP3 -- K10/21 -- K15/21 -- K25/21 -- K30/21 -- K10/25 -- K15/25 -- K25/25 -- K30/25 -- IRR 25/21 -- IRR 25/25 -- IRR 20/21 -- IRR 20/25 -- CAVI 1-BD -- CAVI 2-DC -- CAVI 3-DF -- REDO 1 -- REDO 2 -- REDO 5 -- MAXI MPR	TG30	MD 0401	Ila
Inserts de chirurgie dentaire (intra-orale) stériles et non stériles -- <i>Sterile and non-sterile tips for intraoral surgery</i>			
BS2L II -- BS2R II -- BS4 II -- BS5 II -- BS6 II -- LC1 II -- SL1 II -- SL2 II -- SL3 II -- SL4 II -- SL5 II -- BS1S II -- BS1L II -- TKW1 II -- TKW2 II -- TKW3 II -- TKW4 II -- TKW5 II -- LC2 II -- LC2L II -- LC2R II -- NINJA II -- CE1 II -- CE2 II -- CE3 II -- BS1RD II -- CS1 II -- CS2 II -- CS3 II -- CS4 II -- CS5 II -- CS6 II -- PZ1 II -- PZ2R II -- PZ2L II -- PZ3 II	TE25	MD 0401	Ila
Inserts de chirurgie des os et substituts osseux non stériles -- <i>Non-sterile tips for bone and bone substitute surgery</i>			
Saw 20.9 x 14.1 x 4.0 x 0.6mm for Piezoelectric System (03.000.401) -- Saw 20.1 x 21.4 x 4.0 x 0.6mm for Piezoelectric System (03.000.402) -- Saw Left 17.8 x 15.2 x 4.0 x 0.6mm for Piezoelectric System (03.000.403) -- Saw Right 17.8 x 15.2 x 4.0 x 0.6mm for Piezoelectric System (03.000.404) -- Scalpel round 22.45 x 12.6 x 3.9 dia x 0.7mm for Piezoelectric System (03.000.405) -- Scalpel Flat 21.6 x 11.1 x 3.0 x 0.5mm for Piezoelectric System (03.000.406) -- Scalpel 17.7 x 15.5 x 1.2 x 0.5mm for	TO49	MD 1104	Ilb

**GMED 0459**

On behalf of the President  
**Béatrice LYS**  
 Technical Director

**Addendum au certificat n° 21427 rev. 11** page 3/5  
**Addendum of the certificate n° 21427 rev. 11**  
 Dossiers / Files N° P159650 / P601259

Désignation du dispositif / Accessoires marqués CE Device designation / CE marked accessories	Réf commerciale du dispositif ou code article Device commercial reference or article code	Code NBOG NBOG code	Classe du DM MD class
Piezoelectric System (03.000.407) -- Spade round 24 x 11.05 x 4.0 x 0.65mm for Piezoelectric System (03.000.408) -- Lift 22.9 x 10.1 x 5.0 mm dia x 0.4mm for Piezoelectric System (03.000.411)			
BS1 XXL M+ -- BS1L M+ -- BS1RD M+ -- BS2L XL M+ -- BS2R XL M+ -- BS4 M+ -- BS6 XXL M+ -- RHL5 M+ -- RHS1 M+ -- RHS2Fb M+ -- RHS2Hb M+ -- RHS3L M+ -- RHS3R M+ -- RHS4L M+ -- RHS4R M+ -- RHS5 M+ -- RHS6 M+ -- RHS2F M+ -- RHS2H M+ -- SL1 M+ -- SL2 M+ -- SL3 M+ -- SL4 M+ -- SL5 M+ -- BS1S M+ -- BS2R M+ -- BS2L M+ -- BS5 M+ -- BS6 M+ -- LC1 M+ -- LC1 XL M+ -- LC2 M+ -- LC2 L M+ -- LC2 R M+ -- NINJA M+ -- CE1 M+ -- CE2 M+ -- CE3 M+ -- TKW1 M+ -- TKW2 M+ -- TKW3 M+ -- TKW4 M+ -- TKW5 M+ -- CS1 M+ -- CS2 M+ -- CS3 M+ -- CS4 M+ -- CS5 M+ -- CS6 M+ -- DB1 M+ -- DB2 M+	TO51	MD 1104	IIB
Inserts de chirurgie des os et substituts osseux stériles -- Sterile tips for bone and bone substitute surgery			
Saw 20.9 x 14.1 x 4.0 x 0.6 mm for Piezoelectric System, sterile (03.000.401S) -- Saw 20.1 x 21.4 x 4.0 x 0.6 mm for Piezoelectric System, sterile (03.000.402S) -- Saw left, 17.8 x 15.2 x 4.0 x 0.6 mm for Piezoelectric System, sterile (03.000.403S) -- Saw, right, 17.8 x 15.2 x 4.0 x 0.6 mm for Piezoelectric System, sterile (03.000.404S) -- Scapel round, 22.45 x 12.6 x 3.9 mm dia x 0.7 mm for Piezoelectric System, sterile (03.000.405S) -- Scapel flat, 21.6 x 11.1 x 3.0 x 0.5 mm for Piezoelectric System, sterile (03.000.406S) -- Scapel 17.7 x 15.5 x 1.2 x 0.5 mm for Piezoelectric System, sterile (03.000.407S) -- Spade round, 24 x 11.05 x 4.0 x 0.65 mm for Piezoelectric System, sterile (03.000.408S) -- Diamond, 24.6 x 12.85 x 2.6 x 0.6 mm for Piezoelectric System, sterile (03.000.409S) -- Diamond round, 21.9 x 12.4 x 1.8 mm for Piezoelectric System, sterile (03.000.410S) -- Sinus Lift, 22.9 x 10.1 x 5.0 mm dia x 0.4 mm for Piezoelectric System, sterile (03.000.411S) -- Saw 105.7 x 4.0 x 0.6mm for Piezoelectric System, sterile (03.000.412S) -- Saw left 39.5 x 4.0 x 0.8mm for Piezoelectric System, sterile (03.000.418S) -- Saw right 39.5 x 4.0 x 0.8mm for Piezoelectric System, sterile (03.000.419S) -- Scalpel flat 105.1 x 3.0mm for Piezoelectric System, sterile (03.000.421S) -- Scalpel 32.1 x 1.2mm for Piezoelectric System, sterile (03.000.423S) -- Saw round 38.9 x 5.0mm dia x 0.6mm for Piezoelectric System, sterile (03.000.424S) -- Diamond cylindrical, 30.7 x 4.7 x 2.3mm for Piezoelectric System, sterile (03.000.425S) -- Diamond conical 26.8 x 11.0 x 0.5mm for Piezoelectric System, sterile (03.000.426S)	TO50	MD 1104	IIB

**GMED 0459**

On behalf of the President  
 Beatrice LYS  
 Technical Director

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**GMED 0459**

  
 DocuSigned by:  
  
 On behalf of the President  
 Beatrice LYS  
 Technical Director

Addendum au certificat n° 21427 rev. 11 page 5/5  
 Addendum of the certificate n° 21427 rev. 11  
 Dossiers / Files N° P159650 / P601259

Désignation du dispositif / Accessoires marqués CE <i>Device designation / CE marked accessories</i>	Réf commerciale du dispositif ou code article <i>Device commercial reference or article code</i>	Code NBOG <i>NBOG code</i>	Classe du DM <i>MD class</i>
Nécessaires d'Inserts de chirurgie des os et substituts osseux non stériles – <i>Non-sterile tips for bone and bone substitute surgery Procedure Packs</i>			
KIT CMF NON STERILE (F57803) -- KIT RHS BASIC M+ N STER (F87681) -- KIT RHS EXPERT M+ N STER (F87689) -- KIT PRE-IMPLANTAIRE M+ (F57805) -- KIT SINUS LIFT M+ (F57806) -- KIT INTRALIFT M+ (F57807) -- KIT CREST SPLITTING M+ (F57808) -- KIT EXTRACTION M+ (F57809) -- KIT M+ CROWN EXTENSION N STER (F57810)	TO51	MD 1104	IIb
PIEZOTOME M+ PACK (F57802)	TO53	MD 1104	IIb
Nécessaire d'Inserts de chirurgie des os et substituts osseux stériles – <i>Sterile tips for bone and bone substitute surgery Procedure Packs</i>			
KIT CMF STERILE (F57804) -- RHINOPLASTY ESSENTIAL KIT STERILE (F87999) -- RHINOPLASTY EXPERT KIT STERILE (F88000)	TO52	MD 1104	IIb

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