Popieriaus-plastiko sterilizavimo rulonai STERIDIAMOND

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Gamintojas: E.C.S. Srl (Italija)

TECHNICAL DATA SHEET
STERIDIAMOND - STERIPERFECT
ROLLS

STERIDIAMOND - STERIPERFECT rolls have been designed to be used in the hospitals and dental centres by anyone who needs to have the instruments or other medical devices sterilized by STEAM, EO and/or FORM.

STERIDIAMOND - STERIPERFECT rolls made of moisture-resistant medical 70g/m² paper manufactured with 100% pure pulp and plastic transparent layer made of polyester and polypropylene. The thickness of plastic part is: 12µm polyester, 40µm polypropylene and with the resin, it arrive on a total of 55µm. The medical paper and film has a tension strength of 350kPa. Strength of seals min 140 N/m.

Description:

Sterilization flat and/or gusseted rolls. Preformed Sterile Barrier System according to EN ISO 11607-1

REF No's:

STERIDIAMOND	STERIDIAMOND WITH FO
SRPXXXXX (flat reels)	SRPFXXXX (flat reels)
SRSXXXXX (gusseted reels)	SRSFXXXX (gusseted reels)

Conformity:

CE mark according MDD 93/42/EEC (Class I), 2003/32/CE, EN 868-3, EN 868-4, EN 868-5, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1:2012, EN 1041:2013, EN ISO 14971:2012, process indicators according ISO11140-1, ISO15882:2009

Recommended sealing parameters:

160 – 180°C (best performance in 170 °C)

Storage:

Products should be stored in dry areas away from light and heat. Temperature range 10°C - 30°C. Relative humidity: 30-60 %

Packaging:

Rolls are sealed in a transparent film in order to ensure their protection against dust.

Shelf life:

The product has a 5 years shelf-life from production date.

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STERIDIAMOND® - STERIPERFECT® POUCHES AND ROLLS

PRODUCT DESCRIPTION

- Medical devices belonging to class 1.
- Disposable devices designed for STEAM, EO gas and FORM sterilization.
- Not suitable neither for dry heat sterilization at temperatures higher than 140°C nor gas plasma sterilization.
- Medical paper (60g/mq) and plastic film made of polyester and polypropylene.
- The paper employed is white colour and resistant to humidity. It has been developed in order to grant a high level of protection against bacteria.
- All the printing indications, including the chemical indicators STEAM, EO gas and FORM sterilization, are positioned on the paper side in order to prevent the contamination of the medical devices with the ink and avoiding any release of potentially polluting substances on the product itself (in compliance with EN ISO 11607 and EN 868 standards).
- The product is designed to be welded together.
- The multi-line, hot welding is continuous with rounded corners in order to prevent any deposit of dirt, thus reducing the risk of contamination of the content during stocking and upon package opening.
- STERIDIAMOND® - STERIPERFECT® pouches and rolls are in compliance to European Directive 93/42/CEE.

INSTRUCTIONS FOR USE

- Never use in case the packaging is damaged
- Evaluate the dimension of the product to be sterilized and choose the proper size
- Insert the product to be sterilized, securing sharp parts. Put instrument up to 75% of all spece to allow air flowing.
- Remove most air from the packaging before seeling.
- Ensure that the packaging is well sealed
- Start sterilization process.
- The steam process indicator turns from pink to brown. The EO process indicator turns from green/blue to brown.
- Indications of the process indicator are described in several languages and explained through international symbols
- Minimal release of fiber during opening thanks to the features of the medical paper.

CONDITIONS OF STORAGE

- Avoid direct exposure to sunlight or heat sources and store in environments avoiding any contamination.
- Keep far from heat, light and humidity. Rotate your stock accordingly to arrival order and handle with care.
- Shelf-life: 5 years from date of production
- Batch and production number for a complete traceability.

HEALTH AND SAFETY

When STERIDIAMOND® - STERIPERFECT® pouches and rolls are correctly used, there are no special precautions to adopt during their handling and use.

These products do not contain toxic or irritating agents or material.

REV.4 2015-04

This document is property of E.C.S. s.r.l. and the information herein contained are for industrial purposes only. Any unauthorized copy scanning, different use or publication of this document without prior consent is strictly prohibited. Any improper use will be prosecuted by law

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DECLARATION OF CE CONFORMITY

In compliance with Annex VII, Directive 93/42 CEE and further modifications and integrations, such as those provided to DL n. 37 of 25/01/2010, acknowledging the European Directive 2007/47/CE.

ECS[®] SRL
Via I° Maggio, 18/22
23881 Airuno (LC)
Italy

Declares under its own responsibility that the products, designed as:

"STERIDIAMOND[®] - STERIPERFECT[®] FLAT ROLLS, GUSSETED ROLLS, FLAT POUCHES, GUSSETED POUCHES, SELF-SEALING POUCHES for STEAM, EO gas and FORM sterilization" satisfy all the requirements applicable to European Directive 93/42/CEE and further modifications such as those provided by the D.L. n. 37 of 25/01/2010, acknowledging the European Directive 2007/47/CE, concerning the medical devices and are registered to the Italian Ministry of Health.

For this purpose, E.C.S. s.r.l. grants and declares under its own responsibility that:

- The range of the products referred to subject satisfies the essential requirements as per Annex I of the European Directive 93/42/CEE and further modifications such as those provided by the D.L. n. 37 of 25/01/2010, acknowledging the European Directive 2007/47/CE;
- The range of the products referred to subject belongs to Class I - non sterile;
- The devices referred to subject are not measurement instruments;
- The devices referred to subject are not designed for clinical trials purposes;

- The devices referred to subject comply to the European standards provided by UNI EN ISO 11607-1/-2 and UNI EN 868 - Packaging for medical devices terminally sterilized - part 3: paper to be used in the manufacturing of paper pouches (specified in the EN 868-4) and in the manufacturing of pouches and rolls (specified in the EN 868-5) - requirements and test methods
- ECS srl, manufacturer of these products, is certified according to the rules UNI EN ISO 9001:2008 and UNI EN ISO 13485:2012.

E.C.S. s.r.l. declares that the documentation according to the Annex VII of the European Directive 93/42/CEE and further modifications such as those provided by the D.L. n. 37 of 25/01/2010, acknowledging the European Directive 2007/47/CE, will be held at disposal during a period of five years starting from the last manufacturing batch of the devices referred to subject.

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ECS Srl