HANDLING THE PRODUCT

Once the sterilizer has been taken out of the box, it must be lifted by two persons using appropriate means and possibly handled with a lift truck or similar.

**WARNING**

WE RECOMMEND THAT THE DEVICE BE TRANSPORTED AND STORED AT A TEMPERATURE NO LOWER THAN 5 °C. PROLONGED EXPOSURE TO LOW TEMPERATURE CAN DAMAGE THE PRODUCT.

**NOTE**

KEEP THE ORIGINAL PACKAGING AND USE IT WHENEVER THE DEVICE IS TO BE TRANSPORTED. THE USE OF DIFFERENT PACKAGING COULD DAMAGE THE PRODUCT DURING SHIPMENT.

**DANGER**

BEFORE TRANSPORT, LEAVE THE DEVICE TURNED-OFF FOR ABOUT 30 MINUTES AFTER THE LAST PROGRAM FINISHES AND DRAIN THE DISTILLED WATER AND USED WATER TANKS SO THAT ALL THE HOT INTERNAL PARTS WILL HAVE TIME TO COOL.

1.

The **Classic** series sterilisers are the revolutionary products offered by MOCOM in the field of small water steam sterilisers, equipped with type B (EN 13060) cycles, as well as the new point of reference with respect to safety, performance, flexibility and ease of use.

**PRODUCT INTRO-
DUCTION**
INTRODUCTION

It is a sophisticated but, at the same time, easy to use device that, thanks to its wide range of configuration options and patented operating devices, satisfies every need for sterilizing medical devices, guaranteeing the maximum performance under all conditions.

It also features a better way of relating to users who, rather than having to adapt to the machine and its characteristics, are able to "converse" with it and configure it to meet their own needs.

Thanks to its remarkable ease of use, small size and pleasant appearance, it is the ideal partner for all professional who demand the maximum sterilization safety.

2.
4.

A **Classic** series steriliser is an electronic water steam steriliser that is entirely operated by a micro-processor with a large, printed stainless steel sterilisation chamber.

It is characterized by an advanced fractionated vacuum system for the complete removal of air, even from hollow, porous materials, and an effective final vacuum drying phase capable of eliminating all traces of humidity from any load.

**GENERAL CHA-
RACTERISTICS**

2.

Its exclusive steam generation system, effective plumbing circuit and electronic management (supplemented by high-precision sensors) guarantees high process execution speeds and excellent thermodynamic parameter stability.

Moreover, its Process Evaluation System constantly monitors all the machine's "vital" parameters in real-time, guaranteeing absolute safety and a perfect result.

It offers users 6 sterilization programs (of which one completely PRESETble), all equipped with customizable, optimized drying for the fast, effective sterilization of the various types of loads (instruments and materials) used in a medical environment.

All the cycles can immediately be selected on the clear LCD screen, which also allows extensive configuration of the device according to the user's needs.

Please refer to the chapter, "**Configuration**" for more detail.

Like in the best MOCOM tradition, also the new range of **Classic** autoclaves has the most complete, sophisticated and advanced safety systems available today to guarantee the user against any electrical, mechanical, thermal or biological malfunction.

NOTE	
	FOR THE DESCRIPTION OF THE SAFETY DEVICES, REFER TO THE APPENDIX TECHNICAL CHARACTERISTICS.

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APPENDIX – TECHNICAL CHARACTERISTICS

SUMMARY TABLE

Device	Steam Sterilizer			
Classification (according to the Directive 93/42/EEC and subsequent changes)	B Classic 17	B Classic 22	B Classic 28	
Model	II b			
Manufacturer	Mocom s.r.l. Via Saliceto, 22 40013 Castel Maggiore (BO) - ITALIA			
Power supply voltage	220 V - 240 V~ 50 Hz 220 V - 230 V~ 50 Hz			22.
Mains fuses (6.3 x 32 mm)	F 15A 250V			
On-board fuses (5 x 20 mm)	F1: T3.15A 250V (primary transformer) F2: T 3.15A 250V (vacuum pump)			
External dimensions (HxWxD) (excluding rear connections)	500 x 495 x 610 mm			
Nominal power	2300 W	2300 W	2500 W	16.
Insulation class	Class I			
Installation category	Cat. II			
Environment of use	Internal use			
Sound power level (A weighted)	< 65 db(A)			
Environmental operating conditions	Temperature: +15°C + +35°C Relative humidity: max 80%, non-condensing Altitude: max 3000 m (a.s.l.)			
Net weight: empty empty with trays and supporti empty, with trays and supports and water at MAX level	about 47 kg about 49 kg about 53 kg	about 55 kg about 57 kg about 61 kg	about 58 kg about 60 kg about 64 kg	
Sterilization chamber dimensions (Ø x D)	250 x 350 mm	250 x 450 mm	280 x 450 mm	
Sterilization chamber total volume	about 17 l (0.017 m³)	about 22 l (0.022 m³)	about 28 l (0.028 m³)	3.
Sterilization chamber useful volume (with tray supports inserted)	about 10 l (0.010 m³)	about 13 l (0.013 m³)	about 19 l (0.019 m³)	
Distilled water tank capacity (supply)	about 5,5 l (water at MAX level) about 1 l (water at MIN level)			11., 12.
Sterilization programs	5 standard programs + 1 user-defined program			8.
Test programs	Helix/BD Test Vacuum Test Vacuum Test+Helix/BD Test			8.
Preheating time (from cold)	about 10 min			
USB connection	requires memories formatted with the FAT32 file system			13.
Printer connection	RS232 serial (printer cable length max 2.5 m)			7.
Bacteriological filter (PTFE filtering element)	Porosity: 0,2 mm Connection: male 1/8" NPT connector			

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SAFETY DEVICES

The sterilizer is equipped with the following safety devices for which we provide a brief description of their function:

- **Mains fuses** (see summary table data)
Protection inside the device against a fault in the heating elements.
Action: cuts the electricity.
- **Fuses protecting the electronic circuits** (see summary table data)
Protection against a fault in the primary transformer circuit and low voltage uses.
Action: cuts power to one or more low-voltage circuits.
- **Thermal circuit breakers on the mains voltage windings**
Protection against overheating of the vacuum pump motor and the primary transformer windings.
Action: temporary cut-off (until cooling) of the winding.
- **Safety valve**
Protection against overpressure in the sterilization chamber.
Action: release of the steam and restoration of the safety pressure.
- **Steam generator manual rearm safety thermostat**
Protection against steam generator overheating.
Action: cut-off of the electricity to the steam generator.
- **Heating element manual rearm safety thermostat**
Protection against overheating of the heating elements of the container under pressure.
Action: cut-off of the electricity to the chamber heating element.
- **Door position safety microswitch**
Confirmation of the correct closing position of the door of the container under pressure.
Action: signals wrong door position.
- 6. – **Mechanized door lock mechanism with electromechanical protection (pressure switch)**
Protection against accidental opening of the door (even in a blackout).
Action: prevents accidental opening of the door during a program.
- **Door lock mechanism safety microswitch**
Confirmation of the correct closing of the door lock.
Action: signaling the failure or incorrect operation of the door lock mechanism.
- **Self-leveling plumbing system**
Plumbing system structure for the spontaneous leveling of the pressure in the case of a manual interruption of the cycle, alarm or blackout.
Action: automatic restoration of atmospheric pressure in the sterilization chamber.
- 2. – **Integrated system for evaluating the sterilization process**
Continuous verification of the sterilization process parameters entirely managed by microprocessor.
Action: immediate interruption of the program (in case of anomaly) and generation of alarms.
- **Monitoring of the sterilizer's operation**
Real-time oversight of all significant parameters when the machine is powered.
Action: generation of alarm messages (in the case of anomaly) with possible interruption of the cycle.

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PROGRAM SUMMARY TABLE – B CLASSIC 17

DESCRIPTION CYCLE	NOMINAL VALUES				BASIC CYCLE PARAMETERS					STERILIZABLE MATERIALS			NOTES	
	Temperature (°C)	Pressure (bar)	Retention time (min)	Cycle time (EN 13060)	Pre-vacuum (F=fractionated; S=Single)	Standard drying (L=long; C=short)	Total cycle time (average load + max load)	Average H ₂ O consumption (ml/cycle)	Average energy consumption (kWh/cycle)	TYPE	MAX TOTAL WEIGHT (kg)	MAX WEIGHT FOR TRAY (kg)		MAX WEIGHT FOR ITEM (kg)
134°C UNIVERSAL	134	2.10	4	B	F	L	37+40	500	0.8	Unwrapped porous materials	1.00	0.30	0.30	
										Porous materials in single package	0.75	0.25	0.25	
										Porous materials in double package	0.60	0.20	0.20	
										Solid and hollow materials in single package	3.00	1.00	0.25	
										Unwrapped solid and hollow materials	6.00	1.20	0.50	
										Solid and hollow instruments A" in double package	1.50	0.50	0.25	
134°C PRION	134	2.10	>18	B	F	L	53+56	550	0.9	Unwrapped porous materials	1.00	0.30	0.30	
										Porous materials in single package	0.75	0.25	0.25	
										Porous materials in double package	0.60	0.20	0.20	
										Hollow instruments in single package	3.00	1.00	0.25	
										Unwrapped solid and hollow materials	6.00	1.20	0.50	
										Solid and hollow instruments in double package	1.50	0.50	0.25	
121°C UNIVERSAL	121	1.10	20	B	F	L	52+55	550	0.8	Unwrapped porous materials	1.00	0.30	0.30	
										Porous materials in single package	0.75	0.25	0.25	
										Porous materials in double package	0.60	0.20	0.20	
										Unwrapped solid and hollow materials	6.00	1.20	0.50	
										Hollow instruments in single package	3.00	1.00	0.25	
										Solid and hollow instruments in double package	1.50	0.50	0.25	
134°C HOLLOW UNWRAPPED	134	2.10	4	S	F	C	32+35	500	0.7	Unwrapped hollow instruments	6.00	1.20	0.50	
										Unwrapped solid and hollow materials	6.00	1.20	0.50	
134°C SOLID WRAPPED	134	2.10	4	S	S	L	33+36	300	0.6	Solid and hollow instruments "B" in single package	3.00	1.00	0.25	
										Unwrapped solid and hollow materials	6.00	1.20	0.50	
XXX°C USER (see note)	134 or 121	2.10 or 1.10	>4 or >20	n.a.	F/S	L/C	n.a.	n.a.	n.a.	Unwrapped solid instruments (other load types are possible depending on the user settings)	n.a.	n.a.	n.a.	Variable parameters depending on the settings made
HELIX/BD TEST	134	2.10	3.5	-	F	C	29	-	-	Test device only (without another load)	-	-	-	
VACUUM TEST	-	-0,80	-	-	-	-	27	-	-	Empty chamber	-	-	-	
VACUUM + HELIX/BD TEST (executable in sequence)	-	-	-	-	-	-	64	-	-	-	-	-	-	

5.

For wrapped materials and instruments (single and double package), it is advisable to use the 3-tray configuration

It is advisable to use the 3-tray configuration

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

5.

10.5. DOOR OPENING


To open the autoclave door, press and hold the button shown in the figure.

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


Accessories

A full range of accessories to expand Classic's functions

- 1 **Pure 100**
The Pure 100 resin demineraliser consists of two ion exchange resin cartridges. This system eliminates the ions that contaminate mains water, allowing high quality demineralised water to be produced. The system is fully controlled by the sterilizer, which tells the operator when the cartridges need replacing; replacement is completed in a few short, simple steps.
- 2 **Pure 500**
The Pure 500 reverse osmosis demineraliser is an ecological system that ensures, cycle after cycle, a proper supply of demineralised water to Mocorm autoclaves. Passing the mains water through an osmotic membrane eliminates both organic and inorganic substances, resulting in water that is perfect for sterilization. The autoclave has a visual/acoustic warning system that tells the operator when it is necessary to change filters, thus ensuring performance.
- 3 **Kit Twin Pure 500**
Thanks to this kit, the Pure 500 demineraliser can be used simultaneously on two sterilizers.
- 4 **External printer**
Connected to the Classic sterilizers via an RS232 serial port, it allows executed cycle data to be printed on rolls of thermal paper or bar code labels.
- 5 **Front filling kit**
Allows the operator to fill the main tank frontally through a quick coupling.
- 6 **Automatic filling kit**
Consists of an external pump powered by the sterilizer and allows demineralised water to be withdrawn from a tank or recipient.
- 7 **SV aux kit**
This kit allows the sterilizer to be interfaced with surgery demineralising systems so the sterilizer only draws water from the system, as and when necessary.
- 8 **Modular tray holder**
Modular housing system for filling the sterilizer. Provided with 3 pairs of tray-housing shelves, it allows for the full use of the chamber diameter and the introduction of even cumbersome loads.
- 9 **Bar code reader**
The bar code reader, usable on the most common PCs, makes it possible to complete the traceability process and associate bar code labels to patients.

Pure demineralisers

Technical data	Pure 100	Pure 500
Water input	Mains water	Mains water
Water output (conductivity)	< 15 µS/cm	< 15 µS/cm
Distilled water output (inlet water conductivity approx. 350 µS/cm)	120 litres	800 litres
Min. and max. mains water pressure	1-5 bar	1-5 bar
Min. and max. mains water temperature	+5 °C ÷ +40 °C	+5 °C ÷ +40 °C
Dimensions W x H x D	366 x 365 x 143 mm	461 x 328 x 206 mm
Weight	3,5 kg	5,5 kg


Sterilization cycles 220/230 V

CYCLE	Cycle type	Sterilization time (min.)	B Classic 17	B Classic 22	B Classic 28	Drying *
134 °C UNIVERSAL	B	4	29	31	39	13-17
121 °C UNIVERSAL	B	20	45	48	52	13-17
134 °C UNWRAPPED HOLLOW INSTRUMENTS	S	4	31	34	38	4-6
134 °C WRAPPED SOLID INSTRUMENTS	S	4	20	24	28	13-17
134 °C PRION	B	18	43	45	53	13-17
XXX °C CUSTOM	S	Cycle can be personalised by user at temperatures 134 °C / 121 °C, process times starting from 4' (134 °C) or 20' (121 °C) with adjustable drying				
VACUUM TEST	TEST	TEST	18	18	19	
HELIX / B&D TEST	TEST	TEST	20	24	28	
VACUUM + HELIX / B&D TEST (executed in sequence)	TEST	TEST	42	46	51	


* Drying varies according to machine model and volume. Note: times do not take the max. pre-heating time into account (10 min.). Note: times may vary depending on the load and power supply. Note: 120 V model times may undergo increments up to a max. of 20%.

Technical data	B Classic 17	B Classic 22	B Classic 28
Power supply voltage	220/240V 50Hz 220/230V 60Hz		
Rated power	2300 W		
External dimensions W x H x D	480 x 500 x 600 mm		
Chamber dimensions Diam. x Depth	250 x 350 mm	250 x 450 mm	280 x 450 mm
Total weight	50 kg	51 kg	52 kg
Tank capacity	4,5 litres		
Autonomy (with max. water level)	From 7 to 11 cycles	From 6 to 10 cycles	From 5 to 9 cycles


PT2 PROBE TO GENERAL CARD CONNECTION

- 1 PCB
- 2 PT2 wiring connection

- 1 Water inlet
- 2 Heating element
- 3 Cylindrical outer body
- 4 PT2 probe
- 5 Safety thermostat
- 6 Inner core
- 7 Steam outlet

1 Steam generator complete with heat insulation

9.

FILLING DISTILLED WATER

Manual filling


10.

The first time the sterilizer is used, and later when the MIN water level indicator comes on, you will have to fill, or top-off, the internal distilled water tank.

Open the loading door.

Pour in water taking care not to exceed the maximum level (MAX) indicated inside the tank. Close the door.

Be carefulnot to spill any water on the machine, and if so, immediately dry it off.

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







EN

Ref.	Device description
3	Door
9	Drain tank
10	Filling tank
11	USB port
12	Power switch

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

USB MANAGEMENT

Before carrying out the following operations, insert the USB key.



Selecting USB MANAGEMENT, you can choose whether to download the data or the technical report.

13.

Selecting **CYCLE DATA DOWNLOAD**, you can copy the data relating to the cycles executed (stored in the internal memory of the sterilizer) to the USB key.
 NOTE: The USB key must be formatted using the FAT32 file system.
 You can select the number of cycles to download to the external memory.

If the USB key is not present, you are asked to insert it.

Once the data has been downloaded, you can remove the USB key.

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
Castel Maggiore, 17th September 2015

TO: WHO IT MAY CONCERN

We, MOCOM S.r.l., with its principal offices in Imola, Via Selice Provinciale 23/a, ITALY – and plant in Castel Maggiore (BO), Via Saliceto 22, ITALY hereby declare that our sterilizers model B CLASSIC have menu in Lithuanian and that all information regarding the cycles or other messages are displayed in Lithuanian.

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B Classic sterilizer series

TECHNICAL SPECIFICATIONS

- 6 sterilization programs (custom cycle included)
- 3 test programs (Helix Test, Vacuum Test, Vacuum + Helix combined Test)
- Monochromatic graphic display with "soft touch" keypad
- Fractioned vacuum air removal system (standard vacuum pump)
- Final drying
- Anti-dust filter
- Water conductivity sensor integrated in the water reservoir
- Built-in Wifi
- USB port (for PC data transfer)
- RS232 serial port (for external printer connection)
- Front filling connection (for front filling kit - see Accessories list)
- Provision for automatic filling system (see Accessories list)
- Direct drain system
- Built-in "DataSter" software for automatic download of cycle reports to a PC
- Ready for EasyCheck remote service

Supplied trays:

- B Classic-17: 5 aluminium trays (288 x 188 x 17 mm) **15.**

- B Classic-22: 5 aluminium trays (384 x 188 x 17 mm)

- B Classic-28: 6 stainless steel trays (384 x 208 x 17 mm)

Sterilization chamber dimensions (ØxD):

- B Classic-17: 250 x 350 mm

- B Classic-22: 250 x 450 mm

- B Classic-28: 280 x 450 mm

External dimensions (WxHxD): 480 x 500 x 600 mm

Weight: 49 kg (B Classic-17), 54 kg (B Classic-22), 56 kg (B Classic-28)

Power supply: 220-240Vac 50Hz - 2300W *

* Other power supplies available upon request

Description	Ref.	Price €
B Classic-17 (17 litres) sterilizer	7A210000	
B Classic-22 (22 litres) sterilizer	7A211000	
B Classic-28 (28 litres) sterilizer	7A212000	

GENERAL TERMS OF SUPPLY

CURRENCY

Euro (€)

PAYMENT TERMS

As agreed

VAT

At the purchaser's charge

DELIVERY

Within 30 (thirty) working days from the order confirmation, except unforeseen circumstances and force majeure

PACKING

Included

TRANSPORT

FCA IMOLA (BO) - INCOTERMS 2010

PRICELIST VALIDITY

For orders received from 2020, June 1st

AVAILABILITY OF PRODUCTS IN THE PRICELIST

From 2020, June 1st

WARRANTY TERMS

Sterilizers 36 months or 3,000 cycles **23.**

Tethys H10: 36 months

Tethys T45/T60/D60: 24 months

Highea: 24 months

Thermal sealers: 24 months

Accessories: 12 months

NOTE: For warranty conditions please refer to the relevant Warranty Certificate

Printed in Italy - All rights reserved

CEFLA S.C. reserves the right to make technical and aesthetic changes to products and pieces described herein without further notice.

THIS PRICE LIST CANCELS AND REPLACES THE PREVIOUS ONES