MiniCollect®

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Literatūra

GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture;
Approved Standard-Sixth Edition.

GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens;
Approved Standard-Sixth Edition.

H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based
Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition.

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www.gbo.com/preanalytics
faq.pa@gbo.com
(+43) 07583/6791

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VACUETTE® Blood Collection Set Instructions for Use

Intended Use

The VACUETTE® Blood Collection Set is used in routine venipuncture procedures. The product is to be used by appropriately trained healthcare professionals only in accordance with these instructions.

Product Description

The VACUETTE® Blood Collection Set is a single-use, sterile, winged blood collection needle bonded to a flexible tubing with a female luer adapter. It is available with optional Luer Adapter or Luer Adapter+Holder.

Precautions/Cautions

Precautions

The device will perform as intended when the instructions are followed accordingly.

- Examine individual package for integrity of packaging prior to use. If packaging has been torn, do not use.
- Any used needle is considered contaminated and should be disposed of in an approved “sharps” biohazard container immediately after use.
- Keep hands behind needle at all times during use and disposal.

Caution

- Handle all biological samples and blood collection “sharps” (lancets, needles, luer adapters, and blood collection sets) according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples (for example, through a puncture injury), since they may transmit HIV (AIDS), viral hepatitis, or other infectious disease.
- Discard all blood collection “sharps” (lancets, needles, luer adapters, and blood collection sets) in biohazard containers approved for their disposal.
- Gloves should be worn at all times during venipuncture to minimize exposure hazard.

Storage of VACUETTE® Blood Collection Sets before use

Store the VACUETTE® Blood Collection Set at 4–25°C (40–77°F).

NOTE: Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the VACUETTE® Blood Collection Set quality. **NOTE:** The VACUETTE® Blood Collection Set contains DEHP plasticised polyvinylchloride.

Handling of the VACUETTE® Blood Collection Set

WEAR GLOVES DURING VENIPUNCTURE AND WHEN HANDLING BLOOD COLLECTION TUBES TO MINIMIZE EXPOSURE HAZARD.

1. Select appropriate tube(s) and VACUETTE® Blood Collection Set.
2. Remove VACUETTE® Blood Collection Set from packaging.
NOTE: If individual packaging has been opened or tampered with, please choose another VACUETTE® Blood Collection Set.
3. Select site for venipuncture. Apply tourniquet and prepare site with appropriate antiseptic. **DO NOT PALPATE** site after cleansing. Grasp wings of VACUETTE® Blood Collection Set and carefully remove needle cap from winged needle.
4. Perform venipuncture with patient’s arm in downward position and tube cap upper-most.
5. After completion of venipuncture, gently remove needle from vein, applying pressure to puncture site with dry sterile swab until bleeding stops.
NOTE: Do not apply undue pressure to the puncture site.
6. Promptly dispose of VACUETTE® Blood Collection Set in an approved disposal container in accordance with the procedures of your facility.

Label information

| | |
|--|---|
|  | Item number |
|  | Lot number: Charge number |
|  | Expiry Date: Product can be used until the end of the month indicated |
|  | Symbol for “sterilisation with ethylene oxide gas” |
|  | Symbol for “single use only” |

Reference:

Clinical and Laboratory Standards Institute (CLSI)
H1-A5 “Evacuated Tubes and Additives for Blood Specimen Collection – 5th Edition”;
Approved Standard

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Literatūra

GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition.

GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition.

H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition.

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UAB Diamedica

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Intended Use

To collect, transport, store, and evaluate capillary blood specimens for hematology tests.

Product Description

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes. The interior of the tube wall is coated with dipotassium EDTA (K2EDTA). The EDTA binds calcium ions thus blocking the coagulation cascade.

| Tube Type | Matrix | Cap colour |
|-------------------------------|-----------------|------------|
| MiniCollect® K2E K2EDTA Tubes | Capillary blood | Lavender |

Product versions

MiniCollect® Tubes with optional 13x75mm carrier tubes
MiniCollect® Tubes pre-assembled with 13x75mm carrier tubes

Storage before use

Store tubes at 4–25°C (40–77° F). *Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. evaporation of liquid additives, colouring, etc.) All tubes are able to withstand a freezing temperature of –20°C.*

Specimen Collection and Handling**Equipment required but not provided**

- Labels for positive patient identification of samples
- Gloves and appropriate apparel for protection against exposure to blood borne pathogens
- Alcohol swab
- Dry gauze
- Warming device if required, depending on the volume of blood needed and the test(s) to be performed
- Adhesive plaster or bandage
- Biohazard disposal container
- Safety lancet

Recommended Order of Draw

If multiple specimens are to be collected, including EDTA specimens, the EDTA specimen is drawn first to ensure adequate volume and accurate hematology test results. Specimens with other additives are collected next; serum specimens are collected last.

Cap Removal

The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop.

Specimen Collection

For the collection of blood by skin puncture, please refer to the instructions for use for the lancet device used. Blood flow from the skin puncture site is increased by holding the puncture site in a downward position. The first drop of blood should be wiped away with a gauze pad if appropriate for the test method, as it is likely to contain excess tissue fluids. Apply gentle, intermittent pressure to the surrounding tissue (or proximal to the puncture site when using a finger). Strong, repetitive pressure (milking), must not be applied. It may cause hemolysis or tissue-fluid contamination of the specimen and adversely affect test results.

Drops of blood should be allowed to flow freely into the tube and down the walls of the **MiniCollect®** Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. For correctly filled tubes, observe fill mark. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

Disposal

- The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed.
- Always wear gloves during blood collection and disposal.
- Contaminated or filled blood collection tubes must be disposed of in suitable biohazard disposal containers, which can then be autoclaved and incinerated afterwards.

Sample Stability and Storage

Remix K2EDTA samples immediately prior analysis to avoid result variations.
Refer to the instrument assay instructions for use or literature for information on the correct sample material, correct storage and stability.

Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in platelet clumping, clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

Label Information

| | | | | |
|---|--------------------------|---|------------------------------|--|
|  | Item number |  | Temperature limit | Literature GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved, Standard-Sixth Edition. GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. |
|  | LOT number: Batch number |  | In Vitro Diagnostic Device | |
|  | Expiry Date. |  | Consult Instructions For Use | |
|  | Do Not Reuse |  | Manufacturer | |
| Rx only | Prescription device | | | |

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www.gbo.com/preanalytics
faq.pa@gbo.com
 (+43) 07583/6791

MiniCollect® K3E K3EDTA Tubes

For In Vitro Diagnostic Use
Instructions for use**Intended Use**

To collect, transport, store, and evaluate capillary blood specimens for hematology tests.

Product Description

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes. The interior of the tube wall is coated with tripotassium EDTA (K3EDTA). The EDTA binds calcium ions thus blocking the coagulation cascade.

| Tube Type | Matrix | Cap colour |
|---------------------------|-----------------|------------|
| MiniCollect® K3EDTA Tubes | Capillary blood | Lavender |

Product versions

MiniCollect® Tubes with optional 13x75mm carrier tubes

MiniCollect® Tubes pre-assembled with 13x75mm carrier tubes

Storage before use

Store tubes at 4–25°C (40–77° F). Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. evaporation of liquid additives, colouring, etc.) All tubes are able to withstand a freezing temperature of –20°C.

Specimen Collection and Handling**Equipment required but not provided**

- Labels for positive patient identification of samples
- Gloves and appropriate apparel for protection against exposure to blood borne pathogens
- Alcohol swab
- Dry gauze
- Warming device if required, depending on the volume of blood needed and the test(s) to be performed
- Adhesive plaster or bandage
- Biohazard disposal container
- Safety lancet

Recommended Order of Draw

If multiple specimens are to be collected, including EDTA specimens, the EDTA specimen is drawn first to ensure adequate volume and accurate hematology test results. Specimens with other additives are collected next; serum specimens are collected last.

Cap Removal

The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop.

Specimen Collection

For the collection of blood by skin puncture, please refer to the instructions for use for the lancet device used.

Blood flow from the skin puncture site is increased by holding the puncture site in a downward position. The first drop of blood should be wiped away with a gauze pad if appropriate for the test method, as it is likely to contain excess tissue fluids. Apply gentle, intermittent pressure to the surrounding tissue (or proximal to the puncture site when using a finger). Strong, repetitive pressure (milking), must not be applied. It may cause hemolysis or tissue-fluid contamination of the specimen and adversely affect test results.

Drops of blood should be allowed to flow freely into the tube and down the walls of the MiniCollect® Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. For correctly filled tubes, observe fill mark. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

Disposal

- The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed.
- Always wear gloves during blood collection and disposal.
- Contaminated or filled blood collection tubes must be disposed of in suitable biohazard disposal containers, which can then be autoclaved and incinerated afterwards.

Sample Stability and Storage

Remix K3EDTA samples immediately prior analysis to avoid result variations.

Refer to the instrument assay instructions for use or literature for information on the correct sample material, correct storage and stability.

Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in platelet clumping, clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

Label Information

| | | | | |
|---|--------------------------|---|------------------------------|--|
|  | Item number |  | Temperature limit | Literature GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved. Standard-Sixth Edition. GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. |
|  | LOT number: Batch number |  | In Vitro Diagnostic Device | |
|  | Expiry Date. |  | Consult Instructions For Use | |
|  | Do Not Reuse |  | Manufacturer | |
| Rx only | Prescription device | | | |

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www.gbo.com/preanalytics
faq.pa@gbo.com
 (+43) 07583/6791

Intended Use

To collect, transport, separate and process capillary blood for testing plasma in the clinical laboratory.

Product Description

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes. The interior of the tube wall is coated with lithium heparin. The anticoagulant heparin activates antithrombin, thus blocking the coagulation cascade and producing a whole blood / plasma sample making it ideal for rapid analysis and analysis of blood from patients on anticoagulant therapy. Lithium Heparin Sep Tubes contain a gel on the bottom of the tube which moves upwards during centrifugation to form a stable barrier between the cells and the plasma.

| Tube Type | Matrix | Cap colour |
|---|-----------------|------------|
| MiniCollect® Lithium Heparin Tubes | Capillary blood | Green |
| MiniCollect® Lithium Heparin Sep Tubes (with gel) | Capillary blood | Mint green |

Product versions

MiniCollect® Tubes with optional 13x75mm carrier tubes
MiniCollect® Tubes pre-assembled with 13x75mm carrier tubes

Storage before use

Store tubes at 4–25°C (40–77° F). Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. evaporation of liquid additives, colouring, etc.) All tubes are able to withstand a freezing temperature of –20°C.

Specimen Collection and Handling**Equipment required but not provided**

- Labels for positive patient identification of samples
- For testing bilirubin, **MiniCollect®** tubes must be used with the amber carrier tube or other protective cover.
- Gloves and appropriate apparel for protection against exposure to blood borne pathogens
- Alcohol swab
- Dry gauze
- Warming device if required, depending on the volume of blood needed and the test(s) to be performed
- Adhesive plaster or bandage
- Biohazard disposal container
- Safety lancet

Recommended Order of Draw

If multiple specimens are to be collected, including EDTA specimens, the EDTA specimen is drawn first to ensure adequate volume and accurate hematology test results. Specimens with other additives are collected next; serum specimens are collected last.

Cap Removal

The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop.

Specimen Collection

The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop. For the collection of blood by skin puncture, please refer to the instructions for use for the lancet device used.

Blood flow from the skin puncture site is increased by holding the puncture site in a downward position. The first drop of blood should be wiped away with a gauze pad if appropriate for the test method, as it is likely to contain excess tissue fluids. Apply gentle, intermittent pressure to the surrounding tissue (or proximal to the puncture site when using a finger). Strong, repetitive pressure (milking), must not be applied. It may cause hemolysis or tissue-fluid contamination of the specimen and adversely affect test results.

Drops of blood should be allowed to flow freely into the tube and down the walls of the **MiniCollect®** Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. For correctly filled tubes, observe fill mark. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

Centrifugation

Ensure that tubes are properly seated in the centrifuge carrier. **MiniCollect®** tubes are recommended to be centrifuged at 3000g for a period of 10 minutes. Centrifugation should be done at a temperature of 15°C-24°C (25°C / 77°F). It is not recommended to re-centrifuge tubes once the barrier has been formed. Other centrifugation settings may also provide acceptable separation. Barriers are more stable when tubes are spun in centrifuges with horizontal swing-out rotors rather than those with fixed angle heads.

Disposal

- The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed.
- Always wear gloves during blood collection and disposal.
- Contaminated or filled blood collection tubes must be disposed of in suitable biohazard disposal containers, which can then be autoclaved and incinerated afterwards.

Sample Stability and Storage

Refer to the instrument assay instructions for use or literature for information on the correct sample material, correct storage and stability. Heparin plasma should be separated from cells within 2 hours, either by collection and centrifugation with a gel tube or by transferring plasma into a secondary container if a gel tubes is not used.

Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in platelet clumping, clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

Label Information

| | | | | |
|---|--------------------------|---|------------------------------|---|
|  | Item number |  | Temperature limit | Literature GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition. GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. |
|  | LOT number: Batch number |  | In Vitro Diagnostic Device | |
|  | Expiry Date. |  | Consult Instructions For Use | |
|  | Do Not Reuse |  | Manufacturer | |
| Rx only | Prescription device | | | |

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www.gbo.com/preanalytics
faq.pa@gbo.com
 (+43) 07583/6791

MiniCollect® FX Sodium Fluoride / Potassium Oxalate Tubes

For In Vitro Diagnostic Use
Instructions for use

Intended Use

MiniCollect® FX Sodium Fluoride/ Potassium Oxalate Tubes are used for the determination of glucose and lactate in capillary blood.

Product Description

These MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes. The tubes contain potassium oxalate and sodium fluoride.

| Tube Type | Matrix | Cap colour |
|--|-----------------|------------|
| MiniCollect® Sodium Fluoride/ Potassium Oxalate Tubes - Not available in USA | Capillary blood | Grey |

Product versions

MiniCollect® Tubes with optional 13x75mm carrier tubes

MiniCollect® Tubes pre-assembled with 13x75mm carrier tubes

Storage before use

Store tubes at 4–25°C (40–77° F). Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. evaporation of liquid additives, colouring, etc.) All tubes are able to withstand a freezing temperature of –20°C.

Specimen Collection and Handling

Equipment required but not provided

- Labels for positive patient identification of samples
- Gloves and appropriate apparel for protection against exposure to blood borne pathogens
- Alcohol swab
- Dry gauze
- Warming device if required, depending on the volume of blood needed and the test(s) to be performed
- Adhesive plaster or bandage
- Biohazard disposal container
- Safety lancet

Recommended Order of Draw

If multiple specimens are to be collected, including EDTA specimens, the EDTA specimen is drawn first to ensure adequate volume and accurate hematology test results. Specimens with other additives are collected next; serum specimens are collected last.

Cap Removal

The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop.

Specimen Collection

For the collection of blood by skin puncture, please refer to the instructions for use for the lancet device used.

Blood flow from the skin puncture site is increased by holding the puncture site in a downward position. The first drop of blood should be wiped away with a gauze pad if appropriate for the test method, as it is likely to contain excess tissue fluids. Apply gentle, intermittent pressure to the surrounding tissue (or proximal to the puncture site when using a finger). Strong, repetitive pressure (milking), must not be applied. It may cause hemolysis or tissue-fluid contamination of the specimen and adversely affect test results.

Drops of blood should be allowed to flow freely into the tube and down the walls of the MiniCollect® Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. For correctly filled tubes, observe fill mark. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

Centrifugation

Ensure that tubes are properly seated in the centrifuge carrier. MiniCollect® tubes are recommended to be centrifuged at 3000g for a period of 10 minutes. Centrifugation should be done at a temperature of 15°C-24°C (25°C / 77°F). It is not recommended to re-centrifuge tubes once the barrier has been formed. Other centrifugation settings may also provide acceptable separation. Barriers are more stable when tubes are spun in centrifuges with horizontal swing-out rotors rather than those with fixed angle heads.

Disposal

- The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed.
- Always wear gloves during blood collection and disposal.

Contaminated or filled blood collection tubes must be disposed of in suitable biohazard disposal containers, which can then be autoclaved and incinerated afterwards.

Sample Stability and Storage

Refer to the instrument assay instructions for use or literature for information on the correct sample material, correct storage and stability. Fluoride is known to cause an increase in haemolysis. For further information on substances that may interfere, please consult the assay's instructions for use.

Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in platelet clumping, clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

Label Information

| | | | | |
|---|--------------------------|---|------------------------------|--|
|  | Item number |  | Temperature limit | Literature GP41-A6 Procedures for the Collection of Diagnostic. Blood Specimens by Venipuncture; Approved. Standard-Sixth Edition. GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition |
|  | LOT number: Batch number |  | In Vitro Diagnostic Device | |
|  | Expiry Date |  | Consult Instructions For Use | |
|  | Do Not Reuse |  | Manufacturer | |
| Rx only | Prescription device | | | |

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www.gbo.com/preanalytics
 faq.pa@gbo.com
 (+43) 07583/6791

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Intended Use

To collect, transport, separate and process capillary blood for testing serum in the clinical laboratory.

Product Description

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes.

MiniCollect® Z No Additive Tubes do neither contain any anticoagulant nor clot activator.

| Tube Type | Matrix | Cap colour |
|----------------------------------|-----------------|------------|
| MiniCollect® Z No Additive Tubes | Capillary blood | White |

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Drops of blood should be allowed to flow freely into the tube and down the walls of the **MiniCollect®** Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

Blood in No Additive tubes should be allowed to fully clot prior to centrifugation.

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Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in delayed clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

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www.gbo.com/preanalytics
 faq.pa@gbo.com
 (+43) 07583/6791

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Intended Use

The MiniCollect® 9NC Coagulation Tube is intended for collection of citrate anticoagulated whole blood samples for coagulation assays.

Product Description

These MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes. The tubes are filled with a buffered trisodium citrate solution in a concentration of 0.109mol/L (3.2%).

| Tube Types | Matrix | Cap colour |
|---|--------------|------------|
| MiniCollect® 9NC Coagulation Tubes 3.2% - for venous blood only | Venous blood | Light blue |

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- Labels for positive patient identification of samples
- Gloves and appropriate apparel for protection against exposure to blood borne pathogens
- Alcohol swab
- Dry gauze
- Adhesive plaster or bandage
- Biohazard disposal container
- Needleless transfer device
- Venous blood collection device and accessories

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Specimen Collection

For the collection of venous blood, please refer to your institution's policies. The specimen should be subsequently transferred to the tube by means of a safe, needleless system (such as a plastic syringe).

General Handling

Drops of blood should be allowed to flow freely into the tube and down the walls of the MiniCollect® Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the tube. For correctly filled tubes, observe fill mark. After collection, close the tube with the original cap, an audible click indicates correct closure. Invert 4-5x until the blood completely mixes with the additive. Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

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Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in clotted samples and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

Label Information

| | | | | |
|---|--------------------------|---|------------------------------|---|
|  | Item number |  | Temperature limit | Literature GP41-A6 Procedures for the Collection of Diagnostic. Blood Specimens by Venipuncture; Approved. Standard-Sixth Edition. GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. H21-A5 Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. |
|  | LOT number: Batch number |  | In Vitro Diagnostic Device | |
|  | Expiry Date. |  | Consult Instructions For Use | |
|  | Do Not Reuse |  | Manufacturer | |
| Rx only | Prescription device | | | |

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Intended Use

To collect, transport, separate and process capillary blood for testing serum in the clinical laboratory.

Product Description

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is colour coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes.

The **MiniCollect®** Z Serum Tubes contain a blood clotting activator that functions to initiate the clotting process. Serum tubes should be allowed to clot in an upright position for at least 30 minutes prior to centrifugation. Serum Separator Tubes contain a gel on the bottom of the tube which moves upwards during centrifugation to form a stable barrier between the clot and the serum.

| Tube Type | Matrix | Cap colour |
|--|-----------------|------------|
| MiniCollect® Z Serum Tubes | Capillary blood | Red |
| MiniCollect® Z Serum Sep Tubes (with gel) | Capillary blood | Gold |

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Storage before use

Store tubes at 4–25°C (40–77° F). Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. evaporation of liquid additives, colouring, etc.). All tubes are able to withstand a freezing temperature of –20°C.

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Cap Removal The ribbed area around the tube top indicates the cap opening location. Remove the cap by applying gentle pressure on the cap in an upward direction. The triangle indicates the opposite position of the scoop.

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Drops of blood should be allowed to flow freely into the tube and down the walls of the **MiniCollect®** Tube. If a drop becomes lodged inside the scoop or to mix the content during collection, gently tap the tube on a hard surface. Do not shake and avoid flicking the open tube. For correctly filled tubes, observe fill mark. After collection, close the tubes with the original cap, an audible click indicates correct closure. Invert until the blood completely mixes with the additive (approximately 8-10x). Vigorous shaking may cause foaming and hemolysis. The patient and the patient's blood sample must be positively identified at the time of collection. The specimen should be labelled immediately following collection and mixing.

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Precautions/Cautions

- Insufficient or delayed mixing in tubes with additives may result in delayed clotting and/or incorrect test results.
- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.
- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only

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

Microvette® 100/200



User Guide

1. Collect blood with the assembled End-to-End 100 µl or 200 µl capillary. For collection, best results are achieved if the Microvette® is held in horizontal or slightly inclined position.



2. Collection is complete when the End-to-End capillary is entirely filled with blood.

3. Hold the tube upright to allow the blood to flow from the capillary into the Microvette®.

4. Turn the cap to remove and discard the preassembled capillary as one unit.

5. Remove the cap from the base and seal the Microvette® ('click' position).


6. Mix sample thoroughly by inverting the Microvette®.

All patient blood specimens should be treated with standard precautions.


Wear gloves!

Microvette® 300/500


User Guide


1. Turn the cap to remove it from the Microvette®.


2. Attach the cap to the base of the Microvette®.

3. Collect blood by gravity-flow using any part of the collection rim.

4. Remove the cap from the base and seal the Microvette® ('click' position).







5. Mix sample thoroughly by inverting the Microvette®.

All patient blood specimens should be treated with standard precautions.

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** Label Z Serum will be changed to label CAT Serum according to ISO 6710:2017

Printed copies will not systemically be updated.
Valid at the date of issue. Information is subject to change.

Greiner Bio-One GmbH | Bad Haller Straße 32 | A-4550 Kremsmünster
Tel: (+43) 7583 6791-0 | Fax: (+43) 7583 6318 | E-mail: office@at.gbo.com

www.gbo.com/preanalytics