

# Dynabeads™ HLA CLASS I and CLASS II

Catalog nos. 21001D, 21002D, 21003, 21004D

Store at 2°C to 8°C

Publication No. MAN0009978

Rev. D.0

## Product contents

Cat. No.	Volume	No. Isolations
Dynabeads™ HLA CLASS I		
21001D	2 mL	20
21002D	5 mL	50
Dynabeads™ HLA CLASS II		
21003	2 mL	20
21004D	5 mL	50

Dynabeads™ HLA Class I contains  $1.4 \times 10^8$  beads/mL, and Dynabeads™ HLA Class II contains  $2.8 \times 10^8$  beads/mL. The Dynabeads™ magnetic beads are suspended in phosphate buffered saline (PBS) containing 0.1% bovine serum albumin (BSA) and 0.02% sodium azide.

**Caution:** Sodium azide may react with lead and copper plumbing to form highly explosive metal azides.

## Product description

Dynabeads™ HLA Class I and Dynabeads™ HLA Class II are uniform, superparamagnetic polystyrene beads (4.5 µm diameter). Dynabeads™ HLA Class I is coated with a primary mouse anti-human monoclonal antibody specific for the CD8 membrane antigen on human cells. Dynabeads™ HLA Class II is coated with a primary mouse anti-human monoclonal antibody specific for the HLA Class II membrane antigen on human cells.

**Dynabeads™ HLA Class I** is designed for rapid isolation of CD8<sup>+</sup>T lymphocytes directly from peripheral blood. The isolated cells can be used directly in microcytotoxic HLA Class I (A, B, C) assays.

**Dynabeads™ HLA Class II** is designed for rapid isolation HLA Class II<sup>+</sup> cells (mostly B lymphocytes) directly from peripheral blood. The isolated cells can be used directly in microcytotoxic HLA Class II (DR, DQ, and DP) assays.

The immunomagnetic cell isolation technique using the Dynabeads™ HLA Class I or the Dynabeads™ HLA Class II product enables a rapid and specific isolation of target cells directly from

whole blood in only 15 min. The target cells are bound to the beads after a short incubation, and the bead-bound cells are isolated and washed using a magnet.

Dynabeads™ HLA Class I and Dynabeads™ HLA Class II yield purity up to 99% and viability up to 95% of the isolated cells.

The method is simple, rapid, and specific. When bound to the cell surface, the Dynabeads™ magnetic beads do not interfere with the binding of HLA-typing antibodies and complement to the cells.

**Note:** Viability of Dynabeads™ HLA Class II isolated cells after microcytotoxic incubation is typically lower than its Dynabeads™ HLA Class I counterpart.

## Required materials

- DynaMag™ Magnet (Go to [thermofisher.com/magnets](http://thermofisher.com/magnets) for recommendations).
- Mixer allowing tilting and rotation of tubes.
- Blood collecting tubes containing anti-coagulant (e.g. acid citrate-dextrose [ACD], heparin).
- 5 mL pipettes.
- Container for blood and disposal.
- Cooling device.
- See Table 1 for recommended buffers and solutions.

## General guidelines

- These products are for *in vitro* diagnostic use.
- Because mammalian cells are used in this procedure, follow appropriate laboratory techniques.
- Handle all samples as if they are capable of transmitting disease.
- Wear gloves and appropriate protection at all times.
- Avoid air bubbles (foaming) during pipetting.
- Carefully follow the recommended pipetting volumes and incubation times.

Table 1: Recommended buffers and solutions

PBS, pH 7.4, with sodium citrate*	Add 0.6% v/v sodium citrate.
PBS, pH 7.4, with EDTA*	Add 8% v/v EDTA.
Tris – balanced salt solution (TBSS), pH 7.4*	Add 2% heat inactivated fetal calf serum (FCS), or 2% heat inactivated human serum (HS) before use.

\*Note: Other buffers may be used for resuspension of cells, but they should be supplemented with Ca<sup>2+</sup>, Mg<sup>2+</sup>, and 2% FCS or 2% HS. Use analytical grade reagents.

## General guidelines, continued

- Use a mixer that provides tilting and rotation of the tubes to ensure that Dynabeads™ magnetic beads do not settle in the tube.
- Keep all buffers cold.
- Cells should be isolated within 24 h to achieve the highest yield and viability. However, blood up to 3-days-old can be used. Older samples may yield lower cell viability.
- Blood samples from uremic and leukemic patients should be processed within 24 h after collection.
- Use bead-bound cells immediately (cells can be kept at 2°C to 8°C for up to 24 h if necessary).

## Protocol

### Prepare sample

In some patients, the number of HLA Class I and II cells are low (e.g., hematological disorders, leukemic, uremic), or the serum may contain high quantities of soluble HLA Class I and II molecules, interfering with the cell isolation capacity of the beads. To give a satisfactory yield for HLA Class I and HLA Class II cells for typing, use the following procedure prior to "Isolate cells" for samples known to give low lymphocyte yields in conventional cell isolation procedures.

### For buffy coat isolation:

1. Centrifuge 10 mL of blood at 800–1000 × g for 5 min without brakes.
2. Discard plasma and transfer the buffy coat, and the very upper layer of erythrocytes to a clean 10 mL tube.
3. Add 5 mL cold PBS with 0.6% sodium citrate, pH 7.4 to the buffy coat. Mix gently.
4. Continue to "Isolate cells".

### For whole blood isolation:

1. Collect 5 mL blood (ACD) in a standard blood collecting tube, and cool at 2°C to 8°C to prevent attachment of phagocytic cells to the beads.
2. Add 5 mL cold PBS with 0.6% sodium citrate to each sample.
3. Continue to "Isolate cells".

## Isolate cells

1. Resuspend the Dynabeads™ magnetic beads (i.e., vortex for >30 sec, or tilt and rotate for 5 min). **Note:** keep beads in the rotator until use. **For better performance, wash beads with PBS/citrate buffer** (mix 100 µL of beads with 4–5 mL PBS-citrate solution, place tube in the magnet for 2 min, remove supernatant, resuspend beads in 100 µL of PBS-citrate).
2. Add 100 µL of beads to the blood sample.
3. Mix for 5 min at 2°C to 8°C by gentle tilting and rotating. **Note:** Do not use end-over-end rotation, as this may damage target cells and increase background in the microcytotoxic assay.
4. Add 5 mL cold PBS with 0.6% sodium citrate, pH 7.4, pipet up and down a few times, before placing the tube in the magnet for 2 min.
5. Discard the supernatant without disturbing the pellet containing the bead-bound cells.

6. Remove the tube from the magnet and resuspend the bead-bound cells in 5 mL **cold PBS**, pH 7.4, with 0.6% sodium citrate.
7. Apply to the magnet for 2 min to collect the bead-bound cells, and discard the supernatant while the tube is still in the magnet.
8. Repeat step 6–7 four times, except use 5 mL **cold PBS**, pH 7.4, **without** sodium citrate for each wash.
9. Resuspend the bead-bound cells in 0.4 mL TBSS, pH 7.4, containing Ca<sup>2+</sup> and Mg<sup>2+</sup> with 2% heat inactivated FCS or 2% heat inactivated HS, or equivalent.  
**Note:** Other buffers may be used for resuspension of cells, but they should be supplemented with Ca<sup>2+</sup>, Mg<sup>2+</sup>, and 2% FCS or 2% HS. Use analytical grade reagents.
10. Use bead-bound cells immediately (cells can be kept at 2°C to 8°C for up to 24 h if necessary). To stain cells isolated with Dynabeads™ HLA Class I or Dynabeads™ HLA Class II for use in microcytotoxic assays, follow the instructions that are appropriate for the assay of your choice.

## Related products

Product	Cat. No.
DynaMag™-5 magnet	12303D
DynaMag™-15 magnet	12301D
MPC™-1 magnet	12001D

## Explanation of symbols

	Catalog number
	Contains sufficient for <i>n</i> tests
<b>Read SDS</b>	Read safety data sheet

Open vials have the same shelf life as printed on the lot-specific label, as long as stored at 2°C to 8°C, and contamination is avoided.

References are available on the product specific page on [thermofisher.com](http://thermofisher.com).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

## Limited Use Label License: No. 451: No right to resell

Notice to Purchaser: Resale of this product is expressly prohibited. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

## Valid Edition of the Instructions for use GB, DE, ES, FR, IT

**German:** Gültige Version der Bedienungsanleitung: Siehe Ausgabe der mit diesem Kit gelieferten Bedienungsanleitung. Diese Ausgabe der Bedienungsanleitung ist in Ihrer Landessprache erhältlich und kann kostenlos angefordert werden: 1. Telefonisch unter +1 888-821-4443, Nebenst. 2. Europäische Anwender erhalten auch unter der Rufnummer +44 192 528 2700 Unterstützung. 2. Per E-Mail an [HLATechsupport@lifetech.com](mailto:HLATechsupport@lifetech.com).

**French:** Voir le mode d'emploi révisé fourni dans ce kit. Si vous souhaitez consulter cette révision du mode d'emploi dans la langue du pays, vous pouvez vous la procurer gratuitement: 1. En téléphonant au 888-821-4443, poste 2. Les utilisateurs européens peuvent également appeler le +44 192-528-2700 pour obtenir de l'aide. 2. En envoyant un courriel à [HLATechsupport@lifetech.com](mailto:HLATechsupport@lifetech.com).

**Italian:** Si veda la revisione delle istruzioni per l'uso fornite insieme al kit. Se si necessita di una versione in altra lingua della presente revisione delle Istruzioni per l'uso, è possibile richiederla gratuitamente: 1. Telefonando al numero 888-821-4443, est.2. Gli utenti europei possono anche rivolgersi al numero +44 192-528-2700, per assistenza. 2. Scrivendo un'e-mail a [HLATechsupport@lifetech.com](mailto:HLATechsupport@lifetech.com).

**Spanish:** Consulte la revisión de las instrucciones de uso proporcionadas en este kit. Si necesita que esta revisión de las instrucciones de uso esté en el idioma específico de su país, puede obtenerla de manera gratuita: 1. Llamando al 888-821-4443, ext. 2. Los usuarios europeos también pueden llamar al +44 192-528-2700 para obtener ayuda. 2. Por correo electrónico escribiendo a [HLATechsupport@lifetech.com](mailto:HLATechsupport@lifetech.com).

Manufactured by Thermo Fisher Scientific Baltics UAB, V.A.Graiciuno 8, LT-02241 Vilnius, Lithuania. Thermo Fisher Scientific Baltics complies with the Quality System Standards ISO 9001 and ISO 13485.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

CORPORATE ENTITY: LIFE TECHNOLOGIES | CARLSBAD, CA 92008 USA | TOLL FREE IN USA 1.800.955.6288

© 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Revision history of MAN0009978

Rev	Date	Description
D.0	16 Jun 2021	Removal of fluorescence microscope staining procedures.
C.0	26 Jun 2017	Update of staining procedure
B.0	15 Oct 2015	Document rebranding.
A.0	23 Jan 2014	Baseline for document.

For support visit [www.thermofisher.com/support](http://www.thermofisher.com/support) or email [HLATechSupport@lifetech.com](mailto:HLATechSupport@lifetech.com)

**ThermoFisher**  
SCIENTIFIC

**thermofisher.com**

16 June 2021