



ACUSON Juniper ultrasound system

Quick Reference Guide

VB10 SW Release

[siemens-healthineers.com/juniper](https://www.siemens-healthineers.com/juniper)



At Siemens Healthineers, our mission is to enable healthcare providers to increase value by empowering them on their journey towards expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated five million patients worldwide everyday benefit from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics and molecular medicine as well as digital health and enterprise services.

We're a leading medical technology company with over 120 years of experience and 18,500 patents globally. With over 50,000 employees in more than 70 countries, we'll continue to innovate and shape the future of healthcare.

Table of content

System Overview	3
Getting Started	12
2D and M-mode	22
Color and Spectral Doppler	40
SieScape Panoramic Imaging	58
Measurements and Reports	60
Text, Arrows and Body Markers	70
Data and Image Management	76
Care and Cleaning	88



System Overview


1. Flat Panel Display Monitor
2. 2 USB Ports
3. Power On / Off (QuikStart standby)
4. Pull-out Keyboard (optional)
5. Rear Handle
6. Gel Warmer
7. Articulating Arm (tilt, swivel, raise, lower)
8. Touch Screen
9. Control Panel (height and swivel adjustable)
10. 5 Transducer Ports
11. Auxiliary CW Pencil Port
12. Lock / Steer / Swivel Wheels




System Overview

Power on and QuikStart Stand-by Mode

1. Press the **Power** key to turn the system on or to enter the system shutdown menu
2. Select **Shutdown** to turn the system off
3. Select **Restart** to reboot the system
4. Select **Standby** to enter QuikStart stand-by mode which allows the system to power off/on in less than 5 seconds




 System can remain in Standby Mode for up to 30 minutes before it must be connected to a power source.





Control Panel

1. Soft Key Rotary Controls
2. Patient Study Control
3. Exam Selection
4. Review
5. Report
6. Steer
7. Select
8. UD (user-defined keys)
9. Angle
10. Escape
11. Text, Body Marker, Arrow
12. Pointer
13. Trackball
14. Update
15. Left/Right Set



System Overview


Control Panel

- 16. Pulsed Wave Doppler
- 17. Color Doppler
- 18. 2D Mode
- 19. M-mode
- 20. Continuous Wave Doppler
- 21. 3D/4D
- 22. eSielImage Multi-parametric Optimization
- 23. Dual *Atitikimas į 10.1 reikalavimą*
- 24. Caliper
- 25. Freeze
- 26. Depth / Focus
- 27. DGC Controls
- 28. Zoom
- 29. Print, Store and Clip



Monitor Screen Layout


1. Patient Banner
2. Orientation Marker
3. Imaging Parameters
4. Image Screen Menu Area
5. Measured Results (movable)
6. Image Display
7. User-Defined Key Functions
8. Trackball Functions
9. Thumbnail Display
10. Status Icons




System Overview

Touch Screen Layout

1. Login / Logoff
2. System User Manual
3. System Configuration
4. Digital Keyboard
5. Study Display Options
6. Soft Key functions that correspond to rotary controls
7. Active Modes
8. Advanced Application Tabs
9. Add Exam
10. Blue Highlight = Active Selection
11. Grayed-Out = Unavailable Selection
12. Gray Highlight = Available Selection
13. Page Navigation
14. Transducers





Connecting Transducer

1. Ensure the lever lock is in the left, unlocked position
2. Align the connector and press into the transducer port
3. Place the lever into the right, locked position




System Overview

System Casters/Steering

Easy wheeling casters

- Individual brakes
- Lock/steer/swivel functionality






Getting Started

Patient Registration – Manual

1. Press the **Patient** Key



2. Select the **Patient** Tab on the Touch Screen



3. Enter patient information
4. Select transducer and exam
5. Select **OK**

Tip: Transducer and Exam can also be selected from the Touch Screen after patient registration by selecting Exam on the Control Panel.



Register New Patient

New Patient
Edit Patient Data

Patient Demographics

Last Name: Patient Patient ID: 1234
First Name: A Date of Birth: DD / MM / YYYY Age:
3 MI: Pre: Suf: Gender: Male Female Other
Height: cm BSA: m²
Weight: kg BP: / mmHg

Exam

4 Transducer: 12L3 Accession No.:
Exam: Thyroid Indication:
Protocol: <None>

Medical Information

Additional Info:

Order Information


Performing MD: Referring MD:
Sonographer:

5

Getting Started

Patient Registration – Worklist

1. Press **Patient** Key 
2. Under the **Patient** Tab, select **Worklist** on the Touch Screen



3. Select the desired patient

4. If patient is not displayed, search using the drop-down menu or by dates



5. **Clear Search** will also refresh the overall list
6. Select the scheduled procedure
7. Press **Select** to open the patient registration page



System must be connected to a HIS/RIS server to access patient worklist.

Tip: You can also open the patient registration page by double clicking the scheduled Procedure.



Worklist

New Patient

Edit Patient Data

Worklist

Maps

Search

Patient Last Name **4** from 02/11/2017 US Only This system only Search with Cached List

Patient ID to 03/11/2017 Search Clear Search **5**

Worklist

Patient Name	Patient ID	Study Description	Accession Num	Date/Time
A_FirstName LastNa	39613873		53582877	171102/2043
B_FirstName LastNa	33299687		80032521	171102/2044

Scheduled Procedure(s)

Code Value	Description	State
ScheduledProcedureStepDescription	6	Suspended

Code Value Meaning


JDICOM_IRIS Status : TCP/IP - Ping succeeded.

Ping Echo Select Cancel **7**

Getting Started

Edit Study

1. Press the **Patient** Key 
2. Select the **Patient** Tab on the Touch Screen



3. Select **Edit Patient Data** on the image screen or **Edit Patient** on the Touch Screen
4. Modify the desired data fields
5. Select OK





Getting Started

Transducer / Study

Control Panel and Touch Screen

1. Select Exam 
2. Select the transducer on the Touch Screen
3. Select the desired Exam type






Store, Print, and Clip Store

Control Panel and Touch Screen

1. Press **Freeze**
2. If desired, roll trackball for frame-by-frame CINE review
3. Press **Print Store 1** to save a static image
4. Press **Print Store 2** to save and print¹
5. Press **Clip Store** to acquire and save a clip while scanning




¹ A compatible printer must be connected and assigned to the print function in the Configuration menu before using selection.

Tip: Clip Store length is adjustable on the Touch Screen during live scanning.

Getting Started

Review and End Exam

1. Press **Review**
2. Select the desired layout on the Touch Screen
3. Rotate the Page soft key or rotate **Select** to move between images
4. Select **End Exam** on the Touch Screen to end the study






The screenshot shows a software interface with a top navigation bar containing 'Patient', 'Browser', and 'Review' (highlighted in blue). Below this is a 'General Review' tab. The main area contains a grid of image sizes: 1x1 (highlighted in blue with callout 2), 1X2, 2x2, 3x3, 4x4, and 5x5. To the left of the grid are buttons for 'Select all', 'Deselect All', 'Show Selected', and 'Slide Show'. Below these are 'Teaching File', 'Export', 'Delete', 'PC Printer', 'USB B&W Printer', 'USB Color Printer', and 'AHP'. On the right side, there are buttons for 'Archived Report', 'Stress Echo Review', 'syngo VVI', 'Save', and 'Reset Marker'. A bottom status bar includes 'Page 25 / 25' (with callout 3), 'DICOM Print Color', 'Cine Speed', 'Cine Frame', 'Cine Start Marker', and 'Cine End Marker'. A vertical sidebar on the left contains icons for a user profile, help, settings, keyboard, 'Clear Screen', 'Hide Patient', 'Image Param', and 'End Exam' (with callout 4).

2D and M-mode

2D Control Panel




2D Controls	Function
2D	<i>Enters/exits 2D-mode or adjusts overall 2D gain/brightness</i>
eSiImage (eSiImage multiparametric image optimization)	<i>Automatically and continuously applies TEQ in real-time</i> Press twice to deactivate Note: Available on a live image only
DGC	<i>Manually adjusts gain/brightness at depth</i> Default all slide pods in the center
Zoom	<i>Magnifies the image</i> Rotate to increase or decrease image magnification Note: Zoom is available on a live and frozen image
HD Zoom	<i>Increases image size, detail resolution and frame rate in a region of interest (ROI)/Displays "HD Z" on the image screen</i> Press Zoom key to display the ROI and roll trackball to position Rotate (or, press right or left Set key) to resize the ROI and Press to activate Press again to deactivate Note: Available on a live image only

2D and M-mode

2D Control Panel

2D Controls	Function
# of Focal Zones	<p><i>Adjusts the number of focal zones</i></p> <p>Press the Depth/Focus control to select the desired number of focal zones Note: Availability varies by transducer, exam type and THI activation</p>
Depth/Focus	<p><i>Adjusts the scanning depth and the focal zone</i></p> <p>Press the toggle up/down to adjust the depth Rotate the toggle to adjust the focal zone Note: Direction of depth and focal zone control is customizable in the Configuration Menu/Custom Keys</p>
Dual	<p><i>Displays two separately acquired images side by side</i></p> <p>Press Dual to activate and toggle between images Press the 2D key to exit</p>





Increases the line density/frame rate

25

2D Sector Width

Press the **Set** key and roll trackball to resize the sector


Press the **Set** key and roll trackball to position


Press the **Set** key to lock

Note: Available on a live image only, not available on linear transducers

2D and M-mode

2D Touch Screen





2D Controls

Function

27

Trapezoid

Displays linear or trapezoid image format

Select to activate / deactivate

Note: Available on Linear transducers and on a live image only

Live Dual (Split)*

Activates two simultaneous images in dual format

Select to activate / deactivate

4B

Displays four separately acquired images

Select **4B** to activate and add each image

Select **2D** to exit

Note: Available on a live image only in 2D and Color Doppler

THI

(Tissue Harmonic Imaging)

Increases contrast and lateral resolution, reduces noise and clutter

Select to activate / deactivate

Note: Available on a live image only

Needle Visual

(Enhanced Needle
Visualization)

Activates multiple needle interrogation angles

Select to display available in-plane needle enhancement options

Note: Available on a live image only

*Note: This guide uses the term "Live Dual" to refer to Split imaging displayed as Split on the ultrasound system

2D and M-mode

2D Touch Screen

2D Controls	Function
Focus Width	<p><i>Increases the distance between multiple focal zones</i></p> <p>Rotate to increase / decrease the distance between multiple focal zones</p>
Wide FOV	<p><i>Activates a larger field of view</i></p> <p>Select to display the 220° FOV</p> <p>Note: Available on the 9MC3 transducer only</p>
Full FOV	<p><i>Changes the FOV back to 176°</i></p> <p>Note: Available on the 9MC3 transducer only</p>
FOV Size	<p><i>Activates the FOV size</i></p> <p>Select to increase / decrease the size of the FOV (also can use trackball and select keys)</p> <p>Note: Available in Full FoV on the 9MC3 transducer only</p>





Clarify VE
(Clarify vascular
enhancement [VE]
technology)

*Decreases artifacts, increases contrast resolution and boundary
detection in structures with flow*

Press to activate / deactivate

Note: Available on a live imaging only

Clarify VE Levels

Adjusts levels of sensitivity when using Clarify VE

Select higher level to increase sensitivity

Select lower level to decrease sensitivity

Note: Available on a live image only

U/D and R/L Flip

Adjusts orientation on the image screen

Press to activate / deactivate

Note: Available on a live or frozen image

Clip Settings

Selects the type and length of clip capture

Chronology – prospective / retrospective, capture type, acquisition rate
and trigger type

Note: Available on a live image only

CTI
(Custom Tissue Imaging)

*Optimizes the 2D image by adjusting the speed of sound for fatty
tissue versus dense tissue in Breast imaging*

Press to activate

Adjust the CTI levels based on the type of breast tissue

2D and M-mode

2D Touch Screen

2D Controls

Function


LVO Contrast
(Left Ventricular
Opacification)

Activates contrast imaging for the left ventricle

Press to activate / deactivate

Note: Only compatible with the 5P1 transducer in the Adult Echo preset





2D Soft Keys

2D Controls	Function
Compounding* (Advanced SieClear spatial compounding and SieClear multi-view spatial compounding)	<i>Enhances contrast resolution and boundary detection</i> Increase for more compounding Decrease for less compounding Note: Available on a live image only
TEQ Gain	<i>Adjusts default gain setting for eSieImage multiparametric imaging</i>
R/S (Resolution/Speed)	<i>Balances line density (spatial resolution) and frame rate</i> Increase to improve image detail Decrease to improve frame rate Note: Available on a live image only


*Note: This guide uses the term “Compounding” to refer to Advanced SieClear spatial compounding and SieClear multi-view spatial compounding displayed as Adv. SieClear or SieClear on the ultrasound system

2D and M-mode

2D Soft Keys

2D Controls	Function
DTCE (Dynamic TCE tissue contrast enhancement technology)	<i>Reduces speckle and enhances tissue contrast</i> Settings are Low , Medium and High Note: Available on a live image only
Image Size	<i>Adjusts the size of the image display</i> Increase for a larger image display, decrease for a smaller image display Available options are 0–4 Note: Available on a live image only
Tx Power (Transmit Power)	<i>Adjusts the level of acoustic energy delivered to the patient</i> Rotate clockwise to increase power Rotate counterclockwise to decrease power Note: Available on a live image only





Dyn R
(Dynamic Range)

Adjusts number of gray shades displayed

Rotate clockwise for a softer / grayer image

Rotate counterclockwise for more contrast

Note: Available on a live or frozen image

MultiHertz
**(MultiHertz multiple
frequency imaging)**

Adjusts transmit frequency of the active multi-frequency transducer

Use higher frequency for increased resolution

Use lower frequency to increase penetration

Note: Available on a live image only

Map

Selects a processing curve that assigns echo amplitudes to gray shades

Rotate to select desired Map


Tint

Colorizes the gray scale image

Rotate to select desired Tint

2D and M-mode

M-mode Control Panel



**M-mode Controls****Function****M
(M-mode)**

Enters/exits M-mode or adjusts overall M-mode gain

Press to activate/deactivate

Rotate to adjust overall M-mode gain

Roll trackball to position M-line

Note: Available on a live image only

Update/View

Alternates between 2D/M-mode formats

Press Update to move between modes

Note: Available on a live image only

Color M-mode

Enters/exits color M-mode

With color Doppler activated, press **M** (M-mode)


Position M-line and Press **M**

Press **M** (M-mode) again to deactivate

Note: Only available with cardiac transducers

2D and M-mode

M-mode Touch Screen





M-mode Controls	Function
Full M	<i>Changes the display format to full screen M-mode</i> Press to activate / deactivate
Angle AMM (Anatomical M-mode)	<i>Views the M-mode sweep based on patient anatomy, independent of transducer orientation utilizing angle rotation</i> Press to activate / deactivate
Marker AMM (Anatomical M-mode)	<i>Views the M-mode sweep based on patient anatomy, independent of transducer orientation utilizing freehand drawing of the angle</i> Press to activate / deactivate

2D and M-mode

M-mode Soft Keys

M-mode Controls	Function
Edge	<p><i>Sharpens the contours of a structure in the M-mode sweep</i></p> <p>Increase for more sharpness Decrease for less sharpness To deactivate, select zero</p> <p>Note: Available on a live image only</p>
Map	<p><i>Selects a processing curve that assigns echo amplitudes to gray shades in the M-mode sweep</i></p> <p>Note: Available on a live and frozen image</p>
Tint	<p><i>Colorizes the M-mode sweep</i></p> <p>Note: Available on a live and frozen image</p>
DR (Dynamic Range)	<p><i>Adjusts the number of gray shades displayed in the M-mode sweep</i></p> <p>Rotate clockwise for a softer /grayer sweep Rotate counterclockwise for more contrast</p> <p>Note: Available on a live or frozen image</p>





Sweep

Adjusts the scrolling speed of the M-mode sweep


Increase to display fewer cardiac cycles

Decrease to display more cardiac cycles

Note: Available on a live or frozen image

Color and Spectral Doppler

Color Doppler Control Panel



**Color Doppler Controls****Function****C
(Color Doppler)**

Enters/exits color Doppler or adjusts overall color Doppler gain

Press to activate/deactivate color Doppler mode

Rotate to increase/decrease color Doppler gain

Note: Available on a live image only

Color Box Size/Position

Adjusts the color ROI

Increases frame rate by decreasing the color ROI width or depth

Press **Set** to alternate between Size/Position

Roll the trackball to resize or reposition color ROI

Note: Available on a live image only

Steer


Steers the Color ROI

Toggle the **Steer** key left/right to steer

Note: Steering angle is transducer dependent

Color and Spectral Doppler

Color Doppler Touch screen



Color Doppler Controls	Function
CDE (Color Doppler Energy)	<i>Applies color based on flow amplitude</i> Note: Available on a live image only
Color Display	<i>Turns color display off when in color Doppler</i> Select to activate / deactivate Note: Color ROI remains on the image
Invert	<i>Reverses the colors depicting forward and reverse flow in the color ROI and the color bar</i> Select to activate / deactivate Note: Available on a live or frozen image
4B	<i>Displays four separately acquired images</i> Press to activate and add each image Note: Available on a live image only in 2D and Color Doppler
DTI (Doppler tissue imaging capability)	<i>Activates the Doppler tissue imaging feature for cardiac exams</i> Press to activate / deactivate Note: Available on a live image only

Color and Spectral Doppler

Color Doppler Soft Keys

Color Doppler Controls	Function
MultiHertz (MultiHertz multiple frequency imaging)	<i>Adjusts the color Doppler frequency independent of the 2D frequency</i> Rotate to increase/decrease Note: Available on a live image only
Map	<i>Selects a processing curve that assigns the velocity/velocity variance (Color) or flow amplitudes (Power) to a range of colors</i> Note: Available on a live and frozen image
Tx Power (Transmit Power)	<i>Adjusts the level of acoustic energy delivered to the patient in color Doppler mode</i> Rotate clockwise to increase power Rotate counterclockwise to decrease power Note: Available on a live image only, changes mechanical and thermal indices






Smooth	<i>Adjusts spatial averaging (smoothing) for color flow</i> Increase for flow sensitivity Decrease for color spatial resolution Note: Available on a live image only
Scale	<i>Adjusts the pulse repetition frequency (PRF)</i> Rotate clockwise to increase PRF Rotate counterclockwise to decrease PRF Note: Available on a live image only
Baseline	<i>Adjusts the range of flow velocities displayed above and below the baseline</i> Rotate to adjust the baseline of the active mode
Priority	<i>Prioritizes the display of grayscale or color flow data</i> Increase for greater color fill Decrease to prioritize grayscale and less color fill
Persist (Persistence)	<i>Adjusts the time that the color flow data remains in the ROI before being replaced by another color</i> Increase for greater color fill Decrease to increase temporal resolution Note: Available on a live image only

Color and Spectral Doppler


Color Doppler Soft Keys

Color Doppler Controls	Function
Filter	<i>Adjusts the color wall filter</i> Increase to reduce motion artifacts Decrease to increase sensitivity to low flow states Note: Available on a live image only
R/S (Resolution/Speed)	<i>Balances line density (spatial resolution) and frame rate</i> Increase for image detail Decrease to increase frame rate Note: Available on a live image only





Spectral Doppler Control Panel



Color and Spectral Doppler

Spectral Doppler Soft Keys

Spectral Doppler Controls	Function
PW (pulsed wave) Doppler and CW (continuous wave) Doppler	<p><i>Enters/exits spectral Doppler or adjusts overall spectral Doppler gain</i></p> <p>Press to enter/exit Roll trackball to position Doppler gate Press PW or CW to activate the strip Rotate PW to adjust overall spectral Doppler gain Note: Available on a live image only</p>
Update/View	<p><i>Alternates between the Doppler spectrum and the 2D image</i></p> <p>Press once to alternate between Doppler / 2D Note: Available on a live image only</p>
Steer	<p><i>Adjusts the PW Doppler angle</i></p> <p>Toggle left/right to adjust the steering angle Note: Available on linear transducers and on a live image only</p>






Angle	<i>Adjusts angle correction</i> Press to cycle through 60/60/0° angle selections Rotate to adjust values in 1° increments Note: Available on a live or frozen image
eSiImage (Auto Doppler Optimization)	<i>Optimizes spectral Doppler gain, baseline, scale and dynamic range</i> With Doppler spectrum active, press eSiImage on the Control Panel to activate Press twice to exit Note: Available on a live image only

Color and Spectral Doppler

Spectral Doppler Touch Screen




Spectral Doppler Controls	Function
Full D	<i>Displays Full Screen Doppler</i> Select to activate / deactivate Note: Available on a live image only
Invert	<i>Reverses the Doppler spectrum (negative values display above the baseline)</i> Select to activate / deactivate Note: Available on a live or newly frozen image only
Auto Stat	<i>Automatic tracing and display of PW Doppler measurements</i> Select to activate / deactivate Select settings for data above, below or on both sides of the baseline Note: Available on a live or newly frozen image and in Caliper mode
Sensitivity	<i>Adjusts the sensitivity of the Trace tool when in Auto Stats</i> Increase for higher sensitivity. Choices are 0–7 Note: Available when Auto Stat is activated only
Mean Trace	<i>Displays graphic traces of the mean (yellow) with maximum (blue) velocities or frequencies on the spectrum</i> Note: Available when Auto Stat is activated only

Color and Spectral Doppler

Spectral Doppler Touch Screen

Spectral Doppler Controls	Function
Trace Direction	<p>Select to choose trace direction</p> <p>Choose above, below or both sides of the baseline</p> <p>Note: Available in Auto Stat only</p>
Live/Live* (Simultaneous [Triplex])	<p>Displays both the 2D image and spectral Doppler tracing in real-time</p> <p>Press to activate / deactivate</p> <p>Note: Available on a live image only</p>
DTI (Doppler Tissue imaging capability)	<p>Activates the Doppler Tissue imaging feature for cardiac exams</p> <p>Press to activate / deactivate</p> <p>Note: Available on a live image only</p>
Update Rate	<p>Defines the interval for refreshing the 2D-mode image</p> <p>Select Off, 1 sec, 2 sec, 4, sec or 8 sec</p> <p>Note: Available on a live image only</p>





Mute

Mutes the Doppler audio

Select to activate / deactivate

Note: Available on a live image only


*Note: This guide uses the term "Live/Live" to refer to Simultaneous Triplex imaging displayed as Simult on the ultrasound system

Color and Spectral Doppler

Spectral Doppler Soft Keys

Spectral Doppler Controls	Function
MultiHertz (MultiHertz multiple frequency imaging)	<i>Adjusts spectral Doppler frequency independent of the 2D or color frequency</i> Rotate to adjust Note: Available on a live image only
Map	<i>Selects a processing curve that assigns echo amplitudes to gray shades in the Doppler spectrum</i> Note: Available on a live and frozen image
Tx Power (Transmit Power)	<i>Adjusts the level of acoustic energy delivered to the patient in spectral Doppler mode</i> Rotate clockwise to increase power Rotate counterclockwise to decrease power Note: Available on a live image only, changes mechanical and thermal indices displayed





Scale	<p><i>Adjusts the pulse repetition frequency (PRF)</i></p> <p>Rotate clockwise to increase PRF Rotate counterclockwise to decrease PRF</p> <p>Note: Available on a live image only</p>
Baseline	<p><i>Adjust the range of velocities displayed</i></p> <p>Rotate to adjust the baseline of the active mode</p> <p>Note: Available on a live or frozen image</p>
Dyn R (Dynamic Range)	<p><i>Adjust number of gray shades displayed</i></p> <p>Rotate clockwise for a softer / grayer spectrum Rotate counterclockwise for a blacker / whiter spectrum</p> <p>Note: Available on a live or frozen image</p>
Gate	<p><i>Adjusts sample volume size</i></p> <p>Rotate clockwise to increase size Rotate counterclockwise to decrease size</p> <p>Note: Available on a live image only</p>
Tint	<p><i>Colorizes the Doppler spectrum</i></p> <p>Select desired Tint from 0–11</p> <p>Note: Available on a live and frozen image</p>

Color and Spectral Doppler

Spectral Doppler Soft Keys

Spectral Doppler Controls	Function
Filter	<p><i>Adjusts the spectral Doppler signal filter</i></p> <p>Increase to reduce motion artifacts Decrease to increase sensitivity to low flow states</p> <p>Note: Available on a live image only</p>
Sweep	<p><i>Adjusts the scrolling speed of the Doppler spectrum</i></p> <p>Increase to display fewer cardiac cycles Decrease to display more cardiac cycles</p> <p>Note: Available on a live and frozen image</p>
Volume	<p><i>Adjusts the spectral Doppler audio volume</i></p> <p>Increase to raise audio level Decrease to lower audio level</p> <p>Note: Available on a live and frozen image</p>







Panoramic (SieScape Panoramic Imaging¹)

Perform a Panoramic Acquisition

1. Select the **SieScape*** tab to activate Panoramic imaging
2. Press the right or left **Set** keys to Acquire the image
3. Press **Freeze** or one of the **Set** keys to End acquisition, press **Freeze** again to reactivate acquisition if desired
4. Use the SieScape review controls, zoom rotary control and the trackball and select keys after completing the acquisition to size, rotate, pan, magnify or measure the image
5. Use the Cine function to recall single frames for review
6. Select **Exit** to deactivate Panoramic imaging





¹ *Optional feature*

***Note:** This guide uses the term “Panoramic” to refer to the SieScape Panoramic Imaging displayed as SieScape on the ultrasound system

Measurements and Reports

Unlabeled Measurements

1. Press **Caliper**



2. To change the measurement type, press **More** and select the desired function
3. Roll the trackball to position the first caliper
4. Press **Set** to anchor



5. Roll the trackball to position the second caliper
6. Press **Set** to complete the measurement and / or activate additional caliper sets




Unlabeled measurements do not enter into the exam report.

Tip: Measurements can be performed on frozen, CINE and stored images.





Unlabeled Measurements – Perform Measurements



Measurements and Reports

Labeled Measurements

1. Press **Caliper**



2. Select the desired measurement label
3. Roll the trackball to position the first caliper
4. Press **Set** to anchor



5. Roll the trackball to position the second caliper
6. Press **Set** to complete the measurement, send the labeled value(s) to the patient report and/or activate additional measurements
7. Rotate the Exam type soft key to display additional labeled measurement packages

Tip: The measurement label can be selected before or after calipers are placed.





The screenshot shows a medical software interface for vascular measurements. The top navigation bar includes a profile icon, a red circle with the number '2' next to the 'Carotid' menu, and sub-menus for 'Venous', 'Arterial', and 'TCD'. A 'Style Style1' dropdown is on the far right. Below the navigation bar, a grid of buttons allows selection of measurement sites: 'Right', 'Left', 'None', 'Prox', 'Mid', 'Dist', and 'None'. The 'Prox' button is highlighted in blue. Below this grid, another set of buttons includes 'Rt Prox CCA', 'Rt Prox ICA', 'Rt Prox ECA', 'Rt Prox Vert A', 'Innom A', 'Rt Bulb', and 'Rt Prox Subclav A'. A red circle with the number '2' highlights the 'Rt Prox CCA' button. A 'Result Compare' button is located on the right side. The bottom toolbar contains icons for 'Distance', 'Ellipse', 'Spline', 'Cross Angle', 'Vol 3 Dist', 'ET-%Stenosis', and 'D Flow', along with a 'More' button. The bottom-most bar includes 'Undo/Redo', 'Font Size', 'MDA On', 'Summarize MDA', 'Report', 'Delete', and 'Exam Vascular' (highlighted with a red circle and the number '7').

Measurements and Reports

Edit Measurements (Labeled and Unlabeled)

Edit Completed Measurement

1. Press **Update** to activate calipers and cycle between active calipers
2. Press the **Select Next** soft key to cycle between calipers sets
3. Press **Delete** soft key to remove active measurement
4. Press **Delete All** soft key to delete all measurements
5. Rotate **Undo/Redo** or the **Select** key to edit a trace measurement



Select Next

Delete All

Delete All

Undo/Redo





Reposition Measurement Results

1. Select **Pointer**
2. Click and drag MDA Results to desired position

Tip: Click and drag MDA Results to the bottom of the screen to display in a horizontal format.

Measurements and Reports

Report



Review the Report

1. Press the **Report** key to enter / exit report page 
2. Select **Results** to review additional measurements
3. Select **Print Report** to send to connected external printer
4. Press **Print Store 1** to store report as a thumbnail 
5. To edit the report, roll the trackball to the desired measurement and edit
6. Select **Print Preview** to visualize the revised summary
7. Select **Export PDF** to save the report on a USB



Only labeled measurements appear on the report; an * appears for any edited measurement.






Tip: Press Store Report to capture the report page(s) as an image(s).

Measurements and Reports

Report

Import Images into the Report

1. Press the **Report** key to enter / exit report page
2. Select **Images** on the Touch or Imaging Screen 
3. Select the display format on the imaging screen
4. Roll the trackball to the thumbnail image and press the **Set** key to transfer to the images section
5. To Delete the image from the report, hover the cursor over the image and click the "X" that appears in the upper right corner




Up to 20 images can be stored to the Summary.



Text, Arrows and Body Markers



Text

1. Press the **Text Body Marker** key on the control panel
2. Select the desired annotation on the Touch Screen
3. Select the **Keyboard** icon on the Touch Screen to activate the digital alphanumeric keyboard; press again to return to the text screen



Text Icons

Controls	Function
	Activates / deactivates text cursor and the text Touch Screen menu
	Sets default position for the text cursor
	Positions the cursor at the home location
	Temporarily hides text
	Deletes all text, pictograms or arrows

Controls	Function
	Deletes one line of text at the cursor
	Deletes the last word
	Displays the next or previous page, when available
	Press to activate / deactivate and make all letters typed capitalized / not capitalized

Text, Arrows and Body Markers

Body Markers

1. Press the **Text Body Marker** key on the Control Panel
2. Select the **Body Marker** tab
3. Select the desired Body Marker from the Touch Screen
4. Use the Body Marker widget on the Touch Screen to adjust the transducer orientation / position
5. Press **Set** to anchor
6. Reselect the **Text Body Marker** key on the Control Panel to readjust the Body Marker position
7. Press **Delete** or **Clear Screen** on the Touch Screen to delete the Body Marker



Tip: Body Marker options can be customized in the System Configuration Menu.





The screenshot shows a medical software interface for adding body markers. The interface is divided into several sections:

- Top Navigation:** Three tabs labeled "Text", "Body Marker" (highlighted with a blue bar and a callout circle 2), and "Arrow".
- Left Sidebar:** Contains icons for a user profile, a question mark, a gear, and a keyboard. Below these are buttons for "Clear Screen" (with callout circle 7), "Hide Patient", "Image Param", and "End Exam".
- Marker Grid:** A 3x5 grid of white line-art icons representing different body parts. The top-left icon (a face with a marker) is highlighted in blue. A callout circle 3 is positioned near the bottom-right icon of the grid.
- Central Canvas:** A large circular area showing a white line-art diagram of a human face. A blue dot is at the top, and a green vertical line with a white dot is on the mouth. A callout circle 4 is at the bottom of the diagram. A callout circle 7 is at the bottom of the canvas area.
- Bottom Control Bar:** Contains buttons for "Transducer Marker On", "Delete", "Rotate 0", and "Exam Thyroid".

Text, Arrows and Body Markers

Arrows

1. Press the **Text Body Marker** key on the control panel
2. Select the **Arrow** tab on the Touch Screen
3. Select the desired arrow type / size / color
4. Use the Arrow widget to adjust the arrow orientation
5. Roll the trackball to position
6. Press **Set** to activate additional arrows
7. Press **Delete** to delete the last arrow or **Clear Screen** to delete all markings
8. Press **Update** to cycle through multiple arrows
9. Press **Text Body Marker** to exit Arrow function



Tip: Unlimited arrows can be added.





The screenshot shows a software interface for medical imaging. At the top, there are three tabs: 'Text', 'Body Marker', and 'Arrow', with 'Arrow' selected and marked with a red circle '2'. Below the tabs is a grid of arrow icons in various directions and colors. A red circle '3' highlights a specific arrow icon. To the right of the grid is a circular diagram with a large green arrow pointing towards the top right. A red circle '4' is placed near the tip of the green arrow. On the left side, there is a vertical sidebar with icons for 'Clear Screen' (marked with a red circle '7'), 'Hide Patient', 'Image Param', and 'End Exam'. At the bottom, there is a toolbar with buttons for 'Arrow Color Green' (marked with a red circle '3'), 'Arrow Size Large' (marked with a red circle '3'), 'Delete', and 'Rotate'.



Data and Image Management


Review and Delete Images

Review Images on Current Study

1. Press **Review** on the Control Panel



2. Select the desired layout on the Touch Screen



3. Rotate the **Select** key or **Page** soft key to view more images



4. Press **Review** or **2D** to exit

Tip: Rotate the **Select** key to choose a protocol view and press the **Select** key to pause / resume the protocol.

Review Images from a Closed Study

1. Press **Patient** on the control panel




2. Select **Browser** Tab on the Touch Screen



3. Double click the desired study

4. Select the desired display format on the Touch Screen



5. Rotate the **Select** key or **Page** soft key to view more images




6. Press **Review** or **2D** to exit





Delete Images from a Current Study or in Review


1. Select the **Pointer** from the control panel
2. Select an image on the thumbnail panel to highlight
3. Press the **Delete** icon on the Touch Screen



Data and Image Management

Restart an Exam

1. Press **Patient** Key 
2. Select **Browser** Tab on the Touch Screen



3. Select patient study from list to highlight
4. Select **Restart**
5. Select **OK**



Images, clips and volumes can be added to a closed study within 24 hours.




Data and Image Management

Teaching File

1. Press **Review**



2. Select **Browser** tab



3. Double click the desired exam from **Local Disk** to load study
4. Roll the trackball, highlight an image or clip from the imaging screen (not thumbnail) and press the **Set** key
5. Repeat for all desired images
6. Press **Teaching File** on the Touch Screen

7. Under Teaching File, name the file and select **New** or select an existing file and select **Append**
8. Select **Close**
9. The teaching file is now displayed in the Patient Browser
10. To add images from multiple studies, follow steps 1–5, then select the desired teaching file
11. Press **Append** to add the additional images

Tip: Images and Clips stored in a teaching file removes all patient data.





Review

Select all Deselect All

Show Selected Slide Show

1x1	1X2	2x2
3x3	4x4	5x5

PC Printer USB B&W Printer USB Color Printer

Teaching File Export

Delete

Page 25 / 25

DICOM Print Color

Cine Speed 1

Cine Frame

Teaching File

Select Teaching File

Teaching File Name	Clips	Images
JUNIPER IMAGES 1	14	8

Append

File Name

JUNIPER2 New

Status

Close

Data and Image Management

Export Study

Export Study to Network

1. Select the patient name in the **Patient Browser**
2. Select the desired Output Device
3. Select **Send**

Tip: Images automatically send to PACs if the system is configured to send during the exam or at the end of the exam.





Patient Browser Show Files: 1 year old Source Media: [ID] Local Disk Disk Space: 10.2GB / 293.0GB

Search

Patient Last Name: [] from 2000/3/2018

Patient ID: [] to 2000/3/2018 Search Clear Search

Study List Total Studies: 49(49)

SE	Patient Name	Patient ID	Date/Time	Clips/Images	MB	Archived
<input type="checkbox"/>	18_02_14_MV041364RB165,	165	14/02/2018 17:34	3/5	63.7	HD,USB
<input type="checkbox"/>	18_02_MV073049RR195,	195	13/02/2018 11:00	10/63	267.3	HD,USB
<input checked="" type="checkbox"/>	18_20_01_MV031783JS OB,	TF_20180320194823	20/03/2018 19:48	11/8	203.2	HD,USB
<input type="checkbox"/>	2010_01_18_MV051641TH212,	TF_20180118072845	18/01/2018 07:28	4/19	90.9	HD,USB
<input type="checkbox"/>	2018_01_11_MV060161MS250,	TF_20180111134529	11/01/2018 13:45	0/23	35.8	HD,USB
<input type="checkbox"/>	2018_01_11_MV062938SS131,	TF_20180111135721	11/01/2018 13:57	0/21	54.5	HD,USB
<input type="checkbox"/>	2018_01_11_MV080187AC165,	TF_20180111140236	11/01/2018 14:02	4/9	104.3	HD,USB
<input type="checkbox"/>	2018_01_18_MV042667BR225,	TF_20180118022436	18/01/2018 02:24	5/18	125.5	HD,USB
<input type="checkbox"/>	2018_01_18_MV080189TH180,	TF_20180118073950	18/01/2018 07:39	15/23	340.3	HD,USB
<input type="checkbox"/>			15/02/2018 16:32	8/31	284.7	HD,USB
<input type="checkbox"/>			13/02/2018 06:22	0/13	28.3	HD,USB
<input type="checkbox"/>			12/01/2018 15:30	30/68	536.9	HD,USB
<input type="checkbox"/>			16/01/2018 14:04	3/13	58.5	HD,USB
<input type="checkbox"/>			17/01/2018 17:03	16/20	193.4	HD,USB
<input type="checkbox"/>			15/02/2018 22:04	13/11	309.0	HD,USB
<input type="checkbox"/>			16/02/2018 00:11	6/16	214.3	HD,USB
<input type="checkbox"/>			15/02/2018 22:39	7/5	180.0	HD,USB
<input type="checkbox"/>			25/10/2017 15:16	17/57	180.6	HD,USB

203.2 MB Selected Export Delete **2** Output Device: Network Share: PACS **3** Send

Data and Image Management

Export Study

Export Study to CD-R/DVD in DICOM

1. Insert CD-R/DVD
2. Select the patient name in the **Patient Browser**
3. Select **Export**
4. Select **DICOM** as the export format
5. Select **CD/DVD** as the destination
6. Select **Export**

Tip: Hold **CTRL** to select more than one patient name and select **De-identify** to anonymize a study before exporting.

Export Study to USB in PC Format (TIFF or AVI)

1. Insert USB
2. Select the patient name in the **Patient Browser**
3. Press **Export**
4. Select **Tiff/AVI** as the export format
5. Select a **USB Drive** as the destination
6. Select **Export**

Tip: USB icon is displayed in the top right of the image screen when USB is inserted correctly.





Export to USB or CD-R/DVD

Search

Patient Last Name: [] From: 01/06/2017

Patient ID: [] to: 31/08/2017 Search Clear Search

Study List

SI	Patient Name	Patient ID	Date/Time	Clip/Images	MB	Archived
	Patient, G	31_07_2017_16_59_49	31/07/2017 16:59	205	34.7	HDL/LSB
	Patient, E	345688	31/07/2017 16:44	353	44.4	HD
2	Patient, B	12345	31/07/2017 15:19	200	55.6	HD
	Patient, E	3456	31/07/2017 15:22	05	4.2	HD
	Patient, C	234567	31/07/2017 15:13	11	20.3	HD
	Patient, B	12345	31/07/2017 14:48	04	4.9	HD
	owen,	12234	31/07/2017 14:46	00	0.0	Unarchived
	.A	31_07_2017_14_42_25	31/07/2017 14:42	00	0.9	HD
	26_07_2017_12_32_45	26/07/2017 12:32	01	1.9	HD	
	GROUP D,	2017028-1	28/07/2017 10:39	815	109.8	HDL/LSB
	GYN CASE4,	27_07_2017_11_11	27/07/2017 11:13	85	115.6	HDL/LSB
	GYN CASE3,	27_07_2017_15_22_55	27/07/2017 15:22	718	101.9	HDL/LSB
*	07222017_20	27/07/2017 12:14	233	66.8	HDL/LSB	
*	07222017_10	27/07/2017 11:07	910	488.3	HDL/LSB	
*	07262017_30	26/07/2017 17:03	141	108.4	HDL/LSB	
*	07262017_20	26/07/2017 15:59	2013	370.3	HDL/LSB	
*	26_07_2017_14_54_24	26/07/2017 14:54	33	41.1	HDL/LSB	
*	07262017_10	26/07/2017 10:18	4307	403.6	HDL/LSB	

3 Export Delete Cancel Output Device: Not Active

Export

Destination: [H:] 5 Disk Space: 1.6GB / 7.9GB

Patient Name	Patient ID	Clips	Images	SI
rea	TF_20171025072320	0	1	0.4

Identification

De-Identify Anonymous

Format

DICOM 4

Tiff/Avi

Options

Finalize

Eject

Status

Ready

6 Export Close

Data and Image Management

Import Study

Import Study from CD-R/DVD

1. Insert CD-R/DVD
2. Select **CD/DVD** from Source Media in the Patient Browser
3. Select patient name on the CD/DVD
4. Select **Import**

Import Study from USB to Local Database

1. Insert USB
2. Select USB drive from Source Media in the **Patient Browser**
3. Select patient name on the USB
4. Select **Import**



Import from USB or CD-R/DVD

The screenshot shows the 'Patient Browser' application interface. At the top, there are navigation options like 'Show Files', 'Show All', and 'Source Media'. A search bar is present with fields for 'Patient Last Name' and 'Patient ID', and a date range from '04/09/2018' to '04/09/2018'. Below the search bar is a 'Study List' table. The table has columns: SE, Patient Name, Patient ID, DateTime, Clips/Images, MB, and Archived. One study is listed and selected, with a checkmark in the first column. At the bottom of the interface, there is a status bar showing '8.2 MB Selected' and an 'Import' button. On the right side, there are four preview windows showing image thumbnails.

SE	Patient Name	Patient ID	DateTime	Clips/Images	MB	Archived
✓	jjjj34jjjjrjlejmikm,	04_09_2018_17_05_52	04/09/2018 17:05	0/4	8.2	HD,USB

Care and Cleaning

Cleaning and Disinfecting the System

1. Power off the system and unplug the power supply
2. Use a clean gauze pad or lint-free cloth, lightly moistened with a mild detergent, to wipe the surface of the system
3. Use a clean, lint-free cloth to dry the surface
4. As required, use an approved disinfectant wipe to disinfect the system and accessories
5. Plug the power cord into the power supply

Tip: Ensure that fluids do not seep into any openings on the system.





Cleaning the Trackball

1. Rotate the ring around the trackball counter-clockwise and carefully lift the ring to remove
2. Remove the trackball
3. Clean the ring, trackball and trackball assembly with a cotton swab or lint-free pad moistened with a mild detergent solution
4. As required, use an approved disinfectant wipe to disinfect the ring, trackball and trackball assembly
5. After allowing the trackball components to completely dry place the trackball into the assembly
6. Place the ring over the trackball, aligning the tab with the point of attachment on the assembly and rotate clockwise until the ring snaps into place



Care and Cleaning

Air Filter

Cleaning the Air Filter

1. Power off the system and unplug the power supply
2. Push the air filter tray to release the locking mechanism and then pull the tray from the system
3. Rinse the air filter tray with running water and allow the filter to completely dry
4. Slide the dry air filter back into the system
5. Plug the power cord into the power supply

Tip: Air filter should be checked and cleaned weekly to maintain proper system cooling.







Note: Some functions shown in this material are optional and might not be part of your system.

ACUSON Juniper, SieScape, Clarify, CTI, DTCE, DTI, eSiImage, MultiHertz and SieClear are trademarks of Siemens Medical Solutions USA, Inc.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029, USA
Phone: 1-888-826-9702
siemens-healthineers.com/ultrasound