osteca®

(UAB „OSTECA“)

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

at day is more

10

24. If power supply or power cord is damaged, wires are frayed or exposed, do not use power.
Contact your Smith & Nephew representative for a replacement.
25. Canisters should be changed at least once a week, whenever there is a change in patient or in the event that canister contents reach maximum volume indication (300ml or 800ml fill line).
Do not wait for the Canister Full alarm to sound to change canister.
26. Canisters are single use. Do not reuse.
27. Do not apply SECURA™ No-sting barrier film wipes directly to open wounds. SECURA No-sting barrier film is flammable. Use in a well ventilated area. Avoid using around flames and sources of ignition. Keep out of reach of children. For external use only.
28. As with all adhesive products, apply and remove dressing carefully from sensitive or fragile skin to avoid blistering and skin stripping, especially after frequent dressing changes. Use of skin sealant may assist with protection of periwound skin.
29. If patient must be disconnected, the ends of the dressing tubing and canister tubing should be protected using tethered caps to avoid leakage of fluid and cross contamination.
30. Due to its smaller diameter, the RENASYS™-G 10Fr Round Drain Gauze Kit and Accessory Kit are not recommended for use with RENASYS TOUCH, as reduced pressure in the wound bed may lead to pooling or maceration.

Physician orders

Prior to placement of RENASYS TOUCH, the medical professional treating the wound must assess how to best use the system for an individual wound. It is important to carefully assess wound and patient to ensure clinical indications for Negative Pressure Wound Therapy (NPWT) are met.

All orders should include:

- Wound location, size and type
- Smith & Nephew Wound Dressing Kit
- Pressure settings
- Frequency of dressing changes
- Adjunctive dressings

→ *membrane failure*

FAQs

Nicely talpo ya absorbogian gellis.

RENASYS TOUCH FAQs

What is the flow rate of RENASYS TOUCH pumps?

9.25 Litres/minute free flow without canister and tub set.

What intrinsic mechanism is there in the RENASYS system to prevent backflow if machine is powered off?

Firstly the aeration disc (if Soft Port) or the CLP T-piece provides a flow of air that ensures that the air flow, once the machine is powered off, is away from the wound and towards the canister until the vacuum is dissipated. Secondly the gelling agent helps to retain the liquid in the canister (when the canister includes a solidifier).

How does RENASYS TOUCH monitor pressure?

nicatino netino dekupe
RENASYS TOUCH continuously monitors the pressure delivered and the air flow travelling through the system. If the expected airflow decreases below a pre-defined level, this correlates to a potential blockage and the device will inform the user with an alarm. RENASYS TOUCH also is able to auto correct the system when changes in pressure occur, to overcome a leak state. If a leak can't be compensated by the system, the device will inform the user with an alarm.

Can RENASYS TOUCH deliver and maintain pressure whilst above the wound?

Yes RENASYS TOUCH can deliver and maintain pressure when up to 90cm above the wound.

Can RENASYS TOUCH be taken on an aircraft?

Yes. Certified to all relevant standards that cover this environment.

What performance and safety checks are required between each patient?

- Pump cleaning following the directions provided in the service manual for RENASYS TOUCH.
- Device physical appearance check
- Functionality and alarm checks – complete functionality check, clinician mode, operation check, operation on battery check, blockage alarm check and restore presets as per directions provided in the RENASYS TOUCH service manual

What performance and safety checks are required on an annual basis?

Device vacuum checks as detailed in the RENASYS TOUCH Service manual.
Replace RENASYS TOUCH O-Ring and Odour filter as detailed in the RENASYS TOUCH Service Manual.

- The primary role of the first filter (the one facing towards the liquid in the canister) is to prevent liquid leaving the canister in a direction towards the pump. The secondary purpose of the first filter is to prevent bacteria from the canister leaving the canister in a direction towards the pump.
- The second filter provides odour protection.
- The third filter is at the exit of the canister to the pump. The primary role of the third filter is to prevent any bacteria that have not been stopped by the first filter progressing any further. The secondary purpose of the second filter is to prevent any liquid that might have breached the first filter, (e.g. due to mechanical failure), from progressing any further.

The mechanism by which these filters function is a combination of size exclusion and surface energy.

What materials are the RENASYS TOUCH canisters constructed from?

Body and Lid: Styrolution clearblend

Filters: The bacterial is hydrophobic and ePTFE membrane Solidifier

Clips: Valox 310

Tubing: Polyvinyl chloride (PVC)

Connector plug; nylon; zytel 101 F

Caps for Connector: Santoprene 271- 87

Are RENASYS TOUCH canisters sterile?

No these are not sterile

What is the maximum capacity of the RENASYS 300ml and 800ml canisters with and without solidifier?

The RENASYS canisters has a indication to measure its capacity. Within the design of the canister, there is also an additional space above this which is used to accommodate the impact of the gelling agent before the activation of the canister full alarm. Hence the solidifier does not affect the amount of fluid displaced in the canister. We do not specify the total capacity of the canister as this is dependent on the orientation of the canister/device when in use. The canister is designed to hold either 300ml or 800ml with or without gelling agent before the canister full alarm activation.

Primary Maximum Capacity:

Handwritten note: *Handwritten note: 800ml, ya su pababugjancu pelan mcluyk*

What are the battery's operating conditions when charging?

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




RENASYS TOUCH NPWT System

RENASYS TOUCH Device

Product description	Product code
RENASYS TOUCH - Pump - Power cord - RENASYS TOUCH Clinician User Manual - RENASYS TOUCH Home Healthcare User Manual 	66801260

RENASYS TOUCH Canisters

Product description	Product code
RENASYS TOUCH 300ml canister with solidifier - Sealed canister - Solidifier - Overflow guard with silicone tubing - Canister tubing	66801273
RENASYS TOUCH 300ml canister without solidifier - Sealed canister - Overflow guard with silicone tubing - Canister tubing	66801275
RENASYS TOUCH 800ml canister with solidifier - Sealed canister - Solidifier - Overflow guard with silicone tubing - Canister tubing	66801274

RENASYS TOUCH

Accessories

Product description	Product code
RENASYS TOUCH Carry straps 	66801276
RENASYS TOUCH Carry bag 	66801277
RENASYS TOUCH IV pole/Bed clamp	66801278
RENASYS TOUCH Transit case	66801279
RENASYS TOUCH Odor Filter	66801284
RENASYS TOUCH O-Ring	66801283

RENASYS MAX NPWT System

RENASYS MAX Devices

Product description	Product code
RENASYS EZ MAX (Professional healthcare facility model) - Pump - Canister holder kit - Power cord - RENASYS EZ MAX Clinician User Manual 	66801309
RENASYS EZ Max (Home healthcare model) - Pump - Canister holder kit - Power cord - RENASYS EZ MAX Clinician User Manual - RENASYS EZ MAX Home Healthcare User Manual 	66801310

RENASYS EZ PLUS/MAX Canisters

Product description	Product code
RENASYS EZ 800ml S-canister 10 kits per case - Sealed canister - Overflow guard with silicone tubing - Solidifier - Canister tubing	800ml with solidifier 66800912 800ml without solidifier 66801066
RENASYS EZ 250ml S-canister with solidifier 10 kits per case - Sealed canister - Solidifier - Overflow guard with silicone tubing - Canister tubing 	250ml with solidifier 66800913

RENASYS EZ PLUS/MAX Accessories

Product description	Product code
RENASYS EZ Transit case	66800587
RENASYS EZ Canister holder	66800560

RENASYS GO NPWT System

RENASYS GO Device

Product description	Product code
RENASYS GO - Pump - Power cord - Quick reference guide - Strap - RENASYS GO Clinician User Manual - RENASYS GO Home Healthcare User Manual 	66801496

RENASYS GO Canisters

Product description	Product code
RENASYS GO Large 750ml canister 5 per case - Sealed canister - Solidifier - Canister tubing	750ml with solidifier 66800916
RENASYS GO 300ml canister 5 per case - Sealed canister - Solidifier - Overflow guard with silicone tubing - Canister tubing	300ml with solidifier 66800914 300ml without solidifier 66800915

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
Product Portfolio

Advanced Wound Management

Wound Therapy: Advanced Wound Management

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FILE NR. 1


Device orientation during use

The device is designed to operate in the upright position. Device orientation, rate at which fluid enters the canister and how exudate solidifies can impact filter occlusion. Operation in the upright position optimizes canister volume and alarm functionality. The device should be orientated to face the user's position when in stationary use.

Proper use/correct orientation


The 300ml canister has a kickstand and rubber feet. The device can stand in the upright position with the 300ml canister attached. Open the kickstand for additional stability and to change the viewing angle.

The 800ml canister does not include a kickstand or rubber feet. The IV pole/bed clamp or carry strap accessories can be used to mount the device in the upright position.

Incorrect orientation

Caution: Operating the device in a face-down position could result in damage to the device and inadvertent changes to therapy settings. Operating the device in an inverted position could impact filter occlusion resulting in a blockage alarm and requiring a change of canister.

Confidential – For internal use only by Smith+Nephew employees and sales force. This information may not be publicly disclosed or distributed.

- The primary role of the first filter (the one facing towards the liquid in the canister) is to prevent liquid leaving the canister in a direction towards the pump. The secondary purpose of the first filter is to prevent bacteria from the canister leaving the canister in a direction towards the pump.
- The second filter provides odour protection.
- The third filter is at the exit of the canister to the pump. The primary role of the third filter is to prevent any bacteria that have not been stopped by the first filter progressing any further. The secondary purpose of the second filter is to prevent any liquid that might have breached the first filter, (e.g. due to mechanical failure), from progressing any further.

The mechanism by which these filters function is a combination of size exclusion and surface energy.

What materials are the RENASYS TOUCH canisters constructed from?

Body and Lid: Styrolution clearblend

Filters: The bacterial is hydrophobic and ePTFE membrane Solidifier

Clips: Valox 310

Tubing: Polyvinyl chloride (PVC)

Connector plug: nylon; zytel 101 F

Caps for Connector: Santoprene 271- 87

Are RENASYS TOUCH canisters sterile?

No these are not sterile

What is the maximum capacity of the RENASYS 300ml and 800ml canisters with and without solidifiers?

The RENASYS canisters has a indication to measure its capacity. Within the design of the canister, there is also an additional space above this which is used to accommodate the impact of the gelling agent before the activation of the canister full alarm. Hence the solidifier does not affect the amount of fluid displaced in the canister. We do not specify the total capacity of the canister as this is dependent on the orientation of the canister/device when in use. The canister is designed to hold either 300ml or 800ml with or without gelling agent before the canister full alarm activation.

Battery life and charging

What are the battery's operating conditions when charging?

Operational temperature 5°C to 40°C

Relative humidity 15% to 93% RH

Atmospheric pressure 700mbar to 1060 mbar

Smith & Nephew Medical Limited
 101 Hessle Road
 Hull
 HU3 2BN
 England

T: + 44 (0)1482 225181
 F: + 44 (0)1482 328326
 www.smith-nephew.com


Smith+Nephew

EUROPEAN DECLARATION OF CONFORMITY¹

Declaration confirms that the product listed below meets Regulation 2017/745 and is issued under the sole responsibility of Smith & Nephew Medical Ltd.²

Manufacturer's Name³	Smith & Nephew Medical Limited
Business Address⁴	101 Hessle Road, Hull, HU3 2BN, United Kingdom
Single Registration Number (SRN)⁵	GB-MF-000017580
European Authorised Representative⁶	Smith & Nephew Operations B.V.
Business Address⁷	Bloemlaan 2, 2132 NP Hoofddorp, Netherlands
Product Name: (see attached schedule for product codes/catalogue numbers) ⁸	RENASYS Touch Canisters
Intended Use⁹ See table for other languages	Intended to be used in conjunction with Smith & Nephew Negative Pressure Wound Therapy (NPWT) systems.
Conformity Assessment Procedure (Annex)¹⁰	N/A
Notified Body Name¹¹	N/A
Notified Body Number¹²	N/A
Verification Certificate(s)¹³	N/A

Signed on behalf of Smith & Nephew Medical Limited¹⁴

Signature¹⁵	DocuSigned by: 
Name¹⁶	Ken Fergusson

Smith & Nephew Medical Limited
 101 Hessle Road
 Hull
 HU3 2BN
 England

T: + 44 (0)1482 225181
 F: + 44 (0)1482 328326
www.smith-nephew.com

Smith Nephew

Position¹⁷	Senior Regulatory Affairs Manager
Date¹⁸	07-Apr-2022 17:58:33 BST
Location¹⁹	Hull
Declaration of Conformity Reference²⁰	DOC-WMTF-005/V2

Product Schedule²¹			
Product Code / Catalogue Number²²	Product Description or Product Variant²³	Risk Classification²⁴	Basic UDI²⁵
66801273	RENASYS Touch Canister with Solidifier 300ml	I	5000223SN000095RR
66801274	RENASYS Touch Canister with Solidifier 800ml	I	5000223SN000095RR
66801275	RENASYS Touch Canister without Solidifier 300ml	I	5000223SN000095RR


Standard/Common Specification Number²⁷
EN ISO 13485:2016
EN 60601-1:2005 3.1 edition
EN 60601-1-11:2015 2nd edition
EN ISO 14971:2019
EN 62366-1:2015
EN ISO 15223-1:2021
EN ISO 780:2015
ISO 10993-1:2018
EN ISO 20417:2021
EN ISO 14644-1:2015



Renasy-Touch
puller

RENASYS Accessories

lipnas
infocely

4 bowl laboratory

SKU	Description	Image
66800971	Y-Connector	
66800799	RENASYS Soft Port	
66801021	Foam Filler 10 x 12.5cm Pack 1	
66800394	NPWT LARGE DRAPE (FILM DRSG) 20cm X 30cm	
66800853	RENASYS TRANSPARENT FILM XL 38cm x 60cm	

SKU	Description	Image
66801082	Gel Patch 10cm x 7cm	
66802127	Medicated Gauze Procedure pack	
66801273 66801274 66801275	300ml Canister with Solidifier 800ml Canister with Solidifier 300ml Canister without solidifier	

Deuter
two parallel opening

Minimative
teleg's no how-
claw machine

Menu

14

abidlay co m'eta

smith&nephew

REF 66800971

GTIN (01) 00040565125541

LOT (10) 11419E

2019-04-24

destergetil pagamun
RENASYS° Y-connector

nestenli 7
Non-sterile

- | | |
|----------------------------------|---------------------------|
| (FR) Connecteur en Y Non stérile | (DK) Y-konnektor Usteril |
| (DE) Y-Konnektor Nicht steril | (NO) Y-kobling Usteril |
| (ES) Conector en Y No estéril | (FI) Y-liitin Epästeriili |
| (PT) Conector em Y Não estéril | (SE) Y-koppling Osteril |
| (IT) Connettore a Y Non sterile | (JP) Y コネクタ 未滅菌 |
| (NL) Y-connector Niet-steriel | |

Smith & Nephew Medical Ltd.
101 Hessel Road, Hull HU3 2BN England
www.smith-nephew.com
*Trademark of Smith & Nephew
Made in Mexico

(EC) (REP) Smith & Nephew Orthopaedics GmbH
Alemannerstraße 14, 78532 Tuttlingen, Germany

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3. Indications for use

EN

The RENASYS Y-connector is only intended to be used in conjunction with RENASYS Negative Pressure Wound Therapy (NPWT) system which is indicated for patients who would benefit from a suction pump (Negative

Pressure Wound Therapy), as it may promote wound healing via removal of fluids, including irrigation and body fluids, wound exudates and infectious materials.

Appropriate wound types include:

Wound types	RENASYS Y-connector configuration	
	Single wound	Two wounds of the same type
Chronic	✓	✓
Acute	✓	✓
Traumatic	✓	✓
Acute and dehiscent wounds	✓	✓
Ulcers (such as pressure or diabetic)	✓	✓
Partial-thickness burns	✓	✓
Flaps and grafts	✓	✓
Open abdomen	✓	✗

Indications ✗ Contraindications A Y-connector must not be used to connect dressings on different types of wounds.

4. Contraindications

The use of the RENASYS NPWT System is contraindicated in the presence of:

- Untreated osteomyelitis
- Exposed arteries, veins, organs or nerves

- Necrotic tissue with eschar present
- Malignancy in wound (with exception of palliative care to enhance quality of life)
- Non-enteric and unexplored fistulas
- Exposed Anastomotic sites

5. Warnings

The pump alarms may not sound when a blockage is present in a RENASYS Y-connector branch, so the patient must be monitored more frequently. Monitor the patient more frequently when treating infected wounds. If there are any signs of systemic infection or advancing infections in the area around the wound, contact the treating clinician immediately. Do not use a RENASYS Y-connector with breached or damaged packaging. Do not use the RENASYS Y-connector to treat more than two wounds. Do not use more than one Y-connector per pump.

- Do not connect two wounds of different etiologies (e.g., surgical dehiscence and explored fistula) via the RENASYS Y-connector due to the risk of cross contamination or impaired therapy.
- Do not use the same RENASYS Y-connector to connect two patients to a single pump.
- Care should be taken to ensure that the Y-connector does not lie in a position where it could cause pressure damage to the patient.
- Do not use the Y-connector on two wounds if the pump display refers to the use of the Y-connector on one wound only.

6. Precautions

As a condition of use, the RENASYS NPWT system should only be used by qualified and authorized personnel. The user must have the necessary knowledge of the specific medical application for which the RENASYS Y-connector is being used. The dressing seal may be lost or pooling may occur, without alarm activation, when an occlusion forms on the wound side of the dressing. Viscous, purulent or serosanguinous drainage may contribute to occlusion of the dressing. Regular monitoring of the pump and dressing is required to ensure full delivery of therapy and exudate removal. Ensure the wound dressing is free of air leaks, fully compressed and firm to the wound whenever therapy is active. Monitor the patient for any signs of local or systemic infection. Infected wounds require more frequent monitoring and may require more frequent dressing changes. As NPWT is not intended to directly treat infection, if there are any signs of systemic infection or advancing infection in the area surrounding the wound, contact the treating clinician immediately. When bathing or moving the patient must disconnect from the RENASYS Y-connector leaving it on the pump. Use the connector caps on the Soft Port/Drain adaptor and Y-connector to avoid leakage of wound fluid. Ensure the aeration disc located on the orange quick click connector is free of moisture before reactivation of therapy to ensure proper alarm

- functionality and prevent interruption in therapy.
- RENASYS Y-connector should only be used with Smith & Nephew authorized components. Use of the Y-connector with other products has not been proven safe and effective for use with RENASYS pump.
- NPWT should remain on for duration of treatment. The length of time a patient may be disconnected from the pump is a clinical decision based on individual characteristics of the patient and wound. Factors to consider include location of wound, volume of drainage, integrity of dressing seal, assessment of bacterial burden and patient's risk of infection.
- RENASYS Y-connector is single use only. Use of any part of this system on more than one patient may result in cross contamination that may lead to infection.
- When a RENASYS Y-connector is used to connect two RENASYS Soft Ports and/or drains to the pump in a two-wound configuration, the blockage alarm may not sound when a blockage is present. The patient should be monitored at least every 4 hours.
- When a RENASYS Y-connector is used to connect with RENASYS AB Abdominal Drainage Port, the Y-connector should be used for a single wound only.
- Do not use the RENASYS Y-connector to connect two wounds of different etiologies (e.g., surgical dehiscence and explored fistula) via the RENASYS Y-connector due to the risk of cross contamination or impaired therapy.
- Do not use the same RENASYS Y-connector to connect two patients to a single pump.
- Care should be taken to ensure that the Y-connector does not lie in a position where it could cause pressure damage to the patient.
- Do not use the Y-connector on two wounds if the pump display refers to the use of the Y-connector on one wound only.
- Do not use the RENASYS Y-connector to treat more than two wounds. Do not use more than one Y-connector per pump.
- Do not use the RENASYS Y-connector with breached or damaged packaging.
- Do not use the RENASYS Y-connector to treat more than two wounds. Do not use more than one Y-connector per pump.
- Do not use the RENASYS Y-connector to treat more than two wounds. Do not use more than one Y-connector per pump.

19

Smith & Nephew Medical Limited
101 Hessle Road
Hull
HU3 2BN
England

T: + 44 (0)1482 225181
F: + 44 (0)1482 328326
www.smith-nephew.com


Smith+Nephew

STATEMENT

We confirm that RENASYS TOUCH is suitable for connection to the 230V / 50 Hz electrical network. The kit includes a power cable for the device. The plug is adapted for rosettes used in Lithuania.

Disposable RENASYS tools listed below are suitable for the treatment of wounds under RENASYS TOUCH Negative Pressure Wound Therapy.

66801273	RENASYS TOUCH 300ml Canister with Solidifier
66801274	RENASYS TOUCH 800ml Canister with Solidifier
66800796	RENASYS-F Foam Dressing Kit - Large
66800795	RENASYS-F Foam Dressing Kit - Medium
66800794	RENASYS-F Foam Dressing Kit - Small
66800971	RENASYS Y-Connector
66800394	RENASYS Transparent Film Dressing
66802134	RENASYS TOUCH Non-Connect 4th Ed Device
66800799	RENASYS Soft Port
66801277	RENASYS TOUCH Carry Bag
66801276	RENASYS TOUCH Carry Strap

DocuSigned by:

Signature:

Name: Sam Atkinson, Principal Regulatory Affairs Specialist

Date: 24-MAY-2024

Registered No. 605496 in England and Wales, registered Office: PO Box 81, 101 Hessle Road, Hull, HU3 2BN, England
Smith & Nephew Medical Limited, acting as agent for T.J. Smith and Nephew, Limited.