**Monoclonal Mouse
Anti-Human
Androgen Receptor
Clone AR441**

Code M3562

Intended use

For In Vitro Diagnostic Use.

19 p.d.

Summary and explanation

Androgen receptor is an intracellular protein which belongs to a large family of hormone-induced transcription factors.^{2,3} Androgen receptor mediates the biological actions of physiological androgens such as testosterone and 5 α -dihydrotestosterone, which are essential for differentiation, development and maintenance of the male reproductive organs. Complexes of androgen and receptor regulate growth responses in accessory sex organs by modulating specific gene transcription.³

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Anti-androgen receptor is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide.

Clone: AR441 Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

Anti-human Androgen Receptor may be used at a dilution of 1:50 in the EnVision™+ method determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Synthetic peptide corresponding to amino acids 229–315 of the human androgen receptor coupled to keyhole limpet hemocyanin¹

Specificity

Monoclonal mouse anti-human androgen receptor, clone AR441 (anti-androgen receptor) specifically immunostains the nuclei of positive cells in human tissues. In Western blot analysis, the antibody identifies a 110 and 112 kD doublet in extracts of the metastatic prostate cancer cell line LNCap and in extracts of cells transfected with the gene for androgen receptor.¹ The antibody is unreactive with mouse androgen receptor in immunoblots of prostate extracts. Anti-androgen receptor is capable of recognizing and forming complexes with native androgen bound receptor.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Tissue specimens should be cut in to sections of approximately 4 μ m.

Pre-treatment: Please refer to EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004) for instruction on how to perform pretreatment of slides.

Cryostat Sections and Cell Smears: Anti-androgen receptor can be used for labeling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-androgen is nuclear.

Performance characteristics

Normal Cells

Androgen receptor has been localized by immunohistochemistry (IHC) to the nuclei of a variety of normal cells. In the prostate, nuclear immunoreactivity has been observed in the secretory cells, basal cells, smooth muscle cells and fibroblasts. Moderate immunostaining has been found in spermatogonia, spermatocytes and Sertoli cells of the testes, with weaker staining seen in Leydig and peritubular myoid cells. Strong to moderate expression of androgen receptor has been demonstrated by IHC in uterine myometrial cells, cervical squamous cells, ovarian stromal cells, and syncytiotrophoblasts and cytotrophoblasts of the placenta.⁴ Endometrial glandular and stromal cells also express androgen receptor at all phases of the menstrual cycle.^{4,5} In the breast, the nuclei of acinar cells, myoepithelial cells and stromal cells were reported to stain positively. Androgen receptor was also demonstrated in cardiac muscle, striated muscle and in smooth muscle of arteries, bladder and GI tract. In skin, positive staining has been observed in cells of the hair follicles, sweat glands and squamous cells of the epidermis. Other cells which have been shown to give moderate to strong immunostaining include hepatocytes of the liver and cells of the anterior and posterior pituitary.⁴

Tumor Cells

In adenocarcinomas of the prostate androgen receptor has been demonstrated by IHC to be present in the nuclei of most neoplastic cells with neuroendocrine differentiation.⁶ The majority of breast and ovarian carcinomas evaluated also express androgen receptor, as determined by IHC and biochemical assays.⁷⁻¹⁰ The nuclei of glandular and stromal endometrial cells have been shown to stain positively in adenomyosis and external endometriosis. Malignant cells in both glandular and solid lesions of grade II endometrial adenocarcinomas were found to be immunoreactive.⁵ Androgen receptor expression has also been demonstrated in meningiomas by Northern blot analysis and immunohistochemistry.¹¹

References

1. Report on file.
2. Zhuang YH, et al. Subcellular location of androgen receptor in rat prostate seminal vesicle and human osteosarcoma MG-63 cells J Steroid Biochem Mol Biol 1992; 41(3-8):693
3. Jänne OA, et al. Androgen receptor and mechanism of androgen action. Ann Med 1993; 25:83
4. Kimura N, et al. Immunocytochemical localization of androgen receptor with polyclonal antibody in paraffin-embedded human tissues. J Histochem Cytochem 1993; 41(5):671
5. Horie K, et al. Immunohistochemical localization of androgen receptor in the human endometrium, decidua, placenta and pathological conditions of the endometrium. Hum Reproduct 1992; 7(10):1461
6. Nakada SY, et al. The androgen receptor status of neuroendocrine cells in human benign and malignant prostatic tissue. Canc Res 1993; 53:1967
7. Schulz TJ, et al. Immunocytochemical analysis of androgen receptor in breast cancer and correlation with estrogen and progesterone receptors. Lab Invest 1993; 68(1):19A(#92)
8. Bryan RM, et al. Androgen receptors in breast cancer. Canc 1984; 54(11):2436
9. Chadha S, et al. An immunohistochemical evaluation of androgen and progesterone receptors in ovarian tumors. Hum Pathol 1993; 24(1):90
10. Kühnel R, et al. Androgen receptor predominance in human ovarian carcinoma. J Steroid Biochem 1987; 26(3):393
11. Carroll RS, et al. Androgen receptor expression in meningiomas. J Neurosurg 1995; 82:453

Explanation of symbols

 Catalogue number	 Temperature limitation	 In vitro diagnostic medical device
 Manufacturer	 Batch code	 Contains sufficient for <n> tests
 Use by	 Consult instructions for use	 Authorized representative in the European Community

Agilent Technologies, Inc.
5301 Stevens Creek Blvd.
Santa Clara, CA 95051
United States

Tel. +44 161 492 7050
www.agilent.com

TX02438/01

Revision 2020.07

Monoclonal Mouse
Anti-Human
Prostate-Specific Antigen
Clone ER-PR8
Kodas M0750

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Literatūra

1. Nap M, van der Kwast TM. Immunohistochemical characterization of 53 monoclonal antibodies to prostate-specific antigen. Tumor Biol 1999 (suppl. 1):20:38-42.
2. Papadopoulos I, Sivridis E, Giatromanolaki A, Koukourakis MI. Tumor angiogenesis is associated with MUC1 overexpression and loss of prostate-specific antigen expression in prostate cancer. Clin Cancer Res 2001;7:1533-8.
3. Polascik TJ, Oesterling JE, Partin AW. Prostate specific antigen: a decade of discovery – what we have learned and where we are going [Review Article]. J Urol 1999;162:293-306.
4. Stenman U-H, Leinonen J, Zhang W-M, Finne P. Prostate-specific antigen. Cancer Biol:1999;9:83-93.
5. Wang TJ, Linton HJ, Sokoloff RL, Grauer LS, Rittenhouse HG, Wolfert RL. Western blotting analysis of antibodies to prostate-specific antigen: specificities for prostate-specific antigen and prostate-specific antigen fragments. Tumor Biol 1999;20 (suppl 1):79-85.
6. Stenman U-H, Paus E, Allard WJ, Andersson I, Andrés C, Barnett TR, et al. Summary report of the TD-3 workshop: characterization of 83 antibodies against prostate-specific antigen. Tumor Biol 1999;20 (suppl 1): 1-12.
7. Riesenberger R, Oberneder R, Kriegmair M, Epp M, Bitzer U, Hofstetter A, et al. Immunocytochemical double staining of cytokeratin and prostate specific antigen in individual prostatic tumour cells. Histochemistry 1993;99:61-6.
8. Alanen KA, Kuopio T, Koskinen PJ, Nevalainen TJ. Immunohistochemical labelling for prostate specific antigen in non-prostatic tissues. Path Res Pract 1996;192:233-7.

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Peržiūra 2020.11

Monoclonal Mouse
Anti-Human
CD45, Leucocyte Common Antigen
Clones 2B11 + PD7/26

Code M0701

27 p.d.

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen, Clones 2B11 + PD7/26, is intended for use in immunohistochemistry (IHC). The antibody labels CD45 in both normal and neoplastic cells and is a useful aid for classifying tumor cells of lymphoid origin (1-3). Differential classification is added by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using non-immunological histochemical stains.
Synonyms for antigen	T200, Ly-5.
Summary and explanation	CD45 is a transmembrane glycoprotein expressed on most nucleated cells of hematopoietic origin. CD45, encoded by a single gene mapped to chromosome 1, has various isoforms based on differential splicing of exons 4, 5 and 6. On human leucocytes, five different isoforms of CD45, named ABC, AB, BC, B and 0, have been identified. These isoforms are recognized by CD45RA, CD45RB, CD45RC and CD45R0 antibodies. The isoforms range in Mr from 180 000 to 220 000. All the CD45 isoforms share the same intracellular segment, which has been shown to have tyrosine phosphatase activity. Various leucocytes express characteristic CD45 isoforms, thus T cells express CD45 isoforms corresponding to their development and activation, B cells predominantly express the ABC isoform, and monocytes and dendritic cells predominantly express the B and 0 isoforms. Granulocytes principally express only the B and 0 isoforms (4). Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage, Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN₃. <u>Clone:</u> 2B11 (1) and PD7/26 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
Immunogen	2B11: Isolated neoplastic cells from a case of T-cell lymphoma/leukemia (1). PD7/26: Human peripheral blood lymphocytes maintained in T-cell growth factor (1).
Specificity	Anti-CD45 is a mixture of two monoclonal antibodies, clones 2B11 and PD7/26, directed against different epitopes. Clone 2B11 was clustered as anti-CD45 at the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens, held in Oxford in 1986 and reacts with all the known isotypes of the CD45 family (5). Clone PD7/26 was clustered as anti-CD45RB at the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, held in Boston in 1993 (6).
Precautions	- For in vitro diagnostic use. - For professional users. - This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. - As with any product derived from biological sources, proper handling procedures should be used. - Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. - Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, B5 (2), or Bouin's (1). Pretreatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed, frozen sections (1). The user must validate the staining procedure.
Staining procedure	These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms. <u>Dilution:</u> Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen, Code M0701, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. <u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

Visualization: Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Product-specific limitations	Labeling of the surface membrane of mammary ductal cells in specimens of fibrocystic disease and fibroadenoma has been reported for the antibody. This expression, however, is not considered to present an assessment problem that might lead to misclassification (7).
Staining interpretation	Cells labeled by the antibody predominantly display staining of the cell membrane, but cytoplasmic staining may also occur.
Performance characteristics	<u>Normal tissues:</u> In tonsil, the antibody labels germinal centres, follicular mantle zones, and interfollicular regions (1). In spleen, white pulp and lymphoid cells of red pulp are labeled, as also thymic lymphocytes, bone marrow lymphoid cells, mast cells, cells of probable monocytic derivation, and occasional plasma cells. Variable labeling of immunoblasts, epithelioid histiocytes, sinus histiocytes and plasma cells has been reported. Myeloid cells, erythroid cells, megakaryocytes, Langerhans cells in skin, epithelium, and connective tissue are not labeled by the antibody (2). <u>Abnormal tissues:</u> In non-Hodgkin's lymphoma, the neoplastic cells were labeled by the antibody in 40/40 (100%) of cases (1). In another study (2) the figure was 74/80 (93%). A third study (8) showed that 52/52 (100%) low grade B-cell lymphomas, 99/108 (92%) high grade B-cell lymphomas, and 41/44 (93%) T-cell lymphomas were labeled with the antibody. Altogether 162/162 (100%) non-lymphoid neoplasms were not labeled with the antibody, including small cell anaplastic carcinomas, amelanotic melanomas, alveolar rhabdomyosarcomas, Ewing's sarcoma, and germ cell tumors (1, 2).

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References/ Bibliographie/ Literatur

- Warnke RA, Gatter KC, Falini B, Hildreth P, Woolston R-E, Pulford K, et al. Diagnosis of human lymphoma with monoclonal antileukocyte antibodies. N Engl J Med 1983;309:1275-81.
- Kurtin PJ, Pinkus GS. Leukocyte common antigen - a diagnostic discriminant between hematopoietic and nonhematopoietic neoplasms in paraffin sections using monoclonal antibodies: Correlation with immunologic studies and ultrastructural localization. Hum Pathol 1985;16:353.
- Michie SA, Spagnolo DV, Dunn KA, Warnke RA, Rouse RV. A panel approach to the evaluation of the sensitivity and specificity of antibodies for the diagnosis of routinely processed histologically undifferentiated human neoplasms. Am J Clin Pathol 1987;88:457-62.
- Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502.
- Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leukocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
- Morimoto C. T18. CD45 cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 386-9.
- Herman GE and Elfont E. Aberrant CD45 (leukocyte common antigen) staining of non-malignant breast lesions in zinc formalin fixed tissue. J Histotechnol 1993;16:151-3.
- Hall PA, d'Ardenne AJ, Stansfeld AG. Paraffin section immunohistochemistry. I. Non-Hodgkin's lymphoma. Histopathol 1988;13:149-60.

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Agilent Technologies Singapore (International) Pte Ltd.