**Baxter**

59.1 poz.

# PRISMAFLEX System

## M60/M100/M150

Hemofilter Sets

Making possible personal.

# PRISMAFLEX System M60/M100/M150 Sets

## Preconnected Hemofilter Sets with the AN 69 Dialyzer Membrane for CRRT

The PRISMAFLEX Set is indicated for use only with the PRISMAFLEX Control Unit in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered via the PRISMAFLEX Set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

The PRISMAFLEX M100 & M150 Sets should be restricted to patients with a body weight greater than 30 kg (66 lb). The pediatric use of the PRISMAFLEX M60 Set should be restricted to children with a body weight greater than 11 kg (24 lb) with respect to extracorporeal blood volume.

The M150, M100 and M60 hemofilter sets are designed for prescription flexibility on the PRISMAFLEX System. Both set configurations are created with the AN 69 Dialyzer Membrane and can be used to perform all CRRT therapies (SCUF, CVVH, CVVHD and CVVHDF).

### General data

	PRISMAFLEX M60 Set	PRISMAFLEX M100 Set	PRISMAFLEX M150 Set
Weight	780 g	800 g	860 g
Overall dimensions	27 x 22 x 9 cm	27 x 22 x 9 cm	27 x 22 x 9 cm
Blood volume in set ± 10 %	93 ml	152 ml	189 ml
Minimal patient weight	11 Kg	30 Kg	30 Kg

#### Materials

- AN 69 HF hollow fiber: Acrylonitrile and sodium methallyl sulfonate copolymer
- Filter housing and headers: Polycarbonate
- Filter potting compound: Polyurethane
- Tubing material: plasticized polyvinyl chloride (PVC)
- Cartridge: PETG
- Sterilization mode: EtO (ethylene oxide)
- Not made with natural rubber latex
- Contains DEHP, which may pose risks to infant development

#### Filter operating specifications

- Maximum TMP\* (mmHg/kPa): 450/60
- Maximum blood pressure (mmHg/kPa): 500/66.6
- Minimum blood flow rate:  
PRISMAFLEX M60 Set: 50 ml/min  
PRISMAFLEX M100 Set: 75 ml/min  
PRISMAFLEX M150 Set: 100 ml/min

#### Filter data

- Nominal physical characteristics
  - Effective surface area:  
PRISMAFLEX M60 Set: 0.6 m<sup>2</sup>  
PRISMAFLEX M100 Set: 0.9 m<sup>2</sup>  
PRISMAFLEX M150 Set: 1.5 m<sup>2</sup>
  - Fiber internal diameter (wet): 240 µm
  - Fiber wall thickness: 50 µm

#### Sieving coefficient

- (bovine plasma, Pc 60 g/l, T = 37°C)  
Urea = 1, Creatinine = 1, Vitamin B12 = 1,  
Inulin = 0.95, Myoglobin = 0.55, Albumin <0.01


## In vitro performances

- Blood priming volume in filter  $\pm 10\%$ :  
 PRISMAFLEX M60 Set: 42 ml  
 PRISMAFLEX M100 Set: 66 ml  
 PRISMAFLEX M150 Set: 105 ml
- Blood pressure drop (in post dilution mode)  
 (bovine blood, Hematocrit 32%,  $P_c^{***} = 60$  g/l,  $T = 37^\circ\text{C}$ )

	PRISMAFLEX M60 Set	PRISMAFLEX M100 Set	PRISMAFLEX M150 Set
$Q_B^{**} = 100$ ml/min - $Q_{UF}^{****} = 1$ l/h	47 mmHg	31 mmHg	20 mmHg
$Q_B^{**} = 180$ ml/min - $Q_{UF}^{****} = 2$ l/h	91 mmHg	60 mmHg	—
$Q_B^{**} = 300$ ml/min - $Q_{UF}^{****} = 2$ l/h	—	73 mmHg	51 mmHg
$Q_B^{**} = 400$ ml/min - $Q_{UF}^{****} = 2$ l/h	—	105 mmHg	64 mmHg
$Q_B^{**} = 450$ ml/min - $Q_{UF}^{****} = 2$ l/h	—	—	70 mmHg

## CVWH Performances

- “In vitro” ultrafiltration with blood (in post-dilution)  
 (values  $\pm 15\%$ ) (Continuous venovenous hemofiltration)  
 (Bovine blood at  $37^\circ\text{C}$ , Hematocrit 32%,  $P_c^{***} = 60$  g/l).

## CVWHD Clearances

- Clearances versus inlet dialysate flow rate  
 (Continuous venovenous hemodialysis) (Saline,  $T = 37^\circ\text{C}$ ).

	PRISMAFLEX M60 Set $Q_B^{**} = 100$ ml/min $Q_{UF}^{****} = 0$ ml/min			PRISMAFLEX M100 Set $Q_B^{**} = 150$ ml/min $Q_{UF}^{****} = 0$ ml/min				PRISMAFLEX M150 Set $Q_B^{**} = 200$ ml/min $Q_{UF}^{****} = 0$ ml/min			
QD l/h	1	2.5	4	1	2.5	4	8	1	2.5	4	8
ml/min	17	42	67	17	42	67	133	17	42	67	133
Urea ( $\pm 10\%$ )	17	39	54	17	41	63	95	17	42	66	117
Vitamin B12 ( $\pm 20\%$ )	14	23	28	16	30	37	45	17	37	49	64
Inulin ( $\pm 20\%$ )	12	17	19	14	23	26	30	16	31	37	45

\*Transmembrane pressure.

\*\*Access blood flow rate.

\*\*\*Protein concentration.

\*\*\*\*Ultrafiltration flow rate (1).

(1)The ultrafiltration flow rate is the “patient fluid removal flow rate + replacement flow rate + pre-blood-pump flow rate.”

## Ordering information

	Factory ID	Code No.	No. Units/Box
PRISMAFLEX M60 Set	8353510	106696	4
PRISMAFLEX M100 Set	8353520	106697	4
PRISMAFLEX M150 Set	8353584	109990	4
SP414 - 5-liter bag	6032957	106690	50
SP418 - 9-liter bag	6033765	107650	30

**Contact your local Baxter sales representative to learn more about the product.**  
**For Customer Support, call 800-525-2623.**

---

**Rx Only.** For the safe and proper use of the devices mentioned herein, please refer to the *Instructions for Use*.

Baxter, AN 69 and Prismaflex are trademarks of Baxter International Inc., or its subsidiaries.

[www.baxter.com](http://www.baxter.com)

USMP/MG120/14-0003(2) 1000 01/16

**Baxter**

**Baxter**

59.2 poz.

## Prismaflex sets with the AN 69 ST membrane for CRRT

Prismaflex ST60, ST100, ST150 sets

POWERED BY  
**prismaflex**

Making possible personal.

# Prismaflex ST60 Set / ST100 Set / ST150 Set

The **Prismaflex** set is indicated for use only with the **Prismaflex** Control Unit in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both. This set is intended for use in the following veno-venous therapies: SCUF, CVVH, CVVHD, CVVHDF.

## General Data

	Prismaflex ST60 set	Prismaflex ST100 set	Prismaflex ST150 set
Weight	795 g	830 g	905 g
Overall dimensions	27 x 22 x 9 cm	27 x 22 x 9 cm	27 x 22 x 9 cm
Blood volume in set ± 10 %	93 ml	152 ml	189 ml
Minimal patient weight	11 kg	30 kg	30 kg

## Materials

- AN 69 ST hollow fiber:
  - Acrylonitrile and sodium methallyl sulfonate copolymer
  - Surface treatment agent: Polyethylene imine

Filter housing and headers: Polycarbonate

Filter potting compound: Polyurethane

Tubing material: Plasticized polyvinyl chloride (PVC)

Cartridge: PETG

Sterilization mode: EtO (ethylene oxide)

## Filter operating specifications

Maximum TMP\* (mmHg/kPa) 450/60

Maximum blood pressure (mmHg/kPa) 500/66.6

Range of blood flow rate:

**Prismaflex** ST60 set: 50–180 ml/min

**Prismaflex** ST100 set: 75–400 ml/min

**Prismaflex** ST150 set: 100–450 ml/min

## Filter data

Nominal physical characteristics

Effective surface area:

**Prismaflex** ST60 set: 0.6 m<sup>2</sup>

**Prismaflex** ST100 set: 1 m<sup>2</sup>

**Prismaflex** ST150 set: 1.5 m<sup>2</sup>

Fiber internal diameter (wet) 240 µm

Fiber wall thickness 50 µm

## In Vitro Performances

Blood priming volume in filter ±10%

**Prismaflex** ST60 set: 44 ml

**Prismaflex** ST100 set: 69 ml

**Prismaflex** ST150 set: 105 ml

Blood pressure drop (in post dilution mode)

(Bovine blood, Hematocrit 32%, Cp\*\*\* = 60 g/l, T = 37°C)


	Prismaflex ST60 set	Prismaflex ST100 set	Prismaflex ST150 set
Q <sub>B</sub> ** = 100 ml/min - Q <sub>UF</sub> **** = 1l/h	46 mmHg	31 mmHg	20 mmHg
Q <sub>B</sub> = 180 ml/min - Q <sub>UF</sub> = 2l/h	84 mmHg	55 mmHg	—
Q <sub>B</sub> = 300 ml/min - Q <sub>UF</sub> = 2l/h	—	79 mmHg	51 mmHg
Q <sub>B</sub> = 400 ml/min - Q <sub>UF</sub> = 2l/h	—	100 mmHg	64 mmHg

## Sieving coefficient

Q<sub>B</sub>\*\* = 100 ml/min, Q<sub>UF</sub>\*\*\*\* = 20 ml/min  
(Bovine plasma, Cp 60 g/l, T = 37°C)  
Urea = 1, Creatinine = 1, Vitamin B<sub>12</sub> = 1, Inulin = 0.96, Myoglobin = 0.58, Albumin <0.01

## CVVH performances

"In vitro" ultrafiltration with blood (in post-dilution)  
(values ±20%) (Continuous venovenous hemofiltration)  
(Bovine blood at 37°C, Hematocrit 32%, Cp\*\*\* 60 g/l).

## CVVHD clearances

Clearances versus inlet dialysate flow rate  
(Continuous venovenous hemodialysis) (Saline, T = 37°C).

	Prismaflex ST60 set Q <sub>B</sub> ** = 100 ml/min Q <sub>UF</sub> **** = 0 ml/min			Prismaflex ST100 set Q <sub>B</sub> ** = 150 ml/min Q <sub>UF</sub> **** = 0 ml/min			Prismaflex ST150 set Q <sub>B</sub> ** = 200 ml/min Q <sub>UF</sub> **** = 0 ml/min				
	1	2.5	4	1	2.5	4	1	2.5	4		
Q <sub>D</sub> /h ml/min	17	42	67	17	42	67	133	17	42	67	133
Urea (±10%)	17	40	56	17	41	63	97	17	42	66	117
Vitamin B <sub>12</sub> (±20%)	15	26	30	16	32	41	50	17	38	51	68
Inulin (±20%)	13	19	22	15	26	30	35	16	33	40	49

\*Transmembrane pressure.

\*\*Access blood flow rate.

\*\*\*Protein concentration.

\*\*\*\*Ultrafiltration flow rate<sup>(1)</sup>.

<sup>(1)</sup>The ultrafiltration flow rate is the "patient fluid removal flow rate + replacement flow rate + pre-blood-pump flow rate".

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59.3 poz.

## **OXIRIS Set**

The OXIRIS Set has been Authorized by the FDA for the indication to treat patients with COVID-19 infection.

The OXIRIS Set is indicated for use only with the **PRISMAFLEX** control unit or with the **PRISMAX** control unit.

# OXIRIS SET

## OXIRIS SET GENERAL DATA

Weight	890 g
Overall Dimensions	27 x 22 x 9 cm
Blood volume in set ± 10 %	193 ml
Minimal patient weight	30 kg

### Materials

OXIRIS hollow fiber: Acrylonitrile and sodium methallyl sulfonate copolymer + Polyethylenimine (surface treatment agent) + heparin grafted [4500+/-1500 IU/m<sup>2</sup>]

Filter housing and headers: Polycarbonate

Filter potting compound: Polyurethane

Tubing material: Plasticized polyvinyl chloride (PVC)

Cartridge: PETG

Sterilization mode: EtO (ethylene oxide)

### Filter operating specifications

Maximum TMP* (mmHg/kPa)	450/60
Maximum blood pressure (mmHg/kPa)	500/66.6
Range of blood flow rate	100-450 ml/min

### Filter data

Nominal physical characteristics:

Effective surface area	1.5 m <sup>2</sup>
Fiber internal diameter (wet)	240 µm
Fiber wall thickness	50 µm

## IN VITRO PERFORMANCES

### CVVHD clearances

Clearances versus inlet dialysate flow rate  
(Continuous veno-venous hemodialysis) (Saline, T = 37°C).

QD l/h ml/min	OXIRIS set Q <sub>B</sub> ** = 200 ml/min Q <sub>UF</sub> **** = 0 ml/min			
	1 17	2.5 42	4 67	8 133
Urea (±10%)	17	42	66	117
Vitamin B <sub>12</sub> (±20%)	17	38	51	68
Inulin (±20%)	16	33	40	49

\*Transmembrane pressure.

\*\*Access blood flow rate.

\*\*\*Protein concentration.

\*\*\*\*Ultrafiltration flow rate<sup>(1)</sup>.

<sup>(1)</sup>The ultrafiltration flow rate is the "patient fluid removal flow rate + replacement flow rate + pre-blood-pump flow rate".


## CVVH performances<sup>1</sup>

"In vitro" ultrafiltration with blood (values ±15%)

(Continuous veno-venous hemofiltration)

(Bovine blood at 37°C, Hematocrit 32%, Cp\*\*\* 60 g/l).

Ultrafiltration is controlled by the PRISMAFLEX system and is independent of the ultrafiltration coefficient (KUF)

### Sieving coefficient

(Bovine plasma, Cp 60 g/l, T = 37°C)

Q<sub>B</sub> = 100 ml/min, Q<sub>UF</sub> = 20 ml/min

Urea	1
Vitamin B <sub>12</sub>	1
Inulin	0.96
(Human plasma, Cp 60 g/l, T = 37°C)	
Myoglobin	0.70
Albumin	<0.0045

### Cytokine adsorption

Cytokine adsorption removal rate [%]<sup>(2)</sup>

(human plasma, Cp 60 g/l, 37°C)

Q<sub>B</sub> = 150 ml/min, Q<sub>UF</sub> = 0 ml/min

IL-10 (± 10%)	96
IL-6 (± 10%)	84
HMGB-1 (± 10%)	94
TNF-α (± 30%)	82

<sup>(2)</sup>Removal Rate expressed at t=120 min with a theoretical initial IL-10, IL-6, HMGB-1 and TNF-α respective concentration of 500 pg/ml, 1500 pg/ml, 30 ng/ml and 250 pg/ml.

### Endotoxin adsorption

Lipopolysaccharide adsorption removal rate [%]<sup>(3)</sup>

(human plasma, Cp 60 g/l, 37°C)

Q<sub>B</sub> = 150 ml/min, Q<sub>UF</sub> = 0 ml/min

LPS (± 20%)	75
-------------	----

<sup>(3)</sup>Removal Rate expressed at t=120 min with an initial LPS concentration after stabilization of 50±10 EU/ml

Cp: Protein concentration

RR: removal rate

IL-10: Interleukin-10

IL-6: Interleukin-6

HMGB-1: High-mobility group box 1

TNF-α: Tumor necrosis factor - α

LPS: Lipopolysaccharide

## ORDERING INFORMATION

	Code N°	N° units/box
OXIRIS S set	955503	4

1. Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. Adsorption removal rate obtained in vitro are likely to differ from in vivo results. Adsorption characteristics change with the duration of observation.

## Emergency Use Authorization for the United States

The OXIRIS Set has been Authorized by the FDA for the indication to treat patients with COVID-19 infection.

The OXIRIS Set has been authorized by FDA under EUA200164.

The OXIRIS Set is Authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of the OXIRIS set under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

## Intended Use for Patients with COVID-19

The OXIRIS Set is indicated for use only with the PRISMAFLEX control unit or with the PRISMAX control unit.

It is intended to treat patients 18 years of age or older with confirmed COVID-19 admitted to the ICU with confirmed or imminent respiratory failure in need of blood purification, including continuous renal replacement therapy, to reduce pro-inflammatory cytokine levels, who have any one of the following conditions:

- Early acute lung injury (ALI)/early acute respiratory distress syndrome (ARDS);
- Severe disease, such as:
  - dyspnea,
  - respiratory frequency  $\geq 30$ /min,
  - blood oxygen saturation  $\leq 93\%$ ,
  - partial pressure of arterial oxygen to fraction of inspired oxygen ratio  $< 300$ , and/or
  - lung infiltrates  $>50\%$  within 24 to 48 hours; or
- Life-threatening disease, defined as:
  - respiratory failure,
  - septic shock, and/or
  - multiple organ dysfunction or failure

The use of OXIRIS is contraindicated (as mentioned in the IFU) when:

- Patients present a known allergy to heparin or have type II thrombocytopenia caused by heparin (HIT Syndrome type II)
- A drug used simultaneously with OXIRIS is contraindicated per its Instructions for use

Relative contraindications (individual risk/benefit to be determined by treating physician) for the use of OXIRIS include:

- The inability to establish vascular access to safely perform CRRT/hemoperfusion (SCUF; CVVH; CVVHD; CVVHDF)
- Severe hemodynamic instability
- Known hypersensitivity to any component of the OXIRIS Set

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the OXIRIS Set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

**Rx Only.** For safe and proper use of the devices mentioned herein, please refer to the Instructions for Use.

**For more information, please contact your Baxter representative.**

[www.baxter.com](http://www.baxter.com)

Baxter HealthCare, One Baxter Parkway, Deerfield, IL 60015

Baxter, Oxiris, Prismaflex and Prismax are trademarks of Baxter International Inc. or its subsidiaries.  
USMP/MG146/20-0001 04/20

59.4 poz.

## Prismaflex HF20 set

Prismaflex set for CRRT in low body weight patients.  
The Prismaflex HF20 set should be used only for patients with a body weight above 8kg (18lbs).

POWERED BY  
**prismaflex**

# Prismaflex HF20 set

The **Prismaflex** set is indicated for use only with the **Prismaflex** control unit equipped with software version 4.00 or later in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

## General Data

Prismaflex HF20 Set	
Weight	664 g
Overall dimensions	27 x 22 x 9 cm
Blood volume in set (±10%)	60 ml
Required volume for priming	500 ml

### Materials

PAES hollow fiber: Polyarylethersulfone  
 Housing and headers: Polycarbonate  
 Potting compound: Polyurethane  
 Tubing material: Plasticized polyvinyl chloride (PVC)  
 Cartridge: PETG

### Filter operating specifications

Maximum TMP* (mmHg/kPa)	500/67
Maximum blood pressure (mmHg/kPa)	500/67
Minimum blood flow rate:	20 ml/min
Maximum blood flow rate:	100 ml/min

\*Transmembrane pressure

### Filter data

Nominal physical characteristics	
Effective surface area	0.2m <sup>2</sup>
Fiber internal diameter (wet)	215 µm
Fiber wall thickness	50 µm
Blood priming volume (filter only)	17 ml

## In Vitro Performance Specifications<sup>†</sup>

### Blood pressure drop (post-dilution)

(Bovine blood, Q<sub>UF</sub><sup>\*\*\*</sup> = 0 l/h Hct<sup>\*\*\*\*</sup> 32%, Cp<sup>\*\*\*\*\*</sup> 60 g/l, 37°C)

Q <sub>B</sub> <sup>**</sup> = ml/min	20	50	100
Blood pressure drop (post dilution) mmHg	25	37	59

\*\*Blood flow rate

\*\*\*Ultrafiltration flow rate

\*\*\*\*Hematocrit

\*\*\*\*\*Protein Concentration

<sup>†</sup> Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions.


### Sieving coefficient

(Bovine plasma, Cp 60 g/l, T = 37°C)  
 Q<sub>B</sub><sup>\*\*</sup> = 50 ml/min, Q<sub>UF</sub><sup>\*\*\*</sup> = 10 ml/min  
 Solute = SC, Urea = 1, Vitamin B12 = 1,  
 Inulin = 0.92, Albumin <0.01


## CVVHD clearances

(Continuous veno-venous hemodialysis). Clearances versus inlet dialysate flow rate (37°C) variable QD, variable QB, Q<sub>UF</sub> = 0 ml/min


### Clearance data at QB = 20 ml/min

### Clearance data at QB = 50 ml/min


### Clearance data at QB = 100 ml/min

## CVVH performances

### Ultrafiltration with blood

In vitro ultrafiltration with blood (values ± 20%). (Bovine blood at 37°C, Hct 32%, Protein concentration 60 g/l). Ultrafiltration is controlled by the **Prismaflex** system and is independent of the ultrafiltration coefficient (KUF).

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59.5 poz.

# Prismaflex TPE 1000 and TPE 2000 Sets

Pre-connected sets for  
plasmafiltration

POWERED BY  
POWERED BY  
**prismaflex**

Making possible personal.

# Prismaflex TPE 1000 and TPE 2000 Sets

## General Data

	Prismaflex TPE 1000 set	Prismaflex TPE 2000 set
Weight	680 g	760 g
Blood volume in set ± 10 %	71 ml	125 ml
Blood volume in plasmafilter	23 ml ±10%	41 ml ±10%
Minimal patient weight	9 kg	adult

Materials	
Plasmafilter hollow fiber: Polypropylene	
Filter housing and headers: Polycarbonate	
Filter potting compound: Polyurethane	
Tubing material: Plasticized polyvinyl chloride (PVC)	
Cartridge: PETG	

Sterilization mode: EtO (ethylene oxide)	
Filter operating specifications	
Maximum blood flow rate: Prismaflex TPE 1000 set: 180 ml/min Prismaflex TPE 2000 set: 400 ml/min	
Minimum blood flow rate: Prismaflex TPE 1000 set: 50 ml/min Prismaflex TPE 2000 set: 100 ml/min	
Maximum TMPa* (mmHg/kPa):	

	Prismaflex TPE 1000 set	Prismaflex TPE 2000 set
Q <sub>B</sub> ** = 50 ml/min	100 mmHg - 13.3 kPa	–
Q <sub>B</sub> ** = 100 ml/min	140 mmHg - 18.7 kPa	120 mmHg - 16 kPa
Q <sub>B</sub> ** = 180 ml/min	190 mmHg - 25.3 kPa	–
Q <sub>B</sub> ** = 200 ml/min	–	171 mmHg - 22.8 kPa
Q <sub>B</sub> ** = 250 ml/min	–	193 mmHg - 25.7 kPa
Q <sub>B</sub> ** = 400 ml/min	–	246 mmHg - 32.8 kPa

Filter data	
Nominal physical characteristics	
Effective surface area: Prismaflex TPE 1000 set: 0.15 m <sup>2</sup> Prismaflex TPE 2000 set: 0.35 m <sup>2</sup>	
Fiber internal diameter (wet):	330 µm
Fiber wall thickness:	150 µm

## Ordering Information

	Code N°	N° units/box
Prismaflex TPE 1000 set	107143	4
Prismaflex TPE 1000 set CKT	115314	4
Prismaflex TPE 2000 set	107144	4
Prismaflex TPE 2000 set CKT	115315	4

## In Vitro Performances



Blood pressure drop ±30%  
(Bovine blood, Htc\*\*\*32%, Protein concentration 60 g/l, 37°C)

	Prismaflex TPE 1000 set	Prismaflex TPE 2000 set
Q <sub>B</sub> ** = 50 ml/min	14 mmHg - 1.9 kPa	–
Q <sub>B</sub> ** = 100 ml/min	29 mmHg - 3.9 kPa	25 mmHg - 3.3 kPa
Q <sub>B</sub> ** = 180 ml/min	52 mmHg - 6.9 kPa	–
Q <sub>B</sub> ** = 200 ml/min	–	51 mmHg - 6.8 kPa
Q <sub>B</sub> ** = 250 ml/min	–	64 mmHg - 8.5 kPa
Q <sub>B</sub> ** = 400 ml/min	–	101 mmHg - 13.5 kPa

\*Arterial transmbrane pressure  
\*\*Access blood flow rate  
\*\*\*Hematocrit

Sieving coefficient  
(In vivo data, 19 treatments)  
Albumin = 0.97, IgG = 1, Apolipoprotein B = 0.95, IgM = 0.92

Plasma filtration with blood  
In vitro, values ±20%  
(Bovine blood at 37°C, Htc 32%, Protein concentration 60 g/l)  
Plasma filtration is controlled by the Prismaflex system.

**Sweden:**  
Baxter Medical AB  
Box 63  
164 94 Kista  
T: 08-632 64 00  
www.baxter.se

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# Calcium Line for Prismaflex<sup>®</sup> – CA 250

Prismaflex<sup>®</sup> is a trademark of Gambro Lundia AB

Instruction for Use	GB	3
Brugsanvisning	DK	5
Gebruiksaanwijzing	NL	7
Käyttöohjeet	FI	9
Mode d'emploi	FR	11
Gebrauchsanleitung	DE	13
Istruzioni per l'uso	IT	15
Bruksanvisning	NO	17
Instrukcja użytkowania	PL	19
Instruções de utilização	PT	21
Инструкция по применению	RU	23
Instrucciones de uso	ES	25
Bruksanvisning	SE	27
Kullanma Talimatı	TR	29

## INSTRUCTION FOR USE

### Calcium Line for Prismaflex® – CA 250

Prismaflex® is a trademark of Gambro Lunds AB

**Caution!**  
Read the instructions carefully before using this product.

#### INTENDED USE

The Calcium Line for Prismaflex is a medical device to be used only with the Prismaflex Control Unit software version 5.0 or later for performing the 'Citrate-calcium via Prismaflex syringe pump' anticoagulation method. Carefully read the Prismaflex Operator's Manual and on-line screen instructions. All treatments administered via the Calcium Line for Prismaflex must be prescribed by a physician.

#### CONTRAINDICATIONS AND ADVERSE REACTIONS

There are no known contraindications or adverse reactions to the use of the Calcium Line for Prismaflex if used as indicated in the intended use section above.

#### CAUTIONS AND WARNINGS

-  The lines must be stored in a dry place, between 5°C (41°F) and 30°C (86°F).
-  Do not store in direct sunlight.
-  Expiration date of the product is shown on the packaging.
-  Do not use an individual unit if package is damaged.
- In order to reduce the risk of disease transmission, the use of protective clothing (gloves, glasses, masks, etc.) by nurses, doctors, and other medical staff is recommended when handling the Calcium Line for Prismaflex.
- Before using the Calcium Line for Prismaflex check that all caps are in place.
- The Calcium Line for Prismaflex must be used as soon as the packaging and the protective caps have been removed.
- Before starting, and also during the treatment using Calcium Line for Prismaflex, ensure that the check valve is in place and all connections are secure.
- The minimum temperature for use is 20°C (68°F).
- The Calcium Line for Prismaflex is sterile and non-pyrogenic. Use an aseptic technique when installing the line and throughout the treatment.
- Make sure the Calcium Line for Prismaflex and the patient's vascular access are not kinked.
- Verify safe operation of the system comprising the Prismaflex machine, any accessories and the Calcium Line for Prismaflex during set up.
- The safe connection of the Calcium Line for Prismaflex shall be performed by skilled personnel only or under the supervision of the medical officer in charge.
- Watch carefully for leaks during priming and use.
- All connections must be checked carefully throughout the treatment to prevent any misconnection, leakage or disconnection issues that could result in potential patient injuries such as blood loss, air embolism or fluid imbalance.
- During the treatment, it should be regularly checked that the line is not kinked or obstructed in order to prevent any risk of electrolytic imbalance for the patient.
- Presence of DEHP in the PVC tube should be taken into special consideration in the treatment of children who have not yet reached puberty, pregnant women and nursing mothers because these patient groups exhibit an increased sensitivity.
- Use only drugs compatible with PVC.
-  The Calcium Line for Prismaflex is intended for single use only. Discard it after use, following local laws and regulations for potentially contaminated equipment.
- The Calcium Line for Prismaflex must not be reused in order to avoid bacteriological contamination and possible performance decrease.
- Sterility and performance of this device is guaranteed by the manufacturer only if it is intact and prepared as recommended for single use only.

**OTHER SYMBOLS USED ON THE PRODUCT LABELING**

	Consult Instructions For Use
	Date of Manufacture
	Batch code
	The material used for the manufacture of this specific medical device contains Di-2-ethyl hexyl phthalate (DEHP)
	Keep dry

	Manufacturer
	Catalogue number
	The Calcium Line for Prismaflex is sterilized by irradiation
	The product doesn't contain latex or its derivatives.
	Fragile, handle with care

**INSTRUCTIONS FOR USE AND RECOMMENDATIONS**

Use the Calcium Line for Prismaflex by following the detailed on-line instructions provided by the Prismaflex Control Unit. Additional information is available in the Prismaflex Operator's Manual.

- a) Unpack and remove the protective caps
- b) Follow online-instructions on the Prismaflex screen for the following steps:
  - Connection of the line to the Prismaflex syringe
  - Automatic priming of the line by the Prismaflex system
  - Connection of the line to the patient
- c) Remove eventual air bubbles before use

Please refer to Prismaflex Operator's Manual for more information on the "Citrate-calcium via Prismaflex syringe pump" anticoagulation method and calcium infusion.

**SPECIFICATIONS**

- (1) Line: The Calcium Line for Prismaflex is compliant with the applicable parts of the DIN EN ISO 8536 standard  
**Pump line ISO 8536-9-SPL- P**  
 • Length: 2500 mm  
 • Inner diameter: 0.59 mm  
 • Storage volume at 40°C (104°F): 0.7 mL
- (2) Luer-lock connectors: Both male Luer and female Luer-lock connectors comply with the following standards: EN 20594-1 and EN 1767.

- (3) Check valve is compliant with the applicable parts of the DIN EN ISO 8536 standard.
- (4) Slide clamp
- (5) "Calcium" tags

The materials in direct or indirect contact with blood are: Acrylonitrile Butadiene Styrene (ABS), high and low density polyethylene (HDPE, LDPE), classed Polyvinylchloride (PVC), Styrene Acrylonitrile (SAN), silicone.

**WARRANTY AND LIMITATION OF LIABILITY**


- a) The manufacturer, Medizintechnik Promedt GmbH, warrants that the lines have been manufactured in accordance with their specifications and in compliance with DIN EN ISO 13485, and other regulatory requirements.
- b) Medizintechnik Promedt GmbH Quality system is working in accordance with DIN EN ISO 13485 as it is assessed by the Notified Body MEDCERT
- c) The manufacturer shall not be held responsible for failure to use the lines with the dialysis machines they have been specifically designed for

- a) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the lines, failure or omission to inspect the lines before use in order to insure that they are in proper condition, or any warranty given by independent distributors or dealers.
- b) The manufacturer is Medizintechnik Promedt GmbH, Kleiner Moorweg 4, 25438 Tornesch, Germany.

Sponsor for Australia: Gambro Pty Ltd, Suite 2, Level 4, 62 Norwest Blvd, Baulkham Hills NSW 2153, Australia.

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- e) Gamintojas yra „Medizintechnik Promedt GmbH“, Kleiner Moorweg 4, 25436 Tornesch, Vokietija.
Rėmėjas Australijoje: „Gambro Pty Ltd“, Suite 2, Level 4, 62 Norwest Blvd, Baulkham Hills NSW 2153, Australija.