
ANIOSYME PRIME

ANIOSYME PRIME

Disinfectant
detergent
for instrumentation

- Formulation combining 2 patented concepts of efficiency
- Prevents soil from drying
- Antimicrobial efficacy from 15 min
- Dense foam: instruments are covered evenly
- Wide compatibility

INDICATIONS

Pre-treatment of medical and surgical instruments under foam : during transport in non-humid mode or during storage of soiled instruments before manual or machine cleaning. The product maintains an adequate level of humidity to avoid blood and other organic materials drying.

CHARACTERISTICS

- Tri-enzymatic complex: proteases, lipases, amylases, combined with surfactants
- Non-corrosive for materials
- No production of aerosols, avoiding any risk of emanation from the product
- Ready-to-use, ergonomic foam gun shape
- pH of product : 8 +/- 0.5
- Helps the removal of stain when rinsing
- Enzymatic stability proved
- High rinsability


READY
TO USE

ANIOSYME PRIME

Disinfectant detergent for instrumentation

INSTRUCTIONS FOR USE

1 Beforehand, spray the bottom of the soaking tray with ANIOSYME PRIME and place instruments collected.

2 Recover with foam all the instruments.

3 Cover the tray with the lid.

4 Rinse the instrument before machine or manual cleaning.

COMPOSITION

Quaternary ammonium carbonate, non-ionic surfactants, enzymatic complex (protease, lipase, amylase), sequestering agents, excipients.

PRECAUTIONS FOR USE

Dangerous – respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products).
Storage : between +5°C and +35°C.
Class IIb medical device (Directive 93/42/EEC as amended).

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Page 4

Summary

Regulation and General Information Brochure	Page	7
Composition	Page	9
Stability and storage conditions	Page	11
Studies and expertises		
- Microbiology	Page	13
- Detergency	Page	17
- Corrosion	Page	25
- Compatibilities	Page	27
Regulation	Page	31

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Regulation and general information

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***ANIOSYME PRIME** is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System*

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***ANIOSYME PRIME** has the CE mark for medical devices in the category class IIb, in accordance with the Directive 93/42/CE*

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***ANIOSYME PRIME** formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments*

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***ANIOSYME PRIME** is labelled in accordance with the European regulation related to the classification and labelling of chemical products*

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***ANIOSYME PRIME** meets our voluntary eco-design commitment included in our ANIOSAFE charter*

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Surfactant agents

Complexe enzymatique (protéase, lipase, amylase)

Enzymatic complex (protease, lipase, amylase)

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Pure product not diluted :

- Storage between +5°C and +35°C.
- Stability:
- Not opened bottle: 3 years from the manufacturing date indicated on the label.
- Spray: preserving of the 3 years if the product is stored in a container closed by the sprayer.

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
Present surfactants distinguish themselves by hydrophilic and hydrophobic properties, which ensures performances against proteic and glucidic soilings on one hand, and greasy soiling (lipidic) on another hand, by solubilisation and/or emulsification.

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The tri-enzyme pre-disinfectant detergent ANIOSYME PRIME, ready to use solution, presents a High Degreasing Cleaning Power against an organic soiling fixed on carriers made of stainless steel (hydrophilic character) from 15 minutes at room temperature

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Proteases, lipases and amylases act in the same way as surfactants, degrading, by hydrolysis, organic compounds not much or not soluble in soluble components with lowest molecular weight.

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
The proteasic activity of the ready to use ANIOSYME PRIME solution is demonstrated after storage during 56 days at +30°C, by putting it into contact during 15 minutes with the human albumin substrate at room temperature.

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Lipasic and amylasic activities of ANIOSYME PRIME are demonstrated identically by thin layer chromatography method.

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Humidity retention on surfaces treated (1)

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The humidity retention of a surface after spraying of ANIOSYME PRIME is followed in doing the 3 below techniques:

- *Measure of the loss of weight of the surface during the drying*
- *Checking of the aspect of surfaces : brightness*
- *Determination of humidity after drying by swabbing surface*

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The studies were led in comparison with a version of the product ANIOSYME PRIME without glycerin (negative test) to prove the interest of emollients for this application.

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
Humidity retention on surfaces treated (2)

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TS 17.2.

Aspect of surfaces after a 72 hour drying

ANIOSYME PRIME sans Glycérine
ANIOSYME PRIME without Glycerin

Surface partiellement humide et brillante
Partially wet and weakly shiny surface

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
Humidity retention on surfaces treated (3)

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
It is demonstrated, by electrochemical study of pitting corrosion according to the NF S 94-402-1 standard (May 2004), that ANIOSYME PRIME does not present pitting corrosive character with stainless steel Z30 Cr13 in the conditions of use recommended by Laboratoires ANIOS.

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Incompatibility

- ◆ Aluminium AG3 (ternissement)
Aluminium AG3 (tarnish)
- ◆ Aluminium AU4G (ternissement)
Aluminium AU4G (tarnish)

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Not concerned by biodegradability notion

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26,2 %¹

*Readily biodegradable organic material*²

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Biodegradability within 28 days or less:

- $\geq 60\%$ according to respirometric methods³ (assessment of CO₂ production or oxygen consumption)
- $\geq 70\%$ according to methods measuring dissolved organic carbon³ (DOC) (assessment of DOC die away)

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Substances which not fulfil the readily biodegradability criteria (but possibly displaying inherent biodegradability criteria⁴), or substances without any biodegradability data.

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ANIOSYME PRIME contains at least 99.9 % of inorganic material and of readily biodegradable organic material.

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Serious eye damage, Category 1 (Eye Dam. 1, H318).*

Environment

This mixture does not present an environmental hazard. No known or foreseeable environmental damage under standard conditions of use.*

* Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.

PERSONAL PROTECTION UNDER THE CONDITIONS OF THE INTENDED USE

Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.

Respiratory protection

In case of insufficient ventilation, with the risk of exceeding the TLVs/TWA**, wear suitable breathing apparatus.

Particularly breathing apparatus, A2P2 type.

Hand protection

When handling this product, wear suitable gloves.

Nitrile, latex or vinyl gloves.

Gloves must be replaced immediately if they show signs of wear and tear.

Eye and face protection

Avoid all contact with eyes.

Make safety goggles with side protection available to the personnel.

Provide eye washes in workshops.

If necessary, the nearest water supply.

** please refer to SDS §8 for TLVs (Threshold Limit Values) /TWA (Time-Weighted Average) of concerned substances present in the product

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The Safety Data Sheet (SDS) of this product is available on our Web site www.anios.com. In addition to visualization and/or impression, it is possible for you to download it.

The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.

In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)

* an access is also possible at the site <http://www.quickfds.com>