- For bolt implantation
- "Gold Tip" visible on CT
- For use in brain tissue

## 70 Microdialysis Bolt Catheter

### Intended Purpose

The Brain Microdialysis Bolt Catheter is intended to enable microdialysis of the extracellular (interstitial) fluid of the brain tissue.

### Minimally Invasive

The sterile single use 70 Bolt Microdialysis Catheter is minimally invasive and designed for implantation in brain tissue through the PTO 2L bolt kit from Raumedic or the IM3 Bolt Kit from Integra. The catheter has a 130 mm long shaft and a membrane length of 10 mm. The shaft has a reinforced stainless steel section measuring 16.5 mm in length. This reinforcement is located where the compression screw of the bolt is to be tightened.

The dialysing membrane of the catheter has been especially developed to achieve optimal diffusing characteristics. This allows a high recovery of substances from the extracellular fluid into the catheter.

### The "Gold tip" makes the catheter visible on CT

The tip of the catheter contains a gold thread. The "Gold tip" is visible on CT-scanning and makes it possible to locate the exact position of the catheter.

### Early detection of local tissue chemistry changes in the brain

The 70 Bolt Microdialysis Catheter offers to collect samples manually in Microvials for analysis of Glucose, Lactate, Pyruvate, Glycerol, Glutamate and Urea in ISCUSflex Microdialysis Analyzer and gives the opportunity for early detection of local tissue changes in the brain.

For research purposes the collected samples can be further analyzed by various techniques such as HPLC.

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
label on single 70 Bolt MD cath

**MD pump kit, brain tissue** **μ dialysis**

for Microdialysis Catheter & Microdialysis Pump

1 x Perfusion Fluid CNS  
1 x 106 Syringe  
1 x Battery

**REF** 8003791

**LOT** T12345

2022-11-23

**STERILE**

(01)07332699001062  
(10)T12345  
(11)221123  
(17)221123


2022-11-23  
M Dialysis AB, Hammarby Fabriksv. 43, SE-120 30 Stockholm, Sweden  
US Office: 73 Princeton St, N Chelmsford, MA 01863, USA

8003791

Label on Pumpkit box

# Microdialysis Pump Kit

## MD Pump Kit, Brain Tissue

This Pump Kit includes the accessories needed for single patient microdialysis monitoring in brain tissue. The kit contains a 106 syringe, a battery and a CNS perfusion fluid ampoule necessary for running the 106 and 107 Microdialysis Pumps. The Pump Kit reduces preparation time and allows new users to cost-effectively become familiar with microdialysis.

NOTE: This product may not be approved by your national regulatory authorities. Please contact your local representative for further information.

For USA customers: CAUTION Investigational device. Limited by Federal (or United States) law to investigational use.