



07089627001V1.0

# CARDIAC POC Troponin T Control **cobas**<sup>®</sup>

Roche CARDIAC POC Troponin T Control

[REF] 07089643 190

2 x → 1.0 mL

[SYSTEM] **cobas h 232**

## English

### Intended use

The Roche CARDIAC POC Troponin T Control is used for quality control of the Roche CARDIAC POC Troponin T test on the **cobas h 232** instrument.

### Summary

The Roche CARDIAC POC Troponin T Control is used for accuracy and precision control. Each Roche CARDIAC POC Troponin T Control kit consists of two bottles of lyophilized controls based on human serum containing an elevated Troponin T concentration.

### Reagents - working solutions

- 2 bottles of lyophilized control serum, for 1.0 mL
- Human serum-based solution containing constituents of human origin
- Active ingredients: bovine troponin T

The concentrations of the components are lot-specific.

The exact target values and concentration ranges of the components are given in the code chip provided.

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

2-methyl-2H-isothiazol-3-one hydrochloride.

EUH208 May produce an allergic reaction.

Product safety labeling primarily follows EU GHS guidance.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV.

The testing methods applied were FDA-approved or cleared in compliance with the European Directive 98/79/EC, Annex II, List A. However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.<sup>1,2</sup>

### Handling

Carefully dissolve the contents of one bottle by adding exactly 1.0 mL of distilled or deionized water and allow to stand closed for 15 minutes to reconstitute. Mix carefully, avoiding foam formation.

### Storage and stability

Store at 2-8 °C.

Stability of the lyophilized control serum at 2-8 °C: up to the stated expiration date.

Stability of components in reconstituted control serum:

- at 15-25 °C: up to 4 hours
- at 2-8 °C: 3 days
- at -20 °C: 8 weeks (can be frozen once in the original vial)

Store the controls tightly capped when not in use.

### Materials provided

- Roche CARDIAC POC Troponin T Control
- 1 code chip

### Materials required (but not provided)

- [REF] 07007302190, Roche CARDIAC POC Troponin T
- [REF] 11622889190, Roche CARDIAC Pipettes, 20 disposable syringes (150 µL)
- [REF] 04901126190, **cobas h 232** instrument, software version ≥ 03.00.02, serial number ≥ KQ0120000

- [REF] 04901142190, **cobas h 232** instrument with scanner, software version ≥ 03.00.02, serial number ≥ KS0210000
- General laboratory equipment
- Distilled or deionized water

### Assay

Dispense the required volume into the sample area of the Roche CARDIAC POC Troponin T test and analyze in the same way as a whole blood sample.

Sample volume: 150 µL

Follow the test instructions in the Method Sheet of the Roche CARDIAC POC Troponin T test.

Frozen or refrigerated reconstituted control material must be brought to room temperature prior to use.

### Target values and ranges

Target value determination procedure: Series with a significant number of determinations run on different **cobas h 232** instruments.

The median is calculated as the target value.

The corresponding control range is based on 95 % confidence level.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

### Calibration

Each Roche CARDIAC POC Troponin T Control kit contains a lot-specific code chip.

The code chip must be read into the instrument once for each new lot. The request to insert this code chip is made via the instrument's display. Correct assignment of code chip and control lot can be checked by comparing the lot number in the display with the imprint on the code chip. The code chip provides the **cobas h 232** instrument with all lot-specific information and target value ranges necessary for evaluation. An error flag appears if the wrong code chip is inserted for the control lot in question.

### References

- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.

For further information, please refer to the appropriate operator's manual for the instrument concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

[CONTENT]	Contents of kit
[SYSTEM]	Analyzers/Instruments on which reagents can be used
→	Volume after reconstitution or mixing
[GTIN]	Global Trade Item Number

ROCHE CARDIAC, COBAS and COBAS H are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics



07089627001V1.0

# CARDIAC POC Troponin T Control cobas®

Roche CARDIAC POC Troponin T Control



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

