

# CHOL2

Cholesterol Gen.2

Order information

cobas®

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
03039773 190	Cholesterol Gen.2 400 tests	System-ID 07 6726 3
10759350 190	Calibrator f.a.s. (12 x 3 mL)	Code 401
10759350 360	Calibrator f.a.s. (12 x 3 mL, for USA)	Code 401
12149435 122	Precinorm U plus (10 x 3 mL)	Code 300
12149435 160	Precinorm U plus (10 x 3 mL, for USA)	Code 300
12149443 122	Precipath U plus (10 x 3 mL)	Code 301
12149443 160	Precipath U plus (10 x 3 mL, for USA)	Code 301
10171743 122	Precinorm U (20 x 5 mL)	Code 300
10171735 122	Precinorm U (4 x 5 mL)	Code 300
10171778 122	Precipath U (20 x 5 mL)	Code 301
10171760 122	Precipath U (4 x 5 mL)	Code 301
10781827 122	Precinorm L (4 x 3 mL)	Code 304
11285874 122	Precipath L (4 x 3 mL)	Code 305
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391
05947626 160	PreciControl ClinChem Multi 1 (4 x 5 mL, for USA)	Code 391
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392
05947774 160	PreciControl ClinChem Multi 2 (4 x 5 mL, for USA)	Code 392
04489357 190	Diluent NaCl 9 % (50 mL)	System-ID 07 6869 3

## English

### System information

For **cobas c** 311/501 analyzers:**CHO2I:** ACN 798: ID/MS Standardization**CHO2A:** ACN 433: Abell/Kendall StandardizationFor **cobas c** 502 analyzer:**CHO2I:** ACN 8798: ID/MS Standardization**CHO2A:** ACN 8433: Abell/Kendall Standardization

### Intended use

In vitro test for the quantitative determination of cholesterol in human serum and plasma on Roche/Hitachi **cobas c** systems.

### Summary

Cholesterol is a steroid with a secondary hydroxyl group in the C3 position. It is synthesized in many types of tissue, but particularly in the liver and intestinal wall. Approximately three quarters of cholesterol is newly synthesized and a quarter originates from dietary intake. Cholesterol assays are used for screening for atherosclerotic risk and in the diagnosis and treatment of disorders involving elevated cholesterol levels as well as lipid and lipoprotein metabolic disorders.

Cholesterol analysis was first reported by Liebermann in 1885 followed by Burchard in 1889. In the Liebermann-Burchard reaction, cholesterol forms a blue-green dye from polymeric unsaturated carbohydrates in an acetic acid/acetic anhydride/concentrated sulfuric acid medium. The Abell and Kendall method is specific for cholesterol, but is technically complex and requires the use of corrosive reagents. In 1974, Roeschla and Allain described the first fully enzymatic method. This method is based on the determination of  $\Delta^4$ -cholestenone after enzymatic cleavage of the cholesterol ester by cholesterol esterase, conversion of cholesterol by cholesterol oxidase, and subsequent measurement by the Trinder reaction of the hydrogen peroxide formed. Optimization of ester cleavage (> 99.5 %) allows standardization using primary and secondary standards and a direct comparison with the CDC and NIST reference methods.<sup>1,2,3,4,5,6,7,8,9</sup>

Nonfasting sample results may be slightly lower than fasting results.<sup>10,11,12</sup>


The Roche cholesterol assay meets the 1992 National Institutes of Health (NIH) goal of less than or equal to 3 % for both precision and bias.<sup>12</sup>

The assay is optionally standardized against Abell/Kendall and isotope dilution/mass spectrometry. The performance claims and data presented here are independent of the standardization.

### Test principle

Enzymatic, colorimetric method.

Cholesterol esters are cleaved by the action of cholesterol esterase to yield free cholesterol and fatty acids. Cholesterol oxidase then catalyzes the oxidation of cholesterol to cholest-4-en-3-one and hydrogen peroxide. In the presence of peroxidase, the hydrogen peroxide formed effects the oxidative coupling of phenol and 4-aminophenazone to form a red quinone-imine dye.



The color intensity of the dye formed is directly proportional to the cholesterol concentration. It is determined by measuring the increase in absorbance.

### Reagents – working solutions

**R1** PIPES buffer: 225 mmol/L, pH 6.8;  $\text{Mg}^{2+}$ : 10 mmol/L; sodium cholate: 0.6 mmol/L; 4-aminophenazone:  $\geq 0.45$  mmol/L; phenol:  $\geq 12.6$  mmol/L; fatty alcohol polyglycol ether: 3 %; cholesterol esterase (*Pseudomonas spec.*):  $\geq 25$   $\mu\text{kat/L}$  ( $\geq 1.5$  U/mL); cholesterol oxidase (*E. coli*):  $\geq 7.5$   $\mu\text{kat/L}$  ( $\geq 0.45$  U/mL); peroxidase (horseradish):  $\geq 12.5$   $\mu\text{kat/L}$  ( $\geq 0.75$  U/mL); stabilizers; preservative

R1 is in position B.



# CHOL2

## Cholesterol Gen.2

**cobas®**

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

### Hazardous components:

3,6,9,12,15,18,21,24,27-nonaioxanonatriacantan-1-ol



Warning

H319 Causes serious eye irritation.

### Prevention:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

### Response:

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

Contact phone: all countries: +49-621-7590, USA: +1-800-428-2336

### Reagent handling

Ready for use

### Storage and stability

#### CHOL2

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 4 weeks

*Diluent NaCl 9 %*Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

### Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable. Serum.

Plasma: Li-heparin and K<sub>2</sub>-EDTA plasmaDo not use citrate, oxalate or fluoride.<sup>13</sup>Fasting and nonfasting samples can be used.<sup>11</sup>

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Stability:<sup>14,15</sup>

7 days at 15-25 °C

7 days at 2-8 °C

3 months at (-15)-(-25) °C

### Materials provided

See "Reagents – working solutions" section for reagents.

### Materials required (but not provided)

See "Order information" section

General laboratory equipment

### Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

### Application for serum and plasma

#### cobas c 311 test definition

Assay type	1-Point		
Reaction time / Assay points	10 / 57		
Wavelength (sub/main)	700/505 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R1	47 µL	93 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 µL	–	–
Decreased	2 µL	15 µL	135 µL
Increased	2 µL	–	–

#### cobas c 501 test definition

Assay type	1-Point		
Reaction time / Assay points	10 / 70		
Wavelength (sub/main)	700/505 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		
R1	47 µL	93 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2 µL	–	–
Decreased	2 µL	15 µL	135 µL
Increased	2 µL	–	–

#### cobas c 502 test definition

Assay type	1-Point		
Reaction time / Assay points	10 / 70		
Wavelength (sub/main)	700/505 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting	Diluent (H <sub>2</sub> O)		



# CHOL2

## Cholesterol Gen.2

R1	47 µL	93 µL	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2 µL	–	–
Decreased	2 µL	15 µL	135 µL
Increased	4 µL	–	–

### Calibration

Calibrators	S1: H <sub>2</sub> O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration <ul style="list-style-type: none"> <li>• after reagent lot change</li> <li>• as required following quality control procedures</li> </ul>

Traceability: This method has been standardized according to Abell/Kendall<sup>12</sup> and also by isotope dilution/mass spectrometry.<sup>16</sup>

### Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

### Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factors:	mmol/L x 38.66 = mg/dL
	mmol/L x 0.3866 = g/L
	mg/dL x 0.0259 = mmol/L

### Limitations – interference

Criterion: Recovery within ± 10 % of initial values at a cholesterol concentration of 5.2 mmol/L (200 mg/dL).

Icterus:<sup>17</sup> No significant interference up to an I index of 16 for conjugated bilirubin and 14 for unconjugated bilirubin (approximate conjugated bilirubin concentration 274 µmol/L or 16 mg/dL; approximate unconjugated bilirubin concentration 239 µmol/L or 14 mg/dL).

Hemolysis:<sup>17</sup> No significant interference up to an H index of 700 (approximate hemoglobin concentration: 435 µmol/L or 700 mg/dL).

Lipemia (Intralipid):<sup>17</sup> No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>18,19</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>20</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### ACTION REQUIRED

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/Multiclean/SCCS or the

NaOHD/SMS/SmpCln1+2/SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c** 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

### Limits and ranges

#### Measuring range

0.1-20.7 mmol/L (3.86-800 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 10.

#### Lower limits of measurement

##### Lower detection limit of the test

0.1 mmol/L (3.86 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

#### Expected values

Clinical interpretation according to the recommendations of the European Atherosclerosis Society:<sup>21</sup>

	mmol/L	mg/dL	Lipid metabolic disorder
Cholesterol	< 5.2	(< 200)	No
Triglycerides	< 2.3	(< 200)	No
Cholesterol	5.2-7.8	(200-300)	Yes, if HDL-cholesterol < 0.9 mmol/L (< 35 mg/dL)
Cholesterol	> 7.8	(> 300)	Yes
Triglycerides	> 2.3	(> 200)	Yes

Recommendations of the NCEP Adult Treatment Panel for the following risk-cutoff thresholds for the US American population:<sup>22</sup>

Desirable cholesterol level	< 5.2 mmol/L	(< 200 mg/dL)
Borderline high cholesterol	5.2-6.2 mmol/L	(200-240 mg/dL)
High cholesterol	≥ 6.2 mmol/L	(≥ 240 mg/dL)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

### Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

#### Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

Repeatability	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Precinorm U	2.29 (88.5)	0.02 (0.8)	1.1
Precipath U	4.74 (183)	0.04 (2)	0.9
Human serum 1	2.85 (110)	0.03 (1)	1.1
Human serum 2	7.39 (286)	0.05 (2)	0.7



# CHOL2

## Cholesterol Gen.2



Intermediate precision	Mean mmol/L (mg/dL)	SD mmol/L (mg/dL)	CV %
Precinorm U	2.31 (89.3)	0.04 (1.6)	1.6
Precipath U	4.85 (188)	0.08 (3)	1.6
Human serum 3	1.97 (76.2)	0.03 (1.2)	1.6
Human serum 4	7.13 (276)	0.10 (4)	1.4

### Method comparison

Cholesterol values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c** 501 analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x).

Sample size (n) = 266

Passing/Bablok<sup>23</sup>

$$y = 1.002x + 0.045 \text{ mmol/L}$$

$$r = 0.953$$

Linear regression

$$y = 1.012x - 0.015 \text{ mmol/L}$$

$$r = 0.997$$

The sample concentrations were between 1.53 and 18.5 mmol/L (59.1 and 715 mg/dL).

### References

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer Verlag 1995.
- Liebermann NC. Über das Oxychinerterpen. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell LL, Levy BB, Kendall FE. Cholesterol in serum. In: Seligson D (ed.). Standard Methods of Clinical Chemistry. Vol 2. Academic Press, New York 1958;26-33.
- Allain CC, Poon LS, Chan CS, et al. Enzymatic determination of total serum cholesterol. Clin Chem 1974;20(4):470-475.
- Roeschlau P, Bernt E, Gruber W. Enzymatic determination of total cholesterol in serum. Z Klin Chem Klin Biochem 1974;12(5):226.
- Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. Ann Clin Biochem 1969;6:24-27.
- Siedel J, Hägele EO, Ziegenhorn J, et al. Reagent for the enzymatic determination of serum total cholesterol with improved lipolytic efficiency. Clin Chem 1983;29:1075-1080.
- Wiebe DA, Bernert JT. Influence of incomplete cholesteryl ester hydrolysis on enzymatic measurements of cholesterol. Clin Chem 1984;30:352-356.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995 Dec;119(12):1127-1135.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press p.176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995;130-131.

- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- Siekman L, Hüskes KP, Breuer H. Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 1976;279:145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT



Contents of kit

Volume after reconstitution or mixing

### FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

