

ISE Clean. Soln./Elecsys SysClean cobas®

ISE Cleaning Solution / Elecsys SysClean

Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
11298500 316	ISE Cleaning Solution / Elecsys SysClean (5 x 100 mL)	Roche/Hitachi Roche/Hitachi cobas c Elecsys

English

Intended use

For the cleaning of ISE units on Roche/Hitachi analyzers.
For the cleaning of Elecsys analyzers.

Summary

ISE Cleaning Solution / Elecsys SysClean is an alkaline cleaning solution with antimicrobial properties.

Reagents - working solutions

5 bottles, each containing 100 mL

Components:

Sodium hydroxide 3 mol/L

Sodium hypochlorite solution (< 2 % active chlorine)

Additive

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

Prevention:

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

Response:

P303 + P361 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
+ P310 Immediately call a POISON CENTER or doctor/physician.

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Storage and stability

Store protected from light.

The solution is stable up to the stated expiration date when stored at 2-8 °C.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- General laboratory equipment

- Roche/Hitachi and **cobas c** systems or Elecsys analyzer

Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT

Contents of kit



Volume after reconstitution or mixing

GTIN

Global Trade Item Number

COBAS and ELECSYS are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

