

Insulin

Insulin

cobas®

| REF | | SYSTEM |
|--------------|-----|--|
| 12017547 122 | 100 | Elecsys 2010 MODULAR ANALYTICS E170 cobas e 411 cobas e 601 cobas e 602 |

English

Intended use

Immunoassay for the in vitro quantitative determination of human insulin in human serum and plasma. The determination of insulin is utilized in the diagnosis and therapy of various disorders of carbohydrate metabolism, including diabetes mellitus and hypoglycemia.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

Summary

Insulin is a peptide hormone with a molecular weight of approximately 6000 daltons. It is secreted by the B-cells of the pancreas and passes into circulation via the portal vein and the liver. Insulin is generally released in pulses, with the parallel glucose cycle normally about 2 minutes ahead of the insulin cycle.¹

The insulin molecule consists of two polypeptide chains, the α -chain with 21 and the β -chain with 30 amino acids. Biosynthesis of the hormone takes place in the β -cells of the islets of Langerhans in the form of single-chain proinsulin, which is immediately cleaved to give proinsulin. Specific proteases cleave proinsulin to insulin and C-peptide which pass into the bloodstream simultaneously. About half of the insulin, but virtually none of the C-peptide, is retained in the liver. Circulating insulin has a half-life of 3-5 minutes and is preferentially degraded in the liver, whereas inactivation or excretion of proinsulin and C-peptide mainly takes place in the kidneys.

The amino acid sequence of insulin has remained surprisingly constant during evolution, with the result that prior to the development of genetically engineered human insulin it was possible to successfully use porcine or bovine insulin in the therapy of diabetes mellitus.²

The action of insulin is mediated by specific receptors and primarily consists of facilitation of the uptake of sugar by the cells of the liver, fatty tissue and musculature; this is the basis of its hypoglycemic action.

Serum insulin determinations are mainly performed on patients with symptoms of hypoglycemia. They are used to ascertain the glucose/insulin quotients and for clarification of questions concerning insulin secretion, e.g. in the tolbutamide test and glucagon test or in the evaluation of oral glucose tolerance tests or hunger provocation tests.

Although the adequacy of pancreatic insulin synthesis is frequently assessed via the determination of C-peptide, it is still generally necessary to determine insulin. For example, therapeutic administration of insulins of non-human origin can lead to the formation of anti-insulin antibodies. In this case, measurement of the concentration of serum insulin shows the quantity of free - and hence biologically active - hormone, whereas the determination of C-peptide provides a measure of the patient's total endogenous insulin secretion.^{3,4,5}

A disorder in insulin metabolism leads to massive influencing of a number of metabolic processes. A too low concentration of free, biologically active insulin can lead to the development of diabetes mellitus. Possible causes of this include destruction of the β -cells (type I diabetes), reduced activity of the insulin or reduced pancreatic synthesis (type II), circulating antibodies to insulin, delayed release of insulin or the absence (or inadequacy) of insulin receptors.

On the other hand, autonomous, non-regulated insulin secretion is generally the cause of hypoglycemia. This condition is brought about by inhibition of gluconeogenesis, e.g. as a result of severe hepatic or renal failure, islet cell adenoma, or carcinoma. Hypoglycemia can, however, also be facilitated intentionally or unintentionally (factitious hypoglycemia).

In 3 % of persons with reduced glucose tolerance, the metabolic state deteriorates towards diabetes mellitus over a period of time. Reduced glucose tolerance during pregnancy always requires treatment. The clearly elevated risk of mortality for the fetus necessitates intensive monitoring.

The Elecsys Insulin assay employs two monoclonal antibodies which together are specific for human insulin.

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: Insulin from 20 μ L sample, a biotinylated monoclonal insulin-specific antibody, and a monoclonal insulin-specific antibody labeled with a ruthenium complex^{a)} form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex ($\text{Ru}(\text{bpy})_3^{2+}$)

Reagents - working solutions

The reagent rackpack is labeled as INSULIN.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL:
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-insulin-Ab~biotin (gray cap), 1 bottle, 10 mL:
Biotinylated monoclonal anti-insulin antibody (mouse) 1 mg/L; MES^{b)} buffer 50 mmol/L, pH 6.0; preservative.
- R2 Anti-insulin-Ab~ $\text{Ru}(\text{bpy})_3^{2+}$ (black cap), 1 bottle, 10 mL:
Monoclonal anti-insulin antibody (mouse) labeled with ruthenium complex 1.75 mg/L; MES buffer 50 mmol/L, pH 6.0; preservative.

b) MES = 2-morpholino-ethane sulfonic acid

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

| Stability: | |
|-------------------------|----------------------------------|
| unopened at 2-8 °C | up to the stated expiration date |
| after opening at 2-8 °C | 12 weeks |
| on the analyzers | 4 weeks |

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.



Insulin

Insulin

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K₃-EDTA and sodium citrate plasma.

Hemolysis interferes, as insulin-degrading peptidases are released from erythrocytes.⁶

Criterion: Recovery within 90-110 % of serum value or slope 0.9-1.1 + intercept within $\pm 2x$ analytical sensitivity (LDL) + coefficient of correlation > 0.95 .

Stable for 24 hours at 2-8 °C, 6 months at -20 °C. Freeze only once.⁷

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 12017504122, Insulin CalSet, for 4 x 1 mL
 - [REF] 05341787190, PreciControl Multimarker, for 3 x 2 mL each of PreciControl Multimarker 1 and 2 or [REF] 11731416190, PreciControl Universal, for 2 x 3 mL each of PreciControl Universal 1 and 2
 - [REF] 11731416160, PreciControl Universal, for 2 x 3 mL each of PreciControl Universal 1 and 2 (for USA) or [REF] 05341787160, PreciControl Multimarker, for 3 x 2 mL each of PreciControl Multimarker 1 and 2 (for USA)
 - General laboratory equipment
 - Elecsys 2010, MODULAR ANALYTICS E170 or **cobas e** analyzer
- Accessories for Elecsys 2010 and **cobas e** 411 analyzers:
- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
 - [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
 - [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
 - [REF] 11933159001, Adapter for SysClean
 - [REF] 11706802001, Elecsys 2010 AssayCup, 60 x 60 reaction vessels
 - [REF] 11706799001, Elecsys 2010 AssayTip, 30 x 120 pipette tips

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, Elecsys SysClean, 5 x 100 mL system cleaning solution

- [REF] 11298500160, Elecsys SysClean, 5 x 100 mL system cleaning solution (for USA)

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

Bring the cooled reagents to approx. 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized using the 1st IRP WHO Reference Standard 66/304 (NIBSC).

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 1 month (28 days) when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Multimarker or PreciControl Universal. In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Please note: Commercial controls may contain insulin of animal origin. When assessing results, the corresponding cross-reactivity of this test must be taken into account; see under "Analytical specificity".

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in µU/mL or pmol/L).

Conversion factors: $\mu\text{U/mL} \times 6.945 = \text{pmol/L}$
 $\text{pmol/L} \times 0.144 = \mu\text{U/mL}$

Limitations - interference

The assay is unaffected by icterus (bilirubin $< 1539 \mu\text{mol/L}$ or $< 90 \text{ mg/dL}$), lipemia (Intralipid $< 1800 \text{ mg/dL}$) and biotin ($< 246 \text{ nmol/L}$ or $< 60 \text{ ng/mL}$).

Criterion: Recovery within $\pm 10 \%$ of initial value. Hemolysis interferes.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. $> 5 \text{ mg/day}$) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 18900 IU/mL.

There is no high-dose hook effect at insulin concentrations up to 20000 µU/mL or 138900 pmol/L.

In vitro tests were performed on 20 commonly used pharmaceuticals. No interference with the assay was found.

Samples from patients treated with bovine, porcine or human insulin sometimes contain anti-insulin antibodies which can affect the test results.^{2, 8,9}



Insulin

Insulin

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

0.2-1000 µU/mL or 1.39-6945 pmol/L (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 0.2 µU/mL (< 1.39 pmol/L). Values above the measuring range are reported as > 1000 µU/mL (> 6945 pmol/L).

Lower limits of measurement

Lower detection limit of the test

Lower detection limit: 0.2 µU/mL (1.39 pmol/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

Dilution

Not necessary due to the broad measuring range.

Expected values

Studies with the Elecsys Insulin assay conducted in a clinical center in Germany with samples from 57 healthy, fasting individuals gave the following results (5th-95th percentile range):

2.6-24.9 µU/mL (17.8-173 pmol/L)

Status: Elecsys Insulin MCE, study No.: B99P027 of 29 March 2001.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents and pooled human sera in a modified protocol (EP5-A) of the CLSI (Clinical and Laboratory Standards Institute): 6 times daily for 10 days (n = 60); repeatability on MODULAR ANALYTICS E170 analyzer, n = 21. The following results were obtained:

| Elecsys 2010 and cobas e 411 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Repeatability | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| Human serum 1 | 6.36 | 44.2 | 0.122 | 0.847 | 1.9 |
| Human serum 2 | 20.9 | 145 | 0.391 | 2.71 | 1.9 |
| Human serum 3 | 747 | 5188 | 15.1 | 105 | 2.0 |

| Elecsys 2010 and cobas e 411 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Intermediate precision | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| Human serum 1 | 6.36 | 44.2 | 0.163 | 1.11 | 2.6 |
| Human serum 2 | 20.9 | 145 | 0.593 | 4.10 | 2.8 |
| Human serum 3 | 747 | 5188 | 18.6 | 129 | 2.5 |

cobas[®]

| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Repeatability | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| Human serum 1 | 5.93 | 41.2 | 0.09 | 0.62 | 1.5 |
| Human serum 2 | 14.5 | 101 | 0.13 | 0.92 | 0.9 |
| Human serum 3 | 49.9 | 346 | 0.58 | 4.05 | 1.2 |
| Human serum 4 | 399 | 2768 | 3.32 | 23.1 | 0.8 |

| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Intermediate precision | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| Human serum 1 | 6.85 | 47.6 | 0.336 | 2.33 | 4.9 |
| Human serum 2 | 16.7 | 116 | 0.616 | 4.28 | 3.7 |
| Human serum 3 | 55.1 | 383 | 1.86 | 12.9 | 3.4 |
| Human serum 4 | 425 | 2949 | 10.0 | 69.6 | 2.4 |

Precision was determined using Elecsys reagents and controls in a protocol (EP5-A2) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplication each for 21 days (n = 84). The following results were obtained:

| Elecsys 2010 and cobas e 411 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Repeatability | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| PreciControl MM ^{c)} 1 | 23.7 | 165 | 0.270 | 1.88 | 1.1 |
| PreciControl MM2 | 81.7 | 567 | 1.14 | 7.92 | 1.4 |

c) MM = Multimer

| Elecsys 2010 and cobas e 411 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Intermediate precision | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| PreciControl MM1 | 23.7 | 165 | 0.834 | 5.79 | 3.5 |
| PreciControl MM2 | 81.7 | 567 | 3.04 | 21.1 | 3.7 |

| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Repeatability | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| PreciControl MM1 | 21.9 | 152 | 0.712 | 4.94 | 3.2 |
| PreciControl MM2 | 74.3 | 516 | 2.72 | 18.9 | 3.7 |

| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---|-------|--------|-------|--------|-----|
| Intermediate precision | | | | | |
| Sample | Mean | | SD | | CV |
| | µU/mL | pmol/L | µU/mL | pmol/L | |
| PreciControl MM1 | 21.9 | 152 | 0.926 | 6.43 | 4.2 |
| PreciControl MM2 | 74.3 | 516 | 3.42 | 23.8 | 4.6 |

Method comparison

a) A comparison of the Elecsys Insulin assay (y) with the Enzymun-Test Insulin method (x) using clinical samples gave the following correlations (µU/mL):



Insulin

Insulin

Number of samples measured: 99

Passing/Bablok¹⁰

$$y = 1.00x - 1.16$$

$$r = 0.844$$

Linear regression

$$y = 0.92x + 0.59$$

$$r = 0.958$$

The sample concentrations were between approximately 3.9 and 80 µU/mL (approximately 27 and 550 pmol/L).

b) A comparison of the Elecsys Insulin assay (y) with a commercially available Insulin test (x) using clinical samples gave the following correlations (µU/mL):

Number of samples measured: 99

Passing/Bablok¹⁰

$$y = 0.89x - 0.62$$

$$r = 0.935$$

Linear regression

$$y = 0.93x - 1.02$$

$$r = 0.981$$

The sample concentrations were between approximately 1 and 118 µU/mL (approximately 7 and 820 pmol/L).

Analytical specificity

For the monoclonal antibodies used, the following cross-reactivities were found:

| | Concentration tested | Cross-reactivity % |
|------------------------------|----------------------|--------------------|
| Bovine insulin | 17360 pmol/L | 25.0 |
| Porcine insulin | 8334 pmol/L | 19.2 |
| Human proinsulin | 1000 ng/mL | 0.05 |
| C-peptide | 100 ng/mL | n.d. ^{d)} |
| Glucagon | 1000 pg/mL | n.d. |
| Somatostatin | 100 pg/mL | n.d. |
| Insulin-like growth factor I | 6579 pmol/L | 0.04 |

d) n.d. = not detectable

Results for cross-reactivity with recombinant insulin analogs in a number of insulin methods have been published for example by two groups in France and the USA.^{9,11,12} The following results were published by Owen et al.¹¹ for the Elecsys Insulin assay:

Insulin lispro, insulin aspart, and insulin glargine were each tested in concentrations of 30, 100, 300, and 1000 mIU/L in the absence of insulin. The results obtained were below the detection limit of the Elecsys Insulin assay (< 0.2 µU/mL or < 1.39 pmol/L) at all the concentrations tested.

Moreover, these results also correlate with those published earlier by Sapin et al. for insulin lispro.⁹

References

- Lang DA, Matthews DR, Peto J, et al. Cyclic oscillations of basal plasma glucose and insulin concentrations in human beings. *N Engl J Med* 1979;301:1023-1027.
- Fiedler H. Basiswissen Labordiagnostik: Diabetes mellitus und metabolisches Syndrom. Broschüre Roche Diagnostics 1999;14,67 Best.-Nr. 1951769.
- Arnqvist H, Olsson PO, von Schenck H. Free and Total Insulin as Determined after Precipitation with Polyethylene Glycol: Analytical Characteristics and Effects of Sample Handling and Storage. *Clin Chem* 1987;33(1):93-96.
- Gerbitz KD. Pankreatische B-Zellen Peptide: Kinetik und Konzentration von Proinsulin, Insulin und C-Peptid in Plasma und Urin, Probleme der Meßmethoden, klinische Aussage und Literaturübersicht. *J Clin Chem Clin Biochem* 1980;18(6):313-326.
- Clark PM. Assays for insulin, proinsulin(s) and C-peptide. *Ann Clin Biochem* 1999;36(5):541-564.
- Chevenne D, Letailleur A, Trivin F, et al. Effect of Hemolysis on the Concentration of Insulin in Serum Determined by RIA and IRMA. *Clin Chem* 1998;44(2):354-356.
- DG Klinische Chemie Mitteilungen 1995;26(5):207-224.

- Tietz NW. *Clinical Guide To Laboratory Tests*. 3rd ed. Philadelphia, Pa: WB Saunders Co, 1995:366-367.

- Sapin R, Le Galudec V, Gasser F, et al. Elecsys Insulin Assay: Free Insulin Determination and the Absence of Cross-Reactivity with Insulin Lispro. *Clin Chem* 2001;47:602-605.

- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

- Owen WE, Roberts WL. Letter to the Editor: Cross-Reactivity of Three Recombinant Insulin Analogs with Five Commercial Insulin Immunoassays. *Clin Chem* 2004;50(1):257-259.

- Sapin R. Review: Insulin Assays: Previously Known and New Analytical Features. *Clin Lab* 2003;49(3-4):113-121.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

| | |
|--|---|
| | Contents of kit |
| | Analyzers/Instruments on which reagents can be used |
| | Reagent |
| | Calibrator |
| | Volume after reconstitution or mixing |

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS E, ELECSYS, MODULAR and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

