

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

<p>2D/LAEDVI(A-L)</p> <p>Alias: LAEDV Index (A-L)</p>	<p>(GEU-106-0027, 99GEMS, “Left Atrial End Diastolic Volume Index”)</p>	<p>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125204, DCM, “Area-Length Biplane”)</p>
<p>2D/LAESV(A-L)</p> <p>Alias: LAESV(A-L)</p>	<p>(G-0383, SRT, “Left Atrium Systolic Volume”)</p>	<p>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125204, DCM, “Area-Length Biplane”)</p>
<p>LAEDV(MOD BP)</p> <p>Alias: LAEDV(MOD BP)</p>	<p>(122407, DCM, “Left Atrial End Diastolic Volume”)</p>	<p>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement</p>

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		Method”) = (125207, DCM, “Method of Disks, Biplane”)
<p>LAESV(MOD BP)</p> <p>Alias: LAESV(MOD BP)</p>	<p>(G-0383, SRT, “Left Atrium Systolic Volume”)</p>	<p>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125207, DCM, “Method of Disks, Biplane”)</p>
<p>SD/Laappendix Vmax</p> <p>Alias: Laappendix Vmax</p>	<p>(29486-8, LN, “Left Atrial Appendage Peak Velocity”)</p>	
<p>AFILA/2DLA_ReservoirStrain_R_Wave(A4C)</p> <p>Alias : 2DLA_ReservoirStrain_R_Wave(A4C)</p>	<p>(GEU-106- 0161, 99GEMS, “Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave”)</p>	<p>(G-C036, SRT, “Measurement Method”) = (GEU, GEU- 106-0018, “AFI”) (111031, DCM, “Image View”) = (SRT, G- A19C, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle Point”) = (GEU, GEU-106-0158, “Reservoir phase of the left atrium”)</p>

DIRECTION DOC2652554 REV 3

<p>AFILA/2DLA_ConduitStrain_R_Wave(A4C) Alias : 2DLA_ConduitStrain_R_Wave(A4C)</p>	<p>(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")</p>
<p>AFILA/2DLA_ContractileStrain_R_Wave(A4C) Alias : 2DLA_ContractileStrain_R_Wave(A4C)</p>	<p>(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")</p>
<p>AFILA/2DLA_ReservoirStrain_P_Wave(A4C) Alias : 2DLA_ReservoirStrain_P_Wave(A4C)</p>	<p>(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View")</p>

DIRECTION DOC2652554 REV 3

		= (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_P_Wave(A4C) Alias : 2DLA_ConduitStrain_P_Wave(A4C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_P_Wave(A4C) Alias : 2DLA_ContractileStrain_P_Wave(A4C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		phase of the left atrium")
AFILA/2DLA_ReservoirStrain_R_Wave(A2C) Alias : 2DLA_ReservoirStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_R_Wave(A2C) Alias : 2DLA_ConduitStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_R_Wave(A2C) Alias : 2DLA_ContractileStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	with 0 strain at ECG R-wave"))	(111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium"))
AFILA/2DLA_ReservoirStrain_P_Wave(A2C) Alias : 2DLA_ReservoirStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave"))	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium"))
AFILA/2DLA_ConduitStrain_P_Wave(A2C) Alias : 2DLA_ConduitStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave"))	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_P_Wave(A2C) Alias : 2DLA_ContractileStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")
AFILA/2DLA_ReservoirStrain_R_Wave(BiP) Alias : 2DLA_ReservoirStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_R_Wave(BiP) Alias : 2DLA_ConduitStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_R_Wave(BiP) Alias : 2DLA_ContractileStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")
AFILA/2DLA_ReservoirStrain_P_Wave(BiP) Alias : 2DLA_ReservoirStrain_P_Wave(BiP)	(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_P_Wave(BiP) Alias : 2DLA_ConduitStrain_P_Wave(BiP)	(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

<p>AFILA/2DLA_ContractileStrain_P_Wave(BiP) Alias : 2DLA_ContractileStrain_P_Wave(BiP)</p>	<p>(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")</p>
<p>AFILA/2DLA_EF(A2C) Alias : 2DLA_EF(A2C)</p>	<p>(GEU-106-0165, 99GEMS, "Left Atrium Emptying Fraction by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")</p>
<p>AFILA/2DLA_EV(A2C) Alias : 2DLA_EV(A2C)</p>	<p>(GEU-106-0166, 99GEMS, "Left Atrium Emptying Volume by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")</p>
<p>AFILA/2DLA_Vmax(A2C) Alias : 2DLA_Vmax(A2C)</p>	<p>(GEU-106-0167, 99GEMS, "Left Atrium maximal volume by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-</p>

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		A19B, "Apical two chamber")
AFILA/2DLA_Vmin(A2C) Alias : 2DLA_Vmin(A2C)	(GEU-106-0168, 99GEMS, "Left Atrium minimal volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")
AFILA/2DLA_VpreA(A2C) Alias : 2DLA_VpreA(A2C)	(GEU-106-0169, 99GEMS, "Left Atrium volume at preA time by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")
AFILA/2DLA_EF(A4C) Alias : 2DLA_EF(A4C)	(GEU-106-0165, 99GEMS, "Left Atrium Emptying Fraction by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_EV(A4C) Alias : 2DLA_EV(A4C)	(GEU-106-0166, 99GEMS, "Left Atrium Emptying Volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		= (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_Vmax(A4C) Alias : 2DLA_Vmax(A4C)	(GEU-106-0167, 99GEMS, "Left Atrium maximal volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_Vmin(A4C) Alias : 2DLA_Vmin(A4C)	(GEU-106-0168, 99GEMS, "Left Atrium minimal volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_VpreA(A4C) Alias : 2DLA_VpreA(A4C)	(GEU-106-0169, 99GEMS, "Left Atrium volume at preA time by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_EF(BiP) Alias : 2DLA_EF(BiP)	(GEU-106-0170, 99GEMS, "Left Atrium Emptying Fraction from biplane measurements")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	by speckle tracking")	
AFILA/2DLA_EV(BiP) Alias : 2DLA_EV(BiP)	(GEU-106-0171, 99GEMS, "Left Atrium Emptying Volume from biplane measurements by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")
AFILA/2DLA_Vmax(BiP) Alias : 2DLA_Vmax(BiP)	(GEU-106-0172, 99GEMS, "Left Atrium maximal volume from biplane measurements by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")
AFILA/2DLA_Vmin(BiP) Alias : 2DLA_Vmin(BiP)	(GEU-106-0173, 99GEMS, "Left Atrium minimal volume from biplane measurements by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")
AFILA/2DLA_VpreA(BiP) Alias : 2DLA_VpreA(BiP)	(GEU-106-0174, 99GEMS, "Left Atrium volume at preA time from biplane measurements by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")

Section Right Atrium

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/RA Alias: RA Diam	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/RA Area Alias: RA Area	(G-A166, SRT, “Area”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
2D/RAA diastole Alias: RAA d	(GEU-106-0059, 99GEMS, “Right Atrium Area at Diastole”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32010, SRT, “Diastole”)
2D/RAA systole Alias: RAA s	(GEU-106-0060, 99GEMS, “Right Atrium Area at Systole”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
2D/RAD Major Alias: RA Major	(G-A193, SRT, “Major Axis”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/RAD Minor Alias: RA Minor	(G-A194, SRT, “Minor Axis”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode
RAAs(A4C) Alias: RAAs	(17988-7, LN, “Right Atrium Systolic Area”)	
RALd(A4C) Alias: RALd A4C	(29466-0, LN, Right Atrium Superior- Inferior Dimension, 4- chamber view”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)
RAAd(A4C) Alias: RAAd A4C	(17988-7, LN, “Right Atrium Area on Apical four chamber view”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

RAEDV(A-L A4C) Alias: RAEDV AL A4C	(GEU-106-0104, 99GEMS, “Right Atrium Volume”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area- Length Single Plane”)
RAEDV(MOD A4C) Alias: RAEDV MOD A4C	(GEU-106-0104, 99GEMS, “Right Atrium Volume”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)
RAESV(A-L A4C) Alias: RAESV (A-L A4C)	(GEU-106-106, 99GEMS, “Right Atrium End Systolic Volume”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (109070, DCM, “End Systole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area- Length Single Plane”)
RAESV(MOD A4C) Alias: RAESV(MOD A4C)	(GEU-106-106, 99GEMS, “Right Atrium End Systolic Volume”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (109070, DCM, “End Systole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)
RALs(A4C) Alias: RALs(A4C)	(29466-0, LN, “Right Atrium Superior- Inferior Dimension, 4- chamber view”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (109070, DCM, “End Systole”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)

DIRECTION DOC2652554 REV 3

RAP Alias: RAP	(18070-3, LN, “Right Atrium Systolic Pressure”)	
-------------------------------------	--	--

Section Aortic Valve

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/AVA/AV Diam Alias: AV Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
AV Dec Time Alias: AV Dec Time	(20217-6, LN, “Deceleration Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV Dec Slope Alias: AV Dec Slope	(20216-8, LN, “Deceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PISA/AR/RF Alias: AR RF	(G-0390, SRT, “Regurgitant Fraction”)	(G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
MM/AV Diam Alias: AV Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)
MM/%IVS Thck Alias: %IVS Thck	(18054-7, LN, “Interventricular	(G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	Septum % Thickening")	
2D/AV Diam Alias: AV Diam	(G-038F, SRT, "Cardiovascular Orifice Diameter")	(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode
2D/AV Cusp Alias: AV Cusp	(17996-0, LN, "Aortic Valve Cusp Separation")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode
2D/LAX/Trans AVA diastole Alias: Trans AVA (d)	(G-038E, SRT, "Cardiovascular Orifice Area")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
2D/LAX/Trans AVA systole Alias: Trans AVA (s)	(G-038E, SRT, "Cardiovascular Orifice Area")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole")
2D/SAX/Trans AVA diastole Alias: Trans AVA (d)	(G-038E, SRT, "Cardiovascular Orifice Area")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on

DIRECTION DOC2652554 REV 3

		scan mode (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32010, SRT, “Diastole”)
2D/SAX/Trans AVA systole Alias: Trans AVA (s)	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
2D/AVA Planimetry Alias: AVA Planimetry	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125220, DCM, “Planimetry”)
2D/AV Area Alias: AV Area	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
MM/AV Cusp Alias: AV Cusp	(17996-0, LN, “Aortic Valve Cusp Separation”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

AV Vmax Alias: AV Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV Vmax P Alias: AV Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV maxPG Alias: AV maxPG	(20247-3, LN, “Peak Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV Vmean Alias: AV Vmean	(20352-1, LN, “Mean Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV meanPG Alias: AV meanPG	(20256-4, LN, “Mean Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV Acc Time Alias: AV AccT	(20168-1, LN, “Acceleration Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”)
AV VTI Alias: AV VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AV Env. Ti Alias: AV Env. Ti	(GEU-106-0080, 99GEMS, “Time duration of the VTI trace on Aortic Valve”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AVA (VTI) Alias: AVA (VTI)	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		“Antegrade Flow” (G-C036, SRT, “Measurement Method”) = (125215, DCM, “Continuity Equation by Velocity Time Integral”)
AVA (Vmax) Alias: AVA Vmax	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)
AVA (Vmax)2 Alias: AVA Vmax, Pt	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)
AVA (Vmax)P Alias: AVA Vmax, Pt	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)
AVA (Vmax)P2 Alias: AVA Vmax	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

AV SV Alias: AV SV	(F-32120, SRT, “Stroke Volume”)	(G-C0E3, SRT, «Finding Site») = (T-42000, SRT, «Aorta»)
AV HR Alias: HR	(8867-4, LN, “Heart rate”)	
AV SI Alias: AV SI	(F-00078, SRT, “Stroke Index”)	(G-C0E3, SRT, «Finding Site») = (T-42000, SRT, «Aorta»)
AV Time To Peak Alias: AV Time to Peak	(GEU-106-0006, 99GEMS, “Time to Peak”)	
AV CO Alias: AV CO	(F-32100, SRT, “Cardiac Output”)	(G-C0E3, SRT, «Finding Site») = (T-42000, SRT, «Aorta»)
AV CI Alias: AV CI	(F-32110, SRT, “Cardiac Index”)	(G-C0E3, SRT, «Finding Site») = (T-42000, SRT, «Aorta»)
AV Acc Slope Alias: AV Acc Slope	(20167-3, LN, “Acceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”)
AVET Alias: AVET	(18041-4, LN, “Aortic Valve Ejection Time”)	(G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”)
AV Acc Time/ET Ratio Alias: AV Acc Time/ET	(G-0382, SRT, “Ratio of Aortic Valve Acceleration Time to Ejection Time”)	(G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”)
AV dp/dt Alias: AV dp dt	(59120-6, LN, “Aortic valve antegrade dp/dt [pressure rate] by US”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
AR PHT Alias: AR PHT	(20280-4, LN, “Pressure Half-Time”)	(G-C048, SRT, “Direction of Flow”) =

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		(R-42E61, SRT, “Regurgitant Flow”)
AR Dec Time Alias: AR Dec Time	(20217-6, LN, “Deceleration Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR Dec Slope Alias: AR Dec Slope	(20216-8, LN, “Deceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR Vmax Alias: AR Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR maxPG Alias: AR maxPG	(20247-3, LN, “Peak Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR Env. Ti Alias: AR Env. Ti	(GEU-106-0082, 99GEMS, “Time duration of the VTI trace on Aortic Regurgitant flow”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR HR Alias: AR HR	(8867-4, LN, “Heart rate”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
Arend Vmax Alias: Arend Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

Arend maxPG Alias: Arend PG	(20247-3, LN, “Peak Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”)
AR Vmean Alias: AR Vmean	(20352-1, LN, “Mean Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR meanPG Alias: AR meanPG	(20256-4, LN, “Mean Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR VTI Alias: AR VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
AR dp/dt Alias: AR dp/dt	(59120-6, LN, “Aortic regurgitant dp/dt US pressure by rate”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
PISA/AR/Flow Alias: AR Flow	(34141-2, LN, “Peak Instantaneous Flow Rate”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/AR/Radius Alias: AR Rad	(GEU-106-0004, 99GEMS, “Flow Radius”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		“Proximal Isovelocity Surface Area”)
PISA/AR/Velocity Alias: AR Als.Vel	(GEU-106-0005, 99GEMS, “Alias Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/AR/Vmax Alias: AR Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/AR/VTI Alias: AR VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/AR/ERO Alias: AR ERO	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/AR/RV Alias: AR RV	(33878-0, LN, “Volume Flow”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		“Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
2D/AV Annulus Diam Alias: AV Annulus Diam	(79940-3, LN, “Aortic valve annulus Diameter at end systole by US 2D”)	
AA_DIAMETER(4D) Alias: AA Diameter 4D	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C0E3, SRT, “Finding Site”) = (T-35410, SRT, “Aortic Valve Ring”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0098, 99GEMS, “4DautoAVQ quantification tool”)
AA_AREA(4D) Alias: AA Area 4D	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C036, SRT, “Measurement Method”) = (GEU-106-0098, 99GEMS, “4DautoAVQ quantification tool”)
AA_CIRCUMFERENCE(4D) Alias: AA Circ 4D	(GEU-106-0098, 99GEMS, “Aortic Annulus Circumference on 4D image”)	(G-C036, SRT, “Measurement Method”) = (GEU-106-0098, 99GEMS, “4DautoAVQ quantification tool”)
AA_MAX_DIAMETER(4D) Alias: AA Max Dia 4D	(GEU-106-0099, 99GEMS, “Aortic Annulus major semi- axis from ellipse fit on 4D image”)	(G-C036, SRT, “Measurement Method”) = (GEU-106-0098, 99GEMS, “4DautoAVQ quantification tool”)
AA_MIN_DIAMTER(4D) Alias: AA Min Dia 4D	(GEU-106-0100, 99GEMS, “Aortic Annulus minor semi- axis from ellipse fit on 4D image”)	(G-C036, SRT, “Measurement Method”) = (GEU-106-0098, 99GEMS, “4DautoAVQ quantification tool”)

Section Mitral Valve

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
MV Reg Frac Alias: MV Reg Frac	(G-0390, SRT, “Regurgitant Fraction”)	
MR Acc Slope Alias: MR Acc Slope	(20167-3, LN, “Acceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
MR dp/dt Alias: MR dp/dt	(18035-6, LN, “Mitral Regurgitation dP/dt derived from Mitral Reg. velocity”)	
PISA/MR/RF Alias: MR RF	(G-0390, SRT, “Regurgitant Fraction”)	(G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
2D/MV Annulus Diam Alias: MV Ann Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C0E3, SRT, “Finding Site”) = (T-35313, SRT, “Mitral Annulus”) (G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/MV Annulus Diam AP Alias : MV Annulus Diam AP	(GEU-106-0177, 99GEMS,	(G-C0E3, SRT, "Finding Site") =

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	"Diameter in Anterior to Posterior direction")	(SRT, T-35313, "Mitral Annulus")
<p>2D/SAX/MVA</p> <p>Alias: MVA Planimetry</p>	(G-038E, SRT, "Cardiovascular Orifice Area")	<p>(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")</p> <p>(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (111031, DCM, "Image View") = (G-0397, SRT, "Parasternal short axis")</p>
<p>2D/MVA Planimetry</p> <p>Alias: MVA Planimetry</p>	(G-038E, SRT, "Cardiovascular Orifice Area")	<p>(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")</p> <p>(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (G-C036, SRT, "Measurement Method") = (125220, DCM, "Planimetry")</p>
<p>2D/MV Area</p> <p>Alias: MV Area</p>	(G-038E, SRT, "Cardiovascular Orifice Area")	(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/MV Area (PHN) Alias : MV Area (PHN)	(G-A166, SRT, "Area")	
2D/EPSS Alias: EPSS	(GEU-106-0066, 99GEMS, “E-point Spetal separation in 2D”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
4DAutoMVQ/A-P_Diameter Alias: MV A-P Diam	(GEU-106-0036, 99GEMS, “MV antero-posterior diameter by 4Dauto MV quantification tool”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU- 106-0037, 99GEMS, “4D auto MV quantification tool”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”)
4DAutoMVQ/PM-AL_Diameter Alias: MV PM-AL Diam	(GEU-106-0038, 99GEMS, “Mitral valve Diameter, medLat to AntPost hinge on commissural view”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU- 106-0037, 99GEMS, “4D auto MV quantification tool”) (GEU-106-0039, 99GEMS, “MV commissural view”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

4DAutoMVQ/Annulus_Perimeter Alias: MV Annulus Perimeter	(GEU-106-0040, 99GEMS, “4D mitral annulus perimeter”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU- 106-0037, 99GEMS, “4D auto MV quantification tool”) (T-35313, SRT, “Mitral Annulus”)
4DAutoMVQ/Annulus_Area_3D Alias: MV Annulus Area 3D	(GEU-106-0041, 99GEMS, “4D mitral annulus surface”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU- 106-0037, 99GEMS, “4D auto MV quantification tool”) (T-35313, SRT, “Mitral Annulus”)
4DAutoMVQ/Tenting_Height Alias: MV Tenting Height	(GEU-106-0042, 99GEMS, “MV tenting height from 4D Aplax”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU- 106-0037, 99GEMS, “4D auto MV quantification tool”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”)
4DAutoMVQ/Mitral-Aortic_Angle Alias: Mitral-Aortic Angle	(GEU-106-0043, 99GEMS, “Mitral- Aortic angle from 4D Aplax”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		Method”) = (GEU-106-0037, 99GEMS, “4D auto MV quantification tool”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”)
4DAutoMVQ/Inter_Trigonal_Distance Alias: MV Inter-Trigonal Dist	(GEU-106-0044, 99GEMS, “Distance between mitral trigons”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0037, 99GEMS, “4D auto MV quantification tool”)
4DAutoMVQ/Annulus_Height Alias: MV Annulus Height	(GEU-106-0045, 99GEMS, “MV annulus height”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0037, 99GEMS, “4D auto MV quantification tool”)
4DAutoMVQ/Anterior_Leaflet_Length Alias: MV Ant Leaflet Len	(GEU-106-0046, 99GEMS, “Length of anterior MV leaflet from 4D Aplax”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0037, 99GEMS, “4D auto MV quantification tool”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”) (G-C0E3, SRT, “Finding Site”) =

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		(T-35321, SRT, “Anterior Leaflet of Mitral Valve”)
4DAutoMVQ/Posterior_Leaflet_Length Alias: MV Post Leaflet Len	(GEU-106-0047, 99GEMS, “Length of posterior MV from 4D Aplax”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0037, 99GEMS, “4D auto MV quantification tool”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”) (G-C0E3, SRT, “Finding Site”) = (T-35322, SRT, “Posterior Leaflet of Mitral Valve”)
MV A VTI Alias: MV A VTI	(20354-7, LN, “Velocity Time Integral”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32030, SRT, “Atrial Systole”)
MV Dec Time Alias: MV DecT	(20217-6, LN, “Deceleration Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
MV PHT Alias: MV PHT	(20280-4, LN, “Pressure Half-Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
MV Dec Slope Alias: MV Dec Slope	(20216-8, LN, “Deceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

<p>MVA (PHT) Alias: MVA By PHT</p>	<p>(G-038E, SRT, “Cardiovascular Orifice Area”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125210, DCM, “Area by Pressure Half-Time”)</p>
<p>MVA (VTI) Alias: MVA (VTI)</p>	<p>(G-038E, SRT, “Cardiovascular Orifice Area”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125215, DCM, “Continuity Equation by Velocity Time Integral”)</p>
<p>MV meanPG Alias: MV meanPG</p>	<p>(20256-4, LN, “Mean Gradient”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)</p>
<p>MV Vmax Alias: MV Vmax</p>	<p>(11726-7, LN, “Peak Velocity”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)</p>
<p>MV Vmean Alias: MV Vmean</p>	<p>(20352-1, LN, “Mean Velocity”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)</p>
<p>MV maxPG Alias: MV maxPG</p>	<p>(20247-3, LN, “Peak Gradient”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42047,</p>

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		SRT, “Antegrade Flow”)
MV VTI Alias: MV VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
MV Time To Peak Alias: MV Time to Peak	(GEU-106-0006, 99GEMS, “Time to Peak”)	
MV HR Alias: HR	(8867-4, LN, “Heart rate”)	
MV SV Alias: MV SV	(F-32120, SRT, “Stroke Volume”)	(G-C0E3, SRT, «Finding Site») = (T-35300, SRT, «Mitral Valve»)
MV SI Alias: MV SI	(F-00078, SRT, “Stroke Index”)	(G-C0E3, SRT, «Finding Site») = (T-35300, SRT, «Mitral Valve»)
MV CO Alias: MV CO	(F-32100, SRT, “Cardiac Output”)	(G-C0E3, SRT, «Finding Site») = (T-35300, SRT, «Mitral Valve»)
MV CI Alias: MV CI	(F-32110, SRT, “Cardiac Index”)	(G-C0E3, SRT, «Finding Site») = (T-35300, SRT, «Mitral Valve»)
MV E Env.Ti Alias: MV E Env. Ti	(GEU-106-0076, 99GEMS, “Time duration of the VTI trace on Mitral Valve E-wave”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)
MV A Env. Ti Alias: MV A Env. Ti	(GEU-106-0077, 99GEMS, “Time duration of the VTI trace on Mitral Valve A-wave”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)
MV Env. Ti Alias: MV Env. Ti	(GEU-106-0078, 99GEMS, “Time	(G-0373, SRT, «Image Mode») =

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	duration of the VTI trace on Mitral Valve")	(R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")
MR Env. Ti Alias: MR Env. Ti	(GEU-106-0079, 99GEMS, "Time duration of the VTI trace on Mitral Regurgitant flow")	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow")
MVET Alias: MVET	(GEU-106-0069, 99GEMS, "Mitral Valve Ejection Time")	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») =
MV Eann Velocity Alias: MV Eann Velocity	(18037-2, LN, "Mitral Valve E-Wave Peak Velocity")	(G-C0E3, SRT, "Finding Site") = (T-35313, SRT, "Mitral Annulus")
MV E/A Ratio Alias: MV E/A Ratio	(18038-0, LN, "Mitral Valve E to A Ratio")	
MV Acc Time/MV Dec Time Alias: MV AccT/DecT	(G-0386, SRT, «Mitral Valve AT/DT Ratio»)	
MV dp/dt Alias MV dp dt	(59120-6, LN, "Mitral valve antegrade dP/dt [pressure rate] by US")	(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")
MR Vmax Alias: MR Vmax	(11726-7, LN, "Peak Velocity")	(G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

MR meanPG Alias: MR meanPG	(20256-4, LN, “Mean Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
MR Vmean Alias: MR Vmean	(20352-1, LN, “Mean Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
MR maxPG Alias: MR maxPG	(20247-3, LN, “Peak Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
MR VTI Alias: MR VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)
MCO Alias: MCO	(G-0387, SRT, “Mitral Valve Closure to Opening Time”)	(G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”)
PISA/MR/Flow Alias: MR Flow	(34141-2, LN, “Peak Instantaneous Flow Rate”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/MR/Radius Alias: MR Rad	(GEU-106-0004, 99GEMS, “Flow Radius”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		<p>“Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)</p>
<p>PISA/MR/Velocity</p> <p>Alias: MR Als.Vel</p>	<p>(GEU-106-0005, 99GEMS, “Alias Velocity”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)</p> <p>(G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)</p>
<p>PISA/MR/Vmax</p> <p>Alias: MR Vmax</p>	<p>(11726-7, LN, “Peak Velocity”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)</p> <p>(G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)</p>
<p>PISA/MR/VTI</p> <p>Alias: MR VTI</p>	<p>(20354-7, LN, “Velocity Time Integral”)</p>	<p>(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”)</p> <p>(G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)</p>

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

PISA/MR/ERO Alias: MR ERO	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
PISA/MR/RV Alias: MR RV	(33878-0, LN, “Volume Flow”)	(G-C048, SRT, “Direction of Flow”) = (R-42E61, SRT, “Regurgitant Flow”) (G-C036, SRT, “Measurement Method”) = (125216, DCM, “Proximal Isovelocity Surface Area”)
MV Eprime Velocity Alias: E’	(59133-9, LN, “Peak Tissue Velocity”)	(G-C0E3, SRT, “Finding Site”) = (T-35313, SRT, “Mitral Annulus”) (R-4089A, SRT, “Cardiac Cycle Point”) = (R-40B1B, SRT, “Early Diastole”)
MV E/Eprime Ratio/Calc Alias: E/E’	(59111-5, LN, “E Velocity to Annulus E Velocity Ratio”)	
MV E/A Ratio/Calc Alias: E/A Ratio	(18038-0, LN, “Mitral Valve E to A Ratio”)	
MV Medial Eprime Velocity Alias: E’ Sept	(59133-9, LN, “Peak Tissue Velocity”)	(G-C0E3, SRT, “Finding Site”) = (G-0391, SRT,

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		“Medial Mitral Annulus) (R-4089A, SRT, “Cardiac Cycle Point”) = (R-40B1B, SRT, “Early Diastole”)
MV Medial E/Eprime Ratio/Calc Alias: E/E’ Sept	(59111-5, LN, “E Velocity to Annulus E Velocity Ratio”)	(G-C0E3, SRT, “Finding Site”) = (G-0391, SRT, “Medial Mitral Annulus)
MV Lateral Eprime Velocity Alias: E’ Lat	(59133-9, LN, “Peak Tissue Velocity”)	(G-C0E3, SRT, “Finding Site”) = (G-0392, SRT, “Lateral Mitral Annulus”) (R-4089A, SRT, “Cardiac Cycle Point”) = (R-40B1B, SRT, “Early Diastole”)
MV Lateral E/Eprime Ratio/Calc Alias: E/E’ Lat	(59111-5, LN, “E Velocity to Annulus E Velocity Ratio”)	(G-C0E3, SRT, “Finding Site”) = (G-0392, SRT, “Lateral Mitral Annulus”)
Medial E’/2 + Lateral E’/2 (calc Avg) Alias: E’ Avg	(GEU-106-0031, 99GEMS, “Average Annulus E Velocity”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (R-40B1B, SRT, “Early Diastole”)
E/(Medial E’/2 + Lateral E’/2) Alias: E/E’ Avg	(GEU-106-0032, 99GEMS, “E Velocity to Average Annulus E Velocity”)	
MM/MAPSE Alias: MAPSE	(GEU-106-0035, 99GEMS, “Mitral Annular Plane Systolic Excursion (MAPSE)”)	(G-C0E3, SRT, “Finding Site”) = (T-35313, SRT, “Mitral Annulus”)

DIRECTION DOC2652554 REV 3

MM/MV CE Dist Alias: MV CE Dist	(59122-2, LN, “Valve C-E distance US”)	(G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)
MM/MV D-E Excursion Alias: MV D-E Excursion	(17997-8, LN, “Mitral Valve D-E Excursion”)	(G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)
MM/MV D-E Slope Alias: MV D-E Slope	(59127-1, LN, “Valve D-E slope”)	(G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)

Section Pulmonic Valve

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/PV Annulus Diam Alias: PV Ann Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/PV Annulus Diam (plax) Alias : PV Annulus Diam (plax)	(G-038F, SRT, "Cardiovascular Orifice Diameter")	(111031, DCM, "Image View") = (SRT, G-0396, "Parasternal long axis")
2D/PV Annulus Diam (sax) Alias : PV Annulus Diam (sax)	(G-038F, SRT, "Cardiovascular Orifice Diameter")	(111031, DCM, "Image View") = (SRT, G-0397, "Parasternal short axis")
2D/PV Area Alias: PV Area	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-0373, SRT, “Image Mode”) =

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		(G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
MM/Q-to-PV close Alias: Q-to-PV close	(20295-2, LN, “Time from Q wave to Pulmonic Valve Closes”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
PV Vmax Alias: PV Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV Vmax P Alias: PV Vmax	(11726-7, LN, “Peak Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV maxPG Alias: PV maxPG	(20247-3, LN, “Peak Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV Vmean Alias: PV Vmean	(20352-1, LN, “Mean Velocity”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV meanPG Alias: PV meanPG	(20256-4, LN, “Mean Gradient”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV Acc Time Alias: PV AccT	(20168-1, LN, “Acceleration Time”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV VTI Alias: PV VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV Env.Ti Alias: PV Env.Ti	(GEU-106-0086, “Time duration of the VTI trace on Pulmonic Valve”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PVA (VTI) Alias: PVA (VTI)	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125215, DCM, “Continuity Equation by Velocity Time Integral”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

PVA (Vmax) Alias: PVA (Vmax)	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)
PVA (Vmax)P Alias: PVA (Vmax)	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”) (G-C036, SRT, “Measurement Method”) = (125214, DCM, “Continuity Equation by Peak Velocity”)
PV HR Alias: HR	(8867-4, LN, “Heart rate”)	
PV SV Alias: PV SV	(F-32120, SRT, “Stroke Volume”)	
PV CO Alias: PV CO	(F-32100, SRT, “Cardiac Output”)	
PV SI Alias: PV SI	(F-00078, SRT, “Stroke Index”)	
PV CI Alias: PV CI	(F-32110, SRT, “Cardiac Index”)	
PVO Alias : PVO	(GEU-106-0143, 99GEMS, "Pulmonic Valve Opening Time")	
PVC Alias : PVC	(GEU-106-0144, 99GEMS, "Pulmonic Valve Closing Time")	
PV dp/dt Alias: PV dp/dt	(59120-6, LN, “Pulmonic valve antegrade dp/dt [pressure rate] by US”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)
PV Acc Slope Alias: PV Acc Slope	(20167-3, LN, “Acceleration Slope”)	(G-C048, SRT, “Direction of Flow”) = (R-42047, SRT, “Antegrade Flow”)