Medtronic

OCTOPUS® EVOLUTION AS

Tissue Stabilizer
Stabilisateur tissulaire
Gewebestabilisator
Estabilizador de tejido
Hartspierstabilisator
Stabilizzatore tessutale
Vevsstabilisator
Sydänstabilisator
Vävnadsstabilisator
Σταθεροποιητής ιστού
Vævsstabilisator
Stabilizator pola operacyjnego
Estabilizador de tecido
Doku Stabilizatorü
Stabilizátor tkáně
Szövetstabilizátor
Tkanivový stabilizátor
Koestabilisaator
Стабилизатор на тъкани
Uređaj za stabilizaciju tkiva
Audu stabilizēšanas instruments
Audinių stabilizatorius
Stabilizator de țesuturi
Стабилизатор тканей
Stabilizator tkiva
Pripomoček za stabilizacijo tkiva
TS2500

Instructions for Use • Mode d'emploi
Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso
Bruksanvisning • Käyttöohjeet • Bruksanvisning
Οδηγίες χρήσης • Brugsanvisning
Instrukcja użytkowania • Instruções de utilização
Kullanım Talimatları • Návod k použití
Használati útmutató • Pokyny pre používanie
Kasutusjuhend • Инструкции за употреба
Upute za upotrebu • Lietošanas pamācība
Naudojimo instrukcijos • Instrucțiuni de utilizare
Инструкция по эксплуатации • Uputstva za upotrebu
Navodila za uporabo

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TS2500

Description

The Octopus Evolution AS Tissue Stabilizer is a disposable, retractor-based stabilizer that consists of two tissue stabilizers attached to an articulating arm (arm). The arm fastens to a retractor by use of a mounting clamp. A stopcock provides on/off control of suction (Figure 1). With the tissue stabilizers placed on either side of the anastomotic site, suction is applied to stabilize the tissue. The arm is tightened and loosened by the knob. As the arm tightens, the tissue stabilizers spread in an arc fashion. Sterile, Nonpyrogenic, Disposable, Single use only.

Note: The mounting clamp has been designed to be compatible with most adult median sternotomy retractors. Retractor compatibility should be confirmed prior to beginning the procedure.

Indications

This product is intended to stabilize and minimize the motion of selected areas of the beating heart during cardiac surgery.

Contraindications

This product is not intended for use except as indicated above. Do not position the tissue stabilizers over a coronary artery or newly infarcted or aneurysmal heart tissue.

Warnings

Patient and procedure selection is the responsibility of the medical professional and the outcome is dependent on many variables, including patient anatomy, pathology, and surgical techniques.

This product is not intended for use as a heart positioner.

Precautions

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

The instructions for use are furnished for information purposes only.

This device was designed for single patient use only. Do not reuse, reprocess or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness or death.

The initial spacing between the right and left tissue stabilizers is important for optimal performance. The initial spacing between the tissue stabilizers will affect stabilization and tension on the tissue and should be chosen at the surgeon's discretion.

Canister must be oriented in a vertical position.

Do not fill past the Maximum Capacity line on canister.

Do not exceed (-)400 mm Hg of suction.

Adverse Effects

None known.

Instructions for Use

1. Inspect the packages and products for damage and expiration date. If undamaged and unexpired, open the packages and transfer the tissue stabilizer and canister tubing set onto the sterile field utilizing aseptic technique.
2. Suction circuit hook-up:

- a. Attach tubing (not supplied) from the operating room suction source to the back of the regulator. Turn regulator on and set vacuum to (-)400 mm Hg.
- b. Connect the tubing from the tissue stabilizer to the luer connection of the canister tubing set (Figure 2). Utilizing aseptic technique, pass the suction source connection end of the canister tubing set off of the sterile field and connect it to the vacuum regulator. Utilize the stopcock on the sterile field to control vacuum (on/off) during the procedure.

Note: Canister must be oriented in a vertical position.

Note: Do not fill past the Maximum Capacity line on canister.

3. Attach device to the retractor.
 - a. With clamp lever in the unlocked position, slide open the mounting clamp.
 - b. Place mounting clamp onto the retractor, ensuring clamp contact is flush with the retractor.
 - c. Move the lever to the locked position to secure the clamp to the retractor. Check that the clamp is securely fastened to the retractor (Figure 3).
4. Position tissue stabilizers onto designated anastomotic site. The following steps are recommended:
 - a. Shape the malleable portion of the tissue stabilizers to conform with the heart (Figure 4).

Caution: Do not exceed a 25 degree bend in any axis. Exceeding this angle may occlude the lumen. Repeated bending of the tissue stabilizers may compromise device performance.
 - b. Rotate turret to desired position (Figure 3).
 - c. Turn on suction by turning stopcock to the off-to-air position and gently position tissue stabilizers in the desired location on the surface of the heart (Figure 5a).
 - d. Turn the knob clockwise to immobilize the arm. As the arm tightens, the tissue stabilizers will gradually spread. The tension indicator appears red when the arm reaches maximum tension.
5. To remove the tissue stabilizers, support the heart and terminate suction by turning stopcock to the off-to-vacuum position (Figure 5b).
6. Turn the knob counterclockwise to loosen the arm, and gently lift the tissue stabilizers from the heart.
7. When finished, remove the device from the retractor and discard according to approved hospital procedures.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE OCTOPUS EVOLUTION AS TISSUE STABILIZER, MODEL TS2500, HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE OCTOPUS EVOLUTION AS TISSUE STABILIZER, MODEL TS2500, HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Medtronic® and Octopus® are registered trademarks of Medtronic, Inc.