


FULL SCORE FOR PROFESSIONAL SURFACE DISINFECTION

Incidin® Pro

The aldehyde-free concentrate for surfaces

Because of its SMART-BALANCED formulation, Incidin® Pro brings a new standard to surface cleaning and disinfection.

- ▲ Your Pro for efficacy
- ▲ Your Pro for speed
- ▲ Your Pro for compatibility



Incidin® Pro

Full score for Professional surface disinfection

NEW!

NEW INCIDIN PRO - YOUR PRO FOR EFFICACY

EFFICACY UNDILUTED	CONCENTRATION		CONTACT TIME
	%	ml/L	MIN
Surface disinfection acc. DGHM method 2001 with mechanical action (bactericidal, yeasticidal), clean & dirty	0.5 1.0	5 10	30 15
Bactericidal acc. EN 13727, clean & dirty	0.25	2.5	5
Yeasticidal acc. EN 13624, clean & dirty	0.25	2.5	5
Disinfection of surfaces acc. DGHM method 2001 with mechanical action (tuberculocidal and mycobactericidal)	2.0 4.0	20 40	120 60
Tuberculocidal and mycobactericidal acc. EN 14348, clean & dirty	1.5	15	60
Effective against enveloped viruses acc. RKI 01/2004 (incl. HIV, HBV, HCV)	0.25 0.5	2.5 5	30 15
Adenovirus*	2.0	20	60
Polyomavirus SV 40*	1.5	15	15
Rotavirus*	1.0/1.5	10/15	60/5
Norovirus (MNV)*	1.5	15	30
Norovirus (MNV) acc. EN 14476, clean & dirty	1.5	15	30

*acc. DVV 2008

WHY USE INCIDIN® PRO?

Broad spectrum efficacy with short contact times

Incidin Pro's antimicrobial efficacy is complimented by its short contact times and excellent material compatibility.

Effective against bacteria, yeast and various viruses (including norovirus), Incidin Pro is suitable for high risk areas and a wide range of surfaces, from general wards and public areas to operating rooms and intensive care units.

Excellent cleaning efficacy and material compatibility

It has been extensively tested on all surface types commonly found in the

healthcare settings - from metals, linoleum and PVCs to ceramics, rubbers and plastics - with excellent results.

EN standards compliant

Tests conducted to healthcare relevant EN standards confirmed the broad efficacy spectrum of Incidin Pro, including Norovirus, according to latest DVV and EN standard. See the table above for detailed test results.

Aldehyde-free formulation without perfume and dye

The non-volatile and fragrance free ingredient system ensures a convenient and safe product use.

Further information can be found at

www.ecolab.eu

or speak to your Ecolab representative.

**FAST.
EFFECTIVE.
ADVANCED.**

APPLICATION AREA

Incidin® Pro is a liquid concentrate suitable for cleaning and disinfection of medical devices and all kind of surfaces in healthcare settings.

INSTRUCTIONS FOR USE

Dilute concentrate with water of suitable quality (max 30°C) to the desired concentration. Wet entire surface with a sufficient amount of solution and leave to dry.

To achieve e.g. a 0.5% v/v solution prepare a 1:200 dilution with suitable water e.g. 5ml of Incidin® Pro is added to 995ml of drinking water quality, giving a total volume of 1 litre.

INGREDIENTS

100g contains the following active ingredients: 10.0g 2-Phenoxyethanol, 8.0g N,N-bis-(3-Aminopropyl) dodecylamine, 7.5g Benzalkoniumchloride

PACKAGING AND ORDERING INFORMATION

PRODUCT	OUTER SIZE	ORDER CODE
Incidin® Pro Bottle	3 x 2L	30 766 80
Incidin® Pro Canister	6L	30 767 60
Incidin® Pro Vat	30L	30 767 80
Incidin® Pro Vat	200L	30 767 90

STORAGE

Product should be stored in the original container below 25°C in a dry location out of direct sunlight.

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

ECOLAB HEALTHCARE EUROPE

Richtistrasse 7
8304 Wallisellen
Schweiz
Phone +41 (0)44 877-2000
www.ecolab.eu



ECOLAB DEUTSCHLAND GmbH

Ecolab-Allee 1
40789 Monheim am Rhein
+49 (0) 2173 599 0

