

User Manual

COMBI 200L

COMBI 200

DUO 200

PULSON 200

VACO 200



C € 0344

Manufacturer	GymnaUniphy N.V.
Main office	Pasweg 6A B-3740 BILZEN
Telephone	+(32) (0)89-510.532
Fax	+(32) (0)89-510.541
E-mail	info@gymna-uniphy.com
Website	www.gymna-uniphy.com

Your GymnaUniphy dealer:

All rights reserved. Nothing from this publication may be copied, stored in an automated data file, or made public, in any form or in any way, be it electronically, mechanically, by photocopying, recordings or in any other way, without prior written permission from GymnaUniphy N.V.

Version 2.0.2/May 2016

User Manual 200-series

Devices for electrotherapy, ultrasound therapy, combination therapy and laser therapy

	Electro therapy 	Ultrasound therapy 	Combination therapy 	Laser therapy 	Vacuum (optional) 
Combi 200L	x	x	x	x	x
Combi 200	x	x	x		x
Duo 200	x				x
Pulson 200		x			

Abbreviations

AQ	Accommodation Quotient
CC	Constant Current
CO	Combination therapy
CP	Courte Période
CV	Constant Voltage
DF	Diphassé Fixe
EL	Electrode
EMC	Electromagnetic Compatibility
ESD	Electrostatic Discharge
ET	Electrotherapy
HAC	Hospital Antiseptic Concentrate
LA	Laser therapy
LP	Longue Période
MF	Medium Frequency: with unidirectional and interferential currents Monophasé Fixe: with diadynamic currents
MTP	Myofascial Trigger Point
NMES	Neuro Muscular Electro Stimulation
TENS	Transcutaneous Electrical Nerve Stimulation
US	Ultrasound
VAS	Visual Analogue Scale

Symbols on the equipment



Read the manual



Sensitive to electrostatic discharge



Manufacturer

Symbols on the laser probe



Laser warning sign




Points to laser aperture and shows direction of beam



Type B applied part

Label on the laser probe:



Symbols in the manual



Warning or important information.



ET symbol: only for devices with electrotherapy applications, Combi 200L, Combi 200, Duo 200



US symbol: only for devices with ultrasound applications, Combi 200L, Combi 200, Pulson 200



LA symbol: only for devices with laser applications, Combi 200L



Vaco symbol: only for devices that can work with the vacuum unit, Combi 200L, Combi 200, Duo 200

TABLE OF CONTENTS

1	SAFETY	7
1.1	PURPOSE	7
1.2	SAFETY INSTRUCTIONS	7
1.3	CONTRA INDICATIONS	11
1.4	MEDICAL DEVICES DIRECTIVE	13
1.5	LIABILITY	13
2	INSTALLATION	15
2.1	RECEIPT	15
2.2	PLACING AND CONNECTION	15
2.3	PERFORMING THE FUNCTIONAL TEST	15
2.4	TOUCHSCREEN APPLICATION	15
2.5	SETTING CONTRAST AND LANGUAGE	16
2.6	USE IN COMBINATION WITH ANOTHER DEVICE	16
2.7	TRANSPORT AND STORAGE	16
2.8	RESELLING	16
3	DESCRIPTION OF THE EQUIPMENT	17
3.1	COMBI 200L AND ACCESSORIES	17
3.2	COMBI 200 AND ACCESSORIES	18
3.3	DUO 200 AND ACCESSORIES	19
3.4	PULSON 200 AND ACCESSORIES	20
3.5	VACO 200 AND STANDARD ACCESSORIES	21
3.6	COMPONENTS OF 200-SERIES	21
3.7	COMPONENTS OF VACO 200	23
3.8	DISPLAY	24
3.9	DISPLAY SYMBOLS	25
3.10	SYMBOLS FOR CURRENT SHAPES	25
3.11	PARAMETER SYMBOLS	26
3.12	CURRENT SHAPES	28
4	OPERATION	29
4.1	THERAPY SELECTION	29
4.2	SELECTION BY THE THERAPY MENU	29
4.3	SELECTION BY THE GUIDE MENU	30
4.4	PERFORMING THERAPY	32
4.5	ELECTROTHERAPY	34
4.6	ULTRASOUND THERAPY	43
4.7	COMBINATION THERAPY	46
4.8	LASER THERAPY	47
4.9	DIAGNOSTIC PROGRAMMES	49

4.10	VACUUM TREATMENT	52
4.11	MEMORY	55
4.12	SYSTEM SETTINGS	58
5	INSPECTIONS AND MAINTENANCE	63
5.1	INSPECTIONS	63
5.2	MAINTENANCE	65
6	MALFUNCTIONS, SERVICE AND GUARANTEE	69
6.1	MALFUNCTIONS	69
6.2	SERVICE	70
6.3	GUARANTEE	70
6.4	TECHNICAL LIFE TIME	71
7	TECHNICAL INFORMATION	73
7.1	GENERAL	73
7.2	GENERAL FOR VACO 200	73
7.3	ELECTROTHERAPY	74
7.4	ULTRASOUND THERAPY	77
7.5	LASER THERAPY	79
7.6	ENVIRONMENTAL CONDITIONS	79
7.7	TRANSPORT AND STORAGE	80
7.8	STANDARD ACCESSORIES	80
7.9	OPTIONAL ACCESSORIES ELECTROTHERAPY	83
7.10	OPTIONAL ACCESSORIES ULTRASOUND THERAPY	84
7.11	OPTIONAL ACCESSORIES LASER THERAPY	84
7.12	OPTIONAL ACCESSORIES VACUUM THERAPY	84
8	APPENDICES	85
8.1	AGENTS FOR IONTOPHORESIS	85
8.2	DIAGNOSTIC S-D CURVE	86
8.3	ELECTRODE, US HEAD AND LASER PROBE PLACEMENTS	87
8.4	EMC DIRECTIVE	89
8.5	TECHNICAL SAFETY INSPECTION	94
8.6	DISPOSAL	99
9.2	LITERATURE	101

1 SAFETY

1.1 Purpose

The 200-series is intended solely for medical applications. You can use the 200-series for electrotherapy, ultrasound therapy, combination therapy and laser therapy. The device is suited for continuous use.

1.2 Safety instructions

1.2.1 General



- Only qualified people who are trained in the application of the therapies may use the appliance.
- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- Follow the instructions and directions in these user instructions.
- Place the equipment on a horizontal and stable base.
- Keep the ventilation openings at the bottom and rear of the equipment free.
- Do not place any objects on the equipment.
- Do not place the equipment in the sun or above a heat source.
- Do not use the equipment in a damp area.
- Do not let any liquid flow into the equipment.
- Do not disinfect or sterilise the equipment. Clean the equipment with a dry or moistened cloth. *See §§5.*
- Only treat patients with electrical implants (pacemaker) after obtaining medical advice.
- The 'Directive on Medical Devices' from the European Commission (93/42/EEC) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. *See §§5.1.2.*
- For optimum treatment, a patient investigation must first be performed. On the basis of the findings of the investigation, a treatment plan with objectives will be formulated. Follow the treatment plan during the therapy. This will limit possible risks, related to the treatment, to a minimum.
- Always keep these user instructions with the equipment.

1.2.2 **Electrical safety**



- Only use the equipment in an area with facilities that meet the applicable legal regulations.
- Connect the equipment to an outlet with a protective earth terminal. The outlet must meet the locally applicable requirements for medical areas.

1.2.3 **Prevention of explosion**



- Do not use the equipment in an area where combustible gases or vapours are present.
- Switch off the equipment when it is not used.

1.2.4 **Electro Magnetic Compatibility**



- Medical electrical equipment requires special precautions for Electro Magnetic Compatibility (EMC). Follow the instructions for the installation of the equipment. *See* §2.
- Do not use mobile telephones or other radio, shortwave, or microwave equipment in the vicinity of the equipment. This kind of equipment can cause disturbances.
- Only use the accompanying accessories that are supplied by GymnaUniphy. *See* §7.
Other accessories can lead to an increased emission or a reduced immunity.

1.2.5 **Electrotherapy**



- Do not use the equipment simultaneously with high frequency surgical equipment. This combination can cause burning of the skin under the electrodes.
- Do not use adhesive electrodes with currents that have a galvanic component, such as galvanic, diadynamic, MF rectangular, pulsed rectangular and triangular currents. With these currents, etching of the skin can occur.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Check the electrode cables and the electrodes at least once a month. Check whether the insulation is still intact. *See §5.1.*
- The safety standards for electrical stimulation advise not to exceed the current density of $2.0 \text{ mA}_{\text{rms}}/\text{cm}^2$. However, with iontophoresis treatments, we advise a maximum current density of $0.25 \text{ mA}/\text{cm}^2$, because of using the MF rectangular current. Exceeding this value can result in skin irritation and burns.
- Always use sterilised gauze with iontophoresis treatments.

1.2.6 **Ultrasound therapy**



- Move the US head evenly over the skin during the treatment. This prevents internal burns.
- The US treatment heads are exchangeable. The device detects the characteristics and supplies the right power at the right frequency.
- Handle the US heads carefully. With rough handling, the characteristics can change. Test the US head if it falls on the ground or knocks against something. *See §5.1.1.*
- Check the US head at least once a month. During the check, look for dents, cracks and other damage that could allow liquids to ingress. Check whether the insulation of the cable is still intact. Check whether all pins are present and straight in the connectors. Replace the US head if the head, the cable or the connector is damaged. *See §5.1.*

1.2.7 **Laser therapy**

The laser is a class 3B laser product and has an invisible beam.



- Make sure the laser warning sign is clearly visible outside the entrance to the therapy room.
- The radiation of a laser probe can cause a physiological effect.
- Use the laser therapy only for therapeutic purposes.
- Use of controls or adjustments or performance of procedures other than those specified in this manual may result in hazardous radiation exposure.
- Start a laser therapy only when all persons in the room wear laser goggles for eye protection. If you do not obey this warning, you can cause blindness. Use goggles with at least the characteristics: I 100 - 1000 L2 and with a clear view of the control, the display and the signal lights. *See §7.11.*
- Do not look into the laser beam during a laser therapy.
- Do not point the laser beam into eyes.
- Do not use the laser near flammable materials or liquids.
- Do not use the equipment if any damage shows.
- Regularly check the output of the laser probe with the test facility. *See §4.12.8.*
- Check the laser probe at least once a month. During the check, look for dents, cracks and other damage. Check whether the insulation of the cable is still intact. Check whether all pins are present and straight in the connectors. Replace the laser probe if the laser, the cable or the connector is damaged. *See §5.1.*
- Place the laser probe in the holder when the laser is not used.
- Remove the key of the laser lock when the laser therapy is not used.

1.3 Contra Indications

Treatments should not be given to patients with the listed conditions:

1.3.1 *Electrotherapy*

General

High fever
Severe cardiovascular problems
Psychological problems
Cancer with tumor metastasis
Generalised tuberculosis
Eyes and testis

Specific absolute

On demand pacemakers

Specific relative for monophasic pulses

Skin lesions
Skin infections
Thrombosis, thrombophlebitis
Varices
Epiphyseal disc (children)
Bleeding tissue and increased risk to haemorrhage
Superficially implanted materials
Heart disease, rhythm disorder
Epilepsy (avoid neck electrode positions)
Decreased sensibility (patient is unable to tell what current intensity he feels)
Electrode positions near sinus caroticus
Menses
Pregnancy (not in the vicinity of the foetus, avoid trunk electrode positions)

Specific relative for biphasic pulses

Skin infections
Thrombosis, thrombophlebitis
Epiphyseal disc (children)
Bleeding tissue and increased risk to haemorrhage
Heart disease, rhythm disorder
Epilepsy (avoid neck electrode positions)
Decreased sensibility (patient is unable to tell what current intensity he feels)
Electrode positions near sinus caroticus
Pregnancy (not in the vicinity of the foetus, avoid trunk electrode positions)

Specific relative for vacuum therapy

Internal infections

Hemorrhagic risk in the part of the body where the electrodes are to be placed.

1.3.2 Ultrasound therapy **General**

High fever

Pacemaker (no US head positions in the vicinity of the pacemaker)

Severe cardiovascular problems

Psychological problems

Cancer with tumor metastasis

Generalised tuberculosis

Eyes and testis

Pregnancy (not in the vicinity of the foetus, avoid trunk US head positions)

Specific relative for continuous ultrasound

Infections

Acute inflammations

Thrombosis, thrombophlebitis

Varices

Bleeding tissue and increased increased risk to haemorrhage

Epiphyseal disc (children)

Decreased sensibility

Metal implants

Menses

Cement of endoprosthesis

Diabetes mellitus

1.3.3 Combination therapy 

See contra indications Electrotherapy and US

1.3.4 Laser therapy **General**

High fever

Severe cardiovascular problems

Psychological problems

Cancer with tumor metastasis

Generalised tuberculosis

Specific absolute

Eyes (looking into the laser beam) and testis

Thyroid gland (local applications)

Bleeding tissue and increased increased risk to haemorrhage

Hypertrophic scars

Pregnancy (not in the vicinity of the foetus)

Photo-allergy

1.4 Medical Devices Directive

The device complies with the essential requirements of the Medical Device Directive of the European Committee (93/42/EEC) as most recently changed.

The device contains no human or animal tissue, no medical substances, and no blood or blood products from human or animal origin.

1.5 Liability

The manufacturer cannot be held liable for injury to the therapist, the patient or third parties, or for damage to or by the equipment used, if for example:

- an incorrect diagnosis is made;
- the equipment or the accessories are used incorrectly;
- the user instructions are wrongly interpreted or ignored;
- the equipment is badly maintained;
- maintenance or repairs are performed by people or organisations that are not authorised to do so by GymnaUniphy.

Neither the manufacturer nor the local GymnaUniphy dealer can be held liable, in any way whatsoever, for the transfer of infections via the vaginal, anal and rectal probes and/or other accessories.

2 INSTALLATION

2.1 Receipt

1. Check whether the equipment has been damaged during transport.
2. Check whether the accessories are intact and complete. *See* § 7.
 - Inform your supplier of any damage or defects by no later than within 3 working days after receipt. Report the damage by telephone, fax, e-mail or letter.
 - Do not use the equipment if it is damaged or defective.

2.2 Placing and connection

1. Place the equipment on a horizontal and stable base.
 - Keep the ventilation openings at the bottom and rear of the equipment free.
 - Do not place the equipment in the sun or above a heat source.
 - Do not use the equipment in a wet area.
2. Check whether the mains voltage that is stated on the rear of the equipment corresponds with the voltage of your mains supply. The equipment is suited for a nominal mains voltage from 100 V to 240 VAC / 50-60 Hz.
3. Connect the device to an outlet with protective earth terminal.

2.3 Performing the functional test

1. Switch the equipment on with the switch at the rear of the equipment.
2. When the equipment is switched on, it automatically performs a test. Check whether the indicator lamps next to A and B and light briefly during the test.
3. If the lamps do not light up: *See* § 6.

2.4 Touchscreen application

The 200-series has a touchscreen. Except for the intensity knobs, all settings and treatment possibilities can be selected by touching the applicable option with a fingertip. This includes the functions that follow:

- Direct access to electrotherapy on the left side of the screen;
- The pause, stop, home, memory, return, and enter buttons below the screen;
- Changing values up Δ and down ∇ on the right side of the screen.

2.5 Setting contrast and language

1. Press  for 3 seconds. The **System settings** menu appears. See §4.12.
2. Select **Contrast**.
3. If necessary, change the contrast with Δ and ∇ .
4. Select **Language**.
5. If necessary, change the language with Δ and ∇ .
6. Select  to return to the **Start** menu.

2.6 Use in combination with another device

The Combi 200L, Combi 200, and Duo 200 can be used in combination with the Vaco 200. For information about the use of the 200-series with the Vaco 200, refer to §4.10.

2.7 Transport and storage

Take account of the following matters if the equipment has to be transported or stored:


- Transport or store the equipment in the original packaging.
- The maximum period for transport or storage is: 15 weeks.
- Temperature: -20 °C to +60 °C.
- Relative humidity: 10% to 100%.
- Atmospheric pressure: 200 hPa to 1060 hPa.

2.8 Reselling

This medical equipment must be traceable. The equipment, the US head and some other accessories have a unique serial number. Provide the dealer with the name and address of the new owner.


3 DESCRIPTION OF THE EQUIPMENT

3.1 Combi 200L and accessories




- | | |
|--|--|
| 1. Combi 200L. See §3.6. | 9. Contact gel |
| 2. Power cord | 10. Elastic fixation straps (4 pieces) |
| 3. CD-ROM User manual | 11. EL sponges for rubber electrode (4 pieces) |
| 4. Safety instructions | 12. VAS score card |
| 5. Quickstart guide | 13. Rubber electrodes (4 pieces) |
| 6. Laser goggles. Optional accessory. See §7.11. | 14. Two-ply electrode cable (2 pieces) |
| 7. US head, big | 15. Test connector |
| 8. Laser probe. Optional accessory. See §7.11. | |

3.2 Combi 200 and accessories




- | | |
|---------------------------------------|---|
| 1. Combi 200. See §3.6. | 9. EL sponges for rubber electrode (4 pieces) |
| 2. Power cord | 10. VAS score card |
| 3. CD-ROM User manual | 11. Rubber electrodes (4 pieces) |
| 4. Safety instructions | 12. Two-ply electrode cable (2 pieces) |
| 5. Quickstart guide | 13. Test connector |
| 6. US head, big | |
| 7. Contact gel | |
| 8. Elastic fixation straps (4 pieces) | |

3.3 Duo 200 and accessories




- | | |
|---------------------------------------|---|
| 1. Duo 200. See §3.6. | 7. EL sponges for rubber electrode (4 pieces) |
| 2. Power cord | 8. VAS score card |
| 3. CD-ROM User manual | 9. Rubber electrodes (4 pieces) |
| 4. Safety instructions | 10. Two-ply electrode cable (2 pieces) |
| 5. Quickstart guide | 11. Test connector |
| 6. Elastic fixation straps (4 pieces) | |

3.4 Pulson 200 and accessories



1. Pulson 200. See §3.6.
2. Power cord
3. CD-ROM User manual
4. Safety instructions
5. Quickstart guide
6. US head, big
7. US head, small (optional)
8. Contact gel
9. VAS score card


3.5 Vaco 200 and standard accessories



1. Vaco 200. See §3.7.
2. Power cord
3. Vacuum hose (4 pieces)
4. Connection cable: ET device - Phyaction CL (2 pieces)
5. Sponge for vacuum electrode (4 pieces)
6. Vacuum electrode (4 pieces)

3.6 Components of 200-series


The overview shows Combi 200L, which has the most options. For Combi 200, Duo 200, and Pulson 200, the amount of therapy buttons and rotary buttons differ per device.



- | | |
|-------------------------------|---|
| 1. Display. See §3.8. | 15. Increase button |
| 2. Holder | 16. Decrease button |
| 3. Electrotherapy button | 17. Connector for US head |
| 4. Ultrasound therapy button | 18. Connector for US head |
| 5. Combination therapy button | 19. On/off switch |
| 6. Laser therapy button | 20. Connection to mains supply |
| 7. Pause button | 21. USB port |
| 8. Stop button | 22. Ventilation opening |
| 9. Intensity of channel A | 23. Laser lock |
| 10. Home button | 24. Connector for laser probe |
| 11. Intensity of channel B | 25. Connector for electrotherapy, channel A |
| 12. Memory button | 26. Connector for electrotherapy, channel B |
| 13. Enter button | |
| 14. Return button | |


USB lizdas
12 tech parametrų reikalavimas

3.7 Components of Vaco 200



- | | |
|---|--|
| 1. Vacuum intensity regulator | 9. Output connectors B \rightarrow for vacuum electrodes |
| 2. Pulse regulator | 10. Push button channel B |
| 3. LED water reservoir | 11. Input connector channel A for connection cable |
| 4. Output connectors A \rightarrow for vacuum electrodes | 12. Input connector channel B for connection cable |
| 5. LED: traditional output channel A active | 13. On/off switch |
| 6. Output connector A \rightarrow for two-ply electrode cable | 14. Fuse holder |
| 7. LED traditional output channel B active | 15. Connection to mains supply |
| 8. Output connector B \rightarrow for two-ply electrode cable | 16. Drain hose for water reservoir |
| | 17. Air exhaust |
| | 18. Type plate |

3.8 Display



1. Selected channel
2. Therapy type
3. Current pulseshape
4. Electrotherapy type or subtype
5. Programme number
6. Parameters of the selected channel
7. Scroll through numbers with the up and down arrow keys.
8. Screen for channel A (here, electrotherapy). See §4.5.5.
9. Remaining treatment time
10. Polarity
11. Set intensity
12. Screen for channel B (here, electrotherapy). See §4.6.3.

3.9 Display symbols

3.9.1 General

	Electrotherapy	SEQ	Sequential current shapes
	Ultrasound therapy	A	Channel A
	Combination therapy	B	Channel B
	Laser therapy	A + B	Channel A and B simultaneously
	Treatment time	A ⇄ B	Alternating channels
	Treatment completed		

3.9.2 Current shape groups

	Unidirectional currents		2-pole medium frequency
	Iontophoresis		Dipole vector field
	Diadynamic		Isoplanar vector field
	TENS currents		Diagnostic programmes
	NMES currents		

3.10 Symbols for current shapes

	Medium frequency unidirectional current		Rectangular surge current
	Unidirectional rectangular current		Triangular surge current
	Rectangular pulse		Biphasic surge current
	Unidirectional triangular current		Intrapulse interval surge current

	Triangular pulse		2-pole medium frequency surge current
	Conventional TENS		Russian stimulation
	Low frequency TENS		2-pole medium frequency
	Random TENS		Dipole vector field
	Burst TENS		Isoplanar vector field
CP	CP (diadynamic)		S/d curve rectangular
DF	DF (diadynamic)		S/d curve triangular
LP	LP (diadynamic)		S/d curve rectangular + triangular
MF	MF (diadynamic)		Rheobase and chronaxie
			Rheobase and AQ


3.11 Parameter symbols

3.11.1 *Electrotherapy*

Red+ Red-	Polarity indication	CC	Constant Current
+ ≐ -	Alternating polarity	CV	Constant Voltage
	Biphasic pulse shape, symmetrical	mÂ	mA peak
	Biphasic pulse shape, asymmetrical	∧	Volt peak

Frequency sweep mode

12s/12s



1s/5s - 1s/5s



6s/6s



1s/1s


3.11.2 Ultrasound therapy

10%

US duty cycle 10%

1:10 ms

US on : period time 10%



20%

US duty cycle 20%

2:10 ms

US on : period time 20%



30%

US duty cycle 30%

3:10 ms

US on : period time 30%



40%

US duty cycle 40%

4:10 ms

US on : period time 40%




50%

US duty cycle 50%

5:10 ms

US on : period time 50%



100%

US duty cycle 100%

10:10 ms

US on : period time 100%

 \hat{I}_{set}

Set US intensity

US head, ERA 4 cm² **P_{pk}**

Peak US output power

US head, ERA 1 cm² **W/cm^2**

Unit of the set US intensity

3.11.3 Laser therapy **\bar{P}_{set}**

Set average power

 E_{tot}

Total administered energy

 E_p

Energy per pulse



Wear goggles







Monoprobe





Clusterprobe

3.12 Current shapes



3.12.1 Unidirectional currents







3.12.2 Diadynamic currents




3.12.3 Interferential currents



3.12.4 TENS currents



3.12.5 NMES currents



4 OPERATION

4.1 Therapy selection

You can select a therapy with different keys:


- **Direct access keys:** Select a therapy method. See §4.2.
- **Home symbol:** Gives access to:
 - **Objectives:** Select a therapy on the basis of an objective. See §4.3.1.
 - **Indication:** Select a therapy on the basis of a medical indication. See §4.3.2.
 - **Programme number:** Select a certain programme number. See §4.3.3.
 - **Diagnostic programmes:** Perform a diagnosis, for example to determine the rheobase and the chronaxie, or an S/D curve. See §4.3.4.
 - **Contra indications:** Display an overview with contra indications for the different therapies. See §4.3.5.
- **Memory symbol:** Select a saved therapy. See §4.11.

Besides this, you can change the system settings. See §4.12.

4.2 Selection by the Therapy menu

4.2.1 Electrotherapy

1. Select  to go to the therapy menu.
2. Select the desired current.
3. Select the current shape group.
4. Select the current shape.



4.2.2 Ultrasound therapy

1. Select . The Ultrasound screen appears.

4.2.3 Combination therapy

1. Select . The Combination therapy screen appears.
2. Select the current shape. See §4.7.1.


4.2.4 Laser therapy

1. Select . The Laser therapy screen appears.

4.3 Selection by the Guide menu


4.3.1 Therapy selection via objectives

1. Select  to return to the Start menu.
2. Select **Objectives**.
3. Select **Electrotherapy**, **ET iontophoresis**, **Ultrasound therapy**, **Phonophoresis** or **Laser therapy**.
4. Follow the on-screen options to select the desired treatment.



4.3.2 Therapy selection via indications

1. Select  to return to the Start menu.
2. Select **Indications**.
3. Use Δ and ∇ to select the following indications. See §9.1.2.
4. Select the desired indication.
 -  : Electrotherapy
 -  : Ultrasound therapy
 -  : Combination therapy
 -  : Iontophoresis
 -  : Laser therapy




5. With selection via indication list you can view the placement.
 - Select **Electrode placement (ET, CO)**, **US head placement (US, CO)**, **Treatment method (IO)** or **Laser probe placement (LA)**.
 - If necessary, select the location. You get advice to place the electrodes, US head and laser probe.
 - If available, select a number for the precise anatomic location. See §8.3.

Pasirinkus terapiją per indikacijų sąrašą, galima pasižiūrėti elektrodų išdėstymą.
10 tech. reikalavimų punktus.

4.3.3 Programme number selection


1. Select  to return to the Start menu.
2. Select **Program number**.
3. Select the desired programme with Δ or ∇ . See §9.1.
4. Select 1.



4.3.4 Diagnostic programme selection


With the diagnostic programmes, you can localise / treat pain points, perform an S/D curve and look for stress fractures, etc.

1. Select  to return to the Start menu.
2. Select **Diagnostic programs**.
3. Select the desired diagnosis. See §4.9.



4.3.5 Contra indication selection


1. Select  to return to the Start menu.
2. Select **Contra indications**.
3. Select the therapy for which you want to see the contra indications.
4. Scroll through the text with Δ or ∇ .



4.4 Performing therapy

4.4.1 Set and start therapy

1. Select  to return to the Start menu.
2. Select the desired menu item until the treatment appears.
3. Select the desired parameters.
4. Set the treatment time as follows: select treatment time once to set the minutes, select treatment time twice to set the seconds.
5. Change the value of the parameter with Δ and ∇ . The setting range of the parameter is shown at the bottom of the screen. You can change the parameter as long as the parameter has a black background.
6. Start the therapy:
 - Electrotherapy or ultrasound therapy: Rotate intensity knob A or B to start the treatment and to set the desired intensity. The set intensity is displayed in the screen.
 - Laser therapy: See §4.8.1.



4.4.2 Set channels A and B

The Combi 200L, Combi 200, and Duo 200 have two separated electrotherapy channels A and B. The only restriction is that both are in the CC mode or the CV mode.

The channels A and B can be used independently. You can treat two different indications simultaneously with two different therapies.

1. Press  for 3 seconds. The **System settings** menu appears. See §§4.12.
2. Select ∇ to find copy parameters.
3. If necessary, change the parameter **Copy channel parameters** to **OFF**.
4. The selected channel has a black background. If desired, select A or B to change the channel.
5. Select a direct access key. Select the desired treatment. See §4.1.
6. Set the parameters for the first channel. See §4.4.1.
7. Select A or B to change the other channel.
8. Select the desired treatment for second channel. See §4.1.
9. Set the parameters for the second channel. See §4.4.1.

Both channels are selected simultaneously and automatically in case of:

- 4-pole current shapes
- Alternating channels choice with NMES currents (expert mode)
- Combination therapy

Copy channel

On the second channel, you can set the same parameters for electrotherapy as for the first set channel.


1. Press  for 3 seconds. The **System settings** menu appears. See §§4.12.
2. Select  to find copy parameters.
3. If necessary, change the parameter **Copy channel parameters** to **ON**.
4. Select a direct access key. Select the desired treatment. See §4.1.
5. Set the parameters for the first channel A. See §4.4.1.
6. Select B to select the other channel. The treatment including the settings are copied to the other channel.
7. If desired you can change the parameters or the treatment of the selected channel.

Clear channel

1. Make sure that the intensity is set to zero.
2. Select A or B according to the channel that you want to clear.
3. Select . The channel is cleared.

4.4.3 Opening the intensity screen

1. Set and start the treatment. See §4.4.1.
2. Select A or B. Select . The intensity screen appears. The left part of the screen shows channel A. The right part of the screen shows channel B.
3. Press  to return to the settings menu.



4.4.4 Temporary interruption of treatment

1. If a channel has to pause, select this channel.
2. Select  during the treatment. The treatment time of the selected channel is stopped. **Pause** appears on the screen. The parameter settings are retained.
3. Select  to restart the treatment. The intensity now increases gradually to the set level and the treatment time continues again.

4.4.5 Immediately stop treatment

1. Select . All active treatments are stopped immediately. **Stop** appears on the screen. The parameter settings are retained.
2. Set the intensity of the channel again to continue the treatment.


4.5 Electrotherapy

4.5.1 Performing electrotherapy with electrodes

1. Select the desired electrotherapy programme.
2. Place the electrodes. *See page . Placing rubber electrodes and page . Placing the adhesive electrodes.* With Indication list treatments, the **Electrode placement** parameter becomes available.
3. Rotate intensity knob A or B to start the electrotherapy and to set the desired intensity. *See §4.4.1.*
4. Check the patient's reaction. Repeat this check regularly during the treatment.
5. The equipment stops the treatment and indicates that the treatment is completed. Remove the electrodes.

Placing rubber electrodes

1. Moisten two EL sponges thoroughly.
2. In case of poor conduction, use water with a saline solution to improve the conductivity of the EL sponges.
3. Slide a rubber electrode into each sponge.
4. Place the sponges on the part of the body that must be treated.
5. Fasten the sponges to the part of the body with the elastic fixation straps.
6. Connect the rubber electrode with the red connector to the red connector of the two-ply electrode cable.
7. Connect the rubber electrode with the black connector to the black connector of the two-ply electrode cable.
8. Connect the two-ply cable to connector A or B of the Combi 200L, Combi 200, or Duo 200.




Placing the adhesive electrodes



Do not use adhesive electrodes with currents that have a galvanic component, such as galvanic, diadynamic, MF rectangular, pulsed rectangular and triangular currents. These currents will cause skin burns.

1. If possible, disinfect and, if necessary, shave the parts of the body where the adhesive electrodes are to be placed.
2. Place the electrodes on the part of the body that must be treated.
3. Connect the connectors of the adhesive electrodes to the adapter cables.
4. Connect the adapter cables to the two-ply electrode cable.
5. Connect the two-ply electrode cable to connector A or B of the Combi 200L, Combi 200, or Duo 200.




4.5.2 **Perform electrotherapy with vaginal, anal or rectal stimulation probe**



- Considering the very personal and intimate character of these treatments, a probe may only be used for one patient.
- Never disinfect the probes in an autoclave. The probes can be damaged by the extreme temperature.

1. Clean the probe carefully with soap and water.
2. Select the desired electrotherapy programme.
3. Connect the probe to the Combi 200L, Combi 200, or Duo 200. The vaginal and anal probes are immediately detected by the equipment. To prevent unpleasant stimulations, you can only set alternating currents with a Constant Voltage (CV) setting, such as TENS, NMES, and 2-pole interferential currents.



The rectal stimulation probe is not detected by the equipment. With a rectal stimulation probe, select only alternating currents with a Constant Voltage (CV) setting, such as TENS, NMES, and 2-pole interferential currents. This prevents skin burns and unpleasant stimulations.

4. Apply an antiseptic lubricant to the probe.
5. Place the probe.
6. Rotate intensity knob A or B to start the treatment and to set the desired intensity.
7. Check the patient's reaction. Repeat this check regularly during the treatment.
8. The equipment stops the treatment and indicates that the treatment is completed. Remove the stimulation probe.
9. Clean the stimulation probe. See §5.2.6.

4.5.3 **Electrotherapy with sequential steps**

A treatment with sequential steps consists of a succession of the same current form, but additional with different parameter settings. When the treatment is active, you can set the time and the stimulation beep between steps.

Advantages

Electrotherapy with sequential steps has several advantages:


- In one electrotherapy, you can realise several objectives.
- In a treatment with one objective, you can place different accents in the objective.
- You can distinguish between different phases in a treatment, for example preparation, core effect and cooling.

Set new intensity between sequential steps

The intensity determines the peak value during the treatment. When changing to a following step, the intensity is retained if safety allows. Sometimes, it is necessary to increase the intensity for the following step. If the intensity cannot be maintained for safety reasons, the intensity returns to zero. In this case, the treatment is stopped. You must now set the intensity again.

Setting a treatment with sequential steps

1. Select a treatment whereby you can set sequential steps, for example with **Guide menu**, **Program number**, Select no. **230**.
2. Set the **Step time** and **Stimulation beep** parameters for the start of every individual step. Select **Sequence step number** to select a different step.
3. Rotate intensity knob A or B to start the treatment and to set the desired intensity.



Skip step in treatment

1. Select  to temporarily interrupt the treatment.
2. Select **Sequence step number** and select the desired step.
3. Rotate intensity knob A or B to continue the treatment again and to set the desired intensity.

4.5.4 Performing iontophoresis

With iontophoresis, medicines are administered to the body as electrically charged particles (ions) by means of a direct current. To do this, the **Medium Frequency Rectangular current** (= IO-MF constant) is used.


1. Apply the medicine on a sterile gauze. See §8.1. Care must be taken in administering medicine (allergies, contra indications, ...).
2. Place the gauze on the electrode. Make sure that the polarity corresponds with the medicine used.
3. Place the electrodes.
4. Select an **Iontophoresis** therapy programme.
5. **Set the intensity between 0.1 and 0.25 mA/cm²**. The intensity depends on the surface area of the electrodes. With electrodes of 6 x 8 cm (=48 cm²), the current setting must be between 4.8 and 12 mA. Always discuss the intensity with your patient.



To prevent etching or burns, never exceed 0.25 mA/cm².

4.5.5 Read-out values for electrotherapy

1. Channel
2. Electrotherapy
3. Current shape
4. Remaining treatment time
5. Polarity
6. Present intensity
7. Graphical representation of intensity
8. Progress of current



Progress of current

With NMES currents and 4-pole current shapes, the progress of the current is graphically displayed. This gives a clear insight into the phase in which the current is at that moment. In this way, you can optimally guide the patient during the exercise. With the simultaneous application of two NMES currents, the current is only graphically displayed in the intensity screen.

Select the desired channel to open the intensity screen.

4.5.6 Parameters for electrotherapy

The following parameters are given alphabetically. The setting range or the selection possibilities of the parameters depends on the treatment chosen.

Active rest (s)

The duration of the rest period. During the rest period, a low frequency current is applied to stimulate the recovery process.

Alternating channels

The NMES current alternates between channel A and B.

Burst (Hz)

The frequency of the biphasic pulses. The burst consists of a series of pulses that is repeated several times per second. Each burst consists of a low frequency current with high internal pulse frequency (70 - 100 Hz) and a long pulse duration (100 - 250 msec).

Carrier wave (kHz)

The carrier wave frequency, expressed as the number of cycles per second. The frequency of this medium frequency current corresponds with the cycle duration. A high frequency results in a short pulse duration. A carrier wave frequency of 2 kHz is suited for muscle stimulation.

CC / CV

Constant Current (CC) or Constant Voltage (CV).



- When using a dynamic electrode technique, only use alternating currents with Constant Voltage (CV). This prevents unpleasant stimulations for the patient when the contact is temporarily interrupted during the placement, movement and removal of the electrode.
- With a rectal stimulation probe, select only alternating currents with Constant Voltage (CV), such as TENS, NMES, and 2-pole interferential currents. This prevents skin etching and unpleasant stimulations. The rectal stimulation probe is not detected by the equipment.

Characteristics of Constant Current:

- The voltage increases with an increasing load impedance (a worsening contact).
- Within the stated limits, a variation in the load impedance has hardly any effect on the current.
- Without a load, the voltage will go to the maximum level within a short time. After this, an error message will appear on the screen and the current will be switched off.

Characteristics of Constant Voltage:

- With a decreasing load impedance, the current increases.
- Without a load, the output voltage is equal to the set value.
- With a short circuit, the output current in mA is equal to the set voltage in V.

Electrode placing

Instructions for placing the electrodes. Only available with treatment selection via **Indication list**.

Frequency min./max. (base/top) (Hz)

The minimum and maximum frequency of the current cycles, expressed as the number of cycles per second. Within the set sweep mode, the frequency changes within these limits. During the treatment, frequency modulation is desired to prevent habituation. It is recommended to select a fairly low minimum frequency for this (< 20%).

Isodynamic (on, off), applies with Diadynamic currents only

LP and CP use two phases: MF and DF. The MF phase is more intense than the DF phase. If the patient is very sensitive, this difference in perception can be adjusted with this parameter.

On: Reduce the amplitude of the MF phase by 12.5%.

Off time (off) (s)

The interval between two series of current pulses.

On2 amplitude (%)

The amplitude of the pulses during the **On2** period. This amplitude can be set as a percentage of the set amplitude during the **On** period.

On2 frequency (Hz)

The frequency of the pulses during the **On2** time.

On time (on) (s)

The time that the series of current pulses is switched on.

Polarity

The polarity of the current pulse.

Polarity change (on, off)

Switch polarity between red+ and red- during the treatment.

Pulse pause (ms or s)

The duration between the current pulses.

Pulse shape

The shape of the electrical pulse. See §3.11.

Pulse time (μs , ms or s)

The duration of the current pulse.

Rest amplitude (%)

The amplitude of the pulses that is maintained during the active rest period. The active rest period stimulates recovery, which is otherwise realised by the "Off time". The amplitude during the active rest period is set as a percentage of the amplitude during the "On time".

Rest frequency (Hz)

The frequency that is maintained during the active rest period of the NMES current.

Rotation angle (0 - 355°)

The actual angle between the line with the maximum amplitude and the line between the electrodes of channel B. If **Manual** is selected for **Rotation mode**, you can let this angle rotate step by step. This makes it possible to localise deeper treatment points.

Rotation mode (manual, auto)

The maximum amplitude is present at one line in the rotation field (with 100% modulation depth).

- **Auto:** The line with maximum amplitude and 100% modulation depth automatically rotates 360° through the interference field during the set rotation time.
- **Manual:** Position this line manually in the interference field. You do not need to move the electrodes for this.

Rotation time (0 - 20 s)

The time in which the line with maximum amplitude and 100% modulation depth rotates 360° through the interference field. Use a short rotation time (3 - 5 s) to prevent habituation. Use a long rotation time (10-15 s) to localise deeper treatment points.

Segment angle (0 - 45°)

With the segment angle, a certain segment can be stimulated. The segment angle can be set when the **Rotation angle** is set to **Manual**.

Segment time (s)

The time in which the rotation angle changes within the set segment angle.

Sequence step number (1 - 5)

The number of the sequential step that is activated. See §4.5.3.

Step time (mm:ss)

The time in which the selected sequential step number is performed.

Stimulation beep (on, off)

Switch stimulation beep on or off.

Sweep mode (increase, hold, decrease time)

This parameter is only available if **Frequency min. (base)** deviates from **Frequency max. (top)**. The frequency cycle consists of four steps with variable set values: increase, hold, decrease and hold. During the treatment, frequency modulation is desired to prevent habituation.

Total sequence steps

The maximum number of sequential steps. See §4.5.3.

Treatment method

Treatment method for iontophoresis. Always available with treatment selection via **Indication list**.

Treatment time (mm:ss)

The duration of the treatment.


4.6 Ultrasound therapy

4.6.1 Performing ultrasound therapy



Move the US head evenly over the skin during the treatment. This prevents internal burns.

1. Connect the US head into one of the two connectors (☞ of the Combi 200L, Combi 200, or Pulson 200. You can connect two US heads, but only one US head can be in operation at one time. The device detects which US head is connected to the connector (☞ .
2. Select the desired ultrasound therapy. With **Indication list** treatments, the parameter **Head placement** is available.
3. Set the parameter **ERA** to 1 or 4 cm². The corresponding US head is selected, the green indication led on the US head is on.
4. Apply contact gel to the skin to be treated and to the US head.
5. Place the head on the skin.
6. Rotate intensity knob A or B to start the ultrasound therapy.
7. Move the US head evenly over the skin during the treatment. This prevents internal burns.
8. Check the patient's reaction and the effect of the treatment. Repeat this check regularly during the treatment.
9. The equipment stops the treatment and indicates that the treatment is completed.




4.6.2 Phonophoresis

Phonophoresis is used to enhance transdermal transport of several drugs, especially anti-inflammatory NSAID and local anesthetics.

1. Use the drugs (gel ointment) instead of the US contact gel.
2. Select  to return to the **Start** menu.
3. Select **Objectives**.
4. Select **Phonophoresis**. The frequency is 1 MHz, the duty cycle is 50% and the time is at least 7:30 minutes.

4.6.3 Read-out values for ultrasound therapy

1. Channel
2. Ultrasound therapy
3. Type of US head
4. Remaining treatment time
5. Contact detection
6. \hat{I}_{set}
7. P_{pk}



Test the US head if its conduction is bad. See §5.1.1.

Contact detection

The bargraph represents the contact detection level.

- ▬▬▬▬ No bars filled = no contact
- ▬▬▬▬ Partially or all bars filled = sufficient/very good contact

\hat{I}_{set} (W/cm²)

The power (W) of the US head per cm².

P_{pk} (W)

The peak power of the US head ($\hat{I}_{set} * ERA$). The peak power delivered therefore depends on the size of the US head and the contact with the skin. This value is 0.0 W if the contact with the skin is bad. In this case, the ultrasound treatment of the equipment is stopped to prevent overheating of the transducer.

4.6.4 Parameters for ultrasound therapy

Treatment time (mm:ss)

The duration of the treatment.

Duty cycle (10, 20, 30, 40, 50%, continuous)

Ratio of the pulse duration to the period duration.

- **Continuous:** Continuous ultrasound (100%).
- 10, 20, 30, 40, 50%: Pulsating ultrasound.

Select a high duty cycle for an intensive treatment. Select a low duty cycle for a mild treatment.

ERA (cm²)

The effective radiating area expressed in cm² of the treatment head connected. This area equals the cross-sectional area of the beam at the treatment surface. The ERA depends on the frequency. This parameter remains empty if no US head is connected.

Head placement

Instructions for placing the US head. This is only available with treatment selection via **Indication list**.

US frequency (MHz)

The frequency of the US head. The absorption at a US frequency of 3 MHz is three times higher and the penetration depth is three times less than at a US frequency of 1 MHz. Use 3 MHz for superficial tissue and 1 MHz for deeper tissue.

4.6.5 Indicator light of the US head

The indicator light of the US head provides the following information.

Indication light

Blinking blue
Continuous blue
No light

Situation

Bad contact of the US head with the skin.
The US-emission is in progress.
End of the treatment, the US head is not active


4.7 Combination therapy

4.7.1 Performing combination therapy



- With combination therapy, the US head is always the negative pole. The electrode is the positive pole.
- With combination therapy, a maximum current density of 2.0 mA_{rms}/cm² is advised. Exceeding this current density can result in skin irritation and burns. The intensity depends on the surface area of the US head. For US 214 (4 cm²), the current setting may be a maximum of 8 mA_{rms}; for US 211 (1 cm²), a maximum of 2 mA_{rms}.

1. Select 
2. Select **Combination therapy**.
3. Select the current shape.
4. Connect the two-ply electrode to the connector A and connect the US-head to a US-connector.
5. Place the electrode which is connected to the red plug of A on the patient. See page . *Placing rubber electrodes* and page . *Placing the adhesive electrodes*.
6. Apply contact gel to the skin to be treated and to the US head.
7. Place the head on the skin.
8. Rotate intensity knob A to start the electrotherapy. Set the desired voltage.
9. Rotate intensity knob B to start the ultrasound therapy
10. Check the contact between the US head and the skin. The following indications can indicate a bad contact:
 - The treatment stops.
 - The peak power of the ultrasound treatment goes to 0.0 Watt.
11. Check the patient's reaction and the effect of the treatment. Repeat this check regularly during the treatment.
12. The equipment stops the treatment and indicates that the treatment is completed.




4.8 Laser therapy

4.8.1 Performing laser therapy




Start a laser therapy only when all persons in the room wear laser goggles for eye protection.

1. Make sure all persons wear laser goggles.
2. Plug the connector of the laser probe into the connector  of the Combi 200L.
3. Put the laser lock in the ON position (horizontal).
4. Select the desired laser therapy. The green indicator light on the laser probe lights up. With **Indication list** treatments, the parameter **Laser probe placement** is available.
5. Place the laser probe on the skin.
6. Press the black knob on the laser probe to start the laser therapy. The yellow indicator light on the laser probe lights up. Hold the knob to keep the laser probe active during the treatment.
7. Check the patient's reaction and the effect of the treatment. Repeat this check regularly during the treatment.
8. The equipment stops the treatment and indicates that the treatment is completed.
9. Release the black knob on the laser probe.
10. Put the laser lock in the OFF position (vertical).



4.8.2 Read-out values

1. Channel
2. Laser therapy
3. Type of connected laser probe
4. Remaining treatment time
5. P_{set}
6. E_{tot} (during therapy) or E_p (during laser energy measurement)



E_p (μJ)

The energy per pulse (μJ).

 P_{set} (μW or mW)

The set average power (μW or mW) of the laser probe (E_p * frequency).

 E_{tot} (mJ or J)

The total administered energy (mJ or J) of the current treatment (P_{set} * treatment time).

4.8.3 Parameters**Probetype**

The type of laser probe: monoprobe or clusterprobe.

Frequency (Hz or kHz)

The pulse repetition frequency of the laser beam.

Repeat mode (on, off)

Off: The E_{tot} value is reset to zero after the treatment time ends.

On: The E_{tot} value is kept after the treatment time ends. The E_{tot} value is the sum of the total administered energy from the performed treatments.

Treatment time (mm:ss)

The duration of the treatment.

Laser probe placement

Instructions for placing the laser probe. This is only available with treatment selection via **Indication list**.

4.8.4 Indicator lights on the laser probe

The indicator lights on the laser probe provide the following information.

Indication light

Continuous green

Continuous yellow

Situation

The laser therapy is selected, but the laser probe has no laser emission.

The laser emission is in progress.

4.8.5 Testing laser emission

1. Set a laser therapy. See §4.8.1.
2. Place the laser probe output perpendicular on the laser test eye $\star\text{---}\tau$.
3. Press and hold the black knob on the laser probe during the laser test. The laser test symbol $\star\text{---}$ appears on the read-out screen.
4. Make sure the E_{tot} value increases every second by the P_{set} value.
5. Release the black knob on the laser probe.
6. Select  to return to the **Start** menu.
7. Put the laser lock in the OFF position (vertical).

You can also test the energy per pulse of the laser probe. See §4.12.8.

4.9 Diagnostic programmes

With the diagnostic programmes, you can investigate the state of the electrical sensitivity of the neuro-muscular system:


- Rheobase and chronaxie. See §4.9.1.
- Rheobase and AQ. See §4.9.2.
- Determine a S-D curve. See §4.9.3.

Besides this, there are programmes for localisation:


- Pain points. See §4.9.4.
- Stress fracture search.

4.9.1 Determining rheobase and chronaxie

1. Select  to return to the **Start** menu.
2. Select **Diagnostic programs**.
3. Select **Rheobase and chronaxie**.
4. If desired, change the **Polarity** and **Stimulation beep** settings.
5. Rotate intensity knob A to start the treatment. The set intensity is displayed in the screen.
6. Increase the intensity in steps of 0.1 mA , until you observe a tangible or visible contraction.



7. Confirm pulse amplitude with . The measured rheobase (in mA) is saved.
8. The equipment now doubles the rheobase (mÂ). The pulse time changes to 0.1 ms. Increase the pulse time by Δ , until you observe a tangible or visible contraction.
9. Confirm pulse time with . The chronaxie (in ms) is saved. The results screen appears.
10. If desired, select  to save the data in the memory. See §4.11.1.




4.9.2 Determining Rheobase and Accomodation Quotient (AQ)


1. Select  to return to the Start menu.
2. Select **Diagnostic programs**.
3. Select **Rheobase and AQ**.
4. Determine the rheobase as with **Rheobase and chronaxie**. See §4.9.1.
5. Select **Confirm pulse amplitude**. The measured rheobase is saved.
6. The equipment now selects a triangular pulse. Increase the intensity in steps of 0.1 mÂ, until you observe a tangible or visible contraction.
7. Select **Confirm pulse amplitude**. The measured AQ is saved. The results screen appears.
8. If desired, select  to save the data in the memory. See §4.11.1.

4.9.3 Determine a S-D curve

1. Select  to return to the Start menu.
2. Select **Diagnostic programs**.
3. Select **S-D curve rectangular, S-D curve triangular or S-D curve rect. + tri..**
4. If desired, change the **Recording mode, Polarity** and **Stimulation beep** settings. If **Manual** is selected for the **Recording mode**, you can skip or repeat a measurement with Δ and ∇ .
5. Select **Start recording S-D**.



6. Rotate intensity knob A to start the treatment.
7. Increase the intensity in steps of 0.1 mA , until you observe a tangible and visible contraction.
8. Select **Confirm**.
9. Repeat steps 7 and 8 for all the measurements.
10. When END is marked, the S-D curve is finished. Depending on the measurement the Optimal pulse time (OPT), Rheobase (Rh), Chronaxie (Ch) and Accommodation Quotient (AQ) results are shown.
11. If desired, select  to save the data in the memory. See §4.11.1.



4.9.4 Pain points

1. Select  to return to the **Start** menu.
2. Select **Diagnostic programs**.
3. Select **Pain points**.
4. Select the diagnostic programme for pain points.

4.10 Vacuum treatment

4.10.1 Performing electrotherapy with Vaco 200

1. Connect the Vaco 200 and the electrotherapy device to the mains supply.
2. Connect the input connector A of the Vaco 200 to the output connector A of the electrotherapy device with the connection cable. If necessary, do the same for channel B.
3. Turn the vacuum intensity regulator to the  position.
4. Use the switch on the rear panel to set the Vaco 200 and the electrotherapy device to on. The LEDs from the traditional channels of the Vaco 200 light up. The LEDs show which output channel is connected with the electrotherapy device.
5. Select the desired electrotherapy programme.
6. Connect, adjust and place the electrodes. *See §4.10.2*
7. Turn the intensity knob A or B of the ET device to start the electrotherapy and to set the desired intensity.
8. Check the patient's reaction. Repeat this check regularly during the treatment.
9. The electrotherapy device stops the treatment and indicates that the treatment is completed.
10. If vacuum is used, turn the vacuum intensity regulator to . For safety, the vacuum stops automatically two minutes after the treatment stops and the vacuum channel LEDs flash. Turn the vacuum intensity regulator to  to make a new treatment selection possible.
11. Remove the electrodes. If necessary put your finger under the rim of the vacuum electrode for air inled.


4.10.2 Connecting the electrodes



- Use always demineralized water with vacuum electrodes to avoid lime deposits in the watertank, tubes and sponges. Add a saline solution to improve the electrical conduction.
- Use moist sponges only. Too dry sponges can cause a bad electrical contact and burn the skin.
- Do not use vacuum electrodes with DC current. The DC current causes damage to the vacuum cups by ionization.


Using four vacuum electrodes

1. If present disconnect the vaginal, anal and rectal stimulation probes.
2. Connect the vacuum electrodes to the vacuum hoses.
3. Connect the four vacuum hoses. Select two cables with the same hose colour for each channel.
 1. Connect the red connectors from the vacuum hoses to the output connectors with the red dot.
 2. Connect the black connectors from the vacuum hoses to the output connectors with the black dot.
4. Moisten the round sponges.
5. Put the sponges in the vacuum electrodes.
6. Turn the vacuum intensity regulator to  and regulate the desired suction force.
7. Place the vacuum electrodes on the part of the body that must be treated. The vacuum electrodes stay in place by the underpressure. Too high suction cause patient discomfort.
8. Turn the pulse vacuum suction power regulator to  for patient comfort.
9. Select the stimulation signal on the vacuum output B with the push button. The LED vacuum channel B is lit.




Using two vacuum electrodes

1. If present disconnect the vaginal, anal and rectal stimulation probes.
2. Connect two vacuum electrodes to the vacuum hoses. Select two cables with the same hose colour.
3. Connect the red connector from the vacuum hose to the channel A output connector with the red dot.
4. Connect the black connector from the vacuum hose to the channel A output connector with the black dot.
5. Connect a vacuum hose for a closed vacuum system to make a short-circuit across the channel B output connectors.
6. Moisten the round sponges.
7. Put the sponges in the vacuum electrodes.
8. Turn the vacuum intensity regulator to  and regulate the desired suction force.
9. Place the vacuum electrodes on the part of the body that must be treated. The vacuum electrodes stay in place by the underpressure. Too high suction cause patient discomfort.
10. Turn the pulse vacuum suction power regulator to  for the patient comfort.



Using traditional electrodes


1. Prepare the traditional electrodes. *See the User Manual of the ET device.*
2. Connect the two-ply cable to the corresponding output connector A or B of the Vaco 200.



Using a combination of vacuum and traditional electrodes

Use channel A for the stimulated vacuum electrodes and use channel B for the traditional electrodes.

1. Connect two vacuum electrodes to channel A and short-circuit channel B. See §4.10.2.
2. Prepare the traditional electrodes. See the User Manual of the ET device.
3. Connect the two-ply cable to the output connector B of the Vaco 200.



4.10.3 The water reservoir is full


1. The LED of the water reservoir  light up. Finish the present treatment. After you turn the vacuum off, you can not restart the vacuum.
2. Empty the water reservoir with the drain hose. The start of a vacuum treatment is possible.
3. Clean the water reservoir. See §5.2.9.

4.11 Memory


You can save 50 of your own programmes for later use: programmes 500 up to and including 549. You can modify these programmes for much-used or specific current shapes for a certain patient.

4.11.1 Saving a programme

1. Select a therapy. See §4.1.
2. Change the settings for the patient. See §4.4.
3. Select .
4. Select **Save**.
5. Select a free programme number or overwrite an existing programme number.
If desired, scroll through the list with Δ and ∇ .




6. Enter the name of the programme. Use the name or the number of the patient, for example.
 - Select a character with Δ and ∇ .
 - Select **Move cursor to the left/right** to change the cursor position.
7. Select Ready and save.



4.11.2 Selecting a saved programme


Selecting a programme by the list

1. Select .
2. Select **Open by name list**.
3. Select the desired programme. If necessary, scroll through the list with Δ or ∇ .




Selecting a programme by the number

1. Select .
2. Select **Open by number**.
3. Select the desired programme with Δ or ∇ .
4. Select **Go to selected number**.



4.11.3 Clearing a saved programme


1. Select .
2. Select **Clear program**.
3. Select the desired programme.
If necessary, scroll through the list with Δ or ∇ .
4. Select  to delete the programme.



4.11.4 Editing a standard programme

Standard programmes have a programme number that is lower than 50. You can only edit standard programmes with the therapy keys.

1. Select a programme with the direct access keys.
2. Change the desired parameter settings.
3. Select .
4. Select **Change default program**.
5. Select  to edit the standard programme.



You can also save an edited standard programme under a free programme number. See §4.11.1.


You can reset the standard settings of the standard programmes with Reset Menu. See §4.12.2.

4.12 System settings

With the system settings, you can adapt the Standard settings of the equipment. You cannot change the system settings during a therapy.

4.12.1 Changing the system settings

1. Press  for 3 seconds.
2. Change the desired system setting.



4.12.2 Parameters

Contrast (1 - 20)

The contrast of the display.

Language

The language selection: select the language with which the read-out must work.

Sound settings

Sound settings. See §4.12.3.

Stand-by time (5, 10, 15, 20 minutes, off)

If the device is not used during the stand-by time, the device goes to the stand-by mode. Press any key to reactive the device.

Text for start up screen

The text that appears in the top of the start up screen, after the equipment is switched on. See §4.12.5.

First screen (guide menu, therapy menu)

The first screen you see when activating the device.

Copy channel parameters (on, off)

Choose channel A and B the same or different is set by the copy channel parameter. See §4.4.2.

System information

System information of the equipment

Always have this information available when you contact the technical service department.

Plate electrode test

Test the condition of the rubber electrodes. See §4.12.7.

Cable test

Test the cables. See §4.12.6.

Laser energy measurement

Test the laser probe. See §4.12.8.

Error history

The total number of error reports that the equipment has had and details about the last 10 error reports.

Always have this information available when you contact the technical service department.

Counter working hours (hours, minutes, sec.)

The time that the accessories for electrotherapy or ultrasound therapy have been in use. For this, the output of the channel must have been higher than zero.

Reset menu


- **Reset working hours:** Set the number of working hours of a plate electrode or an US head to zero.
- **Restore program 1-50:** All programme changes will be put back to factory settings.
- **Clear total memory:** Restores the standard settings of the standard programmes and of the edited programmes.

Stop time if bad US (on, off)

When there is a bad US contact, the treatment time counter stops. When the contact is restored, the counting continues.

4.12.3 Setting the sound

1. Press  for 3 seconds.
2. Select **Sound settings**.
3. Change the desired sound setting.



4.12.4 Parameter sound settings

End of treatment

On: A sound signal will be heard at the end of the treatment.

Pressing a key

On: A sound signal will be heard every time a key is pressed.

ET stimulation

On: A sound signal will be heard at each pulse of the electrotherapy.

Beep volume (min.1, standard 5, max.10)

The volume of the sound signals.

ET bad contact

On: A sound signal will be heard if the electrode does not make good contact with the skin.

US bad contact

On: A sound signal will be heard if the US head does not make good contact with the skin.


Laser output

On: A sound signal will be heard at each pulse of the laser therapy.

4.12.5 Set text for start up screen

You can set your own text for the start up screen. For example, you can put your name or address information here.

1. Press  for 3 seconds.
2. Select **Text for start up screen**.
3. Enter the text for the start up screen.
 - Select a character with Δ and ∇ .
 - Select **Move cursor to the left/right** to change the cursor position.
4. Select  to save the changes.



4.12.6 Cable test

1. Press  for 3 seconds.
2. Select **Cable test**.
3. Connect the electrode cable to channel A without the electrodes.
4. Connect the test plug to the connectors of the cable.
5. Set the amplitude to 20 mA with rotary knob A.
6. If the cables function correctly, the following message will appear **Condition cable: OK**.
7. Turn the amplitude back to 0 mA. Select  to exit the window.

4.12.7 Plate electrode test

1. Press  for 3 seconds.
2. Select **Plate electrode test**.
3. Connect the electrode cable to channel A with the electrodes.
4. Place the electrodes on each other, without the sponges. Make sure that the electrodes make contact over the whole surface.
5. Set the amplitude to 20 mA with rotary knob A.
6. If the electrodes function correctly, the following message will appear **Condition electrodes: OK**.
7. Turn the amplitude back to 0 mA.

4.12.8 Laser energy measurement

1. Perform the laser energy measurement on a 'cold' (not recently used) laser probe for a reliable test.
2. Make sure all persons wear laser goggles.
3. Plug the connector of the laser probe into the connector  of the Combi 200L.
4. Put the laser lock in the ON position (horizontal).
5. Press  for 3 seconds. The **System settings** menu appears.
6. Select **Laser energy measurement**.
7. Test the monoprobe or clusterprobe.

Test the monoprobe

1. Place the laser probe output perpendicular on the laser test eye .
2. Press and hold the black knob on the laser probe during the laser test. The laser test symbol  appears on the read-out screen.
3. Move the probe a bit to-and-fro to obtain the maximum value.
4. Release the black knob on the laser probe.
5. Make sure the measured E_p value correspond within $\pm 20\%$ with the E_p value of the supplier control report of the laser probe.
6. Select  to go back to the **System settings** screen.
7. Put the laser lock in the OFF position (vertical).

Test the clusterprobe

1. Place the laser probe output perpendicular over the laser test eye . Place the laser probe against the intensity knob A and B to position the laser diode over the laser test eye.
2. Press and hold the black knob on the laser probe during the laser test. The laser test symbol  appears on the read-out screen.
3. Rotate the probe a bit to-and-fro to obtain the maximum value.
4. Release the black knob on the laser probe.
5. Repeat the measurement for the other laser diodes.
6. Calculate the sum of the four measured E_p value.
7. Make sure the E_p sum value correspond within $\pm 20\%$ with the total E_p value of the supplier control report of the laser probe.
8. Select  to go back to the **System settings** screen.
9. Put the laser lock in the OFF position (vertical).

5 INSPECTIONS AND MAINTENANCE

5.1 Inspections

Component	Visual check or test	Frequency
Electrode cables and electrodes	Damage Insulation intact	At least 1x per month
US head	Dents, cracks or other damage	At least 1x per month
	Test US head. See §5.1.1.	With bad operation or at least 1x per year
Cable of US head	Damage Pins in connector straight	At least 1x per month
Laser probe	Dents, cracks or other damage	At least 1x per month
	Test the laser probe. See §4.8.5. and §4.12.8.	Every day
Cable of laser probe	Damage Insulation intact Pins in connector straight	At least 1x per month
Equipment	Technical safety inspection. See §5.1.2.	At least 1x per year

5.1.1 US head test

Test the US head if its conduction is bad. This is the case when the indication bar for the Ppk value displays .

1. Select an ultrasound therapy.
2. Place the US head in a bowl with water.
3. Rotate intensity knob A or B to start the treatment.
4. Check in the screen whether the Ppk value of the corresponding channel is increasing.
5. If the Ppk value does not increase, contact your local Gymna dealer.

5.1.2 **Technical safety inspection**

The 'Directive on Medical Devices' from the European Commission (93/42/EEG) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. If the legislation in your country or your insurer prescribes a shorter period, you must adhere to this shorter period.



- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- The inspection may only be performed by a suitably qualified person. In some countries this means that the person must be accredited.

Inspection points

The technical safety inspection contains the following tests:

1. Test 1: General: Visual inspection and check on the operating functions
2. Test 2: Electrotherapy
3. Test 3: Ultrasound therapy
4. Test 4: Laser therapy
5. Test 5: Electrical safety inspection: measurement of the earth leakage current and patient leakage current according to DIN/VDE 0751-1 ed. 2.0.

Inspection result

1. A registration must be maintained of the technical safety inspections. Use the inspection report in the appendix for this purpose. See §8.5.
2. Copy this appendix.
3. Complete the copied appendix.
4. Keep the inspection reports for at least 10 years.

The inspection is successful if all inspection items are passed.

Repair all faults on the equipment before the equipment is put back into operation.

By comparing the registered measurement values with previous measurements, a possible slowly-deteriorating deviation can be ascertained.

5.2 Maintenance

Component	Check	Frequency
Main unit	Cleaning. See §5.2.1.	According to requirement
Rubber electrodes	Cleaning. See §5.2.2.	After every treatment
Vacuum electrodes	Cleaning. See §5.2.3.	After every treatment
ET and vacuum unit sponges	Cleaning. See §5.2.4.	After every treatment
Fixation bandages	Cleaning. See §5.2.5.	If necessary
Vaginal, anal and rectal stimulation probe	Cleaning and disinfecting. See §5.2.6.	After each use
US head	Cleaning. See §5.2.7.	After each use
Laser probe	Cleaning. See §5.2.8.	After each use
Vacuum hoses and water reservoir	Cleaning. See §5.2.9.	Weekly



Accessories that come in contact with the body of the patient must be washed with pure water after the disinfection to prevent allergic reactions.

5.2.1 *Cleaning the 200 Series main unit and the vacuum unit*



Do not sterilize the main unit!

1. Remove dust with a dry cloth.
2. If necessary, remove stains or dirt with a damp cloth.
3. If required, clean the device with a non aggressive soap solution, a 70% alcohol solution or any other agent suitable for surface disinfection that does not harm the material (cover: ASA, case: lacquered sheet metal). In case of doubt, please consult an accepted list of disinfectants (e.g. the VAH list of disinfectants).



Do not use chloride based agents as these may affect the plastic parts of the device.

5.2.2 Cleaning the rubber electrodes

1. Clean the electrodes in a non-aggressive soap solution or in a 70% alcohol solution.
2. Rinse the electrodes thoroughly with water.
3. Dry the electrodes.

5.2.3 Cleaning the vacuum electrodes

1. Clean the vacuum electrodes (metal electrodes and rubber suction cups) in a non-aggressive soap solution or in a 70% alcohol solution.
2. Rinse the vacuum electrodes thoroughly with water.
3. If present, remove calcium deposits.
4. Turn the cups inside out.
5. Check for dirt and calcium deposits. Remove it if present.
6. Dry the vacuum electrodes.

5.2.4 Cleaning the EL sponges and vacuum unit sponges

1. Rinse the EL sponges thoroughly with water or clean the EL sponges with a 70% alcohol solution.
2. Rinse the EL sponges thoroughly with water.
3. Let the sponges dry.

5.2.5 Cleaning the fixation bandages

1. Clean the fixation bandages in a 70% alcohol solution or another disinfectant.
2. Rinse the fixation bandages in water.
3. Let the fixation straps dry.

5.2.6 *Cleaning and disinfecting vaginal, anal and rectal stimulation probes*



- Considering the very personal and intimate character of these treatments, a probe may only be used for one patient.
- Never disinfect the probes in an autoclave. The probes can be damaged by the extreme temperature.

Immediately after every treatment

1. Clean the probe carefully with soap and water.
2. Place the probe in an HAC solution of 1% or in a 70% alcohol solution for at least 30 minutes.



- Read the instruction leaflet in the packaging of the HAC.
- Make sure that the probe connector does not get into the HAC solution.

3. Dry the probe with a clean towel.
4. Store the probe in a plastic bag that is provided with the name of the patient.

Before reusing the probe:

1. Clean the probe carefully with soap and water.
2. Apply an antiseptic lubricant to the probe. See §4.5.2.

5.2.7 *Cleaning the US head*

1. Clean the US head with a lightly moistened soft cloth.
2. Disinfect the treatment surface with a cotton bud that is soaked in a 10% HAC solution.
3. Rinse the US head thoroughly with clean water.

5.2.8 *Cleaning the laser probe*



- The laser probe is not waterproof.
- Do not scratch the aperture pane.

1. Clean the laser probe with a lightly moistened soft cloth.
2. Disinfect the treatment surface with a cotton bud that is soaked in a 10% HAC solution.

5.2.9 *Cleaning the vacuum hoses and the water reservoir*

1. Empty the water reservoir with the drain hose.
2. Prepare a recipient with max 150 ml of a 70% alcohol solution.
3. Connect the vacuum hoses.
4. Put the ends of the vacuum hoses in a 70% alcohol solution.
5. Turn the vacuum unit on.
6. Suck the liquid up until the LED of the water reservoir lits.
7. Turn the vacuum unit off.
8. Empty the water reservoir.
9. Repeat steps 4 until 8 with pure water.

6 MALFUNCTIONS, SERVICE AND GUARANTEE

6.1 Malfunctions

Component	Problem	Solution
200-series	Equipment cannot be switched on	See §6.1.1.
	Equipment does not react to commands or a fault report appears	See §6.1.3.
	Foreign language on the screen	Change the language. See §4.12.2.
EL sponges	Furring	Replace the sponges
	Bad conduction	Replace the sponges
Vacuum electrodes	Contamination by ionization	See §6.1.4.

6.1.1 *Equipment cannot be switched on*

1. Check if the mains voltage has failed.
2. Check if the main switch is switched on ("I").
3. Check if the power cord and the fuses are in order. If necessary, replace the fuse. See §6.1.2.
4. Contact your dealer if the equipment still cannot be switched on.

6.1.2 *Replacing a fuse*

1. Switch the main switch off ("O").
2. Unplug the power cord from the equipment.
3. Pull the fuse holder carefully out of the equipment. If necessary, use a screwdriver.
4. Replace the fuse. If necessary, order new fuses from your dealer.
5. Install the fuse holder and plug in the power cord.
6. Switch the main switch on again ("I").

6.1.3 *Equipment does not react to commands or an error message appears*

The safety system of the equipment has ascertained a fault. You cannot continue to work. An instruction usually appears on the screen.

1. Disconnect the connection to the patient.
2. Switch the main switch off ("O").
3. Wait 5 seconds and switch the main switch on again ("I").
4. Contact your dealer if the equipment still does not react to commands.

6.1.4 Remove the contamination from the vacuum electrodes

1. Clean the vacuum electrodes. See §5.2.3.
2. Use steel wool or sandpaper for metal with fine grains ('P 400' or higher) to remove the contamination.
3. Replace the vacuum electrodes if the contamination is still present.

6.2 Service



- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories to perform repairs. The equipment does not contain any components that may be replaced by the user.
- If possible, open the screen with the system settings before you contact the technical service department. See §4.12

Service and guarantee are provided by your local GymnaUniphy dealer. The conditions of delivery of your local GymnaUniphy dealer apply. If you have qualified technical personnel that are authorised by GymnaUniphy to perform repairs, your dealer can provide diagrams, spare parts lists, calibration instructions, spare parts and other information on request, for a fee.

6.3 Guarantee

GymnaUniphy and the local GymnaUniphy dealer declares itself to be solely responsible for the correct operation when:

- all repairs, modifications, extensions or adjustments are performed by authorised people;
- the electrical installation of the relevant area meets the applicable legal regulations;
- the equipment is only used by suitably qualified people, according to these user instructions;
- the equipment is used for the purpose for which it is designed;
- maintenance of the device is regularly performed in the way prescribed. See §5.;
- the technical life time of the equipment and the accessories is not exceeded;
- the legal regulations with regard to the use of the equipment have been observed.

The guarantee period for the equipment is 2 (two) years, beginning on the date of purchase. The date on the purchase invoice acts as proof. This guarantee covers all material and production faults. Consumables, such as sponges, adhesive electrodes and rubber electrodes, do not fall under this guarantee period.

This guarantee does not apply to the repair of defects that are caused:

- by incorrect use of the equipment,
- by an incorrect interpretation or not accurately following these user instructions,
- by carelessness or misuse,
- as a consequence of maintenance or repairs performed by people or organisations that are not authorised to do so by the manufacturer.

6.4 Technical life time

As far as possible GymnaUniphy will supply service, spare parts and accessories for a period of 10 years from the date of manufacture. See the type plate for this information.

7 TECHNICAL INFORMATION

7.1 General

Dimensions of 200-series (w x h x d)	267 x 296 x 125 mm
Weight of 200-series	3,650 kg
Weight including accessories	4,6 kg
Mains voltage	100 - 240 VAC, 50 - 60 Hz
Maximum power, in operation	85 VA
Safety class	Class I (earthed socket required)
Insulation	Type BF (floating patient circuit)
Fuses	2 x T2AL250V
Metal housing	Lacquered sheet metal
Top cover	Asa

7.2 General for Vaco 200

Dimensions Vaco 200 (w x h x d)	267 x 95 x 270 mm
Weight Vaco 200	3,5 kg
Weight including accessories	4,6 kg
Mains voltage	100 - 240 VAC, 50 - 60 Hz
Maximum power consumption	30 VA
Safety class	Class I (earthed socket required)
Insulation	Type BF (floating patient circuit)
Fuses	2 x T2AL250V
Volume water reservoir	± 180 ml
Working pressure continuous vacuum	38 - 320 hPa
Working pressure pulsation vacuum	58 - 480 hPa
Vacuum rhythm	1,5/1,5 - 1,5/4,5 s (on/off time)

7.3 Electrotherapy

7.3.1 General

Gydymo trukmė 0-60 min.
8 tech. reikalavimų punktas

Treatment time	0 - 60 min.
Current limitation	The smallest value: <ul style="list-style-type: none">- 150% of the set value, or:- 110% of the maximum for the selected current shape
Accuracy	Set current value mA at 500 Ω - typically \pm 10%
CC/CV mode	For all current shapes, with the exception of medium frequency rectangular current
Polarity	Red-, red+ and alternating polarity, if applicable

7.3.2 Current shapes

Medium frequency rectangular current

Intensity 0 - 80 mA with 300 to 1000 Ω

Iontophoresis MF rectangular

Electrode area 6 - 300 cm²

Intensity 0 - 80 mA with 300 to 1000 Ω

Rectangular pulsed current, Triangular pulsed current, 2-5 Current (Ultra Reiz)

Pulse time 0,1 ms - 6 s

Pulse pause 1 ms - 6 s

Intensity of CC 0 - 80 mA with 300 to 1000 Ω

Intensity of CV 0 - 80 V_{pk} with I < 80 mA

MF, DF, CP, LP

Intensity of CC 0 - 80 mA with 300 to 1000 Ω

Intensity of CV 0 - 80 V_{pk} with I < 80 mA

Expert parameters:

MF time 1 - 100 s

DF time 1 - 100 s

ISO on / off

Conventional TENS, Low frequency TENS

Pulse time 10 - 650 μ s

Pulse shape symmetrical, asymmetrical

Frequency min. (base) 1 - 150 Hz

Frequency max. (top) 1 - 150 Hz

Freq. increase time 0 - 100 s

Freq. hold time 0 - 100 s

Freq. decrease time 0 - 100 s

Intensity of CC 0 - 120 mA with 300 to 1000 Ω

Intensity of CV 0 - 120 V_{pk} with I < 120 mA

Random frequency TENS

See TENS currents, with the exception of:

Pulse frequency 1 - 150 Hz, with automatic stochastic frequency variation of +/-35% maximum

Burst TENS

See TENS currents, with the exception of:

Pulse frequency 20 - 150 Hz

Burst frequency 1 -10 Hz

Rectangular surge current, Triangular surge current

Pulse time	0,1 - 5 ms
Pulse frequency	1 - 150 Hz
Intensity of CC	0 - 80 mA with 300 to 1000 Ω
Intensity of CV	0 - 80 V _{pk} with I < 80 mA

Biphasic surge current, Biphasic surge intrapulse interval

(with a fixed interval between positive and negative pulses of 100 μ s)

Pulse time	10 - 650 μ s
Pulse frequency	1 - 150 Hz
Pulse shape	symmetrical, asymmetrical (only for Biphasic surge current)
Intensity of CC	0 - 120 mA with 300 to 1000 Ω
Intensity of CV	0 - 120 V _{pk} with I < 120 mA

2-pole medium frequency surge current, Isoplanar vector field surge

Carrier wave frequency	2 - 10 kHz
AM frequency	1 - 200 Hz
Intensity of CC	0 - 100 mA with 300 to 1000 Ω
Intensity of CV	0 - 100 V _{pk} with I < 100 mA

Russian stimulation

Carrier wave frequency	2 - 10 kHz
Burst frequency	20 - 100 Hz
Intensity of CC	0 - 100 mA with 300 to 1000 Ω
Intensity of CV	0 - 100 V _{pk} with I < 100 mA

Expert parameters for NMES currents

On time (ON)	1 - 100 s
Off time (OFF)	0 - 100 s
Rest time	0 - 100 s
Surge time	0 - 100 s
Shrink time	0 - 100 s
Special modes	OFF, REST, ON2, Frequency sweep, Manual stimulation
Alternating channels	ON/OFF (not for Isoplanar vector field surge current)
On2 amplitude	1 - 100%
Rest amplitude	1 - 100%

2-pole medium frequency current, Isoplanar vector field

Carrier wave frequency	2 - 10 kHz
AM frequency min.	0 - 200 Hz
AM frequency max.	0 - 400 Hz
Frequency variation mode	0/1/0, 1/5/1, 6/0/6, 12/0/12
Intensity of CC	0 - 100 mA with 300 to 1000 Ω
Intensity of CV	0 - 100 V _{pk} with I < 100 mA

Dipole vector field

See 2-pole medium frequency current

Rotation time	0 - 20 s
Rotation angle	0 - 355°
Segment angle	0 - $\pm 45^\circ$
Segment time	0 - 10 s

Diagnostic programmes: Rheobase and Chronaxie, Rheobase and AQ, S-D curve rectangular, S-D curve triangular, S-D curve rect. + tri.

Intensity of CC	0 - 80 mA with 300 to 1000 Ω , with Rheobase max. 40 mA
Variable parameter for determination Chronaxie:	
Pulse time	0,1 - 100 ms
Variable parameter for determination S-D curves:	
Pulse time	17 fixed steps between 0,1 - 1000 ms
Recording mode	auto / manual

7.4 Ultrasound therapy**7.4.1 General**

Insulation classification	Type BF
Peak power	0 - 2 W/cm ² , duty cycle = 100% 0 - 3 W/cm ² , duty cycle < 100%
Accuracy of intensity	$\pm 10\%$ of maximum at set values above 10% of this maximum
Treatment time	0 - 30 min.
Deviation of time clock	< 0,5%
Modulation frequency	100 Hz
Modulation type	CW (rectangular on/off)
Repetition period of pulses	10 ms

7.4.2 Modulation and pulse duration

Modulation duty cycle	100	50	40	30	20	10	%
Pulse time	∞	5	4	3	2	1	ms
Ratio of $p_{tm} - p$	1	2	2,50	3,33	5	10	

7.4.3 US heads

US head, model US214			
Acoustic operating frequency	1,0	3,2	MHz
Output power	7,6	8,6	W
Effective intensity of output	2,0	2,0	W/cm ²
Effective Radiating Area (ERA)	3,8	4,3	cm ²
Beam Non-uniform Ratio (BNR)	4,5	7,0	
Maximum intensity of beam	9,0	14,0	W/cm ²
Beam type	Convergent	Collimated	

US head, model US211			
Acoustic operating frequency	1,0	3,2	MHz
Output power	1,6	1,8	W
Effective intensity of output	2,0	2,0	W/cm ²
Effective Radiating Area (ERA)	0,8	0,9	cm ²
Beam Non-uniform Ratio (BNR)	6,9	4,0	
Maximum intensity of beam	13,8	8,0	W/cm ²
Beam type	Divergent	Collimated	

7.5 Laser therapy

7.5.1 General

Safety class 3B laser product

7.5.2 Monoprobe: model LA201

Number of laser diodes	1
Nominal ocular hazard distance	214 mm
Wave length	905 nm
Energy per pulse	2,35 μ J
Peak performance	13,5 W
Maximum average power	70,5 mW
Pulse frequency	2-30000 Hz
Pulse width at 50% of the peak power	155 ns
Beam surface at laser aperture	12,9 mm ²
Beam divergence	Dual mode 10° and 45°

7.5.3 Clusterprobe: model LA204

Number of laser diodes	4
Nominal ocular hazard distance	95 mm
Wave length	904 nm
Total energy per pulse	10,8 μ J
Peak performance	4 x 18 W
Maximum average power	54 mW
Pulse frequency	2-5000 Hz
Pulse width at 50% of the peak power	145 ns
Beam surface at laser aperture	4 x 5,3 mm ²
Composite beam divergence	21°

7.6 Environmental conditions

Temperature	+10 °C to +40 °C
Relative humidity	30% to 75%
Atmospheric pressure	700 hPa to 1060 hPa

7.7 Transport and storage

Transport weight	5,6 kg
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 100%, including condensation
Atmospheric pressure	200 hPa to 1060 hPa
Transport classification	Single piece, by post

The transport and storage specifications apply to equipment in the original packaging.

7.8 Standard accessories

Article numbers can change in the course of time. Check the article numbers in the most recent catalogue or ask your dealer. The drawings are merely indicative, no rights can be derived from them.

7.8.1 General

	Quantity	Description	Art. no.
	1	Power cord ¹	100.689
	1	VAS score card	115.684
	1	CD-ROM User manual Gymna 200-series multi language	EN: 376126
	1	Safety instructions	323.011
	1	Quick start manual Gymna 200 series	376.134

¹ This power cord has a CEE 7/7 type plug. For countries with other outlets, a different power cord with the appropriate plug is available as an option.

7.8.2 *Standard accessories electrotherapy*

	Quantity	Description	Art. no.
	2	2-pole patient cable, 2 mm fiche (2)	340.670
	2	Rubber electrode, 6 x 8 cm, 2 mm fiche (4)	340.468
	1	Chamex bag, 6 x 8 cm (4)	100.658
	4	Fixing strap, elastic, 5 x 60 cm (4)	108.935
	1	Test plug F/F - 2 mm	330.803

7.8.3 *Standard accessories ultrasound therapy*

	Quantity	Description	Art. no.
	1	US head, 1&3 MHz - ERA 4 cm ² incl. holder	320.114
	1	Contact gel, 500 ml	341.088
	1	Power cord ¹	100.689

¹ This power cord has a CEE 7/7 type plug. For countries with other outlets, a different power cord with the appropriate plug is available as an option.

7.8.4 Standard accessories for vacuum therapy

	Quantity	Description	Art. no.
	1	Power cord ¹	100.689
	2	Connection cable: ET device - Vaco	102.032
	1	Vacuum hose dark grey (per 2 pcs: black/red connector)	340.615
	1	Vacuum hose light grey (per 2 pcs: black/red connector)	340.604
	2	Vacuum electrode - Ø 60 mm (per 2 pcs, 2x)	340.626
	1	Sponge for vacuum electrode - Ø 60 mm (per 4 pcs)	340.648

¹ This power cord has a CEE 7/7 type plug. For countries with other outlets, a different power cord with the appropriate plug is available as an option.

7.9 Optional accessories electrotherapy

	Quantity	Description	Art. no.
	1	Vaginal probe Novatys gold	329.978
	1	Anal probe Analys+	330.561
	1	Rectal stimulation probe	112.166
	2	Rubber electrode no. 1 - 4 x 6 cm	340.446
	2	Rubber electrode no. 3 - 8 x 12 cm	340.481
	4	EL sponge no. 1 for electrode 4 x 6 cm	100.657
	4	EL sponge no. 3 for electrode 8 x 12 cm	100.659
	4	Adhesive electrode, 3 cm diameter	326.799
	4	Adhesive electrode, 2,5 x 5 cm	326.810
	4	Adhesive electrode, 5 x 5 cm	326.821
	4	Adhesive electrode, 5 x 10 cm	326.832
	1	Pin electrode 15 mm diameter with grip and sponge	114.142
	10	EL sponges for pin electrode	109.944

Advice: Replace the electrode material at least every 6 months.

7.10 Optional accessories ultrasound therapy

	Quantity	Description	Art. no.
	1	US head, 1&3 MHz - ERA 1 cm ² , incl. holder	320.111
	1	Contact gel, can 5 l	341.088
	1	Pump for can, 5 l	341.121

7.11 Optional accessories laser therapy

	Quantity	Description	Art. no.
	1	Monoprobe, model LA201, incl. holder	320.101
	1	Clusterprobe, model LA204, incl. holder	320.104
	1	Laser goggles EP8-6	339.592
	1	Remote interlock for laser	116.227

7.12 Optional accessories vacuum therapy

	Quantity	Description	Art. no.
	2	Vacuum electrode - 90mm	114.686
	4	Sponge for vacuum electrode - 60mm	114.687

Article numbers can change in the course of time. Check the article numbers in the most recent catalogue or ask your dealer.
The drawings are merely indicative, no rights can be derived from them.


8 APPENDICES

8.1 Agents for iontophoresis

Agent	Property	Application and form
Calcium (+)	Analgeticum and sedative	Application: post-traumatic pain, distorsion, algodystrophic syndromes and neuralgia. Form: 2% calcium chloride solution.
Magnesium (+)	Analgeticum and fibrolyticum	Applications as with calcium. 10% magnesium chloride solution.
Iodine (-)	Sclerolyticum	Application: stubborn scars, cutaneous adhesences, sickness of Dupuytren, stiffness of joints and adhesive capsulitis. Form: 1-2% potassium iodine solution
Salicylate (-)	Anti-inflammation agent	Application: periphlebitis, osteoarthritis, ab-articular rheumatism, articular stiffness and adhesive capsulitis. Form: 2% sodium salicylate solution.
Procaine and lidocaine (+)	Anti-inflammation agent	Application: production of local anaesthesia, in the neuralgia of the trigeminal nerve, e.g. with acute inflammation. Form: 2% solution.
Histamine (+)	Revulsive and vasodilator	Application: degenerative and articular rheumatic pains, such as cramp. Maximum duration of iontophoresis: 3 min. Longer treatment causes allergic reactions and cephalgia. Form: 0,02% bicarbonate solution.
ColtramyI (+)	Myorelaxant	Application: contractures. Form: solutions up to 0,04%. 2 ml coltramyI (4mg/ ampoule), to be dissolved in 8 ml distilled water.
Indocid (-)	A.I.N.S.	Application: inflammatory illnesses. Form: 1% solution. 50 mg freeze-dried powder, to be dissolved in 5 ml distilled water.
Voltaren (-)	A.I.N.S.	Application: inflammatory illnesses. Form: 0,75% solution. 3 ml (75 mg/ampoule), to be dissolved in 7 ml distilled water.
Acetic acid	A.I.N.S.	Application: To dissolve deposition layers caused by ossifying myositis and periarticular ossification. Form: 2% water solution.

8.2 Diagnostic S-D curve

Physiotherapist:		Date of investigation:	
Name of patient:		Date of birth:	M/F
Anamnesis:			
Evaluation (neuro-muscular):		Accommodation Quotient:	
Rheobase:	mA	Chronaxie:	ms
Conclusion:			
Treatment:			



8.3 Electrode, US head and laser probe placements

Select the therapy via indication list to get information about the placement.
See §4.3.2.

8.3.1 Electrotherapy

Select the **Electrode placement** parameter to show the optimal location for the placement of the electrodes.

The numbers in the illustration give information to the precise anatomic location. Select the numbers with the blue keys. The description of the location is often explained with the abbreviations:



pnp	peripheral nerve point	snp	skin nerve point
mnp	motor nerve point	mtp	myofascial trigger point
n	nerve	nn	nervi
m	muscle	mm	musculi
r	ramus	rr	rami

The types of nerves are represented in a different way:


1. Peripheral nerve, deep
2. Fascia
3. Skin nerve

Pasirinkite elektrodų išdėstymo parametra, kad pamatyti optimalias elektrodų pozicijas.
10 tech. reikalavimų punktą

1. Peripheral nerve, superficial
2. Peripheral nerve, deep



1. Skin nerve on the point of the fascia



Other information in the illustrations:

- The electrodes shown on the front of the body are black.
- The electrodes shown on the rear of the body are transparent.
- The type of electrodes to use is not advised.
- The size of the shown electrodes is an indication for the advised size.
- The letters A and B give a recommendation for the choice of channel.
- The symbols + and - recommends the polarity.

8.3.2 *Iontophoresis*

The **electrode placement** parameter is replaced by the parameter **treatment method**. The **treatment method** shows the iontophoresis method on screen.

8.3.3 *Ultrasound therapy*

Select the **US head placement** parameter to show the optimal location for the placement of the US head.

You can select the numbers in the illustration with the blue keys for more information.


- 1 Gives information on the precise anatomic location.
- 2 Numbers with a black background gives specific recommendations.



Relevant bone structures are shown for detailed information on the treated area. The number of points below the US head gives an indication of the dimensions of the treated area. The information in the illustration recommends a treatment technique. This illustration shown an example of the dynamic technique.



If other areas are possible for the US head placement a black area is shown. Select the corresponding number 2 for information on the screen. If the area is on the rear a transparent area is shown.



8.3.4 Combination Therapy

The **US head/electr. placement** parameter for combination therapy shows the US head placement. The electrode is not shown in the illustration. Place the electrode near to the US head.

8.3.5 Laser therapy

Select the **laser probe placement** parameter to show the optimal location for the placement of the laser probe.

8.4 EMC directive

Use only cables, electrodes and US heads that are specified in this manual. See §7. The use of other accessories can have a negative effect on the electromagnetic compatibility of the equipment.

If you use the 200-series in the vicinity of other equipment, you must check that the 200-series is functioning normally.

The following paragraphs contain information about the EMC properties of the equipment. Because this information is intended for technicians, the information is given in English.

8.4.1 Guidance and declarations

Guidance and manufacturer's declaration - electromagnetic emissions		
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The 200-series devices use RF energy only for their internal function. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The 200-series devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
	Class B	
Harmonic emissions IEC 61000-3-3	Class B	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact / ±8 kV air No loss of performance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV power / ±1 kV I/O No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV diff. / ±2 kV comm. No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	U_T - 100% (0,5 period) No loss of performance U_T - 60% (5 periods) No loss of performance U_T - 30% (25 periods) No loss of performance U_T - 100% (5 seconds) Device resets to a safe state. (60601-1 § 49.2)	Mains power quality should be that of a typical commercial or hospital environment. If the user of a 200-series device requires continued operation during power mains interruptions, it is recommended that the 200-series device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level			

Guidance and manufacturer's declaration - electromagnetic immunity			
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} AM 1 kHz 80% 150 kHz to 80 MHz	10 V.....0,15-80 Mhz 51 V.....6,78 Mhz 54 V.....13,56 Mhz 50 V.....27,12 Mhz 45 V.....40,68 Mhz	Portable and mobile RF communications equipment should be used no closer to any part of a 200-series device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,35\sqrt{P}$ $d = 0,07\sqrt{P}$ $d = 0,06\sqrt{P}$ $d = 0,07\sqrt{P}$ $d = 0,08\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m AM 1 kHz 80% 80 MHz to 2,5 GHz	10 V/m..0,08-1,0 Ghz 26 V/m....1,4-2,0 Ghz 30 V/m...433,92 Mhz 30 V/m.....915 Mhz	$d = 0,35\sqrt{P}$ 80 MHz to 800 MHz $d = 0,70\sqrt{P}$ 800 MHz to 2,5 GHz $d = 0,12\sqrt{P}$ $d = 0,23\sqrt{P}$
Radiated RF ENV 50204	3 V/m CW 200 Hz d.c. 50% 895 MHz to 905 MHz	30 V/m.895-905 Mhz	$d = 0,23\sqrt{P}$
<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 			
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 The guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey can be considered. If the measured field strength in the location in which a 200-series device is used exceeds the applicable RF compliance level above, the 200-series devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 200-series device.</p>			
<p>b Over the frequency range 150 kHz to 80 MHz, field strengths must be less than 10 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the 200-series device

The 200-series device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of a 200-series device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 200-series devices as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0,35\sqrt{P}$	80 MHz to 800 MHz $d = 0,35\sqrt{P}$	800 MHz to 2,5 GHz $d = 0,70\sqrt{P}$
0,01	0,04	0,04	0,07
0,1	0,11	0,11	0,22
1	0,35	0,35	0,70
10	1,11	1,11	2,21
100	3,50	3,50	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8.5 Technical safety inspection

200-series with serial number is / is not¹ in good working order		
Location:	Inspection performed by: Name:	Owner: Name:
Date:	Initials:	Initials:

¹ Cross out what does not apply.

If a specific test does not apply to this equipment, place a mark in the NA (not applicable) column.

8.5.1 Test 1: General

	Yes	No	NA
1. The results of earlier safety inspections are available.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2. The logbook is present.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3. The type plate and the supplier's label are legible.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4. The housing, adjusting knobs, keys and display are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5. The power connection and power cord are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6. The output connectors are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7. The electrode connectors and cables are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8. The cables and connectors of the US head(s) are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
9. The US head(s) do not display any cracks or other damage that can endanger the insulation.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10. The automatic self-test at switch-on does not give an error message.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
11. The display does not show any defective points or lines.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8.5.2 Test 2: Electrotherapy

	Yes	No
1. Connect loads of 500 Ω to both normal electrode pairs. Connect an oscilloscope to these pairs (black to ground).		
2. Select channel A, program 4: MF constant.		
3. At maximum intensity, the output currents correspond within 10% with the values on the display.	<input type="checkbox"/>	<input type="checkbox"/>
4. The output signals correspond with figure 1.	<input type="checkbox"/>	<input type="checkbox"/>
5. The polarity changes to negative if "RED(-)" is selected.	<input type="checkbox"/>	<input type="checkbox"/>
6. The warning "Bad contact with the patient" is given if the load is disconnected.	<input type="checkbox"/>	<input type="checkbox"/>
7. Select channel B, program 4: MF constant. Select CC.		
8. At maximum intensity, the output currents correspond within 10% with the values on the display.	<input type="checkbox"/>	<input type="checkbox"/>
9. The output signals correspond with figure 1.	<input type="checkbox"/>	<input type="checkbox"/>
10. The polarity changes to negative if "RED(-)" is selected.	<input type="checkbox"/>	<input type="checkbox"/>
11. The warning "Bad contact with the patient" is given if the load is disconnected.	<input type="checkbox"/>	<input type="checkbox"/>
12. Remove the load, so that the unloaded output voltage can be measured.		
13. Select channel A, program 23: 2-pole medium frequency. Select CV.		
14. At maximum intensity, the output voltage corresponds within 10% with the values on the display.	<input type="checkbox"/>	<input type="checkbox"/>
15. The output signals correspond with figures 2 and 3.	<input type="checkbox"/>	<input type="checkbox"/>
16. The yellow lamp next to the output connectors lights if the intensity is not 0.	<input type="checkbox"/>	<input type="checkbox"/>
17. Select channel B, program 23: 2-pole medium frequency. Select CV		
18. At maximum intensity, the output voltage corresponds within 10% with the values on the display.	<input type="checkbox"/>	<input type="checkbox"/>
19. The output signals correspond with figures 2 and 3.	<input type="checkbox"/>	<input type="checkbox"/>
20. The yellow lamp next to the output connectors lights if the intensity is not 0.	<input type="checkbox"/>	<input type="checkbox"/>

Figure 1




Figure 2





Figure 3



8.5.3 Test 3: Ultrasound

	Yes	No
1. Connect the treatment head and place it in an ultrasound measurement device. Select an ultrasound therapy.		
2. Select 1 MHz, continuous (duty cycle 100%), 2 W/cm ² The measured value is within $\pm 20\%$ of the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
3. Select 1 MHz, duty cycle 50%, 3 W/cm ² The measured value is within $\pm 20\%$ of half the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
4. Select 3 MHz, continuous (duty cycle 100%), 2 W/cm ² The measured value is within $\pm 20\%$ of the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
5. Select 3 MHz, duty cycle 50%, 3 W/cm ² The measured value is $\pm 20\%$ of half the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
6. Select 3 MHz, duty cycle 50%, 0.5 W/cm ² With a dry treatment surface, the Ppk value becomes 0.	<input type="checkbox"/>	<input type="checkbox"/>
7. Select 1 MHz, duty cycle 50%, 0.5 W/cm ² With a dry treatment surface, the Ppk value becomes 0.	<input type="checkbox"/>	<input type="checkbox"/>

The maximum power transfer takes place at the operating frequencies. If the equipment does not function at the correct frequency, this results in a too low output power. It is therefore not necessary to check the operating frequencies.

8.5.4 Test 4: Laser therapy



Start a laser therapy only when all persons in the room wear laser goggles for eye protection.

Use for test A and B a laser radiation measurement device with the following specifications:

- The resolution of the measured energy per pulse value is: $\leq 0,1 \mu\text{J}$.
- The wavelength range is: 900 - 910 nm.
- The capability to measure: 200 ns pulses of $30 W_{pk}$.
- Capable of capturing a divergent beam with a diameter: ≥ 10 mm.
- Tolerance: $\leq 10\%$.

Test A: The monoprobe

	Yes	No
1. Connect the monoprobe to the Combi 200L. See §4.8.1.		
2. Select a laser therapy. The green indicator light lights up.	<input type="checkbox"/>	<input type="checkbox"/>

Test A: The monoprobe**Yes No**

- | | | | |
|----|--|--------------------------|--------------------------|
| 3. | Press the black knob on the laser probe. The yellow indicator light lights up and the green indicator light goes out. | <input type="checkbox"/> | <input type="checkbox"/> |
| 4. | Release the black knob. The green indicator light lights up and the yellow indicator light goes out. | <input type="checkbox"/> | <input type="checkbox"/> |
| 5. | Start the laser therapy to measure the E_p value with the laser radiation measurement device.
The measured E_p value is μ J.
Stop the laser therapy. | | |
| 6. | The measured E_p value corresponds within $\pm 20\%$ with the E_p value of the test protocol of the laser probe. | <input type="checkbox"/> | <input type="checkbox"/> |

Test B: The clusterprobe**Yes No**

- | | | | |
|----|---|--------------------------|--------------------------|
| 1. | Connect the clusterprobe to the Combi 200L. See §4.8.1. | | |
| 2. | Select a laser therapy. The green indicator light lights up. | <input type="checkbox"/> | <input type="checkbox"/> |
| 3. | Press the black knob on the laser probe. The yellow indicator light lights up and the green indicator light goes out. | <input type="checkbox"/> | <input type="checkbox"/> |
| 4. | Release the black knob. The green indicator light lights up and the yellow indicator light goes out. | <input type="checkbox"/> | <input type="checkbox"/> |
| 5. | Start the laser therapy to measure the E_p values with the laser radiation measurement device.
The measured E_p value from laser diode 1 is μ J.
The measured E_p value from laser diode 2 is μ J.
The measured E_p value from laser diode 3 is μ J.
The measured E_p value from laser diode 4 is μ J.
Stop the laser therapy.
The sum of the four measured E_p value is μ J. | | |
| 6. | The sum of the measured E_p values corresponds within $\pm 20\%$ with the total E_p value of the test protocol of the laser probe. | <input type="checkbox"/> | <input type="checkbox"/> |

Test C: Calibration of the laser test eye**Yes No**

- | | | | |
|----|---|--|--|
| 1. | Connect a calibrated monoprobe to the Combi 200L. See §4.8.1. | | |
| 2. | Select Laser energy measurement in the System settings menu. See §4.12.8. | | |

Test C: Calibration of the laser test eye**Yes No**

3. Place the laser probe output perpendicular on the laser test eye . Start the laser energy measurement. Move the probe a bit to-and-fro to obtain the maximum value. The measured E_p value is μJ . Stop the laser energy measurement.
4. The measured E_p value corresponds within $\pm 5\%$ with the E_p value of the calibrated laser probe.
5. If not, contact the service department of your local dealer.

8.5.5 Test 5: Electrical safety test (IEC 62353:2007)**Yes No**

1. The resistance of the safety earth is less than 0.2Ω
2. The housing leakage current is less than $1000 \mu\text{A}$
3. The patient leakage current is less than $5000 \mu\text{A}$

Notes:

8.6 Disposal

Take account of the following environmental aspects when disposing of the equipment and the accessories:

- The basic device, the cables and the electrodes fall under small chemical waste (electrical and electronic equipment waste). These components contain lead, tin, copper, iron, various other metals and various plastics, etc. Dispose according to national regulations.
- Sponges, sponge bags and gels contain only organic material and do not require any special processing.
- Packaging materials and manuals can be recycled. Deliver them to the appropriate collection points or include them with the normal household waste. This depends on the local organisation of the waste processing.

Notify your dealer about the disposal.

9 REFERENCE

9.1 Function overview

9.1.1 *Therapy menu*

This section is listed in the TIO (Therapy menu, Indications, Objectives) document. This is available as a separate document on the DVD.

9.1.2 *Indications list*

This section is listed in the TIO (Therapy menu, Indications, Objectives) document. This is available as a separate document on the DVD.

9.1.3 *Objectives*

This section is listed in the TIO (Therapy menu, Indications, Objectives) document. This is available as a separate document on the DVD.

9.1.4 *System settings*

Press  for 3 seconds.

Contrast

Language

Sound settings

Stand-by time

Text start up screen

First screen

Copy channel parameters

System information

Plate electrode test

Cable test

Laser energy measurement

Error history

Counter working hours

Reset menu

Stop time if bad US

9.2 Literature

A literature list can be sent on request. Please contact GymnaUniphy.

9.3 Terminology

absolute muscle power: The maximum total tension that a muscle can produce.

accomodation: The ability of the nerve tissue to protect itself against stimulations that slowly increase in strength.

Pulse time	Delay in action potential of rectangular pulse: triangular pulse	Accomodation Quotient (AQ)
500 ms	1:1.5 to 1:3	1,5 - 4
1000 ms	1:2 to 1:6	2 - 6

active trigger point: A point that, with stimulation (pressure, stretch or electrical pulse), besides the local pain also generates a projected pain in the area that the patient is complaining about.

antalgic: The pain is reducing.

atrophy: Deterioration in the nourishment state of organs. As a result, the organs become smaller or shrink.

chronaxie: The time threshold that is required for a muscle contraction or a sensory impression, after the occurrence of the necessary minimum required stimulation.

denervation: Switching-off or weakening of the innervation (paralysis).

durability: Being able to frequently repeat a muscle contraction.

epithelisation: Recovery of the epithelium over the bottom of the wound. A unidirectional current can stimulate the epithelisation. Epithelisation can also be activated by an external electrical stimulation.

explosive muscle power: The highest tension that a muscle can produce in the shortest possible time.

hyperalgesia: An increased sensitiveness for pain. Apply a modified dosage in the case of acute hyperalgesia.

injury current: A small unidirectional current between the epidermis and the corium, which occurs after a wound. This current activates the recovery process. With a slow recovery process, an external unidirectional current can be applied to realise the same effect.

innervation: The effect of the nerves on the working of the muscles or glands.

iontophoresis: The flow of ions through a tissue by means of a galvanic current.

isometric contraction: A muscle contraction whereby the length of the muscle remains constant. The external resistance of the muscle must be at least as large as the power that is generated by the contraction. Under isometric circumstances, especially the tension in the muscle increases and muscle cramp is avoided.

loadability: The (maximum) load that can be carried.

loss of muscle tone: The state of tension of muscles reduces.

Myofascial Trigger Point (MTP): A trigger point that is located in the myofascial tissue. The MTP is located in a hard cord of a muscle. The MTPs can be localised with **Pain points** in the **Diagnostics program**.

Neuro Muscular Electro Stimulation (NMES): Contraction of an innervated muscle or muscle group by means of low or medium frequency electrostimulation. The purpose of NMES is to improve or maintain the movement.

pain threshold: The lowest level of stimulation that causes pain.

pain tolerance threshold: The level of stimulation that can just be tolerated by the patient. The pain tolerance threshold is past the pain threshold.

re-innervation: The restoration of the innervation.

responsiveness: The degree to which a tissue or organ reacts to a stimulation. With a high responsiveness, a mild treatment is desired. With a low responsiveness, a more intensive treatment can be desired. Make a good estimate of the responsiveness to determine the correct dosage.

rheobase: The minimum galvanic current strength required with the stimulation of the nerve to cause a muscle contraction.

sclerolysis: The solution of a hardening of the tissue. The tissue can be chemically and electrically softened with a cathode in combination with chlorine or iodine.

skin etching: Electro-chemical reactions that can be threatening for tissues and organs, especially for the skin. With correct application, a desired effect occurs, for example improvement of the circulation. Skin etching occurs with current shapes that have a direct current component.

slow twitch muscle fibre: Muscle fibres with a low contraction speed. The fibres are fairly thin, produce a small amount of power and have a low fatigue level. See also type I muscle tissue.

tetanic contraction: A persistent muscle contraction, on the basis of several contraction waves that are simultaneously in a muscle. You can cause tetanic contractions with an NMES surge current.

tone: The tension state of tissues.

trophic: The state of nourishment.

type I muscle tissue: Muscle tissue with a low contraction speed.

type II muscle tissue: Muscle tissue with a high contraction speed. Set the parameters as follows for stimulation with NMES:

NMES parameter	type I	type II
Pulse time	Long	Short
Pulse frequency	Low	High
Pulse amplitude	-	High
Series duration and series pause	Short	Long
Treatment time	Long	-

VAS score: Score on the Visual Analogue Scale (VAS). Tool for evaluating a clinical complaint from the patient. This usually concerns the degree to which pain is felt. With a high VAS score, a mild treatment is usually adequate. With a lower VAS score, a more intensive treatment is desired.

Index

A

- Abbreviations 3
- Absolute muscle power 102
- Accessories 80, 82, 84
 - optional electrotherapy 83
 - optional laser therapy 84
 - optional ultrasound therapy 84
 - standard 80
- Accommodation 102
- Accommodation Quotient 50
- Active rest 39
- Active trigger point 102
- Agents for iontophoresis 85
- Alternating channels 39
- Anal stimulation probe
 - cleaning 67
 - perform electrotherapy 36
- Antalgic 102
- Atrophy 102

B

- Burst 39

C

- Cable test 59, 61
- Carrier wave 39
- Channel
 - set 32
- Chronaxie 49, 102
- Cleaning 65
- Combination therapy 46
 - perform 46
 - selection 29
- Connection 15
- Constant Current 39
- Constant Voltage 39
- Contra indication
 - selection 31
- Contra Indications 11
- Contra indications 12
- Contrast 16, 58

- Copy channel parameters 58
- Counter working hours 59
- Current shape groups 25
- Current shapes 25, 28

D

- Denervation 102
- Device
 - cleaning 65
- Diadynamic currents 28
- Diagnostic
 - perform 49
 - program selection 31
- Display 24
 - symbols 25
- Disposal 99
- Durability 102
- Duty cycle 44

E

- EL sponges
 - Cleaning 66
- Electrical safety 8
- Electrode
 - cleaning 66
 - test 61
- Electrode placing 40, 42
- Electrodes
 - connection 53
- Electrotherapy 34
 - optional accessories 83
 - parameters 39
 - perform 34
 - read-out values 38
 - safety 9
 - selection 29
 - sequential steps 36
 - technical information 74
- EMC 8
- EMC directive 89
- Environmental conditions 79
- Ep 48
- Epithelisation 102

ERA 44
Error history 59
Etot 48
Explosive muscle power 102

F

First screen 58
Fixation bandage
 cleaning 66
Frequency 40
Function overview 101
Functional test 15

G

Guarantee 70
Guide menu
 selection 30

H

Hyperalgesia 102

I

Indication list 30
Injury current 103
Innervation 103
Inspections 63
Installation 15
Intensity screen 33
Interferential currents 28
Interruption 34
Iontophoresis 38, 103
Îset 44
Isodynamic 40
Isometric contraction 103

L

Language 16, 58
Laser energy measurement 59, 62
Laser frequency 48
Laser probe
 cleaning 67
 placement 48

Laser therapy
 optional accessories 84
 parameters 48
 perform 47
 read-out values 47
 safety 10
 selection 29
 technical information 79

Liability 13

Loadability 103

Loss of muscle tone 103

M

Maintenance 65
Malfunctions 69
Medical Devices Directive 13
Memory 55
Muscle tissue
 type I and type II 104
Myofascial Trigger Point 103

N

NMES 103
 currents 28

O

Objectives 30
Off time (off) 40
On time 40
On2 amplitude 40
On2 frequency 40

P

Pain points 51
Pain threshold 103
Pain tolerance threshold 103
Parameter symbols
 electrotherapy 26
 laser therapy 27
 ultrasound therapy 27
parameters
 sound settings 60

Placing 15
 rubber electrodes 35
Placing the adhesive electrodes 35
Plate electrode test 59
Polarity 40
Polarity change 40
Ppk 44
Prevention of explosion 8
Probetype 48
Program
 clear 57
 number selection 31
 saving 55
 selecting 56
Pset 48
Pulse pause 40
Pulse shape 40
Pulse time 41
Purpose 7

R

Recall
 list 56
 number 56
Rectal stimulation probe
 cleaning 67
 perform electrotherapy 36
Re-innervation 103
Repeat mode 48
Replacing a fuse 69
Reselling 16
Reset menu 59
Responsiveness 103
Rest amplitude 41
Rest frequency 41
Rheobase 49, 50, 104
Rotation angle 41
Rotation mode 41
Rotation time 41

S

Safety 7
 instructions 7
 technical inspection 64, 94
Sclerolysis 104
S-D curve 50
Segment angle 41
Segment time 41
Selection
 therapy 29
Sequence step number 37, 41
Sequential steps
 treatment 37
Service 70
Skin etching 104
Slow twitch muscle fibres 104
Sound
 Set 60
Sound settings 58
Stand-by time 58
Step time 41
Stimulation beep 42
Stop 34
Stop treatment 34
Storage 16
 conditions 80
Sweep mode 42
System information 59
System setting
 parameters 58
System settings 58
 changing 58

T

Technical information 73
Technical life time 71
TENS currents 28
Terminology 102
Tetanic contraction 104
Text for start up screen 58
Text start up screen
 set 61

Therapy
 program selection 31
 selection via indication list 30
 selection via objectives 30
 set 32
 start 32
Therapy menu
 selection 29
Tone 104
Total sequence steps 42
Touchscreen 15
Transport 16
 conditions 80
Treatment
 interruption 34
 sequential steps 37
 stop 34
Treatment time 42, 44, 48
Trophic 104

U

Ultrasound therapy 43
 optional accessories 84
 parameters 44
 Perform 43
 Read-out values 44
 safety 9
 selection 29
 technical information 77

Unidirectional currents 28
US frequency 45
US head
 cleaning 67
 Indicator light 45
 placement 45
 test 63

V

Vacuum electrode
 cleaning 66
Vacuum electrodes
 place 53
Vacuum hoses
 cleaning 68
Vaginal stimulation probe
 cleaning 67
 perform electrotherapy 36
VAS score 104

W

Water reservoir 55
 cleaning 68





Pasweg 6A
B-3740 Bilzen

Tel.: (+32) (0) 89/510.532
Fax: (+32) (0) 89/510.541

www.gymna.com
E-mail: info@gymna.com

Your dealer:

