IMMULITE® 2000 Multi-Diluent 2

English

Intended Use: For the on-board dilution of patient samples exceeding the calibration range of the indicated IMMULITE® 2000 assays. See the IMMULITE 2000 Operator's Manual for an explanation of on-board dilution features.

Catalog Number: **L2M2Z, L2M2Z4**

Test Code: **MD2**

Use with IMMULITE 2000: AFP (L2KAP); Albumin (L2KHA); Beta-2 Microglobulin (L2KBM); BR-MA (L2KBR); Calcitonin (L2KCL); CEA (L2KCE); CK-MB (L2KMB); EPO (L2KEPN); Ferritin (L2KFE); Gastrin (L2KGA); GI-MA (L2KGI); Growth Hormone (L2KGRH); Total IgE (L2KIE); IL2R (L2KIP); IL-6 (L2K6P); Myoglobin (L2KMY); Nicotine Metabolite (L2KNM); PAPP-A (L2KPC); PSA (L2KPS); Free PSA (L2KPF); TSI (L2KSI); IL-8 (L2K8P); IL-10 (L2KXP).

Materials Supplied

One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

Safety data sheets (SDS) available on [siemens-healthineers.com](https://www.siemens-healthineers.com).

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Multi-Diluent 2 is not intended for use in linearity or sensitivity studies, but rather for the dilution of patient samples that exceed the calibration range of the assay.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1%, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

For technical assistance, contact your National Distributor.

[siemens-healthineers.com](https://www.siemens-healthineers.com)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Shading indicates technical content that differs from the previous version.

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IMMULITE is a trademark of Siemens Healthineers.
All rights reserved.