COVID-19 /Flu A&B /RSV /ADENO/M.P Ag combo rapid Test Package Insert

Cat: 3024201

Specimens: **Nasopharyngeal swab**/Nasal swab
For professional in vitro diagnostic use only.

INTENDED USE

COVID-19 /Flu A&B /RSV /ADENO/MP Ag combo rapid Test is a rapid visual immunoassay for the qualitative, presumptive detection of Influenza A and B viral antigens, COVID-19 antigen, RSV antigen, Adenovirus antigen, and Mycoplasma pneumoniae (MP) antigen from nasopharyngeal swab and nasal swab specimens.

INTRODUCTION

Influenza is a highly contagious, acute viral infection of the respiratory tract. The causative agents of the disease are immunologically diverse, single-stranded RNA viruses known as influenza viruses. There are three types of influenza viruses: A, B, and C. Type A viruses are the most prevalent and are associated with the most serious epidemics. Type B viruses generally produce a milder disease compared to type A. Type C viruses have never been associated with a large epidemic of human disease. Both type A and B viruses can circulate simultaneously, but usually, one type is dominant during a given season. Influenza antigens may be detected in clinical specimens by immunoassay.

Coronaviruses are enveloped RNA viruses distributed broadly among humans, other mammals, and birds, causing respiratory, enteric, hepatic, and neurologic diseases. Six coronavirus species are known to cause human disease. Four viruses—229E, OC43, NL63, and HKU1—are prevalent and typically cause common cold symptoms in immunocompetent individuals. The two other strains—severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV)—are zoonotic in origin and have been linked to sometimes fatal illnesses. Given the high prevalence and wide distribution of coronaviruses, their large genetic diversity, frequent genome recombination, and increasing human-animal interface activities, novel coronaviruses are likely to emerge periodically in humans due to frequent cross-species infections and occasional spillover events.

Respiratory syncytial virus (RSV) is an RNA virus belonging to the Paramyxoviridae family. It is transmitted by air droplets, aerosols, and close contact, with an incubation period of 3-7 days. RSV infection is quite common, with over 93% of cord blood specimens testing positive for anti-RSV IgG antibody by immunofluorescence assay. Because maternal antibodies and serum antibodies cannot completely prevent RSV reinfection, RSV frequently causes regional epidemics. RSV is highly infectious, often leading to family cluster infections, predominantly involving upper respiratory tract infections. In the United States, 20%-25% of infant pneumonia and 50%-75% of bronchiolitis cases are caused by RSV infection.

Adenovirus causes respiratory tract infections, such as pharyngitis/tonsillitis, pharyngoconjunctival fever, and pneumonia, but also diarrhea, epidemic keratoconjunctivitis, and other diseases. Pharyngoconjunctival fever is epidemic in summer, while other adenovirus infections occur nearly year-round. The incidence of adenovirus respiratory tract infections is higher in children.

Mycoplasma pneumoniae (MP) is a common cause of respiratory tract infections, particularly in children and young adults. It is transmitted through respiratory droplets and is known to cause mild to moderate community-acquired pneumonia, as well as other respiratory conditions such as bronchitis and pharyngitis. MP infections often present with gradual onset of symptoms, including fever, sore throat, and cough, and are associated with outbreaks in close-contact settings like schools and households.

The Influenza A+B/COVID-19/RSV/Adeno/MP Antigen Rapid Test Device is a lateral-flow immunoassay using highly sensitive monoclonal antibodies specific for influenza types A and B, COVID-19, RSV, adenovirus, and Mycoplasma pneumoniae (MP) antigens. This test is specific to influenza types A and B, COVID-19, RSV, adenovirus, and MP antigens, with no known cross-reactivity to normal flora or other known respiratory pathogens. The test is intended for use with nasopharyngeal swab and nasal swab specimens.

PRINCIPLE

The COVID-19/Flu A&B/RSV/Adeno/MP Antigen Rapid Test Device detects influenza A and B viral antigens, COVID-19 antigen, RSV antigen, Adenovirus

antigen, and Mycoplasma pneumoniae (MP) antigen through visual interpretation of color development on the strip. Anti-influenza A and B antibodies, Anti-COVID-19 antibodies, Anti-RSV antibodies, Anti-Adeno antibodies, and Anti-MP antibodies are immobilized on the test regions of the membrane, respectively. During testing, the extracted specimen reacts with anti-influenza A, B, COVID-19, RSV, Adeno, and MP antibodies conjugated to colored particles, which are pre-coated onto the sample pad of the test. The mixture then migrates through the membrane by capillary action and interacts with reagents on the membrane. If there are sufficient influenza A and B viral antigens, COVID-19 antigen, RSV antigen, Adeno antigen, or MP antigen in the specimen, colored band(s) will form at the corresponding test region of the membrane. The presence of a colored band in these regions indicates a positive result for the particular antigens, while its absence indicates a negative result. The appearance of a colored band at the control region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS

- Individually packed Test Cassette
- Buffer
- Timer
- MATERIALS PROVIDED
 - Sterile nasopharyngeal swab
 - Package insert

MATERIALS NEEDED BUT NOT PROVIDED

PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged. Do not reuse tests.
- The extraction reagent solution contains a salt solution if the solution contacts the skin or eye, flush with copious amounts of water.
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to testing.
- Do not eat, drink or smoke in the area where the specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ is available to receive and culture specimens.
- Do not interchange or mix reagents from different lots.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded in accordance with local regulations.

STORAGE AND STABILITY

- The kit should be stored at 2-30°C until the expiry date printed on the sealed pouch.
- The test must remain in the sealed pouch until use.
- Do not freeze.**
- Care should be taken to protect the components of the kit from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

SPECIMEN COLLECTION AND PREPARATION

Specimen Collection

Nasal swab sample:

- For proper test performance, use the swabs supplied in the kit.
- It is important to obtain as much secretion as possible. Therefore, to collect a nasal swab sample, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril), rotate the swab a few times against the nasal wall.

Nasopharyngeal swab sample:

- It is important to obtain as much secretion as possible. Therefore, to collect a nasopharyngeal swab sample, carefully insert the sterile swab into the nostril that presents the most secretion under visual inspection. Keep the swab near the septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times.


Specimen Transport and Storage:

- Specimens should be tested as soon as possible after collection. If transport of the samples is required, the swab should be stored in the transport tube without media. Alternatively, samples may be stored refrigerated (2-8°C), or at room temperature (15-30°C), in a clean, dry, closed container for up to eight hours before testing. Nasal wash/aspirate specimens may also be stored frozen (-70°C or colder) for up to one month.

PROCEDURE

Bring tests, specimens, and controls to room temperature (15-30°C) before use.


- Remove the test from its sealed pouch, and place it on a clean, level surface. Label the device with patient or control identification. For best results, the assay should be performed within one hour.
- Gently mix Extraction reagent solution. Remove the foil, and place the patient swab specimen into the Extraction Tube with buffer.

- Roll the swab at least 10 times while pressing the swab against the bottom and side of the Extraction Tube. Roll the swab head against the inside of the Extraction Tube as you remove it. Try to release as much liquid as possible. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Put on the tube tip, then add 3 drops of extracted sample into the sample well. Do not handle or move the Test Device until the test is complete and ready for reading.
- As the test begins to work, color will migrate across the membrane. Wait for the colored band(s) to appear. **The result should be read at 10-15 minutes.** Do not interpret the result after 15 minutes.


INTERPRETATION OF RESULTS

POSITIVE RESULT:


***Two lines appear.** One colored line should be in the control region (C), and another apparent colored line adjacent should be in the test region.

NEGATIVE RESULT:

One colored line appears in the control region (C). No line appears in the test region.

INVALID RESULT:

Control band fails to appear. Results from any test which has not produced a control band at the specified reading time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

***NOTE:**

- The intensity of color in the test regions may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test regions should be considered positive. Please note that this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

QUALITY CONTROL

- Internal procedural controls are included in the test. A colored band appearing in the control region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The Influenza A+B/COVID-19/RSV/Adeno/MP Antigen Rapid Test Device is for professional in vitro diagnostic use and should only be used for the qualitative detection of Influenza A and/or B, and/or COVID-19, and/or RSV, and/or Adeno, and/or Mycoplasma pneumoniae (MP).
- The etiology of respiratory infection caused by microorganisms other than Influenza A or B virus, COVID-19, RSV, Adenovirus, or MP will not be established with this test. The Influenza A+B/COVID-19/RSV/Adeno/MP Antigen Rapid Test Device is capable of detecting both viable and non-viable influenza, COVID-19, RSV, Adeno, or MP particles.
- The performance of the Influenza A+B/COVID-19/RSV/Adeno/MP Antigen Rapid Test Device depends on antigen load and may not correlate with cell culture (for influenza) or PCR (for COVID-19, RSV, Adeno, or MP) performed on the same specimen. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time rule out the presence of Influenza A and/or B, and/or COVID-19, and/or RSV, and/or Adenovirus, and/or MP viral antigens in the specimen, as they may be present below the minimum detection level of the test. As with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
- The validity of the Influenza A+B/COVID-19/RSV/Adeno/MP Antigen Rapid Test Device has not been proven for identification or confirmation of cell culture isolates. Inadequate or inappropriate specimen collection, storage, and transport may yield false negative test results.
- Although this test has been shown to detect cultured avian influenza viruses, COVID-19 viruses, RSV, and MP, including avian influenza A subtype H5N1 virus, the performance characteristics of this test with specimens from humans infected with H5N1 or other avian influenza viruses are unknown. Performance characteristics for Influenza A were established when Influenza A/H3 and A/H1 were the predominant Influenza A viruses in circulation.

PERFORMANCE CHARACTERISTICS

Table: COVID-19/Flu A&B/RSV/ADENO/M.P Ag combo rapid Test vs. PCR

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|------------|------------|------------|---------------|
| | | Positive | Negative | |
| Positive | 122 | 0 | 122 | |
| Negative | 2 | 234 | 236 | |
| Total Results | 124 | 234 | 358 | |

Sensitivity = 98.39% (95% CI = 94.31% to 99.49%)*

Specificity = 100.00% (95% CI = 98.38% to 100.00%)*

Accuracy = 99.44% (95% CI = 97.99% to 99.85%)* * 95% Confidence Interval

Flu B

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|------------|------------|------------|---------------|
| | | Positive | Negative | |
| Positive | 105 | 1 | 106 | |
| Negative | 2 | 270 | 272 | |
| Total Results | 107 | 271 | 378 | |

Sensitivity = 98.13% (95% CI = 93.44% to 99.49%)*

Specificity = 99.63% (95% CI = 97.94% to 99.93%)*

Accuracy = 99.21% (95% CI = 97.69% to 99.73%)* * 95% Confidence Interval

COVID-19

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|------------|------------|------------|---------------|
| | | Positive | Negative | |
| Positive | 158 | 0 | 158 | |
| Negative | 3 | 230 | 233 | |
| Total Results | 161 | 230 | 391 | |

Sensitivity = 98.14% (95% CI = 94.67% to 99.36%)*

Specificity = 100.00% (95% CI = 98.36% to 100.00%)*

Accuracy = 99.23% (95% CI = 97.77% to 99.74%)* * 95% Confidence Interval

RSV

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|-----------|-----------|------------|---------------|
| | | Positive | Negative | |
| Positive | 66 | 0 | 66 | |
| Negative | 2 | 96 | 98 | |
| Total Results | 68 | 96 | 164 | |

Sensitivity = 97.06% (95% CI = 89.90% to 99.19%)*

Specificity = 100.00% (95% CI = 96.15% to 100.00%)*

Accuracy = 98.78% (95% CI = 95.66% to 99.66%)* * 95% Confidence Interval

Adeno

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|-----------|------------|------------|---------------|
| | | Positive | Negative | |
| Positive | 76 | 2 | 78 | |
| Negative | 1 | 112 | 113 | |
| Total Results | 77 | 114 | 191 | |

Relative Sensitivity: 98.70% (95% CI = 93.00% to 99.77%)*

Relative Specificity: 99.12% (95% CI = 95.16% to 99.84%)*

Overall Accuracy: 98.95% (95% CI = 96.24% to 99.67%)* * 95% Confidence Interval

MP

| COVID-19/Flu A&B/RSV/Adeno/MP Ag Test | Results | PCR | | Total Results |
|---------------------------------------|------------|------------|------------|---------------|
| | | Positive | Negative | |
| Positive | 142 | 4 | 146 | |
| Negative | 8 | 296 | 304 | |
| Total Results | 150 | 300 | 450 | |

Relative Sensitivity: 94.67% (95% CI = 89.83% to 97.27%)*

Relative Specificity: 98.66% (95% CI = 96.60% to 99.48%)*

Overall Accuracy: 97.32% (95% CI = 95.38% to 98.46%)* * 95% Confidence Interval

ANALYTICAL SPECIFICITY AND CROSS-REACTIVITY

The Influenza A+B/COVID-19/RSV/Adeno antigen Rapid Test Device was evaluated with a total of 30 bacterial and viral isolates. Bacterial isolates were evaluated at a concentration between 107 and 109 org/mL. Viral isolates were evaluated at a concentration of at least 104–108 TCID50/mL. Adenovirus 18 and Parainfluenza virus 3 were tested at 102 TCID50/mL. None of the organisms or viruses listed below gave a positive result in the Influenza A+B/COVID-19/RSV antigen Rapid Test Device.

Bacterial Panel:

- Acinetobacter calcoaceticus
- Neisseria gonorrhoeae
- Pseudomonas aeruginosa
- Streptococcus pneumoniae
- Proteus vulgaris
- Streptococcus sp. Gp. C
- Mycobacterium tuberculosis

- Bacteroides fragilis
- Neisseria meningitidis
- Staphylococcus aureus
- Streptococcus sanguis
- Streptococcus sp. Gp. B
- Streptococcus sp. Gp. G
- Mycoplasma orale

Viral Panel:

- Human Adenovirus B
- Human Adenovirus C
- Human Coronavirus OC43
- Human Coxsackievirus A9
- Coxsackievirus B5
- Human herpesvirus2
- Human Rhinovirus 2
- Human Rhinovirus 14
- Mumps
- Sendai virus
- Parainfluenza virus 2
- Parainfluenza virus 3

INTERFERING SUBSTANCES

Whole blood, and several over-the-counter (OTC) products and common chemicals were evaluated and did not interfere with the Influenza A+B/COVID-19/RSV/Adeno Test at the levels tested: whole blood (2%); three OTC mouthwashes (25%); three OTC throat drops (25%); three OTC nasal sprays (10%); 4-Acetaminophenol (10 mg/mL); Acetylsalicylic Acid (20 mg/mL); Chlorpheniramine (5 mg/mL); Dextromethorphan (10 mg/mL); Diphenhydramine (5 mg/mL); Ephedrine (20 mg/mL); Guaiacol glyceryl ether (20 mg/mL); Oxymetazoline (10 mg/mL); Phenylephrine (100 mg/mL); and Phenylpropanolamine (20 mg/mL).

LITERATURE REFERENCES

- Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. J Autoimmun. 2020; Feb26:102433. doi: 10.1016/j.jaut.2020.102433.
- Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID19): The epidemic and the challenges. Int J Antimicrob Agents. 2020; Mar 55(3): 105924. doi: 10.1016/j.ijantimicag.2020.105924.
- Shaw MW, Arden NH and Massab HF. New aspects of influenza viruses. Clin. Microbiol. Rev. 5: 74-92 (1992).
- WHO recommendations on the use of rapid testing for influenza diagnosis, July 2005.
- Peaper DR, Landry ML. Rapid diagnosis of influenza: state of the art. Clin Lab Med. 2014; 34(2):365-385. doi:10.1016/j.cl.2014.02.009.
- National Committee for Clinical Laboratory Standards. Protection of laboratory workers from infectious disease transmitted by blood, body fluids, and tissue: Tentative guideline. NCCLS Document M29-T. Villanova, PA.: NCCLS, 1989.
- Glezen, WP; Taber, LH; Frank, AL; Kasel, JA (1986). "Risk of primary infection and reinfection with respiratory syncytial virus". American journal of diseases of children (1960). 140(6): 543-6. doi: 10.1001/archpedi.1986.02140200053026. PMID 3706232.
- Schmitz HR, Wigand Heinrich W. Worldwide epidemiology in human adenovirus infection. J. Epidemiol. 1983; 117:455-466.
- Parrott GL, et al. (2016). "A Compendium for Mycoplasma pneumoniae." Frontiers in Microbiology 7, 513.
- Waites KB, Talkington DF. (2004). "Mycoplasma pneumoniae and its role as a human pathogen." Clinical Microbiology Reviews 17(4), 697-728.
- Atkinson TP, et al. (2008). "Epidemiology, clinical manifestations, pathogenesis and laboratory detection of Mycoplasma pneumoniae infections." FEMS Microbiology Reviews 32(6), 956-73.

Index of Symbols

| | | | | | |
|--|----------------------------------|--|---------------|--|---------------------------|
| | Consult Instruction for use | | Tests per kit | | Authorized Representative |
| | For in vitro diagnostic use only | | Use by | | Do not reuse |
| | Store between 2-30°C | | Lot Number | | Catalog # |
| | Do not use if package is damaged | | | | |

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