**Shandon Instant Hematoxylin
Shandon Instant Eosin, Alcoholic
Shandon Instant Eosin, Aqueous**

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

Shandon Instant Hematoxylin and Eosins are intended for use in histological and cytological staining procedures. All are formulated and packaged to provide our customers with the best value. Every reagent and raw material is tested to meet the highest standards. Only those dyes certified by the Biological Stain Commission are used in manufacturing. Each lot of stain is subjected to a series of quantitative and qualitative tests. These tests ensure maximum performance, stability and lot to lot consistency. From scientific testing and spectrophotometry, to staining a simple control sample, we make every effort to guarantee your satisfaction with Thermo Electron stains and reagents.

Shandon Bluing Reagent is a gentle, pH-controlled solution for bluing hematoxylin. It replaces the strong alkaline rinses or long tap water washes which may cause loss of cells or tissue sections.

Nu-Clear™ improves cell clarity by removing excess hematoxylin from the tissue in regressive staining protocols. It may also be used to remove background hematoxylin from the slide. Nu-Clear is designed for hand staining as well as use on automatic stainers. It is available in 2 concentrations: Nu-Clear I is faster acting, with a 10-second differentiating step; Nu-Clear II will accomplish similar results in 20 seconds. The less concentrated Nu-Clear II is recommended for use on Linear stainers and with procedures that use Gill Hematoxylin. The usage and timing depend on the staining intensity preferred by the laboratory.

Cytoplasmic Counterstains are strong stains that produce bright, well-differentiated staining to compliment the nucleus. These stains are recommended in routine H & E or Papanicolaou techniques.

Shandon Xylene Substitute is an odorless Xylene Substitute. A safe and convenient alternative to xylene it is an aliphatic hydrocarbon that is not greasy. It may be used on stainers, automatic coverslippers and tissue processors.

General Comments:

1. Staining times and number of rinses can be adjusted to suit individual preferences in stain intensity.
2. If using a mercuric chloride fixative (i.e. B-5, Zenkers) for histology specimens remove the mercuric chloride crystals with iodine. After removing the wax and before the hematoxylin step, rinse slides in water, clear with a sodium thiosulfate solution then rinse again in water.
3. The coating agent found in Cell-Fixx™ and Cytospin® Collection Fluid must be removed before staining. Generally, soaking in 95% alcohol will remove the coating. Failure to completely remove the coating will result in poor and uneven staining. Users of commercial spray fixatives should follow the manufacturer's recommendation for removal.
4. Appropriate precautions should be taken to minimize cellular contamination from specimens.
5. The pH of tap water varies considerably. Use of tap water in the staining procedures may give unpredictable results. Long-running tap water rinses are not recommended. Chlorine in tap water will bleach out the stain.
6. Control slides should be run daily, or whenever reagents are changed, to verify the quality of the staining process.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet.

Shandon Instant Hematoxylin	Shandon Instant Eosin, Alcoholic	Shandon Instant Eosin, Aqueous
Shandon Instant Hematoxylin is an innovative concept in nuclear staining. It can be used in cytologic and histologic staining procedures, including immunoperoxidase techniques. Just simply substitute Instant Hematoxylin for Gill, Harris, or Mayer in your procedure. Shandon Instant Hematoxylin contains an oxidizing agent which is not mercury-based, so there is minimal hazard to our environment.	Shandon Instant Eosin, Alcoholic, is an innovative concept in cytoplasmic staining for histology. It can be directly substituted for any alcoholic Eosin-Y counterstain.	Shandon Instant Eosin, Aqueous, is an innovative concept in cytoplasmic staining for histology. It can be directly substituted for any aqueous Eosin-Y counterstain.
	Shelf Life	
Shandon Instant Hematoxylin has an indefinite shelf life in the dry state. Once it has been reconstituted, it is good for 6 months if properly stored.	Shandon Instant Eosin, Alcoholic, has an indefinite shelf life in the dry state. Once the stain has been reconstituted, the stain is good for 2 months if properly stored.	Shandon Instant Eosin, Aqueous, has an indefinite shelf life in the dry state. Once the stain has been reconstituted, the stain is good for 2 months if properly stored.
	Storage	
Shandon Instant Hematoxylin (both dry and reconstituted) should be stored at room temperature in tightly closed containers, and kept away from direct sunlight.	Shandon Instant Eosin, Alcoholic (both dry and reconstituted) should be stored at room temperature in tightly closed containers, and kept away from direct sunlight.	Shandon Instant Eosin, Aqueous (both dry and reconstituted) should be stored at room temperature in tightly closed containers, and kept away from direct sunlight.

Shandon Instant Hematoxylin	Shandon Instant Eosin, Alcoholic	Shandon Instant Eosin, Aqueous
	Mixing Instructions	
1. Pour entire contents of one "Part A" and one "Part B" into a container that will hold over 1 liter of solution. 2. Add 1.0 liter of deionized or distilled water. DO NOT USE TAP WATER. 3. Stir or shake the solution for several minutes, or until powders are dissolved. Wait at least 1 hour before using. Best results are obtained after 12 hours. 4. Filter before using. A precipitate and surface scum are normal for this product. 5. A more delicate stain may be desired for cytology procedures. In that case, the amount of water used in mixing can be increased from 1.0 liter to 1.5 or 2.0 liters.	1. Fill flask or other suitable container with 1.0 liter of 70% ethanol or denatured ethanol. To make a 70% solution, combine 700 ml absolute alcohol with 300 ml deionized or distilled water. Freshly prepared 70% alcohol will be warm and will dissolve Shandon Instant Eosin more rapidly than a solution at room temperature. CAUTION: Do not use water alone. 2. Carefully pour contents of one Instant Eosin bottle into the alcohol solution. 3. Mix for 15 – 60 minutes until all components have dissolved. Time required is dependent upon temperature. WARNING: BECAUSE OF THE FLAMMABILITY OF THE SOLVENT, DO NOT HEAT THE SOLUTION TO HASTEN RECONSTITUTION. 4. Use immediately or store in a tightly capped container.	1. Fill a suitable container with 1.0 liter deionized or distilled water. 2. Carefully pour contents of Instant Eosin, Aqueous bottle into the water. 3. Mix for 1 – 2 minutes until contents have dissolved. 4. Use immediately or store in a tightly capped container.

Shandon Instant Hematoxylin	Shandon Instant Eosin, Alcoholic	Shandon Instant Eosin, Aqueous
	Staining Procedures	
- Guidelines for cytology and histology staining procedures can be found at the rear of the booklet. - This stain may be used progressively or regressively, with times varying from 30 seconds (light) to 2 minutes (dark) for progressive staining. Timing will be longer for regressive procedures. - Staining times and number of rinses can be adjusted to suit individual preferences in stain intensity. - The alcohol used in the staining procedure may be either reagent alcohol or ethanol. - When using an automated stainer, follow the procedure recommended in the operator's manual. NOTE: In the cytology staining procedure, the time for Instant Hematoxylin will depend on the dilution used when mixing the powders with water. If 1.5 or 2.0 liters of water was used, the time will be 30 seconds to 4 minutes. If 1.0 liter was used, 30 seconds to 2 minutes is suggested.	- Guidelines for histology staining procedures can be found at the rear of the booklet. - Recommended staining times range from 30 seconds to 4 minutes. If stain is too dark after 30 seconds, dilute the stain 1:1 with 70% ethanol or denatured ethanol (or make the dye up in 2 liters of solvent). - Staining times and number of rinses can be adjusted to suit individual preferences in stain intensity. You may also dilute the reconstituted product with 70% ethanol. - The alcohol used in the staining procedure may be either reagent alcohol or ethanol. - When using an automatic stainer, follow the procedure recommended in the operator's manual. - The formation of a pink precipitate is normal for this product. It may be filtered out with no adverse effect on staining.	- Guidelines for histology staining procedures can be found at the rear of the booklet. - Staining times and number of rinses can be adjusted to suit individual preferences in counterstain intensity. - Recommended staining time is 3 minutes. Wash briefly in water before dehydration. - The alcohol used in the staining procedure may be either reagent alcohol or ethanol. - When using an automatic stainer, follow the procedure recommended in the operator's manual.

**** Thermo has not validated the protocols illustrated on pages 22 and 23 and takes no responsibility for their use. Customers should always validate protocols before placing any reliance on them.**

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HISTOLOGY STAINING GUIDELINES

Instant Hematoxylin and Instant Eosin Alcoholic		
Step	Solution	Time
1	Xylene or Xylene Substitute	3 - 5 minutes
2	Xylene or Xylene Substitute	3 - 5 minutes
3	Xylene or Xylene Substitute	3 - 5 minutes
4	Alcohol, 100%	1 minute
5	Alcohol, 100%	1 minute
6	Alcohol, 95%	1 minute
7	Water Wash	1- 2 minutes
8	Instant Hematoxylin (filter before use) (depending on strength)	2 – 5 minutes
9	Water Wash	1-2 minutes
10	Nu-Clear I	10 seconds
	OR Nu-Clear II	20 seconds
11	Water Wash	1 minute
12	Bluing Reagent	1 minute
13	Water Wash	1-2 minutes
14	Alcohol (70%-95%)	1 minute
15	Instant Eosin, alcoholic	30 seconds to 1 minute
16	Alcohol, 95%	20 seconds to 1 minute
17	Alcohol, 95%	20 seconds to 1 minute
18	Absolute Alcohol	1 minute
19	Absolute Alcohol	1 minute
20	Absolute Alcohol	1 minute
21	Xylene or Shandon Xylene Substitute	1-2 minutes
22	Xylene or Shandon Xylene Substitute	1-2 minutes
23	Xylene or Shandon Xylene Substitute	1-2 minutes
	Mount with Shandon-Mount or Histo-Mount*	

Instant Hematoxylin and Instant Eosin Aqueous		
Step	Solution	Time
1	Xylene or Xylene Substitute	3 - 5 minutes
2	Xylene or Xylene Substitute	3 - 5 minutes
3	Xylene or Xylene Substitute	3 - 5 minutes
4	Alcohol, 100%	1 minute
5	Alcohol, 100%	1 minute
6	Alcohol, 95%	1 minute
7	Water Wash	1-2 minutes
8	Instant Hematoxylin (filter before use) (depending on strength)	2 – 5 minutes
9	Water Wash	1-2 minutes
10	Nu-Clear I	10 seconds
	OR Nu-Clear II	20 seconds
11	Water Wash	1 minute
12	Bluing Reagent	1 minute
13	Water Wash	1-2 minutes
14	Instant Eosin, aqueous	3 minutes
15	Running Water	30 seconds
16	Alcohol, 50% or 95 %	1 minute
16	Alcohol, 95%	1 minute
17	Absolute Alcohol	1 minute
18	Absolute Alcohol	1 minute
20	Absolute Alcohol	1 minute
21	Xylene or Shandon Xylene Substitute	1-2 minutes
22	Xylene or Shandon Xylene Substitute	1-2 minutes
23	Xylene or Shandon Xylene Substitute	1-2 minutes
	Mount with Shandon-Mount or Histo-Mount*	

* Shandon Mount should be used with Xylene and Histo-Mount should be used only with Xylene substitute.

CYTOLOGY STAINING GUIDELINES

Instant Hematoxylin and Instant Eosin Alcoholic		
Step	Solution	Time
	95% Alcohol (if needed)	10 – 15 minute pre-soak
1	Distilled Water Wash	1 minute
2	Instant Hematoxylin (filter before use) (depending on strength)	3 minutes
3	Distilled Water Wash	30 – 60 seconds
4	Nu-Clear I	10 seconds
	OR Nu-Clear II	20 seconds
5	Distilled Water Wash	30 – 60 seconds
6	Bluing Reagent	1 minute
7	Distilled Water Wash	30 – 60 seconds
8	95% Alcohol	10 seconds
9	Orange G-6	1 – 3 minutes
10	95% Alcohol	10 seconds
11	95% Alcohol	10 seconds
12	Instant Eosin, Alcoholic	2 to 4 minutes
13	95% Alcohol	10 seconds
14	95% Alcohol	10 seconds
15	Absolute Alcohol	20 seconds
16	Absolute Alcohol	20 seconds
17	Xylene or Shandon Xylene Substitute	1 minute
18	Xylene or Shandon Xylene Substitute	1 minute
19	Xylene or Shandon Xylene Substitute	5 minutes
	Mount with Histo-Mount or Shandon-Mount*	

STAINING GUIDELINES

Step	Solution	Time
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		
23		

* Shandon Mount should be used with Xylene and Histo-Mount should be used only with Xylene substitute.

ORDER INFORMATION

Product	Unit	REF
Shandon Instant Hematoxylin (one package makes approximately 6 liters of Hematoxylin stain.)	Package of 6 sets	6765015
Shandon Instant Eosin, Alcoholic (one package makes approximately 6 liters of counterstain)	Package of 6 sets	6765040
Shandon Instant Eosin, Aqueous (one package makes approximately 6 liters of counter stain)	Package of 6 sets	6765540
Shandon Bluing Reagent	2 x 1 liter	6769001
Shandon Bluing Reagent	1 x 4 liter	6769002
Nu-Clear I	3.8 Liter	6769008
Nu-Clear II	3.8 Liter	6769009
Shandon Eosin, Alcoholic	2 x 1 liter	6766007
Shandon Eosin, Alcoholic	1 x 4 liter	6766008
Shandon Eosin, Aqueous	2 x 1 liter	6766009
Shandon Eosin, Aqueous	1 x 4 liter	6766010
Shandon Orange G-6	2 x 1 liter	6766005
Shandon Orange G-6	1 x 4 liter	6766006
Shandon EA-50	2 x 1 liter	6766001
Shandon EA-50	1 x 4 liter	6766002
Shandon EA-65	2 x 1 liter	6766003
Shandon EA-65	1 x 4 liter	6766004
Shandon Xylene Substitute	1 x 4 liter	9990505
Shandon Xylene Substitute	4 x 4 liter	6764506
Shandon Xylene Substitute	2 x 10 liter	6764508
Shandon Xylene Substitute	5 Gal. (19L) x 1	9990507
Shandon Xylene Substitute	55 Gal. (209L) x 1	9990509

Refer to Manual

Store Between
15 - 25°

Keep Dry

For in Vitro
Diagnostic use

Irritant

Date of Mfg.

Lot Number

Product number

Thermo Electron Corporation Anatomical Pathology has realized since its beginnings the value of in-house assistance. If you have questions about ordering, shipping, delivery, inventory status or pricing, just phone our Customer Service Department toll-free at 1-800-547-7429 in the U.S., or our direct subsidiaries in the UK, France, Germany and Australia. Our technical support staff consists of certified Cytotechnologists and Histologists.