

ORDERING INFORMATION

THE MYNXGRIP™ VASCULAR CLOSURE DEVICE INCLUDES:

(1) Balloon catheter with integrated sealant

(1) 10ml locking syringe

| SIZE | COLOR | MYNX ORDER NUMBER |
|-------|-------|-------------------|
| 6F/7F | Green | MX6721 |

IN THE UNITED STATES, FAX YOUR MYNXGRIP ORDER TO ACCESSCLOSURE AT (877) 933-0133.

OUTSIDE OF THE UNITED STATES, CONTACT YOUR LOCAL MYNX DISTRIBUTOR.

INDICATIONS FOR USE

The MynxGrip Vascular Closure Device is indicated for use to seal femoral arterial access sites while reducing times to hemostasis and ambulation in patients who have undergone diagnostic or interventional endovascular procedures utilizing a 6F or 7F procedural sheath.

PRECAUTIONS

The MynxGrip device should only be used by a trained licensed physician or healthcare professional.

The MynxGrip device should not be used in patients with a known allergy to PEG.

WARNINGS

Do not use if the sterilization indicator dot on the pouch is yellow/gold.

Do not use if components or packaging appear to be damaged or defective or if any portion of the packaging has been previously opened.

Do not reuse or resterilize. The MynxGrip device is for single use only. The catheter is loaded with a single hydrogel sealant. Reuse of the device would result in no delivery of hydrogel sealant.

Do not use the MynxGrip device if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and/or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma/bleed. Perform a femoral angiogram to verify the location of the puncture site.

Do not use the MynxGrip device if the puncture is through the posterior wall or there are multiple punctures, as such punctures may result in a retroperitoneal hematoma/bleed.

access | closure™

AccessClosure, Inc.
645 Clyde Avenue
Mountain View, CA
94043 USA

PHONE 877-700-MYNX (6969) U.S. only
+1/650-864-5473

FAX 877-933-0133 U.S. only
+1/650-625-8584

MGRIP™

GRIP™
TECHNOLOGY


MGRIP™
MYNXGRIP™ VASCULAR CLOSURE DEVICE

FINISH STRONG

REDEFINING ACTIVE CLOSURE


Built on the proven Mynx™ platform, the MynxGrip™ Vascular Closure Device offers an active, extravascular and patient-friendly vascular closure solution. MynxGrip adds the proprietary Grip Technology™ to the distal end of the original Mynx Sealant. The result is a sealant that actively grips and seals the arteriotomy while expanding and filling the tissue tract, providing effective hemostasis. The MynxGrip Sealant dissolves within thirty days, leaving nothing behind but a healed artery.

MYNXGRIP VASCULAR CLOSURE DEVICE




MYNXGRIP DEPLOYMENT STEPS

DEPLOY BALLOON




Achieve temporary hemostasis at the arteriotomy

PLACE THE SEALANT



The Grip Tip actively seals the arteriotomy while the Mynx Sealant fills the tissue tract

REMOVE THE DEVICE



Platelets and blood cells collect inside the sealant's porous matrix

FINAL RESULT




The sealant dissolves within 30 days leaving nothing behind but a healed artery

SEALANT SCIENCE FEATURING GRIP TECHNOLOGY




THE MYNXGRIP SEALANT

The MynxGrip Sealant consists of Polyethylene Glycol (PEG), a water-soluble, bio-inert, non-thrombogenic polymer. The MynxGrip Sealant is comprised of two configurations of PEG, the Grip Tip and the Mynx Sealant.



GRIP TECHNOLOGY

Once the MynxGrip Sealant enters the tissue tract, the body's temperature and pH level cause the Grip Tip to soften and interlock with the contours of the vessel wall, effectively gripping the artery and providing active closure.



THE MYNX SEALANT

The Mynx Sealant's porous structure absorbs blood and subcutaneous fluids. The Sealant swells three to four times its original size, filling the tissue tract.



THE MYNXGRIP SEALANT

The MynxGrip Sealant actively grips the artery while expanding and filling the tissue tract, providing a durable hemostasis and a platform for natural healing.

Bovine vessel suspended by MynxGrip Sealant

The only **active extravascular** closure device that leaves nothing behind