NAUDOJIMO INSTRUKCIJA

ORTHO CONFIDENCE™ WB

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The following symbols may have been used in the labeling of this product.

 Use by or Expiration Date (Year-Month-Day)	 Manufacturer	 Temperature Limitation
 Lot Number	 In vitro Diagnostic Medical Device	 Consult Instructions for Use
 Catalog Number or Product Code	 Authorised Representative in the European Community	

TRANSFUSION MEDICINE


NAUDOJIMO INSTRUKCIJA

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Peržiūros data	Versija	Sekcija	Techninių pakeitimų aprašymas *
2020-06-03	8.0	Additional Materials Required but not provided for Manual Testing	Pridėtas ORTHO Optix™ Reader
2019-09-19	7.0	Glossary of Symbols	Added EC Representative symbol
		Back Page	Added EC Representative address
2019-05-07	6.0	Back Page	Changed Notified Body number from 0843 (UL) to 1434 (PCBC)
2018-08-02	5.0	Back Page	Updated Manufacturer's address Changed from © Ortho-Clinical Diagnostics, Inc. to © Ortho Clinical Diagnostics.
2017-02-17	4.0	Precautions	First bullet: 'Do not transfer any of these reagent red cells' changed to 'Do not transfer any of the vial contents'.
		Procedure	Addition of ORTHO VISION® Analyzer. Addition of ORTHO VISION® Max Analyzer.
		Additional Materials Required but not Provided for Automated Testing	Addition of ORTHO VISION® Analyzer. Addition of ORTHO VISION® Max Analyzer.
		Additional Materials Required but not provided for Manual Testing	Addition of ORTHO™ Workstation.
		Performance Characteristics	Changed Customer Technical Support Center to Ortho Care™ Technical Solutions Center.
		Back Page	Update to Corporate Logo to remove 'Part of the Johnson and Johnson Family of Companies' . CE Mark updated to include Notified Body number. Updated copyright information.
2013-08-28	3.0	Precautions	Removed information from third bullet: a minimum of Removed ninth bullet: Refer to the Material Safety Data Sheet (MSDS) available from Ortho-Clinical Diagnostics, Inc. for more specific information.
2012-04-03	2.0	Reagents	The Fy ^a antibodies were reclassified from polyclonal origin to monoclonal origin.
		Precautions	Added 3 rd bullet: Once product vials are opened, the reagent's performance characteristics are maintained for a minimum of 14 days.
		Interpretation of Results	Included rows for A ₂ cells and O cells.
		Limitations of the Procedure	Added information to step 5 pertaining to Anti-c specificity and the possible reaction with Affirmagen O cells.
		Back Page	Updated the corporate logo.
2010-06-21	1.0	Throughout	Updated trademarks.
			Initial version of Instructions for Use

* The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.

Quotient Suisse S.A.
Unit B1 Terre Bonne Business Park
Route de Crassier 13
Eysins 1262, Switzerland

Alba Bioscience Limited
James Hamilton Way, Penicuik, EH26 0BF, UK

Distributed by:

Ortho-Clinical Diagnostics, Inc.
1001 US Highway 202
Raritan, NJ 08869 USA

ORTHO CONFIDENCE™ WB Quality Control Kit is a
trademark of Ortho Clinical Diagnostics,
© Ortho Clinical Diagnostics, 2010-2020

Ortho Clinical Diagnostics