

# ISE Clean. Soln./Elecsys SysClean cobas®

ISE Cleaning Solution / Elecsys SysClean

Order information

REF		CONTENT	Analyzer(s) on which kit(s) can be used
11298500316	11298500500	ISE Cleaning Solution / Elecsys SysClean (5 x 100 mL)	cobas c cobas e

## English

### Intended use

For the cleaning of ISE units on Roche/Hitachi analyzers.  
For the cleaning of Elecsys and **cobas e** immunoassay analyzers.

### Summary

ISE Cleaning Solution / Elecsys SysClean is an alkaline cleaning solution with antimicrobial properties.

### Reagents - working solutions

5 bottles, each containing 100 mL

#### Components:

Sodium hydroxide 3 mol/L

Sodium hypochlorite solution (< 2 % active chlorine)

Additive

### Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

- H290 May be corrosive to metals.
- H314 Causes severe skin burns and eye damage.
- H410 Very toxic to aquatic life with long lasting effects.

### Prevention:

- P273 Avoid release to the environment.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

### Response:

- P303 + P361 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.  
+ P353
- P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
+ P310 Immediately call a POISON CENTER/ doctor.
- P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
+ P338 Continue rinsing. Immediately call a POISON CENTER/ doctor.  
+ P310
- P391 Collect spillage.

### Hazardous components:

- sodium hydroxide
- sodium hypochlorite

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

### Reagent handling

Ready for use

### Storage and stability

Store protected from light.

The solution is stable up to the stated expiration date when stored at 2-8 °C.

### Materials provided

See "Reagents – working solutions" section for reagents.

### Materials required (but not provided)

- General laboratory equipment
- cobas c** or **cobas e** analyzers

### Assay

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and method sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, ELECSYS and NAVIFY are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH  
Sandhofer Strasse 116  
68305 Mannheim, Germany  
www.roche.com

+800 5505 6606

