# Complete Early Cleavage Medium® (ECM®) with DSS

Catalog No. 90142

2 x 20 mL

For assisted reproductive procedures.  
 Für assistierte Reproduktionsverfahren.  
 Per tecniche di riproduzione assistita.  
 Para utilización en técnicas de reproducción asistida.  
 Pour les techniques de procréation médicalement assistée.  
 Para técnicas de reprodução assistida.  
 Για διαδικασίες υποβοηθούμενης αναπαραγωγής.  
 Pro postupy asistované reprodukce.  
 Til brug ved assisteret reproduktion.  
 Avusteisiin lisäaöntymismenetelmiin.  
 Paredzēts ar palīgīdzekļiem veicamām reproduktīvām procedūrām.  
 Voor geassisteerde voortplantingsprocedures.  
 Do procedur wspomagane go rozrodu.  
 Pentru proceduri de reproducere asistată.  
 För procedurer för assisterad befruktning.  
 Kasutamiseks reproduktiivse abi protseduurides.  
 Asszisztált reprodukciós eljárásokhoz.  
 Skirta pagalbinio apvaisinimo procedūroms.  
 Yardımla üreme prosedürleri için.  
 Na postupy asistovanej reprodukcie.

**Symbols:**

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Storage Temperature
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged
	Manufacturer: Irvine Scientific®
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

**ENGLISH**

**EU CAUTION:** For Professional Use Only  
**US CAUTION:** Federal (U.S.) law restricts this device to sale by or on the order of a physician.

**INTENDED USE**  
 Complete Early Cleavage Medium® (ECM®) with DSS\* is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of Complete ECM with DSS as a culture medium through day 3 of development.

**PRODUCT DESCRIPTION**  
 Complete ECM with DSS consists of Complete ECM supplemented with 10% Dextran Serum Supplement (DSS), for a total macromolecule concentration of 7 mg/mL. The macromolecule components consist of 5 mg/mL therapeutic grade Human Serum Albumin (HSA) and 2 mg/mL dextran. This medium contains the antibiotic gentamicin sulfate (9 µg/mL).

**COMPOSITION:**

EDTA	Sodium Lactate
Sodium Citrate	Dextran 40
Sodium Chloride	Sodium Bicarbonate
Potassium Chloride	Magnesium Sulfate
Dextrose	Calcium Chloride
L-Alanyl-L-Glutamine	Phenol Red
Taurine	Gentamicin Sulfate
Sodium Pyruvate	Human Serum Albumin

**QUALITY ASSURANCE**  
 Complete ECM with DSS is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of Complete ECM with DSS is tested for:  
 Endotoxin by Limulus Amebocyte Lysate (LAL) methodology  
 Biocompatibility by Mouse Embryo Assay (one-cell)  
 Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

**BUFFER SYSTEM**  
 Complete ECM with DSS uses Sodium Bicarbonate as a buffering system. This is specifically designed for use in a CO<sub>2</sub> incubator.

**DIRECTIONS FOR USE**  
**EQUILIBRATION**  
 Complete ECM with DSS should be warmed to 37°C and equilibrated to the desired pH in a CO<sub>2</sub> incubator prior to use.

Fertilization may be performed by standard in vitro fertilization (IVF) in a culture dish containing pre-equilibrated Complete ECM with DSS, or by ICSI according to the laboratory's standard procedures.

Once fertilization occurs, the reproductive specialist should transfer the zygote/embryo into a new dish with pre-warmed equilibrated Complete ECM with DSS for the growth phase. Allow embryo to grow until desired developmental stage (up to 3 days).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

**STORAGE INSTRUCTIONS AND STABILITY**  
 Store the medium refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

For unopened bottles:  
 The shelf life of the product remains as the expiration date stated on bottle when stored under the recommended conditions of 2° to 8°C.

**PRECAUTIONS AND WARNINGS**  
 This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not reddish-orange in color.

**CAUTION:** Ammonium has been shown to induce retarded embryo development in vitro, specifically in the development of the inner cell mass.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.


**Not for injection use.**  
 Do not use any bottle in which the sterile packaging has been compromised.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used after the initial opening of the container.

**EU:** Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time Irvine Scientific Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

**US:** This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

**CONTRAINDICATION**  
 Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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