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Sterilized using aseptic processing techniques (filtration)

Expiration:  
Year - Month - Day

Caution, consult accompanying documents

Storage Temperature

Do not resterilize

Do not use if package is damaged

Manufacturer:  
Irvine Scientific®

Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner).

CE Mark

Emergo Europe - Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

## ENGLISH

**EU Caution:** For Professional Use Only

### INTENDED USE

Sperm Maintenance Medium with Glycerol is intended for use in assisted reproductive procedures involving the cryopreservation and storage of human sperm.

### PRODUCT DESCRIPTION

The formulation of Sperm Maintenance Medium with Glycerol is a modified version of Human Sperm Preservation Medium (HSPM) that was first described by Mahadevan and Trounson in 1983. The major formulation changes include lowering the osmolality, stabilizing the buffer system and increasing the albumin and glycerol concentrations.

### COMPOSITION

Salts and Ions  
Potassium Chloride  
Calcium Lactate  
Magnesium Chloride  
Calcium Chloride

Protein Source  
Human Serum Albumin

### Amino Acids

L-Glutamine  
Glycine  
Hypotaurine  
Adenosine

### Cryoprotectant

Glycerol  
Sucrose

Copolymer Surfactant  
Pluronic F-68

### Energy Substrates

Dextrose  
Inositol  
Pyruvic Acid, Sodium Salt

Vitamin  
Ascorbic Acid

Buffers  
HEPES  
Potassium Bicarbonate

### QUALITY ASSURANCE

Sperm Maintenance Medium is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of Sperm Maintenance Medium is tested for:  
Endotoxin by Limulus Amebocyte Lysate (LAL) methodology  
Sterility by the current USP Sterility Test <71>  
Sperm Cryo Survival Assay

All results are reported on a lot specific Certificate of Analysis which is available upon request.

### DIRECTIONS FOR USE

1. Semen is collected by masturbation following 2 to 3 days of abstinence.
2. The sample is allowed to liquefy at room temperature or 37°C for 30 minutes.
3. One vial of a previously prepared aliquot of Sperm Maintenance Medium with Glycerol is thawed and brought to room temperature or 37°C. If antibiotics are desired they may be added at this step.
4. The liquefied sample is transferred to a sterile 15 mL conical centrifuge tube. The specimen volume is determined and an appropriate volume of thawed Sperm Maintenance Medium is added drop-wise until a 3:1 sample to medium ratio is reached. For example, for each 1 mL of specimen add 0.33 mL of medium.
5. Aliquot the sample-medium mixture into labelled cryotubes or straws. To allow for expansion, do not overfill cryotubes.
6. Freeze specimens either directly or after optional slow cooling step (see below) using a programmable freezer or vapor freezing procedure.
7. Optional directions for slow cooling prior to freezing: Attach the filled cryovials to an aluminum cane. Submerge the cane(s), sample end down, into an ambient temperature water bath (i.e. a plastic beaker containing 500 mL of water) and then place the water bath in a refrigerator (2 to 8°C) for 60 to 90 minutes before freezing.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

### STORAGE INSTRUCTIONS AND STABILITY

Store Sperm Maintenance Medium with Glycerol at -10°C or colder.

Do not expose medium to repeated freeze-thaw cycles. If smaller aliquots are desired, thaw the product, aliquot into sterile labeled containers and freeze until time of use.

When stored as directed, Sperm Maintenance Medium is stable until the expiration date shown on the bottle label.

### Duration Following Bottle Opening:

Product should be used within (8) weeks from opening when stored under the recommended conditions of 2° to 8°C.

### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures that include the indicated application for which the device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium that shows evidence of particulate matter, cloudiness or is not clear and colorless.

Do not use any bottle in which the sterile packaging has been compromised.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.

**The medium does not contain antibiotics. For procedures requiring antibiotics, these may be optionally added to the medium prior to use. In all cases, antibiotic usage should be determined by appropriate medical personnel to ensure that the patient is not sensitized to these antibiotics.**

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time this product is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

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