



**Thermo Scientific HERASAFE KS and KSP**  
Class II Biological Safety Cabinets

**Setting new Standards for Safety**  
Containment · Comfort · Convenience

# Thermo Scientific Herasafe KS and KSP Setting new standards for safety

Thermo Scientific biological safety cabinets are backed by our global reputation and commitment to provide the safest and most reliable biological safety cabinets available for every lab and customer application.

We deliver decades of market leadership and technology innovations with unmatched quality.

Our philosophy is grounded in compliance with the most stringent regulatory standards in the world that ensure your safety.

**Thermo Scientific – moving science forward.**

Decades of Market  
Leadership



• **EN 12469**

Our complete line of Herasafe® KS and KSP products are independently tested and certified to EN 12469 standard for Class II biological safety cabinets by TÜV Nord, the leading testing agency in Europe.



• **DIN 12980**

Our Herasafe® KSP cabinets with single piece worktops are independently tested and certified to the DIN 12980 standard "Safety Cabinets for Handling Cytotoxic Substances", requiring the highest personnel, product and environmental protection.



• **NSF/ANSI 49**

The Herasafe® KS line is also compliant with NSF/ANSI 49 for Class II biological safety cabinets. NSF/ANSI 49 is recognized as the preeminent safety cabinet standard in North America, Latin America, and parts of Asia Pacific.



**Thermo Scientific biological safety cabinet evolution**

*We launch our first line of clean air products*

*NSF publishes first ANSI 49 standard for Class II biological safety cabinets*

*First NSF-certified Thermo Scientific biological safety cabinet is launched*

*We introduce the world's first TÜV certified biological safety cabinet*

*We launch the world's first safety cabinet with motorized front window and aerosol tight window sealing*

*We launch our first LNE certified safety cabinet*

*We launch the first DIN 12980 certified safety cabinet designed for preparing and handling cytotoxic drugs*

1965

1976

1978

1985

1993

1995

1996

# Our experience brings innovation to you

Herasafe KS and KSP Class II biological safety cabinets deliver maximum safety, comfort and convenience. Featuring an advanced design for easy operation and maximum safety, our Herasafe KS and KSP cabinets may be used wherever protection is of paramount importance in safety levels 1 to 3, including research, pharmaceutical quality control labs, pharmacies, and clinics/hospitals.

- **Superior Personnel and Product Protection** with best-in-class air flow and filter technology
- **Exceptionally Comfortable** ergonomic design prevents user fatigue and promotes safe working habits
- **Maximized Energy Efficiency** for environmental protection and significant cost savings

CONTAINMENT  
COMFORT  
CONVENIENCE

ENT

Europe harmonizes EN 12469, "Performance Criteria for Microbiological Safety Cabinets"  
We launch the **world's first** safety cabinet with innovative brushless DC motor technology

2000

We launch the advanced Thermo Scientific Herasafe KS and KSP biological safety cabinets

2008



Model KS 12 shown with optional electric stand, gas and service valves in side and rear walls.

# Best-in-class containment

The most important factor in selecting a biological safety cabinet is the security of particulate containment. Unique technologies adapted for Thermo Scientific safety cabinets ensure safe air filtration, avoiding sample contamination or environmental hazards.

# CONTAINMENT

*Second to None  
in Containment*

**H14 HEPA filters** guarantee maximum filtration efficiency of 99.995% at the most penetrating particle size (MPPS) for KS and KSP models. The supply and exhaust HEPA filters ensure superior contamination prevention while maximizing filter life.

**Motorized front window** offers comfortable working conditions at a 250 mm working aperture for the KS, and 200 mm for the KSP. The window height can be easily adjusted with only a finger tip.

**Aerosol-tight window sealing** offers additional personnel and sample protection when the KS or KSP biological safety cabinet is set to stand-by mode.

# KS



Model KS 12 shown with optional electric stand, gas and service valves in side and rear walls.

# Exceptional design for personnel, product and environmental protection

The Herasafe KSP meets the latest European and worldwide standards for applications with stringent safety requirements. Robust construction and excellent ergonomic design are combined with significantly more effective air filtration, resulting in exceptional protection for personnel, product and environment.

## 100,000 times better air filtration

An additional segmented H14 HEPA filter system is used directly below the work area of the Herasafe KSP, protecting inner air ducts, plenums and fans from contamination and providing a cumulative filtration efficiency of 99.9999975% to the downflow and exhaust airflow.

**The Herasafe KSP is ideal for pharmacies preparing cytotoxic drugs and in pharmaceutical control labs.**

**SmartFlow™** uses digital technology for KS and KSP models to maintain constant airflow during normal filter loading or temporary airflow obstruction from foreign objects.

**Pressure sensor monitoring** ensures safe airflow across the entire work surface of KS and KSP cabinets. Independent alarm systems immediately alerts to any airflow obstruction.

5 Valdymas

6. Informacinis ekranas

**Interior eye level display** at the rear wall of KS and KSP models allows convenient "at a glance" monitoring of safety parameters, including:

- Performance Factor "PER"\*
- downflow air velocity
- cabinet malfunctions
- time of day
- operating time of optional UV-C lights
- total operating hours
- timer for delayed start
- stop-watch function

**\*Performance factor "PER"** provides important information on the cabinets' current safety performance:

- 60-100, cabinet is safe
- 30-59, cabinet is safe but time for service
- 0-29, cabinet should not be used

**Interlocking H14 HEPA filters** under the KSP work surface provide easy and safe HEPA filter changes in the work area under safe negative pressure.



Model KSP 12 shown with, optional gas and service valves in side and rear walls.



# Unmatched ergonomic design

A comfortable cabinet design promotes user safety. The 10° sloped front allows an ergonomic, fatigue-free working posture. With its frameless front sash and side glass windows, the Herasafe KS and KSP offer an unobstructed view of the complete work area from all angles.

## 3. Priekinis langas

### Adjustable stand options for convenient working positions

Available with electric or manual adjustment, our stands are equipped with levelling feet to securely position the cabinet.

### Hand-held control for ease of use

Our hand-held control is designed for use inside the cabinet and eliminates the need to change your gloves before changing cabinet settings.

## Second to None in Comfort

## 7: Distancinio valdymo pultelis



The optional manual adjustable stand for Herasafe KS can be set from 750 - 950 mm working height in 50 mm increments. The optional foot rest allows a convenient seated working position.



The optional electric stand for Herasafe KS is adjustable with an easy to use controller. You can change the height infinitely from a seated to a standing working position.



The full working height of 950 mm provides a very comfortable standing working position, even for tall workers.

The electric stand with cable remote control provides comfortable height adjustment.



**Bright, open chamber**

The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.



Secure and comfortable armrests enhance your comfort during extended work sessions and prevent carpal-tunnel syndrome. Located just above the intake grill, they also ensure that the air vents are not inadvertently obstructed.

250 mm working aperture on the KS allows most comfortable handling of your samples.

780 mm interior height of the working chamber provides high flexibility for your work inside the cabinet.

ORT

# Innovative design for convenient sample preparation

Thermo Scientific safety cabinets are the perfect choice for labs seeking the highest levels of work efficiency. Our standard cabinets offer extremely low vibration, making them perfect for use in sensitive weighing applications.\* Additional options and accessories provide flexibility to configure the workspace to exactly match your application for enhanced work efficiency.

## *Second to None in Convenience*

### **Customized work surfaces**

Scratch-free, high-quality stainless steel work surfaces are available as single or segmented modules and are easy to clean or autoclave. The indented work modules easily trap spilled liquids.



**Optional vibration damped work surface section with granite plate for maximized stability.**

**Standard worksurface is suitable for sensitive weighing applications.**

\*The standard Herasafe KS and KSP configuration is proven to work with microbalances down to  $10^{-5}$  readings, and depending upon installation environment, may work down to  $10^{-7}$  readings. Optional vibration-dampened worksurfaces are available to provide additional stability.


# CONVENIENCE

## Easy cleaning

The hinged front window can be fully opened for easy and thorough cleaning of all interior surfaces. Rounded edges inside the cabinet facilitate effective disinfection and maximize contamination prevention.

## Efficient disinfection

Powerful cross beam UV-irradiation in both side walls (optional) illuminates the entire working area without causing shadows. This ensures thorough disinfection of the complete chamber for maximum sample protection.



Schematic of Unique Cross Beam UV-irradiation.

Optional UV-irradiator for chamber disinfection can be mounted on the optional hanging bar.

The fully opened front window allows convenient posture for thorough cleaning and disinfection.

The large working aperture provides full access to the work chamber for loading and unloading of large items and accessories like our optional UV light on stand.



# Advanced technology for maximum energy efficiency

## 10. Varikliai

### Superior DC motor technology

We introduced the first innovative DC motor technology to the biological safety cabinet market. DC motors offer significant environmental advantages over traditional AC motors.

### Minimized heat emittance to the lab

DC motors are designed to use energy in an efficient manner, not only consuming less energy to operate the unit, but emitting less heat back into the lab. Lower heat output also means a more comfortable working environment and lower air conditioning costs.

### Maximized cost savings


Combined with the Night-Set-Back mode (reduced blower speed when the window is closed), the Herasafe KS and KSP consume hundreds of kilowatts less energy than competing cabinets. Not only is our Herasafe KS and KSP the preferred cabinet for enhanced containment, comfort, and convenience – but it saves your lab money and reduces your lab’s overall carbon footprint.

*Good for you, good for the environment*

## Less energy consumption delivers

- reduced energy costs
- lower heat output
- reduced air conditioning costs
- comfortable working environment
- environmental protection

## 11. Elektros sunaudojimas



	Herasafe	Alternative A	Alternative B
Operational Energy (W)	210	810	1120
Reduced Flow (W)	40	710	1020
Operational Hours per year (5 days x 8 hours x 52 weeks)	2080	2080	2080
Non-operational hours	6656	6656	6656
Operational kW-hours	437	1685	2330
Reduced Flow kW-hours	266	4726	6789
Total Annual kW-hrs	703	6411	9119


\* Data provided above is based on internal performance testing and has been independently verified in customer applications.

## Options and Accessories

Herasafe KS and Herasafe KSP Options and Accessories				
Order No.	Description	Factory Installed Option	Accessory	Image
<b>Country versions with local sockets and plugs</b>				
51900300	Country version Switzerland	x		
51900303	Country version UK	x		
51900306	Country version Italy	x		
51900448	Country version France	x		
51900449	Country version Australia	x		
51900481	Country version Denmark	x		
51900771	Country version B/CZ/SK/PL	x		
51900900	Country version China	x		
<b>Chamber equipment</b>				
51900648	Powder-coated rear wall instead of stainless steel	x		
51900644	Left side wall stainless steel plate instead of glass	x		
51900645	Right side wall stainless steel plate instead of glass	x		-
51900646	Left side wall polycarbonate instead of glass	x		
51900647	Right side wall polycarbonate instead of glass	x		
51900651	KS/KSP 9 without worksurface (required for installation of non-standard worksurface options)	x		
51900652	KS/KSP 12 without worksurface (required for installation of non-standard worksurface options)	x		
51900653	KS/KSP 15 without worksurface (required for installation of non-standard worksurface options)	x		
51900654	KS/KSP 18 without worksurface (required for installation of non-standard worksurface options)	x		
51900786	Battery backup for motorized front window	x		
<b>Exhaust transitions</b>				
50078791	Direct Duct Exhaust Transition for 0.9 m cabinet. Connection diameter 200 mm.		x	1
50073771	Direct Duct Exhaust Transition for 1.2 and 1.5 m cabinet. Connection diameter 200 mm.		x	1
50078790	Direct Duct Exhaust Transition for 1.8 m cabinet. Connection diameter 200 mm.		x	1
50078796	Thimble Duct Exhaust Transition for 0.9 m cabinet. Connection diameter 200 mm.		x	2
50073770	Thimble Duct Exhaust Transition for 1.2 and 1.5 m cabinet. Connection diameter 200 mm.		x	2
50078795	Thimble Duct Exhaust Transition for 1.8 m cabinet. Connection diameter 200 mm.		x	2
50078873	Exhaust cooling device for 1.2 m and 1.5 m cabinets		x	3
50078874	Exhaust cooling device for 1.8 m cabinet		x	3
<b>UV irradiation</b>				
51900642	Cross-beam UV Light Option in side walls (comes with protective cover)	x		4
50048696	Mobile UV Light on stand (must also order UV-socket 51900643 or 50109996)		x	5
50073962	Mobile UV Light, hanging from rear wall (may not be used with 0.9 m cabinet) (must also order UV-socket 51900643 or 50109996 and hanging bar 50073944)		x	6
51900643	UV socket for mobile UV light option	x		
50109996	UV socket for mobile UV light option, for retro-fit in existing cabinet		x	



1



2




3



4






5



6

## Options and Accessories

Herasafe KS and Herasafe KSP Options and Accessories					
Order No.	Description	Factory Installed	Option	Accessory	Image
<b>Utilities</b>					
50046015	Combustible gas tap for side wall, yellow (we recommend to also order 50074905 solenoid valve)			x	10
50045959	Noncombustible gas tap for side wall, black			x	10
50044678	Vacuum tap for side wall, grey			x	10
50044679	Water tap for side wall, green			x	10
50059017	Hose connection to fit through side wall access port, fits 10-13 mm diameter tubing			x	7
50066668	Hose connection to fit through side wall access port, fits 3 mm diameter tubing			x	8
51900649	Rear wall 3 service taps, left side of cabinet	x			11
51900711	Rear wall 3 service taps, right side of cabinet	x			11
51900722	Rear wall 6 service taps total, 3 on left and 3 on right side of cabinet (may not be used with 0.9m cabinet)	x			11
50068235	Combustible gas tap for rear wall, yellow (we recommend to also order 50074905 solenoid valve)			x	9
50068234	Noncombustible gas tap for rear wall, black			x	9
50068232	Vacuum tap for rear wall, grey			x	9
50068233	Water tap for rear wall, green			x	9
51900650	Magnehelic pressure gauge	x			
50074905	Solenoid valve for combustible gas tap			x	
51900641	Drain pan valve (standard on 120V NSF version)	x			
<b>Filter</b>					
51900657	Additional exhaust HEPA filter for all models, factory installed	x			16
51900658	Additional exhaust Carbon filter for all models, factory installed	x			16
50078857	Additional exhaust HEPA filter for 0.9 m cabinet, field installed			x	14
50073772	Additional exhaust HEPA filter for 1.2 and 1.5 m cabinet, field installed			x	14
50078856	Additional exhaust HEPA filter for 1.8 m cabinet, field installed			x	14
50078783	Additional exhaust Carbon filter for 0.9 m cabinet, field installed			x	
50073773	Additional exhaust Carbon filter for 1.2 and 1.5 m cabinet, field installed			x	
50078784	Additional exhaust Carbon filter for 1.8 m cabinet, field installed			x	
50076846	Grid pre-filter for 0.9 and 1.2 m (Herasafe KS only)			x	
50076921	Grid pre-filter for 1.5 and 1.8 m (Herasafe KS only)			x	
51900781	Activated carbon filter that fits under 0.9 m cabinet worksurface (KSP models only)	x			15
51900782	Activated carbon filter that fits under 1.2 m cabinet worksurface (KSP models only)	x			15
51900783	Activated carbon filter that fits under 1.5 m cabinet worksurface (KSP models only)	x			15
51900784	Activated carbon filter that fits under 1.8 m cabinet worksurface (KSP models only)	x			15
<b>User convenience equipment</b>					
50073943	Hanging bar with set of 6 hooks (for 0.9 m cabinet only)			x	13
50073944	Hanging bar with set of 12 hooks (may not be used with 0.9 m cabinet)			x	like 13
50078900	Utility basket for use with hanging bars 50073943 or 50073944			x	12
50078901	2 Utility baskets including additional hanging bar			x	like 12
50073663	Replacement stainless steel armrest (set of 2)			x	17




## Options and Accessories

### Herasafe KS and Herasafe KSP Options and Accessories

Order No.	Description	Factory Installed Option	Accessory	Image
<b>Work surfaces</b>				
50073666	One-piece flat worksurface module 150 mm wide		x	18
50072685	One-piece flat worksurface module 300 mm wide		x	19
50073667	One-piece flat worksurface module 600 mm wide		x	20
50076398	One-piece flat worksurface module 900 mm wide		x	21
50073668	One-piece indented worksurface module 600 mm wide		x	like 22
50073669	One-piece indented worksurface module 900 mm wide		x	like 22
50073670	One-piece indented worksurface module 1200 mm wide		x	22
50073671	One-piece indented worksurface module 1500 mm wide		x	like 22
50078750	One-piece indented worksurface module 1800 mm wide		x	like 22
50073757	Vibration damped work surface 300 mm width. Constructed in stainless steel with granite plate. For weighing equipment and microscopy (not for use with Herasafe KSP)		x	24
50076162	Vibration damped work surface 600 mm width. Constructed in stainless steel with granite plate. For weighing and microscopy equipment (not for use with Herasafe KSP)		x	like 24
50084339	Vibration damped granite block 300 mm width. For weighing equipment and microscopy equipment for triple filter cabinets only. (To be placed on top of standard worksurface)		x	like 24
50084350	Vibration damped granite block 600 mm width. For weighing and microscopy equipment for triple filter units only. (To be placed on top of standard worksurface)		x	like 24
50078875	Mini sink tabletop insert 300 mm, requires drain pan valve 51900641 to also be fitted. Holds ~3 liters.		x	23
<b>Stand options for Herasafe KS (Manual adjustable floorstand is included with Herasafe KSP)</b>				
50116441	Fixed height floor stand 0.9 m, 750 mm working height		x	25
50109309	Fixed height floor stand 1.2 m, 750 mm working height		x	25
50116442	Fixed height floor stand 1.5 m, 750 mm working height		x	25
50109311	Fixed height floor stand 1.8 m, 750 mm working height		x	25
50116443	Manual adjustable floor stand 0.9 m, 750 to 950 mm working height in 50 mm increments		x	like 25
50109312	Manual adjustable floor stand 1.2 m, 750 to 950 mm working height in 50 mm increments		x	like 25
50116444	Manual adjustable floor stand 1.5 m, 750 to 950 mm working height in 50 mm increments		x	like 25
50109313	Manual adjustable floor stand 1.8 m, 750 to 950 mm working height in 50 mm increments		x	like 25
50117449	Electric adjustable floor stand 0.9 m with cable remote control, 750 to 950 mm working height		x	27
50109314	Electric adjustable floor stand 1.2 m with cable remote control, 750 to 950 mm working height		x	27
50117450	Electric adjustable floor stand 1.5 m with cable remote control, 750 to 950 mm working height		x	27
50109325	Electric adjustable floor stand 1.8 m with cable remote control, 750 to 950 mm working height		x	27
50051982	Footrest bar 0.9 m (not to be used with electric stand)		x	26
50051983	Footrest bar 1.2 m (not to be used with electric stand)		x	26
50051984	Footrest bar 1.5 m (not to be used with electric stand)		x	26
50051985	Footrest bar 1.8 m (not to be used with electric stand)		x	26
50109977	Castors for fixed height floor stands only (not to be used with electric stands)		x	



## Schematic drawing Herasafe KS




		A
seated work position minimum	mm	680
seated work position minimum	inch	26,77
standing work position maximum	mm	880
standing work position maximum	inch	34,65

		B	C	D	E	F
KS 9	mm	81	491	1000	430	900
	inch	3,19	19,33	39,37	16,93	35,43
KS 12	mm	228	644	1300	583	1200
	inch	8,97	25,35	51,18	22,95	47,24
KS 15	mm	478	644	1600	583	1500
	inch	18,82	25,35	62,99	22,95	59,05
KS 18	mm	474,5	951	1900	890	1800
	inch	18,68	37,44	74,80	35,04	70,86

Attention:  
With additional gas valves in the side walls, you must plan with more width (dimension D).

## Schematic drawing Herasafe KSP



		A
seated work position minimum	mm	680
seated work position minimum	inch	26,77
standing work position maximum	mm	880
standing work position maximum	inch	34,65

		B	C	D	E	F
KSP 9	mm	81	491	1000	430	900
	inch	3,19	19,33	39,37	16,93	35,43
KSP 12	mm	228	644	1300	583	1200
	inch	8,97	25,35	51,18	22,95	47,24
KSP 15	mm	478	644	1600	583	1500
	inch	18,82	25,35	62,99	22,95	59,05
KSP 18	mm	474,5	951	1900	890	1800
	inch	18,68	37,44	74,80	35,04	70,86

Attention:  
With additional gas valves in the side walls, you must plan with more width (dimension D).

# SPECIFICATIONS AND ORDERING INFORMATION

## Thermo Scientific Herasafe KS, Class II Biological Safety Cabinet



Model	Unit	KS 9	KS 12	KS 15	KS 18
<b>Order number EN models (230V)</b>		51022481	51022482	51022483	51022484
<b>Order number NSF models (120V)</b>		51023770	51023772	51023774	51023776
<b>Description: Biological Safety Cabinet</b>		Class II (NSF: Type A2)	Class II (NSF: Type A2)	Class II (NSF: Type A2)	Class II (NSF: Type A2)
		0.9 m (3 ft.)	1.2 m (4 ft.)	1.5 m (5 ft.)	1.8 m (6 ft.)
<b>Dimensions</b>					
Exterior Dimensions H/W/D	mm	1586/1000/800	1586/1300/800	1586/1600/800	1586/1900/800
Exterior Dimensions H/W/D	inches	62,2/39,4/31,5	62,2/51,2/31,5	62,2/63/31,5	62,2/74,8/31,5
Usable Work Area/Interior H/W/D	mm	780/900/465	780/1200/465	780/1500/465	780/1800/465
Usable Work Area/Interior H/W/D	inches	30,7/35,4/18,3	30,7/47,2/18,3	30,7/59,1/18,3	30,7/70,9/18,3
Work Surface Height (with adjustable stand)	mm/inches	750-950/30-37	750-950/30-37	750-950/30-37	750-950/30-37
Work Surface Height (with fixed stand)	mm/inches	750/30	750/30	750/30	750/30
Working Height of front window	mm/inches	250/10	250/10	250/10	250/10
Maximum Opening Height of front window	mm/inches	773/30,4	773/30,4	773/30,4	773/30,4
Shipping Dimensions H/W/D <sup>1)</sup>	mm	1760/1110/925	1760/1410/925	1760/1710/925	1760/2010/925
Shipping Dimensions H/W/D <sup>1)</sup>	inches	69,3/43,7/36,4	69,3/55,5/36,4	69,3/67,3/36,4	69,3/79,1/36,4
<b>Weight</b>					
Net Weight	kg/lbs	170/375	200/441	230/507	280/617
Shipping Weight <sup>1)</sup>	kg/lbs	190/419	225/496	260/573	315/694
Max. load capacity of one-piece work surface	kg/lbs	50/110	50/110	50/110	50/110
Max. load capacity of single segmented worksurface per 30 cm segments <sup>4)</sup>	kg/lbs (Total load)	25/55 (50/110)	25/55 (50/110)	25/55 (75/165)	25/55 (75/165)
<b>Ventilation system</b>					
Exhaust/inflow air volume EN models	m <sup>3</sup> /h	366	490	613	737
Exhaust/inflow air volume NSF models	CFM	255	341	427	513
Exhaust volume, thimble connection EN models	m <sup>3</sup> /h	476	637	797	958
Exhaust volume thimble connection NSF models	CFM	322	444	556	668
<b>Heat emission at 25 °C room temperature</b>					
Operation mode (non-vented)	W/BTU/hr	170/580	210/717	275/938	350/1194
Reduced flow mode (non-vented)	W/BTU/hr	40/136	40/136	70/239	70/239
<b>Filter Specification</b>					
Supply/exhaust air filter		H14 HEPA EN 1822, 99.995% at the most penetrating particle size (MPPS)			<b>8. Filtrai</b>
<b>Ergonomics</b>					
Noise pressure level EN/NSF <sup>2)</sup>	dB (A)	58/65	58/65	60/65	60/65
Lighting power	lx/ftc	>800/74	>850/79	>1250/116	>1300/121
<b>Electrical data</b>					
Voltage EN/NSF	V	230/120	230/120	230/120	230/120
Frequency	Hz	50/60	50/60	50/60	50/60
Energy Consumption, EN models, Nominal set point <sup>3)</sup>	W	170	210	275	350
Energy Consumption, NSF models, Nominal set point <sup>3)</sup>	W	210	275	350	435
Energy Consumption, reduced flow mode	W	40	40	70	70
Max Current Consumption	A	7,3	7,3	8,7	8,7
Protection Class		I/IP 20	I/IP 20	I/IP 20	I/IP 20
<b>Features</b>					
Receptacles	One duplex on each side wall. The receptacles have a load capacity of up to 5 A and are protected with T 5 A fuses. When all receptacles are in use simultaneously they must not exceed the maximum total load capacity of 5 A				
Access Ports	4 standard (2 on each side wall)				
Service Valves	Up to 4 may be installed through access ports (2 on each side wall) Up to 6 more may be factory installed on rear wall (optional)				

2: darbinis plotis

4: Apšvietimas

<sup>1)</sup> export packaging dimensions and weight on request

<sup>2)</sup> tested according to EN 12469 and NSF/ANSI 49

<sup>3)</sup> Clean filters, fans operating at nominal set point, interior lighting activated

<sup>4)</sup> The total load on overall working area must not exceed the indicated weight

# SPECIFICATIONS AND ORDERING INFORMATION



## Thermo Scientific Herasafe KSP, Class II Biological Safety Cabinet

Model	Unit	KSP 9	KSP 12	KSP 15	KSP 18
<b>Order No.</b>		51023606	51023607	51023608	51023609
<b>Description: Biological Safety Cabinet</b>		Class II, with triple filter system	Class II, with triple filter system	Class II, with triple filter system	Class II, with triple filter system
<b>Dimensions</b>					
Exterior Dimensions H/W/D	mm	1580/1000/870	1580/1300/870	1580/1600/870	1580/1900/870
Exterior Dimensions H/W/D	inches	62,2/39,4/34,3	62,2/51,2/34,3	62,2/63/34,3	62,2/74,8/34,3
Usable Work Area/Interior H/W/D	mm	780/900/465	780/1200/465	780/1500/465	780/1800/465
Usable Work Area/Interior H/W/D	inches	30,7/35,4/18,3	30,7/47,2/18,3	30,7/59,1/18,3	30,7/70,9/18,3
Work Surface Height (adjustable)	mm/inches	750-950/30-37	750-950/30-37	750-950/30-37	750-950/30-37
Working Height of front window	mm/inches	200/8	200/8	200/8	200/8
Maximum Opening Height of front window	mm/inches	773/30,4	773/30,4	773/30,4	773/30,4
Shipping Dimensions H/W/D <sup>1)</sup>	mm	1760/1110/925	1760/1410/925	1760/1710/925	1760/2010/925
Shipping Dimensions H/W/D <sup>1)</sup>	inches	69,3/43,7/36,4	69,3/55,5/36,4	69,3/67,3/36,4	69,3/79,1/36,4
<b>Weight</b>					
Net Weight	kg/lbs	200/441	240/529	280/617	330/728
Shipping Weight <sup>1)</sup>	kg/lbs	220/485	265/584	310/683	365/805
Max. load capacity of one-piece work surface	kg/lbs	50/110	50/110	50/110	50/110
Max. load capacity of single segmented worksurface per 30 cm segments <sup>4)</sup>	kg/lbs (Total load)	25/55 (50/110)	25/55 (50/110)	25/55 (75/165)	25/55 (75/165)
<b>Ventilation system</b>					
Exhaust/inflow air volume	m <sup>3</sup> /h/CFM	288/170	386/227	483/284	580/341
Exhaust volume, thimble connection	m <sup>3</sup> /h/CFM	375/221	501/295	628/369	754/444
<b>Heat emission at 25 °C room temperature</b>					
Operation mode (non-vented)	W/BTU/hr	200/682	240/819	305/1041	420/1433
Reduced flow mode (non-vented)	W/BTU/hr	40/136	40/136	70/239	70/239
<b>Filter Specification</b>					
Primary/Supply/exhaust air filters		H14 HEPA EN 1822, 99.995% at the most penetrating particle size			
<b>Ergonomics</b>					
Noise pressure level <sup>2)</sup>	dB (A)	56	56	58	58
Lighting power	lx/fc	>800/74	>850/79	>1250/116	>1300/121
<b>Electrical data</b>					
Voltage	V	230	230	230	230
Frequency	Hz	50/60	50/60	50/60	50/60
Energy Consumption, operating set point <sup>3)</sup>	W	200	240	305	420
Energy Consumption, reduced flow mode	W	40	40	70	70
Max Current Consumption	A	7,3	7,3	8,7	8,7
Protection Class		I/IP 20	I/IP 20	I/IP 20	I/IP 20
<b>Features</b>					
Receptacles	One duplex on each side wall. The receptacles have a load capacity of up to 5 A and are protected with T 5 A fuses. When all receptacles are in use simultaneously they must not exceed the maximum total load capacity of 5 A				
Access Ports	4 standard (2 on each side wall)				
Service Valves	Up to 4 may be installed through access ports (2 on each side wall) Up to 6 more may be factory installed on rear wall (optional)				

<sup>1)</sup> export packaging dimensions and weight on request

<sup>2)</sup> tested according to EN 12469

<sup>3)</sup> Clean filters, fans operating at nominal set point, interior lighting activated

<sup>4)</sup> The total load on overall working area must not exceed the indicated weight

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

[www.thermo.com/bsc](http://www.thermo.com/bsc)

**North America:** USA/Canada +1 866 984 3776 (866-9-THERMO)

**Europe:** Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

**Asia:** China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220, Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940 or +33 2 2803 2180

BRO-LECA-HSAFEKS-1208

**Thermo**  
SCIENTIFIC



**Thermo Scientific**  
Class II, Type A2  
Biological Safety Cabinets  
and Clean Benches

# protection

that never takes a day off

**Thermo**  
SCIENTIFIC

# protection

## that never takes a day off

With Thermo Scientific™ biological safety cabinets, the certified performance and protection you receive on Day 1 stays with you every day. Not true with ordinary cabinets.

### 9. Oro srauto greitis

The difference is our design: Thermo Scientific™ SmartFlow™ technology features dual-DC motors to automatically balance the cabinet inflow and downflow air velocities in real time – even as the filters load. That means outstanding user and sample protection you never have to think about. Plus, our Digital Airflow Verification (DAVe) alarm signals any out-of-spec conditions, for added assurance.

With our proven reliability, ergonomics and energy efficiency, the ideal choice is the one you can trust completely. And not just on Day 1.

**5** YEAR

**SmartPromise  
Comprehensive Warranty\***  
delivers peace of mind!

### technology leadership

Technical advancements include:

1965

Our first line of  
clean-air products

1978

Thermo Scientific  
first NSF-certified  
BSC

1985

World's first  
TÜV-certified  
BSC

1993

World's first BSC  
with motorized  
front window and  
aerosol-tight  
sealing

\* 5 year warranty only available on biological safety cabinets in the United States and Canada.

## global market vision

Our design philosophy is grounded in compliance with current worldwide safety standards.

## Globally Certified Biological Safety Cabinet

Not all biological safety cabinets (BSCs) are manufactured to the same safety standards. Our design philosophy is grounded in compliance with the recognized worldwide safety standards.

There are three validated global standards for BSC performance; EN12469 (European Standard), NSF/ANSI 49 (American Standard) and YY0569 (Chinese Standard). Regions without a selected BSC standard are often influenced by one or more of the major global standards.

## quality built, premium features

Our BSCs feature quality engineering for outstanding reliability. Discover why laboratories around the globe rely on Thermo Scientific BSCs every day.



### NSF/ANSI 49

Thermo Scientific 1300 series Class II biological safety cabinets are independently tested and certified to NSF/ANSI 49, by NSF International, a leading testing agency in the USA.



### EN 12469

Thermo Scientific™ Herasafe™ KS 1.2 m cabinet and MSC-Advantage™ 0.9 m and 1.2 m cabinets are independently tested and certified to EN 12469 by LNE, a leading testing agency in France.



### EN 12469

Our Class II biological safety cabinets are independently tested and certified to EN 12469 standard for Class II biological safety cabinets by TÜV Nord, a leading testing agency in Europe.



### DIN 12980

Herasafe KSP and Thermo Scientific™ Maxisafe™ 2020 cabinets with single-piece worktops are independently tested and certified to the DIN 12980 standard "Safety Cabinets for Handling Cytotoxic Substances," for applications requiring the highest personal, sample and environmental protection.



### YY0569

1300 series Class II biological safety cabinets are independently tested and certified to YY0569.

1995

LNE-certified BSC units available

1996

DIN 12980-certified BSC for cytotoxic drugs

2002

World's first BSC with innovative brushless dual-DC motor technology (Herasafe KS BSC)

2007

Launch of the energy efficient 1300 Series A2 BSC

2011

Launch of Thermo Scientific BSC with SmartPort tubing organizer

2013

Launch of Thermo Scientific Heraguard ECO clean bench with DC motor technology

# the right cabinet

## for you



Thermo Scientific Herasafe KS  
biological safety cabinet  
**outstanding protection, comfort and  
performance for advanced applications**



Thermo Scientific 1300 Series A2  
biological safety cabinet  
**exceptional efficiency, safety and  
value for daily applications**

Thermo Scientific™ Heraguard™ ECO\*  
clean bench  
**Ideal when sample protection is required.  
Due to DC motors, Heraguard ECO clean  
benches use less energy to operate so  
there is less heat emitted into the lab\***



\*Based on data from internal testing and The University of Michigan field study, published in American Biotechnology Laboratory.

# outstanding performance

for every application

		Herasafe KS	1300 Series A2	Heraguard ECO			
Types of Protection	Sample	•	•	•			
	Operator	•	•				
	Environmental	•	•				
	Cross-contamination	•	•	•			
Recommended Applications	Tissue culture	•	•	•			
	Biosafety level 1	•	•				
	Biosafety level 2	•	•				
	Biosafety level 3	•	•				
Features	Efficiency	Low-energy DC motors for operational cost savings and low heat emission	•	•	•		
		Automated reduced flow mode for enhanced energy efficiency	•	•			
		Long-life DC motors	•	•	•		
	Containment	Automatic airflow balancing for safety in between certifications (SmartFlow)	•	•	•		
		Real-time monitoring of HEPA filter efficiency (SmartFlow indicator / PER factor)	•	•	•		
		Night Set-Back Mode (maintains sterility)	•	•			
		Timed UV system for long bulb life	•	•	•		
		Real-time airflow alarms to notify user of impaired safety	•	•	•		
		3. Priekinis stiklas, tarpinė	Aerosol-tight sealing front window	•			
			Negative pressure side wall plenums	•	•		
			All metal plenums	•	•	•	
		2. Darbinis paviršius	Convenience	Easy access to interior of front sash for cleaning, without tools	•	•	
				Motorized front window for finger-tip opening and closing	•		
	316 stainless steel interior for enhanced resistance to corrosion			•		•	
	304 stainless steel work surface				•		
	Transparent side windows to maximize light and visibility inside the unit		•		•		
SmartClean window design			•				
SmartClean Plus window design	•						
SmartPort tubing/cable organizer			•	•			
Comfort	Forearm support armrest design to maximize comfort	•	•				
	Angled sash to encourage natural posture	•	•				
	Frameless front sash for maximum visibility	•	•				
	8-inch window option		•				

# innovations

## in containment and efficiency

### Brushless DC motors since 2002

We pioneered the first innovative brushless DC motor technology. Now available in our full line of BSCs and the Heraguard ECO clean benches, DC motor advancements have redefined and transformed BSC energy efficiency, safety, performance and reliability.

#### Use less energy, reduce costs

DC motors use less energy to operate, so there is less heat emitted back into the lab – delivering better comfort and lower air conditioning costs. On average, our BSCs use 75% less energy than traditional AC motor cabinets.\*

#### Exceptional reliability

Our dual-motor BSC design offers an expected lifetime of more than 11 years – at least double the published life of traditional AC motor and single DC motor systems – resulting in significant savings on motor replacements over the life of the cabinet.

#### Long filter life

Larger mini-pleat H14 HEPA filters and balanced airflow, coupled with eight-inch work aperture, results in higher filtration efficiency and longer filter life when compared to previous models.




### check the watts!

**Thermo Scientific biological safety cabinets lead the industry in operational cost efficiency with energy consumption as low as 150W.\***

### optimizing total cost of operation

Thermo Scientific SmartFlow airflow compensation and balanced filters have led to BSCs that are designed to be safer, easier to use and less costly to operate.\*\*



\* Based on data from internal testing and The University of Michigan Field Study of Class II Biological Safety Cabinet Energy Consumption Costs, published in American Biotechnology Laboratory.

\*\* Compared to previous models.


**The more you use it, the more you save!**

## advanced design for peace of mind

A BSC's primary function is to provide you cleanliness and containment. Available on our full range of BSCs, **Thermo Scientific SmartFlow, Digital Airflow Verification and Night Set-Back Mode** are innovations designed to increase safety and minimize sample contamination.

- ▶ H14 HEPA filters increase filtration efficiency at 99.995% of MPPS
- ▶ Solid plenums for long life and lasting containment
- ▶ Negative-pressure, dual side walls protect from leaks to the exterior environment

**SmartFlow:** automated balance of RPM, voltage, and airflow



### SmartFlow designed to ensure safety

#### Unique airflow system takes safety and containment to new levels

- Smart DC motors monitor and control fan speed in real time to maintain constant airflow during filter loading or temporary obstruction.
- Unlike single motor designs, the independent supply and exhaust blowers automatically balance the downflow and inflow velocities to maintain user and sample protection.



### Digital Airflow Verification (DAVe) validates protection

#### Monitoring system alerts user when containment is compromised

- Independent sensors detect pressure changes across the exhaust and downflow plenums.
- Alarm signals when a 20% change in inflow/exhaust or downflow occurs.
- Airflow velocities are displayed on the control panel for monitoring and recording.

5. Valdymas

### Night Set-Back Mode saves energy, maintains sterility

#### When the front sash is closed, our intelligent speed control automatically reduces blower speed

- Continuous airflow helps to maintain a sterile working environment, even when the cabinet is not in use.
- Less air flows through the HEPA filters, extending filter life and decreasing energy consumption.

# designed for comfort and convenience

Thermo Scientific BSCs offer an efficient workspace, ergonomics, and added conveniences that provide you outstanding productivity.

## Ergonomic design

- Sloped front window is designed to reduce user fatigue
- Bright, open chamber improves visibility
- Quiet airflow design offers noticeably reduced noise level

## Ergonomically position arm rest

- Front air-grill arm rest extends user reach inside cabinet
- Provides forearm support to reduce strain
- Ensures clothing does not block front air intake grill
- Supports work midway between rear and front grill

## Thermo Scientific™ Ergolign™ saddle stool

- Exclusive seat design reduces back strain and fatigue
- Clean-room-ready, with 0.5 µm filtration system within the seat, trapping 99.9% of particles exhausted, fulfilling Class 100 requirements
- Fully adjustable seat height and tilt, using a patented icon plate mechanism
- Stores easily under the BSC to save floor space

work smarter, more efficiently

- ▶ Easy to use
- ▶ Easy to clean, without the use of tools
- ▶ Quiet operation

# smarter design



## Thermo Scientific™ SmartPort™ convenience

Clean and safe routing for vacuum tubing and cables through the side of the BSC for improved organization and work efficiency. (1300 Series A2 BSC only)

## Thermo Scientific™ SmartClean™ window design

Our unique window design easily lowers for thorough cleaning of inner surfaces and designed to protect users by maintaining inflow, even when the window is lowered. (1300 Series A2 BSC only)



## Thermo Scientific SmartClean Plus

Fully-open, front hinging window design allows convenient posture for thorough, easy cleaning and disinfection, with full access to BSC interior. (Herasafe KS BSC only)



## Outstanding UV decontamination

Optional UV light is programmable from 30 minutes to 24 hours, in 30-minute increments, extending bulb life and saving energy.

Powerful cross-beam UV irradiation in both side walls illuminates entire work area, designed to ensure thorough disinfection of the complete chamber (Herasafe KS BSC only).

# Thermo Scientific Herasafe KS biological safety cabinet

**outstanding protection, comfort and performance**  
for advanced applications

Featuring an advanced design for protection and ease of use, our Herasafe KS BSC enables quick, simple cleaning and fumigation.

## Containment

- **SmartFlow design** automates airflow compensation and balancing
- **Digital Airflow Verification (DAVe)** alerts user to compromised containment
- **H14 HEPA filters** remove airborne contaminants
- **Performance factor safety indicator** displays performance status
- **Cross beam UV disinfection system** illuminates entire work area without shadows
- **SampleGuard aerosol-tight window sealing** for added sample protection easy sealing fumigation

3:

## Comfort

- **Reduced noise level** is designed to enhance user focus
- **10° sloped front** for outstanding comfort and ergonomics
- **Spacious work area** with comfortable armrests
- **Transparent side-glass windows** maximize light and visibility

2:  
darbinis paviršius

## Convenience

- **SmartClean Plus** hinged-front window design for easy cleaning
- **Intuitive eye-level control panel** with handheld controls to assist with ADA compliance
- **Motorized frameless window sash** designed to ensure safe sash height at the push of a button
- **Scratch-free, 316 stainless steel work surfaces** are available as single or segmented modules and are easy to clean or autoclave
- **Programmable cross-beam UV light** extends bulb life and saves energy

**SampleGuard**  
aerosol-tight window sealing  
provides extra containment  
when window is closed



## Efficiency

- **DC motors** offer low energy consumption and greater than 11 years of use
- **Night Set-Back Mode** reduces airflow to 30% during low-usage periods compared to normal usage, with the window open
- **Longer filter life\***
- **Longer UV bulb life\*** programmable UV light extends bulb life

\*Compared to previous models.



Cross-beam UV lighting

Intuitive eye-level control panel

Transparent side-glass windows

Motorized sash

Optional universal piping/rear wall service taps

316 stainless steel interior

Customizable work surfaces

Optional electric adjustable height stand

## More choices for you.\*

- Choice of 3-, 4-, 5- or 6-foot widths
- Available with 10-inch work aperture
- Choice of electrical or manual adjustable height stands

# Thermo Scientific 1300 Series A2 biological safety cabinet

## optimized efficiency, safety and value

for your daily applications

Thermo Scientific 1300 Series Class II, Type A2 biological safety cabinets combine safety, energy efficiency and value for today's laboratory.

### Containment

- **SmartFlow design** automates airflow compensation and balancing
- **Digital Airflow Verification (DAVe)** alerts user to compromised containment
- **H14 HEPA filters** remove airborne contaminants
- **SmartFlow indicator** displays performance status

### Comfort

- **Low noise level** helps enhance user focus
- **10° sloped front** for outstanding comfort and ergonomics
- **Spacious work area** with comfortable armrests
- **Increased light intensity** for a brighter work environment

### Convenience

- **SmartClean** front window design
- **Intuitive control panel** displays safety and performance data
- **SmartPort** organizes tubing and cables
- **Programmable UV light** extends bulb life and saves energy

**Intuitive Control Panel:** Displays valuable safety and performance data with constant read-out of downflow/inflow velocities



**Airflow Read-Out & Hour Counter:**

Clearly displays realtime performance data, including inflow and downflow velocities and hours of operation

**SmartFlow Indicator:**

Visually demonstrates power and capacity to maintain protection

**Front Window Status:**

Visual and audible alarm indicates whether front window is in correct working position

**Night Set-Back Indicator:**

Displays reduced speed operation when front window is closed

**DAVe Airflow Indicator:**

Audible and visual alarms indicate when airflow is safe or restricted



**Thermo Scientific Ergolign Saddle Stool**  
Designed to relieve back pressure and optimize working posture during long work sessions

## More choices for you.

- Available in 3-, 4-, 5- or 6-foot widths
- Choice of 8- or 10-inch work aperture
- Choice of SmartCoat interior for reduced glare in workspace or seamless one-piece stainless steel
- Simple ordering with package configurations that include accessories (manual adjustable height stand, UV light and armrest)
- Optional electrical adjustable height stands

### 1300 Series A2 BSC

shown with optional manual floor stand, floor-stand hanging shelf, vacuum valve and rear-wall tubing clips (Bottle not included)

**Efficiency**

- **DC motors** offer low energy consumption and designed to offer an expected life of more than 11 years
- **Night Set-Back Mode** reduces airflow to 30% during low-usage periods
- **Longer filter life\***
- **Longer UV bulb life\***
- **Programmable UV light** extends bulb life

**Certifications:** All 1300 Series A2 models have been independently tested and certified to NSF/ANSI 49 for Class II biological safety cabinets

**Manufactured in a ISO 9001 facility**

\*Compared to previous models.

# Thermo Scientific Heraguard ECO\* clean bench

## high performance with optimized efficiency

sterility and functionality

NEW Thermo Scientific™ Heraguard™ ECO clean benches are ideal for laboratory applications where sample protection is required. Due to DC motors, Heraguard ECO clean benches use less energy to operate so there is less heat emitted into the lab.\*

### Containment

- **Thermo Scientific Smartflow design** compensates for filter loading and maintains airflow
- **SmartFlow indicator** displays performance status
- **H14 HEPA filters** remove airborne contaminants

### Comfort

- **Optional transparent hinged wind screen** designed to shield the user against the dry purified air flow.
- **Side windows** increase the available light inside the work area.

### Convenience

- **Prefilter** extends the HEPA filter life and is easily changed by hand without the need of tools
- **Intuitive control panel** displays safety and performance data
- **Easy HEPA filter exchange** from the front of the unit
- **Choice of work tops** in melamine or stainless steel (316L) for different applications and individual work space designs
- **Modular side windows** made of single pane safety glass with an access port for tubing or a media tap. Optional lacquered sheet steel plate with two utility ports and a SmartPort
- **Programmable UV light** extends bulb life and saves energy
- **Variety of sizes** 3,4,5 and 6ft, with the 6ft also available with an increased interior height of 950mm



### Efficiency

#### NEW DC motors

Now available with DC motors which offer over 60% lower energy consumption compared to previous model\*

#### Extended filter life

Smartflow compensates for filter loading and extends filter life

Intuitive control panel displays valuable safety and performance data

\* Based on data from internal testing and The University of Michigan field study, published in American Biotechnology Laboratory.



### UV disinfection (optional)

A UV disinfection light can be fitted into the top of the work space to disinfect the worktop.

### Transparent hinged wind screen (optional)

Designed to shield the user against the dry purified air flow. The wind screen is folded under the light housing and can be flipped down for use.

### Modular side windows

The side windows are made of single pane safety glass with an access port for tubing or a media tap. Optionally, a lacquered sheet steel plate with three holes can be used.



### Prefilter

The prefilter on the intake side can be changed by hand, without the need for tools. Prefilter exchange is made economical and easy.

### NEW Intuitive control panel

Displays safety and performance data.

### Easy HEPA filter exchange

The HEPA filter is easily removed from the front – a simple process that does not involve any disassembly or removal of equipment from the wall. This saves time during service visits.

### Choice of work tops

Standard melamine or optional stainless steel tops are available for different applications and individual work space designs.

# Thermo Scientific Herasafe KS Biological Safety Cabinet

## specifications

Size/Width			3-foot	4-foot	5-foot	6-foot
<b>Dimensions</b>	Exterior dimensions W x H x D	in	39.4 x 62.4 x 31.5	51.2 x 62.4 x 31.5	63.0 x 62.4 x 31.5	74.8 x 62.4 x 31.5
	Interior dimensions W x H x D	in	35.4 x 30.7 x 24.8	47.2 x 30.7 x 24.8	59.1 x 30.7 x 24.8	70.9 x 30.7 x 24.8
	Shipping dimensions <sup>1</sup> W x H x D	in	43.7 x 69.3 x 36.4	55.5 x 69.3 x 36.4	67.3 x 69.3 x 36.4	79.1 x 69.3 x 36.4
	Working height of front window	in	10			
	Maximum opening height of front window	in	30.4			
<b>Weight</b>	Net weight	lbs	375	441	507	617
	Shipping weight <sup>1</sup>	lbs	419	496	573	694
<b>Electrical Data</b>	Voltage	V	120			
	Frequency	Hz	60			
<b>Filter Specification</b>			H14 HEPA EN 1822, 99.995% MPPS			
<b>Certification</b>			NSF/ANSI 49, ETL, CE			
<b>Ergonomics and Utilities</b>	Lighting power	fc	>70	>70	>110	>120
	Receptacles		2 GFI duplex, one on each sidewall			
	Service valves		Up to 4 may be installed through access ports (2 on each side) Up to 6 more may be factory-installed on rear wall (optional)			
<b>10-Inch Work Aperture</b>	Sound pressure level	dB (A)	<65			
	Energy consumption, operating set point <sup>2</sup>	W	210	275	350	435
	Energy consumption, Night Set-Back mode	W	40	40	70	70
	Heat output, operating mode (non-vented)	BTU/hr	717	938	1194	1484
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	136	136	239	239
	Exhaust / inflow air volume, recirculated	cfm	255	342	428	514
	Exhaust volume, thimble ducted	cfm	332	444	556	668

<sup>1</sup> Export packaging dimensions and weight upon request

<sup>2</sup> Clean filters, fans operating at nominal set point, interior lighting activated



# ordering information

The Herasafe KS BSC standard configuration includes cabinet, UV light, armrests, 316 stainless steel indented work surface, and battery back-up for motorized sash. Stands must be ordered separately.

Size/Width		3-foot	4-foot	5-foot	6-foot
<b>Standard Configuration</b>	Includes cabinet, cross-beam UV, one-piece indented worktray and battery back-up for front window. Stands are sold separately.	51027116	51027117	51027118	51027120



Cross-beam UV light



Stainless steel armrests



One-piece indented worktray

Size/Width		3-foot	4-foot	5-foot	6-foot
<b>Stands:</b> <i>(Must be ordered separately)</i>	Manual adjustable floor stand, 30 to 38 inches in 2-inch increments	50116443	50109312	50116444	50109313
	Electric adjustable floor stand, with cable remote control, 30 to 38 inches working height, 120V	50119096	50119126	50119127	50119128
	Castor set (may not be used with electric stands; fixed 30-inch work surface height)	50109977			



Electric adjustable stand



Manual stand



Castors

To build your own Herasafe KS BSC, start with these cabinet-only part numbers (with armrests):

Size/Width		3-foot	4-foot	5-foot	6-foot
<b>Special Order</b>	Cabinet-Only Part Number	51022485	51022486	51022487	51022488

# Thermo Scientific Herasafe KS Biological Safety Cabinet

## field-installed accessories

	Cat. No.	Description
<b>Exhaust Transitions</b>	50078796	Thimble duct exhaust transition for Herasafe KS 9 (3 ft) cabinet
	50073770	Thimble duct exhaust transition for Herasafe KS 12 (4 ft) or KS 15 (5 ft) cabinet
	50078795	Thimble duct exhaust transition for Herasafe KS 18 (6 ft) cabinet
	191780	Alnor exhaust alarm for thimble duct, 120v (independent from BSC, plugs into wall outlet)
<b>User Convenience Equipment</b>	50078900	Utility basket for use with hanging bars 50073943 or 1911312
	50078901	2 Utility baskets including additional hanging bar
	50073943	Hanging bar with set of 6 hooks for KS 9 (3 ft) cabinet only
	1911312	Hanging bar with set of 12 hooks for KS 12, 15, 18 (4, 5, 6 ft.) cabinets
	1911313	Replacement stainless steel armrest (set of 2)
	191127	Adjustable foot rest. For ergonomic foot and posture support
	1910179	Shelf to hang from base stand
	1910083	Interior tubing retainer clips (6 per kit)
	ERG-CR-2330	Thermo Scientific Ergolign Saddle Stool, 23 inch to 30.5 inch adjustable seat height
<b>Sidewall Utilities</b>	1911308	Combustible gas tap for side wall
	1911309	Noncombustible gas tap for side wall
	1911310	Vacuum tap for side wall
	1911311	Water tap for side wall
	50059017	Hose connection to fit through side wall access port, fits 0.5 inch (10-13 mm) inner diameter tubing
	50076408	Hose connection to fit through side wall access port, fits 0.1 inch (3 mm) inner diameter tubing



50073700  
thimble duct exhaust transition



50078900  
utility basket



1911311  
water tap for side wall



Optional rear wall  
3 service tap connections

## factory-installed accessories

	Cat. No.	Description
<b>Customized Models</b>		The Herasafe KS worksurfaces are fully customizable to accommodate any unique work task. Review ordering information with your Thermo Fisher Scientific sales representative. Lead times will be extended for custom orders.
<b>Worksurfaces</b>	51900651	KS 9 (3 ft) without worksurface — required for installation of non-standard worksurface options
	51900652	KS 12 (4 ft) without worksurface — required for installation of non-standard worksurface options
	51900653	KS 15 (5 ft) without worksurface — required for installation of non-standard worksurface options
	51900654	KS 18 (6 ft) without worksurface — required for installation of non-standard worksurface options
	50073666	Flat worksurface module 150 mm wide (316SS)
	50072685	Flat worksurface module 300 mm wide (316SS)
	50073667	Flat worksurface module 600 mm wide (316SS)
	50076398	Flat worksurface module 900 mm wide (316SS)
	50082970	Flat worksurface module 1200 mm wide (316SS)
	50122210	Flat worksurface module 1500 mm wide (316SS)
	50085281	Flat worksurface module 1800 mm wide (316SS)
	50073757	Vibration damped worksurface 300 mm width. Constructed in stainless steel with granite plate. For weighing and microscopy equipment (316SS)
	50076162	Vibration damped worksurface 600 mm width. Constructed in stainless steel with granite plate. For weighing and microscopy equipment (316SS)
	50078875	Minisink tabletop insert 300 mm. Holds ~3 liters
	50073668	Indented worksurface module 600 mm wide (316SS)
	50073669	Indented worksurface module 900 mm wide (316SS)
	50073670	Indented worksurface module 1200 mm wide (316SS)
	50073671	Indented worksurface module 1500 mm wide (316SS)
	50078750	Indented worksurface module 1800 mm wide (316SS)
<b>Universal Piping</b>	51900981	Rear wall 3 service tap connections, NPT 3/8, left side of cabinet
	51900982	Rear wall 3 service tap connections, NPT 3/8, right side of cabinet
	51900984	Rear wall 6 service tap connections, NPT 3/8, 3 on left and 3 on right side of cabinet. (May not be used with KS 9 (3 ft) cabinet.)
	1911304	Combustible gas tap for rear wall (must be ordered with rear wall service tap connections - above)
	1911305	Noncombustible gas tap for rear wall (must be ordered with rear wall service tap connections - above)
	1911306	Vacuum tap for rear wall (must be ordered with rear wall service tap connections - above)
1911307	Water tap for rear wall (must be ordered with rear wall service tap connections - above)	
<b>UV Options</b>	51900642	Cross beam UV-C light option in side walls (comes with protective cover)



50078875  
optional minisink tabletop insert



50076162 optional 600 mm vibration  
dampened granite block



Optional 316 stainless steel  
worksurface configurations

# Thermo Scientific 1300 Series A2 Biological Safety Cabinet

## specifications

Size / Width			3-foot	4-foot	5-foot	6-foot
<b>Dimensions</b>	Exterior dimensions W x H x D	in	39.4 x 61.8 x 31.5	51.2 x 61.8 x 31.5	63 x 61.8 x 31.5	74.8 x 61.8 x 31.5
	Interior dimensions W x H x D	in	35.4 x 30.7 x 24.8	47.2 x 30.7 x 24.8	59.1 x 30.7 x 24.8	70.9 x 30.7 x 24.8
	Shipping dimensions <sup>1</sup> W x H x D	in	43.7 x 68.1 x 36.4	55.5 x 68.1 x 36.4	67.3 x 68.1 x 36.4	79.1 x 68.1 x 36.4
	Working height of front window	in	Choice of 8 or 10			
	Maximum opening height of front window	in	21			
<b>Weight</b>	Net weight	lbs	375	441	507	617
	Shipping weight <sup>1</sup>	lbs	419	496	573	694
<b>Electrical Data</b>	Voltage	V	120			
	Frequency	Hz	60			
<b>Filter Specification</b>			H14 HEPA EN 1822, 99.995% MPPS			
<b>Certification</b>			NSF/ANSI 49, UL, CE			
<b>Ergonomics</b>	Lighting power	fc	>100	>120	>120	>120
	SmartPort		Two 3" plugged cable ports, one on each side wall			
	Receptacles		2 GFI duplex			
	Service valves		Up to 6 total (3 on each side wall)			
<b>8-Inch Work Aperture</b>	Sound pressure level	dB (A)	<62	<62	<64	<64
	Energy consumption, operating set point <sup>2</sup>	W	150	180	280	360
	Energy consumption, Night Set-Back mode	W	70	70	120	120
	Heat output, operating mode (non-vented)	BTU/hr	512	614	955	1228
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	239	239	410	410
	Exhaust/inflow air volume	cfm	207	276	344	413
	Exhaust volume, thimble ducted	cfm	269	358	448	537
<b>10-Inch Work Aperture</b>	Sound pressure level	dB (A)	<63	<63	<65	<65
	Energy consumption, operating set point <sup>2</sup>	W	170	200	310	400
	Energy consumption, Night Set-Back mode	W	70	70	120	120
	Heat output, operating mode (non-vented)	BTU/hr	580	683	1058	1366
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	239	239	410	410
	Exhaust/inflow air volume	cfm	258	344	431	517
	Exhaust volume, thimble ducted	cfm	336	448	560	672

<sup>1</sup> Export packaging dimensions and weight upon request

<sup>2</sup> Clean filters, fans operating at nominal set point, interior lighting activated



## ordering information

Our 1300 Series A2 BSC packages feature everything required to get started, including cabinet, manual adjustable height stand, factory-installed UV light and one set of armrests.

	Size	8-inch work aperture	10-inch work aperture
<b>SmartCoat Interior</b>	3-foot	1340	1321
	4-foot	1395	1385
	5-foot	1351	1369
	6-foot	1397	1387
<b>One-Piece Stainless Steel Interior</b>	3-foot	1342	1323
	4-foot	1335	1375
	5-foot	1353	1371
	6-foot	1337	1377



Manual adjustable stand



UV light



Stainless steel armrests

To build your own 1300 Series A2 BSC, start with these cabinet-only part numbers:

	Size	8-inch work aperture	10-inch work aperture
<b>SmartCoat Interior</b>	3-foot	1339	1320
	4-foot	1325	1355
	5-foot	1350	1368
	6-foot	1327	1357
<b>One-Piece Stainless Steel Interior</b>	3-foot	1341	1322
	4-foot	1305	1345
	5-foot	1352	1370
	6-foot	1307	1347

# Thermo Scientific 1300 Series A2 Biological Safety Cabinet

## field-installed accessories

	Cat. No.	Description
<b>Stands</b>	1911400	Manual adjustable stand for 3 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911302	Manual adjustable stand for 4 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911401	Manual adjustable stand for 5 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911303	Manual adjustable stand for 6 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	50119096	Electric adjustable stand for 3 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119126	Electric adjustable stand for 4 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119127	Electric adjustable stand for 5 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119128	Electric adjustable stand for 6 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	1911409	Rolling castor stand for 3 foot cabinet, provides fixed height of 30 inches
	1911322	Rolling castor stand for 4 foot cabinet, provides fixed height of 30 inches
	1911410	Rolling castor stand for 5 foot cabinet, provides fixed height of 30 inches
	1911346	Rolling castor stand for 6 foot cabinet, provides fixed height of 30 inches
	<b>Exhaust Transitions</b>	1911316
1911317		Thimble duct exhaust transition for 6 foot cabinet
1910185		Alnor exhaust alarm for 1300 Series A2 biological safety cabinets, qualified install only
<b>User Convenience Equipment</b>	1910083	Interior tubing retainer clips (6 per kit)
	1910179	Hanging shelf for base stand
	1911313	Stainless steel armrest (set of two)
	191127	Adjustable footrest
	1910112	Replacement SmartPort grommets (4 per kit)
	ERG-CR-2330	Thermo Scientific Ergolign Saddle Stool, 23 inch to 30.5 inch adjustable seat height
	1911413	IV bag holder kit with 6 hooks (for 3 foot cabinet only)
1911312	IV bag holder kit with 12 hooks (for 4 foot, 5 foot, and 6 foot cabinet)	
<b>Universal Piping Utilities</b>	249098	Stainless steel ball valve (replaces standard non-corrosive PVC drain valve)
<b>Sidewall Utilities</b>	1911308	Combustible gas valve
	1911309	Non-combustible gas valve
	1911310	Vacuum valve
	1911311	Water valve



1910179 hanging shelf for base stand



ERG-CR-2330 Ergolign Saddle Stool



191127 optional adjustable footrest

## factory-installed accessories

	<b>Cat No :</b>	<b>Left or right connection</b>	<b>Top or bottom connection</b>
<b>Universal Piping for 3 Foot Models</b>	1910286	Right	Both
	1911415		Top
	1910287	Left	Both
	1911416		Top
<b>Universal Piping for 4 Foot Models</b>	1911321	Right	Both
	1911417		Top
	1911369	Left	Both
	1911418		Top
<b>Universal Piping for 5 Foot Models</b>	1910288	Right	Both
	1911419		Top
	1910289	Left	Both
	1911420		Top
<b>Universal Piping for 6 Foot Models</b>	1911367	Right	Both
	1911421		Top
	1911370	Left	Both
	1911422		Top
<b>Universal Piping for All Models</b>	1911423	Right	Bottom
	1911424	Left	
	<b>Cat No :</b>	<b>Description</b>	
<b>Universal Piping Utilities</b>	1911304	Combustible gas valve for universal piping option	
	1911305	Non-combustible gas valve for universal piping option	
	1911306	Vacuum valve for universal piping option	
	1911307	Water valve for universal piping option	
<b>UV Light Option</b>	1911319	UV light option for 3 foot or 4 foot cabinet	
	1911368	UV light option for 5 foot or 6 foot cabinet	

# Fumigation Sealing Kits

## field-installed accessories

Available for Herasafe KS and 1300 Series biological safety cabinets. In addition to surface wipe downs within the BSC work area, it may be necessary to fumigate the potentially contaminated interior plenums of the unit. Fumigation to access biologically contaminated BSC plenums is most often done using formalin, paraformaldehyde, or hydrogen peroxide vapor. With these and most other methods, it is necessary to first seal off the unit, then introduce the fumigating agent, allow appropriate time for the agent to be effective, then neutralize and/or release the fumigating agent. The fumigating agent and method selected is usually specified by local regulations or institution guidelines based on the type of contamination. The kits and accessories below facilitate the fumigation process.


	Cat. No.	Description
<b>Complete Kits</b>	50132757	H <sub>2</sub> O <sub>2</sub> adjustable sealing kit (ASK) includes adjustable top box and front window cover. Fits all sizes of BSC. (includes outlet and inlet connections for STERIS™)
<b>Inlets Only</b>	50135768	H <sub>2</sub> O <sub>2</sub> front window sealing plate. This seals the front aperture, and has a STERIS connector
	50137105	STERIS H <sub>2</sub> O <sub>2</sub> Camlock inlet connector. Fits into one access port left or right side of the BSC
	50137107	Bioquell™ H <sub>2</sub> O <sub>2</sub> Camlock inlet connector (includes sample connector). Uses two access ports on side of BSC
<b>Outlets Only</b>	50137106	STERIS or Bioquell Camlock outlet connector (for connection to top box or duct connectors)
	50135767	H <sub>2</sub> O <sub>2</sub> adjustable top box for Thermo Scientific BSC with outlet connector (STERIS or Bioquell)
<b>Others</b>	50114036	Fumigation window sealing kit for 1300 series BSC only
	50109998	Disinfection tent for formaldehyde disinfection

### H<sub>2</sub>O<sub>2</sub> Fumigation Accessories

When fumigating a biological safety cabinet using H<sub>2</sub>O<sub>2</sub>, there are many accessories available. When choosing accessories there are basic questions to consider.

- Is the biological safety cabinet recirculating or ducted? See the chart below for guidance.
- Is the inlet of the H<sub>2</sub>O<sub>2</sub> generator connected via a removable front sealing gasket (STERIS only) or through a fixed side connector (STERIS or Bioquell)? See chart below for guidance.

	Recirculating		Ducted	
	Side Connection	Front Connection	Side Connection	Front Connection
<b>STERIS</b>	50137105 (one per cabinet)	50132757 (one required)	50137105 (one per cabinet)	50135768 (one required)
	50135767 (one required)		50137106 (one per cabinet)	50137106 (one per cabinet)
<b>Bioquell</b>	50137107 (one per cabinet)	Not available	50137107 (one per cabinet)	Not available
	50135767 (one required)		50137106 (one per cabinet)	



50114036 fumigation window sealing kit



H<sub>2</sub>O<sub>2</sub> side wall connections



H<sub>2</sub>O<sub>2</sub> adjustable sealing kit (ASK)

# Thermo Scientific Heraguard ECO Clean Benches

## specifications

Model		ECO 0.9	ECO 1.2	ECO 1.5	ECO 1.8	ECO 1.8/95
<b>Ordering Information</b>						
White melamine worksurface		51029696	51029697	51029698	51029699	51029700
Stainless steel worksurface		51029691	51029692	51029693	51029694	51029695
<b>Dimensions</b>						
Exterior dimensions W x H x D	in	39.4 x 46.1 x 31.9	51.2 x 46.1 x 31.9	63.0 x 46.1 x 31.9	74.8 x 46.1 x 31.9	74.8 x 58.1 x 39.8
Interior dimensions W x H x D	in	36.2 x 25.4 x 23.0	48.0 x 25.4 x 23.0	59.8 x 25.4 x 23.0	71.6 x 25.4 x 23.0	71.6 x 37.4 x 30.9
Shipping dimensions <sup>1</sup> W x H x D	in	44.5 x 52.0 x 38.6	56.3 x 52.0 x 38.6	68.1 x 52.0 x 38.6	79.9 x 52.0 x 38.6	79.9 x 63.8 x 46.5
<b>Weight</b>						
Net weight	lbs	245	275	310	355	375
Shipping weight <sup>1</sup>	lbs	300	330	400	440	485
Max. load capacity of worksurface	lbs	66	66	66	66	66
<b>Electrical Data</b>						
Voltage	V					120
Frequency	Hz					60
Protection class						I
Overvoltage category						II
IP protection system						IP 20
<b>Compliance</b>						
Agency						CSA
Standards						CAN/CSA-C22.2 No. 61010-1; UL Std. No. 61010-1; IEC 61010-1
	prefilter					Type (EN779) G 3 dust filter, filtering efficiency: coarse - medium
<b>Performance Specifications</b>						
Power consumption	W	279	279	574	574	574
Heat output	BTU/hr	951	951	1930	1930	1930
Air velocity	ft/min	89 +- 20%				79 +- 20%
Reduced airflow velocity	ft/min	45 +- 20%				
Lighting intensity	fc	93				
Noise pressure level	dB(A)	< 55		< 60		
Minimum clearance from room ceiling	in	5.9				
Side windows		Tempered safety glass				
Receptacles <sup>2</sup>		2				
Service valves		Up to 2 may be installed through access ports (1 left, 1 right)				

<sup>1</sup> Export packaging dimensions and weight upon request.

<sup>2</sup> The receptacles have a load capacity of up to 5 A and are protected with T 5 A fuses. When all receptacles are in use simultaneously they must not exceed the maximum total load capacity of 5 A.

# Thermo Scientific Heraguard ECO Clean Benches

## ordering information

<b>Heraguard ECO, Stainless Steel (SS)</b>	<b>51029691<sup>1</sup></b>	Heraguard ECO 0.9, SS worksurface (316L, 1.4404), 120V/60Hz, 3ft
	<b>51029692<sup>1</sup></b>	Heraguard ECO 1.2, SS worksurface (316L, 1.4404), 120V/60Hz, 4ft
	<b>51029693<sup>1</sup></b>	Heraguard ECO 1.5, SS worksurface (316L, 1.4404), 120V/60Hz, 5ft
	<b>51029694<sup>1</sup></b>	Heraguard ECO 1.8, SS worksurface (316L, 1.4404), 120V/60Hz, 6ft
	51029695	Heraguard ECO 1.8/95, SS worksurface (316L, 1.4404), 120V/60Hz, 6ft, 37.4" internal height
<b>Heraguard ECO, SS with built-in UV light</b>	<b>51029935<sup>1</sup></b>	Heraguard ECO 0.9, SS worksurface (316L, 1.4404), 120V/60Hz, 3ft, with UV option
	<b>51029937<sup>1</sup></b>	Heraguard ECO 1.2, SS worksurface (316L, 1.4404), 120V/60Hz, 4ft, with UV option
	<b>51029939<sup>1</sup></b>	Heraguard ECO 1.5, SS worksurface (316L, 1.4404), 120V/60Hz, 5ft, with UV option
	<b>51029940<sup>1</sup></b>	Heraguard ECO 1.8, SS worksurface (316L, 1.4404), 120V/60Hz, 6ft, with UV option
<b>Heraguard ECO</b>	51029696	Heraguard ECO 0.9, 120V/60Hz, 3ft
	51029697	Heraguard ECO 1.2, 120V/60Hz, 4ft
	51029698	Heraguard ECO 1.5, 120V/60Hz, 5ft
	51029699	Heraguard ECO 1.8, 120V/60Hz, 6ft
	51029700	Heraguard ECO 1.8/95, 120V/60Hz, 6ft, 37.4" internal height

<sup>1</sup> Model is a stocked configuration.



Left hand glass with metal sheet and two service valves (optional)



UV-resistant night cover



Wind screen protection kit

## field-installed accessories

<b>UV Resistant Night Covers</b>	50139417	UV-resistant night cover for Heraguard ECO 0.9
	50139418	UV-resistant night cover for Heraguard ECO 1.2
	50139419	UV-resistant night cover for Heraguard ECO 1.5
	50139420	UV-resistant night cover for Heraguard ECO 1.8
	50139421	UV-resistant night cover for Heraguard ECO 1.8/95
<b>Stands</b>	50116443	Adjustable floor stand for 3ft. Provides variable work height of 30 to 38 inches in 2-inch increments
	50109312	Adjustable floor stand for 4ft. Provides variable work height of 30 to 38 inches in 2-inch increments
	50116444	Adjustable floor stand for 5ft. Provides variable work height of 30 to 38 inches in 2-inch increments
	50109313	Adjustable floor stand for 6ft. Provides variable work height of 30 to 38 inches in 2-inch increments
	191127	Adjustable foot rest. For ergonomic foot and posture support
<b>Media Valves</b>	1911308	Combustible gas tap for side wall
	1911309	Noncombustible gas tap for side wall
	1911310	Vacuum tap for side wall
	1911311	Water tap for side wall
	50059017	Hose connection to fit through side wall access port, fits 0.5 inch (10-13 mm) inner diameter tubing
	50076408	Hose connection to fit through side wall access port, fits 0.1 inch (3 mm) inner diameter tubing
<b>Hanging Bars</b>	50073943	Hanging bar with 6 hooks for 3ft clean bench
	1911312	Hanging bar with 6 hooks. For use with 4ft, 5ft, 6ft clean bench
<b>Windscreen Protection Kits</b>	50052208	Wind screen protection kit, complete with hinges to fit the control unit for Heraguard ECO 0.9
	50052460	Wind screen protection kit, complete with hinges to fit the control unit for Heraguard ECO 1.2
	50052461	Wind screen protection kit, complete with hinges to fit the control unit for Heraguard ECO 1.5
	50052462	Wind screen protection kit, complete with hinges to fit the control unit for Heraguard ECO 1.8

## factory-installed accessories

<b>Side Window Modifications</b>	51901071	Left hand glass with metal sheet and 2 media valve holes and SmartPort for all models
	51901072	Right hand glass with metal sheet and 2 media valve holes and SmartPort for all models
<b>UV Light Option</b>	51901073	UV option, supplied with UV-resistant night cover

# what you need to meet an application challenge



## Redefining pipetting productivity

Trust the Thermo Scientific S1 Pipet Filler to deliver everything you need to stay productive in the lab: faster performance, advanced speed controls and longer battery life. The portable lightweight S1 Pipet Filler is designed for use with either glass or plastic serological pipets. It offers effective streamlined pipetting performance with maximum pipetting comfort.

[thermoscientific.com/handleanything](http://thermoscientific.com/handleanything)

## Discovery thrives in a culture of confidence

More scientists worldwide trust their valuable cultures to Thermo Scientific CO<sub>2</sub> incubators than any other brand. They depend on proven reliability, exceptional contamination prevention and optimal growing conditions. Delivered with innovative features like HEPA air filtration that surrounds cells with clean room-like air quality and a choice of 100% pure copper or polished stainless steel interior surfaces. Plus a high-temperature decontamination function that eliminates the need for autoclaving and handling components. The inside story is simple: our CO<sub>2</sub> incubators let you culture with complete confidence. Day after day. Year after year.

[thermoscientific.com/co2incubator](http://thermoscientific.com/co2incubator)



To learn more visit [thermoscientific.com/biosafetycabinets](http://thermoscientific.com/biosafetycabinets)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. STERIS is a registered trademark of STERIS. Bioquell is a registered trademark of Bioquell. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Australia** +61 39757 4300  
**Austria** +43 1 801 40 0  
**Belgium** +32 53 73 42 41  
**China** +86 21 6865 4588 or  
+86 10 8419 3588  
**France** +33 2 2803 2180  
**Germany national toll free** 0800 1 536 376  
**Germany international** +49 6184 90 6000

**India toll free** 1800 22 8374  
**India** +91 22 6716 2200  
**Italy** +39 02 95059 552  
**Japan** +81 45 453 9220  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Nordic/Baltic/CIS countries**  
+358 9 329 10200

**Russia** +7 812 703 42 15  
**Spain/Portugal** +34 93 223 09 18  
**Switzerland** +41 44 454 12 22  
**UK/Ireland** +44 870 609 9203  
**USA/Canada** +1 866 984 3766

**Other Asian countries** +852 2885 4613  
**Countries not listed** +49 6184 90 6000

**Thermo**  
S C I E N T I F I C

Part of Thermo Fisher Scientific



Thermo Scientific  
Biological Safety Cabinets

**smart planning**  
for your lab's protection

**Thermo**  
SCIENTIFIC

# Thermo Scientific Biological Safety Cabinets

## protection that never takes a day off

**With our biological safety cabinets (BSCs), the certified performance and protection you get on Day 1 stays with you everyday. Not true with ordinary cabinets.**

### 5. Valdymas

The difference is our design: SmartFlow™ technology features dual-DC motors to automatically balance the cabinet inflow and downflow air velocities in real time – even as the filters load with contaminants. That means exceptional user and sample protection you never have to think about. Plus, our Digital Airflow Verification™ (DAVe) alarm signals any out-of-spec conditions for added assurance. With our proven reliability, ergonomics and energy efficiency, the optimal choice is also the one you can trust completely. And not just on Day 1.



Thermo Scientific Maxisafe 200  
biological safety cabinet



Thermo Scientific 1300 Series B2  
biological safety cabinet

Thermo Scientific 1300 A2  
biological safety cabinet

## Smart planning for your lab's protection

Planning a laboratory or renovating an existing one can be complex and challenging. Critical decisions are required at every stage of the project:

**What are your research requirements – today and in the future?**

**What is your lab's workflow?**

**What are your goals for lab sustainability?**

**How does this impact your budget requirements?**

Whether you are building a new lab or renovating your current one, our team can assist you through the lab planning process with technical and application guidance.

- Application support for class and type of biological safety cabinets
- Specifications for mechanical, electrical and plumbing requirements
- Revit® and CAD drawings available for pre-design planning
- Sustainability and energy savings calculations
- Budgeting and planning support for laboratory equipment



Thermo Scientific Herasafe KS biological safety cabinet

## Thermo Scientific biological safety cabinets

Safe and reliable cabinets for a wide range of laboratory applications.

Thermo Scientific BSC Product Family		North America	Latin America	EMEA		Asia			ANZ (Australia/ New Zealand)
				Europe	Middle East/ Africa	China	Japan	ROA (Rest of Asia)	
Triple Filter	Herasafe KSP	•	•	•	•	•	•	•	•
	Maxisafe 2020			•	•	•		•	•
Class II	Herasafe KS	•	•	•	•	•	•	•	•
	Safe 2020			•	•	•		•	•
	MSC-Advantage			•	•	•		•	•
	1300 Series A2	•	•		•	•	•	•	
	1300 Series B2	•	•		•	•	•	•	
Clean Benches	1800 Series Clean Bench	•	•				•	•	
	Heraguard HPH Clean Bench			•	•	•		•	•

## Outstanding safety and compliance

Not all BSCs are manufactured to the same safety standards. Our design philosophy is grounded in compliance with the recognized worldwide safety standards. We do not compromise on safety – biological safety cabinets available in the following regions are subjected to the following certifications.

### Certification and compliance (including all CE and TÜV certificates)

#### USA / Canada

- IEC 1010-1 / UL 3101-1
- NSF, Standard 49

#### Europe / Middle East / Africa

- IEC 1010-1 / EN61010-1
- EN 12469 / 2000
- DIN 12980
- Low Voltage Directive 73/23 EWG
- EMC Directive 89/336 EWG

#### Australia / China / Asia Pacific

- IEC 1010-1
- YY 0569
- AS 2252.2

EN12469



[thermoscientific.com/bscnewlab](http://thermoscientific.com/bscnewlab)

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Revit is a registered trademark of Autodesk. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Australia** +61 39757 4300  
**Austria** +43 1 801 40 0  
**Belgium** +32 53 73 42 41  
**China** +86 21 6865 4588 or  
+86 10 8419 3588  
**France** +33 2 2803 2180  
**Germany national toll free** 0800 1 536 376  
**Germany international** +49 6184 90 6000

**India toll free** 1800 22 8374  
**India** +91 22 6716 2200  
**Italy** +32 02 95 05 95 52  
**Japan** +81 45 453 9220  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Nordic/Baltic/CIS countries**  
+358 9 329 10200

**Russia** +7 812 703 42 15  
**Spain/Portugal** +34 93 223 09 18  
**Switzerland** +41 44 454 12 22  
**UK/Ireland** +44 870 609 9203  
**USA/Canada** +1 866 984 3766

**Other Asian countries** +852 2885 4613  
**Countries not listed** +49 6184 90 6000

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific