## Serascan Diana 2 Serascan Diana 2P

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## Serascan Diana 2 Serascan Diana 2P

Read carefully before using the product. For use in "in vitro" diagnosis only.

### INTENDED USE

Screening of unexpected antibodies, in gel technique.

### INTRODUCTION

The screening of unexpected antibodies has the aim of detecting the clinically significant antibodies present in the patient's sample. The autcontrol will indicate, in a positive screening of unexpected antibodies, whether it is due to the presence of an autoantibody, an alloantibody or both.

The Serascan Diana 2 and Serascan Diana 2P reagent red blood cells present the most significant antigenic determinants for most of the blood group systems.

### PRINCIPLE

An antibody reacts specifically with the antigen which stimulated its production<sup>1</sup>. According to this, an antibody can be identified depending on its pattern of reactivity when confronted with a panel of the principle of the test is based on the gel technique described by Y. Lapiere<sup>2</sup> for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

### COMPOSITION

Each vial of Serascan Diana 2 (I, II) contains 10 ml of human red blood cells of O blood group, in a 0.8% suspension in a buffered solution with preservatives.

The Serascan Diana 2 reagent is manufactured from a single donor for each vial. It is provided in stoppered glass vials with a white-coloured bulb.

Each vial of Serascan Diana 2P (II) contains 10 ml of papainised human red blood cells of O blood group, in a 0.8% suspension in a buffered solution with preservatives.

The Serascan Diana 2P reagent is manufactured from a single donor for each vial. It is provided in stoppered glass vials with a red-coloured bulb.

The antigenic configuration of the Serascan Diana 2 and Serascan Diana 2P reagent red blood cells is given in the enclosed table (specific for each product lot).

Reagent ready to use. Do not use the reagent if cloudiness and/or large-scale haemolysis are detected.

### STABILITY

Serascan Diana 2 and Serascan Diana 2P are stable up to the expiry date stated on the label, if stored at 2-8 °C in the position indicated on the outer packaging. Do not freeze.

Once the vial has been opened, handle correctly, avoiding contamination of the content, and once the product has been used, keep it at the indicated storage temperature.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Diagnostic Grifols, S.A. gel cards for the screening of unexpected antibodies.

### SAMPLES

Recently extracted blood samples, collected with the usual blood bank anticoagulants (obtaining plasma) or without anticoagulants (obtaining serum).

Do not use haemolysed, cloudy or contaminated samples, or with clot presence.

The procedure for extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives<sup>3</sup>, and following the instructions of the manufacturer of the material used for collecting the sample.

- Screening of unexpected antibodies: use serum or plasma. If necessary, samples stored at 2-8 °C may be used up to 48 hours after their extraction, or frozen samples (-20 °C to -80 °C) up to the time of performing the test.

### TEST METHOD

Serascan Diana 2 and Serascan Diana 2P can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Allow the samples and reagents to reach room temperature (18 - 25 °C).

Inspect the staple of the reagent before use (see section LIMITATIONS "In relation to the product").

Before use, homogenise the red blood cells of the Serascan Diana 2 and Serascan Diana 2P vials, by gentle inversion.

### Manual method

Serascan Diana 2 and Serascan Diana 2P reagent red blood cells may be used with any of the usual methods of screening of unexpected antibodies with Diagnostic Grifols, S.A. gel cards.


Follow the instructions for use of the Diagnostic Grifols, S.A. gel cards used.

## English

An autcontrol should always be performed in parallel with each test, confronting the patient's red blood cells with his/her own serum or plasma.

### RESULTS

#### Result reading

<b>Negative:</b>	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
	+/-	Scas small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations, in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
	DP	Double Population (double band of red blood cells, at the bottom and in the upper part of the column).
	H	Haemolysis (pinkish supernatant and/or gel colour).

Stability of the results: It is recommended an immediate reading of the results after centrifuging the cards.

If necessary, a delayed reading can be made up until 24 hours after processing the cards, if kept in vertical position, refrigerated (2-8 °C) and sealed with parafilm or similar material, to avoid evaporation of the supernatant.

### Interpretation of the results

To know the antigenic configuration of the reagent red blood cells, consult the Serascan Diana 2/ Serascan Diana 2P table enclosed, which is specific for each product lot. Compare the pattern of reactivity obtained with the antigenic profile of the reagent red blood cells used.

Reagent red blood cells	Sample red blood cells (autocontrol)	Interpretation
0	0	Absence of allo- and autoantibodies
+	0	Presence of alloantibodies
0	+	Presence of autoantibodies
+	+	Presence of autoantibodies or autoantibodies + alloantibodies

Consultar, si es necesario, las instrucciones for use of the Diagnostic Grifols, S.A. gel cards used.

The results themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.

### Notes:

- If the autcontrol is positive, investigate the presence of cold antibodies, "rouleaux" and perform a Direct Coombs test.

### QUALITY CONTROL

- It is recommended to include known positive and negative controls in each series of tests.
- In the screening of unexpected antibodies, it is recommended to include weak positive controls.

- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used should be made.

### PERFORMANCE CHARACTERISTICS

Although there is no described procedure or technique capable of detecting all the possible unexpected antibodies present in a sample population with absolute certainty, performance evaluation studies carried out in different hospitals with population samples, back up that Serascan Diana 2 and Serascan Diana 2P own performance characteristics conforming to the intended use of the product. These studies<sup>4</sup> have included tests of screening and identification of unexpected antibodies, being the results obtained similar or better than those obtained with other established products of equivalent intended use.

### Manual method

Serascan Diana 2 and Serascan Diana 2P reagent red blood cells may be used with any of the usual methods of screening of unexpected antibodies with Diagnostic Grifols, S.A. gel cards.

Follow the instructions for use of the Diagnostic Grifols, S.A. gel cards used.

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**LIMITATIONS****In relation to the sample:**

Do not use haemolysed, cloudy or contaminated blood samples, or with clot presence. Sample: Serum/Plasma

- If plasma is used, complement-dependant haemolytic reactions may not be detected.
- If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such. It is recommended to recolt the serum during 10 minutes at 37 °C, centrifuge and repeat the test!
- The formation of "rouleaux", caused by an excess of proteins in the serum, the presence of abnormal proteins, drugs, plasma expanders, etc., may cause false positive reactions!
- If there are antibodies against high-incidence antigens or multiple antibodies in the sample, all the reagent red blood cells may be agglutinated.

**In relation to the product:**

Before use, inspect the state of the reagent:

- Do not use the reagent if cloudiness and/or large-scale hemolysis are observed. These may be caused by microbiological contamination or improper handling/storage.
- Close the vials correctly after use, to prevent contamination. Ensure that the stoppers of the reagent red blood cell vials have not been switched.
- If, due to improper transport or storage, the vial is broken and/or there is a loss of content, do not use the product.
- Improper storage or handling may cause a loss of reactivity of the red blood cells.
- Do not use the product beyond the expiry date.
- The Serascan Diana 2 and Serascan Diana 2P reagent red blood cells do not detect the presence of anti-A or anti-B antibodies.
- Low-incidence antigens may not be represented in the Serascan Diana 2 and Serascan Diana 2P reagent red blood cells, so negative reactions with them do not always indicate absence of an antibody in the sample under study.
- The Serascan Diana 2P (papainised) reagent red blood cells present a reduced or absent reactivity of some antigens, which appear shaded on the table enclosed with the product.

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- Serascan Diana 2 and Serascan Diana 2P are of human origin and are manufactured using a non-reactive material for the HBS antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin do not transmit hepatitis or AIDS. Human blood products and samples should be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

**BIBLIOGRAPHY**

- Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard - 4th edition.
- NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline - 2nd edition.
- Diana System-Reference Studies, Diagnostic Grifols, 2001.

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## LITERATŪRA

1. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
2. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
3. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
4. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard - 4th edition.
5. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline - 2nd edition.
6. Diana System-Reference Studies, Diagnostic Grifols, 2001.

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Read carefully before using the product. For use in “in vitro” diagnosis only

## INTENDED USE

Indirect Coombs and Direct Coombs tests, in gel technique.  
The Indirect Coombs tests include: screening and identification of unexpected antibodies, crossmatch tests, autocontrol and red blood cell typing.

## INTRODUCTION

Coombs introduced the antiglobulin test in clinical medicine in 1945<sup>1</sup>, and this is why it is known as the Coombs test. This test is based on the use of a polyspecific anti-human globulin (AHG) that allows the detection of red blood cells coated with immunoglobulins or complement fractions. The Indirect Coombs test allows the detection of red blood cell antibodies present in the patient's serum or plasma, by “in vitro” sensitisation of red blood cells.

The Direct Coombs test allows the detection of red blood cells sensitised “in vivo” by immunoglobulins or complement fractions.

The screening of unexpected antibodies has the aim of detecting the clinically significant antibodies present in the patient's sample. In a positive screening of unexpected antibodies, the autocontrol will indicate whether it is due to the presence of an autoantibody, an alloantibody or both.

### Crossmatch tests:

- major crossmatch test: the donor's red blood cells confronted with the recipient's serum or plasma will show the presence or absence of unexpected antibodies in the recipient's blood that are specific against the antigens of the donor's red blood cells.
- minor crossmatch test: the recipient's red blood cells confronted with the donor's plasma or serum will show the presence or absence of unexpected antibodies in the donor's blood that are specific against the antigens of the recipient's red blood cells.

## PRINCIPLE

The principle of the test is based on the gel technique described by Y. Lapiere<sup>2</sup> for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

The DG Gel card is a plastic support consisting of 8 microtubes. Each microtube consists of a column and a dispensing/incubation chamber.

Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains anti-human globulin.

The microtubes that contain anti-human globulin act by agglutinating the red blood cells sensitised, “in vivo” or “in vitro”, by IgG antibodies or complement fractions.

During centrifugation and depending on their size, the red blood cell agglutinations are trapped in the surface or along the gel column. The non-agglutinated red blood cells descend to the bottom of the microtube.

## COMPOSITION

Each microtube of the DG Gel Coombs card contains polymerised dextrans in buffered medium with preservatives and mixed with anti-human globulin. The microtubes are identified by the front label of the card:

- microtubes **AHG**: Coombs, buffered low ionic strength solution (LISS) with polyspecific anti-human globulin. Mixture of rabbit polyclonal anti-IgG and monoclonal anti-C3d antibodies, IgM antibodies of murine origin, clone 12011 D10

Reagent ready to use. Use the microtubes immediately once the seal has been opened.

## STABILITY

DG Gel Coombs is stable up to the expiry date stated on the label, with the seal intact and stored at 2-25 °C in the position indicated on the outer packaging. Do not freeze.

## MATERIAL REQUIRED BUT NOT PROVIDED

- 10 µl, 25 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Incubator for DG Gel cards.
- Centrifuge for DG Gel cards.
- Reagent red blood cells for screening/identification of unexpected antibodies from Diagnostic Grifols, S.A.
- Haemoclassification sera.

## SAMPLES

Recently extracted blood samples collected with the usual blood bank anticoagulants (obtaining plasma) or without anticoagulants (obtaining serum).

Do not use haemolysed, cloudy or contaminated samples or with clot presence.

The procedure for extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives<sup>3,4</sup> and following the instructions of the manufacturer of the material used for collecting the sample.

- Screening/identification of unexpected antibodies, crossmatch tests and autocontrol: use serum or plasma. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction; or frozen samples (-20 °C to -80 °C), up to the time of performing the test.
- Direct Coombs, crossmatch tests, autocontrol and red blood cell typing: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in ACD, CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.

If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.

## TEST METHOD

DG Gel Coombs can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Allow the samples and reagents to reach room temperature (18-25 °C).

Inspect the state of the cards before use (see section LIMITATIONS "In relation to the product").

Identify the cards and samples to be used.

Each card is individually identified by a bar code. Unless a bar code reader is used, identify them manually.

If manual method is used: carefully peel off the metal film that covers the microtubes to prevent cross-contaminations among them and carefully dispense the red blood cell suspension, avoiding the pipette tip to come into contact with the wall or the content of the microtubes.

### Manual method:

#### 1. Indirect Coombs:

##### 1.1. Screening/identification of unexpected antibodies:

- Homogenise the vials of reagent red blood cells for screening/identification of unexpected antibodies.
- Dispense 50 µl of reagent red blood cells into the corresponding microtubes.
- Add 25 µl of the patient's serum or plasma.
- Incubate 15 minutes at 37 °C.

##### 1.2. Crossmatch tests (CT):

- Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Use the donor's red blood cells for the major CT and the recipient's red blood cells for the minor CT. Ensure the resuspension of the red blood cells before use.
- Dispense 50 µl of the donor's 1% red blood cell suspension (major CT) or 50 µl of the recipient's 1% red blood cell suspension (minor CT) into the corresponding microtube.
- Add 25 µl of the recipient's serum or plasma (major CT) or 25 µl of the donor's serum or plasma (minor CT).
- Incubate 15 minutes at 37 °C.

##### 1.3. Autocontrol:

- Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol) for the patient. Ensure the resuspension of the red blood cells before use.
- Dispense 50 µl of the patient's 1% red blood cell suspension into the corresponding microtube.
- Add 25 µl of the patient's serum or plasma.
- Incubate 15 minutes at 37 °C.

##### 1.4. Red blood cell typing:

Follow the instructions for use of the haemoclassification serum used.

#### 2. Direct Coombs:

- Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol) for the patient. Ensure the resuspension of the red blood cells before use.
- Add 50 µl of the patient's 1% red blood cell suspension into the corresponding microtube.

#### 3. Centrifuge in centrifuge for DG Gel cards.

#### 4. Read the results.

## RESULTS

### Result reading

<b>Negative:</b>	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column
	+/-	Scarce small-sized agglutinations in the lower half of the column
	1+	Some small-sized agglutinations in the column
	2+	Small or medium-sized agglutinations throughout the column
<b>Positive:</b>	3+	Upper band of medium-sized agglutinations in the upper half of the column
	4+	Band of agglutinated red blood cells in the upper part of the column
<b>DP</b>		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

Stability of the results: it is recommended an immediate reading of the results after centrifuging the cards. Do not leave processed cards in horizontal position.

If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2-8 °C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results:

Indirect Coombs and Direct Coombs tests: tests determined by the result obtained in microtubes AHG.

Red blood cell typing: follow the instructions for use of the haemoclassification serum used.

#### Notes:

1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
2. It is recommended to check/investigate any discrepant results obtained.
3. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
4. Direct Coombs: in newborn's samples, a negative result indicates the absence of detectable antibodies, and a positive result that the red blood cells are sensitised (coated with intrauterine antibodies).

#### QUALITY CONTROL

1. It is recommendable to include known positive and negative controls in each series of tests. It is recommended to include weak positive controls for screening/identification of unexpected antibodies and crossmatch tests, and positive controls with heterozygote antigen expression in the red blood cell typing.
2. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

#### PERFORMANCE CHARACTERISTICS

Although there is no described procedure or technique capable of detecting all the possible unexpected antibodies present in a sample population with absolute certainty, performance evaluation studies<sup>5,6</sup> carried out with DG Gel Coombs card back up that it owns performance characteristics conforming to the intended use of the product. These studies have included Indirect Coombs and Direct Coombs tests, being the results obtained comparable to those obtained with other established products of equivalent intended use.

The precision of the reagents present in the DG Gel Coombs card has been determined in a study<sup>7</sup>, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

#### LIMITATIONS

##### In relation to the sample:

Do not use haemolysed, cloudy or contaminated samples or with clot presence.

Sample: red blood cells

1. The presence of some drugs, dextran solutions or remains of silicone gel from the extraction tube in the sample, can induce a positive result in the Direct Coombs<sup>8</sup>.

Sample: serum/plasma

1. If plasma is used, complement-dependent haemolytic reactions may not be detected.
2. If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such. It is recommended to recloot the serum during 10 minutes at 37 °C, centrifuge and repeat the test<sup>8</sup>.
3. The presence of high concentrations of immunoglobulins or other serum proteins in the sample can neutralise the polyspecific anti-human globulin, even after numerous washes<sup>8</sup>.

##### In relation to the product:

Before use, inspect the microtubes of the cards:

1. Do not use the card if microbiological contamination or change in microtube colour is observed.
2. If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
3. If trapped bubbles, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.
4. In the 24-hour delayed reading of processed cards with weak positive samples, a loss in the agglutination intensity may be observed.

##### In relation to the method:

1. The Indirect Coombs test at 37 °C in gel or glass spheres techniques has been reported to show a level of sensitivity lower than the obtained with the tube technique, in the detection of weak agglutination reactions of the ABO system<sup>9</sup>. To ensure the ABO compatibility between the recipient's and donor's blood, a serological (saline at room temperature with immediate centrifuging) or computer-based (electronic crossmatch) compatibility confirmation, is recommended.
2. Screening and identification of unexpected antibodies/crossmatch tests: the test is acceptable if the volumes of serum or plasma are increased from 25 µl to 50 µl. This variation in the concentration of antibodies brings the antigen/antibody ratio down and can improve the detection of antibodies at very low concentration<sup>8</sup>.

#### WARNINGS AND PRECAUTIONS

The use of "in vitro" diagnostic reagents for professional use requires taking into account the following indications:

- The reagents of the DG Gel Coombs card are not of human origin. HBV, HCV or HIV contamination is excluded. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use red blood cells treated with enzymes in microtubes with anti-human globulin (AHG).
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

#### BIBLIOGRAPHY

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-3), October 2003.
6. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
7. Internal study Diagnostic Grifols, S.A. (DR-103/03), July 2003.
8. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
9. Phillips P et al. An explanation and the clinical significance of the failure of microcolumn tests to detect weak ABO and other antibodies. Transfusion Medicine, 7: 47-53, 1997.

#### PRESENTATION

210342 DG Gel Coombs 50 cards Profile: 8x(AHG)

Revision date: June 2010

This document is available in several languages. The translations have been made from the master document in Spanish. In the event of doubts or discrepancies, the wording in the master document in Spanish shall take precedence.

	"In vitro" diagnostic medical device
	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
	Catalogue number
	Cards

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## LITERATŪRA



1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-3), October 2003.
6. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
8. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
9. Phillips P et al. An explanation and the clinical significance of the failure of microcolumn tests to detect weak ABO and other antibodies. Transfusion Medicine, 7: 47-53, 1997..

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/logotipas/

Read carefully before using the product. For use in "in vitro" diagnosis only

## INTENDED USE

Confirmation of the blood groups of the ABO and Rh (D) system, in gel technique.

## INTRODUCTION

The ABO system was the first human blood group system discovered, by Landsteiner in 1900<sup>1</sup>, and is still the most important in transfusion practice.

The ABO system is defined by the presence or absence of the A (ABO1) and/or B (ABO2) antigens in the red blood cells and by the presence of antibodies in the serum corresponding to the antigen or antigens missing in the red blood cells.

The determination of Rh (D) is defined by the presence or absence of the D antigen (RH1) in the red blood cells.

The anti-A, anti-B and anti-D<sup>0</sup>- reagents are used to perform the ABO and Rh blood group typing.

## PRINCIPLE OF THE TEST

The principle of the test is based on the gel technique described by Y. Lapierre<sup>2</sup> for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

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Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains specific antibodies or a buffer.

The microtubes that contain the specific antibodies into the gel solution act as a reaction medium and the red blood cells agglutinate in contact with the antibodies.

The microtubes without antibodies are used in techniques in which the antibodies react directly with the red blood cells in the incubation chamber and for controls

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## REAGENTS

Each microtube of the DG Gel Confirm P card contains polymerised dextrans in buffered medium, with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtubes A: monoclonal anti-A (IgM antibodies of murine origin, clone Birma-1).
- microtubes B: monoclonal anti-B (IgM antibodies of murine origin, LB-2).
- microtubes D<sup>0</sup>: monoclonal anti-D (IgM antibodies of human origin, clone MS-201).
- microtubes Ctl.: buffered solution without antibodies (control microtubes).

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## STABILITY

DG Gel Confirm P is stable up to the expiry date stated on the label, with the seal intact and stored at 2-25 °C in the position indicated on the outer packaging. Do not freeze.

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- 10 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Centrifuge for DG Gel cards.

## SAMPLES

Blood samples recently collected with the routinely used anticoagulants (EDTA, citrate and heparine). Do not use haemolysed, cloudy or contaminated samples or with clot presence.

The procedure of extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives<sup>3,4</sup> and following the instructions of the manufacturer of the material used for collecting the sample.

- Determination of the antigens of the ABO/Rh system: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.

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## TEST METHOD

DG Gel Confirm P can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

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## Manual method:

1. Confirmation of the ABO/Rh groups (microtubes A/B/D<sup>0</sup>/Ctl.):
  - Prepare a 5% red blood cell suspension in DG Gel Sol (50 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.
  - Add 10 µl of 5% red blood cell suspension to each of the microtubes indicated.
2. Centrifuge in centrifuge for DG Gel cards.
3. Read the results.

## RESULTS

### Result Reading:

<b>Negative:</b>	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column
<b>Positive:</b>	+/-	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
<b>DP</b>		Double Population (double band of red blood cells, at the bottom and in the upper part of the column).

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+ = positive  
0 = negative

## Notes:

- The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
- The microtube Ctl. must be negative. If it is positive, invalidate the test.  
Repeat the determination after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
- ABO/Rh system: in the event of +/- to 3+ reactions, presence of a weak antigen must be investigated.
- Negative reactions with the anti-D<sup>W</sup>- reagent must be verified, using other reagents and techniques which may detect different variants of the D antigen.
- The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
- Occasionally, a red blood cell retention in the incubation chamber may occur with positive 4+ samples, without interfering in result reading.

## QUALITY CONTROL

- It is recommendable to include known positive and negative controls in each series of tests.
- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

## PERFORMANCE CHARACTERISTICS

### Diagnostic sensitivity and specificity:

#### ABO/Rh system

The diagnostic sensitivity and specificity of the antibodies present in the DG Gel Confirm P card<sup>5,6</sup>, for the determination of the antigens of the ABO and Rh systems, have been studied in a representative number of positive and negative samples.

Antibody	No. of Samples	Sensitivity <sup>(a)</sup>	Specificity <sup>(b)</sup>
Anti-A	3041 <sup>5</sup>	100%	100%
Anti-B	3046 <sup>5</sup>	100%	100%
Anti-D <sup>W</sup> -	1329 <sup>(c),6</sup>	99,7%	100%

(a) Sensitivity: (number of true positive results / (number of true positive results + number of false negative results)) x 100

(b) Specificity: (number of true negative results / (number of true negative results + number of false positive results)) x 100

(c) These samples include 49 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D<sup>W</sup>- reagent is 93,9%.

## Precision:

The precision<sup>7,8</sup> of the reagents present in the DG Gel Confirm P card has been determined, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

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Sample: red blood cells

- Red blood cells from individuals with A or B variants may present a weak expression of the antigens. Antigen expression may be weakened in red blood cells of people with leukaemia or other malignant diseases<sup>9</sup>.
- Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test<sup>9</sup>.
- Transfused patients or those subjected to bone marrow transplant may present images of double population<sup>9</sup>.
- Patients with high-potency autoantibodies may coat the red blood cells completely, causing spontaneous agglutination<sup>9</sup>.

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- Do not use the card if microbiological contamination or change in microtube colour is observed.
- If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- If trapped bubbles in the gel, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.
- The anti-A reagent may not react with weak variants of Ax phenotype.

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- The reagents of the DG Gel Confirm P card of human monoclonal origin are manufactured using a non-reactive material for the HB<sub>s</sub> antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

## BIBLIOGRAPHY

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
- NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
- External evaluation Diagnostic Grifols, S.A. (REGD-0002998)M Noviembre 2007.
- External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013
- Internal study Diagnostic Grifols, S.A. (DR-I03/03), Julio 2003.
- Internal study Diagnostic Grifols, S.A. (DG-R1013/07-3), Mayo 2007.
- Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

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## LITERATŪRA



1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
4. NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998)M Noviembre 2007.
6. External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013
7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), Julio 2003.
8. Internal study Diagnostic Grifols, S.A. (DG-RI013/07-3), Mayo 2007.
9. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

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## INTENDED USE

DG Gel Sol is a reagent for preparing red blood cell suspensions used in gel techniques (DG Gel).

## PRINCIPLE

A reduction in the ionic strength of the medium of the red blood cell suspension enhances significantly the association of the antibody with the antigen. DG Gel Sol is a low ionic strength solution that facilitates the attachment of the antibody to the red blood cells by lowering the density of the cation cloud around them, thus improving sensitization.

## COMPOSITION

Each DG Gel Sol kit contains:

### DG Gel Sol:

- 2 vials of 100 ml, liquid. Each vial contains low ionic strength buffered solution, being its main components glycine 1.37% and glucose 0.85%.

## Preparation of the reagents

### DG Gel Sol:

Ready to use.

Stability once opened: at 2-8 °C, it has the same stability as when unopened, provided that it is handled and stored correctly, avoiding contamination of the contents and prolonged exposure to temperatures that are not suitable for storage.

Do not freeze.

## STABILITY

All unopened reagents are stable up to the expiry date stated on the label if stored at 2-8 °C.

## METHOD

DG Gel Sol can be used manually or with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

### Manual method (DG Gel):

Follow the instructions for use of the DG Gel card to be used, particularly the points that indicate the use of DG Gel Sol.

## LIMITATIONS

### Limitations in relation to the reagent:

Before use, inspect the reagent.

1. Do not use if the container is deteriorated and there is a loss of content.
2. Do not use if microbial or fungal contamination is observed.

## WARNINGS AND PRECAUTIONS

The use of “in vitro” diagnosis reagents for professional use requires taking into account the following indications:

- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positives or negatives.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Follow the recommended procedures strictly, periodically checking the instruments used.
- Do not use the product beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have doubts or need further information on the use of this product, consult the authorised distributor in your country.

## PRESENTATION

210354 DG Gel Sol 2x100ml

Revision date: April 2004

This document is available in several languages. The translations have been made from the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

IVD	In vitro diagnostic medical device
LOT	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
REF	Catalogue number

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Diagnostic Grifols, S.A. Passeig Fluvial 24, 08150 Parets del Vallès (Barcelona), Ispanija