NADAL[®] C. difficile Toxins A+B/GDH Test (test cassette)

REF 582016N-10

DE	Gebrauchsanweisung	2	Symbols	23
EN	Instructions for use	6	Our Teams	24
ES	Instrucciones de uso	9		
PL	Sposób użycia	13		
PT	Instruções de Utilização	17		

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14. Referenzen

1. Wren, M.W.D., et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

Rev. 0, 2016-06-15 OM/UJa

1. Intended Use

The NADAL® C. difficile Toxins A+B/GDH Test is a lateral flow, immunochromatographic rapid test for the qualitative detection of *Clostridium difficile* Toxin A and Toxin B (C. difficile Toxins A+B) as well as *Clostridium difficile* glutamate dehydrogenase (C. difficile GDH) in human faecal specimens. The test is intended for use as an aid in the diagnosis of *Clostridium difficile* infection and is designed for professional use only.

2. Introduction and Clinical Significance

Clostridium difficile (C. difficile) is an anaerobic gram-positive, spore-forming bacterium. The ability of C. difficile to form spores is the key feature which enables this bacterium to survive in patients and the physical environment for long periods and thereby facilitates its transmission. C. difficile is transmitted by the faecal-oral route.

C. difficile is the principal pathogen related to antibiotic-associated diarrhoea and/or pseudomembranous colitis in hospitalised patients. Mature colonic bacterial flora in a healthy adult is generally resistant to C. difficile colonisation. However, if the normal enteric flora is altered, resistance to colonisation is lost. Any factor associated with alteration of the normal enteric flora increases the risk of C. difficile colonisation after exposure to antibiotics, especially those with broad-spectrum activity such as penicillins, cephalosporins and clindamycin.

C. difficile releases two high-molecular-weight toxins: toxins A and B, which are responsible for the clinical manifestations ranging from mild, self-limiting, watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death.

Clostridium difficile glutamate dehydrogenase (C. difficile GDH) is an enzyme produced in large quantities by all toxigenic and non-toxicogenic strains, making it an excellent biomarker.

3. Test Principle

The NADAL® C. difficile Toxins A+B/GDH Test is a lateral flow, immunochromatographic rapid test for the qualitative detection of C. difficile Toxin A and Toxin B, as well as C. difficile GDH in human faecal samples.

The test cassette consists of a plastic housing containing three different test strips for the qualitative detection of C. difficile Toxin A, C. difficile Toxin B or C. difficile GDH. Coloured, monoclonal anti-C. difficile Toxin A, anti-C. difficile Toxin B and anti-C. difficile GDH antibodies colloidal gold conjugates are precoated onto the sample pad of the respective membrane.

After sample collection, the faecal sample is extracted in a tube containing buffer and a few drops of this extract are added to each sample well of the test cassette.

As the sample migrates along the membranes, labelled antibody-dye conjugates bind to Toxin A, Toxins B or GDH antigens (if present in the sample), forming antibody-antigen complexes. These complexes bind to the anti-Toxin A, anti-Toxin B and GDH antibodies in the test line regions of the test cassette, producing rose-pink lines.

In the absence of Toxin A, Toxin B or GDH, no lines appear in the corresponding test line region of the test cassette. The reaction mixture continues migrating along the membranes of the test cassette. Unbound conjugates bind to the reagents in

each control line region, producing a rose-pink line, which demonstrates that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® C. difficile Toxins A+B/GDH test cassettes
- 10 disposable plastic pipettes (for liquid faecal samples)
- 10 specimen collection tubes with buffer (2 mL)
- 1 package insert

5. Additional Materials Required

- Timer

6. Storage & Stability

Store unused test cassettes unopened at 4-30°C. The test cassette must remain in the sealed foil pouch until use. The test cassette is stable until the expiration date printed on the sealed foil pouch. Do not freeze test kits.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Humidity and temperature can adversely affect test results.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Used testing materials should be disposed of according to local regulations.
- All reagents and materials coming into contact with potentially infectious specimens should be treated with appropriate disinfectants or autoclaved at 121°C for at least one hour.

8. Specimen Collection and Preparation

Specimen Collection and Storage:

Stool specimens should be collected as soon as possible after the onset of symptoms. Diluted samples can be stored at 2-

8°C for 3 days without interference with assay performance. For long term storage of undiluted specimens, it is recommended to store the specimens at -20°C or colder. Repeated freezing and thawing of samples is not recommended, as it may cause erroneous results.


Note: Do not collect specimens in containers containing media, preservatives, animal serum or detergents, as they may interfere with the test.

Specimen preparation:

Label the specimen collection tube with a patient or control identification as well as the date.

Solid faecal specimens

Unscrew the cap of the tube and insert the applicator stick into three different spots of the faecal specimen to collect the sample volume of a pea.

Place the applicator stick, along with the sample, back into the specimen collection tube and close it properly.


Liquid stool samples

Aspirate the faecal specimen with a pipette, add 200 µL (5 drops) into the specimen collection tube with buffer and close it properly.

9. Test Procedure

Bring the tests and samples with buffer to room temperature (15-30°C) prior to testing. Do not open foil pouches until you are ready to perform the assay.

1. Remove the NADAL® C. difficile Toxins A+B/GDH test cassette from the sealed foil pouch.
2. Extract the sample in the specimen collection tube by swirling the tube vigorously for a long time. Break off the tip of the specimen collection tube.
3. Holding the specimen collection tube vertically, dispense 3 drops (100 µL) of extracted sample into each sample well of the test cassette.
4. Read the result at 10 minutes after adding the sample to the test cassette. Do not interpret results after more than 10 minutes.

10. Result Interpretation

Positive for GDH:

One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'T' for GDH.

Positive for Toxin A (TOX A):


One coloured line develops in each control line region 'C' and one coloured line

develops in the test line region 'T' for Toxin A (TOX A).

Positive for Toxin B (TOX B):


One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'T' for Toxin B (TOX B).

Note: A combination of positive results for various parameters is possible.

Negative:

One coloured line appears in each control line region 'C'. No coloured lines appear in the test line regions 'T' for GDH, Toxin A (TOX A) and Toxin B (TOX B).

Invalid:

At least one control line 'C' fails to appear. Results from any test which has not produced all control lines at the specified reading time must be discarded.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

Internal Control:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

External Control:

The use of external positive and negative controls is recommended to assure the proper performance of the assay.

12. Limitations

- The NADAL® C. difficile Toxins A+B/GDH Test is specifically designed to detect *C. difficile* Toxin A, *C. difficile* Toxin B or *C. difficile* GDH antigens in faecal samples.
- As with any *in-vitro* diagnostic procedure, the physician should confirm the test results with other clinical methods.
- A negative result does not necessarily exclude a *C. difficile* infection. It can be caused by proteolytic digestion of the toxins due to improper specimen storage. If a reasonable suspicion of an infection exists, another stool specimen should be tested.
- A positive result does not exclude the presence of other pathogens.
- The colour of control and test lines may slightly change depending on the faecal sample coloration, e.g. dark green lines (instead of pink lines) have been reported when assaying greenish or darkish faecal samples. This stool

coloration appears in case of treatment of iron deficiency with ferrous fumarate.

13. Performance Characteristics

GDH

Analytical sensitivity

The performance of the NADAL® C. difficile Toxins A+B/GDH Test has been evaluated using a range of dilutions prepared from a commercially available purified GDH antigen. The lower detection limit of the test for GDH is 5 ng/mL.

A competitor GDH test showed an analytical sensitivity of 10 ng/mL when the same GDH antigen was assayed.

Diagnostic sensitivity and specificity

A comparative study between DiaSorin analyser and the NADAL® C. difficile Toxins A+B/GDH Test has been performed using a panel of 59 human faecal samples.

NADAL® C. difficile Toxins A+B/GDH Test (GDH)	DiaSorin analyser			
		+	-	Total
	+	13	1	14
	-	3	39	42
Total		16	40	56

Diagnostic sensitivity: 81% (13/16)

Diagnostic specificity: 97.5% (39/40)

Overall agreement: 93% (52/56)

Two samples were identified as equivocal using the DiaSorin analyser and one sample was found doubtful using the NADAL® C. difficile Toxins A+B/GDH Test. Therefore, they were not included in the calculation.

Cross-reactivity

The following organisms were tested and did not show any cross-reactivity with the NADAL® C. difficile Toxins A+B/GDH Test for GDH:

- *Campylobacter* spp
- *Escherichia coli*
- *Helicobacter pylori*
- *Listeria monocytogenes*
- *Salmonella* spp
- *Shigella* spp
- *Staphylococcus aureus*
- *Yersinia enterocolitica*
- *C. difficile* Toxin A
- *C. difficile* Toxin B
- Rotavirus
- Adenovirus

Intra-assay reproducibility

The intra-assay reproducibility was evaluated by testing 15 replicates of three samples dilution containing different GDH concentrations. All series showed converging results in conformance with the expected results. Intra-lot reproducibility of the NADAL® C. difficile Toxins A+B/GDH Test is 100% for GDH.

C. difficile Toxin A

Analytical sensitivity

The performance of the C. difficile Toxins A+B/GDH Test has been evaluated, using a range of dilutions prepared from a

commercially available purified *C. difficile* Toxin A antigen. The lower detection limit of the test for Toxin A is of 4 ng/mL.

Diagnostic sensitivity and specificity

A study was performed on 63 faecal samples, using the NADAL® C. difficile Toxins A+B/GDH Test in comparison with the cytotoxicity and cell culture methods.

NADAL® C. difficile Toxins A+B/GDH Test (Toxin A)	Cytotoxicity			
		+	-	Total
	+	13	4	17
	-	2	44	46
Total		15	48	63

All samples showed identical results using the cytotoxicity and the cell culture methods.

Compared to the cytotoxicity method, the C. difficile Toxins A+B/GDH Test (Toxin A) showed

Sensitivity: 86.7% (13/15)

Specificity: 91.7% (44/48)

Cross-reactivity

The NADAL® C. difficile Toxins A+B/GDH Test showed consistently negative results for Toxin A with Toxin B concentrations of up to 500 ng/mL.

C. difficile Toxin B

Analytical sensitivity

The sensitivity of the NADAL® C. difficile Toxins A+B/GDH Test has been evaluated using a range of dilutions prepared from a commercially available purified *C. difficile* Toxin B antigen. The lower detection limit of the test for Toxin B is 5 ng/mL.

Specificity

A purified *C. difficile* Toxin A antigen was used to determine the specificity of the NADAL® C. difficile Toxins A+B/GDH Test (Toxin B). The NADAL® C. difficile Toxins A+B/GDH Test showed consistently negative results for Toxin B with Toxin A concentrations of up to 500 ng/mL.

14. References

1. Wren, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

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14. Referencias

- Wren, M.W.D, et al. "Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and

- toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499
 3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

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


14. Referências

1. Wren, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

Rev. 0, 2016-06-15 AL

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeas	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensning
	Código de lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestilingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel – UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH